JAY[®]

J3 BACK

JAY[®] J3 Back

IMPORTANT CONSUMER INFORMATION

NOTICE: This manual contains important instructions that must be passed on to the user of this product. Please do not remove this manual before delivery to the end user.

SUPPLIER: This manual must be given to the user of this product.

USER: Before using this product, read this entire manual and save for future reference.

WARRANTY REGISTRATION: Visit www.SunriseMedical.com to register this product.

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Manuel d'instructions

SUNRISE MEDICAL LISTENS

Thank you for choosing a J3 Back Support System. We want to hear your questions or comments about this manual, the safety and reliability of your product, and the service you receive from your Authorized Sunrise dealer. Please feel free to write or call us at the address and telephone number below:

**Sunrise Medical (US) LLC
Customer Service Department
12002 Volunteer Blvd.
Mount Juliet, TN 37122
(800) 333-4000**

Be sure to return your warranty card, and let us know if you change your address. This will allow us to keep you up to date with information about safety, new products and options to increase your use and enjoyment of this wheelchair. If you lose your warranty card, call or write and we will gladly send you a new one.

FOR ANSWERS TO YOUR QUESTIONS

Your Authorized Dealer knows your product best, and can answer most of your questions about safety, use and maintenance. For future reference, fill in the following:

Supplier: _____

Address: _____

Telephone: _____

Serial #: _____ Date/Purchased: _____

ADDITIONAL INFORMATION YOU SHOULD KNOW

No component of this product was made with Natural Rubber Latex.

DISPOSAL AND RECYCLING INFORMATION

When this product reaches the end of its life, please take it to an approved collection or recycling point designated by your local or state government. This product is manufactured using a variety of materials. Your product should not be disposed of as ordinary household waste. You should dispose of your wheelchair properly, according to local laws and regulations. Most materials that are used in the construction of this product are fully recyclable. The separate collection and recycling of your product at the time of disposal will help conserve natural resources and ensure that it is disposed in a manner that protects the environment.

Ensure you are the legal owner of the product prior to arranging for the product disposal in accordance with the above recommendations

ENGLISH

Introduction.....3
 General Warnings.....4
 Use and Maintenance5
 Fitting the back to the Chair.....6
 Installation & Adjustment7
 Accessories.....10
 Warranty.....11

ESPAÑOL

FRANÇAIS

Sunrise Medical recommends that a clinician such as a doctor or therapist experienced in seating and positioning be consulted to determine if a J3 Back is appropriate. Backs should only be installed by an Authorized Sunrise Medical Dealer.

JAY J3 BACK

The Jay J3 Back is designed to provide the important benefits of a proper back support for a wide variety of users. It is designed to provide simple, comfortable and effective support to help maximize function and increase sitting tolerance. In addition, it is light-weight, durable and easy to use.

Maximum user weight:

300 lbs. (136 kg)

500 lbs. (227 kg) for sizes 22", 24" and 26" widths.


J3 BACK HARDWARE


Locking Pin


10mm Box End Wrench


4mm Hex Key


SH & RH Hardware


AH hardware


RHF Hardware

TH Hardware


FS Hardware

FA Hardware

J3 BACK

Mid Thoracic (example)

A. CENTER OF GRAVITY

⚠ WARNING

Installing a back support on a wheelchair may affect the center of gravity of the wheelchair, and may cause the wheelchair to tip backwards potentially resulting in injury to the rider.

B. ANTI-TIPS

⚠ WARNING

Anti-tip tubes can help keep your chair from tipping over backward in normal conditions.

1. Sunrise recommends the use of anti-tip tubes
2. ALWAYS assess the need for anti-tips or an amputee axle adapter brackets that could be added to the wheelchair to help increase the stability of the wheelchair and seating system.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

C. PUSHING OR LIFTING

⚠ WARNING

Do not use the J3 Back to move or push the wheelchair.

1. Do not use the J3 back handle to lift the wheelchair. It is intended only for Back removal.
2. Improper use of the back by pushing or pulling may cause the back to unexpectedly detach from the wheelchair.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

D. COMPATIBILITY

⚠ WARNING

The J3 Back is designed to be compatible with most wheelchairs with the following exception:

1. Wheelchairs that are angle adjustable with recline or tilt, and result in a back angle greater than 60° from the vertical

J3 back use with these types of wheelchairs may cause the back to accidentally disconnect from the wheelchair. This could result in a fall or potential injury of the rider and/or caregiver.


E. J3 BACK TRANSIT SAFETY 1 2

⚠ WARNING

The J3 back has been dynamically tested for use in a motor vehicle. This system conforms with ANSI-RESNA WC-20. If the J3 back is to be used in a transit situation, it must be mounted to a wheelchair that conforms to the performance requirements of ANSI/RESNA WC-19. For Rider safety, Please follow all installation, use, and maintenance instructions within this manual as well as the transit instructions listed below .

1. If possible and feasible, the rider should transfer to the Original Equipment Manufacturer vehicle seat and use the vehicle restraints.
2. The distance between the top of the user's shoulder and the top of the back should not exceed 6.5" (16.5 cm).
3. The wheelchair must be labeled as appropriate for use as a seat in a motor vehicle, dynamically tested to the performance requirements of ISO Standard 7176-19 and installed, used and maintained according to the manufacturer's instructions.
4. If the original wheelchair seat was replaced, the new seat must be approved for wheelchair transit, installed and used as indicated in the manufacturer's instructions.
5. The Wheelchair Tiedown and Occupant Restraint System (WTORS) should be compatible with the specific wheelchair, used as indicated in the manufacturer's instructions and should comply with the performance requirements of ISO Standard 10542.
6. The Wheelchair must be forward facing during transport.
7. In order to reduce the potential of injury to vehicle occupants wheelchair-mounted accessories such as trays and respiratory equipment should be removed and secured separately.
8. This product is intended for use by adults ranging in weight from 79 to 300lbs or up to 500 lbs. (227 kg) for widths over 22".
9. Postural support devices such as pelvic positioning belts, anterior trunk supports or postural supports such as lateral trunk supports should not be relied on for occupant restraint in a moving vehicle unless they are labeled as conforming to ISO Standard 16840-4.
10. Alterations or substitutions should not be made to the structure of the chair or seating system, parts, and/or components, without consulting the manufacturer.
11. Backshell mounting bracket must be attached in the "leg up" position for transit use. (See below.) Refer to the appropriate section for correct back hardware mounting instructions.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

⚠ WARNING

The owner of this equipment is responsible for making sure that it has been set up and adjusted by a **Authorized Dealer** under the advice of a healthcare professional. It may require periodic safety and function checks or certain in-use adjustments that can be performed by the owner, caregiver or Authorized Dealer if desired.

A. MAINTENANCE AND CLEANING

1. Proper maintenance and cleaning will improve performance and extend the useful life of your equipment.
2. Clean your equipment regularly. This will help you find loose or worn parts and make your equipment easier to use. You will need a mild detergent solution and cleaning rags.
3. If discovered, have loose, worn, bent or damaged parts replaced **before using the back support**.

B. CRITICAL SAFETY CHECKS

1. Fasteners:

Sunrise Medical recommends that all fasteners be checked for wear, such as loose bolts or broken components **Every 6 months**. Loose fasteners should be retightened according to the installation instructions.

2. Replacing components:

Contact your Authorized Dealer immediately to replace any components. Do not continue to use the system after identifying loose or broken components.

C. EASY CLEANING 3

1. Back Cover

- a. Remove the cover from backshell by detaching the Velcro®. Unzip the zipper on the bottom of the cover.
- b. Remove the foam and turn the cover inside out.
- c. Zip the cover closed to avoid snagging.
- d. Machine wash in warm water (50°C).
- e. Drip and tumble dry in on Low Heat.

NOTE – Do not dry clean covers or use industrial washers or dryers to clean covers. Do not steam autoclave. Do not bleach. Laundry care symbol definitions are shown (Fig. 3).

2. Back shell

- a. The backshell should be wiped clean with soap and water on a regular basis to prevent buildup of dirt and grime.

3. Positioning components

- a. Remove the components one at a time, noting the placement for re-assembly.
- b. Positioning components can be wiped clean with soap and water.
- c. **DO NOT** immerse positioning components underwater.

D. REMOVAL AND REPLACEMENT 4

After installation/adjustment, the back can be easily removed from the wheelchair if necessary.

1. Removal:

- a. Push both release levers forward to unlock the hardware (Fig 4).
- b. Once both levers are in the forward position, lift the back straight up, using the handle located on the top of the back.

2. Proper re-attachment after removal


- a. To re-attach the back and hardware receivers, line up the pins to the receivers.
- b. Push the back into place on the receivers.
- c. If properly inserted, the pins will “click” into a locked position.

3. Locking Pins

The locking pins provide a safety feature which prevents the back from being released unintentionally from the wheelchair.

- a. To lock, insert the pins (A) into the mounting receivers (B).
- b. Remove the pins if the quick-release capability is desired.

Laundry Care Symbols		3
	Machine wash in warm water (50° C)	
	Only non-chlorine bleach, when needed	
	Drip dry	
	Do not iron	
	Do not dry clean	

⚠ WARNING

If you are uncertain that you have heard the locking “click”, pull up on the back using the handle to test and make sure it has been securely locked. If you fail to heed this warnings damage to your equipment, a fall, or loss of control may occur and cause severe injury to the rider or others.

⚠ WARNING

The owner of this equipment is responsible for making sure that it has been set up and adjusted by an Authorized Dealer under the advice of a healthcare professional. The back and hardware may require periodic safety and function checks. Adjustments should be performed by your Authorized Dealer. Always use parts and/or accessories that have been recommended or approved by Sunrise Medical.

The J3 back is designed to be compatible with most wheelchairs with the following exception: Wheelchairs that are angle adjustable with recline or tilt that result in a back angle greater than 60° from the vertical should not be used with the J3 back system.

A. WHEELCHAIR DIMENSIONS 5 6

1. The J3 Back is designed to replace the wheelchair's sling upholstery. Prior to installation determine if the wheelchair used has compatible back canes, by measuring the width between your back canes as seen in (Fig. 5).
2. Next check that the back canes are a compatible diameter. The J3 back hardware can be mounted on wheelchair back canes ranging from 3/4" to 1 1/8". Measure your back cane diameter (Fig. 6) and refer to Table 2 to ensure your back canes are compatible.
3. If you cannot determine that the back canes are compatible please contact your authorized Sunrise Medical provider or Sunrise Medical customer service.


TABLE 1 J3 Back Width* (all contours)	Width of back canes from outside to outside of the tubes (A)			
	Minimum width		Maximum width	
	inch	cm	inch	cm
12"	11.7	29.7	14.3	36.3
14"	13.7	34.9	16.3	41.3
16"	15.7	40.0	18.3	46.4
18"	17.7	45.1	20.3	51.5
20"	19.7	50.1	22.3	56.5
22"	21.7	55.1	24.3	61.7
24"	23.7	60.2	26.3	66.8
26"	25.7	65.3	28.3	71.9

* Back Width is identified in your back's part number.
Example: Part number J3PAMTS14SH is a 14" back.


TABLE 2 J3 Back Width (all contours)	Cane diameter (B)			
	Minimum width		Maximum width	
	inch	cm	inch	cm
12" - 26"	3/4	1.9	1 1/8	2.85

A. J3 (SH, RH, RHF, & TH) HARDWARE

Please read the following instructions before beginning the installation. To install the J3 Back, with 2-point hardware, it is best to begin without a user in the wheelchair. Once the back has been installed the user can be seated back in the wheelchair for minor adjustments.

1. Tools required (included with the J3 back):

- a. 4mm hex key
- b. 10mm box end wrench

B. BACK HARDWARE INSTALLATION 7 8 9 10 11

1. Remove existing back components

- a. Remove the existing wheelchair back according to the wheelchair manufacturer's instructions.
- b. Remove any pre-existing back hardware or components.

2. Locating the hardware receivers

- a. The backshell should come with mounting brackets (A) attached.
- b. Visually determine the optimum clamp receiver location (B) on the back canes. The mounting hardware should be at approximately equal heights on each back post (Fig. 7), and parallel to the seat frame (Fig. 8).

NOTE—A higher location on the back canes are generally better.

3. Attaching the Hardware receivers

- a. Using the 4mm hex key, Loosen the clamp screws (Fig. 9) (C)

NOTE—For Large tube diameters, the clamp screws may have to be removed completely.

- b. Attach and align the receiver to the wheelchair back canes (D).
- c. Hand tighten the clamp screws (C).
- d. The receivers may require additional lateral adjustments to ensure proper alignment.

4. Attaching and aligning the J3 Back shell.

- a. Using a 10mm wrench, loosen the bracket nuts (E), and mounting pin bolts (H) until the hardware can move easily in all directions.
- b. Attach the J3 Back to the hardware receivers by inserting left and right mounting pins (F) in the left and right receivers (G).
- c. Locate the back to the approximate location, relative to both the chair and the user.
- d. Tighten bolts on the inside of the mounting pin (F). The recommended torque specification is 100-110 in/lbs. (11.3-12.4 Nm)
- e. Once all components have been properly aligned, tighten the clamp screws (C) to approximately 75-85 in-lbs (8.5-9.3 Nm).
- f. Finish by inserting the locking pin into the receiver, and releasing a few times to ensure smooth operation.






NOTE—Back depth and angle may require additional adjustment to ensure proper fit to the user.

5. Installing RHF hardware to the J3 Back shell.

- a. Using the 4mm hex key, Loosen the clamp screws (Fig. 9) (C)

NOTE—For Large tube diameters, the clamp screws may have to be removed completely.

- b. Attach and align the clamp to the wheelchair back canes (D).
- c. Hand tighten the clamp screws (C).
- d. The clamps may require additional lateral adjustments to ensure proper alignment.
- e. Locate the back to the approximate location, relative to both the chair and the user.
- f. Instead of the pin assembly, the clamp hardware attaches directly to the backshell mounting bracket with hardware (H) shown in (Fig. 11).
- e. The recommended torque specification is 100-110 in/lbs. (11.3-12.4 Nm)
- e. Once all components have been properly aligned, tighten the clamp screws (C) to approximately 75-85 in-lbs (8.5-9.3 Nm).
- f. Finish by inserting the locking pin into the receiver.

C. J3 (FS & FA) HARDWARE 12 13 14 15 16

Please read the following instructions before beginning the installation. To install the J3 Back with 4-point hardware, it is best to begin without a user in the wheelchair. Once the back has been installed and minor adjustments are needed, the user can be seated back in the wheelchair.

1. Tools required (included with the J3 back):

- 4mm hex key
- 10mm box end wrench
- 5/8" box-end wrench

2. Remove existing back components

- Remove the existing wheelchair back according to the wheelchair manufacturer's instructions.
- Remove any pre-existing back hardware or components.

3. Locating the hardware receivers

- The backshell should come with mounting brackets (A) attached.
- Visually determine the optimum clamp receiver location (B) on the back canes. The mounting hardware should be at approximately equal heights on each back post (D) and parallel to the seat frame (Fig. 12).

NOTE—A higher location on the back canes are generally better.

4. Attaching the Hardware receivers

- Using the 4mm hex key, Loosen the clamp screws (C) and attach the clamp receiver (B) to the backpost.

3. Attaching the Hardware receivers 13





- Using the 4mm hex key, Loosen the clamp screws (C)

NOTE—For Large tube diameters, the clamp screws may have to be removed completely.

- Attach and align the receiver to the wheelchair back canes.
- Hand tighten the clamp screws (C).
- Using a 10mm wrench, Loosen the bracket nuts (E) (Fig. 14) and the mounting pin bolts (I) (Fig. 15) until you can move the hardware easily in all directions for adjustment.
- The receivers may require additional lateral adjustments to ensure proper alignment.

4. Attaching and aligning with the upper hardware receivers. 14 15

- Attach the J3 Back to the hardware receivers by inserting left and right mounting pins (F) in the left and right receivers (H).
- Locate the back to the approximate location, relative to both the chair and the user.
- Tighten bolts on the inside of the mounting pin (I) (Fig. 15), to a torque of 100-110 in/lbs. (11.3-12.4 Nm)
- Tighten the bracket nuts (E). Torque to 75-85 in-lbs (8.5-9.3 Nm)
- Test the back release by pressing the release levers (G) forward and removing the back off of the chair. Proper alignment of the mounting hardware and back is correct when it is possible to have a smooth attachment and release of the back assembly.
- If proper alignment has not been achieved, adjust the components until properly aligned.
- Once properly aligned, tighten the hardware receiver clamp screws (C), Torque to 75-85 in-lbs (8.5-9.3 Nm).
- Finish by inserting the locking pin (Fig.15) (J) into the receiver.

5. Attaching and aligning the lower hardware receivers. 16

NOTE - Make sure the upper and lower brackets are aligned and parallel to the Seat frame.

- Test the back release for the upper mounting hardware by pressing the release levers (A) forward and removing the back from the chair. Proper alignment of the mounting hardware and back can be determined by the smooth attachment and release of the back assembly.
- If needed, continue small adjustments until the back can be attached, and/or removed smoothly.
- Loosen screws (B) to install lower receivers on back canes (C)
- Using a 5/8" wrench, loosen the nut (D) that secures the lower mounting pin (E) so that it slides fore and aft.
- Slide lower receiver up (F) until it surrounds the lower mounting pin (E).

NOTE - Leave 1/8" to 1/4" of space between the bottom of the lower mounting pin (G) and the bottom of the lower receiver channel (H).

- Once properly aligned, tighten the hardware receiver clamp screws (B). The recommended torque is 75-85 in-lbs (8.5-9.3 Nm). Then tighten the lower mounting pin nut (D). torque to 100-110 in-lbs 11.3 - 12.4 Nm
- Finish by inserting the locking pin (Fig.15) (H) into the receiver.

6. Function check, and micro adjustments to height, depth, and angle. 13 14 15 17

- Attach the J3 Back to the top and bottom mounting hardware receivers and push down until the assembly is locked. There should be a click when locked properly.
- If needed, continue small adjustments until the back can be attached, and/or removed smoothly, and locks correctly into the mounts.
- To make micro adjustments, loosen bolts (I) and (J). This will allow for minor changes in height, depth and angle. Retighten all bolts to required specs and test back function and quick-release.
- Or, to make macro adjustments, loosen, but do not remove, bolts (A), (C), and (G) (Figs. 13, 14, and 15). Position backrest to desired height, depth and angle.
- Retighten all bolts to required specs and test back function and quick-release.

D. ADJUSTING J3 HEIGHT, DEPTH, AND ANGLE 17 18 19

Once the initial installation is complete, adjustments can be made with the user seated in the chair for a personalized fit. Make sure that the user does not apply full to the back until the hardware has been fully tightened to the specified torque requirements.

1. Height adjustment

- With the user seated in the chair, identify the correct height for the back.
- Have the user slide forward slightly, and using a 10mm wrench Loosen the nuts (I).
- Move back to the correct height, and tighten the nuts to 75-85 in-lbs, or (8.5-9.3 Nm)

2. Depth and angle adjustment

Depth adjustment varies with the type of hardware ordered. Standard hardware including Rigid, Rigid Fixed, and Tilt hardware has approximately 2 inches of depth adjustment. Extended (adjustable) hardware offers up to 3.4 inches of depth adjustment.





- Using a 10mm wrench, loosen the left and right hex screws (J) on the bracket slots.
- Adjust the depth and/or angle settings and securely tighten the screws to lock in the adjustment. Recommended torque specification is 100-110 in-lbs or (11.3 -12.4 Nm).
- If additional depth is needed, remove the J3 back from the receivers. Remove the screws(K) from the receiver using a 4mm hex key and reposition the hardware to the desired depth. Standard hardware (K) has 2 depth settings, Extended (Adjustable) hardware (L) has 4 depth settings.
- After adjusting to the correct depth, tighten the screws to 25-35 in-lbs (2.8 - 4 Nm).

E. CHECK OUT

WARNING

- Ensure that all hardware has been properly tightened according to suggested specifications.
- Make sure the Back Support is locked, and the locking pin is inserted into the receiver.
- Check that the Back Support can be attached smoothly and securely, by attaching and removing the assembly several times.
- Double check fit and comfort for the user before putting the Back Support into use.
- Failure to follow this check out procedure could result in injury and/or damage to your Back Support and/or components.

A. VANITY FLAP INSTALLATION (Option) 20 21 22

Please read the following instructions before beginning the installation. To install the Vanity Flap it is best to begin without a user in the wheelchair. Once the Vanity flap has been installed the user can be seated back in the wheelchair.


1. Upper and lower attachment flaps 20

- a. Lift upper flap (A) of J3 Back Cover
- b. Attach the lower Flap (B) to the Backshell (hook and loop)
- c. Match hook to loop and attach Vanity Flap (C) to the back of the Lower flap of the J3 Back Cover.
- d. Reposition upper flap of cover to sandwich the Vanity Flap between the flaps the upper and lower flaps.

2. Vanity flap and cushion 21 22 23

- a. Remove Seat Cushion from wheelchair and sandwich the “tail” (D) of the Vanity Flap between seat sling and the Seat Cushion.
- b. Adjust the Vanity Flap (D) as necessary by moving the tail in or out from under the Seat Cushion. If necessary the upper part of the vanity flap can be moved up or down on the Seat Back also.
- c. At the rear of the chair, Insert Hook & Loop Wrap ties (Fig. 19) (E) through the appropriate grommets on the Vanity Flap, then wrap around the wheelchair frame and insert the leading end through the slot end of the Wrap. Snug up and close outside gaps in Vanity Flap.

B. SPINE ALIGN POSITIONING COMPONENTS (Option) 24 25 26

Spine Align Positioning Components are designed to provide intimately shaped support within an “off-the-shelf” back. Depending on your market, Spine Align components may be included with your back. Shown below are some typical components that can be used with your J3 Back Spine Align system.

Positioning Components	
Large Lumbar	
Small Lumbar	
Wedge	
Curved Shim	
Flat Shim	
Posterior/Lateral Pelvic	
Lateral Shim	
Large Bean Bag	
Small Bean Bag	



Finished J3 Vanity Flap Installation

I. Spine Align Positioning Component installation 25 26

- a. Have the user remove pressure from the J3 Back by leaning or sliding partially forward if possible.
- b. Release Velcro strip on the top cover flap (A) and lift the cover forward to reveal the inside of the J3 Back.
- c. Place the positioning component(s) (B) on the Velcro inside of the shell and reattach the top cover flap. Pull down making certain that the top seam is aligned with the top of the shell. If seam does not line up, the bottom of the cover may need to be loosened .
- d. Have the user test the current Spine Align setup for comfort and support and repeat steps a-c until adjustment is complete.
- e. Multiple configurations and combinations of Spine Align components can be used. If additional Spine Align shapes or sizes are required, please contact your dealer or Sunrise Medical at 800-333-4000.

NOTE— If you require further information on the J3 Spine Align system, please review our animation at www.sunrisemedical.com.

VIII. WARRANTY

24 MONTH LIMITED WARRANTY

1. Each Jay J3 Back is carefully inspected and tested to provide peak performance. Every Jay J3 Back is guaranteed to be free from defects in materials and workmanship for a period of 24 months from the date of purchase, provided normal use. Should a defect in materials or workmanship occur within 24 months from the original date of purchase, Sunrise Medical will, at its option, repair or replace it without charge. This warranty does not apply to punctures, tears or burns, nor to the back's removable cover.

2. Claims and repairs should be processed through the nearest authorized supplier. Except for express warranties made herein, all other warranties, including implied warranties of merchantability and warranties of fitness for a particular purpose are excluded. There are not warranties which extend beyond the description on the face hereof. Remedies for breach of express warranties herein are limited to repair or replacement of the goods. In no event shall damages for breach of any warranty include any consequential damages or exceed the cost of non-conforming goods sold.

Sunrise Medical S.r.l.
Via Riva, 20 – Montale
29122 Piacenza
Italia
Tel.: +39 0523 573111
www.SunriseMedical.it

Sunrise Medical AG
Erlenauweg 17
CH-3110 Münsingen
Schweiz/Suisse/Svizzera
Tel +41 (0)31 958 3838
www.SunriseMedical.ch

Sunrise Medical AS
Delitoppen 3
1540 Vestby
Norge
Telefon: +47 66 96 38 00
post@sunrisemedical.no
www.SunriseMedical.no

Sunrise Medical AB
Neogatan 5
431 53 Mölndal
Sweden
Tel.: +46 (0)31 748 37 00
post@sunrisemedical.se
www.SunriseMedical.se

MEDICCO s.r.o.
H – Park, Heršpická 1013/11d,
639 00 Brno
Czech Republic
Tel.: (+420) 547 250 955
Fax: (+420) 547 250 956
www.medicco.cz
info@medicco.cz
Bezplatná linka 800 900 809

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Sunrise Medical Australia
11 Daniel Street
Wetherill Park NSW 2164
Australia
Ph: +61 2 9678 6600
Email: enquiries@sunrisemedical.com.au
www.SunriseMedical.com.au

Sunrise Medical (US) LLC
North American Headquarters
12002 Volunteer Blvd.
Mount Juliet, TN 37122, USA
(800) 333-4000
(800) 300-7502
www.SunriseMedical.com

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© 2025 Sunrise Medical (US) LLC
Sunrise Medical B.V.
Sunrise Medical HCM B.V.
Sunrise Medical S.A.S
Sunrise Medical Canada Inc.