Jay Active™ Cushion

User Instruction Manual & Warranty

Jay Active™ pute

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JAY Active™ dynan

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INSTRUCTIONS FOR USE

Your Jay Active™ Cushion has been engineered to provide a high degree of pressure reduction, comfort, stability, and support. To ensure proper use of your cushion, please read these simple directions and save for future reference.

Caution

Installing a cushion on a wheelchair may affect the center of gravity of the wheelchair and may cause the wheelchair to tip backwards potentially resulting in injury. Always assess for the potential need for anti-tippers or amputee axle adapter brackets to be added to the wheelchair to help increase stability

Fitting the cushion

Select a cushion width to match your hip width when in the seated position. The cushion should be long enough to extend to 1" - 2" (3 - 5 cm) from back of knees. Position the cushion back as far as possible in the chair. The curved rear edge and notched corners should fit snugly against your chair back. To achieve the lowest sitting pressures possible, the footrests should be adjusted low enough to allow legs to rest comfortably on the cushion front. (Proper foot rest adjustment can decrease peak pressures.)

When sitting on the Jay Active Cushion, your ischials (seat bones) should be centered on the fluid portion of the cushion. Improper cushion position can cause bottoming out and the possibility of skin breakdown. If you slouch significantly, position the cushion slightly forward in the chair to avoid sitting on the front foam portion of the base. (If you slouch excessively, we recommend the J2® Cushion.) Similarly, if your back upholstery is stretched, it should be replaced. Consider replacing it with a J2 Back or Xtreme™ Back. Solid back support will help position the pelvis and support the trunk to reduce tendency to sacral sit.

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
Check for bottoming out

The Jay Active Cushion, when fit properly, is specifically designed to help prevent bottoming out under the ischials (seat bones) and coccyx (tail bone).


To check for bottoming out, sit on your Jay Active Cushion without its cover for two minutes. Lift yourself off the cushion (or have someone lift you). Examine the points where your ischials and coccyx were positioned. You should be able to push away at least 1/2" (1 cm) of fluid at these points before touching the foam base below. If the cushion is fit properly there should be at least 1/2" (1 cm) of fluid remaining at these points.

Anything less and the cushion is bottoming out and should not be used. (If you are newly injured, a 1" (2 cm) depth of fluid is recommended due to the possible occurrence of post-injury atrophy.) If you are bottoming out, you will need an overfilled pad or fluid supplement pads to add to your cushion. To arrange for an overfilled cushion or fluid supplement pads, call your local supplier or Sunrise Medical Customer Service at (303) 218-4500.

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
When used with solid seats or backs

The Jay Active Cushion is designed to be used on the wheelchair's upholstery sling. If used on a solid seat, the cushion should be modified to cut off the radius built into the bottom of the cushion. To arrange for this modification, call your local authorized supplier or Sunrise Medical Customer Service at (303) 218-4500. For improved posture and back pain relief, use of the J2 or Xtreme Back is highly recommended. Ask your authorized supplier for a free evaluation.

Adjustment of leg separation

To achieve the lowest sitting pressures and maximize leg support, lower the footrests, or raise the seat height, and/or use a longer length cushion. The leg wedges provided in the base cover help decrease leg separation. Also, you may tighten wheelchair upholstery or consider a solid seat insert. If the legs are still too far apart, evaluate a J2 Cushion.

Easy cleaning and assembly

Your wheelchair cushion is a very important part of your life. Keeping your cushion in good condition may ensure many years of use. The fluid in the Jay Active Cushion should require no regular maintenance.

However, monthly cleaning is recommended. During cleaning, it is recommended that you inspect your cushion. Check the cover for tears or excessive wear. Inspect the fluid area for punctures and fluid consistency. If you ever feel the fluid is firmer in one area, simply knead the fluid back to its original consistency and/or call Sunrise Medical.

Proper cleanliness will help extend the life of your cushion. Use the following instructions for cleaning.

Covers

Both the exterior cover and the foam base cover can be unzipped for machine washing and drying.

(Remove leg wedges from base cover.)

Foam base

The internal foam base should not be cleaned, as it will retain moisture and soap. This is why it should be protected by its cover at all times.

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To reassemble after cleaning

1. Insert foam base into incontinence-resistant cover with zipper to rear of base. Position with Velcro® attachment strips over the rear cutaway. Replace leg wedges into base cover, if desired.
2. Attach Jay Flow™ fluid pad to foam base cover at the five Velcro® attachment points. (Start at center rear attachment, then do the rear side attachments, and finally fasten the two front attachments.)
3. Insert the base with the Jay Flow™ fluid pad attached (fluid to the rear) into the cover and zip to close.

Caution

The Jay Active Cushion is designed to help prevent the incidence of pressure sores. However, no cushion can completely eliminate sitting pressure or prevent pressure sores. The Jay Active Cushion is not a substitute for good skin care including: proper diet, cleanliness, and regular pressure reliefs. In addition, prior to prolonged sitting, any cushion should be tried for a few hours at a time while a clinician inspects the skin to ensure that red pressure spots are not developing. Users should regularly check for skin redness. If your skin develops redness, discontinue the use of the cushion immediately and see your doctor or therapist. In addition, avoid sharp objects or exposure to excessive heat or open flame. If the cushion is left outside overnight at temperatures below 40°F (5°C), allow to warm to room temperature before using. Sitting on extremely cold or hot surfaces can cause skin damage.

Use Jay Active Covers only

Jay Active Covers are engineered for pressure relief and perspiration reduction. Use only Jay-designed covers, as the use of an improper cover will adversely affect the performance of the Jay Active Cushion.

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24 MONTH LIMITED WARRANTY

Each Jay Active Cushion is carefully inspected and tested to provide peak performance. Every Jay Active Cushion is guaranteed to be free from defects in materials and workmanship for a period of 24 months from the date of purchase, provided normal use. Should a defect in materials or workmanship occur within 24 months from the original date of purchase, Sunrise Medical will, at its option, repair or replace it without charge. This warranty does not apply to punctures, tears or burns, nor to the cushion's removable cover.

Claims and repairs should be processed through the nearest authorized supplier. Except for express warranties made herein, all other warranties, including implied warranties of merchantability and warranties of fitness for a particular purpose are excluded. There are not warranties which extend beyond the description on the face hereof. Remedies for breach of express warranties herein are limited to repair or replacement of the goods. In no event shall damages for breach of any warranty include any consequential damages or exceed the cost of non-conforming goods sold.

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OTHER JAY PRODUCTS**Xtreme Back™ by Jay**

The first solid back of its kind designed for the very active individual with good trunk control that utilizes a high performance rigid frame wheelchair.

The Xtreme Back is a dynamic high performance solid back that replaces the wheelchair's back upholstery. The ultra-light, low-profile design encourages freedom of movement, helps position and support the body comfortably, and once installed need never be removed from the wheelchair.

J2 Back

The first truly integrated back support for your wheelchair. The J2 Back is comfortable, adaptable, and lightweight. Lateral supports, angle adjustability, and simple one-handed removal and attachment make the J2 Back the most sophisticated back support available today. A back support has never been so important.

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U.K**SUNRISE MEDICAL LTD**

High Street
Wollaston
West Midlands
D48 4PS
England
44-1384-44-66-84

ENGLISH**SUNRISE MEDICAL**

Home Healthcare Group
Mobility Products Division
7477 East Dry Creek Parkway
Longmont, Colorado 80503
USA
(303) 218-4500

NORWEGIAN**SUNRISE MEDICAL**

Sunrise Home Healthcare Group
Mobility Products Division
7477 East Dry Creek Parkway
Longmont, Colorado 80503 USA
+1 (303) 218-4500 eller +1 (800) 648-8282
i Canada (800) 263-3390

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Norsk Rehab

1450 Nesoddtangen
Oslo, Norway 1450

 This product is manufactured to comply with the
'Medical Device Directive' 93/42/EEC.

Manufactured under one or more of the following patents:

U.S. Patent # 5,062,677, 5,524,971, 5,556,168 plus other patents pending.

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SWEDISH**Sunrise Medical AB**