


ORDER FORM




Q500M with Captain's Seat

EIPW30

| Customer | | Ship To | |
|-----------|--|----------|-------------|
| Account # | | Name | |
| Date | | Address | |
| PO# | | City | |
| Buyer | | Province | Postal Code |
| Mark For | | Phone | |

Client Measurements

| | |
|----------------------------|----|
| a) Top of head to seat | in |
| b) Top of shoulder to seat | in |
| c) Underarm to seat | in |
| d) Elbow to seat | in |
| e) Top of leg to seat | in |
| f) Shoulder width | in |
| g) Chest width | in |
| i) Hip width | in |
| j) Actual sitting depth | in |
| m) Foot length | in |
| n) Knee to foot | in |
| o) Seat width | in |



Notes

Chair Specifications

| | |
|--|---|
| Weights starting at: (with Batteries) | NCC: 277 lbs SCC: 301 lbs |
| Power Seating | 50° CG Tilt / 12" seat elevator |
| Seat Type | Captain's Seat w/ 25° Adjustable Recline |
| Seat Size Range | Width 16" to 20" Depth 16" to 22" |
| Seat-to-Floor Height Range | 17" to 19" |

Mid-Wheel Drive Base Specifications

| | |
|------------------------|-------------------------------------|
| Drive Wheel Position | Mid-wheel drive |
| User Weight Capacity | 300 lbs. |
| Base Size | Width 24", Length 33" |
| Minimum Turning Radius | 18.5" (excludes footplates) |
| Drive Control | R-net PM 90A / VR-2 90A |
| Driving Range | Up to 12 miles |
| Drive Motors | 4-pole 6 mph / 10 kph |
| Drive Wheels | 13" w/flat-free inserts (black) |
| Casters | 7" front / 7" rear (black) |
| Charger | 8A Dual Mode (Compact) |
| Battery Size | Group 34 / Group 22 |
| Color | 7 shroud colors and 7 color accents |
| Ground Clearance | 3" |

BASE

| OPTION | DESCRIPTION | RETAIL |
|--------|-------------|--------|
|--------|-------------|--------|

Standard Duty Mid Wheel Drive Base *End Item Number EIPW30 (Captain's Seat)* \$

| | | |
|----------|---|------------|
| 300PC218 | Quickie Q500 M NCC No Power Captain Seat and Back Included in Base Pricing : Q500 M Base, VR-2 90amp Electronics for No Power Seat Functions, 6-wheel Suspension System, 13" Drive Wheels with Dual Density Inserts, 7" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 17". | \$8,274.00 |
| 300PC219 | Quickie Q500 M SCC Single Power Option Captains Seat and Back Included in Base Pricing : Q500 M Base, VR-2 90amp Electronics for No Power Seat Functions, 6-wheel Suspension System, 13" Drive Wheels with Dual Density Inserts, 7" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 17". | \$8,274.00 |
| 300PC220 | Quickie Q500 M NCC No Power Captains Seat Back and Solid Seat Included in Base Pricing : Q500 M Base, VR-2 90amp Electronics for No Power Seat Functions, 6-wheel Suspension System, 13" Drive Wheels with Dual Density Inserts, 7" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 17". | \$8,274.00 |
| 300PC221 | Quickie Q500 M SCC Single Power Option Captains Seat Back and Solid Seat Included in Base Pricing : Q500 M Base, VR-2 90amp Electronics for No Power Seat Functions, 6-wheel Suspension System, 13" Drive Wheels with Dual Density Inserts, 7" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 17". | \$8,274.00 |

User Weight \$

| | | |
|---------|----------------|----|
| 300WG23 | Up to 300 lbs. | NC |
|---------|----------------|----|

Motor Package \$

| | | |
|---------|-----------------------|----|
| 300PF99 | 4-pole 10 kph / 6 mph | NC |
|---------|-----------------------|----|

Drive Wheels  P. 10 \$

| | | |
|---------|---|----|
| 300DW69 | Drive Wheel, 13", Dual Density, Knobby, Black | NC |
|---------|---|----|

Caster Wheels  P. 10 \$

| | | |
|---------|--------------------------------|----|
| 300CW59 | Caster Wheels 7", Solid, Black | NC |
|---------|--------------------------------|----|

Shroud Color \$

| | | |
|---------|-----------------|----|
| 300C119 | Matte Black | NC |
| 300C165 | Candy Apple Red | NC |
| 300C166 | Candy Blue | NC |
| 300C167 | Candy Purple | NC |
| 300C176 | Green Apple | NC |
| 300C178 | Pearl Pink | NC |
| 300C182 | Silver | NC |

| BASE | | | | | |
|---|--|---|-------------------------|----------------------------|------------------|
| OPTION | | DESCRIPTION | | | RETAIL |
| Color Accents | | | | | \$ |
| <i>Note: Accent colors and shroud colors may not match.</i> | | | | | |
| 300EZ20 | | Apple Green | | | NC |
| 300EZ22 | | Black | | | NC |
| 300EZ21 | | Blue | | | NC |
| 300EZ23 | | Orange | | | NC |
| 300EZ24 | | Red | | | NC |
| 300EZ25 | | Silver | | | NC |
| 300EZ19 | | White | | | NC |
| ISO Transit Tie Downs | | | | | \$ |
| 300T050 | | ISO Transit Package | | | NC |
| LED Lights  P. 10 | | | | | \$ |
| 300LI04 | | LED Lights and Indicators | | | \$1,207.00 |
| 300LI08 | | LED Fender Lights <i>(Single Switch Activation)</i> | | | \$636.00 |
| 300LI02 | | Omit LED Lights | | | NC |
| Seat Heights | | | | | \$ |
| | | Seat Height | Power Seating with Tilt | Power Seating with Elevate | No Power Seating |
| 300HT55 | | 17" | N/A | N/A | Group 34 |
| 300HT104 | | 17.5" | N/A | N/A | Group 34 |
| 300HT56 | | 18" | Group 34 | N/A | Group 34 |
| 300HT105 | | 18.5" | Group 34 | Group 34 | Group 34 |
| 300HT57 | | 19" | Group 34 | Group 34 | Group 34 |
| Batteries and Chargers  P. 10 | | | | | \$ |
| Group 34 Battery Options | | | | | |
| 300AG27 | | GP34 Flat Top Pair w/Straps MK Gel | | | \$941.00 |
| 300AG29 | | GP34 Flat Top Pair w/Straps Brightway AGM | | | \$850.00 |
| Group 22 Battery | | | | | |
| 300AG31 | | GP22 Flat Top Pair w/Straps | | | \$668.00 |
| Omit Batteries | | | | | |
| 300AG24 | | Battery Harnesses Only <i>(Omit Batteries)</i> | | | NC |
| Chargers | | | | | |
| 300BJ30 | | 10A Dual Mode (Compact) | | | \$160.00 |
| 300BJ19 | | 8A Dual Mode (Compact) | | | NC |

SEATING

| OPTION | DESCRIPTION | RETAIL |
|--------|-------------|--------|
|--------|-------------|--------|

Power Seating Options  P. 10 \$

Seat Elevator

| | | |
|---------|--------------------------------------|------------|
| 300PI07 | Seat Elevator, Powered, 300 mm / 12" | \$3,426.00 |
|---------|--------------------------------------|------------|

Seat Tilt

| | | |
|---------|--------------------|------------|
| 300PA48 | Posterior 50° Tilt | \$8,535.00 |
|---------|--------------------|------------|

Pre Set Seat Frame Angle (No Power Seating Options, Increases Seat to Floor Height) \$

| | | |
|---------|----------------------|----|
| 300FA10 | 0° (Adds 0" to STFH) | NC |
| 300FA3 | 3° (Adds 1" to STFH) | NC |

Backrest Frame Type \$

| | | |
|---------|--|----|
| 300BF08 | Captain Seat Frame (Includes Back Cushion) | NC |
|---------|--|----|

Seat Cushion Type \$

| | | |
|---------|-----------------------------------|----|
| 300ST77 | Captain's Cushion (See Section A) | NC |
| 300ST80 | JAY® (See Section B) | NC |
| 300ST75 | Omit | NC |

A - Seat Size \$

| Seat Pan Width | | Available Depths | | |
|----------------|-----|------------------|--|----|
| 300R53 | 16" | 16" to 18" | | NC |
| 300R55 | 18" | 16" to 20" | | NC |
| 300R57 | 20" | 18" to 22" | | NC |

Seat Pan Depth (Front of Seat Pan to Back Structure. Back Cushions Can Cause 2" Loss of Effective Depth)

| | | |
|---------|-----|----|
| 300RF51 | 16" | NC |
| 300RF45 | 17" | NC |
| 300RF46 | 18" | NC |
| 300RF47 | 19" | NC |
| 300RF48 | 20" | NC |
| 300RF49 | 21" | NC |
| 300RF50 | 22" | NC |

SEATING

| OPTION | DESCRIPTION | RETAIL |
|--------|-------------|--------|
|--------|-------------|--------|

B - JAY® Seat Cushion  P. 10 \$

| | | |
|----------|---|------------|
| 300CU59 | JAY® Ion® (Sizes 16" - 20" Wide) | \$498.00 |
| 300CU61 | JAY® Union® (Sizes 16" - 20" Wide) | \$753.00 |
| 300CU64 | JAY® Fusion® (Sizes 16" - 20" Wide) | \$1,030.00 |
| 300CU104 | JAY® Fusion® with Cryo™ Technology (Sizes 16" - 20" Wide) | \$1,135.00 |

POSITIONING SUPPORTS & ACCESSORIES

| OPTION | DESCRIPTION | RETAIL |
|---|---|------------|
| Armrest Structure \$ | | |
| 300EG37 | Flip Back <i>(Includes Waterfall Arm Pads) (Height Adjustable 11" to 14")</i> | NC |
| Fixed Centermount Legrest \$ | | |
| 300H183 | SEDEO Fixed Centermount Legrest <i>(3" to 17.5" Leg Length)</i> | \$346.00 |
| 300H184 | Fixed Centermount Legrest, Short <i>(Includes Calf Pads and Single Footplate) (11.5" to 14.5" Leg Length)</i> | NC |
| 300H185 | Fixed Centermount Legrest, Long <i>(Includes Calf Pads and Single Footplate) (13.5" to 18.5" Leg Length)</i> | NC |
| Power Centermount Legrest  P. 10 \$ | | |
| 300H175 | Power Centermount Legrest, 90° - 0° Angle, 7.5" Extension <i>(Mechanical Articulation) (5" to 17.5" Leg Length)</i> | \$2,450.00 |
| 300H109 | Power Centermount Legrest, 90° - 0° Angle, 7.5" Extension <i>(Powered Articulation) (3" to 17.5" Leg Length)</i> | \$3,384.00 |
| Centermount Legrest Accessories \$ | | |
| Footplates | | |
| 300F127 | Dual, Height, Width and Angle Adjustable, Standard (14"W x 10"D) | \$280.00 |
| 300F128 | Dual, Height, Width and Angle Adjustable, Large (16"W x 12"D) | \$280.00 |
| Calf Pad Options | | |
| 300DI34 | Calf Pads, Pair, Standard Integral Foam | NC |
| 300DI35 | Calf Pads, Pair, Padded <i>(N/A w/ Std footplate)</i> | \$351.00 |
| Pelvic Positioning Belts  P. 10 | | |
| 300A837 | Positioning Belt, Auto Style Buckle, 50 mm 1600 mm / 2" | NC |
| 300A839 | Positioning Belt, Aircraft Buckle, 50 mm x 1300 mm / 2" | NC |
| 300A840 | Positioning Belt, Aircraft Buckle, Padded 50 mm x 1300 mm / 2" | \$101.00 |
| 300A845 | Bodypoint – Evoflex Pelvic Stabilizer Belt | \$213.00 |
| 300A846 | Bodypoint – Evoflex Pelvic Stabilizer Belt with Extension | \$267.00 |
| Accessories \$ | | |
| O2 Tank Option | | |
| 300A718 | O2 Tank Holder Mounted to Backrest <i>(D and E Styles)</i> | \$385.00 |
| Captain's Accessories | | |
| 300A608 | Cane / Crutch Holder | \$113.00 |
| 300A610 | Walker Holder | \$188.00 |
| Misc. Accessories | | |
| 3008A52 | Ram Mount Phone | \$76.00 |
| 3008A53 | Ram Mount Tablet | \$303.00 |
| 3008A54 | Ram Mount Cup Holder | \$76.00 |
| 3009A10 | Utility Hooks | \$57.00 |
| Medial Thigh Support | | |
| 3008A80 | Medial Thigh Support | \$267.00 |
| 3008A81 | Pommel Swing-Away Hardware | \$384.00 |

JOYSTICK / ELECTRONICS

| OPTION | DESCRIPTION | RETAIL |
|--------|-------------|--------|
|--------|-------------|--------|

VR-2 Joystick \$

| | | |
|---------|---------------|----|
| 300J127 | VR-2 Joystick | NC |
|---------|---------------|----|

R-net Joystick (Remote)  P. 11 \$

| | | |
|---------|---|------------|
| 300J121 | Standard LED Joystick (w/o Screen Display) | NC |
| 300J116 | Color Joystick (2.3" Color LCD Screen Display, CJSM1) | NC |
| 300J118 | High Resolution Color Joystick, (CJSM2+BT) 3.5" Color LCD Display, Paddle Switches, Programmable Shortcut Keys, ICONs, Stereo Jack, IR and dual (IOS/Android) Bluetooth control | \$1,245.00 |
| 300J119 | Omit Remote Control | NC |

Joystick Mounting \$

| | | |
|---------|--|----------|
| 300JM34 | Fixed on Armrest | NC |
| 300JM35 | Parallel Swing-Away | \$440.00 |
| 300JM36 | Link-It Mid Line Mount on Armrest (Swing-Away) | \$440.00 |

Joystick Mount (Location) \$

| | | |
|---------|--------------------|----|
| 300HM16 | Right Hand Armrest | NC |
| 300HM17 | Left Hand Armrest | NC |

Display Language Selection R-net \$

| | | |
|-------|---------|----|
| 300L1 | English | NC |
| 300L3 | French | NC |
| 300L7 | Spanish | NC |

Joystick Handle (Optional)  P. 11 \$

| | | |
|---------|--|----------|
| 300JH13 | Ball Handle | \$108.00 |
| 300JH15 | T-Handle | \$108.00 |
| 300JH24 | Bodypoint – Rubber Dome Joystick Handle | \$213.00 |
| 300JH25 | Bodypoint – U-Shaped Joystick Handle (4 in) | \$213.00 |
| 300JH26 | Bodypoint – Grooved Mushroom Joystick Handle | \$213.00 |

Electronic Programmers  P. 11 \$

| | | |
|---------|-------------------------------|----------|
| 300PR10 | R-net PC Programmer w/ Dongle | \$696.00 |
|---------|-------------------------------|----------|

R-net Power Modules  P. 11 \$

| | | |
|---------|---|----|
| 300ER12 | Power Module, VR-2 90A, Non-Expandable | NC |
| 300ER22 | Power Module, R-net, 90A, Non-Expandable (Only available with single power) | NC |
| 300ER13 | Power Module, R-net, 90A, Expandable | NC |
| 300EX11 | Expandable 90A Controller - R-net | NC |
| 300EX12 | Harness for 90A Expandable Controller | NC |

JOYSTICK / ELECTRONICS

| OPTION | DESCRIPTION | RETAIL |
|--------|-------------|--------|
|--------|-------------|--------|

Powered Seating Control

| | | | |
|---------|--|--|------------|
| 300E040 | | Omit Thru-Drive <i>(Requires Switch Control)</i> | NC |
| 300E037 | | Thru-drive Control for 1 Actuator | \$1,422.00 |
| 300E038 | | Thru-drive Control for 2+ Actuators | \$2,860.00 |

Switch Control *(Available with Omit Thru-Drive only) (1 actuator for Tilt or Elevate only)*

| | | | |
|---------|--|---|----|
| 300E033 | | Dual Toggle for 1 Actuator | NC |
| 300E034 | | 2.5" Single Button Switch for 1 Actuator <i>(Button switches are not mounted)</i> | NC |
| 300E041 | | Ctrl+5 Switch Box for 2 Actuators or More | NC |

Additional Power Seating Control *(Available with Thru-Drive only)*

| | | | |
|---------|--|--|------------|
| 300QG11 | | Ctrl+5 Switch Box with Free Assignable Functions | \$1,062.00 |
|---------|--|--|------------|

Dual Toggle Mounting Location

| | | | |
|--------|--|---------------------|----|
| 300TN1 | | Toggle Mount, Right | NC |
| 300TN2 | | Toggle Mount, Left | NC |

Ctrl+5 Mounting Location

| | | | |
|---------|--|---------------|----|
| 300EY19 | | Right Armrest | NC |
| 300EY20 | | Left Armrest | NC |

Ctrl+5 Switch Type

| | | | |
|---------|--|-----------------|----------|
| 300EY22 | | Push Buttons | NC |
| 300EY23 | | Toggle Switches | \$351.00 |

R-net Programming**Assignable Buttons**

Assignable Button feature allows for personalization of joystick buttons and/or Ctrl+5 Switches allowing access to various wheelchair and seating functions. When enabled, chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle.

| | | | |
|---------|--|----------------------------|----|
| 300PM12 | | Include Assignable Buttons | NC |
|---------|--|----------------------------|----|

Power Centermount Footplate Extension Option

| | | | |
|---------|--|--|----------|
| 300IE03 | | Include 7.5" Powered Footplate Extension for Transfer Assist <i>(Electronics programmed to allow footplates to be extended with or without legrest being raised.) Only Available with Powered Articulation</i> | \$583.00 |
|---------|--|--|----------|



Q500 M with Captain's Seat

\$

Base Price

Total of additional equipment

Packaging and transport

Total Price

Date

Signature

Month

Day

Year



Customer Service: (800) 263-3390
Fax: (800) 561-5834
E-mail: cscanada@sunmed.com
www.sunrisemedical.ca
Specifications are subject to change without notice.

© 2024 Sunrise Medical Canada Inc.
237 Romina Drive, Unit 3
Concord, ON
L4K 4V3

MK-129645-CAN Rev. 11