

February 2024

**Customer Contact Details**

| | |
|----------------------|--|
| Customer Name: _____ | Ship to Customer Address: <input type="checkbox"/> |
| Advisor: _____ | Shipping Name: _____ |
| Address: _____ | Shipping Address: _____ |
| _____ | _____ |
| _____ | _____ |
| Phone: _____ | Phone: _____ |
| Email: _____ | Email: _____ |

Additional Information

| | |
|--------------|----------------------------|
| Notes: _____ | Current Wheelchair: _____ |
| _____ | Customer Weight: _____ |
| _____ | Customer Height: _____ |
| _____ | Customer Disability: _____ |

| Guide: | | <input type="checkbox"/> Base product | <input checked="" type="checkbox"/> Required option | <input type="checkbox"/> Optional Extra | <input type="checkbox"/> Critical Information | <input type="checkbox"/> Information Needed |
|---|---|---|---|---|---|---|
| <input checked="" type="checkbox"/> SQCXZ010 | CLUB SPORT \$3,282.00 | Max User Weight 125kg (275lbs) | | | <input type="checkbox"/> Yes, I am within the max user weight | |
| Standard Specification Includes: | | | | | | |
| <ul style="list-style-type: none"> <li style="width: 33%;">Basketball or Tennis (option 3 only) <li style="width: 33%;">Removable single sport anti-tip <li style="width: 33%;">72mm caster wheels with extruded forks <li style="width: 33%;">Aluminium screw on sideguards <li style="width: 33%;">Adjustable height backrest 250-300mm <li style="width: 33%;">Adjustable seat and backrest upholstery <li style="width: 33%;">Adjustable position footrest <li style="width: 33%;">Cross laced wire spoke wheels <li style="width: 33%;">5cm medium density foam cushion (*option 4) <li style="width: 33%;">Adjustable centre of gravity (18 degree) <li style="width: 33%;">Aluminium hard anodized push rim/kenda tire <li style="width: 33%;">Calf strap & foot strap | | | | | | |
| Key Product Measurements | | | | | | |
| Frame Option (Choose one) | | | | | | |
| <input type="checkbox"/> XY520 Option 1 <i>(Junior size)</i> | <input type="checkbox"/> XY530 Option 2 <i>(Low Point)</i> | <input type="checkbox"/> XY540 Option 3 <i>(Low/Mid Point)</i> | <input type="checkbox"/> XY550 Option 4 <i>(Max height)</i> | | | |
| Sport | | | | | | |
| <input checked="" type="checkbox"/> Basketball <i>(with wing)</i> | <input checked="" type="checkbox"/> Basketball <i>(with wing)</i> | <input type="checkbox"/> Basketball <i>(with wing)</i> <input type="checkbox"/> Tennis <i>(no wing)</i> <small>(5mm A/T) (50mm A/T)</small> | <input checked="" type="checkbox"/> Basketball <i>(with wing)</i> | | | |
| Frame height | | | | | | |
| <input checked="" type="checkbox"/> Front frame height: 480mm <input checked="" type="checkbox"/> Rear frame height: 420mm | <input checked="" type="checkbox"/> Front frame height: 500mm <input checked="" type="checkbox"/> Rear frame height: 400mm | <input checked="" type="checkbox"/> Front frame height: 500mm <input checked="" type="checkbox"/> Rear frame height: 440mm | <input checked="" type="checkbox"/> Front frame height: 530mm <input checked="" type="checkbox"/> Rear frame height: 530mm | | | |
| | | | *Aluminium seat plate & 5cm Plastazote | | | |
| Wheel size | | | | | | |
| <input type="checkbox"/> 22" XY415 <input type="checkbox"/> 24" XY420 <input type="checkbox"/> 25" XY430 | <input type="checkbox"/> 24" XY420 <input type="checkbox"/> 25" XY430 | <input type="checkbox"/> 24" XY420 <input type="checkbox"/> 25" XY430 | <input type="checkbox"/> 26" XY440 <input type="checkbox"/> 700c XY450 | | | |
| Seat width (inside of sides) | | | | | | |
| <input type="checkbox"/> 340mm <input type="checkbox"/> 370mm | <input type="checkbox"/> 370mm <input type="checkbox"/> 400mm <input type="checkbox"/> 440mm | <input type="checkbox"/> 370mm <input type="checkbox"/> 400mm <input type="checkbox"/> 440mm | <input type="checkbox"/> 370mm <input type="checkbox"/> 400mm <input type="checkbox"/> 440mm | | | |
| Seat length | | | | | | |
| <input type="checkbox"/> 300mm <input type="checkbox"/> 330mm <input type="checkbox"/> 360mm | <input type="checkbox"/> 360mm <input type="checkbox"/> 390mm <input type="checkbox"/> 430mm | <input type="checkbox"/> 360mm <input type="checkbox"/> 390mm <input type="checkbox"/> 430mm | <input type="checkbox"/> 360mm <input type="checkbox"/> 390mm <input type="checkbox"/> 430mm | | | |
| Frame length | | | | | | |
| <input checked="" type="checkbox"/> 520mm | <input checked="" type="checkbox"/> 580mm | <input checked="" type="checkbox"/> 580mm | <input checked="" type="checkbox"/> 580mm (SW = 370/400) <input checked="" type="checkbox"/> 630mm (SW = 440) | | | |
| Options | | | | | | |
| Finish Options | | | Accessories | | | |
| <input checked="" type="checkbox"/> XN010 Powder coated frame | <input type="checkbox"/> Specify Code <input style="width: 100px;" type="text"/> | | <input type="checkbox"/> XS020 Sport ratchet strap (click strap) \$138.00 <i>(not fitted)</i> | <input type="checkbox"/> XS010 100mm velcro adjustable strap \$85.00 | | |
| <input checked="" type="checkbox"/> XN150 Powder coated components <i>(Black)</i> | | | <input type="checkbox"/> XK080 Spinergy SLX 24 spoke \$991.00 <i>(Wheel upgrade)</i> | <input type="checkbox"/> XP170 Individual backrest embroidery \$236.00 <i>(Please provide .jpeg file and positioning)</i> | | |