



MWD Power Wheelchair.

EIPW60



Customer				Ship To			
Account				Name			
Account #		Date		Address			
☐ Quote	☐ Order	PO #		City		Prov	
Buyer				PC		Phone	
Mark For				Email / Fax			

Please submit orders to: cscanada@sunmed.com

Seating Specifications

Typical WC Weight w/ Batteries	Starting at 320 lbs. Weights vary depending on options		
Power Seating	50° of CG Tilt / Power Recline / 12" Seat Elevator	Seat Type	MPS and Rehab Style
Seat Size Range	Width: 13-23" Depth: 13-24"	STF Height Range	17", 18", 19"


Drive Base Specifications

Drive Wheel Position	Mid Wheel Drive	Drive Control	R-Net PM 120A controller
User Weight Capacity	300 lbs.	Drive Motors	4-pole 6mph
V6 Base Size: All Terrain	Width 28" Length 41.4" Turning Radius 23.5"	Drive Wheels	14" All Terrain Drive Wheels 14" Hybrid Drive Wheels
V6 Base Size: Hybrid	Width 27.3" Length 41.5" Turning Radius 23.5"	Tires / Casters	Pneumatic or Solid inserts
		Charger	Dual Mode
		Battery Size	Group 24 or 34
		Color Accents	9 colors
		Ground Clearance	3."



Client Measurements

a)	Top of head to seat	in
b)	Top of shoulder to seat	in
c)	Underarm to seat	in
d)	Elbow to seat	in
e)	Top of leg to seat	in
f)	Shoulder width	in
g)	Chest width	in
i)	Hip width	in
j)	Actual sitting depth	in
m)	Foot length	in
n)	Knee to foot	in
o)	Seat width	in



Notes

✓ OPTION	DESCRIPTION	RETAIL
SEATING OPTIONS		\$
☐ 600SF4	REHAB Seat and Back - Back Canes and adjustable flat seat pan	NC
☐ 600SF5	MPS Seat and Back - Captain style seating	NC
☐ 600SF6	MPS Back and REHAB Seat - Captain style backrest and adjustable flat seat pan	NC
BASE		
Magic 360 Base - Includes Transit and lights		End Item Number: EIMM
		\$
☐ 600MW61	Magic Mobility - 360 No Power Seating Options	\$ 26,350.00
☐ 600MW62	Magic Mobility - 360 Single Power Seating Options	\$ 26,350.00
☐ 600MW63	Magic Mobility - 360 Multi Power Seating Options	\$ 26,350.00
☐ 600MW64	Magic Mobility - 360 No Power Seating Options - MPS (Captain)	\$ 26,350.00
☐ 600MW65	Magic Mobility - 360 Single Power Seating Options - MPS (Captain)	\$ 26,350.00
☐ 600MW66	Magic Mobility - 360 Multi Power Seating Options - MPS (Captain)	\$ 26,350.00
Drive Wheels		\$
☐ 600DW80	All Terrain Drive Wheels and Tires	NC
☐ 600DW81	Hybrid Drive Wheels and Tires - Pneumatic	NC
☐ 600DW82	Hybrid Drive Wheels and Tires - Solid	\$ 62.00
Optional Drive Tire Options		\$
☐ 6008A90	Tire Sealant for 2 Pneumatic Drive Tires	\$ 50.00
Base Accessories		\$
☐ 6008A91	Jack for Drive Tire	\$ 268.00
Frame Color		\$
☐ 600C185	Midnight Black	NC
☐ 600C186	Lucent Lime	NC
☐ 600C187	Sunset Red	NC
☐ 600C188	Horizon Blue	NC
☐ 600C189	Dusk Purple	NC
☐ 600C190	Vivid Pink	NC
☐ 600C191	Desert Orange	NC
☐ 600C192	Camouflage - Woods	NC
☐ 600C193	Camouflage - Desert	NC
☐ 600C171	Carbon Fiber	NC
☐ 600C194	American Flag	NC

Seat To Floor Height \$

360 Seat to Floor Heights start at 17" +/- 0.25"

☐ 600HT44		Lowest Seat To Floor Height	NC
☐ 600HT45		Seat To Floor Height +1"	NC
☐ 600HT46		Seat To Floor Height +2"	NC

Batteries \$

☐ 600AG25		Group 24 Gel Battery (Pair) - Installed	\$ 1,184.00
----------------------------------	--	---	-------------

Battery Charger \$

☐ 6009A14		10A Battery Charger	NC
----------------------------------	--	---------------------	----

Power Retractable Docking Pin \$

☐ 600BJ28		Retractable Docking Pin - Compatible with a Q-Straint QLK 150 docking system	\$ 4,100.00
----------------------------------	--	--	-------------

SEATING Seat Functions

✓ OPTION		DESCRIPTION	RETAIL
----------	--	-------------	--------

No Power Seat Options \$

☐ 600ST87		No Power Seating System with Power Functionality	\$ 9,190.00
----------------------------------	--	--	-------------

Power Tilt \$

Utilized without Power Recline.

☐ 600PA64		50° Power Tilt in Space - 340 lb. weight capacity	\$ 10,085.00
----------------------------------	--	---	--------------

Power Seat Elevate \$

Only available when Power Tilt or Power Tilt and Recline is selected.

☐ 600PI16		12" Adjustable Seat Elevator - 340 lb. weight capacity	\$ 8,785.00
----------------------------------	--	--	-------------

Power Recline \$

Utilized without Power Tilt. Requires No Power Seating with Power Functionality.

☐ 600PA65		150° Power Recline Only - 300 lb. weight capacity	\$ 7,178.00
----------------------------------	--	---	-------------

Power Tilt and Recline \$

☐ 600PA66		50° Tilt & 150° Power Recline - 300 lb. weight capacity	\$ 15,926.00
----------------------------------	--	---	--------------

SEATING: SEAT SIZES

✓ OPTION		DESCRIPTION	RETAIL
----------	--	-------------	--------

MPS (Captain) Seat Cushion (width x depth) \$

MPS Backrest can be paired with the REHAB flat pan. If flat pan is desired, omit the MPS cushion and select the desired REHAB flat pan.

☐ 600SI24		16" x 18" Vinyl Seat	NC
☐ 600SI25		18" x 18" Vinyl Seat	NC
☐ 600SI26		18" x 20" Vinyl Seat	NC
☐ 600SI27		20" x 18" Vinyl Seat	NC
☐ 600SI28		20" x 20" Vinyl Seat	NC
☐ 600SI29		22" x 20" Vinyl Seat	\$ 608.00

REHAB Flat Pan Seat Sizes - Adjustable in width and depth for most sizes \$

REHAB Seat Pan

Standard seat pan available 13-23" wide. Deep seat pan available 13"-23" wide.

☐ 600PN13		Standard Adjustable Flat Seat Pan - 13" to 19" deep	NC
☐ 600PN14		Deep Adjustable Flat Seat Pan - 20" to 22" deep	\$ 612.00
☐ 600PN15		Deep Adjustable Flat Seat Pan - 23" to 24" deep	\$ 612.00

REHAB Seat Width

☐ 13"	☐ 14"	☐ 15"	☐ 16"	☐ 17"	☐ 18"	☐ 19"	☐ 20"	☐ 21"	☐ 22"	☐ 23"
------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------

REHAB Seat Depth

☐ 13"	☐ 14"	☐ 15"	☐ 16"	☐ 17"	☐ 18"	☐ 19"	☐ 20"	☐ 21"	☐ 22"	☐ 23"	☐ 24"
------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------

REHAB Seating Cushion Options - \$

Note cushion size will be determined by Rehab Seat size selected above

☐ 600CU59		JAY® Ion® (Sizes 16" - 21" Wide)	\$ 647.00
☐ 600CU60		JAY® Ion® (Size 22" Wide)	\$ 688.00
☐ 600CU61		JAY® Union® (Sizes 16" - 21" Wide)	\$ 753.00
☐ 600CU62		JAY® Union® (Size 22" Wide)	\$ 794.00
☐ 600CU64		JAY® Fusion® (Sizes 16" - 21" Wide)	\$ 1,030.00
☐ 600CU65		JAY® Fusion® (Size 22" Wide)	\$ 1,170.00

Lateral Thigh Support Pads and Hardware (must pick pad an \$

Lateral Thigh Support Pads

☐ 600KF22		Lateral Thigh Support Pad Short - Right (5"x5")	\$ 74.00
☐ 600KF23		Lateral Thigh Support Pad Short - Left (5"x5")	\$ 74.00

Lateral Thigh Support Mounts

☐ 600KF26		Lateral Thigh Removable Mount - Right	\$ 279.00
☐ 600KF27		Lateral Thigh Removable Mount - Left	\$ 320.00

SEATING: SEATBACK OPTIONS

MPS Backrest Options - Comes standard with Manual Recline \$

Backrest Width

☐ 600BN46		Standard Backrest Width - 16-20" Wide	NC
☐ 600BN47		Wide Backrest Width - 22-23" Wide	\$ 445.00

Backrest Height

☐ 600B80		Standard Backrest Height - 22" Tall	NC
☐ 600B81		Tall Backrest Height - 26" Tall	\$ 733.00

Backrest Material

☐ 600BU106		Cloth Backrest	NC
☐ 600BU107		Vinyl Backrest	\$ 404.00

MPS Back Accessories \$

☐ 600PH33		Push Handle - Needed for O2 holder, backpack, and cane holder	\$ 208.00
☐ 6009A13		Accessory Mount - Needed for aftermarket headrest and/or attendant control	\$ 157.00

REHAB Seating Backrest Options

\$

Backrest will have the width of the MPS seat or flat pan. Please choose your back height. Backrest will come with 2 point hardware with a 500 lb. weight capacity.

JAY J3 Back Contour

☐ 600JP01	J3@ PA Back for 16"-20" Wide - 2.2" of lateral thoracic support	\$ 899.00
☐ 600JP02	J3@ PA Back for 22" Wide - 2.2" of lateral thoracic support	\$ 1,912.00
☐ 600JP11	J3@ PL Back for 16"-20" Wide - 3" of lateral thoracic support	\$ 971.00
☐ 600JP13	J3@ PD Back for 16"-20" Wide (18" and 21" shell height only) - 6" of lateral thoracic support	\$ 1,052.00

JAY® J3® Back Shell Options and Height Ranges

☐ 600JR14	16.5" tall shell (Upper Thoracic Short)	NC
☐ 600JR15	18" tall shell (Upper Thoracic Med)	NC
☐ 600JR16	19.5" tall shell (Upper Thoracic Tall)	NC
☐ 600JR18	18" tall shell (Shoulder Height) for 22" wide	NC
☐ 600JR17	21" tall shell (Shoulder Height)	NC
☐ 600JR27	22.5" tall shell (Shoulder Height)	NC
☐ 600JR28	24" tall shell (Shoulder Height)	NC

REHAB Flip Forward Backrest Option

\$

☐ 600D23	REHAB Flip Forward Backrest Option	#N/A
---------------------------------	------------------------------------	------

JAY® Lateral Thoracic Swing-Away Supports (Requires Pad and Mount) - J3 back only

\$

JAY® Lateral Thoracic Support Pads

☐ 600JW35	Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Right	\$ 74.00
☐ 600JW38	Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Left	\$ 74.00
☐ 600JW37	Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Right	\$ 74.00
☐ 600JW40	Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Left	\$ 74.00

JAY® Lateral Thoracic Swing-Away Mount

☐ 600JW41	Lateral Thoracic Swing Away Bracket - Right	\$ 329.00
☐ 600JW42	Lateral Thoracic Swing Away Bracket - Left	\$ 329.00

JAY® Anterior Thoracic Support

\$

☐ 600ZR10	Bodypoint Monoflex Medium Chest Strap - Pad 3.5" x 18.5" (total length 49.5")	\$ 224.00
☐ 600ZR11	Bodypoint Monoflex Large Chest Strap - Pad 3.5" x 21.5" (total length 52.5")	\$ 224.00
☐ 600K10	JAY® Anterior Thoracic Support - Medium ATS8302	\$ 317.00
☐ 600K11	JAY® Anterior Thoracic Support - Large ATS8303	\$ 317.00

Pelvic Positioning Options - Select One

\$

Pelvic Positioning

☐ 600A848	Positioning Belt, Auto Style Buckle, 50 mm 1600 mm / 2"	NC
☐ 600A849	Positioning Belt, Aircraft Buckle, Padded 50 mm x 1300 mm / 2"	\$ 101.00
☐ 600A850	Bodypoint - Evoflex Pelvic Stabilizer Belt	\$ 157.00
☐ 600A851	Bodypoint - Evoflex Pelvic Stabilizer Belt with Extension	\$ 268.00

Headrest Options - MPS Style			\$
☐ 600NU12		Cloth Headrest - MPS seat back only (12" x 4.5" x 6.25")	\$ 544.00
☐ 600NU13		Vinyl Headrest - MPS seat back only (12" x 4.5" x 6.25")	\$ 544.00
Headrest Options - REHAB Style			\$
Headrest Pad			
☐ 600HE39		Whitmyer® Plush Headrest Pad 254 mm / 10"	\$ 371.00
☐ 600HE46		Whitmyer® Adjust-A-Plush Headrest Pad Medium (Adjustable from 7 to 16" Wide x 4" Tall)	\$ 438.00
Headrest Pad			
☐ 600MH33		Whitmyer® Linx Headrest Mount with Quick Release	\$ 243.00
☐ 600MH34		Whitmyer® Linx-2 Headrest Mount with Quick Release	\$ 294.00
Armrest Style Options			\$
☐ 600EG53		Flip-Up Removable Armrests (Post Style) - Attaches to base. Height adj 3" & angle adj 90°-155°	\$ 1,762.00
☐ 600EG54		Cantilever Armrests (Reclining Flip Back) - Attaches to back. Height adj 4" & angle adj 90°	\$ 2,501.00
☐ 600EG55		Cantilever 2 Point Quick Release Post - Attaches front of armrest to seat frame	\$ 3,120.00
Armrest Pad Options			\$
☐ 600AB80		4" x 16" Tapered Pads - Right	\$ 146.00
☐ 600AE80		4" x 16" Tapered Pads - Left	\$ 146.00
☐ 600AB81		10" Waterfall Armrest Pad - Right	NC
☐ 600AE81		10" Waterfall Armrest Pad - Left	NC
☐ 600AB82		14" Waterfall Armrest Pad - Right	NC
☐ 600AE82		14" Waterfall Armrest Pad - Left	NC
LEGREST OPTIONS			
✓ OPTION		DESCRIPTION	RETAIL
Center Mount Fixed Legrest			\$
Center Mount Fixed Legrest			
☐ 600H210		Center Mount Flip-Up Footplate	\$ 601.00
Center Mount Fixed Legrest - Calfpads			
☐ 600DI59		Fixed Center Mount Calf Pads	\$ 687.00
Center Mount Power Elevating Legrest			\$
Center Mount Power Elevating Legrest			
☐ 600H213		Center Mount Power Elevating - 2 Piece Footplate - Includes calf pads	\$ 5,500.00
Center Mount Power Elevating Legrest Extension			
☐ 600E42		Center Mount Power Elevating Extension - Standard - Footplate range 5"-15"	NC
☐ 600E43		Center Mount Power Elevating Extension - Short - Footplate range 5"-12.5"	NC
Center Mount Power Elevating Legrest Footplates			
☐ 600F117		Center Mount Power Elevating Footplate - Standard	NC
☐ 600F118		Center Mount Power Elevating Footplate - Short	NC

Swing-Away Legrest \$

Swing-Away Legrest

☐ 600H211	Side Mount Swing-Away Legrest - 2 piece flip-up footplate (8.5" x 6"). 13" height adjustable	\$ 864.00
----------------------------------	--	-----------

S/A Non-Power Positioning Options

☐ 600LS15	Heel and Side Support - Attaches to the footplate stem and footplate	\$ 870.00
☐ 600LS16	Foot Strap - 1.5" wide strap	\$ 118.00

Swing-Away Power Elevating Legrest \$

Swing-Away Power Elevating Legrest

☐ 600H214	Side Mount Power Elevating Legrest - Swing away legrest with flip-up 2" height adj footplate	\$ 6,004.00
----------------------------------	--	-------------

S/A Power Legrest Positioning Options

☐ 600DI58	2" Metal Heel Guards for 2 Piece Footplates	\$ 1,134.00
☐ 600HL15	Heel Loops - Attaches to the footplate. Length 11.5-13" and Width 1"	\$ 168.00

ACCESSORIES \$

☐ 6008A93	Waterproof Wheelchair Cover	\$ 219.00
☐ 6008A94	Waterproof Joystick Cover	\$ 68.00
☐ 6008A95	Medical Necessity Bag - Velcro's to armrest (10.5" x 6" x 1.25")	\$ 50.00
☐ 600A221	Deluxe Backpack	\$ 142.00
☐ 600A615	Cup Holder - Self-leveling. Attaches to armrest opposite joystick	\$ 283.00
☐ 600A719	Articulating Oxygen Holder for MPS Seating System	\$ 1,470.00
☐ 600A612	Cane Holder for MPS Seating System	\$ 943.00
☐ 6008A96	Luggage Rack for 360 - Attaches to base behind the chair (10.5" x 13")	\$ 492.00
☐ 6008A97	Charge Mate - 12 Volt/Car outlet and USB/cell phone charger with max draw of 6 amps	\$ 791.00
☐ 600A613	Cell Phone Holder - Adjustable to fit most cell phones	\$ 438.00
☐ 600A614	Fishing Rod Holder	\$ 438.00

Be Adaptive Products

☐ 6008A98	Be Adaptive Limited Mobility Shooting Mount - Gun or crossbow	\$ 2,209.00
☐ 6008A99	Be Adaptive Trigger Mechanism	\$ 1,224.00
☐ 6008A100	Be Adaptive Limited Mobility Fishing Mount	\$ 2,209.00

ELECTRONICS

✓ OPTION		DESCRIPTION	RETAIL
Joystick			\$
☐ 600J121		Standard LED Joystick (w/o Screen Display)	NC
☐ 600J116		Color Joystick (2.3" Color LCD Screen Display, CJSM1)	NC
☐ 600J118		High Resolution Color Joystick,(CJSM2+BT) 3.5" Color LCD Display	\$ 746.00
Joystick Mount			\$
☐ 600JM34		Fixed on Armrest - Depth adjustable 8.25"	NC
☐ 600JM45		Swing-Away Joystick Mount - Includes latch to lock in the forward position	\$ 836.00
☐ 600JM46		Midline Joystick Mount - 2 position mount (Rotation 0° - 180°)	\$ 2,020.00
Joystick Mount - Location			\$
☐ 600HM16		Right Hand Armrest	NC
☐ 600HM17		Left Hand Armrest	NC
Display Language Selection R-net			\$
☐ 600L1		English	NC
☐ 600L3		French	NC
☐ 600L7		Spanish	NC
Joystick Handle Type			\$
Comes standard with PG Cone-Shaped Joystick Handle			
☐ 600JH13		Ball Handle	\$ 108.00
☐ 600JH26		Bodypoint Grooved Mushroom Joystick Handle	\$ 236.00
☐ 600JH24		Bodypoint – Rubber Dome Joystick Handle	\$ 252.00
☐ 600JH25		Bodypoint – U-Shaped Joystick Handle (4 inch)	\$ 252.00
Expandable Power Module			\$
Must select both items			
☐ 600EX05		Expandable 120A Controller - R-net	\$ 1,083.00
☐ 600EX06		Harness for 120A Expandable Controller	\$ 757.00
Power Seat Function			\$
☐ 600EO37		THRU DRIVE CNTRL 1 ACTUATOR- Must be selected if ordering 1 actuator function	\$ 1,422.00
☐ 600EO38		THRU DRIVE CNTRL 2+ ACTUATORS- Must be selected if ordering 2 or more actuator functions	\$ 2,860.00
Power Seat Control Option			\$
☐ 600QG11		Ctrl+5 Switch Box with Free Assignable Functions	\$ 1,002.00
☐ 600QG12		Ctrl+5 Switch Box with Free Assignable Functions (Multi-power only)	NC
Ctrl+5 Switch Type			\$
☐ 600EY22		Push Button	NC
☐ 600EY23		Toggle Switches	\$ 351.00
Power Seat Control Option - Mounting location			\$

☐ 600EY19		CTRL+5 MNT ARMREST RIGHT	NC
☐ 600EY20		CTRL+5 MNT ARMREST LEFT	NC
Attendant Controls			\$
☐ 600A05		Dual Control Module for Attendant	\$ 2,424.00
Attendant Control Mounting Location			\$
☐ 600DN5		Right Back Cane	NC
☐ 600DN6		Left Back Cane	NC
R-net Omni-2 - 3.5" LCD Display w/ Built-In Infrared & Blueto			\$
☐ 600ON14		R-net Omni-2 Only Input Control Module	\$ 3,558.00
☐ 600ON15		R-net Omni2 Bundled Input Control Module Available w/ Switch-It Input Devices	NC
R-net Omni-2 -Mounting Location			\$
☐ 600DS23		Right Hand Armrest	NC
☐ 600DS24		Left Hand Armrest	NC
Head Control Packages			\$
☐ 600EU25		Switch-It DualPro® Proportional Head Array - Facial Pads	\$ 6,953.00
☐ 600EU26		Switch-It DualPro® Proportional Head Array - Spot Pads	\$ 6,953.00