



# Quickie QM-7 Series Power Recline Seating

March 2017  
EIPW11

Account: _____	Ship To: _____
Date: _____	Name: _____
PO#: _____	Address: _____
Buyer: _____	City, Province, Postal Code: _____
Marked For: _____	Phone: _____

✓ Option #	Definition	Price	✓ Option #	Definition	Price
------------	------------	-------	------------	------------	-------

## 1. QUICKIE QM-7 Series Power Recline Seating (Select One)

Tilt AND Recline			
QUICKIE QM-710™ MPC			
<input type="checkbox"/>	Quickie QM-710™ MPC Tilt and Recline		
172PC69	Base Package		\$7,349.00
QUICKIE QM-715™ HD MPC			
<input type="checkbox"/>	Quickie QM-715™ HD MPC Tilt and Recline		
172PC104	Base Package		\$7,885.00
QUICKIE QM-720™ MPD			
<input type="checkbox"/>	Quickie QM-720™ MPD Tilt and Recline		
172PC87	Base Package		\$8,137.00

Recline ONLY			
QUICKIE QM-710™ SC			
<input type="checkbox"/>	Quickie QM-710™ SC for ASAP II Recline ONLY		
172PC68	Base Package		\$7,349.00
QUICKIE QM-715™ HD SC			
<input type="checkbox"/>	Quickie QM-715™ HD SC for ASAP II Recline ONLY		
172PC80	Base Package		\$7,885.00
QUICKIE QM-720™ SD			
<input type="checkbox"/>	Quickie QM-720™ SD for ASAP II Recline ONLY		
172PC86	Base Package		\$8,137.00

### MPC Seating Specifications

<b>Typical Chair Weight</b>	<b>W/battery</b>	375lbs	<b>W/O Battery</b>	285lbs
<b>Transit</b>	WC19 occupied approved			
<b>Seat Type</b>	3.7 Power Recline. NA with QM-715			
<b>Seat size range</b>	<b>Width</b>	14"-22"	<b>Depth</b>	16"-22"
<b>STFH range</b>	<b>MPC</b>	16" - 19.3"		
<b>Power Seating</b>	50° CG tilt and Power recline w/ power sheer reduction / 12" Seat elevator			

### Actuators Required by Power Seating Function

<b>Power Recline with Power Sheer Reduction:</b>	2 actuators
<b>Combination : Recline, Tilt:</b>	3 actuators
<b>Combination : Recline, Seat elevate:</b>	3 actuators
<b>Combination : Recline, Power Legrests:</b>	4 actuators
<b>Combination : Recline, Tilt, Power Legrests:</b>	5 actuators
<b>Combination : Recline, Seat elevate, Power Legrests:</b>	5 actuators
<b>Combination : Recline, Tilt, Power Legrests, Seat elevate:</b>	6 actuators



QM-710™ MPC  
Shown w/ optional Seat Elevator

### QUICKIE QM-7™ Base Specifications

<b>Drive Wheel Position</b>	Midwheel Drive	<b>Drive Wheels</b>	14" Drive Wheels w/ flat free insert (Black or Grey)
<b>User weight Cap</b>	300 lb - 350lbs	<b>Tires/Casters</b>	6" Front / 6" rear (Black or Grey)
<b>Base size</b>	<b>Width</b> 25" <b>Length</b> 36" <b>Turn Radius</b> 20"	<b>Charger</b>	8A Dual Mode Charger (Compact) CTE
<b>Drive Control</b>	R-Net PM120 (Non-Expandable and Expandable)	<b>Battery Size</b>	Group 24 / Group 34.
<b>Driving Range</b>	Up to 13.5 miles	<b>Shroud Colors</b>	12 standard colors, 4 upgrade prints.
<b>Drive Motors</b>	4-pole 6 mph, Optional 5mph and 8mph	<b>Ground Clearance</b>	4"

## 2. BASE CHAIR OPTIONS

Speed Package (4-pole 300lb)			
<input type="checkbox"/>	172PF85	<b>QM-710 Standard</b> 6 mph motor package	N/C
<input type="checkbox"/>	172PF91	QM-710 5 mph High Torque motor package	\$525.00
<input type="checkbox"/>	172PF90	QM-710 8 mph Community motor package	\$1,103.00
<input type="checkbox"/>	172PF84	<b>QM-715 HD Standard</b> 5 mph motor package	N/C
<input type="checkbox"/>	172PF89	<b>QM-720 Standard</b> 8 mph motor package	N/C

Shroud Color			
<input type="checkbox"/>	172EZ01	Candy Apple Red	N/C
<input type="checkbox"/>	172EZ02	High Gloss Black	N/C
<input type="checkbox"/>	172EZ03	Candy Blue	N/C
<input type="checkbox"/>	172EZ04	Silver	N/C
<input type="checkbox"/>	172EZ05	Copper	N/C
<input type="checkbox"/>	172EZ08	Green Apple	N/C
<input type="checkbox"/>	172EZ09	Yellow	N/C
<input type="checkbox"/>	172EZ10	Pearl Pink	N/C
<input type="checkbox"/>	172EZ11	Candy Purple	N/C
<input type="checkbox"/>	172EZ12	Matte Black	N/C
<input type="checkbox"/>	172EZ13	Midnight Blue	N/C
<input type="checkbox"/>	172EZ15	Pearl White	N/C
<input type="checkbox"/>	172EZ06	Digital Camo.	\$525.00
<input type="checkbox"/>	172EZ07	Woodland Camo.	\$525.00
<input type="checkbox"/>	172EZ14	Carbon Fiber Print	\$525.00
<input type="checkbox"/>	172EZ16	American Flag.	\$525.00
<input type="checkbox"/>	172EZ17	Pink Camo	\$525.00
<input type="checkbox"/>	172EZ18	Grunge Skulls	\$525.00

User weight			
<input type="checkbox"/>	172WG06	Less than 150lbs	N/C
<input type="checkbox"/>	172WG13	151-200lb	N/C
<input type="checkbox"/>	172WG14	201-250lb	N/C
<input type="checkbox"/>	172WG03	251-300lb	N/C
<input type="checkbox"/>	172WG15	301-350lb Transit Package NA	N/C

Drive Wheel			
<input type="checkbox"/>	172DW59	14" Cast Wheel Solid Dual Density Inserts <b>Black Knobby Tires</b> and <b>Black Caster</b>	Std

Suretrac Drive Tracking <b>NEW!</b>			
<input type="checkbox"/>	172MI01	Suretrac Module and Motor Encoders	\$1,475.00

Transit Option			
<input type="checkbox"/>	172WN01	4-Point Transit Tie-down Package <i>Meets WC-19 for Occupant Use</i>	<b>NEW!</b> \$373.00
<input type="checkbox"/>	172TO40	Omit Transit Tie-down	N/C


WC19 Transit Lap Belt <b>NEW!</b>			
<input type="checkbox"/>	172WN02	WC19 Transit Lap Belt	\$250.00

Wheel locks			
<input type="checkbox"/>	172WL42	Wheel locks	\$142.00

Ventilator Options			
<input type="checkbox"/>	172VT4	Fixed Vent Tray (Pulmonic LTV Models)	\$525.00

Pre-Set Seat Frame Angle (STFH increases 1" per increment change)			
<input type="checkbox"/>	172FA10	0°	N/C

Seat Height (Use STFH Matrix below)			
<i>(Measured in the front) - Obtain Quote to verify est. STFH</i>			
<input type="checkbox"/>	172HT44	Position "0"	N/C
<input type="checkbox"/>	172HT49	+5" STFH	Requires Group 34 Batteries
<input type="checkbox"/>	172HT45	+1" STFH	
<input type="checkbox"/>	172HT46	+2" STFH	N/C
<input type="checkbox"/>	172HT47	+3" STFH	N/C

LED Light Package			
<input type="checkbox"/>	172LI05	 LED Light Package: Includes headlights, lights, turn indicators and flashers. <b>NEW!</b>	\$1,066.00

STFH Matrix					
Seating Type and Power Option Combinations	Seat-to-Floor Height at 0 tilt				
	0	+5"	+1"	+2"	+3"
SC Power Recline only			16.1	17.1	18.1
SC Power Recline and 9" Lift Seat		16.2	16.7	17.7	18.7
MPC Power Recline and Tilt.			16.3	17.3	18.3
MPC Power Recline with Tilt & 12" Lift combo		16**	16.5	17.5	18.5

STFH measured with batteries but no user. Actual STFH may vary by .5"  
 \*Adding seat angle will increase STFH by 1" for each incremental change.  
 \*\*Ultra low STFH may cause ground clearance of Power Centermount to drop below 2".

## 3. BATTERIES and CHARGERS

Batteries			
<input type="checkbox"/>	172AG04	MK Gr 24 Flat Top Gel Cell (Pair) w/straps <i>Recommended for Max Daily Range</i>	\$830.00
<input type="checkbox"/>	172AG11	MK Gr 34 Flat Top Gel Cell (Pair) w/straps, <i>Required for lower STFH.</i>	\$830.00

Battery Charger			
<input type="checkbox"/>	172BJ19	8A Dual Mode (Compact) CTE	N/C

## 4. ELECTRONICS

### Drive Control Options

#### Expandable Controllers

<input type="checkbox"/>	172ER06	R-Net Expandable Drive Control (120A)	\$1,066.0
	172EX05	Expandable 120A Controller - R-net	\$1,066.00
	172EX06	Harness for 120A Expandable Controller	NC

### R-net Joystick

<input type="checkbox"/>	172J70	Std LED Joystick for Rnet Controls	N/C
<input type="checkbox"/>	172J108	Joystick w/Stereo Jacks & Color Display	N/C
<input type="checkbox"/>	172J88	Joystick W/Toggle, Spd Pot, Stereo Jack & Color Display	\$158.00
<input type="checkbox"/>	172J112	Joystick Advanced: 3.5" Color Display, Paddle Switches, Shortcut Keys and Stereo Jack.	\$1,133.00
<input type="checkbox"/>	172J95	Omit JS. For use with OMNI	N/C

### Joystick Mount(Fixed, Swing-Away, Retractable, Midline)

<input type="checkbox"/>	172JM10	Fixed Mount	N/C
<input type="checkbox"/>	172JM11	Swing Away Mount	\$388.00
<input type="checkbox"/>	172JM15	Retractable Joystick mount w/45° stops	\$721.00
<input type="checkbox"/>	172JM14	Mid Line Mount	\$388.00

### Joystick Mount (Location)

*Height Adjustable.*

<input type="checkbox"/>	172HM1	Right Side Mount	N/C
<input type="checkbox"/>	172HM2	Left Side Mount	N/C

### Joystick Handle (Optional)

<input type="checkbox"/>	172JH1	Ball Handle	\$95.00
<input type="checkbox"/>	172JH2	Foam Ball Handle	\$95.00
<input type="checkbox"/>	172JH4	T-Handle	\$95.00
<input type="checkbox"/>	172JH6	Goalpost	\$95.00

### Attendant Control

<input type="checkbox"/>	172A04	Attendant Control for R-net	\$1,066.00
--------------------------	--------	-----------------------------	------------

### Attendant Mounting location

<input type="checkbox"/>	172DN1	Mounting location: Right Back Cane	N/C
<input type="checkbox"/>	172DN2	Mounting Location: Left Back Cane	N/C

### R-net OMNI and Output Options

*(For expandable controls, choose all that apply)*

<input type="checkbox"/>	172ON6	Rnet OMNI Input Control Module and display with built in Infrared and armrest mount. Choose mount below. (required for Specialty Input Devices.)	\$2,075.00
<input type="checkbox"/>	172ON7	Rnet Output Module (ECM)	\$1,066.00
<input type="checkbox"/>	172ON8	Rnet Bluetooth Mouse-mover for Windows / Android	\$1,066.00
<input type="checkbox"/>	172ON9	Rnet Bluetooth for Apple Devices	\$1,066.00

### OMNI Mount

*(If used with Joystick must be mounted opposite side)*

<input type="checkbox"/>	172DS04	Swing Away Left Side OMNI Mount	NC
<input type="checkbox"/>	172DS03	Swing Away Right Side OMNI Mount	NC

### Switch inputs for Joysticks or OMNI

*Drive or Mode / Profile Control*

<input type="checkbox"/>	172SW22	Stereo to 2 Mono Splitter	\$75.00
<input type="checkbox"/>	172SW6	Ribbon Switch	\$119.00
<input type="checkbox"/>	172SW1	Buddy Button Switch	\$95.00
<input type="checkbox"/>	172SW2	Disc Switch	\$202.00
<input type="checkbox"/>	172SW9	Wobble Switch	\$630.00
<input type="checkbox"/>	172SW4	Micro Light Switch	\$95.00
<input type="checkbox"/>	172SW24	Egg Switch-Black	\$119.00
<input type="checkbox"/>	172SW21	Dual toggle switch and bracket	\$214.00
<input type="checkbox"/>	172SW35	Switch-it Mini 1 button on Mono plug	\$268.00
<input type="checkbox"/>	172SW36	Switch-it Mini 2 button on 2 Mono plugs	\$268.00
<input type="checkbox"/>	172SW37	Switch-it Mini 2 button on 1 Stereo plug	\$268.00
<input type="checkbox"/>	172SW38	Switch-it Mini 4 button on DB9 connector	\$268.00

### Switch Mount Bracket

<input type="checkbox"/>	172EM01	Single Switch Mount Bracket	\$247.00
--------------------------	---------	-----------------------------	----------

### Rnet Electronic Programmer

<input type="checkbox"/>	172PR5	R-net PC Programmer with Dongle	\$615.00
--------------------------	--------	---------------------------------	----------

## 5. Specialty Driver Control Options

*Requires expandable Rnet controls and OMNI unless otherwise noted. Not available with Upholstered backs.*

### 5a. Proportional Control Packages

<b>Head Control Proportional Packages</b>			
<input type="checkbox"/>	<b>Switch-it Dual Pro Proportional Head Control with Whitmyer and Facial Pads (Adult)</b>		<b>\$6,541.00</b>
	172ID72	Switch-it Dual Pro Head Control w/facial pad (Adult)	\$6,139.00
	172IJ09	SA Boom Mount with Egg Switch Control	\$402.00
<input type="checkbox"/>	<b>Switch-it Dual Pro Proportional Head Control with Spot Pads Package</b> <span style="border: 1px solid black; padding: 2px;">NEW!</span>		<b>\$6,541.00</b>
	172ID70	Switch-it Dual Pro Head Control w/Spot pads	\$6,139.00
	172IJ09	SA Boom Mount with Egg Switch Control	\$402.00

<b>Hand Control Proportional Joystick Packages</b>			
<input type="checkbox"/>	<b>Switch-it Versa-Guide Compact Joystick Midline with Gatlin Package</b>		<b>\$2,369.00</b>
	172ID71	Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
	172IJ07	S/A Gatlin Mount	\$330.00
<input type="checkbox"/>	<b>Switch-it Versa-Guide Compact Joystick Midline and prox flex switches with Gatlin Package</b>		<b>\$3,018.00</b>
	172ID71	Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
	172IM08	Dual Flex Proximity Switches (requires no force)	\$649.00
	172IJ07	S/A Gatlin Mount	\$330.00
<input type="checkbox"/>	<b>Switch-it Micro-Pilot Joystick with Space Disc Package</b>		<b>\$4,007.00</b>
	172ID62	Micro-Pilot Proportional Joystick	\$3,534.00
	172IJ08	S/A SpaceDisc mount	\$473.00
<input type="checkbox"/>	<b>Switch-it Micro-Guide Joystick with SpaceDisc Package</b>		<b>\$4,831.00</b>
	172ID63	Micro-Guide Proportional Joystick	\$4,358.00
	172IJ08	S/A SpaceDisc mount	\$473.00
<input type="checkbox"/>	<b>ASL Proportional Micro Extremity Joystick (ASL130) Midline with Gatlin Package</b>		<b>\$4,414.00</b>
	172ID59	Micro Extremity Joystick	\$4,084.00
	172IJ07	S/A Gatlin Mount	\$330.00
<input type="checkbox"/>	<b>ASL Mushroom Joystick (ALS135) Midline with Gatlin Package</b>		<b>\$3,353.00</b>
	172ID60	Mushroom Joystick	\$3,023.00
	172IJ07	S/A Gatlin Mount	\$330.00
<input type="checkbox"/>	<b>ASL Compact Joystick w Integrated buttons (ASL134) Midline with Gatlin Package</b>		<b>\$1,530.00</b>
	172ID61	ASL Compact Joystick	\$1,200.00
	172IJ07	S/A Gatlin Mount	\$330.00
<input type="checkbox"/>	<b>ASL Micro Mini Joystick w Hand Pad (ASL136) Midline with Gatlin Package</b>		<b>\$4,012.00</b>
	172ID68	Micro Mini Joystick	\$3,682.00
	172IJ07	S/A Gatlin Mount	\$330.00

<b>Chin Control Proportional Joystick Packages</b>			
<input type="checkbox"/>	<b>Switch-it Versa-Guide Chin Control Compact Joystick and Egg Switch Package</b>		<b>\$2,406.00</b>
	172ID71	Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
	172IM07	S/A Mount off backrest	\$252.00
	172IJ04	SA Boom Mount with egg switch	\$115.00
<input type="checkbox"/>	<b>Switch-it Versa-Guide Chin Control Compact Joystick and Proximity Flex Switches Package</b>		<b>\$2,940.00</b>
	172ID71	Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
	172IM08	Dual Flex Proximity Switches (requires no force)	\$649.00
	172IM07	S/A Mount off backrest	\$252.00
<input type="checkbox"/>	<b>Switch-it Versa-Guide Bib Chin Control Compact Joystick with Proximity Flex switches Package Std 6"</b>		<b>\$3,445.00</b>
	172ID71	Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
	172IM08	Dual Flex Proximity Switches (requires no force)	\$649.00
	172IM05	Bib mount Small	\$757.00
<input type="checkbox"/>	<b>Switch-It Versa-Guide Bib Chin Control Compact Joystick with Proximity Flex Switches Package Large 8"</b>		<b>\$3,512.00</b>
	172ID71	Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
	172IM08	Dual Flex Proximity Switches (requires no force)	\$649.00
	172IM06	Bib mount Large	\$824.00
<input type="checkbox"/>	<b>Switch-It Versa-Guide Bib Chin Control Compact Joystick with Egg Switch Package Std 6"</b>		<b>\$2,911.00</b>
	172ID71	Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
	172IM05	Bib mount Small	\$757.00
	172IJ04	SA Boom Mount with egg switch	\$115.00
<input type="checkbox"/>	<b>Switch-It Versa-Guide Bib Chin Control Compact Joystick with Egg Switch Package Large 8"</b>		<b>\$2,978.00</b>
	172ID71	Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
	172IM06	Bib mount Large	\$824.00
	172IJ04	SA Boom Mount with egg switch	\$115.00
<input type="checkbox"/>	<b>Switch-It Micro-Pilot Bib Chin Control with Egg Switch Package Std 6"</b>		<b>\$4,406.00</b>
	172ID62	Micro-Pilot Proportional Joystick	\$3,534.00
	172IM05	Bib mount Small	\$757.00
	172IJ04	SA Boom Mount with egg switch	\$115.00
<input type="checkbox"/>	<b>Switch-It Micro-Pilot Bib Chin Control with Egg Switch Package Large 8"</b>		<b>\$4,473.00</b>
	172ID62	Micro-Pilot Proportional Joystick	\$3,534.00
	172IM06	Bib mount Large	\$824.00
	172IJ04	SA Boom Mount with egg switch	\$115.00
<input type="checkbox"/>	<b>Switch-It Micro-Guide Bib Chin Control with Egg Switch Package Std 6"</b>		<b>\$5,230.00</b>
	172ID63	Micro-Guide Proportional Joystick	\$4,358.00
	172IM05	Bib mount Small	\$757.00
	172IJ04	SA Boom Mount with egg switch	\$115.00
<input type="checkbox"/>	<b>Switch-It Micro-Guide Bib Chin Control with Egg Switch Package Large 8"</b>		<b>\$5,297.00</b>
	172ID63	Micro-Guide Proportional Joystick	\$4,358.00
	172IM06	Bib mount Large	\$824.00
	172IJ04	SA Boom Mount with egg switch	\$115.00
<input type="checkbox"/>	<b>ASL Proportional Micro Extremity Joystick Chin / Lip with Egg Switch Control Package</b>		<b>\$4,498.00</b>
	172ID59	Micro Extremity Joystick	\$4,084.00
	172IM07	S/A Mount off backrest	\$299.00
	172IJ04	SA Boom Mount with egg switch	\$115.00

ORDERING INFORMATION: PHONE: (800) 263-3390 FAX: (800) 561-5834 WEBSITE: www.sunrisemedical.ca

Specifications are subject to change without notice. Copyright © 2016

MK-100308 REV A (US) REV 7 (CA)

**5. Specialty Driver Control Options (cont.)**

**5b. Switch Control**

*Requires Headrest and/or Headrest Mount*

**3 Switch Proximity Head Arrays** Choose headrest mount below

<input type="checkbox"/> <b>ASL Proximity Head Array Package (Adult)</b>		<b>\$5,755.00</b>
172II02	ASL Proximity Head Array (Adult)	\$5,353.00
172IJ09	SA Boom Mount with Egg Switch Control	\$402.00
<input type="checkbox"/> <b>ASL Proximity Head Array Package (Ped)</b>		<b>\$5,755.00</b>
172II03	ASL Proximity Head Array (Ped)	\$5,353.00
172IJ09	SA Boom Mount with Egg Switch Control	\$402.00
<input type="checkbox"/> <b>Switch-It Head Array Package 3 Switch Curved (Long)</b>		<b>\$5,755.00</b>
172II13	Switch-it Proximity Head Array Long	\$5,353.00
172IJ09	SA Boom Mount with Egg Switch Control	\$402.00
<input type="checkbox"/> <b>Switch-It Head Array Package 3 switch curved (Stubby)</b>		<b>\$5,691.00</b>
172II14	Switch-it Proximity Head Array Stubby	\$5,289.00
172IJ09	SA Boom Mount with Egg Switch Control	\$402.00

**5 Switch Proximity Head Array (No Egg Switch Needed)** Choose "Headrest Mount for Special Controls"

<input type="checkbox"/> 172II04	Switch-It Head Array 5 switch flat	\$5,353.00
----------------------------------	------------------------------------	------------

**Multiple Switches.**

<input type="checkbox"/> 172SW31	Switch-it 3 Naked proximity switches on DB9 connector. (requires OMNI)	\$3,039.00
<input type="checkbox"/> 172SW32	Switch-it 4 Naked proximity switches on DB9 connector. (requires OMNI)	\$3,718.00
<input type="checkbox"/> 172SW33	Switch-it 5 Naked proximity switches on DB9 connector. (requires OMNI)	\$4,146.00
<input type="checkbox"/> 172SW34	ASL (502) 5-Switch adaptor box on DB9 connector. (for use with 5 mono switches, not included)	\$4,146.00
<input type="checkbox"/> 172SW23	4-way toggle switch button DB9 Connector. Includes bracket	\$670.00

**Sip N Puff Package** Choose headrest and HR mount below

<input type="checkbox"/> <b>Quickie Sip n' Puff System</b>		<b>\$2,635.00</b>
172ID34	Sip N' Puff Kit	\$2,063.00
172IH1	Sip n' Puff Tubing and SA Hardware	\$457.00
172IJ10	SA Boom Mount with egg switch	\$115.00

**Headrest Mount for Specialty Controls**

<input type="checkbox"/> 172MH12	Axys Mount & Bracket for Specialty Controls	\$309.00
<input type="checkbox"/> 172MH13	Cobra Mount & Bracket for Specialty Controls	\$386.00
<input type="checkbox"/> 172MH14	Linx Mount & Bracket for Specialty Controls	\$309.00

**Headrest Pad for Specialty controls**

<input type="checkbox"/> 172HE79	Plush Pad 8" for Specialty controls	\$309.00
<input type="checkbox"/> 172HE80	Plush Pad 10" for Specialty controls	\$309.00
<input type="checkbox"/> 172HE81	Plush Pad 14" for Specialty controls	\$309.00

**Extra Swing Away Head Boom Egg Switch and Mount**

172IJ05	Extra SA Boom Mount with egg switch #1	\$115.00
172IJ06	Extra SA Boom Mount with egg switch #2	\$115.00

**Single Switches**

<input type="checkbox"/> 172IB01	Ribbon Switch	\$119.00
<input type="checkbox"/> 172IB02	Buddy Button Switch	\$108.00
<input type="checkbox"/> 172IB03	Disc Switch	\$202.00
<input type="checkbox"/> 172IB04	Wobble Switch	\$630.00
<input type="checkbox"/> 172IB05	Micro Light Switch	\$95.00
<input type="checkbox"/> 172IB06	Egg Switch-Black	\$134.00
<input type="checkbox"/> 172IB07	Switch-it Mini 1 button on Mono plug	\$268.00

## 6. SEAT TYPE

172ST57 Power Recline NC

### Recline Seat Frame size

<input type="checkbox"/> 172FU26	Power Recline Seat Frame ( G ) 14x16 to 22x20 (Not width Adjustable)	N/C
<input type="checkbox"/> 172FU27	Power Recline Seat Frame ( H ) 18x18 to 22"x22" (Not width Adjustable)	N/C

## 7. SEAT SIZE

### Seat Size Matrix. *Seat type is chosen in section 5.*

G/H = Recline Seat (Width not adjustable)

Note: Power Recline Seat Depth measured from solid surface of back to front of seat pan. Please accommodate for backrest cushion thickness.

Seat Depth (Adjustable)	Recline Seat Width				
	14"	16"	18"	20"	22"
16"	G	G	G	G	G
17"	G	G	G	G	G
18"	G	G	G,H	G,H	G,H
19"	G	G	G,H	G,H	G,H
20"	G	G	G,H	G,H	G,H
21"			H	H	H
22"			H	H	H

### Seat Width

<input type="checkbox"/> 172R23	14"	Seat width	(N/C)
<input type="checkbox"/> 172R5	16"	Seat width	(N/C)
<input type="checkbox"/> 172R7	18"	Seat width	(N/C)
<input type="checkbox"/> 172R9	20"	Seat width	(N/C)
<input type="checkbox"/> 172R25	22"	Seat width	(N/C)

### Seat Depth

<input type="checkbox"/> 172RF5	16"	Seat depth	(N/C)
<input type="checkbox"/> 172RF6	17"	Seat depth	(N/C)
<input type="checkbox"/> 172RF7	18"	Seat depth	(N/C)
<input type="checkbox"/> 172RF8	19"	Seat depth	(N/C)
<input type="checkbox"/> 172RF24	20"	Seat depth	(N/C)
<input type="checkbox"/> 172RF26	21"	Seat depth	(N/C)
<input type="checkbox"/> 172RF25	22"	Seat depth	(N/C)

## 8. SEAT CUSHIONS

### JAY Cushion Bundling

<input type="checkbox"/> 172JC14	Jay Cushion <u>On Chair Order</u> <u>See order Form for codes</u>	See Order Form
----------------------------------	--	----------------

## 9. POWERED SEATING

### Power Seating Options

<input type="checkbox"/> 172PA24	<u>Recline Only</u> w/Power Shear Reduction (2 actuators)	\$8,568.00
<input type="checkbox"/> 172PA25	<u>Tilt and Recline</u> w/Power Shear Reduction (3 actuators)	\$12,324.00
<input type="checkbox"/> 172PA16	QM715 HD <u>Recline Only</u> w/Power Shear Reduction Only (350lbs weight limit) (2 actuators)	\$6,966.00
<input type="checkbox"/> 172PA17	QM-715 HD <u>Tilt and Recline</u> w/Power Shear Reduction (350 weight limit) (3 actuators)	\$12,324.00

### Powered Seat Elevator- (AKA: Power Lift)

<input type="checkbox"/> 172PA26	9" Power Seat Elevate for use without tilt (NA with QM-715 HD)	\$3,257.00
<input type="checkbox"/> 172PA27	12" Power Seat Elevate for use with tilt (Must be ordered with SPOT tilt, NA with QM-715 HD)	\$3,452.00

### Power Seating Control

#### Powered Seating Thru-drive Control Method.

##### Thru-Drive 2+ Actuator

<input type="checkbox"/> 172EO29	Thru Drive 2+ Actuators: (For operating Tilt, recline, lift, power legs) Includes Ctrl+5	\$2,525.00
----------------------------------	--	------------

#### Powered Seating Thru-drive Control Options

##### Assignable Buttons

"Assignable Button" feature allows for customization of joystick buttons and/or Ctrl+5 Switches allowing access to various wheelchair and seating functions. When enabled, chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle.

<input type="checkbox"/> 172PM01	Enable Assignable buttons	N/C
----------------------------------	---------------------------	-----



#### Ctrl+5 Switchbox Options

Ctrl+5 allows direct control to 5 power seating functions (10 total switches). Switches can also be used with assignable buttons for greater flexibility. Box includes 3-Amp USB port for charging portable devices

<input type="checkbox"/> 172QG01	Ctrl+5 Switchbox mounted behind joystick	N/C
<input type="checkbox"/> 172QG02	Ctrl+5 Switchbox mounted opposite joystick	N/C
<input type="checkbox"/> 172QG03	Ctrl+5 Switchbox mounted on back cane	N/C
<input type="checkbox"/> 172QG04	Ctrl+5 Switchbox mounted on seat rail	N/C
<input type="checkbox"/> 172QG05	Omit Ctrl+5	-\$618.00

#### Ctrl+5 Switch Type

<input type="checkbox"/> 172EY03	Push buttons	N/C
<input type="checkbox"/> 172EY04	Toggle switches	\$309.00

ORDERING INFORMATION: PHONE: (800) 263-3390 FAX: (800) 561-5834 WEBSITE: www.sunrisemedical.ca

Specifications are subject to change without notice. Copyright © 2016

MK-100308 REV A (US) REV 7 (CA)

## 11. BACKREST AND BACK OPTIONS for 3.7 POWER RECLINE

Backrest Frame Type			
<input type="checkbox"/>	172BT89	Power Recline (Power Recline Seat only)	N/C

Power Recline Backrest Preset Home Angle N/C					
<input type="checkbox"/>	90° 172QK1	<input type="checkbox"/>	95° 172QK2 default	<input type="checkbox"/>	100° 172QK3

Back Type (Solid Recline)			
<input type="checkbox"/>	172BU83	Standard Curve back with 1.5" Foam cushion. (Reduces Seat Depth)	N/C
<input type="checkbox"/>	172JB03	Jay J3 for Power Recline w/1.5" Foam Choose type below. (Reduces Seat Depth)	N/C

Back Height - Required for Standard Curved Recline Backs and Jay 3 recline backs.			
<i>(Limited combinations of widths) Measured from pan to top of solid back.</i>			
<input type="checkbox"/>	172B42	22"	N/C
<input type="checkbox"/>	172B44	24"	N/C
<input type="checkbox"/>	172B140	26"	N/C
<input type="checkbox"/>	172B141	27"	N/C
<input type="checkbox"/>	172B142	28"	N/C

Push Handle Options			
<input type="checkbox"/>	172BP55	Fixed "Cane" Handles.	N/C
<input type="checkbox"/>	172BP56	Height Adjustable Push Handles. w/2 utility hooks	\$53.00
<input type="checkbox"/>	172BP57	Omit push handles	N/C

### 11a. JAY 3 BACK for 3.7 POWER RECLINE OPTIONS

Jay 3 Contour.			
<input type="checkbox"/>	172JP01	J3 PA Back (2.2" of Contour Depth) Available 14"-20" Wide (J3PAUTZYXX)	\$803.00
<input type="checkbox"/>	172JP02	J3 PA Back (2.2" of Contour Depth) Available 22" Wide (J3PAUTZ2FX)	\$1,708.00
<input type="checkbox"/>	172JP03	J3 PL Back (3.0" of Contour Depth) Available 14"-20" Wide (J3PLUTZYXX)	\$867.00
<input type="checkbox"/>	172JP04	J3 PDL Back (6.0" of Contour Depth) Available 14"-20" Wide (J3PDLUTZYXX)	\$930.00
<input type="checkbox"/>	172JP05	J3 PDC Back (6.0" of Contour Depth) Available 14"-20" Wide (J3PDCUTZYXX)	\$940.00
<input type="checkbox"/>	172JP06	J3 PD Back (6.0" of Contour Depth) Available 14"-20" Wide (J3PDUTZYXX) <b>NEW!</b>	Price List

Jay 3 Back Height. Backs are height adjustable			
<input type="checkbox"/>	172JR06	MT for 22" tall back. (Mid Thoracic Short 15" shell)	N/C
<input type="checkbox"/>	172JR01	UT for 24" tall back. (Upper Thoracic Short 16.5" shell)	N/C
<input type="checkbox"/>	172JR02	UT for 26" tall back. (Upper Thoracic Med 18" shell)	N/C
<input type="checkbox"/>	172JR03	UT for 27" tall back. (Upper Thoracic Tall 19.5" shell)	N/C
<input type="checkbox"/>	172JR04	SH for 28" tall back and widths 14-20" (Shoulder Height 21" shell, )	N/C
<input type="checkbox"/>	172JR05	SH for 26-28" tall back and 22" wide seat (Shoulder Height 18" shell, )	N/C

SPINE-ALIGN MINI FIT KIT 2			
<i>Includes: 5 wedges, 2 lat shims, 2 post/lat pelvic, 1 lg lumbar, 1 flat shim, and 1 curved</i>			
<input type="checkbox"/>	172JS01	SPINE-ALIGN MINI FIT KIT 2	\$190.00




Jay 3 Width			
<i>Back can be ordered 2" narrower than seat pan. Armrest will be adjusted to back width ordered. Future growth can be accomplished by replacing back with wider versions and adjusting arms outward.</i>			
<input type="checkbox"/>	172JQ01	Back width 14" (Use w/ 14-16" wide seat)	N/C
<input type="checkbox"/>	172JQ02	Back width 16" (Use w/ 16-18" wide seat)	N/C
<input type="checkbox"/>	172JQ03	Back width 18" (Use w/ 18-20" wide seat)	N/C
<input type="checkbox"/>	172JQ04	Back width 20" (Use w/ 20-22" wide seat)	N/C
<input type="checkbox"/>	172JQ05	Back width 22" (Use with 22" wide seat)	N/C

## 12. HEADREST OPTIONS

Choose One Each.

Headrest Mount			
<input type="checkbox"/>	172MH10	Axys Mount & Bracket	\$315.00
<input type="checkbox"/>	172MH11	Cobra Mount and Bracket	\$578.00
<input type="checkbox"/>	172MH1	Linx Mount and Bracket	\$237.00

Headrest Pad			
<input type="checkbox"/>	172HE65	Plush Head Support Pad 8"	\$237.00
<input type="checkbox"/>	172HE66	Plush Head Support Pad 10"	\$237.00
<input type="checkbox"/>	172HE78	Plush Head Support Pad 14"	\$237.00
<input type="checkbox"/>	172HE82	Adjusta Plush Head Support Pad Narrow	\$365.00
<input type="checkbox"/>	172HE83	Adjusta Plush Head Support Pad Standard	\$365.00



### 13. POSTURAL SUPPORTS

<b>LATERAL THORACIC SUPPORT (Swing Away)</b>			
<i>Available on all Solid Backs.</i>			
<b>LAT 4X4 S/A MOUNTS W/ 1" OFFSET - PR \$409.00</b>			
<input type="checkbox"/>	172ZV14	LATERAL 4X4 W/ 1" OFFSET - PAIR	\$126.00
<input type="checkbox"/>	172ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$283.00
<b>LAT 4X5 S/A MOUNTS W/ 1" OFFSET - PR \$409.00</b>			
<input type="checkbox"/>	172ZV15	LATERAL 4X5 W/ 1" OFFSET - PAIR	\$126.00
<input type="checkbox"/>	172ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$283.00
<b>LAT 4X6 S/A MOUNTS W/ 1" OFFSET - PR \$409.00</b>			
<input type="checkbox"/>	172ZV16	LATERAL 4X6 W/ 1" OFFSET - PAIR	\$126.00
<input type="checkbox"/>	172ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$283.00
<b>LAT 6X4 S/A MOUNTS W/ 1" OFFSET - PR \$409.00</b>			
<input type="checkbox"/>	172ZV17	LATERAL 6X4 W/ 1" OFFSET - PAIR	\$126.00
<input type="checkbox"/>	172ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$283.00
<b>LAT 6X5 S/A MOUNTS W/ 1" OFFSET - PR \$409.00</b>			
<input type="checkbox"/>	172ZV18	LATERAL 6X5 W/ 1" OFFSET - PAIR	\$126.00
<input type="checkbox"/>	172ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$283.00
<b>LAT 6X6 S/A MOUNTS W/ 1" OFFSET - PR \$409.00</b>			
<input type="checkbox"/>	172ZV19	LATERAL 6X6 W/ 1" OFFSET - PAIR	\$126.00
<input type="checkbox"/>	172ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$283.00

<b>CHEST STRAP</b>			
<input type="checkbox"/>	172ZR03	Padded Chest Strap	\$283.00

<b>Swing Away Knee Adductor Supports</b>			
<i>Available only with Fixed and Power Center Mount Legrests.</i>			
<i>Not available with Dual Post Arm Rests or 2-point Locking Cantilever Arms</i>			
<b>Knee Adductor Supports - LEFT (choose Pad and Arm)</b>			
<b>Pad Size</b>			
<input type="checkbox"/>	172KB5	Knee Adductor Pad SM- Left	\$215.00
<input type="checkbox"/>	172KB7	Knee Adductor Pad LG- Left	\$215.00
<b>S/A Arm Size</b>			
<input type="checkbox"/>	172XX04	Swing Away Short - Left	\$283.00
<input type="checkbox"/>	172XX06	Swing Away Long - Left	\$283.00
<b>Knee Adductor Supports - RIGHT (choose Pad and Arm)</b>			
<b>Pad Size</b>			
<input type="checkbox"/>	172KB6	Knee Adductor Pad SM- Right	\$215.00
<input type="checkbox"/>	172KB8	Knee Adductor Pad LG- Right	\$215.00
<b>S/A Arm Size</b>			
<input type="checkbox"/>	172XX05	Swing Away Short - Right	\$283.00
<input type="checkbox"/>	172XX07	Swing Away Long - Right	\$283.00
<b>Long Thigh Support w/Removable Mount</b>			
<i>Available with Cantilever Arms, NA with Manual Recline Back, Dual Post, Single Post or Cantilever w/ 2 pt Locking arms.</i>			
<b>SPOT Long Thigh Support - Left \$409.00</b>			
<input type="checkbox"/>	172TV2	Long Thigh Support - Left	\$126.00
<input type="checkbox"/>	172TX1	Removable Mount - Left	\$283.00
<b>SPOT Long Thigh Support - Right \$409.00</b>			
<input type="checkbox"/>	172TV1	Long Thigh Support - Right	\$126.00
<input type="checkbox"/>	172TX2	Removable Mount - Right	\$283.00

### 14. ARMREST AND ARMPAD OPTIONS

<b>Armrest for Power Recline Seating.</b>			
<i>(Armrest width can be adjusted 2" narrower than seat)</i>			
<input type="checkbox"/>	172EG26	STD Cantilever/Reclining Arm 10-14" Height Range	N/C
<input type="checkbox"/>	172EG27	Cantilever/Reclining Arm w/ 2 Point Locking Post. (10"-12" Ht Range.)	\$315.00
<input type="checkbox"/>	172EG28	Cantilever/Reclining Arm w/ 2 Point Locking Post. (11"-13" Ht Range.)	\$315.00

<b>Reclining Armpad Options</b>			
<i>(Armpad can be angled inward 15 degrees)</i>			
<b>Standard Reclining Armpad</b>			
<b>Left Armpad (N/C)</b>		<b>Right Armpad (N/C)</b>	
<input type="checkbox"/>	172AE05	St Reclining Arm Pad (16")	<input type="checkbox"/> 172AB05
<b>Comfort Reclining Armpad (Foam)</b>			
<b>Left Armpad (\$98.00)</b>		<b>Right Armpad (98.00)</b>	
<input type="checkbox"/>	172AE22	Comfort 2" wide Reclining Arm Pad	<input type="checkbox"/> 172AB22
<input type="checkbox"/>	172AE23	Comfort 4" wide Reclining Arm Pad	<input type="checkbox"/> 172AB23
<b>Comfort Reclining Armpad (Gel)</b>			
<b>Left Gel Armpad (\$185.00)</b>		<b>Right Gel Armpad (185.00)</b>	
<input type="checkbox"/>	172AE25	Comfort 2" wide Reclining Arm Pad	<input type="checkbox"/> 172AB25
<input type="checkbox"/>	172AE26	Comfort 4" wide Reclining Arm Pad	<input type="checkbox"/> 172AB26
<b>AES Trough/Elbow Block Reclining Armpads (Gel)</b>			
<b>Left Gel Armpad (\$258.00)</b>		<b>Right Gel Armpad (\$258.00)</b>	
<input type="checkbox"/>	172AE27	Trough / Elbow Block Medium (Gel) (3.5" x 13") Reclining Arm Pad	<input type="checkbox"/> 172AB27
<input type="checkbox"/>	172AE28	Trough / Elbow Block Large (Gel) (4" x 14.5") Reclining Arm Pad	<input type="checkbox"/> 172AB28
<b>Otto Bock Reclining Armpad</b>			
<b>Left Armpad (\$158.00)</b>		<b>Right Armpad (\$158.00)</b>	
<input type="checkbox"/>	172AE08	SM Forearm Pads (8")	<input type="checkbox"/> 172AB08
<input type="checkbox"/>	172AE06	MD Forearm Pads (11.5")	<input type="checkbox"/> 172AB06
<input type="checkbox"/>	172AE07	LG Forearm Pads (14")	<input type="checkbox"/> 172AB07
<b>Otto Bock Hand Pad Options</b>			
<b>Left Hand Pad (\$79.00)</b>		<b>Right Hand Pad (\$79.00)</b>	
<input type="checkbox"/>	172OL01	MD Flat Hand Pads (4-1/2" x 7")	<input type="checkbox"/> 172OR01
<input type="checkbox"/>	172OL02	LG Flat Hand Pads (5-1/2" x 9")	<input type="checkbox"/> 172OR02

ORDERING INFORMATION: PHONE: (800) 263-3390 FAX: (800) 561-5834 WEBSITE: www.sunrisemedical.ca

Specifications are subject to change without notice. Copyright © 2016

MK-100308 REV A (US) REV 7 (CA)



## 15. HANGERS AND FOOTRESTS

### Swing-away Hangers

#### Non-Powered Swing Away Legrests

*Sold as Pairs*

<input type="checkbox"/>	172H2	70° Swing Away (B)	N/C
<input type="checkbox"/>	172H138	Manual ELR S/A (B)	\$352.00

#### Power Elevating Legrests (2 Actuators)

<input type="checkbox"/>	172H126	Swing Away Power ELR w/ Articulation Pair (B footplates)	\$2,994.00
--------------------------	---------	--	------------

#### Footplates For Swing Away Hangers

<input type="checkbox"/>	<b>B - Adult Angle Adjust w/Heel Loop for use with 70° hanger and ELR (Pair)</b>		<b>\$205.00</b>
	172F282	Adult Angle Adjustable Footplate	\$205.00
	172HL10	Heel Loop	N/C
<input type="checkbox"/>	<b>B - Adult Angle Adjust w/Heel Loop for use with 90° hanger (Pair)</b>		<b>\$205.00</b>
	172F285	Adult Angle Adjustable Footplate	\$205.00
	172HL10	Heel Loop	N/C
<input type="checkbox"/>	<b>C - Kids Angle Adjust w/Leg Strap for use with 70° hanger (Pair)</b>		<b>\$205.00</b>
	172F283	Kids Angle Adjustable Footplate	\$205.00
	172LS13	Leg Strap	N/C
<input type="checkbox"/>	<b>C - Kids Angle Adjust w/Leg Strap for use with 80° and 90° hanger(Pair)</b>		<b>\$205.00</b>
	172F286	Kids Angle Adjustable Footplate	\$205.00
	172LS13	Leg Strap	N/C
<input type="checkbox"/>	172F279	Omit Footplate and extension tube	N/C

#### Omit Hangers.

<input type="checkbox"/>	172H144	Omit Hangers	N/C
--------------------------	---------	--------------	-----

#### Hanger Extensions. (See range listed in matrix below)

*NOTE: Not all extension tube will be available; depending on the seat height, some extension tube options may not offer enough ground clearance.*

<input type="checkbox"/>	172E16	High Mount Hanger	\$84.00
<input type="checkbox"/>	172E1	Short	N/C
<input type="checkbox"/>	172E2	Medium	N/C
<input type="checkbox"/>	172E5	Long	N/C

#### Footrest Accessories

<input type="checkbox"/>	1725A01	Toe Loop (Pair) <i>Only avail on (opt B).</i>	\$47.00
<input type="checkbox"/>	172LS14	Extra Leg Strap	\$32.00

### Centermount Footplates

#### Non-Powered Center Mount Legrests

<input type="checkbox"/>	172H128	Fixed Center Mount for Power Recline (F)	N/C
--------------------------	---------	--	-----

#### Power Centermount (2 Actuators)

<input type="checkbox"/>	172H127	<b>Power Centermount ELR w/7.5" Power Articulation (F) Centermount footplates plus Custom Options Below</b>	\$2,987.00
--------------------------	---------	---	------------

#### Power Centermount Footplate Extension Option

<input type="checkbox"/>	172IE01	<b>Enable</b> 7.5" Powered footplate extension for transfer assist. (Electronics programmed to allow footplates to be extended <b>with or without</b> legrest being raised)	\$515.00
<input type="checkbox"/>	172IE02	<b>Disable</b> Footplate extension for transfers. (Electronics programmed to <b>ONLY</b> extend footplates when legs are raised.)	N/C

#### Power Centermount Footplate

<input type="checkbox"/>	<b>F - Centermount Dual Footplate (Pair) (Individually Height (3"-19") and Angle Adjustable)</b>		<b>\$247.00</b>
	172F315	CMT Angle Adjustable Footplates	\$247.00

#### Custom Centermount Options

<input type="checkbox"/>	172DI26	Extra Large Single Footplate for Centermount Footrest. Size 16"W x 14"D	\$309.00
<input type="checkbox"/>	172DI27	Gell Foot Pads for XL Single Centermount Footplate.	\$309.00
<input type="checkbox"/>	172DI28	Extra Large Pair Curved Calf Pads for Center Mount Legrests. Fully adjustable	\$309.00

Lower Leg Length Matrix (measured seat pan to the top of the footplate surface)				
Hanger Type	Standard	Short	Medium	Long
<b>70° Swing Away</b>		13.5" - 16"	16" - 18.5"	18" - 20.5"
<b>Manual ELR S/A</b>		12.5" - 17"	14.5" - 19"	17" - 21"
<b>Power ELR S/A</b>		12" - 14"	14" - 17"	16.5" - 19"
<b>Fixed Cntrmnt Recline</b>	3"-19"			
<b>Power Cntrmnt</b>	3"-19"			

**16. ACCESSORIES**

Wheelchair Tray Table			
<input type="checkbox"/>	1727A01	TRAY TABLE SM TUFFAK 12"-14" PWR	\$331.00
<input type="checkbox"/>	1727A02	TRAY TABLE MED TUFFAK 15"-17" PWR	\$331.00
<input type="checkbox"/>	1727A03	TRAY TABLE LRG TUFFAK 18"-20" PWR	\$331.00

Positioning Belt - Select for any style seat			
<input type="checkbox"/>	172EJ01	Auto Style Buckle	N/C
<input type="checkbox"/>	172EJ02	Belt 2" Aircraft Buckle	\$163.00
<input type="checkbox"/>	172EJ03	Belt 2" Aircraft Buckle Padded	\$178.00

Backpack			
<input type="checkbox"/>	172A1	Black Utility	\$97.00

Accessories			
<input type="checkbox"/>	1728A09	Ram Phone Mount on Round JS tube	\$67.00
<input type="checkbox"/>	1728A10	Ram Tablet Mount on Round JS tube	\$268.00
<input type="checkbox"/>	1728A11	Ram Cup holder	\$67.00
<input type="checkbox"/>	1728A12	Switch-it Remote Stop Switch	\$1,200.00

**NEW!**

**17. ACCESSORIES**



Save time! Wheelchair and seating will be delivered together.  
Hardware pre-installed.

Cushions			
<input type="checkbox"/>	Jay Cushion:		
	Width:	Depth:	
	Notes:		

Backs			
<input type="checkbox"/>	Jay Back:		
	Width:	Height:	
	Notes:		

Headrest	
<input type="checkbox"/>	System:
	Pads:
	Notes:
	Hardware:
	Notes:

Jay Positioning Components	
<input type="checkbox"/>	Anterior Trunk Supports:
	Size:
	Hand and Arm Supports:
	Size:
	Pelvic Positioning:
	Size:
	Foot Positioning:
	Size:
	Notes: