



April 1, 2017

**ACCOUNT INFORMATION**

<b>Dealer Contact Name:</b> _____	<b>Bill to (Name):</b> _____
<b>Email Address:</b> _____	<b>Bill to (Location):</b> _____
<b>Phone Number:</b> _____	<b>Ship to (Name):</b> _____
<b>Fax Number:</b> _____	<b>Address:</b> _____
<b>User Name or Number:</b> _____	<b>City:</b> _____
<b>User Weight:</b> _____	<b>Province:</b> _____
<b>User Height:</b> _____	<b>Postal Code:</b> _____
<b>Seat Cushion Thickness:</b> _____	<b>P.O. # or Quote:</b> _____
<b>Knee to Heel Length:</b> Left: _____ " Right: _____ "	
<b>Seat to Top of Scapula:</b> _____	

The order cannot be processed without the following:  
 - A user name or number, the client's measurements and weight  
 - Sunrise Medical is not responsible for configuration or size discrepancies resulting from customer errors on the order form.  
 - Order confirmations will be provided by Amysystems summarizing the items selected. Please review this carefully.  
 - Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

**STANDARD OPTIONS**

The Quickie Xplore Hybrid Power Base, 6-wheel independent wheel suspension, an additional seat suspension system, tie-down hardware for strap systems, offboard charger, 8" front casters, and standard auto style buckle seat positioning strap, 5.9mph motors, Weight Capacity = 300lbs

**QUICKIE XPLORE STD**

- XPLORESTD2 QUICKIE XPLORE STD - 2 Power Wheelchair \$8,741

**USER WEIGHT PACKAGES**

- WCA3 300 lb Maximum user weight capacity STD

**STANDARD FRAME COLORS**

- FCBLK Black STD
- FCSBM Matte Black STD
- FCGMC Gun Metal (Charcoal) STD
- FCSVM Silver Metallic STD
- FCCBU Candy Blue STD
- FCCRE Candy Red STD
- FCMPPL Midnight Purple STD
- FCRZB Razzleberry STD
- FCMWT New Matte White STD

**OPTIONNAL FRAME COLORS (May Increase lead times)**

- FCTXG Toxic Green N/C
- FCBOR Burnt Orange Metallic N/C
- FCBAO Black Opal Metallic N/C
- FCOPM Sparkle Orange N/C
- FCYLW Yellow N/C
- FCPKP Pearl Pink N/C
- FCHDD 1 Desert Digital Camo Hydro Graphics \$320
- FCHWL 1 Woodland Camo Hydro Graphics \$320
- CUSTOM 1 Carbon Fiber Hydro Graphics (Matte Black base) \$320

Note 1 Includes matching accessories with the Hydro Graphics option. The matching accessories option is not required with these colors

**MATCHING ACCESSORIES TO THE FRAME COLOR**

- FC\*\*\*-M Matching Accessories to the Frame Color (Wheels, front access door, headrest mount and back shell when selected) \$242

**TRANSPORT OPTION**

- TTD4 2 4-Point wheelchair transport brackets N/C

Note 2 4-Point wheelchair transport brackets have not yet received WC-19 crash-testing approval.  
 Use these transport brackets only to secure an unoccupied wheelchair during transport.

**MOTOR PACKAGE**

- MP65 6.5 mph 4-Pole High-Speed Motor Package STD
- MP85 8.5 mph 4-Pole High-Speed Motor Package (Max. user weight 250 lbs - Not available with VR2 electronics ECV902A) N/C

**ASSISTIVE DRIVING PACKAGE**

- INTELLIDR 3, 4 Intelli-Drive \$1,475

Note 3 Intelli-Drive technology provides increased tracking capabilities with enhanced traction and steering precision.  
 This is ideal for those with marginal hand function or switch control.  
 Note 4 Requires the R-net PM 120 Amp Expandable Controller (ECR120)

**BATTERY CHARGER**

- BCOB Offboard Battery Charger N/C
- BCLES8OPT Lester 8 amp Offboard Battery Charger \$361

**BATTERY TRAY AND FRAME SIZE**

- |                          |       |   |  |     |
|--------------------------|-------|---|--|-----|
| <input type="checkbox"/> | FS22B | 5 | Battery Tray for Group 22NF Batteries (Total base width = 23 3/8") | STD |
| <input type="checkbox"/> | FS24B |   | Battery Tray for Group 24 Batteries (Total base width = 25 3/8")   | N/C |

Please Select One Option: Front Caster Arms  Flared out (STD)  Straight (Optional)<sup>6</sup>

**Note 5** The 23 3/8" compact frame is not compatible with most of the Center Mount Front Rigging Options due to Mechanical conflicts

**Note 6** The Straight Front Arm option is not available with any of the Center Mount Front Rigging Options due to Mechanical conflicts

**BATTERIES**

- |                          |           |  |  |       |
|--------------------------|-----------|--|--|-------|
| <input type="checkbox"/> | BAT22NFPR |  | Gel Batteries - Group 22NF (Each - On chair) (E2361) | \$620 |
| <input type="checkbox"/> | BAT24PR   |  | Gel Batteries - Group 24 (Each - On chair) (E2363)   | \$830 |

**DRIVE WHEEL**

- |                          |         |   |   |       |
|--------------------------|---------|---|---|-------|
| <input type="checkbox"/> | DW14PT  |   | 14" Cast Wheel Pneumatic Tire                                   | STD   |
| <input type="checkbox"/> | DW14PI  |   | 14" Cast Wheel Pneumatic with Flat-Free inserts                 | \$242 |
| <input type="checkbox"/> | DW144PT | 7 | 14" x 4" Extra Wide Cast Wheel Pneumatic Tire                   | \$263 |
| <input type="checkbox"/> | DW144PI | 7 | 14" x 4" Extra Wide Cast Wheel Pneumatic with Flat-Free inserts | \$263 |
| <input type="checkbox"/> | DW12PT  | 8 | 12 1/2" Cast Wheel Pneumatic Tire                               | N/C   |
| <input type="checkbox"/> | DW12PI  | 8 | 12 1/2" Cast Wheel Pneumatic with Flat-Free inserts             | \$242 |

The 14" x 4" Extra Wide Wheel increases the overall width of the Base by 2 1/2"

**Note 7** The option to Match the Wheels to the Shroud Color (FC...-M) is not available with the 14" x 4" Extra Wide wheels

The 14" x 4" option is not available with VR2 electronics (ECV902A)

**Note 8** 12 1/2" Drive wheels reduce overall seat to floor height by 1/4", lowers ground clearance and speed, but increases motor torque

**FRONT CASTERS**

- |                          |        |   |   |       |
|--------------------------|--------|---|---|-------|
| <input type="checkbox"/> | FC8PTR |   | 8" Pneumatic Casters                        | STD   |
| <input type="checkbox"/> | FC8PIR |   | 8" Pneumatic Casters with Flat-Free inserts | \$242 |
| <input type="checkbox"/> | FC9PT  | 9 | 9" Pneumatic Casters                        | \$163 |
| <input type="checkbox"/> | FC9PI  | 9 | 9" Pneumatic Casters with Flat-Free inserts | \$295 |

9" front casters increases seat to floor height by 1/2".

**Note 9** 9" casters are NOT compatible with the 12 1/2" Drive wheel configuration (DW12PT and DW12PI)

9" casters are NOT compatible with many front rigging options. Please see details in the front rigging sections on page #6

**SEAT TO FLOOR HEIGHT (Approximate measurement to the top of seat pan ± 1/4")**

**Options for a Seat Only**

- |                          |       |          |   |
|--------------------------|-------|----------|---|
| <input type="checkbox"/> | 10,11 | STFH1525 | 15 1/4" (Possible Seat angle of 0° only)    |
| <input type="checkbox"/> | 11    | STFH16   | 16" (Possible Seat angle of 0° only)        |
| <input type="checkbox"/> | 11    | STFH17   | 17" (Possible Seat angle of 0° and 3°)      |
| <input type="checkbox"/> |       | STFH18   | 18" (Possible Seat angle of 0° and 3°)      |
| <input type="checkbox"/> |       | STFH19   | 19" (Possible Seat angle of 0°, 3° and 6°)  |
| <input type="checkbox"/> |       | STFH20   | 20" (Possible Seat angle of 0°, 3°, 6°, 9°) |

**Options for a Power Tilt or a STD Power Elevating Seat**

- |                          |        |          |  |
|--------------------------|--------|----------|--|
| <input type="checkbox"/> | 10, 11 | STFH1575 | 15 3/4" (Possible Seat angle of 0° only) |
| <input type="checkbox"/> | 11     | STFH165  | 16 1/2" (Possible Seat angle of 0° only) |
| <input type="checkbox"/> |        | STFH175  | 17 1/2" (Possible Seat angle of 0° only) |
| <input type="checkbox"/> |        | STFH18   | 18" (Possible Seat angle of 0° only)     |
| <input type="checkbox"/> |        | STFH185  | 18 1/2" (Possible Seat angle of 0° only) |
| <input type="checkbox"/> |        | STFH19   | 19" (Possible Seat angle of 0° only)     |

**Options for a Combo - Power Tilt with a STD Power Elevating Seat**

- |                          |    |          |   |
|--------------------------|----|----------|---|
| <input type="checkbox"/> | 10 | STFH1775 | 17 3/4" (Possible Seat angle of 0° only) <b>New</b> |
| <input type="checkbox"/> |    | STFH185  | 18 1/2" (Possible Seat angle of 0° only) <b>New</b> |

**Note 10** Requires the 12 1/2" Drive wheels (DW12PT or DW12PI). This lowers ground clearance and speed, but increases motor torque.

**Note 11** The Heavy Duty Swing-Away footrests (SOAAF) are not available with this seat to floor height

**SEAT TO FLOOR HEIGHT SPECIFICATION (If not standard)**

- |                          |        |  |     |
|--------------------------|--------|--|-----|
| <input type="checkbox"/> | SOMHSS | Specific seat to floor height requirement other than standard.<br>Custom upcharges may apply. Required at: _____" (± 1/4" Range) | RFQ |
|--------------------------|--------|--|-----|

**PRE-SET SEAT FRAME ANGLE**

- |                          |       |   |     |
|--------------------------|-------|---|-----|
| <input type="checkbox"/> | SFE-0 | 0° Pre-set seat frame seat angle (Must be selected when using power options)                      | STD |
| <input type="checkbox"/> | SFE-3 | 3° Pre-set seat frame seat angle (Only available with a seat configuration without power options) | N/C |
| <input type="checkbox"/> | SFE-6 | 6° Pre-set seat frame seat angle (Only available with a seat configuration without power options) | N/C |
| <input type="checkbox"/> | SFE-9 | 9° Pre-set seat frame seat angle (Only available with a seat configuration without power options) | N/C |

**LIGHTING PACKAGE**

- |                          |      |    |  |         |
|--------------------------|------|----|--|---------|
| <input type="checkbox"/> | EOLP | 12 | LED Light Package (Includes front beams, side flashers, rear tail light with flashers & brake light) | \$1,066 |
|--------------------------|------|----|--|---------|

**Note 12** Requires the upgraded Joystick with Light Operation or the OMNI (ECVJ6KL, ECELJL, ECRJCL, ECRJCTPL or EDIB)  
The Light Package is not compatible with the EPSEL option.

**NON-EXPANDABLE ELECTRONICS - PG DRIVES VR2 CONTROLLER AND JOYSTICKS**

<input type="checkbox"/>	ECV902A	VR2 90 Amp Drive Controller	STD
<input type="checkbox"/>	ECVJ4K	13 VR2 Standard 4-Key Joystick	STD
<input type="checkbox"/>	ECVJ6K	13 VR2 6-Key Joystick (Allows for a single power option to operate through this joystick)	N/C
<input type="checkbox"/>	ECVJ9KL	13 VR2 9-Key Joystick with Lighting Operation (Allows for a single power function to operate through this joystick)	\$158

**Note 13** Requires the VR2 90A Drive Control (ECV902A).

VR2 Electronics can only accommodate a single power seat function. The only choices are Power Tilt or Elevating seat

**EXPANDABLE ELECTRONICS - PG DRIVES RNET CONTROLLERS**

<input type="checkbox"/>	ECEL90	R-net EL 90 Amp Expandable Controller	\$530
<input type="checkbox"/>	ECR120	R-net PM 120 Amp Expandable Controller	\$1,066
<input type="checkbox"/>	ECRNH	14 Harness for Expandable Controller	N/C

**Note 14** Must be ordered with the expandable controller ECEL90 or ECR120 to interface with Joysticks or OMNI screen

**EXPANDABLE ELECTRONICS - PG DRIVES RNET JOYSTICKS**

<input type="checkbox"/>	ECO	Omit Joystick (For use with an Omni Display)	N/C
<input type="checkbox"/>	ECELJ	15 R-net (LED) Joystick	STD
<input type="checkbox"/>	ECRJJC	15 R-net Joystick with Color screen	\$158

**JOYSTICKS WITH LIGHTING OPERATION**

<input type="checkbox"/>	ECELJL	15 R-net (LED) Joystick with Lighting Operation	\$158
<input type="checkbox"/>	ECRJCL	15 R-net Joystick with Color screen and Lighting Operation	\$210
<input type="checkbox"/>	ECRJCTPL	15 R-net Joystick with Color screen w/ On & Off / Mode switch & speed pot and Lighting Operation	\$315

**ADVANCED JOYSTICK**

<input type="checkbox"/>	ECRJ35 New	15 R-net Advanced Joystick (3.5" Color screen). Lighting Operation, paddle switches, shortcut keys and stereo jacks	\$1,133
--------------------------	------------	---	---------

**Note 15** This option requires the expandable controller ECEL90 or ECR120

**JOYSTICK MOUNT**

<input type="checkbox"/>	JSMF	Fixed Joystick Mount	N/C
<input type="checkbox"/>	JSMSA	Swing Away Joystick Mount - Height adjustable, Angle adjustable and Multi-Axis (E1028)	\$388
		<input type="checkbox"/> Set-up with standard installation to swing-away sideways <input type="checkbox"/> Set-up at 35° to swing-away sideways and downwards	
<input type="checkbox"/>	JSMCM	Midline Mounted Swing Away Joystick Mount (E1028)	\$388

**JOYSTICK MOUNT SUPPORT**

<input type="checkbox"/>	JSMFH	Joystick mount support Bracket	STD
--------------------------	-------	--------------------------------	-----

**JOYSTICK HAND MOUNT**

<input type="checkbox"/>	HML	Install Joystick on Left hand side	<input type="checkbox"/>	HMR	Install Joystick on Right hand side
--------------------------	-----	------------------------------------	--------------------------	-----	-------------------------------------

**JOYSTICK HANDLE**

<input type="checkbox"/>	JSH	Standard Joystick Gimble	STD
<input type="checkbox"/>	JBAT	Stick Handle 3" in Length	\$95
<input type="checkbox"/>	JBH	Ball Handle	\$95
<input type="checkbox"/>	JFBH	Foam Ball Handle	\$95
<input type="checkbox"/>	JTH	T-Handle	\$95
<input type="checkbox"/>	JGP	Goalpost	\$95

**ATTENDANT DRIVE CONTROL OPTION**

<input type="checkbox"/>	EADC	Attendant Hand Control for R-net EL & PM Controls	\$1,066
--------------------------	------	---	---------

**ATTENDANT JOYSTICK HAND MOUNT**

<input type="checkbox"/>	EADCL	Install Attendant Joystick on Rear Left hand side	<input type="checkbox"/>	EADCR	Install Attendant Joystick on Rear Right hand side
--------------------------	-------	---	--------------------------	-------	--

**ENHANCED DISPLAY & DIGITAL INTERFACE BOX (Included with specialty controls)**

<input type="checkbox"/>	EDIB	16 Omni (Enhanced Display & Digital Interface Box with Infrared)	\$2,033
--------------------------	------	--	---------

**Note 16** Allows the switch operation of alternative devices through driver controls.

Requires one of the R-net fully Expandable Controllers (ECEL90 or ECR120)

**ENHANCED DISPLAY MOUNT**

<input type="checkbox"/>	EDSAM	User Swing Away Mount	<input type="checkbox"/>	EDML (Install Display on LH side)	<input type="checkbox"/>	EDMR (Install Display on RH side)	N/C
--------------------------	-------	-----------------------	--------------------------	-----------------------------------	--------------------------	-----------------------------------	-----

**ENVIRONMENTAL CONTROL MODULE**

<input type="checkbox"/>	ECU	17	Environment control module - (Output Module)	\$1,045
--------------------------	-----	----	--	---------

1-Adr-17

**Note 17 Requires one of the fully expandable R-net controllers (ECEL90 or ECR120)**

**BLUETOOTH MOUSE EMULATOR (Includes USB Dongle)**

<input type="checkbox"/>	EDBT	18	Bluetooth Mouse Emulator (USB Dongle included)	\$1,066
<input type="checkbox"/>			Windows and MAC OS	
<input type="checkbox"/>			Android Devices	
<input type="checkbox"/>			Iphone, Ipad and Ipod	

**Note 18 Requires one of the fully expandable R-net controllers (ECEL90 or ECR120)  
Not limited to only one selection. Any two can be installed on the same chair, but only one can be used at a time.**

**POWER OPTIONS & SEATING**

<input type="checkbox"/>	SAABA		Amyseat Rehab Seat - Seat Frame only	\$830
<input type="checkbox"/>	SAPTX	19	CG Power Tilt	\$6,961
<input type="checkbox"/>	SAPRX	19, 21	Power Recline with Integrated Shear Reduction (ISR)	\$8,568
<input type="checkbox"/>	SAPTPRX	19, 21	CG Power Tilt and Power Recline with Integrated Shear Reduction (ISR)	\$12,319
<input type="checkbox"/>	SAPESX	19	11" Power Elevating Seat Module for transfer assist	\$3,452
<input type="checkbox"/>	ADDPESX		<b>Module</b> - Add an 11" Power Elevating Seat to either SAPTX, SAPRX or SAPTPRX	\$3,209
<input type="checkbox"/>	ADDMR	20, 21	<b>Module</b> - Manual Recline with Integrated Shear Reduction (ISR)	\$1,066

**Note 19 These options include the fully adjustable Amyseat (SAABA).**

**Note 20 The maximum user weight capacity with this option is 250 lbs**

**Note 21 It is recommended to use the Contour Back Support option with the power recline (STDSB or STDSB/CUSTOM)**

**ELECTRONICS - SINGLE OR DUAL (1 OR 2 ACTUATORS) POWER SEAT FUNCTION CONTROLLERS**

<input type="checkbox"/>	EPSV	22, 23	Tilt or Elevating seat activation through <b>6-key VR2</b> Hand Control (ECVJ6K or ECVJ9KL)	N/C
<input type="checkbox"/>	EPSEL	22, 23	Tilt or Elevating seat activation through <b>RNET 90 Amp Expandable Controller (ECEL90)</b>	N/C
<input type="checkbox"/>	SAC		Single (1 Actuator) Power Seat Function Controller to Interface through R-net Electronics	N/C
<input type="checkbox"/>	DAC		Dual (2 Actuators) Power Seat Function Controller to Interface through R-net Electronics	\$530

**With the SAC or DAC, please select one activation and position from the list below**

<input type="checkbox"/>	SWTJSS	24, 25	Single actuator activation through R-net Joystick or Specialty Controls	\$1,066
<input type="checkbox"/>	SWTJSD	24, 26	Dual actuator activation through R-net Joystick or Specialty Controls	\$1,066
<input type="checkbox"/>	SWDUALACT		Dual Activation ( <b>Allows to operate all seat functions with external switches AND through the joystick</b> )	\$330
<input type="checkbox"/>	SWTOGGLE	25 or 26	Single Toggle Switch (Single Mode: DAC - Dual Mode: SAC) <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation	N/C
<input type="checkbox"/>	SWPUSH1	25	Single Push Button Switch (Single Mode: SAC) <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation	N/C
<input type="checkbox"/>	SWPUSH2	25 or 26	Dual Push Button Switch (Single Mode: DAC - Dual Mode: SAC) <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation	N/C

**Note 22 Allows for the activation of a single power option through the VR2 or R-net 90 Amp Controllers. Only the Power Tilt or Elevating seat options are available with this feature.**

**Note 23 These options do not require the SAC controller as they connect via the chair's power module directly.**

**Note 24 Requires the R-net 90 Amp (ECEL90) or the 120 Amp (ECR120) Fully Expandable Controllers**

**Note 25 Requires the Single (1 Actuator) Power Seat Function Controller (SAC)**

**Note 26 Requires the Dual (2 Actuators) Power Seat Function Controller (DAC)**

**ELECTRONICS - MULTIPLE (3 to 6 ACTUATORS) POWER SEAT FUNCTION CONTROLLER**

<input type="checkbox"/>	MAC		Multiple (3 to 6 Actuators) Power Seat Function Controller to Interface through R-net Electronics	\$2,564
--------------------------	-----	--	---	---------

**With the MAC, please select one activation and position from the list below**

<input type="checkbox"/>	SWTJSM	27, 28	Multiple Actuator activation through R-net Joystick or Specialty Controls	STD
<input type="checkbox"/>	SWDUALACT		Dual Activation (Allows to operate all seat functions with external switches <b>and</b> through the joystick)	\$330
<input type="checkbox"/>	SWMULTI	28	4-Way Toggle Switch (Up to 4 Activations) - <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation	N/C
<input type="checkbox"/>	SWMULTIT	28	Triple Toggle Switch (Up to 5 Activations) - <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation	N/C
<input type="checkbox"/>	SWPUSH4	28	4-Way Push Button Switch - <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation	N/C
<input type="checkbox"/>	SWPUSH5	28	5-Way Push Button Switch - <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation	N/C

**Note 27 Requires the Multiple (3 to 6 Actuators) Power Seat Function Controller (MAC)**

**Note 28 Requires the R-net 90 Amp (ECEL90) or 120 Amp (ECR120) Fully Expandable Controller, and the multiple Actuator Controller (MAC)**

**SWITCH ACTIVATION HAND MOUNT**

Switch Installation:	<input type="checkbox"/>	TSML (Install Switch on LH side)	<input type="checkbox"/>	TSMR (Install Switch on RH side)	N/C
----------------------	--------------------------	----------------------------------	--------------------------	----------------------------------	-----

**PRE-SET SEAT ANGLE (For Non-Reclining Backs)**

SOBAA Please Select one:  86°  90°  94°  98°  102°  106°  110°  114°  118°  122° N/C

**STANDARD ADJUSTABLE SEAT WIDTH RANGE BETWEEN 15" AND 20"** STD

- SOWN15-20-19 (Width adjustment from 15" to 20" with a maximum Depth of 19" )
- SOWN15-20-22 (Width adjustment from 15" to 20" with a maximum Depth of 22" )<sup>29</sup>

---

- 15" (SOWN15)  18" (SOWN18)<sup>29</sup>
- 16" (SOWN16)<sup>29</sup>  19" (SOWN19S)<sup>29</sup>
- 17" (SOWN17)<sup>29</sup>  20" (SOWN20S)<sup>29</sup>

**Note 29** For the Rehab Seat with the Custom Corbeau Sport Back (SAABA-COR3), please choose the long seat option SOWN15-20-22  
 The only possible Seat width options for the CUSTOM Corbeau Sport Seat (SAABA-COR3) are from 16" to 20" wide

**LARGE ADJUSTABLE SEAT WIDTH RANGE BETWEEN 19" AND 24"** \$462

- SOWN19-24-22 (Width adjustment from 19" to 24" and Depth range of 18" to 22" )

---

- 19" (SOWN19L)  22" (SOWN22)
- 20" (SOWN20L)  23" (SOWN23)
- 21" (SOWN21)  24" (SOWN24)

**ADJUSTABLE SEAT DEPTH RANGE BETWEEN 15" AND 22"**

- 15" (SODN15)  16" (SODN16)  17" (SODN17)  18" (SODN18) STD
- 19" (SODN19)  20" (SODN20)  21" (SODN21)  22" (SODN22) \$315

**BACK CANE HEIGHT**<sup>31</sup> (The back cane height measurement is taken from the seat pan to the top of the push handles)

- 16" (SOBCH16)  20" (SOBCH20)  24" (SOBCH24)
- 18" (SOBCH18)  22" (SOBCH22)<sup>30</sup>  26" (SOBCH26)

**Note 30** 22" is the minimum back cane height required with a power or manual recline.

**Note 31** The Back Cane option is not compatible for use with the Rehab Seat with the CUSTOM Corbeau Sport Back (SAABA-COR3)

**BACK SUPPORT**

- STDSB<sup>32</sup> Contour back support with a maximum cushion height of 20" (15" to 20" Wide) \$558
- STDSB/CUSTOM Contour back support (Non-Standard Sizes) Specify: Width: Height: \$1,256

**Note 32** It is recommended to use the Contour Back Support option with the power or manual recline. Not required for Corbeau

**<sup>34</sup> BACK SUPPORT TOTAL HEIGHT (MEASURED FROM THE SEAT PAN TO THE TOP OF THE BACK)**

- 16" (SOBH16)  20" (SOBH20)  24" (SOBH24)  28" (SOBH28)
- 18" (SOBH18)  22" (SOBH22)<sup>33</sup>  26" (SOBH26)

**Note 33** With a power or manual recline, the minimum height is 22".

**Note 34** The space between the bottom of the back and the seat pan is 4" with a non-reclining back and 6" with a reclining back

**SEAT CUSHION**

- SCPSC Contour Seat Cushion with a maximum seat depth of 22" (15" to 20" Wide) \$588
- SCPSC/CUSTOM Contour Seat Cushion (Non-Standard Sizes) Specify: Width: Depth: \$1,281

**CANTILEVER - HEIGHT ADJUSTABLE ARMRESTS - OPENS WITH A 5° ANGLE (For Non-Reclining Backs - Choose 2)**

- SOACORFLH Corbeau Flip-Back Cantilever - Full Length Left Side (Height Adjustable from 9 ½" - 13 ½") STD
- SOACORFRH Corbeau Flip-Back Cantilever - Full Length Right Side (Height Adjustable from 9 ½" - 13 ½") STD
- SOACORDLH Corbeau Flip-Back Cantilever - Desk Length Left Side (Height Adjustable from 9 ½" - 13 ½") N/C
- SOACORDRH Corbeau Flip-Back Cantilever - Desk Length Right Side (Height Adjustable from 9 ½" - 13 ½") N/C
- SOAFBFLH Flip-Back Cantilever - Full Length Left Side (Height Adjustable from 7" to 15" in 1" increments) STD
- SOAFBFRH Flip-Back Cantilever - Full Length Right Side (Height Adjustable from 7" to 15" in 1" increments) STD
- SOAFBDLH Flip-Back Cantilever - Desk Length Left Side (Height Adjustable from 7" to 15" in 1" increments) N/C
- SOAFBDRH Flip-Back Cantilever - Desk Length Right Side (Height Adjustable from 7" to 15" in 1" increments) N/C

**DUAL POST - HEIGHT ADJUSTABLE ARMRESTS - INCLUDES TRANSFER GRIPS (For Non-Reclining Backs - Choose 2)**

- SOADPFLH Dual Post Height Adjustable - Full Length Left Side (Height Adjustable from 9 ½" - 13 ½") \$190
- SOADPFRH Dual Post Height Adjustable - Full Length Right Side (Height Adjustable from 9 ½" - 13 ½") \$190
- SOADPDLH Dual Post Height Adjustable - Desk Length Left Side (Height Adjustable from 9 ½" - 13 ½") \$190
- SOADPDRH Dual Post Height Adjustable - Desk Length Right Side (Height Adjustable from 9 ½" - 13 ½") \$190

**RECLINING FLIP-UP, HEIGHT ADJUSTABLE ARMRESTS (For Reclining Backs - Choose 2)**

- SOANIFLH Reclining armrest - Full Length Left Side (Height Adjustable from 10 ½" to 14" in ½" increments) N/C
- SOANIFRH Reclining armrest - Full Length Right Side (Height Adjustable from 10 ½" to 14" in ½" increments) N/C
- SOANIDLH Reclining armrest - Desk Length Left Side (Height Adjustable from 10 ½" to 14" in ½" increments) N/C
- SOANIDRH Reclining armrest - Desk Length Right Side (Height Adjustable from 10 ½" to 14" in ½" increments) N/C

## ARMREST HEIGHT ADJUSTMENT AT FACTORY

Please set the Armrest Height at: \_\_\_\_\_ " (The Armrest Height is taken from the seat pan to the top of the armpad. Not required for Dual Post)

### ARMPADS

#### Foam Type Selection

<input type="checkbox"/>	SOFAPD	Desk Length Standard (10")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	STD
<input type="checkbox"/>	SOFAPF	Full Length Standard (14")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	STD
<input type="checkbox"/>	SOWAPD	36 Desk Length Waterfall (10")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$32
<input type="checkbox"/>	SOWAPF	36 Full Length Waterfall (14")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$32

#### Hybrid Gel and Foam Type Selection (The Gel always stays soft)

<input type="checkbox"/>	SOGELD new	Desk Length Hybrid Gel and Foam Flat armpad (10" x 2.5")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$32
<input type="checkbox"/>	SOGELF	Full Length Hybrid Gel and Foam Flat armpad (14" x 2.5")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$32
<input type="checkbox"/>	SOWGELD new	36 Desk Length Hybrid Gel and Foam Waterfall armpad (10")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$50
<input type="checkbox"/>	SOWGELF new	36 Full Length Hybrid Gel and Foam Waterfall armpad (14")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$50

#### Specialty Type Selection

<input type="checkbox"/>	SOCAR	37 Arm Trough with Hand Support (18" x 4 1/2")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$210
--------------------------	-------	--	---	--------------	---	--	-------

**Note 36** The Waterfall armpads may interfere with flipping back the armrests when combined with Lateral Trunk supports

**Note 37** Arm Trough Armpads may interfere with Armrests Flipping back.

### FOOTRESTS

<input type="checkbox"/>	SOSAAF	41 Adult Swing-Away Footrests -----	Set at:	<input type="checkbox"/> 60° <sup>42</sup>	<input type="checkbox"/> 70°	\$263
Select Footrest Range: <input type="checkbox"/> SOFLA10 (10 1/2" - 15 1/2") <input type="checkbox"/> SOFLA14 (14 1/2" - 19 1/2")						
<input type="checkbox"/>	SOAAF	41 Adult <b>Heavy Duty</b> Swing-Away footrests -----	Set at:	<input type="checkbox"/> 60° <sup>42</sup>	<input type="checkbox"/> 70°	\$363
Select Footrest Range: <input type="checkbox"/> SOFLHD12 (12 1/2" - 15 1/2") <input type="checkbox"/> SOFLHD14 (14 1/2" - 19 1/2")						
<input type="checkbox"/>	ADDPL	38, 41 Adult Power Elevating & Articulating Legrests	Please Select Function:	<input type="checkbox"/> Individual <sup>39</sup>	<input type="checkbox"/> Combined <sup>40</sup>	\$2,994
Select Footrest Range: <input type="checkbox"/> SOFLPW12 (12 1/2" - 15 1/2") <input type="checkbox"/> SOFLPW14 (14 1/2" - 19 1/2")						

41 Please set the Footrest Length at: \_\_\_\_\_ " (The footrest length is taken from the seat pan to the heel)

**Note 38** This option includes calf pads

**Note 39** Individual legs count as two power functions.

**Note 40** Combined legs count as one power function.

**Note 41** The Footrest Length is taken from the top of the seat pan to the top of the footplate.

The 70° footrest angle is not compatible with the 9" front casters (FC9PT and FC9PI)

**Note 42** The 70° footrest angle may interfere with the casters when selecting the straight front caster arms option

### FOOTPLATES & ACCESSORIES

<input type="checkbox"/>	SOSFAF	Flip-up Footplates	STD
<input type="checkbox"/>	SOAFP	Angle Adjustable & Flip-up Footplates - Each	\$137
<input type="checkbox"/>	SOAHL	Heel Loops - Each	\$41
<input type="checkbox"/>	SOACPFR	Calf Pads for Adult Heavy-Duty Swing-Away footrests ( <b>SOAAF</b> ) with adjustable hardware - Each	\$137

### CENTER MOUNTS

<input type="checkbox"/>	SOAAF90	45, 46, 47, 48, 49, 50 90° Fixed Center Mount Interface Bracket ( <b>Offers Manual adjustability from 95° to 70°</b> )	\$242
Select Foot Platform Range: <input type="checkbox"/> SOFBA12 (12 1/2" - 16 1/2") <input type="checkbox"/> SOFBA16 (16 1/2" - 19 1/2")			
<input type="checkbox"/>	ADDPF	43, 44, 45, 47, 48, 49 90° Power Center Mount Elevating & Articulating Foot Platform - *** <b>Not available with the 12 1/2" Drive Wheel</b>	\$2,994
<input type="checkbox"/> SOFBPW12 (12 1/2" - 16 1/2") <input type="checkbox"/> SOFBPW16 (16 1/2" - 20 1/2")			

48 Please set the Foot Platform Length at: \_\_\_\_\_ " (The foot Platform length is taken from the seat pan to the heel)

**Note 43** The power center mount counts as one power function.

**Note 44** This option includes calf pads

**Note 45** Links and attaches the one piece center mount foot platform (SOFUFP) or the split footplates (SOAFPFB) to the seat

**Note 46** With the split footplates option (SOAFPFB), this bracket offers independent height adjustability for each footplate

**Note 47** Please Select the one piece center mount foot platform (SOFUFP) or the split footplates (SOAFPFB)

**Note 48** The Foot Platform Length is taken from the top of the seat pan to the top of the foot platform or split foot plates.

**Note 49** Center Mount Front Riggings are not available with the Straight Front Arm option due to Mechanical conflicts

**Note 50** This option is not compatible with the 9" Caster option

### FOOT PLATFORMS & ACCESSORIES

<input type="checkbox"/>	SOFUFP	Flip-up one piece Foot Platform	N/C
<input type="checkbox"/>	SOAFPFB	Angle Adjustable & Flip-up Split Footplates for center mount - Each	\$137
<input type="checkbox"/>	SOACPCM	Calf Pads for 90° Fixed Center Mount ( <b>SOAAF90</b> ) with adjustable hardware - Each	\$137

### SEAT POSITIONING

<input type="checkbox"/>	SOAUTO	55" Seat Positioning Strap		STD
--------------------------	--------	----------------------------	--	-----

1-Apr-17

### HEADREST <sup>51</sup>

<input type="checkbox"/>	SONUAHR	Using headrest mounting hardware other than SOMAH located below. Specify: _____		
<input type="checkbox"/>	SOFHRS	Small Curved Headrest Pad 4" x 9"		\$215
<input type="checkbox"/>	SOFHRM	Medium Curved Headrest Pad 5" x 10"		\$215
<input type="checkbox"/>	SOFHRL	Large Curved Headrest Pad 6" x 11"		\$215
<input type="checkbox"/>	SOMAH	Multi-axis Removable Mounting Hardware		\$173
<input type="checkbox"/>	SOHRBC	Bracket to Mount on the Contour Back Support <b>(Without Recline)</b>		\$105
<input type="checkbox"/>	SOHRPR	Bracket to Mount on the Seat Frame <b>(For Reclining Backs, but excludes the Corbeau Back)</b>		\$105

\* Please select one of the headrest options with the mounting hardware and specify where the headrest will be installed

Note 51 The Headrest option is not compatible with the Rehab Seat with the Custom Corbeau Sport Back option (SAABA-COR3)

### LATERAL TRUNK SUPPORTS <sup>52</sup>

<input type="checkbox"/>	SOSLS3	Adjustable Laterals with Small Curved Pads (Each) - 3" x 5"		\$214
<input type="checkbox"/>	SOSLM4	Adjustable Laterals with Medium Curved Pads (Each) - 4" x 6"		\$214
<input type="checkbox"/>	SOSLL5	Adjustable Laterals with Large Curved Pads (Each) - 5" x 7"		\$214
<input type="checkbox"/>	SOSLHM	Adjustable Swing-Away Hardware Mounts - Each		\$316

\*\* Please select one of the lateral pads and the adjustable & removable hardware.

Note 52 The Lateral Trunk Support option is not compatible with the Rehab Seat with the Custom Corbeau Sport Back option (SAABA-COR3)

### HIP & THIGH SUPPORTS

<input type="checkbox"/>	SOHTSS10	Adjustable Hip & Thigh Support Pads (Each) - Small 10" x 4"		\$242
<input type="checkbox"/>	SOHTSM12	Adjustable Hip & Thigh Support Pads (Each) - Medium 12" x 5"		\$242
<input type="checkbox"/>	SOHTSL14	Adjustable Hip & Thigh Support Pads (Each) - Large 14" x 6"		\$242
<input type="checkbox"/>	SOHTSARH	Adjustable and Removable Hip & Thigh Support Hardware - Each		\$457

\*\*\* Please select one of the hip & thigh support pads and the adjustable & removable hardware.

### ACCESSORIES

<input type="checkbox"/>	SOCH	Cup Holder (Self-Leveling)		\$95
<input type="checkbox"/>	SOCPH	Ram Mounts X-Grip Cell Phone Holder		\$163
<input type="checkbox"/>	SOTHS	Ram Mounts X-Grip Small Tablet Holder		\$184
<input type="checkbox"/>	SOTHL	Ram Mounts X-Grip Large Tablet Holder with Swivel		\$278

### REFLECTOR KIT

<input type="checkbox"/>	REFLECT	Red & White Reflector Kit Stickers (They will not be installed - Included in a bag with the owner's manual)		N/C
--------------------------	---------	---	--	-----

### SWITCHES AND OPTIONS

<input type="checkbox"/>	ASCPCP	PGDT RNET Dongle and PC programmer kit		\$615
<input type="checkbox"/>	ASCPVR2	PGDT VR2 Dongle and PC programmer kit		\$615
<input type="checkbox"/>	SSBS	Buddy Button Switch <b>(K0108)</b>		\$115
<input type="checkbox"/>	SSES	Egg Switch-Black <b>(K0108)</b>		\$115
<input type="checkbox"/>	SSML	Micro Light Switch <b>(K0108)</b>		\$158

## SPECIALTY DRIVER CONTROL OPTIONS

### Therafin Sip-N-Puff Kit <sup>60</sup>

<input type="checkbox"/>	SCSNPK	53, 54	Therafin Sip-N-Puff Kit		\$1,603
<input type="checkbox"/>	SCSNPBTk	54	Therafin Sip-N-Puff Breath Tube Kit		\$520
<input type="checkbox"/>	SCSRMSNP	54	Stealth Rock-n-Lock Swingaway Arm for Sip-N-Puff	<input type="checkbox"/> SWL (Swings on Left) <input type="checkbox"/> SWR (Swings on Right)	\$289
<input type="checkbox"/>	SCSRM	54	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/> SWL (Swings Left) <input type="checkbox"/> SWR (Swings Right)	\$520
<input type="checkbox"/>	SSES-STD	54, 56	Egg Switch for Changing Modes - Black	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$115
<input type="checkbox"/>	SSBS-STD	54	Buddy Button Switch (Stop Switch) - Red		\$115
<input type="checkbox"/>	SOMAH	54	Multi-axis Headrest Removable Mounting Hardware		\$173
<input type="checkbox"/>	SOFHRS		Small Curved Headrest Pad 4" x 9"		\$215
<input type="checkbox"/>	SOFHRM	55	Medium Curved Headrest Pad 5" x 10"		\$215
<input type="checkbox"/>	SOFHRL		Large Curved Headrest Pad 6" x 11"		\$215

### Switch-It Proximity Head Array <sup>60</sup>

<input type="checkbox"/>	SCSIHAA	53, 54	Adult - 3 Switch Proximity Head Array with 1/8" Mono Jack for mode/reverse		\$5,353
<input type="checkbox"/>	SCSRM	54	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/> SWL (Swings Left) <input type="checkbox"/> SWR (Swings Right)	\$520
<input type="checkbox"/>	SSES-STD	54, 56	Egg Switch for Changing Modes - Black	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$115
<input type="checkbox"/>	SSBS-STD	54	Buddy Button Switch (Stop Switch) - Red		\$115
<input type="checkbox"/>	SOMAH	54	Multi-axis Headrest Removable Mounting Hardware		\$173

### Switch-It MicroPilot Light Touch Proportional Mini joystick for chin / lip control <sup>60</sup>

<input type="checkbox"/>	SCSIMP	53, 54	MicroPilot Light touch mini Proportional Joystick		\$4,197
<input type="checkbox"/>	SCSRMMP	54	Stealth Rock-n-Lock Swingaway Mount for Chin M-P	<input type="checkbox"/> SWL (Swings Left) <input type="checkbox"/> SWR (Swings Right)	\$520
<input type="checkbox"/>	SCSRM	54	Stealth Rock-n-Lock Swingaway Mount for Switch	<input type="checkbox"/> SWL (Swings Left) <input type="checkbox"/> SWR (Swings Right)	\$520
<input type="checkbox"/>	SSES	54, 56	Egg Switch-Black	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$115
<input type="checkbox"/>	JCC		Chin Cup		\$55
<input type="checkbox"/>	SOMAH	54	Multi-axis Headrest Removable Mounting Hardware		\$173
<input type="checkbox"/>	SOFHRS		Small Curved Headrest Pad 4" x 9"		\$215
<input type="checkbox"/>	SOFHRM	55	Medium Curved Headrest Pad 5" x 10"		\$215
<input type="checkbox"/>	SOFHRL		Large Curved Headrest Pad 6" x 11"		\$215

### Proportional Compact or Mini joystick for Midline Mount

<input type="checkbox"/>	SCPGMJS	53, 54	PG Drives JC200 Proportional Compact Remote Joystick for Midline Mount		\$998
<input type="checkbox"/>	SCGMJSM	57	Gatlin Midline Swing-away Mount for JC200	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$410
<input type="checkbox"/>	SCSIMP	53, 54	Switch-It MicroPilot Light touch mini Proportional Joystick		\$4,197
<input type="checkbox"/>	SCGMPPM	58	Gatlin Midline Swing-away Mount for MicroPilot	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$447
<input type="checkbox"/>	SCBTM	59	Bullet Hand Tray for MicroPilot	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$473

### Proportional Compact or Mini joystick for Armrest Mount

<input type="checkbox"/>	SCPGMJS	55, 56	PG Drives JC200 Proportional Compact Remote Joystick for Armrest Mount		\$998
<input type="checkbox"/>	SCARMJSM	57	Swing-away and Multi-Axis armrest mount	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$410
<input type="checkbox"/>	SCSIMP	53, 54	Switch-It MicroPilot Light touch mini Proportional Joystick		\$4,197
<input type="checkbox"/>	SCARMPPM	58	Swing-away and Multi-Axis armrest mount	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$447
<input type="checkbox"/>	SCBTM	59	Bullet Hand Tray for MicroPilot	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$473

### Proportional RIM Head Control Package <sup>60</sup>

<input type="checkbox"/>	SCPHC	53, 54	Proportional RIM Head Control Package		\$2,621
<input type="checkbox"/>	SCSRM	54	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/> SWL (Swings Left) <input type="checkbox"/> SWR (Swings Right)	\$520
<input type="checkbox"/>	SSES-STD	54, 56	Egg Switch for Changing Modes - Black	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$115
<input type="checkbox"/>	SOMAH	54	Multi-axis Headrest Removable Mounting Hardware		\$173

Note 53 Requires the Omni Enhanced Display & Digital Interface Box (EDIB)

Note 54 This option must be selected for the Kit

Note 55 Choose one (1) Headrest Pad

Note 56 The Mode Switch Type can be changed from this selection on Page #7

Note 57 Must be selected for use with the PG Drives JC200 Mini Joystick

Note 58 Must be selected for use with the Switch-It MicroPilot

Note 59 This option is only compatible with the Switch-It MicroPilot

Note 60 This option is not compatible for use with the Rehab Seat with the Custom Corbeau Sport Back (SAABA-COR3)



**61 CUSTOM SECTION**

**Note 61** Lead times will increase depending on the custom order pipeline at the time the purchase order is received  
Please call for Quote and specifications for these options.

**O<sup>2</sup> Holder**

<input type="checkbox"/>	SUPOXY	62, 63	Custom Non-Articulating Oxygen Tank Holder	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	\$410
<input type="checkbox"/>	SUPOXYA	62, 64	Custom Articulating Oxygen Tank Holder	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	\$1,923
			<input type="checkbox"/> D-Size Oxygen Tank (* 22" Tall x 4.3" Diameter) * The height includes the valve			
			<input type="checkbox"/> E-Size Oxygen Tank (* 31" Tall x 4.3" Diameter) * The height includes the valve			

**Note 62** The O<sup>2</sup> Holder option is not compatible for use with the Rehab Seat with Corbeau Sport Back (SAABA-COR3)

**Note 63** Not available with a power tilt and/or power recline

**Note 64** Required with a power tilt and/or power recline

**VENT TRAY (Please complete in Full)**

<input type="checkbox"/>	SOVENT		Custom Vent Tray			\$2,574
		65, 66	<input type="checkbox"/> Articulating Vent Box (Required with Power Tilt and/or Recline)	Model: _____		
		65, 66	<input type="checkbox"/> Fixed Vent Box	Model: _____		
<input type="checkbox"/>	SOVTBE	65, 66, 67	<input type="checkbox"/> Extended Battery Compartment			N/C
<input type="checkbox"/>	The Vent is being powered by 1 U1 Battery			<input type="checkbox"/> Chair Powered by 2 Group 22NF Batteries	<input type="checkbox"/> Chair Powered by 2 Group 24 Batteries	
<input type="checkbox"/>	The Vent is being powered by 1 Group 22NF Battery			<input type="checkbox"/> Chair Powered by 2 Group 22NF Batteries		
<input type="checkbox"/>	The Vent is being powered by 1 Group 24 Battery			<input type="checkbox"/> Chair Powered by 2 Group 22NF Batteries		

**Note 65** Amysystems does not supply the battery, charger or the wiring for the Vent.  
These items are normally provided by the Vent Manufacturer.

**Note 66** The Vent Tray option is not compatible for use with the Rehab Seat with Corbeau Sport Back (SAABA-COR3)

**Note 67** This option is not available with the Straight Front Caster Arms

**OTHER CUSTOM OPTIONS (Request For Quote)**

<input type="checkbox"/>	PAT	68, 70	Power Anterior Tilt			RFQ	
<input type="checkbox"/>	CUSTOM	68	Power Lateral Tilt 15°	<input type="checkbox"/> Raise the Left Side only	<input type="checkbox"/> Raise the Right Side only	<input type="checkbox"/> Raise the Left & Right Side	RFQ
<input type="checkbox"/>	CUSTOM	68, 70	Power Elevating Stretcher	(Please call CSR for availability with configuration and user weight capacity)		RFQ	
<input type="checkbox"/>	CUSTOM	68, 70	Power Flip-up Foot Platform	(Please call CSR for availability with configuration and user weight capacity)		RFQ	
<input type="checkbox"/>	CUSTOM	68, 72	Power Vertical Lifting Foot Platform	(Please call CSR for availability with configuration and user weight capacity)		RFQ	
<input type="checkbox"/>	CUSTOM	68	Independent Depth Adjustability built into the Footrest Receiver that is attached to the Seat Pan			RFQ	
<input type="checkbox"/>	CUSTOM	71, 72, 73	Rehab Seat with Corbeau Performance Back - Includes a manual Recline Lever	<input type="checkbox"/> Lever LH	<input type="checkbox"/> Lever RH	RFQ	
<input type="checkbox"/>	CUSTOM	68	Endomorphic Back Canes (Pair)	<input type="checkbox"/> 1" per side towards the inside = 2" in	<input type="checkbox"/> 2" per side towards the inside = 4" in	RFQ	
<input type="checkbox"/>	CUSTOM		Transfer Handles (Pair)			RFQ	
<input type="checkbox"/>	CUSTOM		Manual Wheel Locks			RFQ	
<input type="checkbox"/>	CUSTOM	68	Non Standard Frame Colors (Not all custom colors are available, please contact CSR for availability)			RFQ	
<input type="checkbox"/>	CUSTOM # 1	68	Other custom options. Specify in notes section below.			RFQ	
<input type="checkbox"/>	CUSTOM # 2	68	Other custom options. Specify in notes section below.			RFQ	

**Note 68** Please specify details and specifications for these options.

**Note 69** This option is not available with the Rehab Seat and the CUSTOM Corbeau Sport Back (SAABA-COR3)

**Note 70** The maximum user weight capacity with this option is 250 lbs

The Corbeau Back is not available with the Power Recline options, but includes a manual Recline. However, the 50° Power Tilt,

**Note 71** the 11" Power elevating seat module and the power elevating legs or center mount options are compatible.

The Corbeau Sport Back option only comes in one Size, but offers a seat width range from 16" to 20"

**Note 72** This option is not available with the Straight Front Caster Arms

**Note 73** The Custom Corbeau Seat is not available with VR2 Electronics (ECV902A)

**NOTES / SPECIAL REQUESTS / CUSTOM SPECIFICATIONS**