



Quickie QM-7 Series Single Powered Tilt Seating

March 2017
EIPW11

Account: _____
Date: _____
PO#: _____
Buyer: _____
Marked For: _____

Ship To: _____
Name: _____
Address: _____
City, Province, Postal Code: _____
Phone: _____

Option #	Definition	Price	Option #	Definition	Price
1. QUICKIE QM-7 Series Single-Power Tilt Seating (Select One)					

QUICKIE QM-710™ SC

<input type="checkbox"/> Quickie QM-710™ SC for ASAP II Tilt(Fig 1.1)		
172PC68	Base Package	\$7,349.00

QUICKIE QM-715™ HD SC

<input type="checkbox"/> Quickie QM-715™ HD SC for ASAP II Tilt		
172PC80	Base Package	\$7,885.00

QUICKIE QM-720™ SD

<input type="checkbox"/> Quickie QM-720™ SD for ASAP II Tilt (Fig 1.1)		
172PC86	Base Package	\$8,137.00

SC Seating Specifications

Typical Chair Weight	W/battery	350lbs	W/O Battery	260lbs
Transit	ISO occupied approved			
Seat Type	ASAP II			
Seat size range	Width	12"-22"	Depth	12"-22"
STFH range	ASAP II	16.6" -20.6"		
Power Seating Options	<ul style="list-style-type: none"> - 50° CG tilt - Power Swing away Elevating/Articulating Legrest - Power Elevating/Articulating Centermount - 12" Seat elevator 			



Actuators Required by Power Seating Function

50° CG tilt	1 actuator
12" Seat Elevate:	1 actuator
Power Leg Function (Centermount or Swing away):	2 actuators
Combination : Tilt and Seat Elevate:	2 actuators
Combination : Tilt and Power Elevating Legrest. (Centermount or Swing away):	3 actuators
Combination : Tilt, Elevate and Power Legrest (Centermount or Swing away):	4 actuators

QUICKIE QM-710™ Base Specifications

Drive Wheel Position	Midwheel Drive	Drive Wheels	14" Drive Wheels w/ flat free insert (Black or Grey)
User weight Cap	300 lb - 400lbs	Tires/Casters	6" Front / 6" rear (Black or Grey)
Base size	Width 25" Length 36" Turn Radius 20"	Charger	8A Dual Mode Charger (Compact) CTE
Drive Control	R-Net PM120 (Non-Expandable and Expandable)	Battery Size	Group 24 / Group 34.
Driving Range	Up to 13.5 miles	Shroud Colors	12 standard colors, 6 upgrade prints.
Drive Motors	4-pole 6 mph, Optional 5mph and 8mph	Ground Clearance	4"

2. BASE CHAIR OPTIONS

Speed Package (4-pole)			
<input type="checkbox"/>	172PF85	QM-710 Standard 6 mph motor package	N/C
<input type="checkbox"/>	172PF91	QM-710 5 mph High Torque motor package	\$525.00
<input type="checkbox"/>	172PF90	QM-710 8 mph Community motor package	\$1,103.00
<input type="checkbox"/>	172PF84	QM-715 HD Standard 5 mph motor package	N/C
<input type="checkbox"/>	172PF89	QM-720 Standard 8 mph motor package	N/C

Shroud Color			
<input type="checkbox"/>	172EZ01	Candy Apple Red	N/C
<input type="checkbox"/>	172EZ02	High Gloss Black	N/C
<input type="checkbox"/>	172EZ03	Candy Blue	N/C
<input type="checkbox"/>	172EZ04	Silver	N/C
<input type="checkbox"/>	172EZ05	Copper	N/C
<input type="checkbox"/>	172EZ08	Green Apple	N/C
<input type="checkbox"/>	172EZ09	Yellow	N/C
<input type="checkbox"/>	172EZ10	Pearl Pink	N/C
<input type="checkbox"/>	172EZ11	Candy Purple	N/C
<input type="checkbox"/>	172EZ12	Matte Black	N/C
<input type="checkbox"/>	172EZ13	Midnight Blue	N/C
<input type="checkbox"/>	172EZ15	Pearl White	N/C
<input type="checkbox"/>	172EZ06	Digital Camo.	\$525.00
<input type="checkbox"/>	172EZ07	Woodland Camo.	\$525.00
<input type="checkbox"/>	172EZ14	Carbon Fiber Print	\$525.00
<input type="checkbox"/>	172EZ16	American Flag.	\$525.00
<input type="checkbox"/>	172EZ17	Pink Camo	\$525.00
<input type="checkbox"/>	172EZ18	Grunge Skulls	\$525.00

User weight			
<input type="checkbox"/>	172WG06	Less than 150lbs	N/C
<input type="checkbox"/>	172WG13	151-200lb	N/C
<input type="checkbox"/>	172WG14	201-250lb	N/C
<input type="checkbox"/>	172WG03	251-300lb	N/C
<input type="checkbox"/>	172WG15	301-350lb	N/C
<input type="checkbox"/>	172WG16	351-400lb	N/C

Transit Package NA

Drive Wheel			
<input type="checkbox"/>	172DW59	14" Cast Wheel Solid Dual Density Inserts Black Knobby Tires and Black Caster	Std

Suretrac Drive Tracking NEW!			
<input type="checkbox"/>	172MI01	Suretrac Module and Motor Encoders	\$1,475.00


Transit Option			
<input type="checkbox"/>	172TO25	4-Point Transit Tie-down Package Meets ISO-7176 for Occupant Use	\$373.00
<input type="checkbox"/>	172TO40	Omit Transit Tie-down	NC

Wheel locks			
<input type="checkbox"/>	172WL42	Wheel locks	\$142.00

Ventilator Options			
<input type="checkbox"/>	172VT4	Fixed Vent Tray (Pulmonetic LTV Models)	\$525.00

Pre-Set Seat Frame Angle			
<input type="checkbox"/>	172FA10	0°	N/C

Seat Height (Use STFH Matrix below)			
(Measured in the front) - Obtain Quote to verify est. STFH			
<input type="checkbox"/>	172HT44	Position "0"	N/C
<input type="checkbox"/>	172HT49	+5" STFH	Requires Group 34 Batteries
<input type="checkbox"/>	172HT45	+1" STFH	N/C
<input type="checkbox"/>	172HT46	+2" STFH	N/C
<input type="checkbox"/>	172HT47	+3" STFH	N/C

LED Light Package			
<input type="checkbox"/>	172LI05	 ON-CHAIR LED Light Package: Includes headlights, tail lights, turn indicators and flashers.	\$1,066.00

STFH Matrix					
Seating Type and Power Option Combinations	Seat-to-Floor Height at 0 tilt				
	0	+5"	+1"	+2"	+3"
SC ASAP II Tilt only	16.6	17.1	17.6	18.6	19.6
SC ASAP II Tilt and 12" Lift combo.	16.7	17.2	17.7	18.7	
SC ASAP II HD Tilt only	16.8	17.3	17.8	18.8	19.8

STFH measured with batteries but no user. Actual STFH may vary by .5"
 *Adding seat angle will increase STFH by 1" for each incremental change.
 **Ultra low STFH may cause ground clearance of Power Centermount to drop below 2".

3. BATTERIES and CHARGERS

Batteries			
<input type="checkbox"/>	172AG04	MK Gr 24 Flat Top Gel Cell (Pair) w/straps Recommended for Max Daily Range	\$830.00
<input type="checkbox"/>	172AG11	MK Gr 34 Flat Top Gel Cell (Pair) w/straps, Required for lower STFH.	\$830.00

Battery Charger			
<input type="checkbox"/>	172BJ19	8A Dual Mode (Compact) CTE	N/C

4. ELECTRONICS

Drive Control Options			
Non-Expandable Controllers			
<input type="checkbox"/>	172ER05	R-Net Non-Expandable Drive Control 120A (For Use w/ 2 Thru-drive seat Functions)	N/C

Expandable Controllers			
<input type="checkbox"/>	172ER06	R-Net Expandable Drive Control (120A)	\$1,066.00
	172EX05	Expandable 120A Controller - R-net	\$1,066.00
	172EX06	Harness for 120A Expandable Controller	NC

R-net Joystick			
<input type="checkbox"/>	172J70	Std LED Joystick for Rnet Controls	N/C
<input type="checkbox"/>	172J108	Joystick w/Stereo Jacks & Color Display	N/C
<input type="checkbox"/>	172J88	Joystick W/Toggle, Spd Pot, Stereo Jack & Color Display	\$158.00
<input type="checkbox"/>	172J112	Joystick Advanced: 3.5" Color Display, Paddle Switches, Shortcut Keys and Stereo Jack.	\$1,133.00
<input type="checkbox"/>	172J95	Omit JS. For use with OMNI	N/C

Joystick Mount(Fixed, Swing-Away, Retractable, Midline)			
<input type="checkbox"/>	172JM10	Fixed Mount	N/C
<input type="checkbox"/>	172JM11	Swing Away Mount	\$388.00
<input type="checkbox"/>	172JM15	Retractable Joystick mount w/45° stops	\$721.00
<input type="checkbox"/>	172JM14	Mid Line Mount	\$388.00

Joystick Mount (Location)			
<i>Height Adjustable.</i>			
<input type="checkbox"/>	172HM1	Right Side Mount	N/C
<input type="checkbox"/>	172HM2	Left Side Mount	N/C

Joystick Handle (Optional)			
<input type="checkbox"/>	172JH1	Ball Handle	\$95.00
<input type="checkbox"/>	172JH2	Foam Ball Handle	\$95.00
<input type="checkbox"/>	172JH4	T-Handle	\$95.00
<input type="checkbox"/>	172JH6	Goalpost	\$95.00

Attendant Control			
<input type="checkbox"/>	172A04	Attendant Control for R-net	\$1,066.00

Attendant Mounting location			
<input type="checkbox"/>	172DN1	Mounting location: Right Back Cane	N/C
<input type="checkbox"/>	172DN2	Mounting Location: Left Back Cane	N/C

R-net OMNI and Output Options			
<i>(For expandable controls, choose all that apply)</i>			
<input type="checkbox"/>	172ON6	Rnet OMNI Input Control Module and display with built in Infrared and armrest mount. Choose mount below. (required for Specialty Input Devices.)	\$2,075.00
<input type="checkbox"/>	172ON7	Rnet Output Module (ECM)	\$1,066.00
<input type="checkbox"/>	172ON8	Rnet Bluetooth Mouse-mover for Windows / Android	\$1,066.00
<input type="checkbox"/>	172ON9	Rnet Bluetooth for Apple Devices	\$1,066.00

OMNI Mount			
<i>(If used with Joystick must be mounted opposite side)</i>			
<input type="checkbox"/>	172DS04	Swing Away Left Side OMNI Mount	NC
<input type="checkbox"/>	172DS03	Swing Away Right Side OMNI Mount	NC

Switch inputs for Joysticks or OMNI			
<i>Drive or Mode / Profile Control</i>			
<input type="checkbox"/>	172SW22	Stereo to 2 Mono Splitter	\$75.00
<input type="checkbox"/>	172SW6	Ribbon Switch	\$119.00
<input type="checkbox"/>	172SW1	Buddy Button Switch	\$95.00
<input type="checkbox"/>	172SW2	Disc Switch	\$202.00
<input type="checkbox"/>	172SW9	Wobble Switch	\$630.00
<input type="checkbox"/>	172SW4	Micro Light Switch	\$95.00
<input type="checkbox"/>	172SW24	Egg Switch-Black	\$119.00
<input type="checkbox"/>	172SW21	Dual toggle switch and bracket	\$214.00
<input type="checkbox"/>	172SW35	Switch-it Mini 1 button on Mono plug	\$268.00
<input type="checkbox"/>	172SW36	Switch-it Mini 2 button on 2 Mono plugs	\$268.00
<input type="checkbox"/>	172SW37	Switch-it Mini 2 button on 1 Stereo plug	\$268.00
<input type="checkbox"/>	172SW38	Switch-it Mini 4 button on DB9 connector	\$268.00

Switch Mount Bracket			
<input type="checkbox"/>	172EM01	Single Switch Mount Bracket	\$247.00

Rnet Electronic Programmer			
<input type="checkbox"/>	172PR5	R-net PC Programmer with Dongle	\$615.00

5. Specialty Driver Control Options

Requires expandable Rnet controls and OMNI unless otherwise noted. Not available with Upholstered backs.

5a. Proportional Control Packages

Head Control Proportional Packages			
<input type="checkbox"/>	Switch-it Dual Pro Proportional Head Control with Whitmyer and Facial Pads (Adult)		\$6,541.00
	172ID72	Switch-it Dual Pro Head Control w/facial pad (Adult)	\$6,139.00
	172IJ09	SA Boom Mount with Egg Switch Control	\$402.00
<input type="checkbox"/>	Switch-it Dual Pro Proportional Head Control with Spot Pads Package		\$6,541.00
	172ID70	Switch-it Dual Pro Head Control w/Spot pads	\$6,139.00
	172IJ09	SA Boom Mount with Egg Switch Control	\$402.00

Hand Control Proportional Joystick Packages			
<input type="checkbox"/>	Switch-it Versa-Guide Compact Joystick Midline with Gatlin Package		\$2,369.00
	172ID71	Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
	172IJ07	S/A Gatlin Mount	\$330.00
<input type="checkbox"/>	Switch-it Versa-Guide Compact Joystick Midline and prox flex switches with Gatlin Package		\$3,018.00
	172ID71	Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
	172IM08	Dual Flex Proximity Switches (requires no force)	\$649.00
	172IJ07	S/A Gatlin Mount	\$330.00
<input type="checkbox"/>	Switch-it Micro-Pilot Joystick with Space Disc Package		\$4,007.00
	172ID62	Micro-Pilot Proportional Joystick	\$3,534.00
	172IJ08	S/A SpaceDisc mount	\$473.00
<input type="checkbox"/>	Switch-it Micro-Guide Joystick with SpaceDisc Package		\$4,831.00
	172ID63	Micro-Guide Proportional Joystick	\$4,358.00
	172IJ08	S/A SpaceDisc mount	\$473.00
<input type="checkbox"/>	ASL Proportional Micro Extremity Joystick (ASL130) Midline with Gatlin Package		\$4,414.00
	172ID59	Micro Extremity Joystick	\$4,084.00
	172IJ07	S/A Gatlin Mount	\$330.00
<input type="checkbox"/>	ASL Mushroom Joystick (ALS135) Midline with Gatlin Package		\$3,353.00
	172ID60	Mushroom Joystick	\$3,023.00
	172IJ07	S/A Gatlin Mount	\$330.00
<input type="checkbox"/>	ASL Compact Joystick w Integrated buttons (ASL134) Midline with Gatlin Package		\$1,530.00
	172ID61	ASL Compact Joystick	\$1,200.00
	172IJ07	S/A Gatlin Mount	\$330.00
<input type="checkbox"/>	ASL Micro Mini Joystick w Hand Pad (ASL136) Midline with Gatlin Package		\$4,012.00
	172ID68	Micro Mini Joystick	\$3,682.00
	172IJ07	S/A Gatlin Mount	\$330.00

Chin Control Proportional Joystick Packages			
<input type="checkbox"/>	Switch-it Versa-Guide Chin Control Compact Joystick and Egg Switch Package		\$2,406.00
	172ID71	Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
	172IM07	S/A Mount off backrest	\$252.00
	172IJ04	SA Boom Mount with egg switch	\$115.00
<input type="checkbox"/>	Switch-it Versa-Guide Chin Control Compact Joystick and Proximity Flex Switches Package		\$2,940.00
	172ID71	Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
	172IM08	Dual Flex Proximity Switches (requires no force)	\$649.00
	172IM07	S/A Mount off backrest	\$252.00
<input type="checkbox"/>	Switch-it Versa-Guide Bib Chin Control Compact Joystick with Proximity Flex switches Package Std 6"		\$3,445.00
	172ID71	Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
	172IM08	Dual Flex Proximity Switches (requires no force)	\$649.00
	172IM05	Bib mount Small	\$757.00
<input type="checkbox"/>	Switch-it Versa-Guide Bib Chin Control Compact Joystick with Proximity Flex Switches Package Large 8"		\$3,512.00
	172ID71	Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
	172IM08	Dual Flex Proximity Switches (requires no force)	\$649.00
	172IM06	Bib mount Large	\$824.00
<input type="checkbox"/>	Switch-it Versa-Guide Bib Chin Control Compact Joystick with Egg Switch Package Std 6"		\$2,911.00
	172ID71	Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
	172IM05	Bib mount Small	\$757.00
	172IJ04	SA Boom Mount with egg switch	\$115.00
<input type="checkbox"/>	Switch-it Versa-Guide Bib Chin Control Compact Joystick with Egg Switch Package Large 8"		\$2,978.00
	172ID71	Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
	172IM06	Bib mount Large	\$824.00
	172IJ04	SA Boom Mount with egg switch	\$115.00
<input type="checkbox"/>	Switch-it Micro-Pilot Bib Chin Control with Egg Switch Package Std 6"		\$4,406.00
	172ID62	Micro-Pilot Proportional Joystick	\$3,534.00
	172IM05	Bib mount Small	\$757.00
	172IJ04	SA Boom Mount with egg switch	\$115.00
<input type="checkbox"/>	Switch-it Micro-Pilot Bib Chin Control with Egg Switch Package Large 8"		\$4,473.00
	172ID62	Micro-Pilot Proportional Joystick	\$3,534.00
	172IM06	Bib mount Large	\$824.00
	172IJ04	SA Boom Mount with egg switch	\$115.00
<input type="checkbox"/>	Switch-it Micro-Guide Bib Chin Control with Egg Switch Package Std 6"		\$5,230.00
	172ID63	Micro-Guide Proportional Joystick	\$4,358.00
	172IM05	Bib mount Small	\$757.00
	172IJ04	SA Boom Mount with egg switch	\$115.00
<input type="checkbox"/>	Switch-it Micro-Guide Bib Chin Control with Egg Switch Package Large 8"		\$5,297.00
	172ID63	Micro-Guide Proportional Joystick	\$4,358.00
	172IM06	Bib mount Large	\$824.00
	172IJ04	SA Boom Mount with egg switch	\$115.00
<input type="checkbox"/>	ASL Proportional Micro Extremity Joystick Chin / Lip with Egg Switch Control Package		\$4,498.00
	172ID59	Micro Extremity Joystick	\$4,084.00
	172IM07	S/A Mount off backrest	\$299.00
	172IJ04	SA Boom Mount with egg switch	\$115.00

5. Specialty Driver Control Options (cont.)

5b. Switch Control

Requires Headrest and/or Headrest Mount

3 Switch Proximity Head Arrays *Choose headrest mount below*

<input type="checkbox"/>	ASL Proximity Head Array Package (Adult)		\$5,755.00
	172II02	ASL Proximity Head Array (Adult)	\$5,353.00
	172IJ09	SA Boom Mount with Egg Switch Control	\$402.00
<input type="checkbox"/>	ASL Proximity Head Array Package (Ped)		\$5,755.00
	172II03	ASL Proximity Head Array (Ped)	\$5,353.00
	172IJ09	SA Boom Mount with Egg Switch Control	\$402.00
<input type="checkbox"/>	Switch-It Head Array Package 3 Switch Curved (Long)		\$5,755.00
	172II13	Switch-it Proximity Head Array Long	\$5,353.00
	172IJ09	SA Boom Mount with Egg Switch Control	\$402.00
<input type="checkbox"/>	Switch-It Head Array Package 3 switch curved (Stubby)		\$5,691.00
	172II14	Switch-it Proximity Head Array Stubby	\$5,289.00
	172IJ09	SA Boom Mount with Egg Switch Control	\$402.00

5 Switch Proximity Head Array (No Egg Switch Needed) *Choose "Headrest Mount for Special Controls"*

<input type="checkbox"/>	172II04	Switch-It Head Array 5 switch flat	\$5,353.00
--------------------------	---------	------------------------------------	------------

Multiple Switches.

<input type="checkbox"/>	172SW31	Switch-it 3 Naked proximity switches on DB9 connector. (requires OMNI)	\$3,039.00
<input type="checkbox"/>	172SW32	Switch-it 4 Naked proximity switches on DB9 connector. (requires OMNI)	\$3,718.00
<input type="checkbox"/>	172SW33	Switch-it 5 Naked proximity switches on DB9 connector. (requires OMNI)	\$4,146.00
<input type="checkbox"/>	172SW34	ASL (502) 5-Switch adaptor box on DB9 connector. (for use with 5 mono switches, not included)	\$4,146.00
<input type="checkbox"/>	172SW23	4-way toggle switch button DB9 Connector. Includes bracket	\$670.00

Sip N Puff Package *Choose headrest and HR mount below*

<input type="checkbox"/>	Quickie Sip n' Puff System		\$2,635.00
	172ID34	Sip N' Puff Kit	\$2,063.00
	172IH1	Sip n' Puff Tubing and SA Hardware	\$457.00
	172IJ10	SA Boom Mount with egg switch	\$115.00

Headrest Mount for Specialty Controls

<input type="checkbox"/>	172MH12	Axys Mount & Bracket for Specialty Controls	\$309.00
<input type="checkbox"/>	172MH13	Cobra Mount & Bracket for Specialty Controls	\$386.00
<input type="checkbox"/>	172MH14	Linx Mount & Bracket for Specialty Controls	\$309.00

Headrest Pad for Specialty controls

<input type="checkbox"/>	172HE79	Plush Pad 8" for Specialty controls	\$309.00
<input type="checkbox"/>	172HE80	Plush Pad 10" for Specialty controls	\$309.00
<input type="checkbox"/>	172HE81	Plush Pad 14" for Specialty controls	\$309.00

Extra Swing Away Head Boom Egg Switch and Mount

	172IJ05	Extra SA Boom Mount with egg switch #1	\$115.00
	172IJ06	Extra SA Boom Mount with egg switch #2	\$115.00

Single Switches

<input type="checkbox"/>	172IB01	Ribbon Switch	\$119.00
<input type="checkbox"/>	172IB02	Buddy Button Switch	\$108.00
<input type="checkbox"/>	172IB03	Disc Switch	\$202.00
<input type="checkbox"/>	172IB04	Wobble Switch	\$630.00
<input type="checkbox"/>	172IB05	Micro Light Switch	\$95.00
<input type="checkbox"/>	172IB06	Egg Switch-Black	\$134.00
<input type="checkbox"/>	172IB07	Switch-it Mini 1 button on Mono plug	\$268.00

6. SEAT TYPE

<input type="checkbox"/> 172ST56	Quickie A.S.A.P II (Adjustable Seating and Positioning) Choose Seat Frame Size Platform.	\$830.00
----------------------------------	--	----------

ASAP II Frame Size Platform

(each frame is width and depth adjustable within given range.)

<input type="checkbox"/> 172FU17	Adjustable 12x12 to 16X19 (A)	N/C
<input type="checkbox"/> 172FU18	Adjustable 16x12 to 20x20 (B)	N/C
<input type="checkbox"/> 172FU19	Adjustable 18x16 to 22x22 (C)	N/C

7. SEAT SIZE

Seat Size Matrix. *Seat type is chosen in section 5.*

A = ASAP II Seat Frame
(Adjustable from 12x12 to 16X19)

B = ASAP II Seat Frame
(Adjustable from 16X12 to 20X20)

C = ASAP II Seat Frame
(Adjustable from 18X16 to 22X22)

		ASAP II Seat Width									
		12"	13"	14"	15"	16"	17"	18"	19"	20"	22"
Seat Depth	12"	A	A	A	A	A,B	B	B	B	B	
	13"	A	A	A	A	A,B	B	B	B	B	
	14"	A	A	A	A	A,B	B	B	B	B	
	15"	A	A	A	A	A,B	B	B	B	B	
	16"	A	A	A	A	A,B	B	B,C	B,C	B,C	C
	17"	A	A	A	A	A,B	B	B,C	B,C	B,C	C
	18"	A	A	A	A	A,B	B	B,C	B,C	B,C	C
	19"	A	A	A	A	A,B	B	B,C	B,C	B,C	C
	20"					B	B	B,C	B,C	B,C	C
	21"							C	C	C	C
	22"							C	C	C	C

Seat Width			
<input type="checkbox"/> 172R21	12"	Seat width	(N/C)
<input type="checkbox"/> 172R22	13"	Seat width	(N/C)
<input type="checkbox"/> 172R23	14"	Seat width	(N/C)
<input type="checkbox"/> 172R24	15"	Seat width	(N/C)
<input type="checkbox"/> 172R5	16"	Seat width	(N/C)
<input type="checkbox"/> 172R6	17"	Seat width	(N/C)
<input type="checkbox"/> 172R7	18"	Seat width	(N/C)
<input type="checkbox"/> 172R8	19"	Seat width	(N/C)
<input type="checkbox"/> 172R9	20"	Seat width	(N/C)
<input type="checkbox"/> 172R25	22"	Seat width	(N/C)

Seat Depth			
<input type="checkbox"/> 172RF20	12"	Seat depth	(N/C)
<input type="checkbox"/> 172RF21	13"	Seat depth	(N/C)
<input type="checkbox"/> 172RF22	14"	Seat depth	(N/C)
<input type="checkbox"/> 172RF23	15"	Seat depth	(N/C)
<input type="checkbox"/> 172RF5	16"	Seat depth	(N/C)
<input type="checkbox"/> 172RF6	17"	Seat depth	(N/C)
<input type="checkbox"/> 172RF7	18"	Seat depth	(N/C)
<input type="checkbox"/> 172RF8	19"	Seat depth	(N/C)
<input type="checkbox"/> 172RF24	20"	Seat depth	(N/C)
<input type="checkbox"/> 172RF26	21"	Seat depth	(N/C)
<input type="checkbox"/> 172RF25	22"	Seat depth	(N/C)

24" Seat Width			
<input type="checkbox"/> 172VU16	24"	Seat width (fixed width, depth is adjustable 16" to 22" deep) <i>N/A with Power or manual recline back or sling back upholstery.</i>	\$525.00

8. SEAT CUSHIONS

JAY Cushion Bundling			
<input type="checkbox"/> 172JC14		Jay Cushion <u>On Chair Order</u> <u>See order Form for codes</u>	See Order Form

9. POWERED SEATING

Power Seating Options

<input type="checkbox"/>	172PJ7	Quickie SPOT 50° CG Tilt (1 Actuator)	\$6,961.00
<input type="checkbox"/>	172PJ5	Quickie HD SPOT 50° CG Tilt (400lbs weight limit) (1 Actuator)	\$6,961.00

Powered Seat Elevator- (AKA: Power Lift)

<input type="checkbox"/>	172PA27	12" Power Seat Elevate for use with tilt (Must be ordered with SPOT tilt, NA with QM-715 HD)(1 Actuator)	\$3,452.00
--------------------------	---------	--	------------

Power Seating Control

Powered Seating Switch Control Method

(For use with power seating actuators. Switch activation is stand-alone and does not require expandable or thru-drive electronics.)

Single Actuator Switch Control (Tilt)

<input type="checkbox"/>	172EO10	Switch Control 1 actuator with <u>Dual Toggle</u>	N/C
<input type="checkbox"/>	172EO11	Switch Control 1 actuator with <u>Buddy Button</u>	N/C
<input type="checkbox"/>	172EO19	Switch Control 1 actuator with <u>Switch-it Mini 1 button</u>	N/C
<input type="checkbox"/>	172EO20	Switch Control 1 actuator with <u>Switch-it Mini 2 button</u>	N/C

2-3 Actuator Switch Control (Tilt/Legs, Lift/Tilt)

<input type="checkbox"/>	172EO21	Switch Control 2-3 actuator with <u>Dual Toggle</u>	\$750.00
<input type="checkbox"/>	172EO12	Switch Control 2-3 actuator with <u>4-way toggle</u>	\$750.00
<input type="checkbox"/>	172EO23	Switch Control 2-3 actuator with <u>Switch-it Mini 2 button</u>	\$750.00
<input type="checkbox"/>	172EO22	Switch Control 2-3 actuator with <u>Switch-it Mini 4 button</u>	\$750.00

Switch Mount

(for use with switch control powered seating)

<input type="checkbox"/>	172EY01	Switch Mount Right	N/C
<input type="checkbox"/>	172EY02	Switch Mount Left	N/C

Powered Seating Thru-drive Control Method.

Thru-Drive 1 Seat Function

<input type="checkbox"/>	172EO28	Thru Drive 1 seat Function (For operating Tilt thru Joystick.) Includes Ctrl+5	\$1,751.00
--------------------------	---------	--	------------

Thru-Drive 2+ Seat Function

<input type="checkbox"/>	172EO29	Thru Drive 2+ Seat Functions: (For operating Tilt, lift, power legs) Includes Ctrl+5	\$2,525.00
--------------------------	---------	--	------------

Powered Seating Thru-drive Control Options

Assignable Buttons

"Assignable Button" feature allows for customization of joystick buttons and/or Ctrl+5 Switches allowing access to various wheelchair and seating functions. When enabled, chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle.

<input type="checkbox"/>	172PM01	Enable Assignable buttons	N/C
--------------------------	---------	---------------------------	-----

Ctrl+5 Switchbox Options

Ctrl+5 allows direct control to 5 power seating functions (10 total switches). Switches can also be used with assignable buttons for greater flexibility. Box includes 3-Amp USB port for charging portable devices

<input type="checkbox"/>	172QG01	Ctrl+5 Switchbox mounted behind joystick	N/C
<input type="checkbox"/>	172QG02	Ctrl+5 Switchbox mounted opposite joystick	N/C
<input type="checkbox"/>	172QG03	Ctrl+5 Switchbox mounted on back cane	N/C
<input type="checkbox"/>	172QG05	Omit Ctrl+5	-\$618.00

Ctrl+5 Switch Type

<input type="checkbox"/>	172EY03	Push buttons	N/C
<input type="checkbox"/>	172EY04	Toggle switches	\$309.00

10. BACKREST AND BACK OPTIONS

Backrest Frame Type

<input type="checkbox"/>	172BT8	Angle Adjustable (ASAP II Seat only)	N/C
<input type="checkbox"/>	172BT98	Manual Recline 40° (ASAP II Seat only) <i>not width adjustable</i>	\$536.00

Manual Recline Adjustment Method Option

Back Cane Height 23".

<input type="checkbox"/>	172MN01	Tool	N/C
<input type="checkbox"/>	172MN02	Knob (NA with Cantilever or dual post arms.)	N/C

Angle Adjustable Backrest Options

(NA with Manual recline)

<input type="checkbox"/>	172OD8	Non-Folding Backrest	N/C
<input type="checkbox"/>	172OD9	Folding Backrest	\$137.00

Cane Angle

<input type="checkbox"/>	172BP40	Straight Back Post	N/C
<input type="checkbox"/>	172BP41	8° Bend Back Post (Avail. only w/ Upholstery or J3 Back)	\$126.00

Towel Bar Depth

<input type="checkbox"/>	172QT03	Shallow (NA with Upholstery backs.)	N/C
<input type="checkbox"/>	172QT04	Deep (Extends 4" from back canes for deep backs)	N/C

Back Cane

(Upholstery Height Measured pan to top of upholstery.)

<input type="checkbox"/>	Low (15.5") 172B1	<input type="checkbox"/>	Med (17.5") 172B5	<input type="checkbox"/>	Tall (19.5") 172B9
--------------------------	-------------------	--------------------------	-------------------	--------------------------	--------------------

Back Type

<input type="checkbox"/>	172BU76	S.P.O.T. Curved	N/C
<input type="checkbox"/>	172BU77	Jay SPO Posterior Positioning (ASAP II only)	\$598.00
<input type="checkbox"/>	172BU78	ART SPO Custom Back (ASAP II only)	Order form
<input type="checkbox"/>	172BU14	Omit Back.	NC

Jay 3 Option (select Omit above)

<input type="checkbox"/>	172JC08	Jay J3 <u>On Chair Order</u> <i>Choose back, headrest, laterals</i>	See List
--------------------------	---------	--	----------

Back Height - Required for SPOT Curved, SPO Backs.




(Limited combinations of widths) Measured from pan to top of solid back.

<input type="checkbox"/>	172B132	14"	N/C
<input type="checkbox"/>	172B36	16"	N/C
<input type="checkbox"/>	172B38	18"	N/C
<input type="checkbox"/>	172B40	20"	N/C
<input type="checkbox"/>	172B42	22"	N/C
<input type="checkbox"/>	172B44	24"	N/C
<input type="checkbox"/>	172B140	26"	N/C


12. HEADREST OPTIONS

Choose One Each.

Headrest Mount			
<input type="checkbox"/>	172MH10	Axys Mount & Bracket	\$315.00
<input type="checkbox"/>	172MH11	Cobra Mount and Bracket	\$578.00
<input type="checkbox"/>	172MH1	Linx Mount and Bracket	\$237.00



Headrest Pad			
<input type="checkbox"/>	172HE65	Plush Head Support Pad 8"	\$237.00
<input type="checkbox"/>	172HE66	Plush Head Support Pad 10"	\$237.00
<input type="checkbox"/>	172HE78	Plush Head Support Pad 14"	\$237.00
<input type="checkbox"/>	172HE82	Adjusta Plush Head Support Pad Narrow	\$365.00
<input type="checkbox"/>	172HE83	Adjusta Plush Head Support Pad std	\$365.00



13. POSTURAL SUPPORTS

LATERAL THORACIC SUPPORT (Swing Away)

Available on all Solid Backs.

<input type="checkbox"/>	LAT 4X4 S/A MOUNTS W/ 1" OFFSET - PR		\$409.00
	172ZV14	LATERAL 4X4 W/ 1" OFFSET - PAIR	\$126.00
	172ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$283.00
<input type="checkbox"/>	LAT 4X5 S/A MOUNTS W/ 1" OFFSET - PR		\$409.00
	172ZV15	LATERAL 4X5 W/ 1" OFFSET - PAIR	\$126.00
	172ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$283.00
<input type="checkbox"/>	LAT 4X6 S/A MOUNTS W/ 1" OFFSET - PR		\$409.00
	172ZV16	LATERAL 4X6 W/ 1" OFFSET - PAIR	\$126.00
	172ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$283.00
<input type="checkbox"/>	LAT 6X4 S/A MOUNTS W/ 1" OFFSET - PR		\$409.00
	172ZV17	LATERAL 6X4 W/ 1" OFFSET - PAIR	\$126.00
	172ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$283.00
<input type="checkbox"/>	LAT 6X5 S/A MOUNTS W/ 1" OFFSET - PR		\$409.00
	172ZV18	LATERAL 6X5 W/ 1" OFFSET - PAIR	\$126.00
	172ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$283.00
<input type="checkbox"/>	LAT 6X6 S/A MOUNTS W/ 1" OFFSET - PR		\$409.00
	172ZV19	LATERAL 6X6 W/ 1" OFFSET - PAIR	\$126.00
	172ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$283.00

CHEST STRAP

<input type="checkbox"/>	172ZR03	Padded Chest Strap	\$283.00
--------------------------	---------	--------------------	----------

Swing Away Knee Adductor Supports

Available only with Fixed and Power Center Mount Legrests.

Not available with Dual Post Arm Rests or 2-point Locking Cantilever Arms

Knee Adductor Supports - LEFT (choose Pad and Arm)			
Pad Size			
<input type="checkbox"/>	172KB5	Knee Adductor Pad SM- Left	\$215.00
<input type="checkbox"/>	172KB7	Knee Adductor Pad LG- Left	\$215.00
S/A Arm Size			
<input type="checkbox"/>	172XX04	Swing Away Short - Left	\$283.00
<input type="checkbox"/>	172XX06	Swing Away Long - Left	\$283.00
Knee Adductor Supports - RIGHT (choose Pad and Arm)			
Pad Size			
<input type="checkbox"/>	172KB6	Knee Adductor Pad SM- Right	\$215.00
<input type="checkbox"/>	172KB8	Knee Adductor Pad LG- Right	\$215.00
S/A Arm Size			
<input type="checkbox"/>	172XX05	Swing Away Short - Right	\$283.00
<input type="checkbox"/>	172XX07	Swing Away Long - Right	\$283.00

Long Thigh Support w/Removable Mount

Available with Cantilever Arms, NA with Manual Recline Back, Dual Post, Single Post or Cantilever w/ 2 pt Locking arms.

<input type="checkbox"/>	SPOT Long Thigh Support - Left		\$409.00
	172TV2	Long Thigh Support - Left	\$126.00
	172TX1	Removable Mount - Left	\$283.00
<input type="checkbox"/>	SPOT Long Thigh Support - Right		\$409.00
	172TV1	Long Thigh Support - Right	\$126.00
	172TX2	Removable Mount - Right	\$283.00

14. ARMREST AND ARMPAD OPTIONS

Cantilever Armrest Mounting

Cantilever Armrests
(Adjustable length arm pad included (choose Pad setting below, except for Low Arm Option) NA with Manual Recline Back

<input type="checkbox"/>	172EG12	STD Cantilever Arm 9"-12" Ht Range	Std
<input type="checkbox"/>	172EG13	Low Cantilever Arm 7"-10" Ht Range (Choose ASAPII pad type and size below)	NC
<input type="checkbox"/>	172EG10	Standard Cantilever Arm w/ 2 Point Locking Post. (9"-10" Ht Range.) NA w/ seat depths below 16"	\$289.00
<input type="checkbox"/>	172EG11	Standard Cantilever Arm w/ 2 Point Locking Post. (11"-12" Ht Range.) NA w/ seat depths below 16"	\$289.00

Cantilever Armpad

Pad Setting Options for **Standard Cantilever Arm**

Left Arm Pad (Std)		Right Arm Pad (Std)		
<input type="checkbox"/>	172AE11	Full Length Cantilever (14")	<input type="checkbox"/>	172AB11
<input type="checkbox"/>	172AE10	Desk Length Cantilever (10")	<input type="checkbox"/>	172AB10

Single - Dual Post Armrest Mounting

Single Post Armrests			
<input type="checkbox"/>	172EG14	Single Post Height Adjustable - Std.	\$289.00
<input type="checkbox"/>	172EG15	Single Post Height Adjustable - Low	\$289.00

Dual Post Armrests (NA with Infinite Seat Depth)

<input type="checkbox"/>	172EG09	Shrt Dual Post Flip Back HT Adj 10"-15"	\$378.00
<input type="checkbox"/>	172EG23	Lng Dual Post Flip Back HT Adj 10"-15"	\$378.00

Pad Type / Size Options for Dual Post / Single Post / Low Cantilever

Left Arm Pad (N/C)		Right Arm Pad (N/C)		
<input type="checkbox"/>	172AE01	Desk Length Waterfall (10")	<input type="checkbox"/>	172AB01
<input type="checkbox"/>	172AE02	Full Length Waterfall (14") NA w/Low Cantilever arm	<input type="checkbox"/>	172AB02
<input type="checkbox"/>	172AE03	Desk Length Standard (10")	<input type="checkbox"/>	172AB03
<input type="checkbox"/>	172AE04	Full Length Standard (14") NA w/Low Cantilever arm	<input type="checkbox"/>	172AB04

15. HANGERS AND FOOTRESTS

Swing-away Hangers

Non-Powered Swing Away Legrests

Sold as Pairs

<input type="checkbox"/>	172H2	70° Swing Away (B)	N/C
<input type="checkbox"/>	172H6	90° Swing Away (B) (5.5-9")	N/C
<input type="checkbox"/>	172H138	Manual ELR S/A (B) (Avail. with Fixed Centermount, choose below)	\$352.00

Power Elevating Legrests (2 Actuators)

<input type="checkbox"/>	172H126	Swing Away Power ELR w/ Articulation Pair (B footplates) (N/A for 12" or 13" wide)	\$2,994.00
--------------------------	---------	--	------------

Footplates For Swing Away Hangers

<input type="checkbox"/>		B - Adult Angle Adjust w/Heel Loop for use with 70° hanger and ELR (Pair)	\$205.00
	172F282	Adult Angle Adjustable Footplate	\$205.00
	172HL10	Heel Loop	N/C
<input type="checkbox"/>		B - Adult Angle Adjust w/Heel Loop for use with 90° hanger (Pair)	\$205.00
	172F285	Adult Angle Adjustable Footplate	\$205.00
	172HL10	Heel Loop	N/C
<input type="checkbox"/>		C - Kids Angle Adjust w/Leg Strap for use with 70° hanger (Pair)	\$205.00
	172f283	Kids Angle Adjustable Footplate	\$205.00
	172LS13	Leg Strap	N/C
<input type="checkbox"/>		C - Kids Angle Adjust w/Leg Strap for use with 80° and 90° hanger(Pair)	\$205.00
	172f286	Kids Angle Adjustable Footplate	\$205.00
	172LS13	Leg Strap	N/C
<input type="checkbox"/>	172F279	Omit Footplate and extension tube	N/C

Omit Hangers.

<input type="checkbox"/>	172H144	Omit Hangers	N/C
--------------------------	---------	--------------	-----

Hanger Extensions. (See range listed in matrix below)

NOTE: Not all extension tube will be available; depending on the seat height, some extension tube options may not offer enough ground clearance.

<input type="checkbox"/>	172E16	High Mount Hanger	\$84.00
<input type="checkbox"/>	172E1	Short	N/C
<input type="checkbox"/>	172E2	Medium	N/C
<input type="checkbox"/>	172E5	Long	N/C

Footrest Accessories

<input type="checkbox"/>	1725A01	Toe Loop (Pair) <i>Only avail on (opt B).</i>	\$47.00
<input type="checkbox"/>	172LS14	Extra Leg Strap	\$32.00

Centermount Legrests

Non-Powered Center Mount Legrests

Fixed Center Mount Legrests (includes footplate) Plus options below

<input type="checkbox"/>	172H134	Fixed Center ASAP II Mount Short	N/C
<input type="checkbox"/>	172H135	Fixed Center ASAP II Mount Long	N/C

Angle Adjustable Center Mount Legrests for ASAP II

<input type="checkbox"/>	172H104	Angle adjustable Center mount Med (I) <i>Independent leg length and angle adjustable</i>	<i>Built-4-Me Only</i>
<input type="checkbox"/>	172H108	Angle adjustable Center mount Large (I) <i>Independent leg length and angle adjustable</i>	<i>Built-4-Me Only</i>

Power Centermount (2 Actuators)

<input type="checkbox"/>	172H127	Power Centermount ELR w/7.5" Power Articulation (F) Centermount footplates plus Custom Options Below	\$2,987.00
--------------------------	---------	--	------------

Power Centermount Footplate Extension Option (For Thru-Drive Operation)

<input type="checkbox"/>	172IE01	Enable 7.5" Powered footplate extension for transfer assist. (Electronics programmed to allow footplates to be extended with or without legrest being raised)	\$515.00
<input type="checkbox"/>	172IE02	Disable Footplate extension for transfers. (Electronics programmed to ONLY extend footplates when legs are raised.)	N/C

Centermount Footplates

<input type="checkbox"/>		F - Centermount Dual Footplate (Pair) (Individually Height (3"-19") and Angle Adjustable)	\$247.00
	172F315	CMT Angle Adjustable Footplates	\$247.00
<input type="checkbox"/>		I - Angle adj Centermount Footplate pair	<i>Built-4-Me Only</i>
	172F125	Angle Adjustable Footplates	<i>Built-4-Me Only</i>

Custom Centermount Options

<input type="checkbox"/>	172DI26	Extra Large Single Footplate for Centermount Footrest. Size 16"W x 14"D	\$309.00
<input type="checkbox"/>	172DI27	Gell Foot Pads for XL Single Centermount Footplate.	\$309.00
<input type="checkbox"/>	172DI28	Extra Large Pair Curved Calf Pads for Center Mount Legrests. Fully adjustable	\$309.00

Hanger Type	Lower Leg Length Matrix (measured seat pan to the top of the footplate surface)					Footplate Matrix			
	Standard	High Mount	Short	Medium	Long	B	F	I	J
70° Swing Away		7"-14"	13.5"- 16"	16"-18.5"	18"-20.5"	X			
90° Swing Away		5.5-9" Footplate mounts to front of hanger.				X			
Manual ELR S/A			12.5"-17"	14.5"-19"	17"-21"	X			
Power ELR S/A			12"-14"	14"-17"	16.5"-19"	X			
Fixed Cntrmnt ASAP			12.5"-15.5"		14.5"-19.5"				
Angle Adj Cntrmnt				3"-13"	3"-15"			X	X
Power Cntrmnt	3"-19"						X		

16. ACCESSORIES

Wheelchair Tray Table			
<input type="checkbox"/>	1727A01	TRAY TABLE SM TUFFAK 12"-14" PWR	\$331.00
<input type="checkbox"/>	1727A02	TRAY TABLE MED TUFFAK 15"-17" PWR	\$331.00
<input type="checkbox"/>	1727A03	TRAY TABLE LRG TUFFAK 18"-20" PWR	\$331.00

Positioning Belt - Select for any style seat			
<input type="checkbox"/>	172EJ01	Auto Style Buckle	N/C
<input type="checkbox"/>	172EJ02	Belt 2" Aircraft Buckle	\$163.00
<input type="checkbox"/>	172EJ03	Belt 2" Aircraft Buckle Padded	\$178.00

Backpack			
<input type="checkbox"/>	172A1	Black Utility	\$97.00

Accessories			
<input type="checkbox"/>	1728A09	Ram Phone Mount on Round JS tube	\$67.00
<input type="checkbox"/>	1728A10	Ram Tablet Mount on Round JS tube	\$268.00
<input type="checkbox"/>	1728A11	Ram Cup holder	\$67.00
<input type="checkbox"/>	1728A12	Switch-it Remote Stop Switch	\$1,200.00

17. ACCESSORIES



*Save time! Wheelchair and seating will be delivered together.
Hardware pre-installed.*

Cushions			
<input type="checkbox"/>	Jay Cushion:		
	Width:	Depth:	
	Notes:		

Backs			
<input type="checkbox"/>	Jay Back:		
	Width:	Height:	
	Notes:		

Headrest	
<input type="checkbox"/>	System:
	Pads:
	Notes:
	Hardware:
	Notes:

Jay Positioning Components	
<input type="checkbox"/>	Anterior Trunk Supports:
	Size:
	Hand and Arm Supports:
	Size:
	Pelvic Positioning:
	Size:
	Foot Positioning:
	Size:
	Notes: