



Zippie Kidz
July 2018
EIKZ



Date: _____ **SUBMITTING FOR:** Quote Order

Dealer Acct #: _____

Dealer: _____ PO#: _____

Dealer Contact: _____ Ship To: _____

Dealer Address: _____ Attention: _____

Dealer City: _____ PV: _____ PC: _____ Address: _____

Dealer Phone: () _____ Fax: () _____ Address: _____

Ship To City: _____ PV: _____ PC: _____

Ship To Phone: () _____ Fax: () _____

Confirm Via: Fax Email

Mark For: _____

✓	Option #	Definition	Price	✓	Option #	Definition	Price
---	----------	------------	-------	---	----------	------------	-------

ZIPPIE KIDZ

<input type="checkbox"/>	013WM36	Zippie KIDZ	\$3,070.00
--------------------------	---------	-------------	------------

1. FRAME STYLE

<input type="checkbox"/>	013FC3	Standard Configuration	Std	<input type="checkbox"/>	013CV1	Reverse Configuration	\$97.00
<input type="checkbox"/>	013FC4	Reverse Configuration (Lg. wheels in front & casters in back. Not avail. w/ 20" rear wheels or 5" casters.	NC			Conversion Kit (Not avail. w/ 20" rear wheels or 5" casters. Used to reverse configure chair w/Lg. wheels in front & casters in back.	

2. FRAME COLOR

Standard Colors

<input type="checkbox"/>	013C140	Aztec Gold	NC
<input type="checkbox"/>	013C2	Black (Not avail as base w/Zebra)	NC
<input type="checkbox"/>	013C141	Black Cherry	NC
<input type="checkbox"/>	013C10	Black Opal	NC
<input type="checkbox"/>	013C38	Blue Opal	NC
<input type="checkbox"/>	013C36	Candy Blue	NC
<input type="checkbox"/>	013C14	Candy Purple	NC
<input type="checkbox"/>	013C35	Candy Red	NC
<input type="checkbox"/>	013C142	Electric Blue	NC
<input type="checkbox"/>	013C143	Evergreen	NC
<input type="checkbox"/>	013C126	Glow	NC
<input type="checkbox"/>	013C144	Green Apple	NC
<input type="checkbox"/>	013C147	Hot Sparkle Pink	NC
<input type="checkbox"/>	013C148	Mauve Pink	NC
<input type="checkbox"/>	013C103	Rootbeer	NC
<input type="checkbox"/>	013C114	Sparkle Silver	NC
<input type="checkbox"/>	013C145	Sunrise Orange	NC
<input type="checkbox"/>	013C113	Titanium Color Paint	NC
<input type="checkbox"/>	013C8	Yellow	NC

Standard Matte Colors

<input type="checkbox"/>	013C116	Matte Black	NC
<input type="checkbox"/>	013C151	Matte Black Cherry	NC
<input type="checkbox"/>	013C152	Matte Purple	NC
<input type="checkbox"/>	013C150	Matte Electric Blue	NC
<input type="checkbox"/>	013C149	Matte Evergreen	NC

Custom Colors

<input type="checkbox"/>	013C43	Mossy Oak Camouflage	\$340.00
<input type="checkbox"/>	013C42	Stars and Stripes	\$340.00
<input type="checkbox"/>	013C108	Zebra (Choose base color above: White Std, n/a black)	\$340.00

3. FRAME WIDTH

(Seat Width 1/2" Narrower)

<input type="checkbox"/>	013W1	10"	STD	<input type="checkbox"/>	013W4	13"	NC
<input type="checkbox"/>	013W2	11"	NC	<input type="checkbox"/>	013W5	14"	NC
<input type="checkbox"/>	013W3	12"	NC				

Mark For: _____

✓	Option #	Definition	Price	✓	Option #	Definition	Price
---	----------	------------	-------	---	----------	------------	-------

4. SEAT UPHOLSTERY DEPTH

12" rear wheels on regular frame are not available with high-pull wheel locks.

Regular Frame			
<input type="checkbox"/>	013SD1	8" Seat Depth	STD
<input type="checkbox"/>	013SD2	9" Seat Depth	NC
<input type="checkbox"/>	013SD3	10" Seat Depth	NC

Long Frame			
<input type="checkbox"/>	013SD4	11" Seat Depth	STD
<input type="checkbox"/>	013SD5	12" Seat Depth	NC
<input type="checkbox"/>	013SD6	13" Seat Depth	NC
<input type="checkbox"/>	013SD7	14" Seat Depth	NC

5. SEATING

<input type="checkbox"/>	013SL1	Seat Upholstery	STD
<input type="checkbox"/>	013SL6	Omit Seat Upholstery	NC
<input type="checkbox"/>	013CU1	2" Seat Cushion	\$173.00

Forward-Folding			
<input type="checkbox"/>	013B101	Regular Height 8"-10"	\$111.00
<input type="checkbox"/>	013B102	Tall Height 10"-12"	\$111.00

Non-Folding			
<input type="checkbox"/>	013B103	Regular Height 8"-10"	STD
<input type="checkbox"/>	013B104	Tall Height 10"-12"	NC

JAY Cushions (Optional)			
<i>To order any JAY cushions on-chair, please see the JAY price list. Cushion name or part number must be filled out with desired width and depth</i>			
<input type="checkbox"/>	013JC20	JAY Bundled Cushion	Price List
		Cushion Name or Part #	
		Width: Depth:	

6. BACKREST

<input type="checkbox"/>	013SR1	Stroller Handle	\$150.00
<input type="checkbox"/>	013BU14	Omit Back Upholstery	(\$21.00)
<input type="checkbox"/>	013BP8	Omit Backrest & Back Upholstery	(\$47.00)

Jay Backs (Optional)			
<i>Ordering JAY backs on-chair, See the JAY Price List. Requires Part #</i>			
<input type="checkbox"/>	013JB02	JAY Bundled Backrest	Price List
		Backrest Name or Part #	
		Width: Height:	

7. FOOTREST

<input type="checkbox"/>	013F54	Regular 8 1/2" to 10 3/4"	STD
<input type="checkbox"/>	013F46	Long 11" to 13 1/4"	NC

<input type="checkbox"/>	013LS5	Leg Strap	NC
<input type="checkbox"/>	013F38	High Mount	\$90.00

8. ARMRESTS (T-POST)

Adjust in 1/4" increments

<input type="checkbox"/>	013AA12	Kidz Std 4 1/2" to 5 3/4"	STD
<input type="checkbox"/>	013AM14	Foam Pad	NC
<input type="checkbox"/>	013AM1	Desk Pad	NC

<input type="checkbox"/>	013AA13	Kidz Std 5 1/2" to 6 3/4"	NC
<input type="checkbox"/>	013AM14	Foam Pad	STD
<input type="checkbox"/>	013AM1	Desk Pad	NC
<input type="checkbox"/>	013AA14	Omit Armrest	NC

9. CASTERS

Caster Wheel			
<input type="checkbox"/>	013CW1	3" Micro	NC
<input type="checkbox"/>	013CW30	3" Micro Light up	\$67.00
<input type="checkbox"/>	013CW5	5" Low-Profile Polyurethane	STD

Rear Wheel Size			
<input type="checkbox"/>	013RW13	12" Mag (Avail. w/ 3" casters only.) 12" rear wheels on reg. frame are not avail. With high pull wheel locks.	NC
<input type="checkbox"/>	013RW14	16" Mag (Not available with full poly tires)	STD
<input type="checkbox"/>	013RW3	20" Mag Wheel	NC
<input type="checkbox"/>	013RW4	20" Spoke	\$295.00

Caster Options			
<input type="checkbox"/>	013CF4	3/4" Ext. Fork Stem Bolt	NC
<input type="checkbox"/>	013CW12	Omit Caster Wheels	NC

Rear Wheel Tire			
<input type="checkbox"/>	013RT1	Pneumatic	NC
<input type="checkbox"/>	013RT2	Pneumatic Airless Insert	\$228.00
<input type="checkbox"/>	013RT3	Full Polyurethane	NC

Mark For: _____

--	--

✓	Option #	Definition	Price	✓	Option #	Definition	Price
---	----------	------------	-------	---	----------	------------	-------

10. AXLES

Axles			
<input type="checkbox"/>	013AX1	Quick Release Axles	STD
<input type="checkbox"/>	013AX11	Quad Quick Release Axles (Only avail. w/20" rear wheels. Not avail. w/16" rear wheel rev. config.)	\$168.00

UPDATED

Extension Handle			
<input type="checkbox"/>	013EH1	6" Ext. Handles (Avail. w/ high-push wheel locks only.)	\$67.00

Handrim			
<i>*Not available with 12" or 16" rear wheels.</i>			
<input type="checkbox"/>	013HR1	Aluminum Anodized	NC
<input type="checkbox"/>	013HR4	Plastic Coated	\$134.00
<input type="checkbox"/>	013HR7	Omit Handrims	NC

Anti-Tip (Standard Configuration)			
<input type="checkbox"/>	013AT1	Rear Anti-Tip Tubes	\$111.00
<input type="checkbox"/>	013AO1	Front Anti-Tip Tubes	\$145.00

UPDATED

Wheel Lock			
<input type="checkbox"/>	013WL9	High Push	STD
<input type="checkbox"/>	013WL8	High Pull (12" rear wheels on reg. frame are not avail. w/ high pull wheel locks. Not avail. w/16" rear wheel rev. config.)	NC
<input type="checkbox"/>	013WL10	Omit	NC

Anti-Tip (Reverse Configuration)			
<input type="checkbox"/>	013AT19	Rear Anti-Tip Tubes	\$111.00
<input type="checkbox"/>	013AO1	Front Anti-Tip Tubes	\$145.00

11. ACCESSORIES

Upper Extremity Support Tray			
<input type="checkbox"/>	013A29	X-Sm (Suggested 10"-12" Width)	\$351.00
<input type="checkbox"/>	013A30	Sm (Suggested 13"-14" Width)	\$351.00

Miscellaneous			
<input type="checkbox"/>	013A83	Tool Kit	\$39.00

Spoke Guard			
<i>(Not available with 12" & 16" rear wheels)</i>			
<input type="checkbox"/>	013A23	Black	\$145.00
<input type="checkbox"/>	013A303	Clear Zippie	\$145.00