

X′CAPE™

The Zippie X'CAPE weighs as little as 19 pounds and is the first ever wheelchair with XLOCK[™] Technology for the performance of a rigid frame and convenience of a folding frame.



Lily explores her world in the Zippie X'CAPE

FEATURES AND BENEFITS

Patented XLOCK[™] Technology

XLOCK lockout crossbrace and fewer moving parts reduce weight and frame play for greater efficiency and less fatigue.

Built-in Growth

2" of width growth in the crossbrace and 3" of depth growth in the seat rail keep your child positioned properly as they grow.

WC-19 Transit Approved

WC-19 transit approved for use as an occupant's seat in a motor vehicle when used with an approved Wheelchair Tiedown and Occupant Restraint System in accordance with SAE (J2249).

$Z\text{-finity}^{^{\rm TM}} System$

Lower leg length setting can be adjusted between -.5" and 21" from the seat and footrest angle can be adjusted infinitely to achieve the best positioning as your child grows.

Swing-in/Swing-out & Fixed Hanger Options

Swing-in/Swing-out hanger option aids in transfers and pulling up close to cabinets and counters. Fixed front frame option offers the lightest weight configuration and an active style.

VersaRail[™] Seat Rail

VersaRail provides your child access to all of the components and accessories needed – custom seating, armrests, wheel locks, side guards, etc. – even on the shortest frames.











One-step Push-to-Lock One-step push-to-lock tab locks out the frame.

One-step Pull-to-Fold One-step pull-to-fold strap collapses the frame effortlessly

Folding Push Handle Tucks away or raises to your height.

Color Options

27 trame colors and 24 crossbrace colors can create over 600 combinations!

Refer to the Quickie®/Zippie® Color Brochure or the Sunrise Medical website to view color choices.

HCPCS Code: E1238

Chair Weight: 19 lbs 8x8" fixed front frame, 20" lite spoke wheels. Excludes armrests, push handles & anti-tips. Actual weight may vary according to chair configuration & accessories.

Overall Width 14.8" to 26.5"

Folded Width Starting at 5" (without wheel hangers, and footplates)

Weight Capacity 165 lbs (w/o seating)

Seat Widths 8" to 16" (grows to 18")

Seat Depths 8" to 18" (grows to 20")

Front Seat-to-Floor Heights 13" to 19"

Rear Seat-to-Floor Heights 11" to 18.5"

Footrest/Frame Angle

80° (Zfinity™ adjusts to a wide range of footrest angles) Swing-in/Swing-out 80°, 90°

Transit Option: Yes

Conforms with ANSI/RESNA WC-19 standards. Refer to Quickie® & Zippie® WC-19 Transit Securement Supplement (MK-100072) for additional information. Supplement is available for download at www.sunrisemedical.com.

Aztec Gold	Black	Black Cherry	Black Opal
Blue Opal	Candy Blue	Candy Purple	Candy Red
Electric Blue	Evergreen	Glow	Green Apple
Hot Sparkle Pink	Mauve Pink	Rootbeer	Sparkle Silver
Sunrise Orange	Yellow	Titanium Color Paint	Matte Black
Matte Black Cherry	Matte Purple	Matte Electric Blue	Matte Evergreen
Mossy Oak Camo	Stars and Stripes	Zebra	 Kolorfusion not avail. on the crossbrace



©01.2017 Sunrise Medical (US) LLC MK-100118 Rev. B Click here for printer-friendly version

