



# JAY GO Cushion Price List



July 2016


Skin Protection Wheelchair Cushion  
Skin Protection Wheelchair Cushion 22" and over in width

## How to Order:

Part #	JAY GO Cushions	Price
<input checked="" type="checkbox"/> JGO1818	JAY GO Cushion 18" x 18"	\$263

Each JAY GO part number is simply based upon the size you are ordering. You may also reference the size matrix on page 1 of this price list.

Part # Example



### JAY GO Cushion Cover

The JAY GO cover is designed using dartex™ fabric which is known to be both stretchable and water resistant. This cover's sizing and sewing pattern is designed to reduce shear and peak pressure buildup. The material itself provides stretch which is an uncommon property in most incontinent products.



Figure A

Please visit [www.dartexcoatings.com](http://www.dartexcoatings.com) for additional information on the use of dartex coatings



### JAY GO Cushion Foam

The JAY GO foam base is designed to contain an accommodating seating surface in a thicker, softer dual layered foam base. The optimal combination of two different foam densities provides skin protection.



Cross section view

		JAY GO Size Matrix											
		WIDTH											
Depth		14	15	16	17	18	19	20	21	22	23	24	
	14												
	15												
	16												
	17												
	18												
	19												
	20												
	21												
	22												
	23												
	24												

**NO Modified sizes are available**

Size is available     Size is not available

### JAY GO Cushion Sealed Zipper

The JAY GO zipper and cover provide protection to the foam base against odor because the location is low and away from fluid. The water resistant material, combined with the Aquaguard™ zipper and anti wicking thread provides a thorough protection system for a cushion base.



Figure B

Optimized Aquaguard Zipper location

*The JAY GO Cushion is an excellent compliment to the JAY GO Back.*

Please see the **JAY GO Back** order form (part # 119284) for additional information.



**Each JAY GO cushion includes an outer cover, dual layered foam base, plastic bag with label and owner's manual**

Ordering Information PHONE (800) 333-4000, Fax (800) 300-7502, WEB: [www.sunrisemedical.com](http://www.sunrisemedical.com)

Specifications Subject to change without notice. Copyright © 2011 Sunrise Medical (8.11) Jay GO Cushion 119288 Rev A Page 1 of 2