



Quickie GP/GPV November 2016

EIGP70

Mark For: Date: _____ Dealer Acct #: _____ Dealer: _____ Dealer Contact: _____ Dealer Address: _____ Dealer City: _____ PR: _____ PC: _____ Dealer Phone: () Fax: () Confirmation Email: _____ Confirm Via: <input type="checkbox"/> Fax <input type="checkbox"/> Email	Submitting for: <input type="checkbox"/> Quote <input type="checkbox"/> Order PO#: _____ <hr/> <p style="text-align: center;">ADDITIONAL SHIPPING INFORMATION</p> Ship To: _____ Attention: _____ Address: _____ Address: _____ Ship To City: _____ PR: _____ PC: _____ Ship To Phone: () Fax: ()
--	---

	Option #	Definition	Price		Option #	Definition	Price
--	----------	------------	-------	--	----------	------------	-------

1. FRAME

Base Model			
<input type="checkbox"/>	024WM23	Quickie GP-70°	\$2,459.00
<input type="checkbox"/>	024WM23	Quickie GPV-70° (NA 12"-13")	\$2,459.00

Frame Angle			
<input type="checkbox"/>	024FA6	60° (NA 18 1/4" and 19 1/4" frame length)	NC
<input type="checkbox"/>	024FA7	70°	NC
<input type="checkbox"/>	024FA1	75° (NA with rigid footrest)	NC
<input type="checkbox"/>	024FA8	80° (NA with rigid footrest)	NC

2. SIZE / SEATING

Seat Width			
<i>Seat Width 1/2" Narrower</i>			
<input type="checkbox"/>	024W3	12"	NC
<input type="checkbox"/>	024W4	13"	NC
<input type="checkbox"/>	024W5	14"	NC
<input type="checkbox"/>	024W6	15"	NC
<input type="checkbox"/>	024W7	16"	NC
<input type="checkbox"/>	024W8	17"	NC
<input type="checkbox"/>	024W9	18"	NC
<input type="checkbox"/>	024W10	19"	\$268.00
<input type="checkbox"/>	024W11	20"	\$268.00
<input type="checkbox"/>	024W12	21"	\$268.00
<input type="checkbox"/>	024W13	22"	\$268.00

Frame Length			
<input type="checkbox"/>	024FL30	18-1/4" (Avail. on 12" Sling Depth)	NC
<input type="checkbox"/>	024FL31	19-1/4" (Avail. on 13" Sling Depth)	NC
<input type="checkbox"/>	024FL32	20-1/4" (Avail. on 14", 15" Sling Depth)	NC
<input type="checkbox"/>	024FL33	21-1/4" (Avail. on 14", 15", 16" Sling Depth)	NC
<input type="checkbox"/>	024FL34	22-1/4" (Avail. on 16", 17" Sling Depth)	NC
<input type="checkbox"/>	024FL35	23-1/4" (Avail. on 17", 18" Sling Depth)	NC
<input type="checkbox"/>	024FL36	24-1/4" (Avail. on 18" Sling Depth)	NC
<input type="checkbox"/>	024FL37	25-1/4" (Avail. on 19" Sling Depth)	NC
<input type="checkbox"/>	024FL38	26-1/4" (Avail. on 20" Sling Depth)	NC
<input type="checkbox"/>	024FL39	27-1/4" (Avail. on 21" Sling Depth)	NC
<input type="checkbox"/>	024FL40	28-1/4" (Avail. on 22" Sling Depth)	NC

Seat Depth			
<i>21 & 22" seat depths only available with 16" or wider chairs</i>			
<input type="checkbox"/>	024SD5	12"	NC
<input type="checkbox"/>	024SD6	13"	NC
<input type="checkbox"/>	024SD7	14"	NC
<input type="checkbox"/>	024SD8	15"	NC
<input type="checkbox"/>	024SD9	16"	NC
<input type="checkbox"/>	024SD10	17"	NC
<input type="checkbox"/>	024SD11	18"	NC
<input type="checkbox"/>	024SD72	19"	NC
<input type="checkbox"/>	024SD13	20"	\$268.00
<input type="checkbox"/>	024SD14	21"	\$268.00
<input type="checkbox"/>	024SD33	22"	\$268.00

Frame Height			
<i>See Appendix 11 for height chart. Short frame NA w/Xtender</i>			
<input type="checkbox"/>	024FH2	Short (front 9.25"-rear 8.25") NA with Xtender Power Assist	\$200.00
<input type="checkbox"/>	024FH5	Standard (front 10" - rear 9")	NC
<input type="checkbox"/>	024FH7	Tall (front 11" - rear 10")	\$200.00

Cushion			
<input type="checkbox"/>	024CU9	Reinforced Seat Sling	\$47.00
<input type="checkbox"/>	024CU1	2" Cushion	\$163.00
<input type="checkbox"/>	024CU2	Omit Cushion	NC
<input type="checkbox"/>	024CU4	Omit Cushion and Seat Sling	-\$19.00

JAY Bundled Cushions - Optional

To order any JAY cushions on-chair, please see JAY Price List. Cushion name or part number must be filled out with desired width & depth.

<input type="checkbox"/>	024JC20	JAY Bundled Cushion Cushion Name or Part #: _____ Width: _____ Depth: _____	Price List
--------------------------	---------	---	------------

3. CASTERS

Caster Wheels

<input type="checkbox"/>	024CW1	3" Micro	NC
<input type="checkbox"/>	024CW3	4" Low Profile Polyurethane	NC
<input type="checkbox"/>	024CW31	4" Micro Lighted	\$95.00
<input type="checkbox"/>	024CW6	5" Soft Roll Semi-Pneumatic	\$110.00
<input type="checkbox"/>	024CW5	5" Low Profile Polyurethane	NC
<input type="checkbox"/>	024CW7	5" x 1.5" Semi-Pneumatic <i>only avail. For use in bottom fork holes</i>	\$95.00
<input type="checkbox"/>	024CW33	6" Pneumatic	\$95.00
<input type="checkbox"/>	024CW34	6" Polyurethane	NC
<input type="checkbox"/>	024CW35	6" x 1.5" Semi-Pneumatic <i>only avail. For use in bottom fork holes</i>	\$95.00
<input type="checkbox"/>	024CW21	8" Pneumatic <i>Only available on 60° frame option</i>	\$126.00

Caster Wheels - Continued

<input type="checkbox"/>	024CW32	8" Polyurethane <i>Only available on 60° frame option</i>	NC
<input type="checkbox"/>	024CW36	8" x 2" Pneumatic <i>Only available on 60° frame option</i>	\$215.00
<input type="checkbox"/>	024CW37	8" x 2" Pneumatic w/Airless Insert <i>Only available on 60° frame option</i>	\$342.00

Caster Fork

<input type="checkbox"/>	024CS5	3/4" Longer Fork Stem Bolt	NC
<input type="checkbox"/>	024CS9	1-1/2" Longer Fork Stem Bolt	NC
<input type="checkbox"/>	024CS3	Caster Pin Locks	\$132.00

4. REAR WHEELS

Xtender Power Assist

NEW!

RSH 15"-19" with 22" wheels, RSH 16"-20" with 24" wheels. Spare batteries must match set up. Anti-Tips included. Dual Drive mode available w/B4ME

<input type="checkbox"/>	024XT6	22" Xtender Power Assist JWX-2 NICKEL	\$8,695.00
<input type="checkbox"/>	024XT7	22" Xtender Power Assist JWX-2 LITHIUM	\$10,095.00
<input type="checkbox"/>	024XT4	24" Xtender Power Assist JWX-2 NICKEL	\$8,695.00
<input type="checkbox"/>	024XT5	24" Xtender Power Assist JWX-2 LITHIUM	\$10,095.00

Please select one of the following handrim options

<input type="checkbox"/>	024XH5	Plastic Coated	NC
<input type="checkbox"/>	024XH6	Stainless Steel	NC

Xtender - Other Options

<input type="checkbox"/>	024XA10	Spare NiMH Power Pack - NICKEL	\$1,675.00
<input type="checkbox"/>	024XA11	Spare Li-ion Power Pack - LITHIUM	\$2,155.00

Rear Wheel Size

Spinergy NA w/ Xtender, spoke guards. Contact Built 4 Me for Spinergy Color Options

<input type="checkbox"/>	024RW8	22" Mag	NC
<input type="checkbox"/>	024RW9	22" Spoke	\$278.00
<input type="checkbox"/>	024RW1	24" Mag	NC
<input type="checkbox"/>	024RW2	24" Spoke	NC
<input type="checkbox"/>	024RW11	26" Spoke (NA with Xtender Power Assist)	\$79.00
<input type="checkbox"/>	024RW42	24" Spinergy (540) SPOX Black	\$751.00
<input type="checkbox"/>	024RR25	Omit Rear Wheel and Axle Intended Size for Set Up :	-\$70.00

Rear Wheel Tire

<input type="checkbox"/>	024RT1	Pneumatic	NC
<input type="checkbox"/>	024RT2	Pneumatic Airless Insert <i>NA with 26" wheels & Xtender</i>	\$215.00
<input type="checkbox"/>	024RT3	Full Poly <i>NA with 26" wheels</i>	NC
<input type="checkbox"/>	024RT6	High-Pressure Clincher	\$90.00

Hand Rim

<input type="checkbox"/>	024HR1	Aluminum Anodized	NC
<input type="checkbox"/>	024HR4	Plastic Coated - <i>NA on all rims</i>	\$126.00
<input type="checkbox"/>	024HR7	Omit Hand Rim	NC

Projection Hand Rim

Projections NA w/ plastic coated handrims. NA w/ Spoke guards

<input type="checkbox"/>	024HR8	8 Oblique (24"/26")	\$126.00
--------------------------	--------	---------------------	----------

Axle

<input type="checkbox"/>	024AX1	Quick Release Steel	NC
<input type="checkbox"/>	024AX11	Quad Quick Release	\$158.00

Axle Plate

<input type="checkbox"/>	024AL2	Standard Adjustable	NC
<input type="checkbox"/>	024AL13	Double Slotted Adjustable	NC
<input type="checkbox"/>	024AL6	Amputee Adjustable <i>NA with padded swing-away armrests or stroller handles</i>	\$126.00

5. WHEEL LOCKS

Wheel Lock

<input type="checkbox"/>	024WL1	High Mount Push	NC
<input type="checkbox"/>	024WL2	High Mount Pull	NC
<input type="checkbox"/>	024WL5	High Mount Scissor <i>NA w/ extension handles, grade aides</i>	NC
<input type="checkbox"/>	024WL12	Do Not Mount Wheel Locks	NC
<input type="checkbox"/>	024WL10	Omit	NC

Extension Handle

<input type="checkbox"/>	024EH1	6" Pair	\$63.00
--------------------------	--------	---------	---------

Grade Aid

<input type="checkbox"/>	024GA1	Grade Aid <i>NA w/ poly, high-pressure clincher tires, or Xtender</i>	\$137.00
--------------------------	--------	---	----------

6. FOOTREST

Footplates

Heel loop or leg strap not included, See Footplate Options Section

1. GP FOOTREST OPTIONS

<input type="checkbox"/>	024F56	Rigid NA w/Ext. Tubes; NA with 75° or 80° frame angle	NC
<input type="checkbox"/>	024F50	Platform NA w/Ext. Tubes	\$132.00
<input type="checkbox"/>	024F1	Composite Footplate	NC
<input type="checkbox"/>	024F3	Adult Angle Adj.	\$148.00
<input type="checkbox"/>	024F48	Platform Flip Up	\$200.00

2. GPV FOOTREST OPTIONS

<input type="checkbox"/>	024F56	Rigid NA w/Ext. Tubes; NA with 75° or 80° frame angle	NC
<input type="checkbox"/>	024F50	Platform NA w/Ext. Tubes	\$132.00
<input type="checkbox"/>	024F48	Platform Flip Up	\$200.00

Footrest Extension Tubes

<input type="checkbox"/>	024E1	Short (14" to 16-1/2")	NC
<input type="checkbox"/>	024E2	Medium (16-1/2" to 19")	NC
<input type="checkbox"/>	024E5	Long (19" to 21-1/2")	NC
<input type="checkbox"/>	024E16	High Mount Only avail. with adult angle adjustable footplates	\$84.00

Footplate Options (Ext Tube & Direct Mount)

<input type="checkbox"/>	024HL1	Heel Loop	NC
<input type="checkbox"/>	024HL5	Leg Strap	NC
<input type="checkbox"/>	024HL4	Omit Heel Loop	NC
<input type="checkbox"/>	024HL6	Omit Leg Strap	NC
<input type="checkbox"/>	024F7	Omit Footplate & Ext Tube	-\$16.00
<input type="checkbox"/>	024IG2	Impact Guard Neoprene	\$53.00

7. BACKREST

Back Type

<input type="checkbox"/>	024BT3	Non-Folding Backrest (Fixed) <i>NA with push handles or adj. upholstery</i>	
<input type="checkbox"/>	024BB1	8" Back Height	NC
<input type="checkbox"/>	024BB2	9" Back Height	NC
<input type="checkbox"/>	024BB3	10" Back Height	NC
<input type="checkbox"/>	024BB4	11" Back Height	NC
<input type="checkbox"/>	024BB5	12" Back Height	NC
<input type="checkbox"/>	024BB6	13" Back Height	NC
<input type="checkbox"/>	024BB7	14" Back Height	NC
<input type="checkbox"/>	024BB8	15" Back Height	NC
<input type="checkbox"/>	024BB9	16" Back Height	NC
<input type="checkbox"/>	024BT1	Folding Backrest (Adjustable)	
<input type="checkbox"/>	024BA1	11" - 14" Back Height	\$105.00
<input type="checkbox"/>	024BA2	12" - 16" Back Height	\$105.00
<input type="checkbox"/>	024BA3	14" - 18" Back Height	\$105.00
<input type="checkbox"/>	024BA4	16" - 20" Back Height	\$105.00

Back Upholstery

<input type="checkbox"/>	024BU1	Standard	NC
<input type="checkbox"/>	024BU8	Adjustable	\$178.00
<input type="checkbox"/>	024BU16	Reinforced (NA w/ adj. upholstery)	\$47.00
<input type="checkbox"/>	024BU14	Omit Back Upholstery	-\$19.00

Backrest Options

<input type="checkbox"/>	024BT33	Folding w/ Integral Push Handles	NC
<input type="checkbox"/>	024A85	Bolt-On Push Handle (Folding only)	\$121.00
<input type="checkbox"/>	024SR2	Regular Stroller Handle <i>NA with amputee axle plate</i>	\$168.00
<input type="checkbox"/>	024SR3	Tall Stroller Handle <i>NA with amputee axle plate</i>	\$168.00

JAY Back (Optional)

<input type="checkbox"/>	024JB01	J3 Back (See J3 Order Form)	Price List
--------------------------	---------	-----------------------------	------------

8. ARMREST / SIDEGUARDS

Single Post

<input type="checkbox"/>	024AA9	Single Post Ht. Adjust. Low <i>Must be removed to fold back down</i>	\$295.00
<input type="checkbox"/>	024AM1	Desk Length - Classic (10")	NC
<input type="checkbox"/>	024AM2	Full Length - Contour (14")	NC
<input type="checkbox"/>	024AA8	Single Post Ht. Adjust. Std <i>Must be removed to fold back down</i>	\$295.00
<input type="checkbox"/>	024AM1	Desk Length - Classic (10")	NC
<input type="checkbox"/>	024AM2	Full Length - Contour (14")	NC

Padded Swing Away Armrest

<i>All Padded Swing-Away NA with amputee axle plate or Xtender Power Assist</i>			
<input type="checkbox"/>	024AA19	Padded S-A Medium	\$273.00
<input type="checkbox"/>	024AA16	Padded S-A Adult	\$273.00

Side Guard

<input type="checkbox"/>	024SG3	Plastic Regular	\$163.00
<input type="checkbox"/>	024SG4	Plastic Kids	\$163.00

9. KEY OPTIONS AND ACCESSORIES

Rear Anti-Tip

<input type="checkbox"/>	024AT1	Rear Anti-Tip	\$105.00
--------------------------	--------	---------------	----------

Backpack

<input type="checkbox"/>	024A1	Quickie Backpack Black	\$97.00
<input type="checkbox"/>	024A2	Quickie Kids Backpack Black	\$90.00

Seat Pouch

<input type="checkbox"/>	024A8	Quickie Seat Pouch Black	\$73.00
--------------------------	-------	--------------------------	---------

Positioning Belts

<input type="checkbox"/>	024A142	Belt 1.5" Auto Buckle w/Clamps	\$84.00
<input type="checkbox"/>	024A12	Belt 2" Aircraft Buckle w/Clamps	\$121.00
<input type="checkbox"/>	024A13	Belt 2" Aircraft Buckle Pad w/Clamp	\$178.00

Spoke Guards

<i>NA w/ Xtender Power Assist; Only available w/ Mag and Spoke wheels</i>			
<input type="checkbox"/>	024A23	Black Spoke Guards	\$137.00
<input type="checkbox"/>	024A24	Clear Spoke Guards	\$137.00

Crutch Holder

<input type="checkbox"/>	024A53	Left Side	\$100.00
<input type="checkbox"/>	024A52	Right Side	\$100.00

Extra Leg Strap

<input type="checkbox"/>	024A55	Extra Leg Strap Black	\$53.00
--------------------------	--------	-----------------------	---------

Miscellaneous

<input type="checkbox"/>	024A83	Tool Kit	\$37.00
--------------------------	--------	----------	---------

Tray Table

<input type="checkbox"/>	024A31	Tray Table	\$331.00
--------------------------	--------	------------	----------

10. COLOR

Standard Frame Colors

<input type="checkbox"/>	024C140	Aztec Gold	NC
<input type="checkbox"/>	024C2	Black <i>(Not avail as base w/Zebra)</i>	NC
<input type="checkbox"/>	024C141	Black Cherry	NC
<input type="checkbox"/>	024C10	Black Opal	NC
<input type="checkbox"/>	024C38	Blue Opal	NC
<input type="checkbox"/>	024C36	Candy Blue	NC
<input type="checkbox"/>	024C14	Candy Purple	NC
<input type="checkbox"/>	024C35	Candy Red	NC
<input type="checkbox"/>	024C142	Electric Blue	NC
<input type="checkbox"/>	024C143	Evergreen	NC
<input type="checkbox"/>	024C126	Glow	NC
<input type="checkbox"/>	024C144	Green Apple	NC
<input type="checkbox"/>	024C147	Hot Sparkle Pink	NC
<input type="checkbox"/>	024C148	Mauve Pink	NC
<input type="checkbox"/>	024C103	Rootbeer	NC
<input type="checkbox"/>	024C114	Sparkle Silver	NC
<input type="checkbox"/>	024C145	Sunrise Orange	NC
<input type="checkbox"/>	024C113	Titanium Color Paint	NC
<input type="checkbox"/>	024C8	Yellow	NC

Matte Frame Colors

<input type="checkbox"/>	024C116	Matte Black	NC
<input type="checkbox"/>	024C151	Matte Black Cherry	NC
<input type="checkbox"/>	024C150	Matte Electric Blue	NC
<input type="checkbox"/>	024C149	Matte Evergreen	NC
<input type="checkbox"/>	024C152	Matte Purple	NC

Custom Frame Color

<input type="checkbox"/>	024C112	Desert Camouflage	\$320.00
<input type="checkbox"/>	024C43	Mossy Oak Camouflage	\$320.00
<input type="checkbox"/>	024C42	Stars & Stripes	\$320.00
<input type="checkbox"/>	024C108	<i>Zebra (Choose above base color) Std with White unless above base paint is chosen NA w/ Black, Matte Black, Black Opal or Titanium.</i>	\$320.00

11. BUNDLING



*Save time! Wheelchair and seating will be delivered together.
Hardware pre-installed.*

Cushions

<input type="checkbox"/>	Jay Cushion:		
Width:		Depth:	
Notes:			

Backs

<input type="checkbox"/>	Jay Back:		
Width:		Height:	
Notes:			

Headrest

<input type="checkbox"/>	System:	
Pads:		
Notes:		
Hardware:		
Notes:		

Jay Positioning Components

<input type="checkbox"/>	Anterior Trunk Supports:	
Size:		
Hand and Arm Supports:		
Size:		
Pelvic Positioning:		
Size:		
Foot Positioning:		
Size:		
Notes:		

12. APPENDIX

GP Series Rear STF Chart

**All measurements given are an estimate and have been rounded to the nearest .5"*

SHORT Frame (9 1/4" Front / 8 1/4" Rear)			Rear Axle Plate Hole Position					
Wheel	Rim	Tire	Standard			Inverted		
			Top	Middle	Bottom	Top	Middle	Bottom
22"	Mag / Spoke	Full Poly	14.5	15	15.5	16	16.5	17
22"	Mag / Spoke	Pneum/ Insert	14.5	15	16	16	17	17.5
24"	Mag / Spoke	Full Poly	15	15.5	16.5	16.5	17	17.5
24"	Spoke / Spinerger	Clincher	15	16	16.5	16.5	17	18
24"	Mag / Spoke / Spinerger	Pneum/ Insert	15.5	16	17	17	17.5	18
26"	Spoke / Spinerger	Clincher	16	17	17.5	17.5	18	19
26"	Spoke / Spinerger	Pneumatic	16.5	17	17.5	18	18.5	19

TALL Frame (11" Front / 10" Rear)			Rear Axle Plate Hole Position					
Wheel	Rim	Tire	Standard			Inverted		
			Top	Middle	Bottom	Top	Middle	Bottom
22"	Mag / Spoke	Full Poly	16	16.5	17.5	17.5	18	19
22"	Mag / Spoke	Pneum/ Insert	16.5	17	17.5	18	18.5	19
24"	Mag / Spoke	Full Poly	17	17.5	18	18.5	18.5	19.5
24"	Spoke / Spinerger	Clincher	17	17.5	18	18.5	19	19.5
24"	Mag / Spoke / Spinerger	Pneum/ Insert	17	17.5	18.5	19	19	20
26"	Spoke	Clincher	18	18.5	19	19.5	20	20.5
26"	Spoke	Pneumatic	18	19	19.5	19.5	20	21

STANDARD Frame (10" Front / 9" Rear)			Rear Axle Plate Hole Position					
Wheel	Rim	Tire	Standard			Inverted		
			Top	Middle	Bottom	Top	Middle	Bottom
22"	Mag / Spoke	Full Poly	15	15.5	16.5	16.5	17.5	18
22"	Mag / Spoke	Pneum/ Insert	15.5	16	17	17	17.5	18
24"	Mag / Spoke	Full Poly	16	16.5	17	17	18	18.5
24"	Spoke / Spinerger	Clincher	16	16.5	17	17.5	18	18.5
24"	Mag / Spoke / Spinerger	Pneum/ Insert	16	16.5	17.5	17.5	18	18.5
26"	Spoke / Spinerger	Clincher	17	17.5	18	18.5	19	19.5
26"	Spoke / Spinerger	Pneumatic	17	18	18.5	18.5	19.5	20

Maximum Dump 2 3/4"

*** Custom Frame Height Through Configurator - REAR 6 1/2 - 13" (Max frame squeeze = 3 1/2") Calculate Accordingly

GP Series Front STF Chart

**All measurements given are an estimate and have been rounded to the nearest .5"*

SHORT Frame (9 1/4" Front / 8 1/4" Rear)

Fork	Spacer	Fork Hole	3"		4" Poly		5" Poly 5" x 1.5		6" Poly		6" Pneu		8" Poly		8" Pneu		8x2	
			4" Lighted	5" Soft Roll	6" x 1.5													
4 1/4	STD	Bottom	15.50	16.00	16.50	17.00	17.00					x	x	x	x			
4 1/4	3/4	Bottom	16.50	17.00	17.00	17.50	18.00					x	x	x	x			
4 1/4	1 1/2	Bottom	17.00	17.50	18.00	18.50	18.50					x	x	x	x			
5 1/4	STD	Top	x	x	16.50	17.00	17.00				x	x	x	x				
5 1/4	STD	Bottom	16.50	17.00	17.5	18.00	18.00				19.00	20.5	x					
5 1/4	3/4	Top	x	x	17.50	18.00	18.00				x	x	x	x				
5 1/4	3/4	Bottom	17.50	18.00	18.50	18.50	19.00				19.50	20.00	20.00	x				
5 1/4	1 1/2	Top	x	x	18.00	18.50	18.50				x	x	x	x				
5 1/4	1 1/2	Bottom	18.00	18.50	19.00	19.50	19.50				20.50	20.50	20.50	x				
6	STD	Top	x	x	17.50	18.00	18.00				18.50	19.00	19.00	x				
6	STD	Bottom	17.50	18.00	18.50	19.00	19.00				19.50	20.00	20.00	x				
6	3/4	Top	x	x	18.00	18.50	18.50				19.50	19.50	19.50	x				
6	3/4	Bottom	18.00	18.50	19.00	19.50	19.50				20.50	20.50	20.50	x				
6	1 1/2	Top	x	x	19.00	19.50	19.50				20.00	20.50	20.50	x				
6	1 1/2	Bottom	19.00	19.50	20.00	20.50	20.50				21.00	21.00	21.00	21.50	x			
7	STD	Top	x	x	18.50	19.00	19.00				19.50	19.50	20.00	x				
7	STD	Bottom	18.50	19.00	19.50	20.00	20.00				20.50	21.00	21.00	x				
7	3/4	Top	x	x	19.00	19.50	19.50				20.50	20.50	20.50	x				
7	3/4	Bottom	19.00	19.50	20.00	20.50	20.50				21.50	21.50	21.50	x				
7	1 1/2	Top	x	x	20.00	20.50	20.50				21.00	21.50	21.50	x				
7	1 1/2	Bottom	20.00	20.50	21.00	21.00	21.00				22.00	22.00	22.50	x				
8x2	STD	Top	x	x	x	x	x				x	x	x	x				
8x2	STD	Bottom	x	x	x	x	x				x	x	x	x				
8x2	3/4	Top	x	x	x	x	x				x	x	x	x				
8x2	3/4	Bottom	x	x	x	x	x				x	x	x	x				
8x2	1 1/2	Top	x	x	x	x	x				x	x	x	x				
8x2	1 1/2	Bottom	x	x	x	x	x				x	x	x	x				

STANDARD Frame (10" Front / 9" Rear) ADD 3/4" to each measurement above

TALL Frame (11" Front / 10" Rear) ADD 1 3/4" to each measurement above

*** Custom Frame Height through Configurator - Front 9 - 13" - Calculate Accordingly (Minimum Front STF = 15 1/2" / Max = 26")

Frame Dimensions

