



# Quickie Shark

July 2017



Account: \_\_\_\_\_  
 Date: \_\_\_\_\_  
 PO#: \_\_\_\_\_  
 Buyer: \_\_\_\_\_  
 Marked For: \_\_\_\_\_

Ship To: \_\_\_\_\_  
 Name: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 City, Province, Postal Code: \_\_\_\_\_  
 Phone: \_\_\_\_\_

### QUICKIE SHARK

- 157WM38 Quickie Shark \$4,204.00  
*Only available with X-Link Fork*

**X-Link Fork Option**  
 Easier access-  
 Raises seat 3.94" (10 cm) higher for transfers  
 Moves the crank 5.91" (15 cm) forward for transfers

### FRAMES

- FRAME COLOR**
- 157C57 Black NC
  - 157C94 Tutti Frutti Orange NC
  - 157C95 White RAL NC
  - 157C90 Silver NC

- BRILLIANT COLOR**
- 157C73 Brilliant Lime NC
  - 157C72 Brilliant Marine NC
  - 157C74 Brilliant Red NC
  - 157C76 Brilliant Silver NC

- MATTE COLOR**
- 157C68 Matte Lime NC
  - 157C67 Matte Marine NC
  - 157C63 Matte Red NC
  - 157C62 Matte Silver NC


- FRAME LENGTH**
- 157FL70 Frame Length Short NC

**Frame Length Notes:**  
 Overall with:  
 20" Wheels 77.56" (197 cm)  
 26" Wheels 79.53" (202 cm)  
 without Rear Wheels 55.91" (142 cm)  
 without steering fork & wheels 41.34" (105 cm)  
 Transfer Seat Height 8.66" (22 cm)

### POSITION FORK MOUNT

- 157PZ01 15" Setting (88 cm) NC
- 157PZ02 15.75" Setting (40 cm) NC
- 157PZ03 16.5" Setting (42 cm) w/ adj. ball bearing NC

**Position Fork Mount Diagram:**  
 High Mount Fork Preset  
 @ 15.75" (40 cm)



### CONTROL OPTIONS

- HAND CRANK**
- 157CK4 Crank Standard V Crank NC
- HANDLE TYPE**
- 157HY11 Vertical Handle Small NC
  - 157HY12 Vertical Handle Std NC

- HAND MOUNT**
- 157HM6 Handle Right Side NC

**Control Notes:**  
**Standard V Crank**  
 6.89" (17.5 cm) long  
 19.29" (49 cm) wide (center to center of Handle)  
**Handle**  
 Vertical Small and Std @ 15 Degree Diagonal  
**Brake**  
 Parking Brake Std  
 Brake lever can be locked when brake is on

- SPEED**
- 157SA4 3 X 9 Gear (27 speed) NC

**Speed Notes:**  
 27 Speed with rim brake on left handle

## SEATING

### SEAT SLING

<input type="checkbox"/>	157SL1	Seat Sling Standard	NC
--------------------------	--------	---------------------	----

### SEAT WIDTH

<input type="checkbox"/>	157S34	15" Seat Width (36 cm)	NC
<input type="checkbox"/>	157S35	16.5" Seat Width (40 cm)	NC
<input type="checkbox"/>	157S36	18" Seat Width (44 cm)	NC


## BACKREST

### BACK TYPE

<input type="checkbox"/>	157BT63	Standard Backrest <i>Adj. height &amp; angle adj.</i>	NC
<input type="checkbox"/>	157BT64	Folding Backrest	\$157.00

#### Backrest Notes:


Standard Backrest (shown on right)  
22.83" (58 cm) long  
Height and Angle Adjustable  
(30 to 90 degree)  
Fold Down Backrest  
22.83" (58 cm) long  
Uses Quick Release Pin and  
Retainer Ring for ease of adjustment



### BACKREST POSITION

<input type="checkbox"/>	157BQ07	Adj (10"-21.5") Preset @ 12"	NC
--------------------------	---------	------------------------------	----


Backrest Position Diagram:  
Preset @ 12" (30.48 cm)



## NECK SUPPORT

<input type="checkbox"/>	157NT01	Neck Support	\$157.00
--------------------------	---------	--------------	----------

Neck Support Diagram:  
Height and Depth Adjustable



## FOOTREST

### FOOTREST

<input type="checkbox"/>	157F236	Adj Short (16"-31.5") (40-80 cm)	NC
<input type="checkbox"/>	157F237	Adj Long (27.5"-37.5") (70-95 cm)	NC

#### Footrest Notes:

Adjustable Short (16"-31.5")  
Adjustable Long (27.5"-37.5")

## WHEELS

### DRIVE WHEEL

<input type="checkbox"/>	157DW34	26 X 1" Hi Perf	\$75.00
<input type="checkbox"/>	157DW35	Omit Drive Wheel	-\$70.00

### REAR WHEEL SIZE

<input type="checkbox"/>	157WS4	20" Wheel	NC
<input type="checkbox"/>	157WS18	26" Wheel Dt Rader	\$304.00
<input type="checkbox"/>	157WS20	Omit Rear Wheels	-\$98.00

### REAR WHEEL TIRE

<input type="checkbox"/>	157RT20	1" Hi Performance Reifen	\$130.00
--------------------------	---------	--------------------------	----------

## AXLE

<input type="checkbox"/>	157AX37	Threaded Axle	NC
<input type="checkbox"/>	157AX38	Quick Release Axle	NC

## ACCESSORIES

### ACCESSORIES SPORTS

<input type="checkbox"/>	157AZ01	Front Fender	\$48.00
<input type="checkbox"/>	157AZ02	Rear Fender (20")	\$48.00
<input type="checkbox"/>		Rear Fender (26")	\$48.00
<input type="checkbox"/>	157AZ03	Front Headlight	\$48.00
<input type="checkbox"/>	157AZ04	Rear Light	\$48.00
<input type="checkbox"/>	157AZ05	Cushion W/Composite Foam	\$151.00
<input type="checkbox"/>	157AZ06	Waist Belt	\$151.00
<input type="checkbox"/>	157AZ07	Signal Flag	NC
<input type="checkbox"/>	157AZ08	Rear Bumper 20" Rear Wheels	\$211.00
<input type="checkbox"/>	157AZ09	Rear Bumper 26" Rear Wheels	\$211.00