



Quickie Shark

July 2016



Account: _____
 Date: _____
 PO#: _____
 Buyer: _____
 Marked For: _____

Ship To: _____
 Name: _____
 Address: _____
 City, Province, Postal Code: _____
 Phone: _____

QUICKIE SHARK		
<input type="checkbox"/>	157WM38 Quickie Shark <i>Only available with X-Link Fork</i>	\$4,082.00

X-Link Fork Option
 Easier access-
 Raises seat 3.94" (10 cm) higher for transfers
 Moves the crank 5.91" (15 cm) forward for transfers

FRAMES		
FRAME COLOR		
<input type="checkbox"/>	157C57 Black	NC
<input type="checkbox"/>	157C94 Tutti Frutti Orange	NC
<input type="checkbox"/>	157C95 White RAL	NC
<input type="checkbox"/>	157C90 Silver	NC

BRILLIANT COLOR		
<input type="checkbox"/>	157C73 Brilliant Lime	NC
<input type="checkbox"/>	157C72 Brilliant Marine	NC
<input type="checkbox"/>	157C74 Brilliant Red	NC
<input type="checkbox"/>	157C76 Brilliant Silver	NC


MATTE COLOR		
<input type="checkbox"/>	157C68 Matte Lime	NC
<input type="checkbox"/>	157C67 Matte Marine	NC
<input type="checkbox"/>	157C63 Matte Red	NC
<input type="checkbox"/>	157C62 Matte Silver	NC

FRAME LENGTH		
<input type="checkbox"/>	157FL70 Frame Length Short	NC

Frame Length Notes:
 Overall with:
 20" Wheels 77.56" (197 cm)
 26" Wheels 79.53" (202 cm)
 without Rear Wheels 55.91" (142 cm)
 without steering fork & wheels 41.34" (105 cm)
 Transfer Seat Height 8.66" (22 cm)

POSITION FORK MOUNT		
<input type="checkbox"/>	157PZ01 15" Setting (88 cm)	NC
<input type="checkbox"/>	157PZ02 15.75" Setting (40 cm)	NC
<input type="checkbox"/>	157PZ03 16.5" Setting (42 cm) w/ adj. ball bearing	NC

Position Fork Mount Diagram:
 High Mount Fork Preset
 @ 15.75" (40 cm)



CONTROL OPTIONS		
HAND CRANK		
<input type="checkbox"/>	157CK4 Crank Standard V Crank	NC
HANDLE TYPE		
<input type="checkbox"/>	157HY11 Vertical Handle Small	NC
<input type="checkbox"/>	157HY12 Vertical Handle Std	NC

HAND MOUNT		
<input type="checkbox"/>	157HM6 Handle Right Side	NC

Control Notes:
Standard V Crank
 6.89" (17.5 cm) long
 19.29" (49 cm) wide (center to center of Handle)
Handle
 Vertical Small and Std @ 15 Degree Diagonal
Brake
 Parking Brake Std
 Brake lever can be locked when brake is on

SPEED		
<input type="checkbox"/>	157SA4 3 X 9 Gear (27 speed)	NC

Speed Notes:
 27 Speed with rim brake on left handle

SEATING**SEAT SLING**

<input type="checkbox"/>	157SL1	Seat Sling Standard	NC
--------------------------	--------	---------------------	----

SEAT WIDTH


<input type="checkbox"/>	157S34	15" Seat Width (36 cm)	NC
<input type="checkbox"/>	157S35	16.5" Seat Width (40 cm)	NC
<input type="checkbox"/>	157S36	18" Seat Width (44 cm)	NC

BACKREST**BACK TYPE**

<input type="checkbox"/>	157BT63	Standard Backrest <i>Adj. height & angle adj.</i>	NC
<input type="checkbox"/>	157BT64	Folding Backrest	\$152.00


Backrest Notes:

Standard Backrest (shown on right)
 22.83" (58 cm) long
 Height and Angle Adjustable
 (30 to 90 degree)
 Fold Down Backrest
 22.83" (58 cm) long
 Uses Quick Release Pin and
 Retainer Ring for ease of adjustment

**BACKREST POSITION**


<input type="checkbox"/>	157BQ07	Adj (10"-21.5") Preset @ 12"	NC
--------------------------	---------	------------------------------	----

Backrest Position Diagram:
 Preset @ 12" (30.48 cm)

**NECK SUPPORT**

<input type="checkbox"/>	157NT01	Neck Support	\$152.00
--------------------------	---------	--------------	----------

Neck Support Diagram:
 Height and Depth Adjustable

**FOOTREST****FOOTREST**

<input type="checkbox"/>	157F236	Adj Short (16"-31.5") (40-80 cm)	NC
<input type="checkbox"/>	157F237	Adj Long (27.5"-37.5") (70-95 cm)	NC

Footrest Notes:

Adjustable Short (16"-31.5")
 Adjustable Long (27.5"-37.5")

WHEELS**DRIVE WHEEL**

<input type="checkbox"/>	157DW34	26 X 1" Hi Perf	\$73.00
<input type="checkbox"/>	157DW35	Omit Drive Wheel	-\$68.00

REAR WHEEL SIZE

<input type="checkbox"/>	157WS4	20" Wheel	NC
<input type="checkbox"/>	157WS18	26" Wheel Dt Rader	\$295.00
<input type="checkbox"/>	157WS20	Omit Rear Wheels	-\$95.00

REAR WHEEL TIRE

<input type="checkbox"/>	157RT20	1" Hi Performance Reifen	\$126.00
--------------------------	---------	--------------------------	----------

AXLE

<input type="checkbox"/>	157AX37	Threaded Axle	NC
<input type="checkbox"/>	157AX38	Quick Release Axle	NC

ACCESSORIES**ACCESSORIES SPORTS**

<input type="checkbox"/>	157AZ01	Front Fender	\$47.00
<input type="checkbox"/>	157AZ02	Rear Fender (20")	\$47.00
<input type="checkbox"/>		Rear Fender (26")	\$47.00
<input type="checkbox"/>	157AZ03	Front Headlight	\$47.00
<input type="checkbox"/>	157AZ04	Rear Light	\$47.00
<input type="checkbox"/>	157AZ05	Cushion W/Composite Foam	\$147.00
<input type="checkbox"/>	157AZ06	Waist Belt	\$147.00
<input type="checkbox"/>	157AZ07	Signal Flag	NC
<input type="checkbox"/>	157AZ08	Rear Bumper 20" Rear Wheels	\$205.00
<input type="checkbox"/>	157AZ09	Rear Bumper 26" Rear Wheels	\$205.00