



# Quickie TS EIQT

For Alberta Aids to Daily Living only  
AADL #W132



July 01, 2015

ACCOUNT #: \_\_\_\_\_ SHIP TO: \_\_\_\_\_  
 DATE: \_\_\_\_\_ NAME: \_\_\_\_\_  
 P.O. #: \_\_\_\_\_ ADDRESS: \_\_\_\_\_  
 CITY: \_\_\_\_\_ PROVINCE: \_\_\_\_\_  
 POSTAL CODE: \_\_\_\_\_ PHONE #: \_\_\_\_\_  
 BUYER: \_\_\_\_\_ MARKED FOR: \_\_\_\_\_



### QUICKIE TS

- 043WM35 **WMK5QU005** Quickie TS \$4,136.00  
 The Quickie TS is NA w/ seat or back upholstery nor does it include seating.

### BASE FRAME (Required for Contracture Footrest System)

- 043FS4 Growing Base NC

### SEAT FRAME

- 043FT49 Growing Seat Frame STD

### HEAVY DUTY OPTION (Maximum weight 350 lbs.)

- 43WT2 Heavy-Duty Option \$469.00  
 Use 24" Pneumatic rear tires and 8" x 2" pneumatic front tires only

### SEAT FRAME POSITION

- 043SP1 Std. Seat Frame Position NC  
 043SP2 -1" Seat Frame Position NC  
 043SP3 -2" Seat Frame Position NC

### TILT OPTION (NA w/ reclining back)

- 043TT1 Gas Cylinders NC  
 043TT2 Dual Mech Locks NC

### FRAME COLOUR

- 043C2 Black (Not available as base w/ Zebra) NC  
 043C36 Candy Blue NC

### FRAME WIDTH

- 043W5 14" NC  
 043W6 15" NC  
 043W7 16" NC  
 043W8 17" NC  
 043W9 18" NC  
 043W10 **WAMA** 19" \$260.00  
 043W11 **WAMA** 20" \$260.00

Product Info ▶

### FRAME DEPTH & BACK POST LOCATION

Dual post armrests NA w/ 14" short, 16" regular, or 19" long options

- 043FL5 Short**
- 043BL1 14" NC  
 043BL2 15" NC  
 043BL3 16" NC  
 043BL4 17" NC
- 043FL11 Regular**
- 043BL3 16" NC  
 043BL4 17" NC  
 043BL5 18" NC  
 043BL6 19" NC
- 043FL21 Long**
- 043BL6 **WAMB** 19"<sup>5</sup> NC  
 043BL7 **WAMB** 20" \$260.00  
 043BL8 **WAMB** 21" \$260.00  
 043BL9 **WAMB** 22" \$260.00

### SEATING

Attach Order Form for Seating For:

- Jay Fit Seating
  - Standard Conventional Seating
  - Custom Conventional Seating
- 043PN1 Aluminum Solid Seat Base \$112.00

### BACK TYPE

- 043BT8 Angle Adjustable Backrest NC  
 043BT5 Reclining Backrest \$587.00  
 NA w/ low, SA adj stroller handles, or dual mech locks

**STROLLER HANDLE BACK POSTS** (For Adaptive Seating)

<input type="checkbox"/>	043BP9	Low (25-3/4" Ht.)	NC
<input type="checkbox"/>	043BP10	Medium (28-3/4" Ht.)	NC
<input type="checkbox"/>	043BP11	Tall (31-3/4" Ht.)	NC
<input type="checkbox"/>	043BP12	WAMV Removable Stroller Handle Ext.	\$148.00
<input type="checkbox"/>	043BP15	Swing Away Adjustable Stroller	\$163.00

**ARMREST**

Dual post armrests NA w/ 14" short, 16" regular, or 19" long frame depth/backpost options

**ARMREST**

<input type="checkbox"/>	043AA26	Dual Post	NC
<input type="checkbox"/>	043AA39	Dual Post Height Adjustable	\$321.00
<input type="checkbox"/>	043AA43	Dual Post Flip Back	\$143.00
<input type="checkbox"/>	043AA48	Dual Post Flip Back Height Adjust.	\$286.00
<input type="checkbox"/>	043AA8	Height Adjustable Adult	\$286.00
<input type="checkbox"/>	043AA9	Height Adjustable Kids	\$286.00
<input type="checkbox"/>	043AA1	Adjustable Locking Flip-Up NA w/ Armpads	\$286.00

**ARMPAD**

<input type="checkbox"/>	043AM1	Desk Length (10")	NC
<input type="checkbox"/>	043AM2	Full Length (14")	NC

**FOOTREST** (Heel loop or leg strap not included, see footplate options section)

**1. HANGER MOUNTED FOOTREST**

(Footrests Mount to Hanger Tubing)

<input type="checkbox"/>	043H7	70° Swing-Away Front Mount	NC
<input type="checkbox"/>	043H8	80° Swing-Away Front Mount	NC
<input type="checkbox"/>	043H6	90° Swing-Away Front Mount	NC
<input type="checkbox"/>	043H9	70° Heavy-Duty Front Mount Lift-Off	\$179.00
<input type="checkbox"/>	043H10	80° Heavy-Duty Front Mount Lift-Off	\$179.00
<input type="checkbox"/>	043H23	WAMK 90° Elevating (NA w/ widths over 16")	\$286.00

**Extensions**

<input type="checkbox"/>	043E23	2"	\$46.00
<input type="checkbox"/>	043E16	High Mount (NA w/ ELRs)	\$82.00

**Footplate**

<input type="checkbox"/>	043F4	WAMJ Kids Angle Adjustable w/Leg Strap Only available on 14"-16" width frames	\$143.00
<input type="checkbox"/>	043F3	WAMJ Adult Angle Adjustable w/Heel Loops	\$143.00

**2. EXTENSION TUBE MOUNTED FOOTRESTS**

(Footrests Mount to Ext. Tubes)

<input type="checkbox"/>	043H1	60o Swing-Away	NC
<input type="checkbox"/>	043H2	70o Swing-Away	STD
<input type="checkbox"/>	043H14	WAMK Elevating Only available w/ aluminum and angle adj footplates	\$286.00
<input type="checkbox"/>	043H24	Articulating Adult NA w/ high-push wheel locks and short footrest ext tubes or 14" width.	\$444.00

**Extension Tube**

<input type="checkbox"/>	043E1	Short (11.5" - 14") Only available w/ aluminum and angle adj footplates	NC
<input type="checkbox"/>	043E2	Medium (14" - 17")	NC
<input type="checkbox"/>	043E5	Long (17" - 19")	NC

**Footplate**

<input type="checkbox"/>	043F1	Composite w/Heel Loops Only available with 60° or 70° hangers and ALRs	STD
<input type="checkbox"/>	043F5	Aluminum (ELRs only)	NC
<input type="checkbox"/>	043F3	WAMJ Adult Angle Adjust. w/Heel Loops	\$148.00
<input type="checkbox"/>	043F48	Platform Flip-Up w/Leg Strap	\$194.00

**FOOTREST OPTIONS**

<input type="checkbox"/>	043HL1	Heel Loop	NC
<input type="checkbox"/>	043HL5	Leg Strap NA w/ ELR, ALR, 90o ELR or contracture footrest	NC
<input type="checkbox"/>	043IG2	Impact Guard Neoprene* Available w/ 60o or 70o hangers with ext tubes only	\$51.00
<input type="checkbox"/>	043H32	Omit Hanger	-\$16.00
<input type="checkbox"/>	043F7	Omit Hanger & Extension Tube	-\$16.00
<input type="checkbox"/>	043HL4	Omit Heel Loop	NC
<input type="checkbox"/>	043HL6	Omit Single Leg Strap	NC

**3. CONTRACTURE FOOTREST SYSTEM\* (CFS)**

Multi-Adjustable Padded Hangers & Angle Adjustable Aluminum Footrests  
Best when used with Flared Front End  
Length measurement taken from top of seat rail to top of footrest plate when CFS is set at 90°. Footrest plate width determined by seat frame width. CFS may interfere with front caster rotation when Quickie TS is in non-tilt position. \*Client's responsibility

<input type="checkbox"/>	043H26	Contracture Footrest System	\$1,010.00
<b>043E25 Short</b>			
<input type="checkbox"/>	043EP1	10.5" Extension Height	NC
<input type="checkbox"/>	043EP2	11.0" Extension Height	NC
<input type="checkbox"/>	043EP3	11.5" Extension Height	NC
<input type="checkbox"/>	043EP4	12.0" Extension Height	NC
<input type="checkbox"/>	043EP5	12.5" Extension Height	NC
<input type="checkbox"/>	043EP6	13.0" Extension Height	NC
<b>043E26 Regular</b>			
<input type="checkbox"/>	043EP5	12.5" Extension Height	NC
<input type="checkbox"/>	043EP6	13.0" Extension Height	NC
<input type="checkbox"/>	043EP7	13.5" Extension Height	NC
<input type="checkbox"/>	043EP8	14.0" Extension Height	NC
<input type="checkbox"/>	043EP9	14.5" Extension Height	NC
<input type="checkbox"/>	043EP10	15.0" Extension Height	NC
<b>043E27 Long</b>			
<input type="checkbox"/>	043EP9	14.5" Extension Height	NC
<input type="checkbox"/>	043EP10	15.0" Extension Height	NC
<input type="checkbox"/>	043EP11	15.5" Extension Height	NC
<input type="checkbox"/>	043EP12	16.0" Extension Height	NC
<input type="checkbox"/>	043EP13	16.5" Extension Height	NC
<input type="checkbox"/>	043EP14	17.0" Extension Height	NC

**CASTERS**

<input type="checkbox"/>	043CW5	5" Low-Profile Polyurethane	STD
<input type="checkbox"/>	043CW7	5" x 1.5" Semi-Pneumatic Caster wheel interferes with caster fork in the upper hole position.	\$92.00
<input type="checkbox"/>	043CW33	6" Pneumatic	\$92.00
<input type="checkbox"/>	043CW34	6" Polyurethane	NC
<input type="checkbox"/>	043CW35	6" x 1.5" Semi-Pneumatic Caster wheel interferes with caster fork in the upper hole position.	\$92.00
<input type="checkbox"/>	043CW21	8" Pneumatic NA w/ 90° hangers and 90° ELR	\$122.00
<input type="checkbox"/>	043CW32	8" Polyurethane NA w/ 90° hangers and 90° ELR	NC
<input type="checkbox"/>	043CW39	8" x 1.5" Semi-Pneumatic Caster wheel interferes with caster fork in the upper hole position. NA w/ 90° hangers and 90° ELR	\$122.00
<input type="checkbox"/>	043CW36	8" x 2" Pneumatic NA w/ 90° hangers and 90° ELR	\$209.00
<input type="checkbox"/>	043CW37	8" x 2" Pneumatic w/Airless Ins.* NA w/ 90° hangers and 90° ELR	\$332.00

\* Not covered by AADL



**CASTER OPTIONS**

<input type="checkbox"/>	043CS5	3/4" Longer Fork Stem Bolt	NC
<input type="checkbox"/>	043CS9	1-1/2" Longer Fork Stem Bolt	NC
<input type="checkbox"/>	043CW12	Omit caster wheel and tire	NC
<input type="checkbox"/>	043CF30	Omit Caster Wheels and Forks	NC
<input type="checkbox"/>	043CS4	Quick-Release Caster Stems *	\$235.00
NA w/ 3/4" or 1-1/2" fork stem bolts or 8" x 2" front casters			

**AXLE AND WHEEL**

**AXLE**

<input type="checkbox"/>	043AX1	WAMU Quick Release Axles	NC
<input type="checkbox"/>	043RR25	Omit Rear Wheel and Axle	-\$67.00

**REAR WHEEL**

12", 16", and 18" wheels only available w/ Pneumatic and Pneumatic w/ Airless Inserts. 20", 22" only available w/ Pneumatic, Pneumatic w/ Airless Inserts and Full Poly 12" and 16" Mags NA w/ handrims

<input type="checkbox"/>	043RW13	12" Mag Wheel	\$71.00
<input type="checkbox"/>	043RW14	16" Mag Wheel	\$71.00
<input type="checkbox"/>	043RW3	20" Mag Wheel	NC
<input type="checkbox"/>	043RW8	22" Mag Wheel	NC
<input type="checkbox"/>	043RW1	24" Mag Wheel	NC

**REAR WHEEL TIRE**

<input type="checkbox"/>	043RT3	Full Polyurethane (NA w/ 12", 16", and 18")	NC
<input type="checkbox"/>	043RT1	Pneumatic	NC
<input type="checkbox"/>	043RT2	Pneumatic Airless Insert*	\$209.00

**HANDRIM**

NA w/ 12", 16", and 18"

<input type="checkbox"/>	043HR1	Aluminum Anodized	NC
<input type="checkbox"/>	043HR4	WAMO Plastic Coated	\$66.00
<input type="checkbox"/>	043HR7	Omit Handrim	-\$12.00
Does not apply to 12" & 16" wheels			

**PROJECTIONS (not available on plastic coated)**

<input type="checkbox"/>	043HR8	WAMQ 8 Oblique	\$128.00
--------------------------	--------	----------------	----------

**WHEEL LOCK**

<input type="checkbox"/>	043WL1	High Mount Push Wheel Lock	NC
<input type="checkbox"/>	043WL2	High Mount Pull Wheel Lock	NC
<input type="checkbox"/>	043WL11	Foot Wheel Lock (Only 12" whl)	\$87.00
<input type="checkbox"/>	043WL10	Omit Wheel Lock	NC

**GRADE AID**

<input type="checkbox"/>	043GA1	WAMR Grade Aid	\$133.00
NA w/ poly or high-pressure clincher			

**EXTENSION HANDLE**

<input type="checkbox"/>	043EH1	6" Pair	\$61.00
--------------------------	--------	---------	---------

**ANTI-TIP**

<input type="checkbox"/>	043AO1	Front Anti-Tip	\$133.00
<input type="checkbox"/>	043AT5	Rear Anti-Tip	NC

**ACCESSORIES**

**BACKPACK**

<input type="checkbox"/>	043A1	Black Utility Backpack*	\$94.00
<input type="checkbox"/>	043A2	Black/Charcoal Kids*	\$87.00

**SEAT POUCH**

<input type="checkbox"/>	043A10	Black	\$71.00
--------------------------	--------	-------	---------

**EXTRA LEG STRAP**

<input type="checkbox"/>	043A54	Black*	\$51.00
--------------------------	--------	--------	---------

**MISCELLANEOUS**

<input type="checkbox"/>	043A82	Caddy*	\$77.00
Only available w/ footrests with ext tubes			

**POSITIONING BELT**

<input type="checkbox"/>	043A15	1.5" Auto Buckle w/Clamps	\$82.00
--------------------------	--------	---------------------------	---------

**SIDE GUARD**

<input type="checkbox"/>	043SG3	WAMV Plastic Regular	\$158.00
<input type="checkbox"/>	043SG4	WAMV Plastic Kids	\$158.00

**JAY CUSHIONS (OPTIONAL)**

To order any JAY cushions on-chair, please see the JAY price list. Cushion name or part number must be filled out with desired width and depth

<input type="checkbox"/>	043JC20	JAY Bundled Cushion	Price List
		Cushion Name or Part #	
		Width: _____	Depth: _____

\* Not covered by AADL, clients responsibility.

