



# Zippie ZM-310™ Group 3

July 2016

Account: \_\_\_\_\_  
 Date: \_\_\_\_\_  
 PO#: \_\_\_\_\_  
 Buyer: \_\_\_\_\_  
 Marked For: \_\_\_\_\_

Ship To: \_\_\_\_\_  
 Name: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 City, Province, Postal Code: \_\_\_\_\_  
 Phone: \_\_\_\_\_

✓ Option #	Definition	Price	✓ Option #	Definition	Price
------------	------------	-------	------------	------------	-------

## 1. Zippie ZM-310™ GROUP 3 CATEGORY (Select One)

### Zippie ZM-310™ BC (Group 3)

Zippie ZM-310™ BC w/ASAPII Seat **\$6,850.00**

195PC106	Base Package	
----------	--------------	--

### Zippie ZM-310™ SC (Group 3)

Zippie ZM-310™ SC for SPOT Tilt **\$6,850.00**

195PC107	Base Package	
----------	--------------	--

#### BC Seating Specifications

<u>Typical Chair Weight</u>	<u>W/battery</u>	235lbs	<u>W/O Battery</u>	157lbs
<u>Transit</u>	WC19 w/optional crash tested belt.			
<u>Seat Type</u>	ASAP II / Optional Manual tilt and/or Manual Recline			
<u>Seat size range</u>	<u>Width</u>	10"-16"	<u>Depth</u>	10"-19"
<u>STFH range</u>	15.5" - 21.6"			
<u>Power Seating</u>	9" Seat elevator			

#### SC Seating Specifications

<u>Typical Chair Weight</u>	<u>W/battery</u>	260lbs	<u>W/O Battery</u>	182lbs
<u>Transit</u>	WC19 w/optional crash tested belt.			
<u>Seat Type</u>	ASAP II / Optional Manual Recline.			
<u>Seat size range</u>	<u>Width</u>	10"-16"	<u>Depth</u>	10"-19"
<u>STFH range</u>	17.6" - 21.8"			
<u>Power Seating</u>	50° CG tilt / optional 12" seat elevator			



#### Zippie ZM-310 Base Specifications

<u>Drive Wheel Position</u>	Midwheel Drive					
<u>User weight Cap</u>	165 lb					
<u>Base size</u>	<u>Width</u>	24"	<u>Length</u>	34"	<u>Turn Radius</u>	18"
<u>Drive Control</u>	VR-2 90A (non-expand)		R-net EL90 (Expandable)			
<u>Driving Range</u>	Up to 13.5 miles					
<u>Drive Motors</u>	4-pole 6 mph					

<u>Drive Wheels</u>	13" Drive Wheels w/ flat free insert (Black or Grey)		
<u>Tires/Casters</u>	7" Front / 6" rear (Black or Grey)		
<u>Charger</u>	8A Dual Mode Charger (Compact) CTE		
<u>Battery Size</u>	Group 22		
<u>Shroud Colors</u>	12 standard colors, 3 upgrade prints.		
<u>Ground Clearance</u>	3"		

## 2. BASE CHAIR OPTIONS

Shroud Color			
<input type="checkbox"/>	195EZ01	Candy Apple Red	N/C
<input type="checkbox"/>	195EZ02	High Gloss Black	N/C
<input type="checkbox"/>	195EZ03	Candy Blue	N/C
<input type="checkbox"/>	195EZ04	Silver	N/C
<input type="checkbox"/>	195EZ05	Copper	N/C
<input type="checkbox"/>	195EZ08	Green Apple	N/C
<input type="checkbox"/>	195EZ09	Yellow	N/C
<input type="checkbox"/>	195EZ10	Pearl Pink	N/C
<input type="checkbox"/>	195EZ11	Candy Purple	N/C
<input type="checkbox"/>	195EZ12	Matte Black	N/C
<input type="checkbox"/>	195EZ13	Midnight Blue	N/C
<input type="checkbox"/>	195EZ15	Pearl White	N/C
<input type="checkbox"/>	195EZ17	Pink Camo	\$515.00
<input type="checkbox"/>	195EZ18	Grunge Skulls	\$515.00

Drive wheel/tire package			
<input type="checkbox"/>	195DW60	Cast Wheel Solid Dual Density Inserts Black Knobby Tires and Black Casters	N/C

User weight			
<input type="checkbox"/>	195WG17	< 75	N/C
<input type="checkbox"/>	195WG18	76-100lbs	N/C
<input type="checkbox"/>	195WG19	101-125lbs	N/C
<input type="checkbox"/>	195WG20	126-165lbs	N/C

Transit Option (Required)			
<input type="checkbox"/>	195WN03	WC19 Transit Package (order belt separately below.	N/C

WC19 Transit Lap Belt			
<input type="checkbox"/>	195WN02	WC19 Transit Lap Belt	\$210.00

Pre-Set Seat Frame Angle			
<input type="checkbox"/>	195FA10	0°	N/C
<input type="checkbox"/>	195FA3	3° NA with +3" or +4" STFH or lift	N/C
<input type="checkbox"/>	195FA4	6° NA with +3" or +4" STFH or Lift/Tilt	N/C

Seat Height			
<i>(Measured in the front) - Obtain Quote to verify est. STFH</i>			
<input type="checkbox"/>	195HT44	Lowest STFH	NC
<input type="checkbox"/>	195HT45	+1" Seat to Floor Height	NC
<input type="checkbox"/>	195HT46	+2" Seat to Floor Height	NC
<input type="checkbox"/>	195HT47	+3" STFH	NC
<input type="checkbox"/>	195HT48	+4" STFH	NC

Motor Package			
<input type="checkbox"/>	195PF96	165lbs 6mph 4-Pole	Std

Ventilator Options			
<i>Not available w/ Transit</i>			
<input type="checkbox"/>	195VT4	Fixed Vent Tray (Pulmonetic LTV Models)	\$530.00
<input type="checkbox"/>	195VT9	Lift off Vent Tray (Pulmonetic LTV Models)	\$530.00

STFH Matrix						
Seating Type and Power Option Combinations	Seat-to-Floor Height at 0 tilt					
	0	+1"	+2"	+3"	+4"	
BC ASAP II No power seating.*	15.5	16.5	17.5	18.5	19.5	*For seat widths 12-16" add .5"
BC ASAP II 9" Lift only.	18.4	19.4	20.4	21.4	22.4	
BC ASAP II Manual Tilt	17.6	18.6	19.6	20.6	21.6	
SC ASAP II Tilt only	17.6	18.6	19.6	20.6	21.6	
SC ASAP II Tilt and 12" Lift combo.	17.8	18.8	19.8	20.8	21.8	

## 3. BATTERIES and CHARGERS

Batteries			
<input type="checkbox"/>	195AG01	22NF Std post Batteries (Pair) w/straps	\$620.00
<input type="checkbox"/>	195AG05	Omit batteries (Includes harness)	N/C

Battery Charger			
<input type="checkbox"/>	195BJ19	8A Dual Mode Charger (Compact) CTE	N/C

## 4. ELECTRONICS

Drive Control Options			
Non-Expandable Controllers			
<input type="checkbox"/>	195ER01	Non-expandable VR2-90 controller (NA with High Speed Motor package)	N/C
Expandable Controllers			
<input type="checkbox"/>	195ER02	R-Net Expandable Drive Control Upgrade (90A)	N/C
	195EX01	Expandable Controller - R-net	\$530.00
	195EX02	Harnessing for Expandable Controller	N/C

VR2 Joystick (non-expandable)			
<input type="checkbox"/>	195J94	Non Expandable VR2 Joystick	N/C

***R-net Joystick			
<input type="checkbox"/>	195J70	Std LED Joystick for Rnet Ctrl's	N/C
<input type="checkbox"/>	195J108	Joystick w/Stereo Jacks & Color Display	N/C
<input type="checkbox"/>	195J88	Joystick W/Toggle, Spd Pot, Stereo Jack & Color Display	\$158.00
<input type="checkbox"/>	195J112	Joystick Advanced: 3.5" Color Display, Paddle Switches, Shortcut Keys and Stereo Jack.	\$1,040.00
<input type="checkbox"/>	195J95	Omit JS. For use with OMNI	N/C

*Joystick Mount (Fixed, Swing Away, Retractable, Mid)			
<input type="checkbox"/>	195JM10	Fixed Mount	N/C
<input type="checkbox"/>	195JM11	Swing Away Mount	\$388.00
<input type="checkbox"/>	195JM15	Retractable Joystick mount w/45° stops	\$721.00
<input type="checkbox"/>	195JM14	Mid Line Mount	\$388.00

Joystick Mount (Location)			
<i>Height Adjustable.</i>			
<input type="checkbox"/>	195HM1	Right Side Mount	N/C
<input type="checkbox"/>	195HM2	Left Side Mount	N/C

Joystick Handle (Optional)			
<input type="checkbox"/>	195JH1	Ball Handle	\$95.00
<input type="checkbox"/>	195JH2	Foam Ball Handle	\$95.00
<input type="checkbox"/>	195JH4	T-Handle	\$95.00
<input type="checkbox"/>	195JH6	Goalpost	\$95.00

Attendant Control			
<input type="checkbox"/>	1950A01	Attendant Control for VR-2	\$1,066.00
<input type="checkbox"/>	1950A02	Attendant Control for R-net	\$1,066.00

Attendant Mounting location			
<input type="checkbox"/>	195DN1	Mounting location: Right Back Cane	N/C
<input type="checkbox"/>	195DN2	Mounting Location: Left Back Cane	N/C

R-net OMNI and Output Options			
<i>(For expandable controls, choose all that apply)</i>			
<input type="checkbox"/>	195ON1	Rnet OMNI Input Control Module and display with built in Infrared and armrest mount. Choose mount below. (required for Specialty Input Devices.)	\$2,075.00
<input type="checkbox"/>	195ON2	Rnet Output Module (ECM)	\$1,066.00

OMNI Mount			
<i>(If used with Joystick must be mounted opposite side)</i>			
<input type="checkbox"/>	195DS04	Swing Away Left Side OMNI Mount	N/C
<input type="checkbox"/>	195DS03	Swing Away Right Side OMNI Mount	N/C

Switch inputs for R-Net Joysticks or OMNI			
<i>Drive or Mode / Profile Control</i>			
<input type="checkbox"/>	195SW22	Stereo to 2 Mono Splitter	\$75.00
<input type="checkbox"/>	195SW6	Ribbon Switch	\$119.00
<input type="checkbox"/>	195SW1	Buddy Button Switch	\$95.00
<input type="checkbox"/>	195SW2	Disc Switch	\$75.00
<input type="checkbox"/>	195SW9	Wobble Switch	\$321.00
<input type="checkbox"/>	195SW4	Micro Light Switch	\$95.00
<input type="checkbox"/>	195SW24	Egg Switch-Black	\$119.00
<input type="checkbox"/>	195SW21	Dual toggle switch and bracket	\$214.00
<input type="checkbox"/>	195SW35	Switch-it Mini 1 button on Mono plug	\$268.00
<input type="checkbox"/>	195SW36	Switch-it Mini 2 button on 2 Mono plugs	\$268.00
<input type="checkbox"/>	195SW37	Switch-it Mini 2 button on 1 Stereo plug	\$268.00
<input type="checkbox"/>	195SW38	Switch-it Mini 4 button on DB9 connector	\$268.00

Single Switch Mount Bracket			
<input type="checkbox"/>	195EM01	Single Switch Mount Bracket	\$247.00

Electronic Programmers			
<input type="checkbox"/>	195PR2	VR-2 Hand Held Programmer	\$530.00
<input type="checkbox"/>	195PR5	R-net PC Programmer with Dongle	\$615.00

## 5. Specialty Driver Control Options

*Requires expandable Rnet controls and OMNI unless otherwise noted. Not available with Upholstered backs.*

### 5a. Proportional Control Packages

Head Control Proportional Packages			
<input type="checkbox"/>	Switch-it Dual Pro Proportional Head Control with Whitmyer and Facial Pads (Adult)		<b>\$6,953.00</b>
	195ID72	Switch-it Dual Pro Head Control w/facial pad (Adult)	\$6,551.00
	195J09	SA Boom Mount with Egg Switch Control	\$402.00
<input type="checkbox"/>	Switch-it Dual Pro Proportional Head Control with Spot Pads Package <b>NEW!</b>		<b>\$6,953.00</b>
	195ID70	Switch-it Dual Pro Head Control w/Spot pads	\$6,551.00
	195J09	SA Boom Mount with Egg Switch Control	\$402.00

Hand Control Proportional Joystick Packages			
<input type="checkbox"/>	Switch-it Versa-Guide Compact Joystick Midline with Gatlin Package		<b>\$2,369.00</b>
	195ID71	Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
	195J07	S/A Gatlin Mount	\$330.00
<input type="checkbox"/>	Switch-it Versa-Guide Compact Joystick Midline and prox flex switches with Gatlin Package		<b>\$3,018.00</b>
	195ID71	Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
	195IM08	Dual Flex Proximity Switches (requires no force)	\$649.00
	195J07	S/A Gatlin Mount	\$330.00
<input type="checkbox"/>	Switch-it Micro-Pilot Joystick with SpaceDisc Package		<b>\$4,007.00</b>
	195ID62	Micro-Pilot Proportional Joystick	\$3,534.00
	195J08	S/A SpaceDisc mount	\$473.00
<input type="checkbox"/>	Switch-it Micro-Guide Joystick with SpaceDisc Package		<b>\$4,831.00</b>
	195ID63	Micro-Guide Proportional Joystick	\$4,358.00
	195J08	S/A SpaceDisc mount	\$473.00
<input type="checkbox"/>	ASL Proportional Micro Extremity Joystick (ASL130) Midline with Gatlin Package		<b>\$4,414.00</b>
	195ID59	Micro Extremity Joystick	\$4,084.00
	195J07	S/A Gatlin Mount	\$330.00
<input type="checkbox"/>	ASL Mushroom Joystick (ALS135) Midline with Gatlin Package		<b>\$3,353.00</b>
	195ID60	Mushroom Joystick	\$3,023.00
	195J07	S/A Gatlin Mount	\$330.00
<input type="checkbox"/>	ASL Compact Joystick w Integrated buttons (ASL134) Midline with Gatlin Package		<b>\$1,530.00</b>
	195ID61	ASL Compact Joystick	\$1,200.00
	195J07	S/A Gatlin Mount	\$330.00
<input type="checkbox"/>	ASL Micro Mini Joystick w Hand Pad (ASL136) Midline with Gatlin Package		<b>\$4,012.00</b>
	195ID68	Micro Mini Joystick	\$3,682.00
	195J07	S/A Gatlin Mount	\$330.00

Chin Control Proportional Joystick Packages			
<input type="checkbox"/>	Switch-it Versa-Guide Chin Control Compact Joystick and Egg Switch Package		<b>\$2,154.00</b>
	195ID71	Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
	195IM07	S/A Mount off backrest	N/C
	195J04	SA Boom Mount with egg switch	\$115.00
<input type="checkbox"/>	Switch-it Versa-Guide Chin Control Compact Joystick and Proximity Flex Switches Package		<b>\$2,688.00</b>
	195ID71	Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
	195IM08	Dual Flex Proximity Switches (requires no force)	\$649.00
	195IM07	S/A Mount off backrest	N/C
<input type="checkbox"/>	Switch-It Versa-Guide Bib Chin Control Compact Joystick with Proximity Flex switches Package Std 6"		<b>\$3,445.00</b>
	195ID71	Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
	195IM08	Dual Flex Proximity Switches (requires no force)	\$649.00
	195IM05	Bib mount Small	\$757.00
<input type="checkbox"/>	Switch-It Versa-Guide Bib Chin Control Compact Joystick with Proximity Flex Switches Package Large 8"		<b>\$3,512.00</b>
	195ID71	Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
	195IM08	Dual Flex Proximity Switches (requires no force)	\$649.00
	195IM06	Bib mount Large	\$824.00
<input type="checkbox"/>	Switch-It Versa-Guide Bib Chin Control Compact Joystick with Egg Switch Package Std 6"		<b>\$2,911.00</b>
	195ID71	Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
	195IM05	Bib mount Small	\$757.00
	195J04	SA Boom Mount with egg switch	\$115.00
<input type="checkbox"/>	Switch-It Versa-Guide Bib Chin Control Compact Joystick with Egg Switch Package Large 8"		<b>\$2,978.00</b>
	195ID71	Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
	195IM06	Bib mount Large	\$824.00
	195J04	SA Boom Mount with egg switch	\$115.00
<input type="checkbox"/>	Switch-It Micro-Pilot Bib Chin Control with Egg Switch Package Std 6"		<b>\$4,406.00</b>
	195ID62	Micro-Pilot Proportional Joystick	\$3,534.00
	195IM05	Bib mount Small	\$757.00
	195J04	SA Boom Mount with egg switch	\$115.00
<input type="checkbox"/>	Switch-It Micro-Pilot Bib Chin Control with Egg Switch Package Large 8"		<b>\$4,473.00</b>
	195ID62	Micro-Pilot Proportional Joystick	\$3,534.00
	195IM06	Bib mount Large	\$824.00
	195J04	SA Boom Mount with egg switch	\$115.00
<input type="checkbox"/>	Switch-It Micro-Guide Bib Chin Control with Egg Switch Package Std 6"		<b>\$5,230.00</b>
	195ID63	Micro-Guide Proportional Joystick	\$4,358.00
	195IM05	Bib mount Small	\$757.00
	195J04	SA Boom Mount with egg switch	\$115.00
<input type="checkbox"/>	Switch-It Micro-Guide Bib Chin Control with Egg Switch Package Large 8"		<b>\$5,297.00</b>
	195ID63	Micro-Guide Proportional Joystick	\$4,358.00
	195IM06	Bib mount Large	\$824.00
	195J04	SA Boom Mount with egg switch	\$115.00
<input type="checkbox"/>	ASL Proportional Micro Extremity Joystick Chin / Lip with Egg Switch Control Package		<b>\$4,199.00</b>
	195ID59	Micro Extremity Joystick	\$4,084.00
	195IM07	S/A Mount off backrest	N/C
	195J04	SA Boom Mount with egg switch	\$115.00

## 5. Specialty Driver Control Options (cont.)

### 5b. Switch Control

*Requires Headrest and/or Headrest Mount*

#### 3 Switch Proximity Head Arrays *Choose headrest mount below*

<input type="checkbox"/>	<b>ASL Proximity Head Array Package (Adult)</b>		<b>\$5,755.00</b>
	195II02	ASL Proximity Head Array (Adult)	\$5,353.00
	195IJ09	SA Boom Mount with Egg Switch Control	\$402.00
<input type="checkbox"/>	<b>ASL Proximity Head Array Package (Ped)</b>		<b>\$5,353.00</b>
	195II03	ASL Proximity Head Array (Ped)	\$4,951.00
	195IJ09	SA Boom Mount with Egg Switch Control	\$402.00
<input type="checkbox"/>	<b>Switch-It Head Array Package 3 Switch Curved (Long)</b>		<b>\$5,353.00</b>
	195II13	Switch-it Proximity Head Array Long	\$4,951.00
	195IJ09	SA Boom Mount with Egg Switch Control	\$402.00
<input type="checkbox"/>	<b>Switch-It Head Array Package 3 Switch Curved (Stubby)</b>		<b>\$5,691.00</b>
	195II14	Switch-it Proximity Head Array Stubby	\$5,289.00
	195IJ09	SA Boom Mount with Egg Switch Control	\$402.00

#### 5 Switch Proximity Head Array (No Egg Switch Needed) *Choose "Headrest Mount for Special Controls"*

<input type="checkbox"/>	195II04	Switch-It Head Array 5 switch flat	\$5,353.00
--------------------------	---------	------------------------------------	------------

#### Multiple Switches.

<input type="checkbox"/>	195SW31	Switch-it 3 Naked proximity switches on DB9 connector. (requires OMNI)	\$3,574.00
<input type="checkbox"/>	195SW32	Switch-it 4 Naked proximity switches on DB9 connector. (requires OMNI)	\$4,398.00
<input type="checkbox"/>	195SW33	Switch-it 5 Naked proximity switches on DB9 connector. (requires OMNI)	\$4,954.00
<input type="checkbox"/>	195SW34	ASL (502) 5-Switch adaptor box on DB9 connector. (for use with 5 mono switches, not included)	\$4,954.00
<input type="checkbox"/>	195SW23	4-way toggle switch button DB9 Connector. Includes bracket	\$258.00

#### Sip N Puff Package *Choose headrest and HR mount below*

<input type="checkbox"/>	<b>Quickie Sip n' Puff System</b>		<b>\$840.00</b>
	195ID34	Sip N' Puff Kit	\$268.00
	195IH1	Sip n' Puff Tubing and SA Hardware	\$457.00
	195IJ10	SA Boom Mount with egg switch	\$115.00

#### Headrest Mount for Specialty Controls

<input type="checkbox"/>	195MH12	Axys Mount & Bracket for Specialty Controls	\$309.00
<input type="checkbox"/>	195MH13	Cobra Mount & Bracket for Specialty Controls	\$386.00
<input type="checkbox"/>	195MH14	Linx Mount & Bracket for Specialty Controls	\$309.00

#### Headrest Pad for Specialty controls

<input type="checkbox"/>	195HE79	Plush Pad 8" for Specialty controls	\$309.00
<input type="checkbox"/>	195HE80	Plush Pad 10" for Specialty controls	\$309.00
<input type="checkbox"/>	195HE81	Plush Pad 14" for Specialty controls	\$309.00

#### Extra Swing Away Head Boom Egg Switch and Mount

	195IJ05	Extra SA Boom Mount with egg switch #1	\$115.00
	195IJ06	Extra SA Boom Mount with egg switch #2	\$115.00

#### Single Switches

<input type="checkbox"/>	195IB01	Ribbon Switch	\$119.00
<input type="checkbox"/>	195IB02	Buddy Button Switch	\$95.00
<input type="checkbox"/>	195IB03	Disc Switch	\$202.00
<input type="checkbox"/>	195IB04	Wobble Switch	\$630.00
<input type="checkbox"/>	195IB05	Micro Light Switch	\$95.00
<input type="checkbox"/>	195IB06	Egg Switch-Black	\$119.00
<input type="checkbox"/>	195IB07	Switch-it Mini 1 button on Mono plug	\$268.00

## 6. SEAT TYPE

<input type="checkbox"/>	195ST56	A.S.A.P II (Adjustable Seating and Positioning) Choose Seat Frame Size Platform.	N/C
--------------------------	---------	--	-----

### ASAP II Frame size

*(each frame is width and depth adjustable within given range. )*

<input type="checkbox"/>	195FU16	Adjustable 10x10 to 13x16 ( Z )	N/C
<input type="checkbox"/>	195FU17	Adjustable 12x12 to 16X19 ( A )	N/C

## 7. SEAT SIZE

### Seat Size Matrix. *Seat type is chosen in section 5.*

Z = ASAP II Seat Frame Adjustable from 10x10 to 13x16

A = ASAP II Seat Frame Adjustable from 12x12 to 16X19

		ASAPII Seat Width							
		10	11	12	13	14	15	16	
Seat Depth	10	Z	Z	Z	Z				
	11	Z	Z	Z	Z				
	12	Z	Z	Z,A	Z,A	A	A	A	A
	13	Z	Z	Z,A	Z,A	A	A	A	A
	14	Z	Z	Z,A	Z,A	A	A	A	A
	15	Z	Z	Z,A	Z,A	A	A	A	A
	16	Z	Z	Z,A	Z,A	A	A	A	A
	17				A	A	A	A	A
	18				A	A	A	A	A
19				A	A	A	A	A	

Seat Width			
<input type="checkbox"/>	195R45	10" Seat width	N/C
<input type="checkbox"/>	195R46	11" Seat width	N/C
<input type="checkbox"/>	195R21	12" Seat width	N/C
<input type="checkbox"/>	195R22	13" Seat width	N/C
<input type="checkbox"/>	195R23	14" Seat width	N/C
<input type="checkbox"/>	195R24	15" Seat width	N/C
<input type="checkbox"/>	195R5	16" Seat width	N/C

Seat Depth			
<input type="checkbox"/>	195RF40	10" Seat depth	N/C
<input type="checkbox"/>	195RF41	11" Seat depth	N/C
<input type="checkbox"/>	195RF20	12" Seat depth	N/C
<input type="checkbox"/>	195RF21	13" Seat depth	N/C
<input type="checkbox"/>	195RF22	14" Seat depth	N/C
<input type="checkbox"/>	195RF23	15" Seat depth	N/C
<input type="checkbox"/>	195RF5	16" Seat depth	N/C
<input type="checkbox"/>	195RF6	17" Seat depth	N/C
<input type="checkbox"/>	195RF7	18" Seat depth	N/C
<input type="checkbox"/>	195RF8	19" Seat depth	N/C

## 8. SEAT CUSHIONS

### JAY/Zip Cushions (Optional)

*To order any JAY/Zip cushions on-chair, please see the JAY/Zip price list.*

*Cushion name or part number must be filled out w/ desired width & depth*

<input type="checkbox"/>	195JC20	JAY/Zip Bundled Cushion		Price List
		Cushion Name or Part #		
		Width:	Depth:	

## 9. SEATING SYSTEMS

### Non-Powered Seating

<input type="checkbox"/>	195PJ8	Manual Tilt 45° CG Tilt <i>Must be ordered with Stroller handle NA with height adjustable back cane</i>	\$876.00
--------------------------	--------	--	----------

### Power Seating Options (SC Models)

<input type="checkbox"/>	195PJ2	SPOT 50° CG Tilt	\$4,523.00
--------------------------	--------	------------------	------------

### Power Seating Control

#### Powered Seating Switch Control Method

*(For use with power seating actuators. Switch activation is stand-alone and does not require expandable or thru-drive electronics.)*

#### Single Actuator Switch Control (Lift or Tilt)

<input type="checkbox"/>	195EO10	Switch Control 1 actuator with <u>Dual Toggle</u>	N/C
<input type="checkbox"/>	195EO11	Switch Control 1 actuator with <u>Buddy Button</u>	N/C
<input type="checkbox"/>	195EO19	Switch Control 1 actuator with <u>Switch-it Mini 1 button</u>	N/C
<input type="checkbox"/>	195EO20	Switch Control 1 actuator with <u>Switch-it Mini 2 button</u>	N/C

#### 2-3 Actuator Switch Control (Lift and Tilt)

<input type="checkbox"/>	195EO21	Switch Control 2-3 actuator with <u>Dual Toggle</u>	\$750.00
<input type="checkbox"/>	195EO12	Switch Control 2-3 actuator with <u>4-way toggle</u>	\$750.00
<input type="checkbox"/>	195EO23	Switch Control 2-3 actuator with <u>Switch-it Mini 2 button</u>	\$750.00
<input type="checkbox"/>	195EO22	Switch Control 2-3 actuator with <u>Switch-it Mini 4 button</u>	\$750.00

#### Switch Mount

*(for use with switch control powered seating)*

<input type="checkbox"/>	195EY01	Switch Mount Right	N/C
<input type="checkbox"/>	195EY02	Switch Mount Left	N/C

### Powered Seat Elevator- (AKA: Power Lift)

<input type="checkbox"/>	195PA8	Power Seat Elevate (9")	\$3,414.00
<input type="checkbox"/>	195PA15	Power Seat Elevate for use with tilt (12") <i>(Must be ordered with SPOT tilt)</i>	\$3,257.00

#### Powered Seating Thru-drive Control Method.

##### Thru-Drive 1 or 2 Actuator

<input type="checkbox"/>	195EO26	Thru Drive 1 or 2 Actuator(s) For operating Tilt and/or Lift thru Joystick. Includes Ctrl+5 when ordered with Expandable Controls (Rnet)	\$1,751.00
--------------------------	---------	--	------------

#### Powered Seating Thru-drive Control Options

Only available with Expandable Controls (Rnet)

##### Assignable Buttons

"Assignable Button" feature allows for customization of joystick buttons and/or Ctrl+5 Switches allowing access to various wheelchair and seating functions. When enabled, chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle.

<input type="checkbox"/>	195PM01	Enable Assignable buttons	\$618.00
--------------------------	---------	---------------------------	----------

##### Ctrl+5 Switchbox Options



Ctrl+5 allows direct control to 5 power seating functions (10 total switches). Switches can also be used with assignable buttons for greater flexibility. Box includes 3-Amp USB port for charging portable devices

<input type="checkbox"/>	195QG01	Ctrl+5 Switchbox mounted behind joystick	\$618.00
<input type="checkbox"/>	195QG02	Ctrl+5 Switchbox mounted opposite joystick	\$618.00
<input type="checkbox"/>	195QG03	Ctrl+5 Switchbox mounted on back cane	\$618.00
<input type="checkbox"/>	195QG05	Omit Ctrl+5	-\$618.00

#### Ctrl+5 Switch Type

<input type="checkbox"/>	195EY03	Push buttons	N/C
<input type="checkbox"/>	195EY04	Toggle switches	\$309.00

## 10. BACKREST AND BACK OPTIONS FOR ASAP II SEAT

### Backrest Frame Type

<input type="checkbox"/>	195BT8	Angle Adjustable	N/C
<input type="checkbox"/>	195BT98	Manual Recline 40° (not width adjust.)	\$536.00

### Manual Recline Adjustment Method Option

*Back Cane Height 23". Fits all seat types. NA with Cantilever Arms*

<input type="checkbox"/>	195MN01	Tool	N/C
<input type="checkbox"/>	195MN02	Knob (NA with dual post arms.)	N/C

### Angle Adjustable Backrest Options

*(NA with Manual recline)*

<input type="checkbox"/>	195OD8	Non-Folding Backrest	N/C
<input type="checkbox"/>	195OD9	Folding Backrest	\$137.00

### Cane Angle

<input type="checkbox"/>	195BP40	Straight Back Post	N/C
<input type="checkbox"/>	195BP41	8° Bend Back Post (Avail. only on Zip or J3 Back) NA with Manual recline	\$126.00

### Cane Adjustable height

<input type="checkbox"/>	195B1	Fixed Low (15.5")	N/C
<input type="checkbox"/>	195B5	Fixed Med (17.5")	N/C
<input type="checkbox"/>	195B9	Fixed Tall (19.5")	N/C
<input type="checkbox"/>	195B243	Adjustable Low 14" - 18"	N/C
<input type="checkbox"/>	195B244	Adjustable High 20" - 26"	N/C

### Towel Bar Depth

<input type="checkbox"/>	195QT03	Shallow	N/C
<input type="checkbox"/>	195QT04	Deep (Ext. 4" from back canes for deep backs)	N/C

### Back Type (Solid)

<input type="checkbox"/>	195BU86	S.P.O.T. Curved back with rail for mounting laterals	N/C
<input type="checkbox"/>	195BU87	Jay SPO Posterior (Curved back with T-Foam inserts and mounting rails)	\$600.00
<input type="checkbox"/>	195BU14	Omit Backrest	-\$15.00

### Jay 3 / ZIP Option

<input type="checkbox"/>	195JB01	Jay J3/Zip (select OMIT above)	See J3 OF
--------------------------	---------	--------------------------------	-----------

### Back Height - Required for SPOT Curved and SPO Backs.

*(Limited combinations of widths) Measured from pan to top of solid back.*

<input type="checkbox"/>	195B118	12"	Available on 10" - 16" wide seats	N/C
<input type="checkbox"/>	195B132	14"		N/C
<input type="checkbox"/>	195B36	16"	Available on 12" - 16" wide seats	N/C
<input type="checkbox"/>	195B38	18"		N/C
<input type="checkbox"/>	195B40	20"		N/C

### Push Handle Options



<input type="checkbox"/>	195PH14	Stroller Handle: Single height adjustable Bar Push Handle. w/2 utility hooks. Required with Manual Tilt (NA with adjustable height cane.)	\$57.00
--------------------------	---------	---	---------



## 11. HEADREST OPTIONS

Choose One Each.

Headrest Mount			
<input type="checkbox"/>	195MH10	Axys Mount & Bracket	\$315.00
<input type="checkbox"/>	195MH11	Cobra Mount and Bracket	\$578.00
<input type="checkbox"/>	195MH1	Linx Mount and Bracket	\$237.00



Headrest Pad			
<input type="checkbox"/>	195HE65	Plush Head Support Pad 8"	\$237.00
<input type="checkbox"/>	195HE66	Plush Head Support Pad 10"	\$237.00
<input type="checkbox"/>	195HE82	Adjust-a-plush pad Narrow	\$365.00
<input type="checkbox"/>	195HE83	Adjust-a-plush pad Standard	\$365.00
<input type="checkbox"/>	195HE84	Contour Cradle Small	\$365.00
<input type="checkbox"/>	195HE85	Contour Cradle Medium	\$365.00
<input type="checkbox"/>	195HE11	Omit headrest pad only.	N/C

## 12. POSTURAL SUPPORTS

LATERAL THORACIC SUPPORT (Swing Away)			
<i>Available on all Solid Backs.</i>			
<b>LAT 3x3 S/A FLUSH MOUNTS - PR \$578.00</b>			
<input type="checkbox"/>	195ZV20	LATERAL 3x3 Flush Mount - PAIR	\$295.00
<input type="checkbox"/>	195ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$283.00
<b>LAT 4X4 S/A MOUNTS W/ 1" OFFSET - PR \$409.00</b>			
<input type="checkbox"/>	195ZV14	LATERAL 4X4 W/ 1" OFFSET - PAIR	\$126.00
<input type="checkbox"/>	195ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$283.00
<b>LAT 4X5 S/A MOUNTS W/ 1" OFFSET - PR \$409.00</b>			
<input type="checkbox"/>	195ZV15	LATERAL 4X5 W/ 1" OFFSET - PAIR	\$126.00
<input type="checkbox"/>	195ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$283.00
<b>LAT 4X6 S/A MOUNTS W/ 1" OFFSET - PR \$409.00</b>			
<input type="checkbox"/>	195ZV16	LATERAL 4X6 W/ 1" OFFSET - PAIR	\$126.00
<input type="checkbox"/>	195ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$283.00

Thigh Support w/Removable Mount			
<i>Available with Cantilever Arms, NA with Manual Recline Back, Dual Post, Single Post or Cantilever w/ 2 pt Locking arms.</i>			
<b>Jay Rectangular Thigh Support - Left</b>			
<input type="checkbox"/>	195TV3	4" x 5" Thigh Support - Left	\$118.00
<input type="checkbox"/>	195TX1	Removable Mount - Left	\$283.00
<b>Jay Rectangular Thigh Support - Right</b>			
<input type="checkbox"/>	195TV4	4" x 5" Thigh Support - Right	\$118.00
<input type="checkbox"/>	195TX2	Removable Mount - Right	\$283.00

Jay Anterior Trunk Support Style (Choose Size below)			
<input type="checkbox"/>	195JJ12	Classic Dynamic	\$210.00
<input type="checkbox"/>	195JJ13	Classic Structured	\$210.00
<input type="checkbox"/>	195JJ14	Contoured Dynamic	\$210.00
<input type="checkbox"/>	195JJ15	Contoured Structured	\$210.00
<input type="checkbox"/>	195ZR03	Padded Chest Strap	\$283.00

Jay Anterior Trunk Support Size			
<input type="checkbox"/>	195JK01	Small	NC
<input type="checkbox"/>	195JK02	Medium	NC

JAY Pelvic Positioning Belts			
<input type="checkbox"/>	195ZZ19	2 Point Pelvic Belt 1' Push Button Padded	\$127.00
<input type="checkbox"/>	195ZZ20	2 Point Pelvic Belt 1' Side Release Padded	\$77.00
<input type="checkbox"/>	195ZZ21	2 Point Pelvic Belt 1.5" Push Button Padded	\$121.00
<input type="checkbox"/>	195ZZ22	2 Point Pelvic Belt 1.5" Side Release Padded	\$77.00
<input type="checkbox"/>	195ZZ23	4 Pnt V Contour Pelvic Belt 1.5" Push Btn	\$187.00

Swing Away Knee Adductor Supports			
<i>Available only with Fixed and Power Center Mount Legrests.</i>			
<i>Not available with Dual Post Arm Rests or 2-point Locking Cantilever Arms</i>			

Knee Adductor Supports - LEFT (choose Pad and Arm)			
<b>Pad Size</b>			
<input type="checkbox"/>	195KB5	Knee Adductor Pad SM 3"- Left	\$215.00
<b>S/A Arm</b>			
<input type="checkbox"/>	195XX04	Swing Away Short - Left	\$283.00

Knee Adductor Supports - RIGHT (choose Pad and Arm)			
<b>Pad Size</b>			
<input type="checkbox"/>	195KB6	Knee Adductor Pad SM 3"- Right	\$215.00
<b>S/A Arm</b>			
<input type="checkbox"/>	195XX05	Swing Away Short - Right	\$283.00

Medial Thigh Support			
<input type="checkbox"/>	<b>Small Jay Style Wedge Medial Thigh Support (Height adjustable)</b>		<b>\$500.00</b>
<input type="checkbox"/>	195KH12	Small Wedge Pad (3" x 4")	\$201.00
<input type="checkbox"/>	195LM3	Push Button Flip Down Bracket	\$299.00
<input type="checkbox"/>	<b>Medium Jay Style Wedge Medial Thigh Support (Height adjustable)</b>		<b>\$500.00</b>
<input type="checkbox"/>	195KH13	Medium Wedge Pad (4.5" x 5")	\$201.00
<input type="checkbox"/>	195LM3	Push Button Flip Down Bracket	\$299.00
<input type="checkbox"/>	<b>Small Jay Style OVAL Medial Thigh Support (Height adjustable)</b>		<b>\$500.00</b>
<input type="checkbox"/>	195KH14	Small Wedge Pad (3" x 4")	\$201.00
<input type="checkbox"/>	195LM3	Push Button Flip Down Bracket	\$299.00
<input type="checkbox"/>	<b>Medium Jay Style Oval Medial Thigh Support (Height adjustable)</b>		<b>\$500.00</b>
<input type="checkbox"/>	195KH15	Medium Wedge Pad (4.5" x 5")	\$201.00
<input type="checkbox"/>	195LM3	Push Button Flip Down Bracket	\$299.00

### 13. ARMREST AND ARMPAD OPTIONS

Armrest Mounting for K0848 Models with ASAP II			
<b>Dual Post Armrests</b>			
<input type="checkbox"/>	195EG01	Short Dual Post Flip Back 10" Height.	N/C
<input type="checkbox"/>	195EG21	Long Dual Post Flip Back 10" Height.	N/C
<input type="checkbox"/>	195EG02	Shrt Dual Post Flip Back HT Adj 10"-15"	\$378.00
<input type="checkbox"/>	195EG22	Lng Dual Post Flip Back HT Adj 10"-15"	\$378.00
<b>Cantilever Armrests</b>			
<i>(Adjustable length arm pad included (choose Pad setting below, except for Low Arm Option) NA with Manual Recline Back</i>			
<input type="checkbox"/>	195EG03	<b>Standard Cantilever Arm</b> w/ 2 Point Locking Post. (9"-10" Ht Range.) NA w/ seat depths below 16"	\$289.00
<input type="checkbox"/>	195EG04	<b>Standard Cantilever Arm</b> w/ 2 Point Locking Post. (11"-12" Ht Range.) NA w/ seat depths below 16"	\$289.00
<input type="checkbox"/>	195EG05	<b>STD Cantilever Arm</b> 9"-12" Ht Range	\$289.00
<input type="checkbox"/>	195EG06	<b>Low Cantilever Arm</b> 7"-10" Ht Range (Choose ASAPII pad type and size below)	\$289.00
<b>Single Post Armrests</b>			
<input type="checkbox"/>	195EG07	Single Post Height Adjustable - Std.	\$289.00
<input type="checkbox"/>	195EG08	Single Post Height Adjustable - Low	\$289.00

Armrest Mounting for K0856 SC Models with Tilt.			
<b>Cantilever Armrests</b>			
<i>(Adjustable length arm pad included (choose Pad setting below, except for Low Arm Option) NA with Manual Recline Back</i>			
<input type="checkbox"/>	195EG12	<b>STD Cantilever Arm</b> 9"-12" Ht Range	N/C
<input type="checkbox"/>	195EG13	<b>Low Cantilever Arm</b> 7"-10" Ht Range (Choose ASAPII pad type and size below)	N/C
<input type="checkbox"/>	195EG10	<b>Standard Cantilever Arm</b> w/ 2 Point Locking Post. (9"-10" Ht Range.) NA w/ seat depths below 16"	N/C
<input type="checkbox"/>	195EG11	<b>Standard Cantilever Arm</b> w/ 2 Point Locking Post. (11"-12" Ht Range.) NA w/ seat depths below 16"	N/C
<b>Single Post Armrests</b>			
<input type="checkbox"/>	195EG14	Single Post Height Adjustable - Std.	N/C
<input type="checkbox"/>	195EG15	Single Post Height Adjustable - Low	N/C
<b>Dual Post Armrests</b>			
<input type="checkbox"/>	195EG09	Shrt Dual Post Flip Back HT Adj 10"-15"	N/C
<input type="checkbox"/>	195EG23	Lng Dual Post Flip Back HT Adj 10"-15"	N/C

Arm Pad Options for ASAPII				
Pad Type / Size Options for Dual Post / Single Post / Low Cantilever				
<b>Left Arm Pad (N/C)</b>		<b>Right Arm Pad (N/C)</b>		
<input type="checkbox"/>	195AE01	Desk Length Waterfall (10")	<input type="checkbox"/>	195AB01
<input type="checkbox"/>	195AE02	Full Waterfall (14") NA w/Low Cant	<input type="checkbox"/>	195AB02
<input type="checkbox"/>	195AE03	Desk Length Standard (10")	<input type="checkbox"/>	195AB03
<input type="checkbox"/>	195AE04	Full Standard (14") NA w/Low Cant	<input type="checkbox"/>	195AB04
Pad Setting Options for <b>Standard Cantilever Arm</b>				
<b>Left Arm Pad (N/C)</b>		<b>Right Arm Pad (N/C)</b>		
<input type="checkbox"/>	195AE11	Full Length Cantilever (14")	<input type="checkbox"/>	195AB11
<input type="checkbox"/>	195AE10	Desk Length Cantilever (10")	<input type="checkbox"/>	195AB10

## 14. HANGERS AND FOOTRESTS FOR ASAP II

*Hanger and Extension Tubes Options (see footplate compatibility in the brackets)*

### Swing-away Hangers

#### Non-Powered Swing Away Legrests

*Sold as Pairs*

<input type="checkbox"/>	195H2	70° Swing Away (B,C)	N/C
<input type="checkbox"/>	195H106	80° Swing In-Out (B,C)	N/C
<input type="checkbox"/>	195H107	90° Swing In-out (B,C)	N/C
<input type="checkbox"/>	195H138	Manual ELR S/A (B,C) (Avail. with Fixed Centermount, choose below)	\$352.00

### Centermount Legrests

#### Center Mount Legrests

<input type="checkbox"/>	195H134	Fixed Center Mount Short w/ footplt	N/C
<input type="checkbox"/>	195H135	Fixed Center Mount Long w/footplt	N/C
<input type="checkbox"/>	195H103	Angle adjustable Center mount Short (I) <i>Independent leg length and angle adjustable</i>	\$381.00
<input type="checkbox"/>	195H104	Angle adjustable Center mount Medium (I) <i>Independent leg length and angle adjustable</i>	\$381.00

#### Swing Away Footplates

<input type="checkbox"/>	<b>B - Adult Angle Adjust w/Heel Loop for use with 65° and 70° hanger (Pair)</b>		<b>\$205.00</b>
	195F282	Adult Angle Adjustable Footplate	\$205.00
	195HL10	Heel Loop	N/C
<input type="checkbox"/>	<b>C - Kids Angle Adjust w/Leg Strap for use with 70° hanger (Pair)</b>		<b>\$205.00</b>
	195F283	Kids Angle Adjustable Footplate	\$205.00
	195LS13	Leg Strap	N/C
<input type="checkbox"/>	<b>C - Kids Angle Adjust w/Leg Strap for use with 80° and 90° hanger(Pair)</b>		<b>\$205.00</b>
	195F286	Kids Angle Adjustable Footplate	\$205.00
	195LS13	Leg Strap	N/C
<input type="checkbox"/>	195F279	Omit Footplate and extension tube	N/C

#### Centermount Footplates

<input type="checkbox"/>	<b>F - Centermount Dual Footplate (Pair) (Individually Height (3"-19") and Angle Adjustable)</b>		<b>\$247.00</b>
	195F315	CMT Angle Adjustable Footplates	\$247.00
<input type="checkbox"/>	<b>I - Angle adj Centermount Footplate pair</b>		<b>\$206.00</b>
	195F125	Angle Adjustable Footplates 5.75" x 6"	\$206.00

#### Omit Hangers.

<input type="checkbox"/>	195H144	Omit Hangers	N/C
--------------------------	---------	--------------	-----

#### Hanger Extensions. (See range listed in matrix below)

*NOTE: Not all extension tube will be available; depending on the seat height, some extension tube options may not offer enough ground clearance.*

<input type="checkbox"/>	195E16	High Mount Hanger	\$84.00
<input type="checkbox"/>	195E1	Short	N/C
<input type="checkbox"/>	195E2	Medium	N/C
<input type="checkbox"/>	195E5	Long	N/C
<input type="checkbox"/>	195E23	2" Extension for swing in/out hanger	N/C
<input type="checkbox"/>	195E24	4" Extension for swing in/out hanger	N/C
<input type="checkbox"/>	195E21	Omit Extension	N/C

#### Footrest Accessories

<input type="checkbox"/>	1955A01	Toe Loop (Pair) <i>Only avail on (opt B).</i>	\$37.00
<input type="checkbox"/>	195LS14	Extra Leg Strap	\$53.00
<input type="checkbox"/>	1952A01	Shoe holders Small Pair w/ padded strap	\$100.00
<input type="checkbox"/>	1952A02	Shoe holders Medium Pair w/ padded strap	\$100.00

**Lower Leg Length Matrix**  
(measured seat pan to the top of the footplate surface)

Hanger Type	High Mount	Short	Medium	Long
70° Swing Away	7-14"	13.5- 16"	16-18.5"	18-20.5"
80° Swing in/out	7-14"	13.5- 16"	14.5-19"	15-18"
90° Swing in/out	5.5-9" Footplate mounts to front of hanger.			
Manual ELR S/A		12.5-17"	14.5-19"	17-21"
Angle adj Cntmt		3"-7"	3"-13"	
Fixed Cntmnt		12.5"-15.5"		14.5"-19.5"

**Footplate Matrix**

Hanger Type	B	C	I
70° Swing Away	X	X	
80° Swing in/out	X	X	
90° Swing in/out	X	X	
Manual ELR S/A	X	X	
Ang Adj Cntmt			X
Fixed Cntmnt			

## 15. ACCESSORIES

Wheelchair Tray Table			
<input type="checkbox"/>	1957A01	TRAY TABLE SM TUFFAK 12"-14" PWR	\$331.00
<input type="checkbox"/>	1957A02	TRAY TABLE MED TUFFAK 15"-17" PWR	\$331.00

Backpack			
<input type="checkbox"/>	195A1	Black Utility	\$97.00
<input type="checkbox"/>	195A2	Black/Charcoal Kids	\$90.00



Accessories			
<input type="checkbox"/>	1958A09	Ram Phone Mount on Round JS tube	\$67.00
<input type="checkbox"/>	1958A10	Ram Tablet Mount on Round JS tube	\$268.00
<input type="checkbox"/>	1958A11	Ram Cup holder	\$67.00
<input type="checkbox"/>	1958A12	Switch-it Remote Stop Switch	\$1,200.00

Mounting Accessories (available with All Seats)			
<i>Only one accessory can be mounted at a time.</i>			
<input type="checkbox"/>	195A603	Cane/Crutch Holder	\$100.00
<input type="checkbox"/>	195A604	O2 holder	\$100.00
<input type="checkbox"/>	195A605	Walker holder	\$100.00