



Quickie QM-710™ Quickship



July 2016

Account: _____
 Date: _____
 PO#: _____
 Buyer: _____
 Marked For: _____

Ship To: _____
 Name: _____
 Address: _____
 City, Province, Postal Code: _____
 Phone: _____

QUICKIE QM-710™

172PC67 Quickie QM-710™ \$7,349.00
 Adult: WP0000686 Base Package
 Kids: WP0000687 R-Net 120A Non-Expandable Drive Control
 300 lb 6 mph 4-pole Motor Package
 14" Caster Wheel Pneumatics
 8A Dual Mode (Compact) CTE
 172TO25 4-Point Transit Package **(\$362.00)**
 172JM11 Swing Away Mount **(\$377.00)**
 172ST56 Quickie A.S.A.P II **(\$806.00)**

Speed Package (4-pole 300lb)

172PF85 6 mph Standard package N/C

Shroud Color

172EZ02 High Gloss Black N/C

User Weight

172WG03 WEPF 251-300lb N/C

Pre-Set Seat Frame Angle (STFH increases 1" per increment change)

172FA10 0° N/C

Batteries (Gel Cells)

172AG04 MK Group 34 Top Gel Cell (Pair) w/ Straps \$830.00
 Recommended for Max Daily Range
 Omit

Battery Charger

172BJ19 8A Dual Mode (Compact) CTE N/C

Drive Control Options

R-Net Non-Expandable Drive Control (120A)
 172EX02 Harnessing for Expandable Controller N/C
 R-Net Joystick
 172J108 Joystick w/ Stereo Jacks and Colour Display N/C

Armrests & Armpads

172EG22 Dual Post Flip Back HT Adj 10"-15" \$378.00
 172AE02 Full Length Waterfall (14") - Left Arm Pad N/C
 172AB02 Full Length Waterfall (14") - Right Arm Pad N/C

Positioning Belt - Select for any style seat

172EJ01 Auto Style Buckle N/C

Seat Width

172R5 16" Seat width N/C
 172R6 17" Seat width N/C
 172R7 18" Seat width N/C
 172R8 WEPA 19" Seat width N/C
 172R9 WEPA 20" Seat width N/C

Seat Depth

172RF5 16" Seat depth N/C
 172RF6 17" Seat depth N/C
 172RF7 18" Seat depth N/C
 172RF8 WEPB 19" Seat depth N/C
 172RF24 WEPB 20" Seat depth N/C

Backrest Frame Type

172BT8 Angle Adjustable (ASAP II Seat only) N/C

Angle Adjustable Backrest Options (NA with Manual recline)

172OD8 Non-Folding Backrest N/C

Cane Angle

172BP40 Medium Straight Back Post N/C

Towel Bar Depth

172QT04 Deep (Extends 4" from back canes for deep backs) N/C

Back Type

172BU1 Upholstery N/C

Non-Powered Swing Away Legrests

172H2 70° Swing Away (A,B,C) \$263.00
 172F285 Adult Angle Adjustable Footplate \$137.00
 172HL10 Heel Loop \$43.00

Fixed Center Mount Footboard

172H134 Fixed Center Mount Short N/C

*Seat Height (Use STFH Matrix below)

(Measured in the front) - Obtain Quote to verify e st. STFH

		Position "0"	Requires Group 34 Batteries	
<input type="checkbox"/>	172HT44	-		N/C
<input type="checkbox"/>	172HT49	+5" STFH		N/C
<input checked="" type="checkbox"/>	172HT45	+1" STFH		N/C
<input type="checkbox"/>	172HT46	+2" STFH	N/C	
<input type="checkbox"/>	172HT47	+3" STFH	N/C	
<input type="checkbox"/>	172HT48	+4" STFH	N/C	

Seating Type and Power Option Combinations	STFH Matrix					
	Seat-to-Floor Height					
	0	+5"	+1"	+2"	+3"	+4"
BC ASAP II No power seating.*	16.3	16.8	17.3	18.3	19.3	20.3

STFH measured with batteries but no user. Actual ST FH may vary by .5"
 *Adding seat angle will increase STFH by 1" for each incremental change.
 Yellow Highlight indicates what is shipped on chair.