



# Zippie XPLORE - XPERIENCE2

December 2017



### ACCOUNT INFORMATION

Dealer Contact Name:	_____	Bill to (Name):	_____
Email Address:	_____	Bill to (Location):	_____
Phone Number:	_____	Ship to (Name):	_____
Fax Number:	_____	Address:	_____
User Name or Number:	_____	City:	_____
User Weight:	_____	Province:	_____
User Height:	_____	Postal Code:	_____
Seat Cushion Thickness:	_____	P.O. # or Quote:	_____
Knee to Heel Length:	Left: _____" Right: _____"		
Seat to Top of Scapula:	_____		

The order cannot be processed without the following:  
 - A user name or number, the client's measurements and weight  
 - Sunrise Medical is not responsible for configuration or size discrepancies resulting from customer errors on the order form.  
 - Order confirmations will be provided by Amysystems summarizing the items selected. Please review this carefully.  
 - Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

### STANDARD OPTIONS

Zippie Power Bases, 6-wheel independent wheel suspension, an additional seat suspension system, tie-down hardware for strap systems, offboard charger, 6" front casters on MWD or 8" on RWD, and standard auto style buckle seat positioning strap, 5.9mph motors, Weight Capacity = 300lbs

#### (MID-WHEEL CONFIGURATION) ZIPPIE XPERIENCE POWER WHEELCHAIR

- XPERPED2 QUICKIE XPERIENCE PED - 2 Power Wheelchair \$8,295.00

#### (HYBRID-WHEEL CONFIGURATION) ZIPPIE XPLORE POWER WHEELCHAIR

- XPLOREPED2 QUICKIE XPLORE PED - 2 Power Wheelchair \$9,003.00

### USER WEIGHT PACKAGES

- WCA3 300 lb Maximum user weight capacity STD

#### STANDARD FRAME COLORS

- |                          |                  |                      |     |
|--------------------------|------------------|----------------------|-----|
| <input type="checkbox"/> | FCBLK            | Black                | STD |
| <input type="checkbox"/> | FCSBM            | Matte Black          | STD |
| <input type="checkbox"/> | FCGMC            | Gun Metal (Charcoal) | STD |
| <input type="checkbox"/> | FCSVM            | Silver Metallic      | STD |
| <input type="checkbox"/> | FCCBU            | Candy Blue           | STD |
| <input type="checkbox"/> | FCCRE            | Candy Red            | STD |
| <input type="checkbox"/> | FCMPL            | Midnight Purple      | STD |
| <input type="checkbox"/> | FCRZB            | Razzleberry          | STD |
| <input type="checkbox"/> | FCMWT <b>New</b> | Matte White          | STD |

#### OPTIONNAL FRAME COLORS (May increase lead times)

- |                          |                 |  |                  |
|--------------------------|-----------------|--|------------------|
| <input type="checkbox"/> | FCTXG           | Toxic Green                                    | N/C              |
| <input type="checkbox"/> | FCBOR           | Burnt Orange Metallic                          | N/C              |
| <input type="checkbox"/> | FCBAO           | Black Opal Metallic                            | N/C              |
| <input type="checkbox"/> | FCOPM           | Sparkle Orange                                 | N/C              |
| <input type="checkbox"/> | FCYLW           | Yellow   | N/C              |
| <input type="checkbox"/> | FCPKP           | Pearl Pink                                     | N/C              |
| <input type="checkbox"/> | FCHDD <b>†</b>  | Desert Digital Camo Hydro Graphics             | <b>New</b> \$330 |
| <input type="checkbox"/> | FCHWL <b>†</b>  | Woodland Camo Hydro Graphics                   | <b>New</b> \$330 |
| <input type="checkbox"/> | CUSTOM <b>†</b> | Carbon Fiber Hydro Graphics (Matte Black base) | <b>New</b> \$330 |

Note 1 Includes matching accessories with the Hydro Graphics option. The matching accessories option is not required with these colors

### MATCHING ACCESSORIES TO THE FRAME COLOR

- FC\*\*\*-M Matching Accessories to the Frame Color (Wheels, front access door, headrest mount and back shell when selected) \$249

### TRANSPORT OPTION

- TTD4 <sup>2</sup> 4-Point wheelchair transport brackets N/C

Note 2 **4-Point wheelchair transport brackets have not yet received WC-19 crash-testing approval. Use these transport brackets only to secure an unoccupied wheelchair during transport.**

### MOTOR PACKAGE

- MP65 6.5MPH 4 - Pole High-Speed Motor Package STD
- MP85 8.5MPH 4 - Pole High-Speed Motor Package (For Rear-Wheel Configuration only - Not available with VR2 electronics ECV902A) N/C

### ASSISTIVE DRIVING PACKAGE

- INTELLIDR <sup>3,4</sup> Intelli-Drive \$1,644

Note 3 **Intelli-Drive technology provides increased tracking capabilities with enhanced traction and steering precision. This is ideal for those with marginal hand function or switch control.**

Note 4 **Requires the R-net PM 120 Amp Expandable Controller (ECR120)**

**BATTERY CHARGER**

<input type="checkbox"/>	BCOB	Offboard Battery Charger	N/C
<input type="checkbox"/>	BCLES8OPT	Lester 8 amp Offboard Battery Charger	\$366

**BATTERY TRAY AND FRAME SIZE (TYPE)**

<input type="checkbox"/>	FS22B	5	Battery Tray for Group 22NF Batteries (Total base width = 23 3/8")	STD
<input type="checkbox"/>	FS2424B	6	Custom Battery Tray for Group 24 Batteries with a Compact 24" frame width ( <b>Mid-Wheel Drive configuration only</b> )	N/C

Please Select One Option for the Rear Hybrid Configuration: Front Caster Arms  Flared out (STD)  Straight (**Optional**)<sup>7</sup>

**Note 5** The 23 3/8" compact frame is not compatible with most of the Center Mount Front Rigging Options due to Mechanical conflicts

**Note 6** The 90° Power Center Mount Elevating & Articulating Foot Platform (ADDPF) is not available with this option

**Note 7** The Straight Front Arm option is not available with any of the Center Mount Front Rigging Options due to Mechanical conflicts

**BATTERIES**

<input type="checkbox"/>	BAT22NFPR	Gel Batteries - Group 22NF (Each - On chair) ( <b>E2361</b> )	\$319.00
<input type="checkbox"/>	BAT24PR	Gel Batteries - Group 24 (Each - On chair) ( <b>E2363</b> )	\$428.00

**DRIVE WHEEL**

<input type="checkbox"/>	DW14PT	14" Cast Wheel Pneumatic Tire	STD
<input type="checkbox"/>	DW14PI	14" Cast Wheel Pneumatic with Flat-Free inserts	\$249.00
<input type="checkbox"/>	DW144PT	8 14" x 4" Extra Wide Cast Wheel Pneumatic Tire	\$271.00
<input type="checkbox"/>	DW144PI	8 14" x 4" Extra Wide Cast Wheel Pneumatic with Flat-Free inserts	\$271.00
<input type="checkbox"/>	DW12PT	9 12 1/2" Cast Wheel Pneumatic Tire	N/C
<input type="checkbox"/>	DW12PI	9 12 1/2" Cast Wheel Pneumatic with Flat-Free inserts	\$249.00

The 14" x 4" Extra Wide Wheel increases the overall width of the Base by 1.5" to 2"

**Note 8** The option to Match the Wheels to the Shroud Color (FC...-M) is not available with the 14" x 4" Extra Wide wheels

The 14" x 4" option is not available with VR2 electronics (ECV902A)

**Note 9** 12 1/2" Drive wheels reduce overall seat to floor height by 3/4", lowers ground clearance and speed, but increases motor torque

**FRONT CASTERS (Mid-Wheel Configuration)**

<input type="checkbox"/>	FC6S	6" Solid Casters	STD
<input type="checkbox"/>	FC8PTM	10 8" Pneumatic Casters	\$168.00
<input type="checkbox"/>	FC8PIM	10 8" Pneumatic Casters with Flat-Free inserts	\$249.00

**Note 10** The 8" casters are NOT compatible with the 12 1/2" Drive wheel configuration (DW12PT and DW12PI)

**FRONT CASTERS (Rear-Wheel Configuration)**

<input type="checkbox"/>	FC8PTR	8" Pneumatic Casters	STD
<input type="checkbox"/>	FC8PIR	8" Pneumatic Casters with Flat-Free inserts	\$249.00
<input type="checkbox"/>	FC9PT	11 9" Pneumatic Casters	\$168.00
<input type="checkbox"/>	FC9PI	11 9" Pneumatic Casters with Flat-Free inserts	\$303.00

9" front casters increases seat to floor height by 1/2".

**Note 11** 9" casters are NOT compatible with the 12 1/2" Drive wheel configuration (DW12PT and DW12PI)

9" casters are NOT compatible with many front rigging options. Please see details in the front rigging sections on page #6

**SEAT TO FLOOR HEIGHT (Approximate measurement to the top of seat pan ± 1/4")**

**Options for a Seat Only**

<input type="checkbox"/>	12,13,15	STFH155	15 1/2"	(Possible Seat angle of 0° only)
<input type="checkbox"/>	13,15	STFH1625	16 1/4"	(Possible Seat angle of 0° only)
<input type="checkbox"/>	13,15	STFH1725	17 1/4"	(Possible Seat angle of 0° and 3°)
<input type="checkbox"/>	13	STFH1825	18 1/4"	(Possible Seat angle of 0° and 3°)
<input type="checkbox"/>	13	STFH1925	19 1/4"	(Possible Seat angle of 0°, 3° and 6°)
<input type="checkbox"/>	13	STFH2025	20 1/4"	(Possible Seat angle of 0°, 3°, 6°, 9°)

**Options for a Power Tilt or a STD Power Elevating Seat**

<input type="checkbox"/>	12,14,15	STFH155	15 1/2"	(Possible Seat angle of 0° only)
<input type="checkbox"/>	14,15	STFH1625	16 1/4"	(Possible Seat angle of 0° only)
<input type="checkbox"/>	14	STFH1725	17 1/4"	(Possible Seat angle of 0° only)
<input type="checkbox"/>	14	STFH1775	17 3/4"	(Possible Seat angle of 0° only)
<input type="checkbox"/>	14	STFH1825	18 1/4"	(Possible Seat angle of 0° only)
<input type="checkbox"/>	14	STFH1875	18 3/4"	(Possible Seat angle of 0° only)

**Options for a Combo - Power Tilt with a STD Power Elevating Seat**

<input type="checkbox"/>	12,14	STFH175	17 1/2"	(Possible Seat angle of 0° only) <b>New</b>
<input type="checkbox"/>	14	STFH1825	18 1/4"	(Possible Seat angle of 0° only) <b>New</b>

**Note 12** Requires the 12 1/2" Drive wheels (DW12PT or DW12PI). This lowers ground clearance and speed, but increases motor torque.

**Note 13** Seat to Floor Heights are taken with the Mid-Wheel Drive Configurations. Please remove 1/4" for the Rear-Wheel Drive Height

**Note 14** Seat to Floor Heights are taken with the Mid-Wheel Drive Configurations. Please Add 1/4" for the Rear-Wheel Drive Height

**Note 15** The Adult Heavy Duty Swing-Away footrests (SOAAF) are not available with this seat to floor height on the RWD configuration only

**PRE-SET SEAT FRAME ANGLE**

<input type="checkbox"/>	SFE-0	16	0° Pre-set seat frame seat angle ( <b>Must be selected when using power options</b> )	STD
<input type="checkbox"/>	SFE-3	16	3° Pre-set seat frame seat angle (Only available with a seat configuration without power options)	N/C
<input type="checkbox"/>	SFE-6	16	6° Pre-set seat frame seat angle (Only available with a seat configuration without power options)	N/C
<input type="checkbox"/>	SFE-9	16	9° Pre-set seat frame seat angle (Only available with a seat configuration without power options)	N/C

**Note 16** Please refer to the special notes in the above seat to floor height section to confirm seat angle possibilities.

The Light Package is not compatible with the EPSEL option.

### LIGHTING PACKAGE

<input type="checkbox"/>	<u>EOLP</u>	17	LED Light Package (Includes front beams, side flashers, rear tail light with flashers & brake light)	\$1,098.00
<b>Note 17 Requires the upgraded Joystick with Light Operation or the OMNI (ECVJ6KL, ECELJL, ECRJCL, ECRJCTPL or EDIB)</b>				

### NON-EXPANDABLE ELECTRONICS - PG DRIVES VR2 CONTROLLER AND JOYSTICKS

<input type="checkbox"/>	<u>ECV902A</u>		VR2 90 Amp Drive Controller	STD
<input type="checkbox"/>	<u>ECVJ4K</u>	18	VR2 Standard 4-Key Joystick	STD
<input type="checkbox"/>	<u>ECVJ6K</u>	18	VR2 6-Key Joystick (Allows for a single power option to operate through this joystick)	N/C
<input type="checkbox"/>	<u>ECVJ9KL</u>	18	VR2 9-Key Joystick with Lighting Operation (Allows for a single power function to operate through this joystick)	\$162.00
<b>Note 18 Requires the VR2 90A Drive Control (ECV902A).</b>				
<b>Note 18 VR2 Electronics can only accommodate a single power seat function. The only choices are Power Tilt or Elevating seat</b>				

### EXPANDABLE ELECTRONICS - PG DRIVES RNET CONTROLLERS

<input type="checkbox"/>	<u>ECEL90</u>		R-net EL 90 Amp Expandable Controller	\$546.00
<input type="checkbox"/>	<u>ECR120</u>		R-net PM 120 Amp Expandable Controller	\$1,098.00
<input type="checkbox"/>	<u>ECRNH</u>	19	Harness for Expandable Controller	N/C
<b>Note 19 Must be ordered with the expandable controller ECEL90 or ECR120 to interface with Joysticks or OMNI screen</b>				

### EXPANDABLE ELECTRONICS - PG DRIVES RNET JOYSTICKS

<input type="checkbox"/>	<u>ECO</u>		Omit Joystick (For use with an Omni Display)	N/C
<input type="checkbox"/>	<u>ECELJ</u>	20	R-net (LED) Joystick	STD
<input type="checkbox"/>	<u>ECRJCL</u>	20	R-net Joystick with Color screen	\$162.00

#### JOYSTICKS WITH LIGHTING OPERATION

<input type="checkbox"/>	<u>ECELJL</u>	20	R-net (LED) Joystick with Lighting Operation	\$157.00
<input type="checkbox"/>	<u>ECRJCL</u>	20	R-net Joystick with Color screen and Lighting Operation	\$271.00
<input type="checkbox"/>	<u>ECRJCTPL</u>	20	R-net Joystick with Color screen w/ On & Off / Mode switch & speed pot and Lighting Operation	\$325.00

#### ADVANCED JOYSTICK

<input type="checkbox"/>	<u>ECRJ35 New</u>	20	R-net Advanced Joystick (3.5" Color screen). Lighting Operation, paddle switches, shortcut keys and stereo jacks	\$1,167
<b>Note 20 This option requires the expandable controller ECEL90 or ECR120</b>				

### JOYSTICK MOUNT

<input type="checkbox"/>	<u>JSMF</u>		Fixed Joystick Mount	N/C
<input type="checkbox"/>	<u>JSMSA</u>		Swing Away Joystick Mount - Height adjustable, Angle adjustable and Multi-Axis	\$400.00
<input type="checkbox"/> Set-up with standard installation to swing-away sideways <input type="checkbox"/> Set-up at 35° to swing-away sideways and downwards				
<input type="checkbox"/>	<u>JSMCM</u>		Midline Mounted Swing Away Joystick Mount	\$400.00

### JOYSTICK MOUNT SUPPORT

<input type="checkbox"/>	<u>JSMFH</u>		Joystick mount support Bracket	STD
--------------------------	--------------	--	--------------------------------	-----

### JOYSTICK HAND MOUNT

<input type="checkbox"/>	<u>HML</u>		Install Joystick on Left hand side	
<input type="checkbox"/>	<u>HMR</u>		Install Joystick on Right hand side	

### JOYSTICK HANDLE

<input type="checkbox"/>	<u>JSH</u>		Standard Joystick Gimble	STD
<input type="checkbox"/>	<u>JBAT</u>		Stick Handle 3" in Length	\$93.00
<input type="checkbox"/>	<u>JBH</u>		Ball Handle	\$98.00
<input type="checkbox"/>	<u>JFBH</u>		Foam Ball Handle	\$98.00
<input type="checkbox"/>	<u>JTH</u>		T-Handle	\$98.00
<input type="checkbox"/>	<u>JGP</u>		Goalpost	\$98.00

### ATTENDANT DRIVE CONTROL OPTION

<input type="checkbox"/>	<u>EADC</u>		Attendant Hand Control for R-net EL & PM Controls	\$1,098.00
--------------------------	-------------	--	---	------------

### ATTENDANT JOYSTICK HAND MOUNT

<input type="checkbox"/>	<u>EADCL</u>		Install Attendant Joystick on Rear Left hand side	
<input type="checkbox"/>	<u>EADCR</u>		Install Attendant Joystick on Rear Right hand side	

### ENHANCED DISPLAY & DIGITAL INTERFACE BOX (Included with specialty controls)

<input type="checkbox"/>	<u>EDIB</u>	21	Omni (Enhanced Display & Digital Interface Box with Infrared)	\$2,094.00
<b>Note 21 Allows the switch operation of alternative devices through driver controls. Requires one of the R-net fully Expandable Controllers (ECEL90 or ECR120)</b>				

### ENHANCED DISPLAY MOUNT

<input type="checkbox"/>	<u>EDSAM</u>		User Swing Away Mount	
<input type="checkbox"/>	<u>EDML</u>		(Install Display on LH side)	
<input type="checkbox"/>	<u>EDMR</u>		(Install Display on RH side)	N/C

### ENVIRONMENTAL CONTROL MODULE

<input type="checkbox"/>	ECU	22	Environment control module - (Output Module)	\$1,098.00
--------------------------	-----	----	--	------------

Note 22 Requires one of the fully expandable R-net controllers (ECEL90 or ECR120)

### BLUETOOTH MOUSE EMULATOR (Includes USB Dongle)

<input type="checkbox"/>	EDBT	23	Bluetooth Mouse Emulator (USB Dongle included)	\$1,098.00
	<input type="checkbox"/>		Windows and MAC OS	
	<input type="checkbox"/>		Android Devices	
	<input type="checkbox"/>		Iphone, Ipad and Ipod	

Note 23 Requires one of the fully expandable R-net controllers (ECEL90 or ECR120)  
Not limited to only one selection. Any two can be installed on the same chair, but only one can be used at a time.

### POWER OPTIONS & SEATING

<input type="checkbox"/>	SAABA		Amyseat Rehab Seat - Seat Frame only	\$855.00
<input type="checkbox"/>	SAPTX	24	CG Power Tilt	\$7,170.00
<input type="checkbox"/>	SAPRX	24, 25, 26	Power Recline with Integrated Shear Reduction (ISR)	\$8,825.00
<input type="checkbox"/>	SAPTRX	24, 25, 26	CG Power Tilt and Power Recline with Integrated Shear Reduction (ISR)	\$12,688.00
<input type="checkbox"/>	SAPESX	24	11" Power Elevating Seat Module for transfer assist	\$3,555.00
<input type="checkbox"/>	ADDPESX		<b>Module</b> - Add an 11" Power Elevating Seat to either SAPTX, SAPRX or SAPTRX	\$3,306.00
<input type="checkbox"/>	ADDMR	25, 26	<b>Module</b> - Manual Recline with Integrated Shear Reduction (ISR)	\$1,098.00

- Note 24 These options include the fully adjustable Amyseat (SAABA).  
 Note 25 This option is not available with the Pediatric Seat Frame (SOWN11-16)  
 Note 26 It is recommended to use the Contour Back Support option with the power recline (STDSB or STDSB/CUSTOM)

### ELECTRONICS - SINGLE OR DUAL (1 OR 2 ACTUATORS) POWER SEAT FUNCTION CONTROLLERS

<input type="checkbox"/>	EPSV	27, 28	Tilt or Elevating seat activation through <b>6-key VR2</b> Hand Control (ECVJ6K or ECVJ9KL)	\$9,003
<input type="checkbox"/>	EPSEL	27, 28	Tilt or Elevating seat activation through <b>RNET 90 Amp Expandable Controller (ECEL90)</b>	N/C
<input type="checkbox"/>	SAC		Single (1 Actuator) Power Seat Function Controller to Interface through R-net Electronics	N/C
<input type="checkbox"/>	DAC		Dual (2 Actuators) Power Seat Function Controller to Interface through R-net Electronics	\$546.00

**With the SAC or DAC, please select one activation and position from the list below**

<input type="checkbox"/>	SWTJSS	29, 30	Single actuator activation through R-net Joystick or Specialty Controls	\$1,098
<input type="checkbox"/>	SWTJSD	29, 31	Dual actuator activation through R-net Joystick or Specialty Controls	\$1,098
<input type="checkbox"/>	SWDUALACT		Dual Activation ( <b>Allows to operate all seat functions with external switches AND through the joystick</b> )	\$340
<input type="checkbox"/>	SWTOGGLE	30 or 31	Single Toggle Switch (Single Mode: DAC - Dual Mode: SAC) <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation	N/C
<input type="checkbox"/>	SWPUSH1	30	Single Push Button Switch (Single Mode: SAC) <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation	N/C
<input type="checkbox"/>	SWPUSH2	30 or 31	Dual Push Button Switch (Single Mode: DAC - Dual Mode: SAC) <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation	N/C

- Note 27 Allows for the activation of a single power option through the VR2 or R-net 90 Amp Controllers. Only the Power Tilt or Elevating seat options are available with this feature.  
 Note 28 These options do not require the SAC controller as they connect via the chair's power module directly.  
 Note 29 Requires the R-net 90 Amp (ECEL90) or the 120 Amp (ECR120) Fully Expandable Controllers  
 Note 30 Requires the Single (1 Actuator) Power Seat Function Controller (SAC)  
 Note 31 Requires the Dual (2 Actuators) Power Seat Function Controller (DAC)

### ELECTRONICS - MULTIPLE (3 to 6 ACTUATORS) POWER SEAT FUNCTION CONTROLLER

<input type="checkbox"/>	MAC		Multiple (3 to 6 Actuators) Power Seat Function Controller to Interface through R-net Electronics	\$2,641.00
--------------------------	-----	--	---	------------

**With the MAC, please select one activation and position from the list below**

<input type="checkbox"/>	SWTJSM	32, 33	Multiple Actuator activation through R-net Joystick or Specialty Controls	STD
<input type="checkbox"/>	SWDUALACT		Dual Activation ( <b>Allows to operate all seat functions with external switches AND through the joystick</b> )	\$340
<input type="checkbox"/>	SWMULTI	32	4-Way Toggle Switch (Up to 4 Activations) - <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation	N/C
<input type="checkbox"/>	SWMULTIT	32	Triple Toggle Switch (Up to 5 Activations) - <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation	N/C
<input type="checkbox"/>	SWPUSH4	32	4-Way Push Button Switch - <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation	N/C
<input type="checkbox"/>	SWPUSH5	32	5-Way Push Button Switch - <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation	N/C

- Note 32 Requires the Multiple (3 to 6 Actuators) Power Seat Function Controller (MAC)  
 Note 33 Requires the R-net 90 Amp (ECEL90) or 120 Amp (ECR120) Fully Expandable Controller, and the multiple Actuator Controller (MAC)

### SWITCH ACTIVATION HAND MOUNT

Switch Installation:	<input type="checkbox"/> TSML (Install Switch on LH side)	<input type="checkbox"/> TSMR (Install Switch on RH side)	N/C
----------------------	---	---	-----

**PRE-SET SEAT ANGLE (For Non-Reclining Backs)**

SOBAA      Please Select one:  86°  90°  94°  98°  102°  106°  110°  114°  118°  122°      N/C

**PEDIATRIC ADJUSTABLE SEAT WIDTH RANGE BETWEEN 11" AND 16"** STD

SOWN11-16 (Seat to Floor Height increases by 1 ¼")

<input type="checkbox"/> 11" (SOWN11)	<input type="checkbox"/> 14" (SOWN14)
<input type="checkbox"/> 12" (SOWN12)	<input type="checkbox"/> 15" (SOWN15)
<input type="checkbox"/> 13" (SOWN13)	<input type="checkbox"/> 16" (SOWN16)

**PEDIATRIC ADJUSTABLE SEAT DEPTH RANGE BETWEEN 11" AND 16"** STD

<input type="checkbox"/> 11" (SODN11)	<input type="checkbox"/> 14" (SODN14)
<input type="checkbox"/> 12" (SODN12)	<input type="checkbox"/> 15" (SODN15)
<input type="checkbox"/> 13" (SODN13)	<input type="checkbox"/> 16" (SODN16)

**STANDARD ADJUSTABLE SEAT WIDTH RANGE BETWEEN 15" AND 20"** N/C

SOWN15-20-19 (Width adjustment from 15" to 20" with a maximum Depth of 19" )

<input type="checkbox"/> 15" (SOWN15)	<input type="checkbox"/> 17" (SOWN17)	<input type="checkbox"/> 19" (SOWN19S)
<input type="checkbox"/> 16" (SOWN16)	<input type="checkbox"/> 18" (SOWN18)	<input type="checkbox"/> 20" (SOWN20S)

**ADJUSTABLE SEAT DEPTH RANGE BETWEEN 15" AND 19"** N/C

15" (SODN15)       16" (SODN16)       17" (SODN17)       18" (SODN18)       19" (SODN19)

**BACK CANE HEIGHT** (The back cane height measurement is taken from the seat pan to the top of the push handles)

<input type="checkbox"/> 16" (SOBCH16)	<input type="checkbox"/> 20" (SOBCH20)	<input type="checkbox"/> 24" (SOBCH24)
<input type="checkbox"/> 18" (SOBCH18)	<input type="checkbox"/> 22" (SOBCH22) <sup>34</sup>	<input type="checkbox"/> 26" (SOBCH26)

**Note 34** 22" is the minimum back cane height required with a power or manual recline.

**BACK SUPPORT**

<input type="checkbox"/> STDSB	<sup>35</sup> Contour back support with a maximum cushion height of 20" (11" to 20" Wide)	\$606
<input type="checkbox"/> STDSB/CUSTOM	Contour back support (Non-Standard Sizes) Specify:      Width:      Height:	\$1,294

**Note 35** It is recommended to use the Contour Back Support option with the power or manual recline.

**37 BACK SUPPORT TOTAL HEIGHT (MEASURED FROM THE SEAT PAN TO THE TOP OF THE BACK)**

<input type="checkbox"/> 16" (SOBH16)	<input type="checkbox"/> 20" (SOBH20)	<input type="checkbox"/> 24" (SOBH24)	<input type="checkbox"/> 28" (SOBH28)
<input type="checkbox"/> 18" (SOBH18)	<input type="checkbox"/> 22" (SOBH22) <sup>36</sup>	<input type="checkbox"/> 26" (SOBH26)	

**Note 36** With a power or manual recline, the minimum height is 22".

**Note 37** The space between the bottom of the back and the seat pan is 4" with a non-reclining back and 6" with a reclining back

**SEAT CUSHION**

<input type="checkbox"/> SCPSC	Contour Seat Cushion with a maximum seat depth of 19" (11" to 20" Wide)	\$606
<input type="checkbox"/> SCPSC/CUSTOM	Contour Seat Cushion (Non-Standard Sizes) Specify:      Width:      Depth:	\$1,320

**ARMRESTS**

**CANTILEVER - HEIGHT ADJUSTABLE ARMRESTS - OPENS WITH A 5° ANGLE (For Non-Reclining Backs - Choose 2)**

<input type="checkbox"/> SOAFBDLH	Pediatric Flip-Back Cantilever - Left Side (Height Adjustable from 7" to 15" in 1" increments)	STD
<input type="checkbox"/> SOAFBDRH	Pediatric Flip-Back Cantilever - Right Side (Height Adjustable from 7" to 15" in 1" increments)	STD
<input type="checkbox"/> SOAFBFLH	Adult Flip-Back Cantilever - Full Length Left Side (Height Adjustable from 7" to 15" in 1" increments)	N/C
<input type="checkbox"/> SOAFBFRH	Adult Flip-Back Cantilever - Full Length Right Side (Height Adjustable from 7" to 15" in 1" increments)	N/C
<input type="checkbox"/> SOAFBDLH	Adult Flip-Back Cantilever - Desk Length Left Side (Height Adjustable from 7" to 15" in 1" increments)	N/C
<input type="checkbox"/> SOAFBDRH	Adult Flip-Back Cantilever - Desk Length Right Side (Height Adjustable from 7" to 15" in 1" increments)	N/C

**DUAL POST - HEIGHT ADJUSTABLE ARMRESTS - INCLUDES TRANSFER GRIPS (For Non-Reclining Backs - Choose 2)**

<input type="checkbox"/> SOADPPLH	Pediatric Adjustable Height Dual Post - Left Side (Height Adjustable from 7" - 10")	\$195.00
<input type="checkbox"/> SOADPPRH	Pediatric Adjustable Height Dual Post - Right Side (Height Adjustable from 7" - 10")	\$195.00
<input type="checkbox"/> SOADPFLH	Dual Post Height Adjustable - Full Length Left Side (Height Adjustable from 9 ½" - 13 ½")	\$195.00
<input type="checkbox"/> SOADPFRH	Dual Post Height Adjustable - Full Length Right Side (Height Adjustable from 9 ½" - 13 ½")	\$195.00
<input type="checkbox"/> SOADPDLH	Dual Post Height Adjustable - Desk Length Left Side (Height Adjustable from 9 ½" - 13 ½")	\$195.00
<input type="checkbox"/> SOADPDRH	Dual Post Height Adjustable - Desk Length Right Side (Height Adjustable from 9 ½" - 13 ½")	\$195.00

**RECLINING FLIP-UP, HEIGHT ADJUSTABLE ARMRESTS (For Reclining Backs - Choose 2)**

<input type="checkbox"/> SOANIFLH	Adult Reclining armrest - Full Length Left Side (Height Adjustable from 10 ½" to 14" in ½" increments)	N/C
<input type="checkbox"/> SOANIFRH	Adult Reclining armrest - Full Length Right Side (Height Adjustable from 10 ½" to 14" in ½" increments)	N/C
<input type="checkbox"/> SOANIDLH	Adult Reclining armrest - Desk Length Left Side (Height Adjustable from 10 ½" to 14" in ½" increments)	N/C
<input type="checkbox"/> SOANIDRH	Adult Reclining armrest - Desk Length Right Side (Height Adjustable from 10 ½" to 14" in ½" increments)	N/C

## ARMREST HEIGHT ADJUSTMENT AT FACTORY

Please set the Armrest Height at: \_\_\_\_\_ " (The Armrest Height is taken from the seat pan to the top of the armpad. Not required for Dual Post)

### ARMPADS

#### Foam Type Selection

<input type="checkbox"/>	SOFAPD	Desk Length Standard (10")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	STD
<input type="checkbox"/>	SOFAPF	Full Length Standard (14")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	STD
<input type="checkbox"/>	SOWAPD	39 Desk Length Waterfall (10")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$33.00
<input type="checkbox"/>	SOWAPF	39 Full Length Waterfall (14")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$33.00

#### Hybrid Gel and Foam Type Selection (The Gel always stays soft)

<input type="checkbox"/>	SOGELD <span style="color: red;">new</span>	Desk Length Hybrid Gel and Foam Flat armpad (10" x 2.5")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$33
<input type="checkbox"/>	SOGELF	Full Length Hybrid Gel and Foam Flat armpad (14" x 2.5")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$33
<input type="checkbox"/>	SOWGELD <span style="color: red;">new</span>	39 Desk Length Hybrid Gel and Foam Waterfall armpad (10")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$50
<input type="checkbox"/>	SOWGELF <span style="color: red;">new</span>	39 Full Length Hybrid Gel and Foam Waterfall armpad (14")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$50

#### Specialty Type Selection

<input type="checkbox"/>	SOCAR	40 Arm Trough with Hand Support (18" x 4 1/2")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$216.00
--------------------------	-------	--	---	--------------	---	--	----------

**Note 39** The Waterfall armpads may interfere with flipping back the armrests when combined with Lateral Trunk supports

**Note 40** Arm Trough Armpads may interfere with Armrests Flipping back.

### FOOTRESTS

<input type="checkbox"/>	SOPA AF	46 Pediatric Swing-Away Footrests -----	Set at:	<input type="checkbox"/> 60°	<sup>47</sup> <input type="checkbox"/> 70°	<sup>41,47</sup> <input type="checkbox"/> 80°	<sup>41,47</sup> <input type="checkbox"/> 90°	\$324.00
		Select Footrest Range:	<input type="checkbox"/> SOFLP5 (5 1/2" - 7 1/2")	<input type="checkbox"/> SOFLP8 (8 1/2" - 10 1/2")	<input type="checkbox"/> SOFLP10 (10 1/2" - 15 1/2")			
<input type="checkbox"/>	SOSA AF	45, 46 Adult Swing-Away Footrests -----	Set at:	<input type="checkbox"/> 60°	<sup>47</sup> <input type="checkbox"/> 70°			\$271.00
		Select Footrest Range:	<input type="checkbox"/> SOFLA10 (10 1/2" - 15 1/2")	<input type="checkbox"/> SOFLA14 (14 1/2" - 19 1/2")				
<input type="checkbox"/>	SOAAF	45, 46 Adult <b>Heavy-Duty</b> Swing-Away footrests -----	Set at:	<input type="checkbox"/> 60°	<sup>47</sup> <input type="checkbox"/> 70°			\$379.00
		Select Footrest Range:	<input type="checkbox"/> SOFLHD12 (12 1/2" - 15 1/2")	<input type="checkbox"/> SOFLHD14 (14 1/2" - 19 1/2")				
<input type="checkbox"/>	ADDPL	42, 45, 46 Adult Power Elevating & Articulating Legrests	Please Select Function:	<input type="checkbox"/> Individual <sup>43</sup>	<input type="checkbox"/> Combined <sup>44</sup>		\$9,003.00	
		Select Footrest Range:	<input type="checkbox"/> SOFLPW12 (12 1/2" - 15 1/2")	<input type="checkbox"/> SOFLPW14 (14 1/2" - 19 1/2")				

<sup>46</sup> Please set the Footrest Length at: \_\_\_\_\_ " (The footrest length is taken from the seat pan to the heel)

**Note 41** Depending on the seat configuration, this option may not be available.

**Note 42** This option includes calf pads

**Note 43** Individual legs count as two power functions.

**Note 44** Combined legs count as one power function.

**Note 45** This option is not available with the pediatric seat frame (SOWN11-16)

**Note 46** The Footrest Length is taken from the top of the seat pan to the top of the footplate.

This footrest angle is not compatible with the 9" front casters (FC9PT and FC9PI)

**Note 47**

This footrest angle may interfere with the casters when selecting the straight front caster arms option

### FOOTPLATES & ACCESSORIES

<input type="checkbox"/>	SOSFAF	Flip-up Footplates	STD
<input type="checkbox"/>	SOAFP	Angle Adjustable & Flip-up Footplates - Each	\$141.00
<input type="checkbox"/>	SOAHL	Heel Loops - Each	\$43.00
<input type="checkbox"/>	SOACPFR	Calf Pads for Adult Heavy-Duty Swing-Away footrests ( <b>SOAAF</b> ) with adjustable hardware - Each	\$141.00

### CENTER MOUNTS

<input type="checkbox"/>	SOAAF90	<sup>50, 51, 52, 55, 56</sup> 90° Fixed Center Mount Interface Bracket ( <b>Offers Manual adjustability from 95° to 70°</b> )	\$249.00			
		Select Foot Platform Range:	<input type="checkbox"/> SOFBP7 (8 1/2" - 12 1/2")	<input type="checkbox"/> SOFBA12 (12 1/2" - 16 1/2")	<input type="checkbox"/> SOFBA16 (16 1/2" - 19 1/2")	
<input type="checkbox"/>	ADDPF	<sup>48, 49, 50, 52, 53, 54, 55</sup> 90° Power Center Mount Elevating & Articulating Foot Platform - *** <b>Not available with the 12 1/2" Drive Wheel</b>	\$3,084.00			
		<input type="checkbox"/> SOFBPW12 (12 1/2" - 16 1/2")	<input type="checkbox"/> SOFBPW16 (16 1/2" - 20 1/2")			

<sup>55</sup> Please set the Footrest Length at: \_\_\_\_\_ " (The footrest length is taken from the seat pan to the heel)

**Note 48** The power center mount counts as one power function.

**Note 49** This option includes calf pads

**Note 50** Links and attaches the one piece center mount foot platform (SOFUFP) or the split footplates (SOAFPFB) to the seat

**Note 51** With the split footplates option (SOAFPFB), this bracket offers independent height adjustability for each footplate

**Note 52** Please Select the one piece center mount foot platform (SOFUFP) or the split footplates (SOAFPFB)

**Note 53** This option is not available with the pediatric seat frame (SOWN11-16)

**Note 54** The power center mount option is not available with the straight front caster arms (RWD) and Group 24/24 Frame (FS2424B - MWD)

**Note 55** The Foot Platform Length is taken from the top of the seat pan to the top of the foot platform or split foot plates.

**Note 56** This option is not compatible with the 9" Caster option on Rear Hybrid Configuration



### FOOT PLATFORMS & ACCESSORIES

<input type="checkbox"/>	SOFUFP	Flip-up one piece Foot Platform	N/C
<input type="checkbox"/>	SOAFPFB	Angle Adjustable & Flip-up Split Footplates for center mount - Each	\$141.00
<input type="checkbox"/>	SOACPCM 57	Calf Pads for 90° Fixed Center Mount ( <b>SOAAF90</b> ) with adjustable hardware - Each	\$141.00

Note 57 Not available with the pediatric length (SOFBP7 - 7 ½" - 12 ½")

### SEAT POSITIONING

<input type="checkbox"/>	SOAUTO	55" Seat Positioning Strap	STD
--------------------------	--------	----------------------------	-----

### HEADREST

<input type="checkbox"/>	SONUAHR	Using headrest mounting hardware other than SOMAH located bellow. Specify: _____	
<input type="checkbox"/>	SOFHRS	Small Curved Headrest Pad 4" x 9"	\$220.00
<input type="checkbox"/>	SOFHRM	Medium Curved Headrest Pad 5" x 10"	\$220.00
<input type="checkbox"/>	SOFHRL	Large Curved Headrest Pad 6" x 11"	\$220.00
<input type="checkbox"/>	SOMAH	Multi-axis Removable Mounting Hardware	\$178.00
<input type="checkbox"/>	SOHRBC	Bracket to Mount on the Contour Back Support ( <b>Without Recline</b> )	\$108
<input type="checkbox"/>	SOHRPR	Bracket to Mount on the Seat Frame ( <b>For Reclining Backs</b> )	\$108

\* Please select one of the headrest options with the mounting hardware and specify where the headrest will be installed

### LATERAL TRUNK SUPPORTS

<input type="checkbox"/>	SOSLS3	Adjustable Laterals with Small Curved Pads (Each) - 3" x 5"	\$655
<input type="checkbox"/>	SOSLM4	Adjustable Laterals with Medium Curved Pads (Each) - 4" x 6"	\$9,003
<input type="checkbox"/>	SOSLL5	Adjustable Laterals with Large Curved Pads (Each) - 5" x 7"	\$655
<input type="checkbox"/>	SOSLHM	Adjustable Swing-Away Hardware Mounts - Each	\$325

\*\* Please select one of the lateral pads and the adjustable & removable hardware.

\$9,003

### HIP & THIGH SUPPORTS

<input type="checkbox"/>	SOHTSS10	Adjustable Hip & Thigh Support Pads (Each) - Small 10" x 4"	\$249
<input type="checkbox"/>	SOHTSM12	Adjustable Hip & Thigh Support Pads (Each) - Medium 12" x 5"	\$249
<input type="checkbox"/>	SOHTSL14	Adjustable Hip & Thigh Support Pads (Each) - Large 14" x 6"	\$249
<input type="checkbox"/>	SOHTSARH	Adjustable and Removable Hip & Thigh Support Hardware - Each	\$471

\*\*\* Please select one of the hip & thigh support pads and the adjustable & removable hardware.

### ACCESSORIES

<input type="checkbox"/>	SOCH	Cup Holder (Self-Leveling)	\$98
<input type="checkbox"/>	SOCPH	Ram Mounts X-Grip Cell Phone Holder	\$168
<input type="checkbox"/>	SOTHS	Ram Mounts X-Grip Small Tablet Holder	\$190
<input type="checkbox"/>	SOTHL	Ram Mounts X-Grip Large Tablet Holder with Swivel	\$286

### REFLECTOR KIT

<input type="checkbox"/>	REFLECT	Red & White Reflector Kit Stickers (They will not be installed - Included in a bag with the owner's manual)	N/C
--------------------------	---------	---	-----

### SWITCHES AND OPTIONS

<input type="checkbox"/>	ASCPCP	PGDT RNET Dongle and PC programmer kit	\$633.00
<input type="checkbox"/>	ASCPVR2	PGDT VR2 Dongle and PC programmer kit	\$633.00
<input type="checkbox"/>	SSBS	Buddy Button Switch	\$119.00
<input type="checkbox"/>	SSES	Egg Switch-Black	\$119.00
<input type="checkbox"/>	SSML	Micro Light Switch	\$162.00

## SPECIALTY DRIVER CONTROL OPTIONS

### Therafin Sip-N-Puff Kit

<input type="checkbox"/>	SCSNPK	58, 59	Therafin Sip-N-Puff Kit		\$1,651
<input type="checkbox"/>	SCSNPBTK	59	Therafin Sip-N-Puff Breath Tube Kit		\$536
<input type="checkbox"/>	SCSRMSNP	59	Stealth Rock-n-Lock Swingaway Arm for Sip-N-Puff	<input type="checkbox"/> SWL (Swings on Left) <input type="checkbox"/> SWR (Swings on Right)	\$298
<input type="checkbox"/>	SCSRM	59	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/> SWL (Swings Left) <input type="checkbox"/> SWR (Swings Right)	\$487
<input type="checkbox"/>	SSES-STD	59, 61	Egg Switch for Changing Modes - Black	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$119
<input type="checkbox"/>	SSBS-STD	59	Buddy Button Switch (Stop Switch) - Red		\$119
<input type="checkbox"/>	SOMAH	59	Multi-axis Headrest Removable Mounting Hardware		\$178
<input type="checkbox"/>	SOFHRS		Small Curved Headrest Pad 4" x 9"		\$222
<input type="checkbox"/>	SOFHRM	60	Medium Curved Headrest Pad 5" x 10"		\$222
<input type="checkbox"/>	SOFHRL		Large Curved Headrest Pad 6" x 11"		\$222

### Switch-It Proximity Head Array

<input type="checkbox"/>	SCSIHAP	58, 59	Pediatric - 3 Switch Proximity Head Array with 1/8" Mono Jack for mode/reverse		\$5,513
<input type="checkbox"/>	SCSRM	59	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/> SWL (Swings Left) <input type="checkbox"/> SWR (Swings Right)	\$487
<input type="checkbox"/>	SSES-STD	59, 61	Egg Switch for Changing Modes - Black	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$119
<input type="checkbox"/>	SSBS-STD	59	Buddy Button Switch (Stop Switch) - Red		\$119
<input type="checkbox"/>	SOMAH	59	Multi-axis Headrest Removable Mounting Hardware		\$178

### Switch-It MicroPilot Light Touch Proportional Mini joystick for chin / lip control

<input type="checkbox"/>	SCSIMP	58, 59	MicroPilot Light touch mini Proportional Joystick		\$4,323
<input type="checkbox"/>	SCSRMMP	59	Stealth Rock-n-Lock Swingaway Mount for Chin M-P (E1028)	<input type="checkbox"/> SWL (Swings Left) <input type="checkbox"/> SWR (Swings Right)	\$536
<input type="checkbox"/>	SCSRM	59	Stealth Rock-n-Lock Swingaway Mount for Switch	<input type="checkbox"/> SWL (Swings Left) <input type="checkbox"/> SWR (Swings Right)	\$487
<input type="checkbox"/>	SSES	59, 61	Egg Switch-Black	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$119
<input type="checkbox"/>	JCC		Chin Cup		\$9,003
<input type="checkbox"/>	SOMAH	59	Multi-axis Headrest Removable Mounting Hardware		\$178
<input type="checkbox"/>	SOFHRS		Small Curved Headrest Pad 4" x 9"		\$222
<input type="checkbox"/>	SOFHRM	60	Medium Curved Headrest Pad 5" x 10"		\$222
<input type="checkbox"/>	SOFHRL		Large Curved Headrest Pad 6" x 11"		\$9,003

### Proportional Compact or Mini joystick for Midline Mount

<input type="checkbox"/>	SCPGMJS	58, 59	PG Drives JC200 Proportional Compact Remote joystick for Midline Mount		\$1,028
<input type="checkbox"/>	SCGMJSM	62	Gatlin Midline Swing-away Mount for JC200	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$422
<input type="checkbox"/>	SCSIMP	58, 59	Switch-It MicroPilot Light touch mini Proportional Joystick		\$4,323
<input type="checkbox"/>	SCGMPPM	63	Gatlin Midline Swing-away Mount for MicroPilot	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$460
<input type="checkbox"/>	SCBTM	64	Bullet Hand Tray for MicroPilot	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$487

### Proportional Compact or Mini joystick for Armrest Mount

<input type="checkbox"/>	SCPGMJS	58, 59	PG Drives JC200 Proportional Compact Remote joystick for Armrest Mount		\$1,028
<input type="checkbox"/>	SCARMJSM	62	Swing-away Armrest Mount for JC200	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$422
<input type="checkbox"/>	SCSIMP	58, 59	Switch-It MicroPilot Light touch mini Proportional Joystick		\$4,323
<input type="checkbox"/>	SCARMPPM	63	Swing-away and Multi-Axis armrest mount	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$460
<input type="checkbox"/>	SCBTM	64	Bullet Hand Tray for MicroPilot	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$487

### Proportional RIM Head Control Package

<input type="checkbox"/>	SCPHC	58, 59	Proportional RIM Head Control Package		\$2,700
<input type="checkbox"/>	SCSRM	59	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/> SWL (Swings Left) <input type="checkbox"/> SWR (Swings Right)	\$487
<input type="checkbox"/>	SSES-STD	59, 61	Egg Switch for Changing Modes - Black	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$119
<input type="checkbox"/>	SOMAH	59	Multi-axis Headrest Removable Mounting Hardware		\$178

**Note 58** Requires the Omni Enhanced Display & Digital Interface Box (EDIB)

**Note 59** This option must be selected for the Kit

**Note 60** Choose one (1) Headrest Pad

**Note 61** The Mode Switch Type can be changed from this selection on Page #7

**Note 62** Must be selected for use with the PG Drives JC200 Mini Joystick

**Note 63** Must be selected for use with the Switch-It MicroPilot

**Note 64** This option is only compatible with the Switch-It MicroPilot



**65 CUSTOM SECTION**

**Note 65** Lead times will increase depending on the custom order pipeline at the time the purchase order is received  
Please call for Quote and specifications for these options.

**O<sup>2</sup> Holder**

<input type="checkbox"/>	SUPOXY	66	Custom Non-Articulating Oxygen Tank Holder	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	\$422
<input type="checkbox"/>	SUPOXYA	67	Custom Articulating Oxygen Tank Holder	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	\$1,981
			<input type="checkbox"/> D-Size Oxygen Tank (* 22" Tall x 4.3" Diameter) * The height includes the valve			
			<input type="checkbox"/> E-Size Oxygen Tank (* 31" Tall x 4.3" Diameter) * The height includes the valve			

**Note 66** Not available with a power tilt and/or power recline

**Note 67** Required with a power tilt and/or power recline

**VENT TRAY (Please complete in Full)**

<input type="checkbox"/>	SOVENT		Custom Vent Tray			\$2,651
		68	<input type="checkbox"/> Articulating Vent Box (Required with Power Tilt and/or Recline)	Model: _____		\$1,273
		68	<input type="checkbox"/> Fixed Vent Box	Model: _____		\$477
<input type="checkbox"/>	SOVTBE	68, 69, 70	<input type="checkbox"/> Extended Battery Compartment			\$583
<input type="checkbox"/>			The Vent is being powered by 1 U1 Battery	<input type="checkbox"/> Chair Powered by 2 Group 22NF Batteries	<input type="checkbox"/> Chair Powered by 2 Group 24 Batteries	
<input type="checkbox"/>			The Vent is being powered by 1 Group 22NF Battery	<input type="checkbox"/> Chair Powered by 2 Group 22NF Batteries		
<input type="checkbox"/>			The Vent is being powered by 1 Group 24 Battery	<input type="checkbox"/> Chair Powered by 2 Group 22NF Batteries		

**Note 68** Amysystems does not supply the battery, charger or the wiring for the Vent.  
These items are normally provided by the Vent Manufacturer.

**Note 69** This option is not available with the Straight Front Caster Arms

**Note 70** The Vent Tray option is not available with the CUSTOM Battery Tray for Group 24 Batteries with a 24" frame width (FS2424B)

**OTHER CUSTOM OPTIONS (Request For Quote)**

<input type="checkbox"/>	PAT	71	Power Anterior Tilt			RFQ	
<input type="checkbox"/>	CUSTOM	71	Power Lateral Tilt 15°	<input type="checkbox"/> Raise the Left Side only	<input type="checkbox"/> Raise the Right Side only	<input type="checkbox"/> Raise the Left & Right Side	RFQ
<input type="checkbox"/>	CUSTOM	71, 73	Power Elevating Stretcher (Please call CSR for availability with configuration and user weight capacity)			RFQ	
<input type="checkbox"/>	CUSTOM	71, 73	Power Flip-up Foot Platform (Please call CSR for availability with configuration and user weight capacity)			RFQ	
<input type="checkbox"/>	CUSTOM	71, 73	Power Vertical Lifting Foot Platform (Please call CSR for availability with configuration and user weight capacity)			RFQ	
<input type="checkbox"/>	CUSTOM	71, 72, 74	Power Stander (Requires Power Tilt, Power Recline and 90° Power elevating & articulating center mount foot platform)			RFQ	
<input type="checkbox"/>	CUSTOM		Transfer Handles (Pair)			RFQ	
<input type="checkbox"/>	CUSTOM	71	Independent Depth Adjustability built into the Footrest Receiver that is attached to the Seat Pan			RFQ	
<input type="checkbox"/>	CUSTOM		Manual Wheel Locks			\$406.00	
<input type="checkbox"/>	CUSTOM		Non Standard Frame Colors (Not all custom colors are available, please contact CSR for availability)			\$812.00	
<input type="checkbox"/>	CUSTOM # 1	71	Other custom options. Specify in notes section below.			RFQ	
<input type="checkbox"/>	CUSTOM # 2	71	Other custom options. Specify in notes section below.			RFQ	

**Note 71** Please specify details and specifications for these options.

**Note 72** Please call CSR to verify if the power standing option is available with the required configuration for your user.

**Note 73** This option is not available with the Straight Front Caster Arms on Rear Hybrid Configuration

**Note 74** The power Stander is only available with the mid-wheel drive base configuration

**NOTES / SPECIAL REQUESTS / CUSTOM SPECIFICATIONS**

---



---



---



---