



April 1, 2017

ACCOUNT INFORMATION

Dealer Contact Name, Email Address, Phone Number, Fax Number, User Name or Number, User Weight, User Height, Seat Cushion Thickness, Knee to Heel Length, Bill to (Name), Bill to (Location), Ship to (Name), Address, City, Province, Postal Code, P.O. # or Quote, Left, Right, Seat to Top of Scapula

The order cannot be processed without the following:
- A user name or number, the client's measurements and weight
- Sunrise Medical is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Amysystems summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

STANDARD OPTIONS

The Quickie Xcel Power Base, 6-wheel independent wheel suspension, an additional seat suspension system, tie-down hardware for strap systems, offboard charger, 8" front casters, and standard auto style buckle seat positioning strap, 5.4mph motors, Weight Capacity = 550lbs

QUICKIE XCEL 2 POWER WHEELCHAIR

- XCEL2 QUICKIE XCEL -2 Power Wheelchair \$10,679

USER WEIGHT PACKAGES

- WCA550 550 lb Maximum user weight capacity STD

STANDARD FRAME COLORS

- FCBLK Black STD
FCSBM Matte Black STD
FCGMC Gun Metal (Charcoal) STD
FCSVM Silver Metallic STD
FCCBU Candy Blue STD
FCCRE Candy Red STD
FCMPL Midnight Purple STD
FCRZB Razzleberry STD
FCMWT New Matte White STD

OPTIONNAL FRAME COLORS (May Increase lead times)

- FCTXG Toxic Green N/C
FCBOR Burnt Orange Metallic N/C
FCBAO Black Opal Metallic N/C
FCOPM Sparkle Orange N/C
FCYLW Yellow N/C
FCPKP Pearl Pink N/C
FCHDD 1 Desert Digital Camo Hydro Graphics \$320
FCHWL 1 Woodland Camo Hydro Graphics \$320
CUSTOM 1 Carbon Fiber Hydro Graphics (Matte Black base) \$320

Note 1 Includes matching accessories with the Hydro Graphics option. The matching accessories option is not required with these colors

MATCHING ACCESSORIES TO THE FRAME COLOR

- FC\*\*-M Matching Accessories to the Frame Color (Wheels, front access door, headrest mount and back shell when selected) \$242

TRANSPORT OPTION

- TTD4 2 4-Point wheelchair transport brackets N/C

Note 2 4-Point wheelchair transport brackets have not yet received WC-19 crash-testing approval. Use these transport brackets only to secure an unoccupied wheelchair during transport.

MOTOR PACKAGE

- MP544 5.4MPH 4-Pole Heavy-Duty Motor Package N/C

ASSISTIVE DRIVING PACKAGE

- INTELLIDR 3, 4 Intelli-Drive \$1,597

Note 3 Intelli-Drive technology provides increased tracking capabilities with enhanced traction and steering precision. This is ideal for those with marginal hand function or switch control.
Note 4 Requires the R-net PM 120 Amp Expandable Controller (ECR120)

BATTERY CHARGER

- BCOB Offboard Battery Charger N/C
BCLES8OPT Lester 8 amp Offboard Battery Charger \$361

**BATTERY TRAY AND FRAME SIZE**

FS24B Battery Tray for Group 24 Batteries (Total frame width = 27 3/8") N/C

**BATTERIES**

BAT24PR Gel Batteries - Group 24 (Each - On chair) \$830

**DRIVE WHEEL**

DW14PT 14" Cast Wheel with Pneumatic Tire STD  
 DW14PI 14" Cast Wheel with Pneumatic Tire with Flat-Free inserts \$242

**FRONT CASTERS**

FC8PTM-STD 8" Pneumatic Casters STD  
 FC8PIM 8" Pneumatic Casters with Flat-Free inserts \$242

**SEAT TO FLOOR HEIGHT (Approximate measurement to the top of seat pan ± ¼")****Options for a Seat Only**

STFH1625 16 ¼" (Possible Seat angle of 0° only)  
 STFH1725 17 ¼" (Possible Seat angle of 0° and 3°)  
 STFH1825 18 ¼" (Possible Seat angle of 0° and 3°)  
 STFH1925 19 ¼" (Possible Seat angle of 0°, 3° and 6°)  
 STFH2025 20 ¼" (Possible Seat angle of 0°, 3°, 6°, 9°)

**Options for the 45° HD Power Tilt**

STFH17 17" (Possible Seat angle of 0° only)  
 STFH175 17 ½" (Possible Seat angle of 0° only)  
 STFH18 18" (Possible Seat angle of 0° only)  
 STFH185 18 ½" (Possible Seat angle of 0° only)  
 STFH19 19" (Possible Seat angle of 0° only)

**Options for the 11" HD Power Elevating Seat Lift Only**

STFH175 17 ½" (Possible Seat angle of 0° only)

**Options for a Combo - Power Tilt with an HD Power Elevating Seat**

STFH1925 19 ¼" (Possible Seat angle of 0° only)

**SEAT TO FLOOR HEIGHT SPECIFICATION (If not standard)**

SOMHSS Specific seat to floor height requirement other than standard. RFQ  
 Custom upcharges may apply. Required at: \_\_\_\_\_" (± ¼" Range)

**PRE-SET SEAT FRAME ANGLE**

SFE-0 0° Pre-set seat frame seat angle (**Must be selected when using power options**) STD  
 SFE-3 3° Pre-set seat frame seat angle (Only available with a seat configuration without power options) N/C  
 SFE-6 6° Pre-set seat frame seat angle (Only available with a seat configuration without power options) N/C  
 SFE-9 9° Pre-set seat frame seat angle (Only available with a seat configuration without power options) N/C

**LIGHTING PACKAGE**

EOLP 5 LED Light Package (Includes front beams, side flashers, rear tail light with flashers & brake light) \$1,066

**Note 5** Requires the upgraded Joystick with Light Operation or the OMNI (ECVJ6KL, ECELJL, ECRJCL, ECRJCTPL or EDIB)  
 The Light Package is not compatible with the EPSEL option.

**NON-EXPANDABLE ELECTRONICS - PG DRIVES VR2 CONTROLLER AND JOYSTICKS**

<input type="checkbox"/>	ECV902A		VR2 90 Amp Drive Controller		STD
<input type="checkbox"/>	ECVJ4K	6	VR2 Standard 4-Key Joystick		STD
<input type="checkbox"/>	ECVJ6K	6	VR2 6-Key Joystick (Allows for a single power option to operate through this joystick)		N/C
<input type="checkbox"/>	ECVJ9KL	6	VR2 9-Key Joystick with Lighting Operation (Allows for a single power function to operate through this joystick)		\$158.00

**Note 6** Requires the VR2 90A Drive Control (ECV902A).

VR2 Electronics can only accommodate a single power seat function. The only choices are Power Tilt or Elevating seat

**EXPANDABLE ELECTRONICS - PG DRIVES RNET CONTROLLERS**

<input type="checkbox"/>	ECEL90		R-net EL 90 Amp Expandable Controller		\$530.00
<input type="checkbox"/>	ECR120		R-net PM 120 Amp Expandable Controller		\$1,066.00
<input type="checkbox"/>	ECRNH	7	Harness for Expandable Controller		N/C

**Note 7** Must be ordered with the expandable controller ECEL90 or ECR120 to interface with Joysticks or OMNI screen

**EXPANDABLE ELECTRONICS - PG DRIVES RNET JOYSTICKS**

<input type="checkbox"/>	ECO		Omit Joystick (For use with an Omni Display)		N/C
<input type="checkbox"/>	ECELJ	8	R-net (LED) Joystick		STD
<input type="checkbox"/>	ECRJJC	8	R-net Joystick with Color screen		\$158

**JOYSTICKS WITH LIGHTING OPERATION**

<input type="checkbox"/>	ECELJL	8	R-net (LED) Joystick with Lighting Operation		\$158
<input type="checkbox"/>	ECRJCL	8	R-net Joystick with Color screen and Lighting Operation		\$210
<input type="checkbox"/>	ECRJCTPL	8	R-net Joystick with Color screen w/ On & Off / Mode switch & speed pot and Lighting Operation		\$315

**ADVANCED JOYSTICK**

<input type="checkbox"/>	ECRJJC35 New	8	R-net Advanced Joystick (3.5" Color screen). Lighting Operation, paddle switches, shortcut keys and stereo jacks		\$1,133
--------------------------	--------------	---	--	--	---------

**Note 8** This option requires the expandable controller ECEL90 or ECR120

**JOYSTICK MOUNT**

<input type="checkbox"/>	JSMF		Fixed Joystick Mount		N/C
<input type="checkbox"/>	JMSA		Swing Away Joystick Mount - Height adjustable, Angle adjustable and Multi-Axis		\$388
		<input type="checkbox"/>	Set-up with standard installation to swing-away sideways	<input type="checkbox"/>	Set-up at 35° to swing-away sideways and downwards
<input type="checkbox"/>	JSMCM		Midline Mounted Swing Away Joystick Mount (E1028)		\$388

**JOYSTICK MOUNT SUPPORT**

<input type="checkbox"/>	JSMFH		Joystick mount support Bracket		STD
--------------------------	-------	--	--------------------------------	--	-----

**JOYSTICK HAND MOUNT**

<input type="checkbox"/>	HML		Install Joystick on Left hand side	<input type="checkbox"/>	HMR		Install Joystick on Right hand side
--------------------------	-----	--	------------------------------------	--------------------------	-----	--	-------------------------------------

**JOYSTICK HANDLE**

<input type="checkbox"/>	JSH		Standard Joystick Gimble		STD
<input type="checkbox"/>	JBAT		Stick Handle 3" in Length		\$95
<input type="checkbox"/>	JBH		Ball Handle		\$95
<input type="checkbox"/>	JFBH		Foam Ball Handle		\$95
<input type="checkbox"/>	JTH		T-Handle		\$95
<input type="checkbox"/>	JGP		Goalpost		\$95

**ATTENDANT DRIVE CONTROL OPTION**

<input type="checkbox"/>	EADC		Attendant Hand Control for R-net EL & PM Controls		\$1,066
--------------------------	------	--	---	--	---------

**ATTENDANT JOYSTICK HAND MOUNT**

<input type="checkbox"/>	EADCL		Install Attendant Joystick on Rear Left hand side	<input type="checkbox"/>	EADCR		Install Attendant Joystick on Rear Right hand side
--------------------------	-------	--	---	--------------------------	-------	--	--

**ENHANCED DISPLAY & DIGITAL INTERFACE BOX (Included with specialty controls)**

<input type="checkbox"/>	EDIB	9	Omni (Enhanced Display & Digital Interface Box with Infrared)		\$2,033
--------------------------	------	---	---	--	---------

**Note 9** Allows the switch operation of alternative devices through driver controls.

Requires one of the R-net fully Expandable Controllers (ECEL90 or ECR120)

**ENHANCED DISPLAY MOUNT**

<input type="checkbox"/>	EDSAM		User Swing Away Mount (E1028)	<input type="checkbox"/>	EDML		Install Display on LH side	<input type="checkbox"/>	EDMR		Install Display on RH side	N/C
--------------------------	-------	--	-------------------------------	--------------------------	------	--	----------------------------	--------------------------	------	--	----------------------------	-----

**ENVIRONMENTAL CONTROL MODULE**

ECU 10 Environment control module - (Output Module) \$1,045

**Note 10 Requires one of the fully expandable R-net controllers (ECEL90 or ECR120)**

**BLUETOOTH MOUSE EMULATOR (Includes USB Dongle)**

EDBT 11 Bluetooth Mouse Emulator (USB Dongle included) \$1,066  
 Windows and MAC OS       Android Devices       Iphone, Ipad and Ipod

**Note 11 Requires one of the fully expandable R-net controllers (ECEL90 or ECR120)  
 Not limited to only one selection. Any two can be installed on the same chair, but only one can be used at a time.**

**POWER OPTIONS & SEATING**

SAABAHD Bariatric Amyseat - Seat Frame only \$1,355  
 SAPTXVHD 12 Bariatric CG Power Tilt \$8,568  
 SAPRXVHD 12, 14 Bariatric Power Recline with Integrated Shear Reduction (ISR) \$9,318  
 SAPTPRXVHD 12, 14 Bariatric CG Power Tilt and Power Recline with Integrated Shear Reduction (ISR) \$13,926  
 SAPEXVHD 12, 13 11" Heavy-Duty Power Elevating Seat Module for transfer assist \$5,353  
 ADDPESXHD 13 **Module** - Add a Heavy-Duty 11" Power Elevating Seat to either SAPTXVHD, SAPRXVHD or SAPTPRXVHD \$4,282

**Note 12 These options include the fully adjustable Amyseat (SAABAHD).  
 Note 13 The maximum user weight capacity with this option is 400 lbs  
 Note 14 It is recommended to order the Contour Back Support option with the power recline (STDSB)**

**ELECTRONICS - SINGLE OR DUAL (1 OR 2 ACTUATORS) POWER SEAT FUNCTION CONTROLLERS**

EPSV 15, 16 Tilt or Elevating seat activation through **6-key VR2** Hand Control (ECVJ6K or ECVJ9KL) N/C  
 EPSEL 15, 16 Tilt or Elevating seat activation through **RNET 90 Amp Expandable Controller (ECEL90)** N/C  
 SAC Single (1 Actuator) Power Seat Function Controller to Interface through R-net Electronics N/C  
 DAC Dual (2 Actuators) Power Seat Function Controller to Interface through R-net Electronics \$530

**With the SAC or DAC, please select one activation and position from the list below**

SWTJSS 17, 18 Single actuator activation through R-net Joystick or Specialty Controls \$1,066  
 SWTJSD 17, 19 Dual actuator activation through R-net Joystick or Specialty Controls \$1,066  
 SWDUALACT Dual Activation (**Allows to operate all seat functions with external switches AND through the joystick**) \$330  
 SWTOGGLE 18 or 19 Single Toggle Switch (Single Mode: DAC - Dual Mode: SAC)     User Activation     Attendant Activation N/C  
 SWPUSH1 18 Single Push Button Switch (Single Mode: SAC)     User Activation     Attendant Activation N/C  
 SWPUSH2 18 or 19 Dual Push Button Switch (Single Mode: DAC - Dual Mode: SAC)     User Activation     Attendant Activation N/C

**Note 15 Allows for the activation of a single power option through the VR2 or R-net 90 Amp Controllers. Only the Power Tilt or Elevating seat options are available with this feature.**  
**Note 16 These options do not require the SAC controller as they connect via the chair's power module directly.**  
**Note 17 Requires the R-net 90 Amp (ECEL90) or the 120 Amp (ECR120) Fully Expandable Controllers**  
**Note 18 Requires the Single (1 Actuator) Power Seat Function Controller (SAC)**  
**Note 19 Requires the Dual (2 Actuators) Power Seat Function Controller (DAC)**

**ELECTRONICS - MULTIPLE (3 to 6 ACTUATORS) POWER SEAT FUNCTION CONTROLLER**

MAC Multiple (3 to 6 Actuators) Power Seat Function Controller to Interface through R-net Electronics \$2,354

**With the MAC, please select one activation and position from the list below**

SWTJSM 20, 21 Multiple Actuator activation through R-net Joystick or Specialty Controls STD  
 SWDUALACT Dual Activation (Allows to operate all seat functions with external switches **and** through the joystick) \$330  
 SWMULTI 20 4-Way Toggle Switch (Up to 4 Activations)    -     User Activation     Attendant Activation N/C  
 SWMULTIT 20 Triple Toggle Switch (Up to 5 Activations)    -     User Activation     Attendant Activation N/C  
 SWPUSH4 20 4-Way Push Button Switch    -     User Activation     Attendant Activation N/C  
 SWPUSH5 20 5-Way Push Button Switch    -     User Activation     Attendant Activation N/C

**Note 20 Requires the Multiple (3 to 6 Actuators) Power Seat Function Controller (MAC)**  
**Note 21 Requires the R-net 90 Amp (ECEL90) or 120 Amp (ECR120) Fully Expandable Controller, and the multiple Actuator Controller (MAC)**

**SWITCH ACTIVATION HAND MOUNT**

Switch Installation:     TSML (Install Switch on LH side)     TSMR (Install Switch on RH side)    N/C

**PRE-SET SEAT ANGLE (For Non-Reclining Backs)**

SOBAA Please Select one:  82°  90°  98°  106°

N/C

**STANDARD ADJUSTABLE SEAT WIDTH RANGE BETWEEN 19" and 28"** STD

- SOWN19-24-22 (Width adjustment from 19" to 24" with a maximum Depth of 22" )
- SOWN24-28-24 (Width adjustment from 24" to 28" with a maximum Depth of 24" )

- |  |  |
|--|--|
| <input type="checkbox"/> 19" (SOWN19L) | <input type="checkbox"/> 24" (SOWN24HDL) |
| <input type="checkbox"/> 20" (SOWN20L) | <input type="checkbox"/> 25" (SOWN25HDL) |
| <input type="checkbox"/> 21" (SOWN21)  | <input type="checkbox"/> 26" (SOWN26HDL) |
| <input type="checkbox"/> 22" (SOWN22)  | <input type="checkbox"/> 27" (SOWN27HDL) |
| <input type="checkbox"/> 23" (SOWN23)  | <input type="checkbox"/> 28" (SOWN28HDL) |
| <input type="checkbox"/> 24" (SOWN24)  |  |

**ADJUSTABLE SEAT DEPTH RANGE BETWEEN 17" and 24"**

- |                                       |                                       |                                       |                                       |
|---------------------------------------|---------------------------------------|---------------------------------------|---------------------------------------|
| <input type="checkbox"/> 17" (SODN17) | <input type="checkbox"/> 18" (SODN18) | <input type="checkbox"/> 19" (SODN19) | <input type="checkbox"/> 20" (SODN20) |
| <input type="checkbox"/> 21" (SODN21) | <input type="checkbox"/> 22" (SODN22) | <input type="checkbox"/> 23" (SODN23) | <input type="checkbox"/> 24" (SODN24) |

**BACK CANE HEIGHT (The back cane height measurement is taken from the seat pan to the top of the push handles)**

- |  |  |  |
|--|--|--|
| <input type="checkbox"/> 16" (SOBCH16) | <input type="checkbox"/> 20" (SOBCH20)               | <input type="checkbox"/> 24" (SOBCH24) |
| <input type="checkbox"/> 18" (SOBCH18) | <input type="checkbox"/> 22" (SOBCH22) <sup>22</sup> | <input type="checkbox"/> 26" (SOBCH26) |

Note 22 22" is the minimum back cane height required with a power recline.

**BACK SUPPORT**

- |  |  |       |
|--|--|-------|
| <input type="checkbox"/> STDSB <sup>23</sup> | Contour back support with a maximum cushion height of 20" (19" and 20" Wide)           | \$588 |
| <input type="checkbox"/> STDSB/CUSTOM        | Contour back support (Non-Standard Sizes) <b>Specify:</b> <b>Width:</b> <b>Height:</b> | \$698 |

Note 23 It is recommended to use the Contour Back Support option with the power or manual recline.

<sup>25</sup> **BACK SUPPORT TOTAL HEIGHT (MEASURED FROM THE SEAT PAN TO THE TOP OF THE BACK)**

- |                                       |   |                                       |                                       |
|---------------------------------------|---|---------------------------------------|---------------------------------------|
| <input type="checkbox"/> 16" (SOBH16) | <input type="checkbox"/> 20" (SOBH20)               | <input type="checkbox"/> 24" (SOBH24) | <input type="checkbox"/> 28" (SOBH28) |
| <input type="checkbox"/> 18" (SOBH18) | <input type="checkbox"/> 22" (SOBH22) <sup>24</sup> | <input type="checkbox"/> 26" (SOBH26) |                                       |

Note 24 With a power recline, the minimum height is 22".

Note 25 The space between the bottom of the back and the seat pan is 4" with a non-reclining back and 6" with a reclining back

**SEAT CUSHION**

- |                                       |   |         |
|---------------------------------------|---|---------|
| <input type="checkbox"/> SCPSC        | Contour Seat Cushion with a maximum seat depth of 22" (19" and 20" Wide)              | \$588   |
| <input type="checkbox"/> SCPSC/CUSTOM | Contour Seat Cushion (Non-Standard Sizes) <b>Specify:</b> <b>Width:</b> <b>Depth:</b> | \$1,281 |

**CANTILEVER - HEIGHT ADJUSTABLE ARMRESTS - OPENS WITH A 5° ANGLE (For Non-Reclining Backs - Choose 2)**

- |                                   |   |     |
|-----------------------------------|---|-----|
| <input type="checkbox"/> SOAFBFLH | Flip-Back Cantilever - Full Length Left Side (Height Adjustable from 7" to 15" in 1" increments)  | STD |
| <input type="checkbox"/> SOAFBFRH | Flip-Back Cantilever - Full Length Right Side (Height Adjustable from 7" to 15" in 1" increments) | STD |
| <input type="checkbox"/> SOAFBDLH | Flip-Back Cantilever - Desk Length Left Side (Height Adjustable from 7" to 15" in 1" increments)  | N/C |
| <input type="checkbox"/> SOAFBDRH | Flip-Back Cantilever - Desk Length Right Side (Height Adjustable from 7" to 15" in 1" increments) | N/C |

**DUAL POST - HEIGHT ADJUSTABLE ARMRESTS - INCLUDES TRANSFER GRIPS (For Non-Reclining Backs - Choose 2)**

- |                                   |  |       |
|-----------------------------------|--|-------|
| <input type="checkbox"/> SOADPFLH | Dual Post Height Adjustable - Full Length Left Side (Height Adjustable from 9 ½" - 13 ½")  | \$190 |
| <input type="checkbox"/> SOADPFRH | Dual Post Height Adjustable - Full Length Right Side (Height Adjustable from 9 ½" - 13 ½") | \$190 |
| <input type="checkbox"/> SOADPDLH | Dual Post Height Adjustable - Desk Length Left Side (Height Adjustable from 9 ½" - 13 ½")  | \$190 |
| <input type="checkbox"/> SOADPDRH | Dual Post Height Adjustable - Desk Length Right Side (Height Adjustable from 9 ½" - 13 ½") | \$190 |

**RECLINING FLIP-UP, HEIGHT ADJUSTABLE ARMRESTS (For Reclining Backs - Choose 2)**

- |                                   |   |     |
|-----------------------------------|---|-----|
| <input type="checkbox"/> SOANIFLH | Reclining armrest - Full Length Left Side (Height Adjustable from 10 ½" to 14" in ½" increments)  | N/C |
| <input type="checkbox"/> SOANIFRH | Reclining armrest - Full Length Right Side (Height Adjustable from 10 ½" to 14" in ½" increments) | N/C |
| <input type="checkbox"/> SOANIDLH | Reclining armrest - Desk Length Left Side (Height Adjustable from 10 ½" to 14" in ½" increments)  | N/C |
| <input type="checkbox"/> SOANIDRH | Reclining armrest - Desk Length Right Side (Height Adjustable from 10 ½" to 14" in ½" increments) | N/C |

**ARMREST HEIGHT ADJUSTMENT AT FACTORY**

Please set the Armrest Height at: \_\_\_\_\_ " (The Armrest Height is taken from the seat pan to the top of the armpad. Not required for Dual Post)

**ARMPADS**

**Foam Type Selection**

<input type="checkbox"/>	SOFAFD	Desk Length Standard (10")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	STD
<input type="checkbox"/>	SOFAFP	Full Length Standard (14")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	STD
<input type="checkbox"/>	SOWAPD	34 Desk Length Waterfall (10")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$32 ea.
<input type="checkbox"/>	SOWAFP	34 Full Length Waterfall (14")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$32 ea.

**Hybrid Gel and Foam Type Selection (The Gel always stays soft)**

<input type="checkbox"/>	SOGELD new	Desk Length Hybrid Gel and Foam Flat armpad (10" x 2.5")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$32 ea.
<input type="checkbox"/>	SOGELF	Full Length Hybrid Gel and Foam Flat armpad (14" x 2.5")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$32 ea.
<input type="checkbox"/>	SOWGELD new	34 Desk Length Hybrid Gel and Foam Waterfall armpad (10")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$50 ea.
<input type="checkbox"/>	SOWGELF new	34 Full Length Hybrid Gel and Foam Waterfall armpad (14")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$50 ea.

**Specialty Type Selection**

<input type="checkbox"/>	SOCAR	35 Arm Trough with Hand Support (18" x 4 1/2")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$210
--------------------------	-------	--	---	--------------	---	--	-------

**Note 27** The Waterfall armpads may interfere with flipping back the armrests when combined with Lateral Trunk supports

**Note 28** Arm Trough Armpads may interfere with Armrests Flipping back.

**FOOTRESTS**

<input type="checkbox"/>	SOAAF**-STD	32 <b>Heavy Duty</b> Swing-Away footrests ----- Set at: <input type="checkbox"/> 60° <input type="checkbox"/> 70°					\$370
		Select Footrest Range:	<input type="checkbox"/> SOFLHD12 (12 1/2" - 15 1/2")	<input type="checkbox"/> SOFLHD14 (14 1/2" - 19 1/2")			
<input type="checkbox"/>	ADDPLHD	29, 32 <b>Heavy Duty</b> Power Elevating & Articulating Legrests	Select Function:	<input type="checkbox"/> Individual <sup>30</sup>	<input type="checkbox"/> Combined <sup>31</sup>		\$2,994
		Select Footrest Range:	<input type="checkbox"/> SOFLPW12 (12 1/2" - 15 1/2")	<input type="checkbox"/> SOFLPW14 (14 1/2" - 19 1/2")			

32 Please set the Footrest Length at: \_\_\_\_\_ " (The footrest length is taken from the seat pan to the heel)

**Note 29** This option includes calf pads

**Note 30** Individual legs count as two power functions.

**Note 31** Combined legs count as one power function.

**Note 32** The Footrest Length is taken from the top of the seat pan to the top of the footplate.

**FOOTPLATES & ACCESSORIES**

<input type="checkbox"/>	SOSFAFHD	Flip-up Footplates					STD
<input type="checkbox"/>	SOAFPHD	Angle Adjustable & Flip-up Footplates - Each					\$137
<input type="checkbox"/>	SOAHL	Heel Loops - Each					\$42
<input type="checkbox"/>	SOACPFR	Calf Pads for Adult Heavy-Duty Swing-Away footrests (SOAAF) with adjustable hardware - Each					\$137

**CENTER MOUNTS**

<input type="checkbox"/>	SOAAF90HD	36, 37, 38, 40 90° <b>Heavy Duty</b> Fixed Center Mount Interface Bracket ( <b>Offers Manual adjustability from 95° to 70°</b> )					\$242
		Select Foot Platform Range:	<input type="checkbox"/> SOFBP7 (7 1/2" - 12 1/2")	<input type="checkbox"/> SOFBA12 (12 1/2" - 16 1/2")	<input type="checkbox"/> SOFBA16 (16 1/2" - 19 1/2")		
<input type="checkbox"/>	ADDPFHHD	33, 35, 36, 38, 39, 40 90° <b>Heavy Duty</b> Power Center Mount Elevating & Articulating Foot Platform - (Uses 1 center mounted actuator)					\$3,204
		<input type="checkbox"/> SOFBPW12 (12 1/2" - 16 1/2")	<input type="checkbox"/> SOFBPW16 (16 1/2" - 20 1/2")				
<input type="checkbox"/>	ADDPFVHD	34, 35, 36, 38, 40 90° <b>Very Heavy Duty</b> Power Center Mount Elevating & Articulating Foot Platform - (Uses 2 side mounted actuators)					\$3,500
		<input type="checkbox"/> SOFBPW15 (15" - 19")	<input type="checkbox"/> SOFBPW19 (19" - 23")				

40 Please set the Foot Platform Length at: \_\_\_\_\_ "

**Note 33** ADDPFHD counts as one power function.

**Note 34** ADDPFVHD counts as two power functions. Please select (SOFUFP) on the next page

**Note 35** This option includes calf pads

**Note 36** Links and attaches the one piece center mount foot platform (SOFUFP) or the split footplates (SOAFPFB) to the seat

**Note 37** With the split footplates option (SOAFPFB), this bracket offers independent height adjustability for each footplate

**Note 38** Please Select the one piece center mount foot platform (SOFUFP) or the split footplates (SOAFPFB) if compatible (See footnotes)

**Note 39** The Standard user weight capacity is 450 Lbs.

**Note 40** The Foot Platform Length is taken from the top of the seat pan to the top of the foot platform or split foot plates.

**FOOT PLATFORMS & ACCESSORIES**

<input type="checkbox"/>	SOFUFP	Flip-up one piece Foot Platform					N/C
<input type="checkbox"/>	SOAFPFB	Angle Adjustable & Flip-up Split Footplates for center mount - Each					\$137
<input type="checkbox"/>	SOACPCM	Calf Pads for 90° Fixed Center Mount (SOAAF90) with adjustable hardware - Each					\$137

### SEAT POSITIONING

<input type="checkbox"/>	SOAUTO80-STD	80" Seat Positioning Strap	\$79
--------------------------	--------------	----------------------------	------

### HEADREST

<input type="checkbox"/>	SONUAHR	Using headrest mounting hardware other than SOMAH located below. Specify: _____	
<input type="checkbox"/>	SOFHRS	Small Curved Headrest Pad 4" x 9"	\$215
<input type="checkbox"/>	SOFHRM	Medium Curved Headrest Pad 5" x 10"	\$215
<input type="checkbox"/>	SOFHRL	Large Curved Headrest Pad 6" x 11"	\$215
<input type="checkbox"/>	SOMAH	Multi-axis Removable Mounting Hardware	\$173
<input type="checkbox"/>	SOHRBC	Bracket to Mount on the Contour Back Support <b>(Without Recline)</b>	\$105
<input type="checkbox"/>	SOHRPR	Bracket to Mount on the Seat Frame <b>(For Reclining backs)</b>	\$105

\* Please select one of the headrest options with the mounting hardware and specify where the headrest will be installed

### LATERAL TRUNK SUPPORTS

<input type="checkbox"/>	SOSLS3	Adjustable Laterals with Small Curved Pads (Each) - 3" x 5"	\$214
<input type="checkbox"/>	SOSLM4	Adjustable Laterals with Medium Curved Pads (Each) - 4" x 6"	\$214
<input type="checkbox"/>	SOSLL5	Adjustable Laterals with Large Curved Pads (Each) - 5" x 7"	\$214
<input type="checkbox"/>	SOSLHM	Adjustable Swing-Away Hardware Mounts - Each	\$316

\*\* Please select one of the lateral pads and the adjustable & removable hardware.

### HIP & THIGH SUPPORTS

<input type="checkbox"/>	SOHTSS10	Adjustable Hip & Thigh Support Pads (Each) - Small 10" x 4"	\$242
<input type="checkbox"/>	SOHTSM12	Adjustable Hip & Thigh Support Pads (Each) - Medium 12" x 5"	\$242
<input type="checkbox"/>	SOHTSL14	Adjustable Hip & Thigh Support Pads (Each) - Large 14" x 6"	\$242
<input type="checkbox"/>	SOHTSARH	Adjustable and Removable Hip & Thigh Support Hardware - Each	\$457

\*\*\* Please select one of the hip & thigh support pads and the adjustable & removable hardware.

### ACCESSORIES

<input type="checkbox"/>	SOCH	Cup Holder (Self-Leveling)	\$95
<input type="checkbox"/>	SOPHP	Ram Mounts X-Grip Cell Phone Holder	\$163
<input type="checkbox"/>	SOTHS	Ram Mounts X-Grip Small Tablet Holder	\$184
<input type="checkbox"/>	SOTHL	Ram Mounts X-Grip Large Tablet Holder with Swivel	\$278

### REFLECTOR KIT

<input type="checkbox"/>	REFLECT	Red & White Reflector Kit Stickers (They will not be installed - Included in a bag with the owner's manual)	N/C
--------------------------	---------	---	-----

### SWITCHES AND OPTIONS

<input type="checkbox"/>	ASCPCP	PGDT RNET Dongle and PC programmer kit	\$615
<input type="checkbox"/>	ASCPVR2	PGDT VR2 Dongle and PC programmer kit	\$615
<input type="checkbox"/>	SSBS	Buddy Button Switch	\$95
<input type="checkbox"/>	SSES	Egg Switch-Black	\$95
<input type="checkbox"/>	SSML	Micro Light Switch	\$158

## SPECIALTY DRIVER CONTROL OPTIONS

### Therafin Sip-N-Puff Kit

<input type="checkbox"/>	SCSNPK	41, 42	Therafin Sip-N-Puff Kit		\$1,603
<input type="checkbox"/>	SCSNPBTk	42	Therafin Sip-N-Puff Breath Tube Kit		\$520
<input type="checkbox"/>	SCSRMSNP	42	Stealth Rock-n-Lock Swingaway Arm for Sip-N-Puff	<input type="checkbox"/> SWL (Swings on Left) <input type="checkbox"/> SWR (Swings on Right)	\$289
<input type="checkbox"/>	SCSRM	42	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/> SWL (Swings Left) <input type="checkbox"/> SWR (Swings Right)	\$520
<input type="checkbox"/>	SSES-STD	42, 44	Egg Switch for Changing Modes - Black	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$115
<input type="checkbox"/>	SSBS-STD	42	Buddy Button Switch (Stop Switch) - Red		\$115
<input type="checkbox"/>	SOMAH	42	Multi-axis Headrest Removable Mounting Hardware		\$173
<input type="checkbox"/>	SOFHRS		Small Curved Headrest Pad 4" x 9"		\$215
<input type="checkbox"/>	SOFHRM	43	Medium Curved Headrest Pad 5" x 10"		\$215
<input type="checkbox"/>	SOFHRL		Large Curved Headrest Pad 6" x 11"		\$215

### Switch-It Proximity Head Array

<input type="checkbox"/>	SCSIHAA	41, 42	Adult - 3 Switch Proximity Head Array with 1/8" Mono Jack for mode/reverse		\$5,353
<input type="checkbox"/>	SCSRM	42	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/> SWL (Swings Left) <input type="checkbox"/> SWR (Swings Right)	\$520
<input type="checkbox"/>	SSES-STD	42, 44	Egg Switch for Changing Modes - Black	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$115
<input type="checkbox"/>	SSBS-STD	42	Buddy Button Switch (Stop Switch) - Red		\$115
<input type="checkbox"/>	SOMAH	42	Multi-axis Headrest Removable Mounting Hardware		\$173

### Switch-It MicroPilot Light Touch Proportional Mini joystick for chin / lip control

<input type="checkbox"/>	SCSIMP	41, 42	MicroPilot Light touch mini Proportional Joystick		\$4,197
<input type="checkbox"/>	SCSRMMP	42	Stealth Rock-n-Lock Swingaway Mount for Chin M-P	<input type="checkbox"/> SWL (Swings Left) <input type="checkbox"/> SWR (Swings Right)	\$520
<input type="checkbox"/>	SCSRM	42	Stealth Rock-n-Lock Swingaway Mount for Switch	<input type="checkbox"/> SWL (Swings Left) <input type="checkbox"/> SWR (Swings Right)	\$520
<input type="checkbox"/>	SSES	42, 44	Egg Switch-Black	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$95
<input type="checkbox"/>	JCC		Chin Cup		\$53
<input type="checkbox"/>	SOMAH	42	Multi-axis Headrest Removable Mounting Hardware		\$173
<input type="checkbox"/>	SOFHRS		Small Curved Headrest Pad 4" x 9"		\$215
<input type="checkbox"/>	SOFHRM	43	Medium Curved Headrest Pad 5" x 10"		\$215
<input type="checkbox"/>	SOFHRL		Large Curved Headrest Pad 6" x 11"		\$215

### Proportional Compact or Mini joystick for Midline Mount

<input type="checkbox"/>	SCPGMJS	41, 42	PG Drives JC200 Proportional Compact Remote Joystick for Midline Mount		\$998
<input type="checkbox"/>	SCGMJSM	45	Gatlin Midline Swing-away Mount for JC200	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$410
<input type="checkbox"/>	SCSIMP	41, 42	Switch-It MicroPilot Light touch mini Proportional Joystick		\$4,197
<input type="checkbox"/>	SCGMPPM	46	Gatlin Midline Swing-away Mount for MicroPilot	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$447
<input type="checkbox"/>	SCBTM	47	Bullet Hand Tray for MicroPilot	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$473

### Proportional Compact or Mini joystick for Armrest Mount

<input type="checkbox"/>	SCPGMJS	41, 42	PG Drives JC200 Proportional Compact Remote Joystick for Armrest Mount		\$998
<input type="checkbox"/>	SCARMJSM	45	Swing-away and Multi-Axis armrest mount	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$410
<input type="checkbox"/>	SCSIMP	41, 42	Switch-It MicroPilot Light touch mini Proportional Joystick		\$4,197
<input type="checkbox"/>	SCARMPM	46	Swing-away and Multi-Axis armrest mount	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$447
<input type="checkbox"/>	SCBTM	47	Bullet Hand Tray for MicroPilot	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$473

### Proportional RIM Head Control Package

<input type="checkbox"/>	SCPHC	41, 42	Proportional RIM Head Control Package		\$2,621
<input type="checkbox"/>	SCSRM	42	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/> SWL (Swings Left) <input type="checkbox"/> SWR (Swings Right)	\$520
<input type="checkbox"/>	SSES-STD	42, 44	Egg Switch for Changing Modes - Black	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$115
<input type="checkbox"/>	SOMAH	42	Multi-axis Headrest Removable Mounting Hardware		\$173

Note 41 Requires the Omni Enhanced Display & Digital Interface Box (EDIB)

Note 42 This option must be selected for the Kit

Note 43 Choose one (1) Headrest Pad

Note 44 The Mode Switch Type can be changed from this selection on Page #7

Note 45 Must be selected for use with the PG Drives JC200 Mini Joystick

Note 46 Must be selected for use with the Switch-It MicroPilot

Note 47 This option is only compatible with the Switch-It MicroPilot



**48 CUSTOM SECTION**

**Note 48** Lead times will increase depending on the custom order pipeline at the time the purchase order is received  
Please call for Quote and specifications for these options.

**O<sup>2</sup> Holder**

<input type="checkbox"/>	SUPOXY	49	Custom Non-Articulating Oxygen Tank Holder	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	\$410
<input type="checkbox"/>	SUPOXYA	50	Custom Articulating Oxygen Tank Holder	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	\$1,923
			<input type="checkbox"/> D-Size Oxygen Tank (* 22" Tall x 4.3" Diameter) * The height includes the valve			
			<input type="checkbox"/> E-Size Oxygen Tank (* 31" Tall x 4.3" Diameter) * The height includes the valve			

**Note 49** Not available with a power tilt and/or power recline  
**Note 50** Required with a power tilt and/or power recline

**VENT TRAY (Please complete in Full)**

<input type="checkbox"/>	SOVENT	51	Custom Vent Tray	<input type="checkbox"/> Articulating Vent Box (Required with Power Tilt and/or Recline)	Model: _____	\$2,574
		51	<input type="checkbox"/> Fixed Vent Box	Model: _____		
<input type="checkbox"/>	SOVTBE	51	<input type="checkbox"/> Extended Battery Compartment			N/C
<input type="checkbox"/>	The Vent is being powered by 1 U1 Battery <input type="checkbox"/> Chair Powered by 2 Group 22NF Batteries <input type="checkbox"/> Chair Powered by 2 Group 24 Batteries					
<input type="checkbox"/>	The Vent is being powered by 1 Group 22NF Battery <input type="checkbox"/> Chair Powered by 2 Group 22NF Batteries					
<input type="checkbox"/>	The Vent is being powered by 1 Group 24 Battery <input type="checkbox"/> Chair Powered by 2 Group 22NF Batteries					

**Note 51** Amysystems does not supply the battery, charger or the wiring for the Vent.  
These items are normally provided by the Vent Manufacturer.

**OTHER CUSTOM OPTIONS (Request For Quote)**

<input type="checkbox"/>	CUSTOM		Manual Wheel Locks			RFQ
<input type="checkbox"/>	CUSTOM		Non Standard Frame Colors (Not all custom colors are available, please contact CSR for availability)			RFQ
<input type="checkbox"/>	CUSTOM		Transfer Handles (Pair)			RFQ
<input type="checkbox"/>	CUSTOM	52	Endomorphic Back Canes (Pair)	<input type="checkbox"/> 1" per side towards the inside = 2" in	<input type="checkbox"/> 2" per side towards the inside = 4" in	RFQ
<input type="checkbox"/>	CUSTOM # 1	52	Other custom options. Specify in notes section bellow.			RFQ
<input type="checkbox"/>	CUSTOM # 2	52	Other custom options. Specify in notes section bellow.			RFQ
<input type="checkbox"/>	CUSTOM # 3	52	Other custom options. Specify in notes section bellow.			RFQ

**Note 52** Please specify details and specifications for these options.

**NOTES / SPECIAL REQUESTS / CUSTOM SPECIFICATIONS**

---

---

---

---

---

---

---

---

---

---