



Quickie Xcel 2

October 1, 2017



ACCOUNT INFORMATION

Dealer Contact Name:	_____	Bill to (Name):	_____
Email Address:	_____	Bill to (Location):	_____
Phone Number:	_____	Ship to (Name):	_____
Fax Number:	_____	Address:	_____
User Name or Number:	_____	City:	_____
User Weight:	_____	Province:	_____
User Height:	_____	Postal Code:	_____
Seat Cushion Thickness:	_____	P.O. # or Quote:	_____
Knee to Heel Length:	Left: _____" Right: _____"		
Seat to Top of Scapula:	_____		

The order cannot be processed without the following:

- A user name or number, the client's measurements and weight
- Sunrise Medical is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Sunrise Medical summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

STANDARD OPTIONS

Quickie Xcel Power Base, 6-wheel independent wheel suspension, an additional seat suspension system, tie-down hardware for strap systems, offboard charger, 8" front casters, and standard auto style buckle seat positioning strap, 5.4mph motors, Weight Capacity = 550lbs

QUICKIE XCEL POWER WHEELCHAIR

- | | | | |
|--------------------------|-------|-----------------------------------|-------------|
| <input type="checkbox"/> | XCEL2 | QUICKIE XCEL - 2 Power Wheelchair | \$10,999.00 |
|--------------------------|-------|-----------------------------------|-------------|

USER WEIGHT PACKAGES

- | | | | |
|--------------------------|--------|-------------------------------------|-----|
| <input type="checkbox"/> | WCA550 | 550 lb Maximum user weight capacity | STD |
|--------------------------|--------|-------------------------------------|-----|

STANDARD FRAME COLORS

- | | | | | | | | |
|--------------------------|------------------|----------------------|-----|--------------------------|----------|--|-------|
| <input type="checkbox"/> | FCBLK | Black | STD | <input type="checkbox"/> | FCTXG | Toxic Green | N/C |
| <input type="checkbox"/> | FCSBM | Matte Black | STD | <input type="checkbox"/> | FCBOR | Burnt Orange Metallic | N/C |
| <input type="checkbox"/> | FCGMC | Gun Metal (Charcoal) | STD | <input type="checkbox"/> | FCBAO | Black Opal Metallic | N/C |
| <input type="checkbox"/> | FCSVM | Silver Metallic | STD | <input type="checkbox"/> | FCOPM | Sparkle Orange | N/C |
| <input type="checkbox"/> | FCCBU | Candy Blue | STD | <input type="checkbox"/> | FCYLW | Yellow | N/C |
| <input type="checkbox"/> | FCCRE | Candy Red | STD | <input type="checkbox"/> | FCPKP | Pearl Pink | N/C |
| <input type="checkbox"/> | FCMPL | Midnight Purple | STD | <input type="checkbox"/> | FCHDD 1 | Desert Digital Camo Hydro Graphics | \$330 |
| <input type="checkbox"/> | FCRZB | Razzleberry | STD | <input type="checkbox"/> | FCHWL 1 | Woodland Camo Hydro Graphics | \$330 |
| <input type="checkbox"/> | FCMWT New | Matte White | STD | <input type="checkbox"/> | CUSTOM 1 | Carbon Fiber Hydro Graphics (Matte Black base) | \$330 |

Note 1 Includes matching accessories with the Hydro Graphics option. The matching accessories option is not required with these colors

MATCHING ACCESSORIES TO THE SHROUD COLOR

- | | | | |
|--------------------------|---------|---|-------|
| <input type="checkbox"/> | FC***-M | Matching Accessories to the Shroud Color (Wheels, front access door & headrest mount when selected) | \$249 |
|--------------------------|---------|---|-------|

TRANSPORT OPTION

- | | | | | |
|--------------------------|------|---|---------------------------------------|-----|
| <input type="checkbox"/> | TTD4 | 1 | 4-Point wheelchair transport brackets | N/C |
|--------------------------|------|---|---------------------------------------|-----|

Note 1 4-Point wheelchair transport brackets have not yet received WC-19 crash-testing approval. Use these transport brackets only to secure an unoccupied wheelchair during transport.

MOTOR PACKAGE

- | | | | |
|--------------------------|-----------|--|-----|
| <input type="checkbox"/> | MP544-STD | 5.4MPH 4-Pole Heavy-Duty Motor Package | N/C |
|--------------------------|-----------|--|-----|

ASSISTIVE DRIVING PACKAGE

- | | | | | |
|--------------------------|-----------|------|---------------|---------|
| <input type="checkbox"/> | INTELLIDR | 2, 3 | Intelli-Drive | \$1,644 |
|--------------------------|-----------|------|---------------|---------|

Note 2 Intelli-Drive technology provides increased tracking capabilities with enhanced traction and steering precision.

This is ideal for those with marginal hand function or switch control.

Note 3 Requires the R-net PM 120 Amp Expandable Controller (ECR120)

BATTERY CHARGER

- | | | | |
|--------------------------|-----------|---------------------------------------|-------|
| <input type="checkbox"/> | BCOB | Offboard Battery Charger | N/C |
| <input type="checkbox"/> | BCLES8OPT | Lester 8 amp Offboard Battery Charger | \$366 |

BATTERY TRAY AND FRAME SIZE

FS24B Battery Tray for Group 24 Batteries (Total frame width = 27 3/8") N/C

BATTERIES

BAT24PR Gel Batteries - Group 24 (Each - On chair) \$428.00

DRIVE WHEEL

DW14PT 14" Cast Wheel with Pneumatic Tire STD
 DW14PI 14" Cast Wheel with Pneumatic Tire with Solid inserts \$249.00

FRONT CASTERS

FC8PT-STD 8" Pneumatic Casters STD
 FC8PI 8" Pneumatic Casters with Solid inserts \$249.00

SEAT TO FLOOR HEIGHT (Approximate measurement to the top of seat pan ± ¼")**Options for a Seat Only**

STFH1625 16 ¼" (Possible Seat angle of 0° only)
 STFH1725 17 ¼" (Possible Seat angle of 0° and 3°)
 STFH1825 18 ¼" (Possible Seat angle of 0° and 3°)
 STFH1925 19 ¼" (Possible Seat angle of 0°, 3° and 6°)
 STFH2025 20 ¼" (Possible Seat angle of 0°, 3°, 6°, 9°)

Options for the 45° HD Power Tilt

STFH17 17" (Possible Seat angle of 0° only)
 STFH175 17 ½" (Possible Seat angle of 0° only)
 STFH18 18" (Possible Seat angle of 0° only)
 STFH185 18 ½" (Possible Seat angle of 0° only)
 STFH19 19" (Possible Seat angle of 0° only)

Options for the 11" HD Power Elevating Seat Lift Only

STFH175 17 ½" (Possible Seat angle of 0° only)

Options for a Combo - Power Tilt with an HD Power Elevating Seat

STFH1925 19 ¼" (Possible Seat angle of 0° only)

SEAT TO FLOOR HEIGHT SPECIFICATION (If not standard)

SOMHSS Specific seat to floor height requirement other than standard. RFQ
 Custom upcharges may apply. Required at: _____" (± ¼" Range)

PRE-SET SEAT FRAME ANGLE

SFE-0 0° Pre-set seat frame seat angle (**Must be selected when using power options**) STD
 SFE-3 3° Pre-set seat frame seat angle (Only available with a seat configuration without power options) N/C
 SFE-6 6° Pre-set seat frame seat angle (Only available with a seat configuration without power options) N/C
 SFE-9 9° Pre-set seat frame seat angle (Only available with a seat configuration without power options) N/C

LIGHTING PACKAGE

EOLP 4 LED Light Package (Includes front beams, side flashers, rear tail light with flashers & brake light) \$1,098.00
 Note 4 Requires the upgraded Joystick with Light Operation or the OMNI (ECVJ6KL, ECELJL, ECRJCL, ECRJCTPL or EDIB)

NON-EXPANDABLE ELECTRONICS - PG DRIVES VR2 CONTROLLER AND JOYSTICKS

<input type="checkbox"/>	ECV902A	VR2 90 Amp Drive Controller	STD
<input type="checkbox"/>	ECVJ4K	5 VR2 Standard 4-Key Joystick	STD
<input type="checkbox"/>	ECVJ6K	5 VR2 6-Key Joystick (Allows for a single power option to operate through this joystick)	\$103
<input type="checkbox"/>	ECVJ9KL	5 VR2 9-Key Joystick with Lighting Operation (Allows for a single power function to operate through this joystick)	\$162.00

Note 5 Requires the VR2 90A Drive Control (ECV902A).
VR2 Electronics can only accommodate a single power seat function. The only choices are Power Tilt or Elevating seat

EXPANDABLE ELECTRONICS - PG DRIVES RNET CONTROLLERS

<input type="checkbox"/>	ECEL90	R-net EL 90 Amp Expandable Controller	\$546.00
<input type="checkbox"/>	ECR120	R-net PM 120 Amp Expandable Controller	\$1,098.00
<input type="checkbox"/>	ECRNH	6 Harness for Expandable Controller	N/C

Note 6 Must be ordered with the expandable controller ECEL90 or ECR120 to interface with Joysticks or OMNI screen

EXPANDABLE ELECTRONICS - PG DRIVES RNET JOYSTICKS

<input type="checkbox"/>	ECO	Omit Joystick (For use with an Omni Display)	N/C
<input type="checkbox"/>	ECELJ	8 R-net (LED) Joystick	STD
<input type="checkbox"/>	ECRJC	8 R-net Joystick with Color screen	\$162.00

JOYSTICKS WITH LIGHTING OPERATION

<input type="checkbox"/>	ECELJL	8 R-net (LED) Joystick with Lighting Operation	\$157.00
<input type="checkbox"/>	ECRJCL	8 R-net Joystick with Color screen and Lighting Operation	\$271.00
<input type="checkbox"/>	ECRJCTPL	8 R-net Joystick with Color screen w/ On & Off / Mode switch & speed pot and Lighting Operation	\$325.00

ADVANCED JOYSTICK

<input type="checkbox"/>	ECRJC35 New	8 R-net Advanced Joystick (3.5" Color screen). Lighting Operation, paddle switches, shortcut keys and stereo jacks	\$1,167.00
--------------------------	-------------	--	------------

Note 8 This option requires the expandable controller ECEL90 or ECR120

JOYSTICK MOUNT

<input type="checkbox"/>	JSMF	Fixed Joystick Mount	N/C
<input type="checkbox"/>	JSMSA	Swing Away Joystick Mount - Height adjustable, Angle adjustable and Multi-Axis	\$400.00
		<input type="checkbox"/> Set-up with standard installation to swing-away sideways <input type="checkbox"/> Set-up at 35° to swing-away sideways and downwards	
<input type="checkbox"/>	JSMCM	Midline Mounted Swing Away Joystick Mount (E1028)	\$400.00

JOYSTICK MOUNT SUPPORT

<input type="checkbox"/>	JSMFH	Joystick mount support Bracket	STD
--------------------------	-------	--------------------------------	-----

JOYSTICK HAND MOUNT

<input type="checkbox"/>	HML	Install Joystick on Left hand side	<input type="checkbox"/>	HMR	Install Joystick on Right hand side
--------------------------	-----	------------------------------------	--------------------------	-----	-------------------------------------

JOYSTICK HANDLE

<input type="checkbox"/>	JSH	Standard Joystick Gimble	STD
<input type="checkbox"/>	JBAT	Stick Handle 3" in Length	\$93.00
<input type="checkbox"/>	JBH	Ball Handle	\$98.00
<input type="checkbox"/>	JFBH	Foam Ball Handle	\$98.00
<input type="checkbox"/>	JTH	T-Handle	\$98.00
<input type="checkbox"/>	JGP	Goalpost	\$98.00

ATTENDANT DRIVE CONTROL OPTION

<input type="checkbox"/>	EADC	Attendant Hand Control for R-net EL & PM Controls	\$1,035.00
--------------------------	------	---	------------

ATTENDANT JOYSTICK HAND MOUNT

<input type="checkbox"/>	EADCL	Install Attendant Joystick on Rear Left hand side	<input type="checkbox"/>	EADCR	Install Attendant Joystick on Rear Right hand side
--------------------------	-------	---	--------------------------	-------	--

ENHANCED DISPLAY & DIGITAL INTERFACE BOX (Included with specialty controls)

<input type="checkbox"/>	EDIB	8 Omni (Enhanced Display & Digital Interface Box with Infrared)	\$2,094.00
--------------------------	------	---	------------

Note 8 Allows the switch operation of alternative devices through driver controls.
Requires one of the R-net fully Expandable Controllers (ECEL90 or ECR120)

ENHANCED DISPLAY MOUNT

<input type="checkbox"/>	EDSAM	User Swing Away Mount	<input type="checkbox"/>	EDML (Install Display on LH side)	<input type="checkbox"/>	EDMR (Install Display on RH side)	N/C
--------------------------	-------	-----------------------	--------------------------	-----------------------------------	--------------------------	-----------------------------------	-----

ENVIRONMENTAL CONTROL MODULE

<input type="checkbox"/>	ECU	9	Environment control module - (Output Module)	\$1,098.00
Note 9 Requires one of the fully expandable R-net controllers (ECEL90 or ECR120)				

BLUETOOTH MOUSE EMULATOR (Includes USB Dongle)

<input type="checkbox"/>	EDBT	10	Bluetooth Mouse Emulator (USB Dongle included)	\$1,098.00
<input type="checkbox"/> Windows and MAC OS <input type="checkbox"/> Android Devices <input type="checkbox"/> Iphone, Ipad and Ipod				
Note 10 Requires one of the fully expandable R-net controllers (ECEL90 or ECR120)				
Note 10 Not limited to only one selection. Any two can be installed on the same chair, but only one can be used at a time.				

POWER OPTIONS & SEATING

<input type="checkbox"/>	SAABAHD	14	Bariatric Amyseat - Seat Frame only	\$1,396.00
<input type="checkbox"/>	SAPT XVHD	11	Bariatric CG Power Tilt	\$8,825.00
<input type="checkbox"/>	SAPR XVHD	11, 13	Bariatric Power Recline with Integrated Shear Reduction (ISR)	\$9,598.00
<input type="checkbox"/>	SAPTR XVHD	11, 13	Bariatric CG Power Tilt and Power Recline with Integrated Shear Reduction (ISR)	\$14,343.00
<input type="checkbox"/>	SAPESXHD	11, 12	11" Power Elevating Seat Module for transfer assist	\$5,513.00
<input type="checkbox"/>	ADDPESXHD	12	Module - Add a Heavy-Duty 11" Power Elevating Seat to either SAPTXVHD, SAPRXVHD or SAPTRXVHD	\$4,410.00

Note 11 These options include the fully adjustable Amyseat (SAABAHD).
Note 12 The maximum user weight capacity with this option is 400 lbs
Note 13 It is recommended to order the Contour Back Support option with the power recline (STDSB)

ELECTRONICS - SINGLE OR DUAL (1 OR 2 ACTUATORS) POWER SEAT FUNCTION CONTROLLERS

<input type="checkbox"/>	EPSV	14, 15	Tilt or Elevating seat activation through 6-key VR2 Hand Control (ECVJ6K or ECVJ9KL)	N/C
<input type="checkbox"/>	EPSEL	14, 15	Tilt or Elevating seat activation through RNET 90 Amp Expandable Controller (ECEL90)	N/C
<input type="checkbox"/>	SAC		Single (1 Actuator) Power Seat Function Controller to Interface through R-net Electronics	N/C
<input type="checkbox"/>	DAC		Dual (2 Actuators) Power Seat Function Controller to Interface through R-net Electronics	\$546.00
With the SAC or DAC, please select one activation and position from the list below				
<input type="checkbox"/>	SWTJSS	16, 17	Single actuator activation through R-net Joystick or Specialty Controls	\$1,098.00
<input type="checkbox"/>	SWTJSD	16, 18	Dual actuator activation through R-net Joystick or Specialty Controls	\$1,098.00
<input type="checkbox"/>	SWDUALACT		Dual Activation (Allows to operate all seat functions with external switches AND through the joystick)	\$340
<input type="checkbox"/>	SWTOGGLE	17 or 18	Single Toggle Switch (Single Mode: DAC - Dual Mode: SAC) <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation	N/C
<input type="checkbox"/>	SWPUSH1	17	Single Push Button Switch (Single Mode: SAC) <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation	N/C
<input type="checkbox"/>	SWPUSH2	17 or 18	Dual Push Button Switch (Single Mode: DAC - Dual Mode: SAC) <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation	N/C

Note 14 Allows for the activation of a single power option through the VR2 or R-net 90 Amp Controllers
Only the Power Tilt or Elevating seat options are available with this feature.
Note 15 These options do not require the SAC controller as they connect via the chair's power module directly.
Note 16 Requires the R-net 90 Amp (ECEL90) or the 120 Amp (ECR120) Fully Expandable Controllers
Note 17 Requires the Single (1 Actuator) Power Seat Function Controller (SAC)
Note 18 Requires the Dual (2 Actuators) Power Seat Function Controller (DAC)

ELECTRONICS - MULTIPLE (3 to 6 ACTUATORS) POWER SEAT FUNCTION CONTROLLER

<input type="checkbox"/>	MAC		Multiple (3 to 6 Actuators) Power Seat Function Controller to Interface through R-net Electronics	\$2,424.00
With the MAC, please select one activation and position from the list below				
<input type="checkbox"/>	SWTJSM	19, 20	Multiple Actuator activation through R-net Joystick or Specialty Controls	N/C
<input type="checkbox"/>	SWDUALACT		Dual Activation (Allows to operate all seat functions with external switches and through the joystick)	\$340
<input type="checkbox"/>	SWMULTI	19	4-Way Toggle Switch (Up to 4 Activations) - <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation	N/C
<input type="checkbox"/>	SWMULTIT	19	Triple Toggle Switch (Up to 5 Activations) - <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation	N/C
<input type="checkbox"/>	SWPUSH4	19	4-Way Push Button Switch - <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation	N/C
<input type="checkbox"/>	SWPUSH5	19	5-Way Push Button Switch - <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation	N/C

Note 19 Requires the Multiple (3 to 6 Actuators) Power Seat Function Controller (MAC)
Note 20 Requires the R-net 90 Amp (ECEL90) or 120 Amp (ECR120) Fully Expandable Controller, and the multiple Actuator Controller (MAC)

SWITCH ACTIVATION HAND MOUNT

Switch Installation:	<input type="checkbox"/> TSML (Install Switch on LH side)	<input type="checkbox"/> TSMR (Install Switch on RH side)	N/C
----------------------	---	---	-----

PRE-SET SEAT ANGLE (For Non-Reclining Backs)

SOBAA Please Select one: 82° 90° 98° 106°

N/C

STANDARD ADJUSTABLE SEAT WIDTH RANGE BETWEEN 19" and 28"

STD

- SOWN19-24-22 (Width adjustment from 19" to 24" with a maximum Depth of 22")
 SOWN24-28-24 (Width adjustment from 24" to 28" with a maximum Depth of 24")

- | | |
|--|--|
| <input type="checkbox"/> 19" (SOWN19L) | <input type="checkbox"/> 24" (SOWN24HDL) |
| <input type="checkbox"/> 20" (SOWN20L) | <input type="checkbox"/> 25" (SOWN25HDL) |
| <input type="checkbox"/> 21" (SOWN21) | <input type="checkbox"/> 26" (SOWN26HDL) |
| <input type="checkbox"/> 22" (SOWN22) | <input type="checkbox"/> 27" (SOWN27HDL) |
| <input type="checkbox"/> 23" (SOWN23) | <input type="checkbox"/> 28" (SOWN28HDL) |
| <input type="checkbox"/> 24" (SOWN24) | |

ADJUSTABLE SEAT DEPTH RANGE BETWEEN 17" and 24"

- | | | | |
|---------------------------------------|---------------------------------------|---------------------------------------|---------------------------------------|
| <input type="checkbox"/> 17" (SODN17) | <input type="checkbox"/> 18" (SODN18) | <input type="checkbox"/> 19" (SODN19) | <input type="checkbox"/> 20" (SODN20) |
| <input type="checkbox"/> 21" (SODN21) | <input type="checkbox"/> 22" (SODN22) | <input type="checkbox"/> 23" (SODN23) | <input type="checkbox"/> 24" (SODN24) |

BACK CANE HEIGHT (The back cane height measurement is taken from the seat pan to the top of the push handles)

- | | | |
|--|--|--|
| <input type="checkbox"/> 16" (SOBCH16) | <input type="checkbox"/> 20" (SOBCH20) | <input type="checkbox"/> 24" (SOBCH24) |
| <input type="checkbox"/> 18" (SOBCH18) | <input type="checkbox"/> 22" (SOBCH22) ²¹ | <input type="checkbox"/> 26" (SOBCH26) |

Note 21 22" is the minimum back cane height required with a power or manual recline.

BACK SUPPORT

- | | | |
|---------------------------------------|--|---------|
| <input type="checkbox"/> STDSB | ²² Contour back support with a maximum cushion height of 20" (19" and 20" Wide) | \$606 |
| <input type="checkbox"/> STDSB/CUSTOM | Contour back support (Non-Standard Sizes up to 28" Wide) Specify: Width: Height: | \$1,320 |

Note 22 It is recommended to use the Contour Back Support option with the power or manual recline.

24 BACK SUPPORT TOTAL HEIGHT (MEASURED FROM THE SEAT PAN TO THE TOP OF THE BACK)

- | | | | |
|---------------------------------------|---|---------------------------------------|---------------------------------------|
| <input type="checkbox"/> 16" (SOBH16) | <input type="checkbox"/> 20" (SOBH20) | <input type="checkbox"/> 24" (SOBH24) | <input type="checkbox"/> 28" (SOBH28) |
| <input type="checkbox"/> 18" (SOBH18) | <input type="checkbox"/> 22" (SOBH22) ²³ | <input type="checkbox"/> 26" (SOBH26) | |

Note 23 With a power or manual recline, the minimum height is 22".

Note 24 The space between the bottom of the back and the seat pan is 4" with a non-reclining back and 6" with a reclining back

SEAT CUSHION

- | | | |
|---------------------------------------|--|---------|
| <input type="checkbox"/> SCPSC | Contour Seat Cushion with a maximum seat depth of 22" (19" and 20" Wide) | \$606 |
| <input type="checkbox"/> SCPSC/CUSTOM | Contour Seat Cushion (Non-Standard Sizes) Specify: Width: Depth: | \$1,320 |

CANTILEVER - HEIGHT ADJUSTABLE ARMRESTS (For Non-Reclining Backs - Choose 2)

- | | | |
|-----------------------------------|---|-----|
| <input type="checkbox"/> SOAFBFLH | Flip-Back Cantilever - Full Length Left Side (Height Adjustable from 7" to 15" in 1" increments) | STD |
| <input type="checkbox"/> SOAFBFRH | Flip-Back Cantilever - Full Length Right Side (Height Adjustable from 7" to 15" in 1" increments) | STD |
| <input type="checkbox"/> SOAFBDLH | Flip-Back Cantilever - Desk Length Left Side (Height Adjustable from 7" to 15" in 1" increments) | N/C |
| <input type="checkbox"/> SOAFBDRH | Flip-Back Cantilever - Desk Length Right Side (Height Adjustable from 7" to 15" in 1" increments) | N/C |

DUAL POST - HEIGHT ADJUSTABLE ARMRESTS - INCLUDES TRANSFER GRIPS (For Non-Reclining Backs - Choose 2)

- | | | |
|-----------------------------------|--|----------|
| <input type="checkbox"/> SOADPFLH | Dual Post Height Adjustable - Full Length Left Side (Height Adjustable from 9 ½" - 13 ½") | \$195.00 |
| <input type="checkbox"/> SOADPFRH | Dual Post Height Adjustable - Full Length Right Side (Height Adjustable from 9 ½" - 13 ½") | \$195.00 |
| <input type="checkbox"/> SOADPDLH | Dual Post Height Adjustable - Desk Length Left Side (Height Adjustable from 9 ½" - 13 ½") | \$195.00 |
| <input type="checkbox"/> SOADPDRH | Dual Post Height Adjustable - Desk Length Right Side (Height Adjustable from 9 ½" - 13 ½") | \$195.00 |

RECLINING FLIP-UP, HEIGHT ADJUSTABLE ARMRESTS (For Reclining Backs - Choose 2)

- | | | |
|-----------------------------------|---|-----|
| <input type="checkbox"/> SOANIFLH | Reclining armrest - Full Length Left Side (Height Adjustable from 10 ½" to 14" in ½" increments) | N/C |
| <input type="checkbox"/> SOANIFRH | Reclining armrest - Full Length Right Side (Height Adjustable from 10 ½" to 14" in ½" increments) | N/C |
| <input type="checkbox"/> SOANIDLH | Reclining armrest - Desk Length Left Side (Height Adjustable from 10 ½" to 14" in ½" increments) | N/C |
| <input type="checkbox"/> SOANIDRH | Reclining armrest - Desk Length Right Side (Height Adjustable from 10 ½" to 14" in ½" increments) | N/C |

ARMREST HEIGHT ADJUSTMENT AT FACTORY

Please set the Armrest Height at: _____ " (The Armrest Height is taken from the seat pan to the top of the armpad. Not required for Dual Post)

ARMPADS

<input type="checkbox"/>	SOFAPD	Desk Length Standard (10")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	STD
<input type="checkbox"/>	SOFAPF	Full Length Standard (14")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	STD
<input type="checkbox"/>	SOWAPD	34 Desk Length Waterfall (10")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$33.00
<input type="checkbox"/>	SOWAPF	34 Full Length Waterfall (14")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$33.00

Hybrid Gel and Foam Type Selection (The Gel always stays soft)

<input type="checkbox"/>	SOGELD <i>new</i>	Desk Length Hybrid Gel and Foam Flat armpad (10" x 2.5")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$33
<input type="checkbox"/>	SOGELF	Full Length Hybrid Gel and Foam Flat armpad (14" x 2.5")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$33
<input type="checkbox"/>	SOWGELD <i>new</i>	34 Desk Length Hybrid Gel and Foam Waterfall armpad (10")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$50
<input type="checkbox"/>	SOWGELF <i>new</i>	34 Full Length Hybrid Gel and Foam Waterfall armpad (14")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$50

Specialty Type Selection

<input type="checkbox"/>	SOCAR	35 Arm Trough with Hand Support (18" x 4 1/2")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$216.00
--------------------------	-------	--	---	--------------	---	--	----------

Note 27 The Waterfall armpads may interfere with flipping back the armrests when combined with Lateral Trunk supports

Note 28 Arm Trough Armpads may interfere with Armrests Flipping back.

FOOTRESTS

<input type="checkbox"/>	SOAAF**-STD	29 Heavy Duty Swing-Away footrests ----- Set at: <input type="checkbox"/> 60° <input type="checkbox"/> 70°	\$379.00
		Select Footrest Range: <input type="checkbox"/> SOFLHD12 (12 1/2" - 15 1/2") <input type="checkbox"/> SOFLHD14 (14 1/2" - 19 1/2")	
<input type="checkbox"/>	ADDPLHD	26, 29 Heavy Duty Power Elevating & Articulating Legrests Select Function: <input type="checkbox"/> Individual ²⁷ <input type="checkbox"/> Combined ²⁸	\$3,084.00
		Select Footrest Range: <input type="checkbox"/> SOFLPW12 (12 1/2" - 15 1/2") <input type="checkbox"/> SOFLPW14 (14 1/2" - 19 1/2")	

29 Please set the Footrest Length at: _____ " (The footrest length is taken from the seat pan to the heel)

Note 26 This option includes calf pads

Note 27 Individual legs count as two power functions.

Note 28 Combined legs count as one power function.

Note 29 The Footrest Length is taken from the top of the seat pan to the top of the footplate.

FOOTPLATES & ACCESSORIES

<input type="checkbox"/>	SOSFAFHD	Flip-up Footplates	STD
<input type="checkbox"/>	SOAFPHD	Angle Adjustable & Flip-up Footplates - Each	\$141.00
<input type="checkbox"/>	SOAHL	Heel Loops - Each	\$43.00
<input type="checkbox"/>	SOACPFH	Calf Pads for Adult Heavy-Duty Swing-Away footrests (SOAAF) with adjustable hardware - Each	\$141.00

CENTER MOUNTS

<input type="checkbox"/>	SOAAF90HD	33, 34, 36 90° Heavy Duty Fixed Center Mount Interface Bracket	\$249.00
		Select Foot Platform Range: <input type="checkbox"/> SOFBP7 (7 1/2" - 12 1/2") <input type="checkbox"/> SOFBA12 (12 1/2" - 16 1/2") <input type="checkbox"/> SOFBA16 (16 1/2" - 20 1/2")	
<input type="checkbox"/>	ADDPFH	30, 32, 33, 34, 35, 36 90° Heavy Duty Power Center Mount Elevating & Articulating Foot Platform - (Uses 1 center mounted actuator)	\$3,084.00
		<input type="checkbox"/> SOFBPW12 (12 1/2" - 16 1/2") <input type="checkbox"/> SOFBPW16 (16 1/2" - 20 1/2")	
<input type="checkbox"/>	ADDPFVHD	31, 32, 33, 34 90° Very Heavy Duty Power Center Mount Elevating & Articulating Foot Platform - (Uses 2 side mounted actuators)	\$3,713.00

36 Please set the Foot Platform Length at: _____ "

Note 30 ADDPFH counts as one power function.

Note 31 ADDPFVHD counts as two power functions. Please select (SOFUFP) on the next page

Note 32 This option includes calf pads

Note 33 Links and attaches the one piece center mount foot platform (SOFUFP) or the split footplates (SOAFPFH) to the seat

Note 34 Please Select the one piece center mount foot platform (SOFUFP) or the split footplates (SOAFPFH) if compatible (See footnotes)

Note 35 The Standard user weight capacity is 450 Lbs.

Note 36 The Foot Platform Length is taken from the top of the seat pan to the top of the foot platform or split foot plates.

FOOT PLATFORMS & ACCESSORIES

<input type="checkbox"/>	SOFUFP	Flip-up one piece Foot Platform	N/C
<input type="checkbox"/>	SOAFPFH	Angle Adjustable & Flip-up Split Footplates for center mount - Each	\$141.00
<input type="checkbox"/>	SOACPCM	Calf Pads for 90° Fixed Center Mount (SOAAF90) with adjustable hardware - Each	\$141.00

ORDER INFORMATION: Customer Service: 800-263-3390 Fax: 800-561-5834 www.sunrisemedical.ca

Specifications are subject to change without notice. Copyright © 2017 Sunrise Medical Canada Inc. QUICKIE XCEL 2, Rev 1

SEAT POSITIONING

<input type="checkbox"/>	SOAUTO80	80" Seat Positioning Strap	\$82.00
--------------------------	----------	----------------------------	---------

HEADREST

<input type="checkbox"/>	SONUAHR	Using headrest mounting hardware other than SOSAH or SOMAH located below. Specify: _____	
<input type="checkbox"/>	SOFHRS	Small Curved Headrest Pad 4" x 9"	\$222.00
<input type="checkbox"/>	SOFHRM	Medium Curved Headrest Pad 5" x 10"	\$222.00
<input type="checkbox"/>	SOFHRL	Large Curved Headrest Pad 6" x 11"	\$222.00
<input type="checkbox"/>	SOMAH-A	Multi-axis Removable Mounting Hardware	\$178.00
<input type="checkbox"/>	SOHRBC	Bracket to Mount on the Contour Back Support (Without Recline)	\$108
<input type="checkbox"/>	SOHRSB	Bracket to Mount on the Seat Frame Crossbar (Without Recline)	\$108
<input type="checkbox"/>	SOHRPR	Bracket to Mount on the Seat Frame (For Reclining backs, but excludes the Corbeau Back)	\$108

* Please select one of the headrest options with the mounting hardware and specify where the headrest will be installed (The frame or back cushion)

LATERAL TRUNK SUPPORTS

<input type="checkbox"/>	SOSLS3	Adjustable Laterals with Small Curved Pads (Each) - 3" x 5"	\$221.00
<input type="checkbox"/>	SOSLM4	Adjustable Laterals with Medium Curved Pads (Each) - 4" x 6"	\$221.00
<input type="checkbox"/>	SOSLL5	Adjustable Laterals with Large Curved Pads (Each) - 5" x 7"	\$221.00
<input type="checkbox"/>	SOSLHM	Adjustable Swing-Away Hardware Mounts - Each	\$326.00

** Please select one of the lateral pads and the adjustable & removable hardware.

HIP & THIGH SUPPORTS

<input type="checkbox"/>	SOHTSS10	Adjustable Hip & Thigh Support Pads (Each) - Small 10" x 4"	\$221.00
<input type="checkbox"/>	SOHTSM12	Adjustable Hip & Thigh Support Pads (Each) - Medium 12" x 5"	\$221.00
<input type="checkbox"/>	SOHTSL14	Adjustable Hip & Thigh Support Pads (Each) - Large 14" x 6"	\$221.00
<input type="checkbox"/>	SOHTSARH	Adjustable and Removable Hip & Thigh Support Hardware - Each	\$471.00

*** Please select one of the hip & thigh support pads and the adjustable & removable hardware.

ACCESSORIES

<input type="checkbox"/>	SOCH	Cup Holder (Self-Leveling)	\$98.00
<input type="checkbox"/>	SOPHP	Ram Mounts X-Grip Cell Phone Holder	\$168.00
<input type="checkbox"/>	SOTHS	Ram Mounts X-Grip Small Tablet Holder	\$190.00
<input type="checkbox"/>	SOTHL	Ram Mounts X-Grip Large Tablet Holder with Swivel	\$286.00

REFLECTOR KIT

<input type="checkbox"/>	REFLECT	Red & White Reflector Kit Stickers (They will not be installed - Included in a bag with the owner's manual)	N/C
--------------------------	---------	---	-----

SWITCHES AND OPTIONS

<input type="checkbox"/>	ASCPCP	PGDT RNET Dongle and PC programmer kit	\$633.00
<input type="checkbox"/>	ASCPVR2	PGDT VR2 Dongle and PC programmer kit	\$633.00
<input type="checkbox"/>	SSBS	Buddy Button Switch	\$119.00
<input type="checkbox"/>	SSES	Egg Switch-Black	\$119.00
<input type="checkbox"/>	SSML	Micro Light Switch	\$162.00

SPECIALTY DRIVER CONTROL OPTIONS

Therafin Sip-N-Puff Kit

<input type="checkbox"/>	SCSNPK	37, 38	Therafin Sip-N-Puff Kit		\$1,651
<input type="checkbox"/>	SCSNPBTK	38	Therafin Sip-N-Puff Breath Tube Kit		\$536
<input type="checkbox"/>	SCSRMSNP	38	Stealth Rock-n-Lock Swingaway Arm for Sip-N-Puff	<input type="checkbox"/> SWL (Swings on Left) <input type="checkbox"/> SWR (Swings on Right)	\$298
<input type="checkbox"/>	SCSRM	38	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/> SWL (Swings Left) <input type="checkbox"/> SWR (Swings Right)	\$536
<input type="checkbox"/>	SSES-STD	38, 40	Egg Switch for Changing Modes - Black	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$119
<input type="checkbox"/>	SSBS-STD	38	Buddy Button Switch (Stop Switch) - Red		\$119
<input type="checkbox"/>	SOMAH	38	Multi-axis Headrest Removable Mounting Hardware		\$178
<input type="checkbox"/>	SOFHRS		Small Curved Headrest Pad 4" x 9"		\$222
<input type="checkbox"/>	SOFHRM	39	Medium Curved Headrest Pad 5" x 10"		\$222
<input type="checkbox"/>	SOFHRL		Large Curved Headrest Pad 6" x 11"		\$222

Switch-It Proximity Head Array

<input type="checkbox"/>	SCSIHAA	37, 38	Adult - 3 Switch Proximity Head Array with 1/8" Mono Jack for mode/reverse		\$5,513
<input type="checkbox"/>	SCSRM	38	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/> SWL (Swings Left) <input type="checkbox"/> SWR (Swings Right)	\$536
<input type="checkbox"/>	SSES-STD	38, 40	Egg Switch for Changing Modes - Black	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$119
<input type="checkbox"/>	SSBS-STD	38	Buddy Button Switch (Stop Switch) - Red		\$119
<input type="checkbox"/>	SOMAH	38	Multi-axis Headrest Removable Mounting Hardware		\$178

Switch-It Proximity Head Array & Sip-N-Puff Combo

<input type="checkbox"/>	SCSICMBOA	37, 38	Adult - 4 Switch Combo Head Array with two proximity sensors for Left & Right and Low pressure Sip-N-Puff for forward and reverse and 1/8" jack for mode		\$5,513
<input type="checkbox"/>	SCSNPBTK	38	Therafin Sip-N-Puff Breath Tube Kit		\$536
<input type="checkbox"/>	SCSRMSNP	38	Stealth Rock-n-Lock Swingaway Arm for Sip-N-Puff	<input type="checkbox"/> SWL (Swings on Left) <input type="checkbox"/> SWR (Swings on Right)	\$298
<input type="checkbox"/>	SCSRM	38	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/> SWL (Swings Left) <input type="checkbox"/> SWR (Swings Right)	\$536
<input type="checkbox"/>	SSES-STD	38, 40	Egg Switch for Changing Modes - Black	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$119
<input type="checkbox"/>	SSBS-STD	38	Buddy Button Switch (Stop Switch) - Red		\$119
<input type="checkbox"/>	SOMAH	38	Multi-axis Headrest Removable Mounting Hardware		\$178

Switch-It MicroPilot Light Touch Proportional Mini joystick for chin / lip control

<input type="checkbox"/>	SCSIMP	37, 38	MicroPilot Light touch mini Proportional Joystick		\$4,323
<input type="checkbox"/>	SCSRMMP	38	Stealth Rock-n-Lock Swingaway Mount for Chin M-P	<input type="checkbox"/> SWL (Swings Left) <input type="checkbox"/> SWR (Swings Right)	\$536
<input type="checkbox"/>	SCSRM	38	Stealth Rock-n-Lock Swingaway Mount for Switch	<input type="checkbox"/> SWL (Swings Left) <input type="checkbox"/> SWR (Swings Right)	\$536
<input type="checkbox"/>	SSES	38, 40	Egg Switch-Black	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$98
<input type="checkbox"/>	JCC		Chin Cup		\$54
<input type="checkbox"/>	SOMAH	38	Multi-axis Headrest Removable Mounting Hardware		\$178
<input type="checkbox"/>	SOFHRS		Small Curved Headrest Pad 4" x 9"		\$222
<input type="checkbox"/>	SOFHRM	39	Medium Curved Headrest Pad 5" x 10"		\$222
<input type="checkbox"/>	SOFHRL		Large Curved Headrest Pad 6" x 11"		\$222

Proportional Compact or Mini joystick for Midline Mount

<input type="checkbox"/>	SCPGMJS	37, 38	PG Drives JC200 Proportional Compact Remote Joystick for Midline Mount		\$1,028
<input type="checkbox"/>	SCGMJSM	41	Gatlin Midline Swing-away Mount for JC200	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$422
<input type="checkbox"/>	SCSIMP	37, 38	Switch-It MicroPilot Light touch mini Proportional Joystick		\$4,323
<input type="checkbox"/>	SCGMPM	42	Gatlin Midline Swing-away Mount for MicroPilot	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$460
<input type="checkbox"/>	SCBTM	43	Bullet Hand Tray for MicroPilot	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$487

Proportional Compact or Mini joystick for Armrest Mount

<input type="checkbox"/>	SCPGMJS	37, 38	PG Drives JC200 Proportional Compact Remote Joystick for Armrest Mount		\$1,028
<input type="checkbox"/>	SCARMJSM	41	Swing-away and Multi-Axis armrest mount	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$129
<input type="checkbox"/>	SCSIMP	37, 38	Switch-It MicroPilot Light touch mini Proportional Joystick		\$4,323
<input type="checkbox"/>	SCARMMPM	42	Swing-away and Multi-Axis armrest mount	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$460
<input type="checkbox"/>	SCBTM	43	Bullet Hand Tray for MicroPilot	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$487

Proportional RIM Head Control Package

<input type="checkbox"/>	SCPHC	37, 38	Proportional RIM Head Control Package		\$2,700
<input type="checkbox"/>	SCSRM	38	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/> SWL (Swings Left) <input type="checkbox"/> SWR (Swings Right)	\$536
<input type="checkbox"/>	SSES-STD	38, 40	Egg Switch for Changing Modes - Black	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$119
<input type="checkbox"/>	SOMAH	38	Multi-axis Headrest Removable Mounting Hardware		\$178

- Note 37** Requires the Omni Enhanced Display & Digital Interface Box (EDIB)
- Note 38** This option must be selected for the Kit
- Note 39** Choose one (1) Headrest Pad
- Note 40** The Mode Switch Type can be changed from this selection on Page #6
- Note 41** Must be selected for use with the PG Drives JC200 Mini Joystick
- Note 42** Must be selected for use with the Switch-It MicroPilot
- Note 43** This option is only compatible with the Switch-It MicroPilot

