



April 1, 2017

ACCOUNT INFORMATION

| | |
|--|----------------------------------|
| Dealer Contact Name: _____ | Bill to (Name): _____ |
| Email Address: _____ | Bill to (Location): _____ |
| Phone Number: _____ | Ship to (Name): _____ |
| Fax Number: _____ | Address: _____ |
| User Name or Number: _____ | City: _____ |
| User Weight: _____ | Province: _____ |
| User Height: _____ | Postal Code: _____ |
| Seat Cushion Thickness: _____ | P.O. # or Quote: _____ |
| Knee to Heel Length: Left: _____ " Right: _____ " | |
| Seat to Top of Scapula: _____ | |

The order cannot be processed without the following:
 - A user name or number, the client's measurements and weight
 - Sunrise is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
 - Order confirmations will be provided by Amysystems summarizing the items selected. Please review this carefully.
 - Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

STANDARD OPTIONS

The Quickie Xperience Power Base, 6-wheel independent wheel suspension, an additional seat suspension system, tie-down hardware for strap systems, offboard charger, 6" front casters, and standard auto style buckle seat positioning strap, 5.9mph motors, Weight Capacity = 300lbs

QUICKIE XPERIENCE STD POWER WHEELCHAIR

- XPERSTD2 QUICKIE XPERIENCE STD - 2 Power Wheelchair \$8,053

USER WEIGHT PACKAGES

- WCA3 300 lb Maximum user weight capacity STD

| STANDARD FRAME COLORS | | | OPTIONNAL FRAME COLORS (May Increase lead times) | | | | |
|--------------------------|------------------|----------------------|--|--------------------------|----------|--|-------|
| <input type="checkbox"/> | FCBLK | Black | STD | <input type="checkbox"/> | FCTXG | Toxic Green | N/C |
| <input type="checkbox"/> | FCSBM | Matte Black | STD | <input type="checkbox"/> | FCBOR | Burnt Orange Metallic | N/C |
| <input type="checkbox"/> | FCGMC | Gun Metal (Charcoal) | STD | <input type="checkbox"/> | FCBAO | Black Opal Metallic | N/C |
| <input type="checkbox"/> | FCSVM | Silver Metallic | STD | <input type="checkbox"/> | FCOPM | Sparkle Orange | N/C |
| <input type="checkbox"/> | FCCBU | Candy Blue | STD | <input type="checkbox"/> | FCYLW | Yellow | N/C |
| <input type="checkbox"/> | FCCRE | Candy Red | STD | <input type="checkbox"/> | FCPKP | Pearl Pink | N/C |
| <input type="checkbox"/> | FCMPL | Midnight Purple | STD | <input type="checkbox"/> | FCHDD 1 | Desert Digital Camo Hydro Graphics | \$320 |
| <input type="checkbox"/> | FCRZB | Razzleberry | STD | <input type="checkbox"/> | FCHWL 1 | Woodland Camo Hydro Graphics | \$320 |
| <input type="checkbox"/> | FCMWT New | Matte White | STD | <input type="checkbox"/> | CUSTOM 1 | Carbon Fiber Hydro Graphics (Matte Black base) | \$320 |

Note 1 Includes matching accessories with the Hydro Graphics option. The matching accessories option is not required with these colors

MATCHING ACCESSORIES TO THE FRAME COLOR

- FC***-M Matching Accessories to the Frame Color (Wheels, front access door, headrest mount and back shell when selected) \$242

TRANSPORT OPTION

- TTD4 2 4-Point wheelchair transport brackets N/C

Note 2 4-Point wheelchair transport brackets have not yet received WC-19 crash-testing approval.
 Use these transport brackets only to secure an unoccupied wheelchair during transport.

MOTOR PACKAGE

- MP65 6.5MPH 4-Pole High Speed Motor Package STD

ASSISTIVE DRIVING PACKAGE

- INTELLIDR 3, 4 Intelli-Drive \$1,475

Note 3 Intelli-Drive technology provides increased tracking capabilities with enhanced traction and steering precision.
 This is ideal for those with marginal hand function or switch control.
 Note 4 Requires the R-net PM 120 Amp Expandable Controller (ECR120)

BATTERY CHARGER

- BCOB Offboard Battery Charger N/C
- BCLES8OPT Lester 8 amp Offboard Battery Charger \$361

BATTERY TRAY AND FRAME SIZE

| | | | |
|--------------------------|---------|---|-----|
| <input type="checkbox"/> | FS22B | Battery Tray for Group 22NF Batteries (Total frame width = 23 3/8") | STD |
| <input type="checkbox"/> | FS24B | Battery Tray for Group 24 Batteries (Total frame width = 25 3/8") | N/C |
| <input type="checkbox"/> | FS2424B | 5 CUSTOM Battery Tray for Group 24 Batteries with a Compact 24" frame width | N/C |

Note 5 The 90° Power Center Mount Elevating & Articulating Foot Platform (ADDPF) is not available with this option

BATTERIES

| | | | |
|--------------------------|-----------|--|-------|
| <input type="checkbox"/> | BAT22NFPR | Gel Batteries - Group 22NF (Each - On chair) | \$620 |
| <input type="checkbox"/> | BAT24PR | Gel Batteries - Group 24 (Each - On chair) | \$830 |

DRIVE WHEEL

| | | | |
|--------------------------|---------|---|-------|
| <input type="checkbox"/> | DW14PT | 14" Cast Wheel Pneumatic Tire | STD |
| <input type="checkbox"/> | DW14PI | 14" Cast Wheel Pneumatic with Flat-Free inserts | \$242 |
| <input type="checkbox"/> | DW144PT | 6 14" x 4" Extra Wide Cast Wheel Pneumatic Tire | \$263 |
| <input type="checkbox"/> | DW144PI | 6 14" x 4" Extra Wide Cast Wheel Pneumatic with Flat-Free inserts | \$263 |
| <input type="checkbox"/> | DW12PT | 7 12 1/2" Cast Wheel Pneumatic Tire | N/C |
| <input type="checkbox"/> | DW12PI | 7 12 1/2" Cast Wheel Pneumatic with Flat-Free inserts | \$242 |

The 14" x 4" Extra Wide Wheel increases the overall width of the Base by 2 1/2"

Note 6 The option to Match the Wheels to the Shroud Color (FC...-M) is not available with the 14" x 4" Extra Wide wheels

The 14" x 4" option is not available with VR2 electronics (ECV902A)

Note 7 12 1/2" Drive wheels reduce overall seat to floor height by 3/4", lowers ground clearance and speed, but increases motor torque

FRONT CASTERS

| | | | |
|--------------------------|--------|---|-----|
| <input type="checkbox"/> | FC6S | 6" Solid Casters (Not available with the 14" x 4" Extra Wide Cast Wheel options DW144PT or DW144PI) | STD |
| <input type="checkbox"/> | FC8PTM | 8" Pneumatic Casters | N/C |
| <input type="checkbox"/> | FC8PIM | 8" Pneumatic Casters with Flat-Free inserts | N/C |

Note 8 The 8" casters are NOT compatible with the 12 1/2" Drive wheel configuration (DW12PT and DW12PI)

SEAT TO FLOOR HEIGHT (Approximate measurement to the top of seat pan ± 1/4")

Options for a Power Tilt or a STD Power Elevating Seat

| Options for a Seat Only | | | |
|--------------------------|----------|---|--------------------------|
| <input type="checkbox"/> | 9 STF155 | 15 1/2" (Possible Seat angle of 0° only) | <input type="checkbox"/> |
| <input type="checkbox"/> | STFH1625 | 16 1/4" (Possible Seat angle of 0° only) | <input type="checkbox"/> |
| <input type="checkbox"/> | STFH1725 | 17 1/4" (Possible Seat angle of 0° and 3°) | <input type="checkbox"/> |
| <input type="checkbox"/> | STFH1825 | 18 1/4" (Possible Seat angle of 0° and 3°) | <input type="checkbox"/> |
| <input type="checkbox"/> | STFH1925 | 19 1/4" (Possible Seat angle of 0°, 3° and 6°) | <input type="checkbox"/> |
| <input type="checkbox"/> | STFH2025 | 20 1/4" (Possible Seat angle of 0°, 3°, 6°, 9°) | <input type="checkbox"/> |
| | | | <input type="checkbox"/> |
| | | | <input type="checkbox"/> |
| | | | <input type="checkbox"/> |
| | | | <input type="checkbox"/> |
| | | | <input type="checkbox"/> |
| | | | <input type="checkbox"/> |

Options for a Combo - Power Tilt with a STD Power Elevating Seat

| | | |
|--------------------------|----------|---|
| <input type="checkbox"/> | 9 STF175 | 17 1/2" (Possible Seat angle of 0° only) New |
| <input type="checkbox"/> | STFH1825 | 18 1/4" (Possible Seat angle of 0° only) New |

Note 9 Requires the 12 1/2" Drive wheels (DW12PT or DW12PI). This lowers ground clearance and speed, but increases motor torque.

SEAT TO FLOOR HEIGHT SPECIFICATION (If not standard)

| | | | |
|--------------------------|--------|--|-----|
| <input type="checkbox"/> | SOMHSS | Specific seat to floor height requirement other than standard. Custom upcharges may apply. Required at: _____" (± 1/4" Range) | RFQ |
|--------------------------|--------|--|-----|

PRE-SET SEAT FRAME ANGLE

| | | | |
|--------------------------|-------|---|-----|
| <input type="checkbox"/> | SFE-0 | 0° Pre-set seat frame seat angle (Must be selected when using power options) | STD |
| <input type="checkbox"/> | SFE-3 | 3° Pre-set seat frame seat angle (Only available with a seat configuration without power options) | N/C |
| <input type="checkbox"/> | SFE-6 | 6° Pre-set seat frame seat angle (Only available with a seat configuration without power options) | N/C |
| <input type="checkbox"/> | SFE-9 | 9° Pre-set seat frame seat angle (Only available with a seat configuration without power options) | N/C |

LIGHTING PACKAGE

| | | | |
|--------------------------|------|---|---------|
| <input type="checkbox"/> | EOLP | 10 LED Light Package (Includes front beams, side flashers, rear tail light with flashers & brake light) | \$1,066 |
|--------------------------|------|---|---------|

Note 10 Requires the upgraded Joystick with Light Operation or the OMNI (ECVJ6KL, ECELJL, ECRJCL, ECRJCTPL or EDIB)
The Light Package is not compatible with the EPSEL option.

NON-EXPANDABLE ELECTRONICS - PG DRIVES VR2 CONTROLLER AND JOYSTICKS

| | | | | | |
|--------------------------|---------|----------------|---|--|-------|
| <input type="checkbox"/> | ECV902A | | VR2 90 Amp Drive Controller | | STD |
| <input type="checkbox"/> | ECVJ4K | 11 | VR2 Standard 4-Key Joystick | | STD |
| <input type="checkbox"/> | ECVJ6K | 11 | VR2 6-Key Joystick (Allows for a single power option to operate through this joystick) | | N/C |
| <input type="checkbox"/> | ECVJ9KL | 11 | VR2 9-Key Joystick with Lighting Operation (Allows for a single power function to operate through this joystick) | | \$158 |
| | | Note 11 | Requires the VR2 90A Drive Control (ECV902A). VR2 Electronics can only accommodate a single power seat function. The only choices are Power Tilt or Elevating seat | | |

EXPANDABLE ELECTRONICS - PG DRIVES RNET CONTROLLERS

| | | | | | |
|--------------------------|--------|----------------|---|--|---------|
| <input type="checkbox"/> | ECEL90 | | R-net EL 90 Amp Expandable Controller | | \$530 |
| <input type="checkbox"/> | ECR120 | | R-net PM 120 Amp Expandable Controller | | \$1,066 |
| <input type="checkbox"/> | ECRNH | 12 | Harness for Expandable Controller | | N/C |
| | | Note 12 | Must be ordered with the expandable controller ECEL90 or ECR120 to interface with Joysticks or OMNI screen | | |

EXPANDABLE ELECTRONICS - PG DRIVES RNET JOYSTICKS

| | | | | | |
|--|------------|----------------|--|--|---------|
| <input type="checkbox"/> | ECO | | Omit Joystick (For use with an Omni Display) | | N/C |
| <input type="checkbox"/> | ECELJ | 13 | R-net (LED) Joystick | | STD |
| <input type="checkbox"/> | ECRJJC | \$13 | R-net Joystick with Color screen | | \$158 |
| JOYSTICKS WITH LIGHTING OPERATION | | | | | |
| <input type="checkbox"/> | ECELJL | \$13 | R-net (LED) Joystick with Lighting Operation | | \$158 |
| <input type="checkbox"/> | ECRJCL | \$13 | R-net Joystick with Color screen and Lighting Operation | | \$210 |
| <input type="checkbox"/> | ECRJCTPL | \$13 | R-net Joystick with Color screen w/ On & Off / Mode switch & speed pot and Lighting Operation | | \$315 |
| ADVANCED JOYSTICK | | | | | |
| <input type="checkbox"/> | ECRJ35 New | 13 | R-net Advanced Joystick (3.5" Color screen). Lighting Operation, paddle switches, shortcut keys and stereo jacks | | \$1,133 |
| | | Note 13 | This option requires the expandable controller ECEL90 or ECR120 | | |

JOYSTICK MOUNT

| | | | | | |
|--------------------------|-------|--|---|--|-------|
| <input type="checkbox"/> | JSMF | | Fixed Joystick Mount | | N/C |
| <input type="checkbox"/> | JSMSA | | Swing Away Joystick Mount - Height adjustable, Angle adjustable and Multi-Axis | | \$388 |
| | | | <input type="checkbox"/> Set-up with standard installation to swing-away sideways <input type="checkbox"/> Set-up at 35° to swing-away sideways and downwards | | |
| <input type="checkbox"/> | JSMCM | | Midline Mounted Swing Away Joystick Mount | | \$388 |

JOYSTICK MOUNT SUPPORT

| | | | | | |
|--------------------------|-------|--|--------------------------------|--|-----|
| <input type="checkbox"/> | JSMFH | | Joystick mount support Bracket | | STD |
|--------------------------|-------|--|--------------------------------|--|-----|

JOYSTICK HAND MOUNT

| | | | | | | | |
|--------------------------|-----|--|------------------------------------|--------------------------|-----|--|-------------------------------------|
| <input type="checkbox"/> | HML | | Install Joystick on Left hand side | <input type="checkbox"/> | HMR | | Install Joystick on Right hand side |
|--------------------------|-----|--|------------------------------------|--------------------------|-----|--|-------------------------------------|

JOYSTICK HANDLE

| | | | | | |
|--------------------------|------|--|---------------------------|--|------|
| <input type="checkbox"/> | JSH | | Standard Joystick Gimble | | STD |
| <input type="checkbox"/> | JBAT | | Stick Handle 3" in Length | | \$95 |
| <input type="checkbox"/> | JBH | | Ball Handle | | \$95 |
| <input type="checkbox"/> | JFBH | | Foam Ball Handle | | \$95 |
| <input type="checkbox"/> | JTH | | T-Handle | | \$95 |
| <input type="checkbox"/> | JGP | | Goalpost | | \$95 |

ATTENDANT DRIVE CONTROL OPTION

| | | | | | |
|--------------------------|------|--|---|--|---------|
| <input type="checkbox"/> | EADC | | Attendant Hand Control for R-net EL & PM Controls | | \$1,066 |
|--------------------------|------|--|---|--|---------|

ATTENDANT JOYSTICK HAND MOUNT

| | | | | | | | |
|--------------------------|-------|--|---|--------------------------|-------|--|--|
| <input type="checkbox"/> | EADCL | | Install Attendant Joystick on Rear Left hand side | <input type="checkbox"/> | EADCR | | Install Attendant Joystick on Rear Right hand side |
|--------------------------|-------|--|---|--------------------------|-------|--|--|

ENHANCED DISPLAY & DIGITAL INTERFACE BOX (Included with specialty controls)

| | | | | | |
|--------------------------|------|----------------|--|--|---------|
| <input type="checkbox"/> | EDIB | 14 | Omni (Enhanced Display & Digital Interface Box with Infrared) | | \$2,033 |
| | | Note 14 | Allows the switch operation of alternative devices through driver controls. Requires one of the R-net fully Expandable Controllers (ECEL90 or ECR120) | | |

ENHANCED DISPLAY MOUNT

| | | | | | | | | | |
|--------------------------|-------|--|-----------------------|--------------------------|-----------------------------------|--------------------------|-----------------------------------|--|-----|
| <input type="checkbox"/> | EDSAM | | User Swing Away Mount | <input type="checkbox"/> | EDML (Install Display on LH side) | <input type="checkbox"/> | EDMR (Install Display on RH side) | | N/C |
|--------------------------|-------|--|-----------------------|--------------------------|-----------------------------------|--------------------------|-----------------------------------|--|-----|

ENVIRONMENTAL CONTROL MODULE

| | | | | |
|--------------------------|-----|----|--|---------|
| <input type="checkbox"/> | ECU | 15 | Environment control module - (Output Module) | \$1,045 |
|--------------------------|-----|----|--|---------|

Note 15 Requires one of the fully expandable R-net controllers (ECEL90 or ECR120)

BLUETOOTH MOUSE EMULATOR (Includes USB Dongle)

| | | | | |
|--------------------------|------|----|--|---------|
| <input type="checkbox"/> | EDBT | 16 | Bluetooth Mouse Emulator (USB Dongle included) | \$1,066 |
| <input type="checkbox"/> | | | Windows and MAC OS | |
| <input type="checkbox"/> | | | Android Devices | |
| <input type="checkbox"/> | | | Iphone, Ipad and Ipod | |

Note 16 Requires one of the fully expandable R-net controllers (ECEL90 or ECR120)
 Not limited to only one selection. Any two can be installed on the same chair, but only one can be used at a time.

POWER OPTIONS & SEATING

| | | | | |
|--------------------------|---------|--------|---|----------|
| <input type="checkbox"/> | SAABA | | Amyseat Rehab Seat - Seat Frame only | \$830 |
| <input type="checkbox"/> | SAPTX | 17 | CG Power Tilt | \$6,961 |
| <input type="checkbox"/> | SAPRX | 17, 19 | Power Recline with Integrated Shear Reduction (ISR) | \$8,568 |
| <input type="checkbox"/> | SAPTPRX | 17, 19 | CG Power Tilt and Power Recline with Integrated Shear Reduction (ISR) | \$12,319 |
| <input type="checkbox"/> | SAPESX | 17 | 11" Power Elevating Seat Module for transfer assist | \$3,452 |
| <input type="checkbox"/> | ADDPESX | | Module - Add an 11" Power Elevating Seat to either SAPTX, SAPRX or SAPTPRX | \$3,209 |
| <input type="checkbox"/> | ADDMR | 18, 19 | Module - Manual Recline with Integrated Shear Reduction (ISR) | \$1,066 |

Note 17 These options include the fully adjustable Amyseat (SAABA).

Note 18 The maximum user weight capacity with this option is 250 lbs

Note 19 It is recommended to use the Contour Back Support option with the power recline (STDSB or STDSB/CUSTOM)

ELECTRONICS - SINGLE OR DUAL (1 OR 2 ACTUATORS) POWER SEAT FUNCTION CONTROLLERS

| | | | | |
|--------------------------|-------|--------|---|-------|
| <input type="checkbox"/> | EPSV | 20, 21 | Tilt or Elevating seat activation through 6-key VR2 Hand Control (ECVJ6K or ECVJ9KL) | N/C |
| <input type="checkbox"/> | EPSEL | 20, 21 | Tilt or Elevating seat activation through RNET 90 Amp Expandable Controller (ECEL90) | N/C |
| <input type="checkbox"/> | SAC | | Single (1 Actuator) Power Seat Function Controller to Interface through R-net Electronics | N/C |
| <input type="checkbox"/> | DAC | | Dual (2 Actuators) Power Seat Function Controller to Interface through R-net Electronics | \$530 |

With the SAC or DAC, please select one activation and position from the list below

| | | | | |
|--------------------------|-----------|----------|--|---------|
| <input type="checkbox"/> | SWTJSS | 22, 23 | Single actuator activation through R-net Joystick or Specialty Controls | \$1,066 |
| <input type="checkbox"/> | SWTJSD | 22, 24 | Dual actuator activation through R-net Joystick or Specialty Controls | \$1,066 |
| <input type="checkbox"/> | SWDUALACT | | Dual Activation (Allows to operate all seat functions with external switches AND through the joystick) | \$330 |
| <input type="checkbox"/> | SWTOGGLE | 23 or 24 | Single Toggle Switch (Single Mode: DAC - Dual Mode: SAC) <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation | N/C |
| <input type="checkbox"/> | SWPUSH1 | 23 | Single Push Button Switch (Single Mode: SAC) <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation | N/C |
| <input type="checkbox"/> | SWPUSH2 | 23 or 24 | Dual Push Button Switch (Single Mode: DAC - Dual Mode: SAC) <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation | N/C |

Note 20 Allows for the activation of a single power option through the VR2 or R-net 90 Amp Controllers. **Only**
 the Power Tilt or Elevating seat options are available with this feature.

Note 21 These options do not require the SAC controller as they connect via the chair's power module directly.

Note 22 Requires the R-net 90 Amp (ECEL90) or the 120 Amp (ECR120) Fully Expandable Controllers

Note 23 Requires the Single (1 Actuator) Power Seat Function Controller (SAC)

Note 24 Requires the Dual (2 Actuators) Power Seat Function Controller (DAC)

ELECTRONICS - MULTIPLE (3 to 6 ACTUATORS) POWER SEAT FUNCTION CONTROLLER

| | | | | |
|--------------------------|-----|--|---|---------|
| <input type="checkbox"/> | MAC | | Multiple (3 to 6 Actuators) Power Seat Function Controller to Interface through R-net Electronics | \$2,564 |
|--------------------------|-----|--|---|---------|

With the MAC, please select one activation and position from the list below

| | | | | |
|--------------------------|-----------|--------|---|-------|
| <input type="checkbox"/> | SWTJSM | 25, 26 | Multiple Actuator activation through R-net Joystick or Specialty Controls | STD |
| <input type="checkbox"/> | SWDUALACT | | Dual Activation (Allows to operate all seat functions with external switches and through the joystick) | \$330 |
| <input type="checkbox"/> | SWMULTI | 26 | 4-Way Toggle Switch (Up to 4 Activations) - <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation | N/C |
| <input type="checkbox"/> | SWMULTIT | 26 | Triple Toggle Switch (Up to 5 Activations) - <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation | N/C |
| <input type="checkbox"/> | SWPUSH4 | 26 | 4-Way Push Button Switch - <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation | N/C |
| <input type="checkbox"/> | SWPUSH5 | 26 | 5-Way Push Button Switch - <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation | N/C |

Note 25 Requires the Multiple (3 to 6 Actuators) Power Seat Function Controller (MAC)

Note 26 Requires the R-net 90 Amp (ECEL90) or 120 Amp (ECR120) Fully Expandable Controller, and the multiple Actuator Controller (MAC)

SWITCH ACTIVATION HAND MOUNT

| | | | | | |
|----------------------|--------------------------|----------------------------------|--------------------------|----------------------------------|-----|
| Switch Installation: | <input type="checkbox"/> | TSML (Install Switch on LH side) | <input type="checkbox"/> | TSMR (Install Switch on RH side) | N/C |
|----------------------|--------------------------|----------------------------------|--------------------------|----------------------------------|-----|

PRE-SET SEAT ANGLE (For Non-Reclining Backs)

SOBAA Please Select one: 86° 90° 94° 98° 102° 106° 110° 114° 118° 122° N/C

STANDARD ADJUSTABLE SEAT WIDTH RANGE BETWEEN 15" AND 20" STD

- SOWN15-20-19 (Width adjustment from 15" to 20" with a maximum Depth of 19")
- SOWN15-20-22 (Width adjustment from 15" to 20" with a maximum Depth of 22") ²⁷

- 15" (SOWN15) 18" (SOWN18) ²⁷
- 16" (SOWN16) ²⁷ 19" (SOWN19S) ²⁷
- 17" (SOWN17) ²⁷ 20" (SOWN20S) ²⁷

Note 27 For the Rehab Seat with the Custom Corbeau Sport Back (SAABA-COR3), please choose the long seat option SOWN15-20-22
The only possible Seat width options for the CUSTOM Corbeau Sport Seat (SAABA-COR3) are from 16" to 20" wide

LARGE ADJUSTABLE SEAT WIDTH RANGE BETWEEN 19" AND 24" \$462

- SOWN19-24-22 (Width adjustment from 19" to 24" and Depth range of 18" to 22")

- 19" (SOWN19L) 22" (SOWN22)
- 20" (SOWN20L) 23" (SOWN23)
- 21" (SOWN21) 24" (SOWN24)

ADJUSTABLE SEAT DEPTH RANGE BETWEEN 15" AND 22"

- 15" (SODN15) 16" (SODN16) 17" (SODN17) 18" (SODN18) N/C
- 19" (SODN19) 20" (SODN20) 21" (SODN21) 22" (SODN22) \$315

BACK CANE HEIGHT ²⁹ (The back cane height measurement is taken from the seat pan to the top of the push handles)

- 16" (SOBCH16) 20" (SOBCH20) 24" (SOBCH24)
- 18" (SOBCH18) 22" (SOBCH22) ²⁸ 26" (SOBCH26)

Note 28 22" is the minimum back cane height required with a power or manual recline.

Note 29 The Back Cane option is not compatible for use with the Rehab Seat with the CUSTOM Corbeau Sport Back (SAABA-COR3)

BACK SUPPORT

- STDSB ³⁰ Contour back support with a maximum cushion height of 20" (15" to 20" Wide) \$588
- STDSB/CUSTOM Contour back support (Non-Standard Sizes) Specify: Width: Height: \$1,256

Note 30 It is recommended to use the Contour Back Support option with the power or manual recline. Not required for Corbeau

³² BACK SUPPORT TOTAL HEIGHT (MEASURED FROM THE SEAT PAN TO THE TOP OF THE BACK)

- 16" (SOBH16) 20" (SOBH20) 24" (SOBH24)
- 18" (SOBH18) 22" (SOBH22) ³¹ 26" (SOBH26)

Note 31 With a power or manual recline, the minimum height is 22".

Note 32 The space between the bottom of the back and the seat pan is 4" with a non-reclining back and 6" with a reclining back

SEAT CUSHION

- SCPSC Contour Seat Cushion with a maximum seat depth of 22" (15" to 20" Wide) \$588
- SCPSC/CUSTOM Contour Seat Cushion (Non-Standard Sizes) Specify: Width: Depth: \$1,281

CANTILEVER - HEIGHT ADJUSTABLE ARMRESTS - OPENS WITH A 5° ANGLE (For Non-Reclining Backs - Choose 2)

- SOAFBFLH Flip-Back Cantilever - Full Length Left Side (Height Adjustable from 7" to 15" in 1" increments) STD
- SOAFBFRH Flip-Back Cantilever - Full Length Right Side (Height Adjustable from 7" to 15" in 1" increments) STD
- SOAFBDLH Flip-Back Cantilever - Desk Length Left Side (Height Adjustable from 7" to 15" in 1" increments) N/C
- SOAFBDRH Flip-Back Cantilever - Desk Length Right Side (Height Adjustable from 7" to 15" in 1" increments) N/C

DUAL POST - HEIGHT ADJUSTABLE ARMRESTS - INCLUDES TRANSFER GRIPS (For Non-Reclining Backs - Choose 2)

- SOADPFLH Dual Post Height Adjustable - Full Length Left Side (Height Adjustable from 9 ½" - 13 ½") \$190
- SOADPFRH Dual Post Height Adjustable - Full Length Right Side (Height Adjustable from 9 ½" - 13 ½") \$190
- SOADPDLH Dual Post Height Adjustable - Desk Length Left Side (Height Adjustable from 9 ½" - 13 ½") \$190
- SOADPDRH Dual Post Height Adjustable - Desk Length Right Side (Height Adjustable from 9 ½" - 13 ½") \$190

RECLINING FLIP-UP, HEIGHT ADJUSTABLE ARMRESTS (For Reclining Backs - Choose 2)

- SOANIFLH Reclining armrest - Full Length Left Side (Height Adjustable from 10 ½" to 14" in ½" increments) N/C
- SOANIFRH Reclining armrest - Full Length Right Side (Height Adjustable from 10 ½" to 14" in ½" increments) N/C
- SOANIDLH Reclining armrest - Desk Length Left Side (Height Adjustable from 10 ½" to 14" in ½" increments) N/C
- SOANIDRH Reclining armrest - Desk Length Right Side (Height Adjustable from 10 ½" to 14" in ½" increments) N/C

ARMREST HEIGHT ADJUSTMENT AT FACTORY

Please set the Armrest Height at: _____ " (The Armrest Height is taken from the seat pan to the top of the armpad. Not required for Dual Post)

ARMPADS

Foam Type Selection

| | | | | | | | |
|--------------------------|--------|--------------------------------|---|--------------|---|--|------|
| <input type="checkbox"/> | SOFAPD | Desk Length Standard (10") | - | Select Side: | <input type="checkbox"/> Left Hand Side | <input type="checkbox"/> Right Hand Side | STD |
| <input type="checkbox"/> | SOFAPF | Full Length Standard (14") | - | Select Side: | <input type="checkbox"/> Left Hand Side | <input type="checkbox"/> Right Hand Side | STD |
| <input type="checkbox"/> | SOWAPD | 34 Desk Length Waterfall (10") | - | Select Side: | <input type="checkbox"/> Left Hand Side | <input type="checkbox"/> Right Hand Side | \$32 |
| <input type="checkbox"/> | SOWAPF | 34 Full Length Waterfall (14") | - | Select Side: | <input type="checkbox"/> Left Hand Side | <input type="checkbox"/> Right Hand Side | \$32 |

Hybrid Gel and Foam Type Selection (The Gel always stays soft)

| | | | | | | | |
|--------------------------|--|---|---|--------------|---|--|------|
| <input type="checkbox"/> | SOGELD new | Desk Length Hybrid Gel and Foam Flat armpad (10" x 2.5") | - | Select Side: | <input type="checkbox"/> Left Hand Side | <input type="checkbox"/> Right Hand Side | \$32 |
| <input type="checkbox"/> | SOGELF | Full Length Hybrid Gel and Foam Flat armpad (14" x 2.5") | - | Select Side: | <input type="checkbox"/> Left Hand Side | <input type="checkbox"/> Right Hand Side | \$32 |
| <input type="checkbox"/> | SOWGELD new | 34 Desk Length Hybrid Gel and Foam Waterfall armpad (10") | - | Select Side: | <input type="checkbox"/> Left Hand Side | <input type="checkbox"/> Right Hand Side | \$50 |
| <input type="checkbox"/> | SOWGELF new | 34 Full Length Hybrid Gel and Foam Waterfall armpad (14") | - | Select Side: | <input type="checkbox"/> Left Hand Side | <input type="checkbox"/> Right Hand Side | \$50 |

Specialty Type Selection

| | | | | | | | |
|--------------------------|-------|--|---|--------------|---|--|-------|
| <input type="checkbox"/> | SOCAR | 35 Arm Trough with Hand Support (18" x 4 ½") | - | Select Side: | <input type="checkbox"/> Left Hand Side | <input type="checkbox"/> Right Hand Side | \$210 |
|--------------------------|-------|--|---|--------------|---|--|-------|

Note 34 The Waterfall armpads may interfere with flipping back the armrests when combined with Lateral Trunk supports

Note 35 Arm Trough Armpads may interfere with Armrests Flipping back.

FOOTRESTS

| | | | | | | |
|--------------------------|---------|---|---|---|---|---------|
| <input type="checkbox"/> | SOSA AF | 39 Adult Swing-Away Footrests ----- | Set at: | <input type="checkbox"/> 60° | <input type="checkbox"/> 70° | \$263 |
| | | Select Footrest Range: | <input type="checkbox"/> SOFLA10 (10 ½" - 15 ½") | <input type="checkbox"/> SOFLA14 (14 ½" - 19 ½") | | |
| <input type="checkbox"/> | SOAAF | 39 Adult Heavy Duty Swing-Away footrests ----- | Set at: | <input type="checkbox"/> 60° | <input type="checkbox"/> 70° | \$370 |
| | | Select Footrest Range: | <input type="checkbox"/> SOFLHD12 (12 ½" - 15 ½") | <input type="checkbox"/> SOFLHD14 (14 ½" - 19 ½") | | |
| <input type="checkbox"/> | ADDPL | 36, 39 Adult Power Elevating & Articulating Legrests | Please Select Function: | <input type="checkbox"/> Individual ³⁷ | <input type="checkbox"/> Combined ³⁸ | \$2,994 |
| | | Select Footrest Range: | <input type="checkbox"/> SOFLPW12 (12 ½" - 15 ½") | <input type="checkbox"/> SOFLPW14 (14 ½" - 19 ½") | | |

39 Please set the Footrest Length at: _____ " (The footrest length is taken from the seat pan to the heel)

Note 36 This option includes calf pads

Note 37 Individual legs count as two power functions.

Note 38 Combined legs count as one power function.

Note 39 The Footrest Length is taken from the top of the seat pan to the top of the footplate.

FOOTPLATES & ACCESSORIES

| | | | |
|--------------------------|---------|--|-------|
| <input type="checkbox"/> | SOSFAF | Flip-up Footplates | STD |
| <input type="checkbox"/> | SOAFP | Angle Adjustable & Flip-up Footplates - Each | \$137 |
| <input type="checkbox"/> | SOAHL | Heel Loops - Each | \$42 |
| <input type="checkbox"/> | SOACPFR | Calf Pads for Adult Heavy-Duty Swing-Away footrests (SOAAF) with adjustable hardware - Each | \$137 |

CENTER MOUNTS

| | | | | |
|--------------------------|---------|--|---|--|
| <input type="checkbox"/> | SOAAF90 | 42, 43, 44, 45 90° Fixed Center Mount Interface Bracket (Offers Manual adjustability from 95° to 70°) | \$242 | |
| | | Select Foot Platform Range: | <input type="checkbox"/> SOFBA12 (12 ½" - 16 ½") | <input type="checkbox"/> SOFBA16 (16 ½" - 19 ½") |
| <input type="checkbox"/> | ADDPF | 40, 41, 42, 44, 45 90° Power Center Mount Elevating & Articulating Foot Platform - *** Not available with the 12 ½" Drive Wheel | \$2,994 | |
| | | <input type="checkbox"/> SOFBPW12 (12 ½" - 16 ½") | <input type="checkbox"/> SOFBPW16 (16 ½" - 20 ½") | |

45 Please set the Foot Platform Length at: _____ " (The foot Platform length is taken from the seat pan to the heel)

Note 40 The power center mount counts as one power function.

Note 41 This option includes calf pads

Note 42 Links and attaches the one piece center mount foot platform (SOFUFP) or the split footplates (SOAFPFB) to the seat

Note 43 With the split footplates option (SOAFPFB), this bracket offers independent height adjustability for each footplate

Note 44 Please Select the one piece center mount foot platform (SOFUFP) or the split footplates (SOAFPFB)

Note 45 The Foot Platform Length is taken from the top of the seat pan to the top of the foot platform or split foot plates.

FOOT PLATFORMS & ACCESSORIES

| | | | |
|--------------------------|---------|---|-------|
| <input type="checkbox"/> | SOFUFP | Flip-up one piece Foot Platform | N/C |
| <input type="checkbox"/> | SOAFPFB | Angle Adjustable & Flip-up Split Footplates for center mount - Each | \$137 |
| <input type="checkbox"/> | SOACPCM | Calf Pads for 90° Fixed Center Mount (SOAAF90) with adjustable hardware - Each | \$137 |

SEAT POSITIONING

| | | | |
|--------------------------|--------|----------------------------|-----|
| <input type="checkbox"/> | SOAUTO | 55" Seat Positioning Strap | STD |
|--------------------------|--------|----------------------------|-----|

HEADREST ⁴⁶

| | | | |
|--------------------------|---------|--|-------|
| <input type="checkbox"/> | SONUAHR | Using headrest mounting hardware other than SOMAH located below. Specify: _____ | |
| <input type="checkbox"/> | SOFHRS | Small Curved Headrest Pad 4" x 9" | \$215 |
| <input type="checkbox"/> | SOFHRM | Medium Curved Headrest Pad 5" x 10" | \$215 |
| <input type="checkbox"/> | SOFHRL | Large Curved Headrest Pad 6" x 11" | \$215 |
| <input type="checkbox"/> | SOMAH | Multi-axis Removable Mounting Hardware | \$173 |
| <input type="checkbox"/> | SOHRBC | Bracket to Mount on the Contour Back Support (Without Recline) | \$105 |
| <input type="checkbox"/> | SOHRPR | Bracket to Mount on the Seat Frame (For Reclining Backs, but excludes the Corbeau Back) | \$105 |

* Please select one of the headrest options with the mounting hardware and specify where the headrest will be installed

Note 46 The Headrest option is not compatible with the Rehab Seat with the Custom Corbeau Sport Back option (SAABA-COR3)

LATERAL TRUNK SUPPORTS ⁴⁷

| | | | |
|--------------------------|--------|--|-------|
| <input type="checkbox"/> | SOSLS3 | Adjustable Laterals with Small Curved Pads (Each) - 3" x 5" | \$214 |
| <input type="checkbox"/> | SOSLM4 | Adjustable Laterals with Medium Curved Pads (Each) - 4" x 6" | \$214 |
| <input type="checkbox"/> | SOSLL5 | Adjustable Laterals with Large Curved Pads (Each) - 5" x 7" | \$214 |
| <input type="checkbox"/> | SOSLHM | Adjustable Swing-Away Hardware Mounts - Each | \$316 |

** Please select one of the lateral pads and the adjustable & removable hardware.

Note 47 The Lateral Trunk Support option is not compatible with the Rehab Seat with the CUSTOM Corbeau Sport Back option (SAABA-COR3)

HIP & THIGH SUPPORTS

| | | | |
|--------------------------|----------|--|-------|
| <input type="checkbox"/> | SOHTSS10 | Adjustable Hip & Thigh Support Pads (Each) - Small 10" x 4" | \$242 |
| <input type="checkbox"/> | SOHTSM12 | Adjustable Hip & Thigh Support Pads (Each) - Medium 12" x 5" | \$242 |
| <input type="checkbox"/> | SOHTSL14 | Adjustable Hip & Thigh Support Pads (Each) - Large 14" x 6" | \$242 |
| <input type="checkbox"/> | SOHTSARH | Adjustable and Removable Hip & Thigh Support Hardware - Each | \$457 |

*** Please select one of the hip & thigh support pads and the adjustable & removable hardware.

ACCESSORIES

| | | | |
|--------------------------|-------|---|-------|
| <input type="checkbox"/> | SOCH | Cup Holder (Self-Leveling) | \$95 |
| <input type="checkbox"/> | SOCPH | Ram Mounts X-Grip Cell Phone Holder | \$163 |
| <input type="checkbox"/> | SOTHS | Ram Mounts X-Grip Small Tablet Holder | \$184 |
| <input type="checkbox"/> | SOTHL | Ram Mounts X-Grip Large Tablet Holder with Swivel | \$278 |

REFLECTOR KIT

| | | | |
|--------------------------|---------|---|-----|
| <input type="checkbox"/> | REFLECT | Red & White Reflector Kit Stickers (They will not be installed - Included in a bag with the owner's manual) | N/C |
|--------------------------|---------|---|-----|

SWITCHES AND OPTIONS

| | | | |
|--------------------------|---------|--|-------|
| <input type="checkbox"/> | ASCPCP | PGDT RNET Dongle and PC programmer kit | \$615 |
| <input type="checkbox"/> | ASCPVR2 | PGDT VR2 Dongle and PC programmer kit | \$615 |
| <input type="checkbox"/> | SSBS | Buddy Button Switch | \$115 |
| <input type="checkbox"/> | SSES | Egg Switch-Black | \$115 |
| <input type="checkbox"/> | SSML | Micro Light Switch | \$158 |

SPECIALTY DRIVER CONTROL OPTIONS

| Therafin Sip-N-Puff Kit ⁵⁵ | | | |
|---|----------|--------|---|
| <input type="checkbox"/> | SCSNPK | 48, 49 | Therafin Sip-N-Puff Kit \$1,603 |
| <input type="checkbox"/> | SCSNPBTk | 49 | Therafin Sip-N-Puff Breath Tube Kit \$520 |
| <input type="checkbox"/> | SCSRMSNP | 49 | Stealth Rock-n-Lock Swingaway Arm for Sip-N-Puff <input type="checkbox"/> SWL (Swings on Left) <input type="checkbox"/> SWR (Swings on Right) \$289 |
| <input type="checkbox"/> | SCSRM | 49 | Stealth Rock-n-Lock Swingaway Mount for Mode Switch <input type="checkbox"/> SWL (Swings Left) <input type="checkbox"/> SWR (Swings Right) \$520 |
| <input type="checkbox"/> | SSES-STD | 49, 51 | Egg Switch for Changing Modes - Black <input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right) \$115 |
| <input type="checkbox"/> | SSBS-STD | 49 | Buddy Button Switch (Stop Switch) - Red \$115 |
| <input type="checkbox"/> | SOMAH | 49 | Multi-axis Headrest Removable Mounting Hardware \$173 |
| <input type="checkbox"/> | SOFHRS | | Small Curved Headrest Pad 4" x 9" \$215 |
| <input type="checkbox"/> | SOFHRM | 50 | Medium Curved Headrest Pad 5" x 10" \$215 |
| <input type="checkbox"/> | SOFHRL | | Large Curved Headrest Pad 6" x 11" \$215 |
| Switch-It Proximity Head Array ⁵⁵ | | | |
| <input type="checkbox"/> | SCSIHAA | 48, 49 | Adult - 3 Switch Proximity Head Array with 1/8" Mono Jack for mode/reverse \$5,353 |
| <input type="checkbox"/> | SCSRM | 49 | Stealth Rock-n-Lock Swingaway Mount for Mode Switch <input type="checkbox"/> SWL (Swings Left) <input type="checkbox"/> SWR (Swings Right) \$520 |
| <input type="checkbox"/> | SSES-STD | 49, 51 | Egg Switch for Changing Modes - Black <input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right) \$115 |
| <input type="checkbox"/> | SSBS-STD | 49 | Buddy Button Switch (Stop Switch) - Red \$115 |
| <input type="checkbox"/> | SOMAH | 49 | Multi-axis Headrest Removable Mounting Hardware \$173 |
| Switch-It MicroPilot Light Touch Proportional Mini joystick for chin / lip control ⁵⁵ | | | |
| <input type="checkbox"/> | SCSIMP | 48, 49 | MicroPilot Light touch mini Proportional Joystick \$4,197 |
| <input type="checkbox"/> | SCSRMMP | 49 | Stealth Rock-n-Lock Swingaway Mount for Chin M-P <input type="checkbox"/> SWL (Swings Left) <input type="checkbox"/> SWR (Swings Right) \$520 |
| <input type="checkbox"/> | SCSRM | 49 | Stealth Rock-n-Lock Swingaway Mount for Switch <input type="checkbox"/> SWL (Swings Left) <input type="checkbox"/> SWR (Swings Right) \$520 |
| <input type="checkbox"/> | SSES | 49, 51 | Egg Switch-Black <input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right) \$115 |
| <input type="checkbox"/> | JCC | | Chin Cup \$53 |
| <input type="checkbox"/> | SOMAH | 49 | Multi-axis Headrest Removable Mounting Hardware \$173 |
| <input type="checkbox"/> | SOFHRS | | Small Curved Headrest Pad 4" x 9" \$215 |
| <input type="checkbox"/> | SOFHRM | 50 | Medium Curved Headrest Pad 5" x 10" \$215 |
| <input type="checkbox"/> | SOFHRL | | Large Curved Headrest Pad 6" x 11" \$215 |
| Proportional Compact or Mini joystick for Midline Mount | | | |
| <input type="checkbox"/> | SCPGMJS | 48, 49 | PG Drives JC200 Proportional Compact Remote Joystick for Midline Mount \$998 |
| <input type="checkbox"/> | SCGMJSM | 52 | Gatlin Midline Swing-away Mount for JC200 <input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right) \$410 |
| <input type="checkbox"/> | SCSIMP | 48, 49 | Switch-It MicroPilot Light touch mini Proportional Joystick \$4,197 |
| <input type="checkbox"/> | SCGMPM | 53 | Gatlin Midline Swing-away Mount for MicroPilot <input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right) \$447 |
| <input type="checkbox"/> | SCBTM | 54 | Bullet Hand Tray for MicroPilot <input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right) \$473 |
| Proportional Compact or Mini joystick for Armrest Mount | | | |
| <input type="checkbox"/> | SCPGMJS | 48, 49 | PG Drives JC200 Proportional Compact Remote Joystick for Armrest Mount \$998 |
| <input type="checkbox"/> | SCARMJSM | 52 | Swing-away and Multi-Axis armrest mount <input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right) \$410 |
| <input type="checkbox"/> | SCSIMP | 48, 49 | Switch-It MicroPilot Light touch mini Proportional Joystick \$4,197 |
| <input type="checkbox"/> | SCARMPM | 53 | Swing-away and Multi-Axis armrest mount <input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right) \$447 |
| <input type="checkbox"/> | SCBTM | 54 | Bullet Hand Tray for MicroPilot <input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right) \$473 |
| Proportional RIM Head Control Package ⁵⁵ | | | |
| <input type="checkbox"/> | SCPHC | 48, 49 | Proportional RIM Head Control Package \$2,621 |
| <input type="checkbox"/> | SCSRM | 49 | Stealth Rock-n-Lock Swingaway Mount for Mode Switch <input type="checkbox"/> SWL (Swings Left) <input type="checkbox"/> SWR (Swings Right) \$520 |
| <input type="checkbox"/> | SSES-STD | 49, 51 | Egg Switch for Changing Modes - Black <input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right) \$115 |
| <input type="checkbox"/> | SOMAH | 49 | Multi-axis Headrest Removable Mounting Hardware \$173 |

Note 48 Requires the Omni Enhanced Display & Digital Interface Box (EDIB)

Note 49 This option must be selected for the Kit

Note 50 Choose one (1) Headrest Pad

Note 51 The Mode Switch Type can be changed from this selection on Page #7

Note 52 Must be selected for use with the PG Drives JC200 Mini Joystick

Note 53 Must be selected for use with the Switch-It MicroPilot

Note 54 This option is only compatible with the Switch-It MicroPilot

Note 55 This option is not compatible for use with the Rehab Seat with the Custom Corbeau Sport Back (SAABA-COR3)

56 CUSTOM SECTION

Note 56 Lead times will increase depending on the custom order pipeline at the time the purchase order is received
Please call for Quote and specifications for these options.

O² Holder

| | | | | | | |
|--------------------------|---------|--------|--|--|---|---------|
| <input type="checkbox"/> | SUPOXY | 57, 58 | Custom Non-Articulating Oxygen Tank Holder | <input type="checkbox"/> MLH (Mount on Left) | <input type="checkbox"/> MRH (Mount on Right) | \$410 |
| <input type="checkbox"/> | SUPOXYA | 57, 59 | Custom Articulating Oxygen Tank Holder | <input type="checkbox"/> MLH (Mount on Left) | <input type="checkbox"/> MRH (Mount on Right) | \$1,923 |
| | | | <input type="checkbox"/> D-Size Oxygen Tank (* 22" Tall x 4.3" Diameter) * The height includes the valve | | | |
| | | | <input type="checkbox"/> E-Size Oxygen Tank (* 31" Tall x 4.3" Diameter) * The height includes the valve | | | |

Note 57 The O² Holder option is not compatible for use with the Rehab Seat with Corbeau Sport Back (SAABA-COR3)
Note 58 Not available with a power tilt and/or power recline
Note 59 Required with a power tilt and/or power recline

VENT TRAY (Please complete in Full)

| | | | | | | |
|--------------------------|---|------------|--|--------------|--|---------|
| <input type="checkbox"/> | SOVENT | | Custom Vent Tray | | | \$2,574 |
| | | 60, 61, 62 | <input type="checkbox"/> Articulating Vent Box (Required with Power Tilt and/or Recline) | Model: _____ | | |
| | | 60, 61, 62 | <input type="checkbox"/> Fixed Vent Box | Model: _____ | | |
| <input type="checkbox"/> | SOVTBE | 60, 61, 62 | <input type="checkbox"/> Extended Battery Compartment | | | N/C |
| <input type="checkbox"/> | The Vent is being powered by 1 U1 Battery <input type="checkbox"/> Chair Powered by 2 Group 22NF Batteries <input type="checkbox"/> Chair Powered by 2 Group 24 Batteries | | | | | |
| <input type="checkbox"/> | The Vent is being powered by 1 Group 22NF Battery <input type="checkbox"/> Chair Powered by 2 Group 22NF Batteries | | | | | |
| <input type="checkbox"/> | The Vent is being powered by 1 Group 24 Battery <input type="checkbox"/> Chair Powered by 2 Group 22NF Batteries | | | | | |

Note 60 Amysystems does not supply the battery, charger or the wiring for the Vent.
These items are normally provided by the Vent Manufacturer.
Note 61 The Vent Tray option is not compatible for use with the Rehab Seat with Corbeau Sport Back (SAABA-COR3)
Note 62 The Vent Tray option is not available with the CUSTOM Battery Tray for Group 24 Batteries with a 24" frame width (FS2424B)

OTHER CUSTOM OPTIONS (Request For Quote)

| | | | | | | | |
|--------------------------|------------|------------|--|---|---|--|-----|
| <input type="checkbox"/> | PAT | 63, 65 | Power Anterior Tilt | | | RFQ | |
| <input type="checkbox"/> | CUSTOM | 63 | Power Lateral Tilt 15° | <input type="checkbox"/> Raise the Left Side only | <input type="checkbox"/> Raise the Right Side only | <input type="checkbox"/> Raise the Left & Right Side | RFQ |
| <input type="checkbox"/> | CUSTOM | 63, 65 | Power Latitude System (Features 14" of forward travel and a 5" seat to floor height from the pan to the ground) | | | RFQ | |
| <input type="checkbox"/> | CUSTOM | 63, 65 | Power Elevating Stretcher (Please call CSR for availability with configuration and user weight capacity) | | | RFQ | |
| <input type="checkbox"/> | CUSTOM | 63 | Power Flip-up Foot Platform (Please call CSR for availability with configuration and user weight capacity) | | | RFQ | |
| <input type="checkbox"/> | CUSTOM | 63 | Power Vertical Lifting Foot Platform (Please call CSR for availability with configuration and user weight capacity) | | | RFQ | |
| <input type="checkbox"/> | CUSTOM | 63, 64, 65 | Power Stander (Requires Power Tilt, Power Recline and 90° Power elevating & articulating center mount foot platform) | | | RFQ | |
| <input type="checkbox"/> | CUSTOM | 63 | Independent Depth Adjustability built into the Footrest Receiver that is attached to the Seat Pan | | | RFQ | |
| <input type="checkbox"/> | CUSTOM | 66, 67, 68 | Rehab Seat with Corbeau Performance Back - Includes a manual Recline Lever | <input type="checkbox"/> Lever LH | <input type="checkbox"/> Lever RH | RFQ | |
| <input type="checkbox"/> | CUSTOM | 63 | Endomorphic Back Canes (Pair) | <input type="checkbox"/> 1" per side towards the inside = 2" in | <input type="checkbox"/> 2" per side towards the inside = 4" in | RFQ | |
| <input type="checkbox"/> | CUSTOM | | Transfer Handles (Pair) | | | RFQ | |
| <input type="checkbox"/> | CUSTOM | | Manual Wheel Locks | | | RFQ | |
| <input type="checkbox"/> | CUSTOM | 63 | Non Standard Frame Colors (Not all custom colors are available, please contact CSR for availability) | | | RFQ | |
| <input type="checkbox"/> | CUSTOM # 1 | 63 | Other custom options. Specify in notes section bellow. | | | RFQ | |
| <input type="checkbox"/> | CUSTOM # 2 | 63 | Other custom options. Specify in notes section bellow. | | | RFQ | |

Note 63 Please specify details and specifications for these options.
Note 64 Please call CSR to verify if the power standing option is available with the required configuration for your user.
Note 65 This option is not available with the Rehab Seat and CUSTOM Corbeau Sport Back (SAABA-COR3)
Note 66 The maximum user weight capacity with this option is 250 lbs
The Corbeau Back is not available with the Power Recline options, but includes a manual Recline. However, the 50° Power Tilt,
Note 67 the 11" Power elevating seat module and the power elevating legs or center mount options are compatible.
The Corbeau Sport Back option only comes in one Size, but offers a seat width range from 16" to 20"
Note 68 The Custom Corbeau Seat is not available with VR2 Electronics (ECV902A)

NOTES / SPECIAL REQUESTS / CUSTOM SPECIFICATIONS