



Quickie Xplore2 STD July 2018

EIPW24



▼ ACCOUNT INFORMATION

Dealer Contact Name:	_____	P.O. # or Quote:	_____
Email Address:	_____	User Name or Number:	_____
Phone Number:	_____ Fax Number: _____	User Weight:	_____
Bill to (Name):	_____	User Height:	_____
Bill to (Location):	_____	Seat Cushion Thickness:	_____
Ship to (Name):	_____	Knee to Heel Length: Left:	_____ Right:
Address:	_____		
City:	_____		
Province:	_____ Postal Code:		

The order cannot be processed without the following:
 - A user name or number, the client's measurements and weight
 - Sunrise Medical is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
 - Order confirmations will be provided by Sunrise Medical summarizing the items selected. Please review this carefully.
 - Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

▼ STANDARD OPTIONS

Quickie Xplore 2 STD Power Base, 6-wheel independent wheel suspension, an additional seat suspension system, tie-down hardware for strap systems, 8A offboard charger, 8" front casters, and standard auto style buckle seat positioning strap, 6.5 mph motors, Weight Capacity = **300 lbs.**

▼ QUICKIE Xplore 2 STD POWER WHEELCHAIR

XPLORESTD2 QUICKIE Xplore 2 STD Power Wheelchair \$9,275.00

▼ USER WEIGHT PACKAGES

WCA3 300 lbs. Maximum user weight capacity STD

▼ STANDARD FRAME COLORS

- FCBLK Black STD
- FCSBM Matte Black STD
- FCGMC Gun Metal (Charcoal) STD
- FCSVM Silver Metallic STD
- FCCBU Candy Blue STD
- FCCRE Candy Red STD
- FCMPL Midnight Purple STD
- FCRZB Razzleberry STD
- FCMWT Matte White STD

▼ OPTIONAL FRAME COLORS (May Increase lead times)

- FCMYG Military Green (New) N/C
- FCBOR Burnt Orange Metallic N/C
- FCBAO Black Opal Metallic N/C
- FCOPM Sparkle Orange N/C
- FCMYT Military Tan (New) N/C
- FCPKP Pearl Pink N/C
- FCHDD Desert Digital Camo Hydro Graphics \$340.00
- FCHWL Woodland Camo Hydro Graphics \$340.00
- FCHCF Carbon Fiber Hydro Graphics (Matte Black base) \$340.00

Includes matching accessories with the Hydro Graphics Option. The matching accessories option is not required with these colours.

▼ MATCHING ACCESSORIES TO THE FRAME COLOR

FC***-M Matching Accessories to the Frame Color (Wheels, front access door, headrest mount and back shell when selected) \$256.00
 This option is included at no charge with the Hydro Graphics Frame Color option.

▼ TRANSPORT OPTION

- TTD4 1 4-Point wheelchair transport brackets STD
- WC19-ISO New WC-19 and ISO 7176-19 Occupied Transit Package N/C

To confirm if your overall wheelchair configuration is WC-19 and ISO 7176-19 compliant, please refer to this website: <http://wc-transportation-safety.umtri.umich.edu/crash-tested-product-lists/wheelchairs>
 Other restrictions may apply, please contact CSR for details. This option is only available on new power wheelchairs as an "on-chair" option. It is not available to install on a delivered wheelchair

Note 1 4-Point wheelchair transport brackets have not yet received WC-19 crash-testing approval.
 Use these transport brackets only to secure an unoccupied wheelchair during transport.

▼ MOTOR PACKAGE

- MP65 6.5 MPH 4-Pole Motor Package STD
- MP85 8.5 MPH 4-Pole High Speed Motor Package (Max. user weight 250 lbs--Not available with VR2 electronics ECV902A) N/C

▼ ASSISTIVE DRIVING PACKAGE

INTELLIDR 2,3 Intelli-Drive \$1,694.00

Note 2 Intelli-Drive technology provides increased tracking capabilities with enhanced traction and steering precision. This is ideal for those with marginal hand function or switch control.

Note 3 Requires the R-net PM 120 Amp Expandable Controller (ECR120)

▼ BATTERY CHARGER

<input type="checkbox"/>	BCOB	8A Offboard Battery Charger	N/C
<input type="checkbox"/>	BCLES8OPT	Lester 8 amp Offboard Charger	\$383.00

▼ BATTERY TRAY AND FRAME SIZE

<input type="checkbox"/>	FS22B	4	Battery Tray for Group 22NF Batteries (Total frame width = 23 3/8")	STD
<input type="checkbox"/>	FS24B		Battery Tray for Group 24 Batteries (Total frame width = 25 3/8")	N/C
Please Select One Option: <input type="checkbox"/> Front Caster Arms <input type="checkbox"/> Flared out (STD) <input type="checkbox"/> Straight (Optional)				
Note 4 The 23 3/8" compact frame is not compatible with most of the Center Mount Front Rigging Options due to Mechanical conflicts				
Note 5 The Straight Front Arm option is not available with any of the Center Mount Front Rigging Options due to Mechanical conflicts				

▼ BATTERIES

<input type="checkbox"/>	BAT22NFPR	MK Gel Batteries-Group 22NF (Each-On chair)	\$658.00
<input type="checkbox"/>	BAT24PR	MK Gel Batteries-Group 24 (Each-On Chair)	\$881.00

▼ CHOOSE TIRE COLOR (APPLIES TO ALL TIRES OR CASTERS)

- BLACK (New) GREY

▼ DRIVE WHEEL

<input type="checkbox"/>	DW14PT		14" Cast Wheel Pneumatic Tire	STD
<input type="checkbox"/>	DW14PI		14" Cast Wheel Pneumatic with Flat-Free Inserts	\$256.00
<input type="checkbox"/>	DW144PT	6	14" x 4" Extra Wide Cast Wheel Pneumatic Tire	\$279.00
<input type="checkbox"/>	DW144PI	6	14" x 4" Extra Wide Cast Wheel Pneumatic with Flat-Free Inserts	\$279.00
<input type="checkbox"/>	DW12PT	7	12 1/2" Cast Wheel Pneumatic Tire	N/C
<input type="checkbox"/>	DW12PI	7	12 1/2" Cast Wheel Pneumatic with Flat-Free Inserts	\$256.00
Note 6 The 14" x 4" Extra Wide Wheel increases the overall width of the Base by 2 1/2"				
The option to match the wheels to the shroud colour (FC...M) is not available with the 14" x 4" Extra Wide Wheels.				
The 14" x 4" option is not available with the VR2 electronics (ECV902A)				
Note 7 12 1/2" Drive wheels reduce overall seat to floor height by 1/4", lowers ground clearance and speed, but increases motor torque				

▼ FRONT CASTERS

<input type="checkbox"/>	FC8PTR		8" Pneumatic Casters	STD
<input type="checkbox"/>	FC8PIR		8" Pneumatic Casters with Flat-Free Inserts	\$256.00
<input type="checkbox"/>	FC9PT	8	9" Pneumatic Casters	\$173.00
<input type="checkbox"/>	FC9PI	8	9" Pneumatic Casters with Flat-Free Inserts	\$204.00
Note 8 The 9" casters will increase the minimum seat to floor height by 1/2" and are NOT compatible with the 12 1/2" Drive wheel configuration (DW12PT and DW12PI)				
The 9" casters are NOT compatible with many front rigging options. Please see details in the front rigging section.				

▼ SEAT TO FLOOR HEIGHT (Approximate measurement to the top of seat pan ± 1/4")

► Options for a Seat Only

<input type="checkbox"/>	STFH1525	9, 10	15 1/4" (Possible Seat angle of 0° only)	<input type="checkbox"/>	STFH18	18"	(Possible Seat angle of 0° and 3°)
<input type="checkbox"/>	STFH16	10	16" (Possible Seat angle of 0° only)	<input type="checkbox"/>	STFH19	19"	(Possible Seat angle of 0°, 3° and 6°)
<input type="checkbox"/>	STFH17	10	17" (Possible Seat angle of 0° and 3°)	<input type="checkbox"/>	STFH20	20"	(Possible Seat angle of 0°, 3°, 6°, 9°)

► Options for a Power Tilt or a STD Power Elevating Seat

<input type="checkbox"/>	STFH1575	9, 10	15 3/4" (Possible Seat angle of 0° only)	<input type="checkbox"/>	STFH18	18"	(Possible Seat angle of 0° only)
<input type="checkbox"/>	STFH165	10	16 1/2" (Possible Seat angle of 0° only)	<input type="checkbox"/>	STFH185	18 1/2"	(Possible Seat angle of 0° only)
<input type="checkbox"/>	STFH175		17 1/2" (Possible Seat angle of 0° only)	<input type="checkbox"/>	STFH19	19"	(Possible Seat angle of 0° only)

► Options for a Combo - Power Tilt with a STD Power Elevating Seat

<input type="checkbox"/>	STFH1775	9	17 3/4" (Possible Seat angle of 0° only)
<input type="checkbox"/>	STFH185		18 1/2" (Possible Seat angle of 0° only)

Note 9 Requires the 12 1/2" Drive wheels (DW12PT or DW12PI). This lowers ground clearance and speed, but increases motor torque.
 Note 10 The Heavy Duty Swing Away footrests (SOSAFAF) are not available with this seat to floor height.

▼ SEAT TO FLOOR HEIGHT SPECIFICATION (If not standard)

<input type="checkbox"/>	SOMHSS	Specific seat to floor height requirement other than standard. Custom upcharges may apply. Required at: _____" (± 1/4" Range)	RFQ
--------------------------	--------	----------------------------------------------------------------------------------------------------------------------------------	-----

▼ PRE-SET SEAT FRAME ANGLE

<input type="checkbox"/>	SFE-0	0° Pre-set seat frame seat angle (Must be selected when using power options)	STD
<input type="checkbox"/>	SFE-3	3° Pre-set seat frame seat angle (Only available with a seat configuration without power options)	N/C
<input type="checkbox"/>	SFE-6	6° Pre-set seat frame seat angle (Only available with a seat configuration without power options)	N/C
<input type="checkbox"/>	SFE-9	9° Pre-set seat frame seat angle (Only available with a seat configuration without power options)	N/C

▼ LIGHTING PACKAGE

<input type="checkbox"/>	EOLP	11	LED Light Package (Includes front beams, side flashers, rear tail light with flashers & brake light)	\$1,131.00
Note 11 Requires the upgraded Joystick with Light Operation or the OMNI (ECVJ6KL, ECELJL, ECRJCL, ECRJCTPL, ECRJC35B or EDIB). The Light Package is not compatible with the EPSEL option.				

▼ NON-EXPANDABLE ELECTRONICS - PG DRIVES VR2 CONTROLLER AND JOYSTICKS

<input type="checkbox"/>	ECV902A	VR2 90 Amp Drive Controller	STD
<input type="checkbox"/>	ECVJ4K	12 VR2 Standard 4-Key Joystick	STD
<input type="checkbox"/>	ECVJ6K	12 VR2 6-Key Joystick (Allows for a single power option to operate through this joystick)	N/C
<input type="checkbox"/>	ECVJ9KL	12 VR2 9-Key Joystick with Lighting Operation (Allows for a single power function to operate through this joystick)	\$168.00

Note 12 Requires the VR2 90A Drive Control (ECV902A)

VR2 Electronics can only accommodate a single power seat function. The only choices are Power Tilt or Elevating seat

▼ EXPANDABLE ELECTRONICS - PG DRIVES RNET CONTROLLERS

<input type="checkbox"/>	ECEL90	R-net EL 90 Amp Expandable Controller	\$562.00
<input type="checkbox"/>	ECR120	R-net PM 120 Amp Expandable Controller	\$1,131.00
<input type="checkbox"/>	ECRNH	13 Harness for Expandable Controller	N/C

Note 13 Must be ordered with the expandable controller ECEL90 or ECR120 to interface with Joysticks or OMNI screen

▼ EXPANDABLE ELECTRONICS - PG DRIVES RNET JOYSTICKS

<input type="checkbox"/>	ECO	Omit Joystick (For use with an Omni Display)	N/C
<input type="checkbox"/>	ECELJ	14 R-net (LED) Joystick	STD
<input type="checkbox"/>	ECRJJC	14 R-net Joystick with Color screen	\$168.00

► JOYSTICKS WITH LIGHTING OPERATION

<input type="checkbox"/>	ECELJL	14 R-net (LED) Joystick with Lighting Operation	\$168.00
<input type="checkbox"/>	ECRJCL	14 R-net Joystick with Color screen and Lighting Operation	\$222.00
<input type="checkbox"/>	ECRJCTPL	14 R-net Joystick with Color screen w/ On & Off/ Mode switch & speed pt and Lighting Operation.	\$334.00

► ADVANCED JOYSTICK

<input type="checkbox"/>	ECRJJC35	14 R-net Advanced Joystick (3.5" Screen). Lighting Operation, paddle switches, shortcut keys and stereo jacks	\$1,202.00
<input type="checkbox"/>	ECRJJC35B	14 R-net Advanced Bluetooth Joystick (3.5" Screen). Lighting Operation, paddle switches, shortcut keys and stereo jacks	\$1,202.00

Note 14 This option requires the expandable controller ECEL90 or ECR120

▼ JOYSTICK MOUNT

<input type="checkbox"/>	JSMF	Fixed Joystick Mount	N/C
<input type="checkbox"/>	JSMSA	Swing Away Joystick Mount - Height adjustable, Angle adjustable and Multi-Axis <input type="checkbox"/> Set-up with standard installation to swing-away sideways <input type="checkbox"/> Set-up at 35° to swing-away sideways and downwards	\$412.00
<input type="checkbox"/>	JSMCM	Midline Mounted Swing Away Joystick Mount	\$412.00

▼ JOYSTICK MOUNT SUPPORT

<input type="checkbox"/>	JSMFH	Joystick mount support Bracket	STD
--------------------------	-------	--------------------------------	-----

▼ JOYSTICK HAND MOUNT

<input type="checkbox"/>	HML	Install Joystick on Left hand side	<input type="checkbox"/>	HMR	Install Joystick on Right hand side
--------------------------	-----	------------------------------------	--------------------------	-----	-------------------------------------

▼ JOYSTICK HANDLE

<input type="checkbox"/>	JSH	Standard Joystick Gimble	STD
<input type="checkbox"/>	JBAT	Stick Handle 3" in Length	\$101.00
<input type="checkbox"/>	JBH	Ball Handle	\$101.00
<input type="checkbox"/>	JFBH	Foam Ball Handle	\$101.00
<input type="checkbox"/>	JTH	T-Handle	\$101.00
<input type="checkbox"/>	JMUSH	New 2.5" Mushroom Handle	\$122.00
<input type="checkbox"/>	JGP	Goalpost Handle	\$101.00

▼ ATTENDANT DRIVE CONTROL OPTION

<input type="checkbox"/>	EADC	Attendant Hand Control for R-net EL & PM Controls	\$1,131.00
--------------------------	------	---------------------------------------------------	------------

▼ ATTENDANT JOYSTICK HAND MOUNT

<input type="checkbox"/>	EADCL	Install Attendant Joystick on Rear Left hand side	<input type="checkbox"/>	EADCR	Install Attendant Joystick on Rear Right hand side
--------------------------	-------	---------------------------------------------------	--------------------------	-------	----------------------------------------------------

▼ ENHANCED DISPLAY & DIGITAL INTERFACE BOX (Included with specialty controls)

<input type="checkbox"/>	EDIB	15 Omni (Enhanced Display & Digital Interface Box with Infrared)	\$2,157.00
--------------------------	------	------------------------------------------------------------------	------------

Note 15 Allows the switch operation of alternative devices through driver controls. Requires one of the R-net fully Expandable Controllers (ECEL90 or ECR120)

▼ ENHANCED DISPLAY MOUNT

<input type="checkbox"/>	EDSAM	User Swing Away Mount	<input type="checkbox"/>	EDML (Install Display on LH side)	<input type="checkbox"/>	EDMR (Install Display on RH side)	N/C
--------------------------	-------	-----------------------	--------------------------	-----------------------------------	--------------------------	-----------------------------------	-----

▼ ENVIRONMENTAL CONTROL MODULE

<input type="checkbox"/> ECU	16	Environment control module - (Output Module)	\$1,131.00
Note 16		Requires one of the fully expandable R-net controllers (ECEL90 or ECR120)	

▼ BLUETOOTH MOUSE EMULATOR (Includes USB Dongle)

<input type="checkbox"/> EDBT	17	Bluetooth Mouse Emulator (USB Dongle included)	\$1,131.00
Note 17		Requires one of the fully expandable R-net controllers (ECEL90 or ECR120) . Not required if the ECRJC35B joystick has been selected Not limited to only one selection. Any two can be installed on the same chair, but only one can be used at a time.	
		<input type="checkbox"/> Windows and MAC OS	<input type="checkbox"/> Android Devices
			<input type="checkbox"/> iPhone, iPad and iPod

▼ POWER OPTIONS & SEATING

<input type="checkbox"/> SAABA		Amyseat Rehab Seat - Seat Frame only	\$881.00
<input type="checkbox"/> SAPTX	18	CG Power Tilt	\$7,385.00
<input type="checkbox"/> SAPRX	18, 20	Power Recline with Integrated Shear Reduction (ISR)	\$9,090.00
<input type="checkbox"/> SAPTPRX	18, 20	CG Power Tilt and Power Recline with Integrated Shear Reduction (ISR)	\$13,070.00
<input type="checkbox"/> SAPESX	18	11" Power Elevating Seat Module for transfer assist	\$3,663.00
<input type="checkbox"/> ADDPESX		Module - Add an 11" Power Elevating Seat to either SAPTX, SAPRX or SAPTPRX	\$3,404.00
<input type="checkbox"/> ADDMR	19, 20	Module - Manual Recline with Integrated Shear Reduction (ISR)	\$1,131.00
Note 18		These options include the fully adjustable Amyseat (SAABA)	
Note 19		The maximum user weight capacity with this option is 250 lbs	
Note 20		It is recommended to use the Contour Back Support option with the power or manual reclines (STDSB or STDSB/CUSTOM)	

▼ ELECTRONICS - SINGLE OR DUAL (1 OR 2 ACTUATORS) POWER SEAT FUNCTION CONTROLLERS

<input type="checkbox"/> EPSV	21, 22	Tilt or Elevating seat activation through 6-key VR2 Hand Control (ECVJ6K or ECVJ9KL)	N/C
<input type="checkbox"/> EPSEL	21, 22	Tilt or Elevating seat activation through RNET 90 Amp Expandable Controller (ECEL90)	N/C
<input type="checkbox"/> SAC		Single (1) Power Seat Function Controller to Interface through R-net Electronics	N/C
<input type="checkbox"/> DAC		Dual (2) Power Seat Function Controller to Interface through R-net Electronics	\$562.00
With the SAC or DAC, please select one activation and position from the list below			
<input type="checkbox"/> SWTJSS	23, 24	Single actuator activation through R-net Joystick or Specialty Controls EPSR2A	\$1,131.00
<input type="checkbox"/> SWTJSD	23, 25	Dual actuator activation through R-net Joystick or Specialty Controls EPSRA	\$1,131.00
<input type="checkbox"/> SWDUALACT		Dual Activation (Allows to operate all seat functions with external switches and through the joystick)	\$350.00
<input type="checkbox"/> SWTOGGLE	24 or 25	Single Toggle Switch (Single Mode: DAC - Dual Mode: SAC)	<input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation N/C
<input type="checkbox"/> SWPUSH1	24	Single Push Button Switch (Single Mode: SAC)	<input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation N/C
<input type="checkbox"/> SWPUSH2	24 or 25	Dual Push Button Switch (Single Mode: DAC - Dual Mode: SAC)	<input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation N/C
Note 21		Allows for the activation of a single power option through the VR2 or R-net 90 Amp Controllers. Only the Power Tilt or Elevating seat options are available with this feature.	
Note 22		These options do not require the SAC controller as they connect via the chair's power module directly	
Note 23		Requires the R-net 90 Amp (ECEL90) or the 120 Amp (ECR120) Fully Expandable Controllers	
Note 24		Requires the Single (1 Actuator) Power Seat Function (SAC)	
Note 25		Requires the Dual (2 Actuators) Power Seat Function (DAC)	

▼ ELECTRONICS - MULTIPLE (3 to 6 ACTUATORS) POWER SEAT FUNCTION CONTROLLER

<input type="checkbox"/> MAC		Multiple (3 to 6 Actuators) Power Seat Function Controller to Interface through R-net Electronics	\$2,720.00
With the MAC, please select one activation and position from the list below			
<input type="checkbox"/> SWTJSM	26, 27	Multiple Actuator activation through R-net Joystick or Specialty Controls MAC-EPSR	STD
<input type="checkbox"/> SWDUALACT		Dual Activation (Allows to operate all seat functions with external switches and through the joystick)	\$350.00
<input type="checkbox"/> SWMULTI	27	4-Way Toggle Switch (Up to 4 Activations)	<input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation N/C
<input type="checkbox"/> SWMULTIT	27	Triple Toggle Switch (Up to 5 Activations)	<input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation N/C
<input type="checkbox"/> SWPUSH4	27	4-Way Push Button Switch	<input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation N/C
<input type="checkbox"/> SWPUSH5	27	5-Way Push Button Switch	<input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation N/C
Note 26		Requires the Multiple (3 to 6 Actuators) Power Seat Function Controller (MAC)	
Note 27		Requires the R-net 90 Amp (ECEL90) or 120 Amp (ECR120) Fully Expandable Controller, and the multiple Actuator Controller (MAC)	

▼ SWITCH ACTIVATION HAND MOUNT

Switch Installation:	<input type="checkbox"/> TSMC (Install Switch on LH side)	<input type="checkbox"/> TSMR (Install Switch on RH side)	N/C
----------------------	-----------------------------------------------------------	-----------------------------------------------------------	-----

▼ STANDARD ADJUSTABLE SEAT WIDTH RANGE BETWEEN 15" AND 20"		
<input type="checkbox"/> SOWN15-20-19	(Width adjustment from 15" to 20" with a maximum Depth of 19")	N/C
<input type="checkbox"/> SOWN15-20-22	(Width adjustment from 15" to 20" with a maximum Depth of 22")	N/C

<input type="checkbox"/> 15" (SOWN15)	<input type="checkbox"/> 17" (SOWN17)	<input type="checkbox"/> 19" (SOWN19S)
<input type="checkbox"/> 16" (SOWN16)	<input type="checkbox"/> 18" (SOWN18)	<input type="checkbox"/> 20" (SOWN20S)

▼ LARGE ADJUSTABLE SEAT WIDTH RANGE BETWEEN 19" AND 24"		
<input type="checkbox"/> SOWN19-24-22	(Width adjustment from 19" to 24" and Depth range of 18" to 22")	\$490.00
<input type="checkbox"/> 19" (SOWN19L)	<input type="checkbox"/> 21" (SOWN21)	<input type="checkbox"/> 23" (SOWN23)
<input type="checkbox"/> 20" (SOWN20L)	<input type="checkbox"/> 22" (SOWN22)	<input type="checkbox"/> 24" (SOWN24)

▼ ADJUSTABLE SEAT DEPTH RANGE BETWEEN 15" AND 22"				
<input type="checkbox"/> 15" (SODN15)	<input type="checkbox"/> 16" (SODN16)	<input type="checkbox"/> 17" (SODN17)	<input type="checkbox"/> 18" (SODN18)	STD
<input type="checkbox"/> 19" (SODN19)	<input type="checkbox"/> 20" (SODN20)	<input type="checkbox"/> 21" (SODN21)	<input type="checkbox"/> 22" (SODN22)	\$334.00

▼ BACK CANE HEIGHT (The back cane height measurement is taken from the seat pan to the top of the push handles)		
<input type="checkbox"/> 16" (SOBCH16)	<input type="checkbox"/> 22" (SOBCH22)	
<input type="checkbox"/> 18" (SOBCH18)	<input type="checkbox"/> 24" (SOBCH24)	
<input type="checkbox"/> 20" (SOBCH20)	<input type="checkbox"/> 26" (SOBCH26)	

Note 28 22" is the minimum back cane height required with a power or manual recline.

▼ PRE-SET SEAT ANGLE (For Non-Reclining Backs)		
<input type="checkbox"/> SOBAA	Please Select one: <input type="checkbox"/> 86° <input type="checkbox"/> 90° <input type="checkbox"/> 94° <input type="checkbox"/> 98° <input type="checkbox"/> 102° <input type="checkbox"/> 106°	N/C

▼ BACK SUPPORT		
<input type="checkbox"/> STDSB	29 Contour back support 15" to 20" Wide with a maximum cushion height of 20" (15" to 20" Wide)	\$624.00
<input type="checkbox"/> STDSB/CUSTOM	Contour back support (Non Standard Sizes) Specify: _____ Width: _____ Depth: _____	\$1,333.00

Note 29 It is recommended to use the Contour Back Support option with the power or manual recline.

▼ BACK SUPPORT TOTAL HEIGHT (MEASURED FROM THE SEAT PAN TO THE TOP OF THE BACK)		
<input type="checkbox"/> 16" (SOBH16)	<input type="checkbox"/> 24" (SOBH24)	
<input type="checkbox"/> 18" (SOBH18)	<input type="checkbox"/> 26" (SOBH26)	
<input type="checkbox"/> 20" (SOBH20)	<input type="checkbox"/> 28" (SOBH28)	
<input type="checkbox"/> 22" (SOBH22)		

Note 30 22" is the minimum back height required with a power or manual recline.

Note 31 The space between the bottom of the back and the seat pan is 4" with a non-reclining back and 6" with a reclining back.

▼ SEAT CUSHIONS		
► Contour Seat Cushion		
<input type="checkbox"/> SCPSC	Contour Seat Cushion (15" to 20" Wide, with maximum 22" seat depth)	\$624.00
<input type="checkbox"/> SCPSC/CUSTOM	Contour Seat Cushion (Non Standard Sizes) Width: _____ Depth: _____	\$1,359.00

▼ CANTILEVER - HEIGHT ADJUSTABLE ARMRESTS - OPENS WITH A 5° ANGLE (For Non-Reclining Backs - Choose 2)		
<input type="checkbox"/> SOAFBFLH	Flip-Back Cantilever - Full Length Left Side (Height Adjustable from 7" to 15" with 1" increments)	STD
<input type="checkbox"/> SOAFBFRH	Flip-Back Cantilever - Full Length Right Side (Height Adjustable from 7" to 15" with 1" increments)	STD
<input type="checkbox"/> SOAFBDLH	Flip-Back Cantilever - Desk Length Left Side (Height Adjustable from 7" to 15" with 1" increments)	N/C
<input type="checkbox"/> SOAFBDRH	Flip-Back Cantilever - Desk Length Right Side (Height Adjustable from 7" to 15" with 1" increments)	N/C

▼ DUAL POST - HEIGHT ADJUSTABLE ARMRESTS - INCLUDES TRANSFER GRIPS (For Non-Reclining Backs - Choose 2)		
<input type="checkbox"/> SOADPFLH	Dual Post Height Adjustable - Full Length Left Side (Height Adjustable from 10 ¼" - 15 ¼")	\$202.00
<input type="checkbox"/> SOADPFRH	Dual Post Height Adjustable - Full Length Right Side (Height Adjustable from 10 ¼" - 15 ¼")	\$202.00
<input type="checkbox"/> SOADPDLH	Dual Post Height Adjustable - Desk Length Left Side (Height Adjustable from 10 ¼" - 15 ¼")	\$202.00
<input type="checkbox"/> SOADPDRH	Dual Post Height Adjustable - Desk Length Right Side (Height Adjustable from 10 ¼" - 15 ¼")	\$202.00

▼ RECLINING FLIP-UP, HEIGHT ADJUSTABLE ARMRESTS (For Reclining Backs - Choose 2)		
<input type="checkbox"/> SOANIFLH	Reclining armrest - Full Length Left Side (Height Adjustable from 10 ½" to 14" with ½" increments)	N/C
<input type="checkbox"/> SOANIFRH	Reclining armrest - Full Length Right Side (Height Adjustable from 10 ½" to 14" with ½" increments)	N/C
<input type="checkbox"/> SOANIDLH	Reclining armrest - Desk Length Left Side (Height Adjustable from 10 ½" to 14" with ½" increments)	N/C
<input type="checkbox"/> SOANIDRH	Reclining armrest - Desk Length Right Side (Height Adjustable from 10 ½" to 14" with ½" increments)	N/C

▼ ARMREST HEIGHT ADJUSTMENT AT FACTORY

Please set the Armrest Height at: _____ " (The Armrest Height is taken from the seat pan to the top of the armpad. Not required for Dual Post)

▼ ARMPADS

► Foam Type Selection

<input type="checkbox"/> SOFAPD	Desk Length Standard (10")	Select Side: <input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	STD
<input type="checkbox"/> SOFAPF	Full Length Standard (14")	Select Side: <input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	STD
<input type="checkbox"/> SOWAPD	32 Desk Length Waterfall (10")	Select Side: <input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$ 34.00
<input type="checkbox"/> SOWAPF	32 Full Length Waterfall (14")	Select Side: <input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$ 34.00

► Hybrid Gel and Foam Type Selection (The Gel always stays soft)

<input type="checkbox"/> SOGELD	Desk Length Hybrid Gel and Foam Flat armpad (10" x 2.5")	Select Side: <input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$ 34.00
<input type="checkbox"/> SOGELF	Full Length Hybrid Gel and Foam Flat armpad (14" x 2.5")	Select Side: <input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$ 34.00
<input type="checkbox"/> SOWGELD	32 Desk Length Hybrid Gel and Foam Waterfall armpad (10")	Select Side: <input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$ 54.00
<input type="checkbox"/> SOWGELF	32 Full Length Hybrid Gel and Foam Waterfall armpad (14")	Select Side: <input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$ 54.00

► Specialty Type Selection

<input type="checkbox"/> SOCAR	33 Arm Trough with Hand Support (18" x 4 1/2")	Select Side: <input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$ 222.00
--------------------------------	------------------------------------------------	------------------------------------------------------	------------------------------------------	-----------

Note 32 The Waterfall armpads may interfere with flipping back the armrests when combined with Lateral Trunk supports
 Note 33 Arm Trough Armpads may interfere with Armrests Flipping back.

▼ FOOTRESTS

<input type="checkbox"/> SOSAAF	37 Adult Swing-Away Footrests	Set at: <input type="checkbox"/> 60°	<input type="checkbox"/> 38 70°	\$ 279.00
<input type="checkbox"/> SOAAF	37 Adult Heavy Duty Swing-Away Footrests	Select Footrest Range: <input type="checkbox"/> SOFLA10 (10 1/2" - 15 1/2")	<input type="checkbox"/> SOFLA14 (14 1/2" - 19 1/2")	\$ 390.00
<input type="checkbox"/> ADDPL	34, 37 Adult Power Elevating & Articulating Legrests	Set at: <input type="checkbox"/> 60°	<input type="checkbox"/> 38 70°	\$ 3,177.00
		Please Select Function: <input type="checkbox"/> Individual 35	<input type="checkbox"/> Combined 36	
		Select Footrest Range: <input type="checkbox"/> SOFLPW12 (12 1/2" - 15 1/2")	<input type="checkbox"/> SOFLPW14 (14 1/2" - 19 1/2")	

Note 34 This options includes calf pads
 Note 35 Individual legs count as two power functions.
 Note 36 Combined legs count as one power function.
 Note 37 The Footrest Length is taken from the top of the seat pan to the top of the footplate.
 Note 38 The 70° footrest angle is not compatible with the 9° front casters (FC9PT and FC9PI) and may interfere with the casters when selecting the straight front caster arm option

▼ FOOTPLATES & ACCESSORIES

<input type="checkbox"/> SOSFAF	Flip-up Footplates	STD
<input type="checkbox"/> SOAFP	Angle Adjustable & Flip-up Footplates - Each	\$ 145.00
<input type="checkbox"/> SOAHL	Heel Loops - Each	\$ 44.00
<input type="checkbox"/> SOACPFPR	Calf Pads for Adult Heavy-Duty Swing-Away footrests (SOAAF) with adjustable hardware - Each	\$ 145.00

▼ CENTER MOUNTS

<input type="checkbox"/> SOAAF90	41,42,43,44, 45,46 90° Fixed Center Mount Interface Bracket (Offers Manual adjustability from 95° to 70°)	\$ 256.00
	Select Foot Platform Range: <input type="checkbox"/> SOFBA12 (12 1/2" - 16 1/2")	<input type="checkbox"/> SOFBA16 (16 1/2" - 19 1/2")
<input type="checkbox"/> ADDPF	39,40,41,43, 44,45 90° Power Center Mount Elevating & Articulating Foot Platform	\$3,177.00
	Select Foot Platform Range: <input type="checkbox"/> SOFBPW12 (12 1/2" - 16 1/2")	<input type="checkbox"/> SOFBPW16 (16 1/2" - 20 1/2")

Note 39 The power centermount counts as one power function
 Note 40 This option includes calf pads
 Note 41 Links and attaches the one piece center mount foot platform (SOFUPF) or the split footplates (SOAFPFB) to the seat
 Note 42 With the split footplates option (SOAFPFB), this bracket offers independent height adjustability for each footplate.
 Note 43 Please select the one piece center mount foot platform (SOFUPF) or the split footplates (SOAFPFB)
 Note 44 The Foot Platform length is taken from the top of the seat pan to the top of the foot platform or split footplates.
 Note 45 Centermount Front Riggings are not available with the Straight Front Arm option due to Mechanical conflicts.
 Note 46 This option is not compatible with the 9° Caster option.

▼ FOOT PLATFORMS & ACCESSORIES

<input type="checkbox"/> SOFUPF	Flip-up one piece Foot Platform	N/C
<input type="checkbox"/> SOAFPFB	Angle Adjustable & Flip-up Split Footplates for center mount - Each	\$ 145.00
<input type="checkbox"/> SOACPCM	Calf Pads for 90° Fixed Centermount (SOAAF90) with adjustable hardware - Each	\$ 145.00

▼ SEAT POSITIONING (NOT COMPATIBLE WITH THE WC19-ISO OPTION)

<input type="checkbox"/>	SOAUTO	55" Seat Positioning Strap	STD
--------------------------	--------	----------------------------	-----

▼ HEADREST

Please select one of the headrest options with the mounting hardware and specify where the headrest will be installed

<input type="checkbox"/>	SONUAHR	Using headrest mounting hardware other than SOMAH located below. Specify: _____	N/C
<input type="checkbox"/>	SOFHRS	Small Curved Headrest Pad 4" x 9"	\$228.00
<input type="checkbox"/>	SOFHRM	Medium Curved Headrest Pad 5" x 10"	\$228.00
<input type="checkbox"/>	SOFHRL	Large Curved Headrest Pad 6" x 11"	\$228.00
<input type="checkbox"/>	SOMAH	Multi-axis Removable Mounting Hardware	\$183.00
<input type="checkbox"/>	SOHRBC	Bracket to Mount on the Contour Back Support (Without Recline)	\$111.00
<input type="checkbox"/>	SOHRPR	Bracket to Mount on the Seat Frame (For Reclining Backs)	\$111.00

Please select one of the headrest options with the mounting hardware and specify where the headrest will be installed (The frame or back cushion)

▼ LATERAL TRUNK SUPPORTS

<input type="checkbox"/>	SOSLS3-C	Adjustable Laterals with Small Curved Pads (Each) - 3" x 5"	\$675.00
<input type="checkbox"/>	SOSLM4-C	Adjustable Laterals with Medium Curved Pads (Each) - 4" x 6"	\$675.00
<input type="checkbox"/>	SOSLL5-C	Adjustable Laterals with Large Curved Pads (Each) - 5" x 7"	\$675.00

▼ HIP & THIGH SUPPORTS

<input type="checkbox"/>	SOHTSS10	Adjustable Hip & Thigh Support Pads (Each) - Small 10" x 4"	\$256.00
<input type="checkbox"/>	SOHTSM12	Adjustable Hip & Thigh Support Pads (Each) - Medium 12" x 5"	\$256.00
<input type="checkbox"/>	SOHTSL14	Adjustable Hip & Thigh Support Pads (Each) - Large 14" x 6"	\$256.00
<input type="checkbox"/>	SOHTSARH	Adjustable and Removable Hip & Thigh Support Hardware - Each	\$485.00

Please select one of the hip & thigh support pads and the adjustable & removable hardware.

▼ ACCESSORIES

<input type="checkbox"/>	SOCH	Cup Holder (Self-Leveling)	\$101.00
<input type="checkbox"/>	SOCPPH	Ram Mounts X-Grip Cell Phone Holder	\$173.00
<input type="checkbox"/>	SOTHS	Ram Mounts X-Grip Small Tablet Holder	\$196.00
<input type="checkbox"/>	SOTHL	Ram Mounts X-Grip Large Tablet Holder with Swivel	\$295.00

▼ REFLECTOR KIT

<input type="checkbox"/>	REFLECT	Red & White Reflector Kit Stickers (They will not be installed - Included in a bag with the owner's manual)	N/C
--------------------------	---------	-------------------------------------------------------------------------------------------------------------	-----

▼ SWITCHES AND OPTIONS

<input type="checkbox"/>	ASCPCP	PGDT R-NET Dongle and PC programmer kit	\$652
<input type="checkbox"/>	ASCPVR2	PGDT VR2 Dongle and PC programmer kit	\$652.00
<input type="checkbox"/>	SSBS-STD	Buddy Button Switch	\$122.00
<input type="checkbox"/>	SSES-STD	Egg Switch-Black	\$122.00
<input type="checkbox"/>	SSML	Ultra Light Switch	\$168.00

▼ SPECIALTY DRIVER CONTROL OPTIONS

▼ Therafin Sip n' Puff Kit

<input type="checkbox"/>	SCSNPK	47, 48	Therafin Sip n' Puff Kit				\$1,701.00	
<input type="checkbox"/>	SCSNPBTk	48	Therafin Sip n' Puff Breath Tube Kit				\$552.00	
<input type="checkbox"/>	SCSRMSNP	48	Stealth Rock-n-Lock Swingaway Arm for Sip-N-Puff	<input type="checkbox"/>	SWL (Swings on Left)	<input type="checkbox"/>	SWR (Swings on Right)	\$307.00
<input type="checkbox"/>	SCSRM	48	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/>	SWL (Swings on Left)	<input type="checkbox"/>	SWR (Swings on Right)	\$552.00
<input type="checkbox"/>	SSES-STD	48, 50	Egg Switch for Changing Modes - Black	<input type="checkbox"/>	MLH (Mount on Left)	<input type="checkbox"/>	MRH (Mount on Right)	\$122.00
<input type="checkbox"/>	SSBS-STD	48	Buddy Button Switch (Stop Switch) - Red				\$122.00	
<input type="checkbox"/>	SOMAH	48	Multi-axis Headrest Removable Mounting Hardware				\$183.00	
<input type="checkbox"/>	SOFHRS	49	Small Curved Headrest Pad 4" x 9"				\$228.00	
<input type="checkbox"/>	SOFHRM	49	Medium Curved Headrest Pad 5" x 10"				\$228.00	
<input type="checkbox"/>	SOFHRL	49	Large Curved Headrest Pad 6" x 11"				\$228.00	

▼ Switch-It Proximity Head Array

<input type="checkbox"/>	SCSIHAA	47, 48	Adult - 3 Switch Proximity Head Array with 1/8" Mono Jack for mode/reverse (E2330)				\$5,679.00	
<input type="checkbox"/>	SCSRM	48	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/>	SWL (Swings on Left)	<input type="checkbox"/>	SWR (Swings on Right)	\$552.00
<input type="checkbox"/>	SSES-STD	48, 50	Egg Switch for Changing Modes - Black	<input type="checkbox"/>	MLH (Mount on Left)	<input type="checkbox"/>	MRH (Mount on Right)	\$122.00
<input type="checkbox"/>	SSBS-STD	48	Buddy Button Switch (Stop Switch) - Red				\$122.00	
<input type="checkbox"/>	SOMAH	48	Multi-axis Headrest Removable Mounting Hardware				\$183.00	

▼ Proportional RIM Head Control Package

<input type="checkbox"/>	SCPHC	47, 48	Proportional RIM Head Control Package				\$2,781.00	
<input type="checkbox"/>	SCSRM	48	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/>	SWL (Swings on Left)	<input type="checkbox"/>	SWR (Swings on Right)	\$552.00
<input type="checkbox"/>	SSES-STD	48, 50	Egg Switch for Changing Modes - Black	<input type="checkbox"/>	MLH (Mount on Left)	<input type="checkbox"/>	MRH (Mount on Right)	\$122.00
<input type="checkbox"/>	SOMAH	48	Multi-axis Headrest Removable Mounting Hardware				\$183.00	

▼ Switch-It MicroPilot Light Touch Proportional Mini joystick for chin / lip control

<input type="checkbox"/>	SCSIMP	47, 48	MicroPilot Light touch mini Proportional Joystick				\$4,453.00	
<input type="checkbox"/>	SCSRMMP	48	Stealth Rock-n-Lock Swingaway Arm for Chin MP	<input type="checkbox"/>	SWL (Swings on Left)	<input type="checkbox"/>	SWR (Swings on Right)	\$552.00
<input type="checkbox"/>	SCSRM	48	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/>	SWL (Swings on Left)	<input type="checkbox"/>	SWR (Swings on Right)	\$552.00
<input type="checkbox"/>	SSES-STD	48, 50	Egg Switch for Changing Modes - Black	<input type="checkbox"/>	MLH (Mount on Left)	<input type="checkbox"/>	MRH (Mount on Right)	\$122.00
<input type="checkbox"/>	JCC		Chin Cup				\$57.00	
<input type="checkbox"/>	SOMAH	48	Multi-axis Headrest Removable Mounting Hardware				\$183.00	
<input type="checkbox"/>	SOFHRS	49	Small Curved Headrest Pad 4" x 9"				\$228.00	
<input type="checkbox"/>	SOFHRM	49	Medium Curved Headrest Pad 5" x 10"				\$228.00	
<input type="checkbox"/>	SOFHRL	49	Large Curved Headrest Pad 6" x 11"				\$228.00	

▼ Proportional Compact or Mini joystick for Midline Mount

<input type="checkbox"/>	SCPGMJS	47, 48	PG Drives JC200 Proportional Compact Remote Joystick for Midline Mount				\$1,059.00	
<input type="checkbox"/>	SCGMJSM	51	Gatlin Midline Swing-away Mount for JC200	<input type="checkbox"/>	MLH (Mount on Left)	<input type="checkbox"/>	MRH (Mount on Right)	\$435.00
<input type="checkbox"/>	SCSIMP	47, 48	Switch-It MicroPilot Light touch mini Proportional Joystick				\$4,453.00	
<input type="checkbox"/>	SCGMPPM	52	Gatlin Midline Swing-away Mount for MicroPilot	<input type="checkbox"/>	MLH (Mount on Left)	<input type="checkbox"/>	MRH (Mount on Right)	\$474.00
<input type="checkbox"/>	SCBTM	53	Bullet Hand Tray for MicroPilot	<input type="checkbox"/>	MLH (Mount on Left)	<input type="checkbox"/>	MRH (Mount on Right)	\$502.00

▼ Proportional Compact or Mini joystick for Armrest Mount

<input type="checkbox"/>	SCPGMJS	49, 50	PG Drives JC200 Proportional Compact Remote Joystick for Armrest Mount				\$1,059.00	
<input type="checkbox"/>	SCARMJSM	51	Swing-away and Multi-Axis armrest mount	<input type="checkbox"/>	MLH (Mount on Left)	<input type="checkbox"/>	MRH (Mount on Right)	\$435.00
<input type="checkbox"/>	SCSIMP	47, 48	Switch-It MicroPilot Light touch mini Proportional Joystick				\$4,453.00	
<input type="checkbox"/>	SCARMPM	52	Swing-away and Multi-Axis armrest mount	<input type="checkbox"/>	MLH (Mount on Left)	<input type="checkbox"/>	MRH (Mount on Right)	\$474.00
<input type="checkbox"/>	SCBTM	53	Bullet Hand Tray for MicroPilot	<input type="checkbox"/>	MLH (Mount on Left)	<input type="checkbox"/>	MRH (Mount on Right)	\$502.00

Note 47 Requires the OMNI Enhanced Display and Digital Interface Box (EDIB)

Note 48 This option must be selected for the kit

Note 49 Choose one (1) Headrest Pad

Note 50 The mode switch type can be changed from this selection on Page #6

Note 51 Must be selected for use with the PG Drives JC200 Mini Joystick

Note 52 Must be selected for use with the Switch-It MicroPilot

Note 53 This option is only compatible with the Switch-It MicroPilot

▼ CUSTOM SECTION

Note 54 Lead times will increase depending on the custom order pipeline at the time the purchase order is received. Please call for Quote and specifications for these options.

▼ O² Holder

<input type="checkbox"/>	SUPOXY	55	Custom Non-Articulating Oxygen Tank Holder	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	\$435.00
<input type="checkbox"/>	SUPOXYA	56	Custom Articulating Oxygen Tank Holder	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	\$2,040.00
			<input type="checkbox"/> D-Size Oxygen Tank (* 22" Tall x 4.3" Diameter) * The height includes the valve			
			<input type="checkbox"/> E-Size Oxygen Tank (* 31" Tall x 4.3" Diameter) * The height includes the valve			
		Note 55	Not available with a power tilt and/or power recline			
		Note 56	Required with a power tilt and/or power recline			

▼ VENT TRAY (Please complete in Full)

<input type="checkbox"/>	SOVENT		Custom Vent Tray			
<input type="checkbox"/>	SOVENT-A	57	Articulating Vent Box (Required with Power Tilt and/or Recline)	Model: _____		\$1,957.00
<input type="checkbox"/>	SOVENT-F	57	Fixed Vent Box	Model: _____		\$491.00
<input type="checkbox"/>	SOVTBE	57, 58	Extended Battery Compartment			\$600.00
			<input type="checkbox"/> The Vent is being powered by 1 U1 Battery and the chair by 2 Group 22NF batteries			
			<input type="checkbox"/> The Vent is being powered by 1 U1 Battery and the chair by 2 Group 24 batteries			
			<input type="checkbox"/> The Vent is being powered by 1 Group 22NF Battery and the chair Powered by 2 Group 22NF Batteries			
			<input type="checkbox"/> The Vent is being powered by 1 Group 24 Battery and the chair Powered by 2 Group 22NF Batteries			
		Note 57	Sunrise Medical does not supply the battery, charger or the wiring for the Vent. These items are normally provided by the Vent Manufacturer.			
		Note 58	This option is not available with the Straight Front Caster Arms			

▼ OTHER CUSTOM OPTIONS (Request For Quote)

<input type="checkbox"/>	PAT	59, 60	Power Anterior Tilt			RFQ
<input type="checkbox"/>	CUSTOM	59	Power Lateral Tilt 15° <input type="checkbox"/> Raise the Left Side only <input type="checkbox"/> Raise the Right Side only <input type="checkbox"/> Raise the Left & Right Side			RFQ
<input type="checkbox"/>	CUSTOM	59, 60	Power Elevating Stretcher (Please call CSR for availability with configuration and user weight capacity)			RFQ
<input type="checkbox"/>	CUSTOM	59, 60	Power Flip-up Foot Platform (Please call CSR for availability with configuration and user weight capacity)			RFQ
<input type="checkbox"/>	CUSTOM	59, 62	Power Vertical Lifting Foot Platform (Please call CSR for availability with configuration and user weight capacity)			RFQ
<input type="checkbox"/>	CUSTOM	59	Independent Depth Adjustability built into the Footrest Receiver that is attached to the Seat Pan			RFQ
<input type="checkbox"/>	CUSTOM	59	Endomorphic Back Canes (Pair) <input type="checkbox"/> 1" per side towards the inside = 2" in <input type="checkbox"/> 2" per side towards the inside = 4" in			RFQ
<input type="checkbox"/>	CUSTOM		Transfer Handles (Pair)			RFQ
<input type="checkbox"/>	CUSTOM		Manual Wheel Locks			RFQ
<input type="checkbox"/>	CUSTOM	59	Non Standard Frame Colors (Not all custom colors are available, please contact CSR for availability)			RFQ
<input type="checkbox"/>	CUSTOM # 1	59	Other custom options. Specify in notes section below.			RFQ
<input type="checkbox"/>	CUSTOM # 2	59	Other custom options. Specify in notes section below.			RFQ
		Note 59	Please specify details and specifications for these options.			
		Note 60	The maximum user weight capacity with this option is 250 lbs			
		Note 61	The 11" Power Elevating Seat Module and the Power Elevating Legs or Centermount options are compatible.			
		Note 62	This option is not available with the Straight Front Caster Arms.			

▼ NOTES / SPECIAL REQUESTS / CUSTOM SPECIFICATIONS
