



Distributor Name: _____

Contact Name: _____

Address: _____

Postcode: _____

Country: _____

Telephone: _____

Email: _____

Dealers Name: _____

Client's Name: _____

Client's Weight: _____

Client's Height: _____

Client's Disability: _____

Purchase Order No: _____

Delivery Address: _____

Postcode: _____

Contact Name: _____

Telephone: _____

Email: _____

Notes/Comments:

Do Not Use for Specification Information

FOR RGK OFFICE USE ONLY

Distributor: _____

A/C No: _____

Net Value: \$ _____

Freight: \$ _____

Invoice Value: \$ _____

Payment Received: \$ _____


	Who	Date
OK to Quote:		
OK to Confirm:		
Document No:	DOC	

Serial Number: _____


Manufacturing Date: _____ /20

Base product

Standard configuration

 Mandatory information

Option with no extra cost

 Information required for option

Option with extra cost


















SQCXZ010 Club Sport Exclusive - Chrome Moly
\$ 2,575.00

Maximum user weight 125kg


Standard Specification Includes:

- Aluminium sideguards
- Height, angle and length adjustable footrest
- Quick release single anti-tip basketball style - height adjustable
- 5cm medium density foam cushion*
- Cross laced wire spoke wheels complete with aluminium hard anodised hand rim & kenda tyre
- 72mm castor wheels in extruded forks
- Adjustable seat and backrest upholstery (black)*
- Calf strap & foot strap
- Adjustable height backrest 250-300mm (25mm increments)
- 18° camber angle
- 130mm centre of gravity (adjustable +/- 50mm) (Maximum 180mm)

Key Product Measurements

<p><input type="checkbox"/> XY520 Option 1 (Kids)</p> <p>Front frame height = 480mm Rear frame height = 420mm</p> 	<p><input type="checkbox"/> XY530 Option 2</p> <p>Front frame height = 500mm Rear frame height = 400mm</p> 	<p><input type="checkbox"/> XY540 Option 3</p> <p>Front frame height = 500mm Rear frame height = 440mm</p> 	<p><input type="checkbox"/> XY550 Option 4</p> <p>Front frame height = 530mm Rear frame height = 530mm</p>  <p><small>*Aluminium seat plate & 5cm Plastazote cushion</small></p>
Wheels	Wheels	Wheels	Wheels
 <input type="checkbox"/> 22" XY415 <input type="checkbox"/> 24" XY420 <input type="checkbox"/> 25" XY430	 <input type="checkbox"/> 24" XY420 <input type="checkbox"/> 25" XY430	 <input type="checkbox"/> 24" XY420 <input type="checkbox"/> 25" XY430	 <input type="checkbox"/> 26" XY440 <input type="checkbox"/> 700c XY450
Seat width (inside of sides)	Seat width (inside of sides)	Seat width (inside of sides)	Seat width (inside of sides)
 <input type="checkbox"/> 340 mm <input type="checkbox"/> 370 mm	 <input type="checkbox"/> 370 mm <input type="checkbox"/> 400 mm <input type="checkbox"/> 440 mm	 <input type="checkbox"/> 370 mm <input type="checkbox"/> 400 mm <input type="checkbox"/> 440 mm	 <input type="checkbox"/> 370 mm <input type="checkbox"/> 400 mm <input type="checkbox"/> 440 mm
Seat length	Seat length	Seat length	Seat length
 <input type="checkbox"/> 300 mm <input type="checkbox"/> 330 mm <input type="checkbox"/> 360 mm	 <input type="checkbox"/> 360 mm <input type="checkbox"/> 390 mm <input type="checkbox"/> 430 mm	 <input type="checkbox"/> 360 mm <input type="checkbox"/> 390 mm <input type="checkbox"/> 430 mm	 <input type="checkbox"/> 360 mm <input type="checkbox"/> 390 mm <input type="checkbox"/> 430 mm
Frame length	Frame length	Frame length	Frame length
<input checked="" type="checkbox"/> 520 mm	<input checked="" type="checkbox"/> 580 mm	<input checked="" type="checkbox"/> 580 mm	 <input type="checkbox"/> 580 mm <input type="checkbox"/> 630 mm

OPTIONS

<p><input type="checkbox"/> XN010 Finish frame powder coated</p> <p> Specify Code <input style="width: 150px; height: 20px;" type="text"/></p> <p><input checked="" type="checkbox"/> XN150 Finish components powder coated (Standard finish Black)</p>	<p><input type="checkbox"/> XS020 Strapping click basic sport \$ 72.00 (not fitted)</p> <p><input type="checkbox"/> XS010 Strapping 100mm velcro adjustable sport \$ 67.00</p> <p><input type="checkbox"/> XS080 Strapping 50mm buckle fasten \$ 25.00</p> <p><input type="checkbox"/> XK080 Wheel Spinergy SLX 24 (Black only) \$ 778.00 (Upgrade from standard wheels)</p> <p><input type="checkbox"/> XP170 Upholstery with individual embroidery \$ 115.00 (Please provide .jpeg file and</p>