



# Breezy Relax<sup>2</sup>



July 2018

Account: \_\_\_\_\_  
 Date: \_\_\_\_\_  
 PO#: \_\_\_\_\_  
 Buyer: \_\_\_\_\_  
 Marked For: \_\_\_\_\_

Ship To: \_\_\_\_\_  
 Name: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 City, Province, Postal Code: \_\_\_\_\_  
 Phone: \_\_\_\_\_

### Relax<sup>2</sup> Features (Included With All Chairs):

- Adjustable Seat Depth (Adjustable 16.5-20")
- Length and Height Adjustable Armrests
- Comfort Cushion
- Comfort Back (Integrated Lateral Support)
- Push Handle Height Adjustable
- ELR - Footplate Angle Adjustable
- Silver Colour
- Relax Mounting Headrest
- Relax Headrest Bracket
- Relax Headrest Pad
- Anti Tips
- Wheel Locks - Pull To Lock
- Max. User Weight - 285 lbs

### BREEZY RELAX<sup>2</sup> SPECIFICATIONS

<input type="checkbox"/>	Breezy Relax <sup>2</sup>	\$2,224.00
<b>STF Range: 17.5" - 19.5"</b>		
<input type="checkbox"/>	069001-001 16" wide 24" spoke/solid tire/aluminum handrim. 8" caster	NC
<input type="checkbox"/>	069001-002 18" wide 24" spoke/solid tire/aluminum handrim. 8" caster	NC
<input type="checkbox"/>	069001-003 20" wide 24" spoke/solid tire/aluminum handrim. 8" caster	NC
<b>STF Range: 16" - 18.5"</b>		
<input type="checkbox"/>	069002-001 16" wide 22" spoke/solid tire/aluminum handrim. 6" caster	NC
<input type="checkbox"/>	069002-002 18" wide 22" spoke/solid tire/aluminum handrim. 6" caster	NC
<input type="checkbox"/>	069002-003 20" wide 22" spoke/solid tire/aluminum handrim. 6" caster	NC

### OPTIONS (SHIPPED SEPARATELY)

<input type="checkbox"/>	069701 Tray Slide On	\$162.00
<input type="checkbox"/>	069744 Option Kit (Trunk Support, Tray and Belt)	\$323.00
<input type="checkbox"/>	069327 Abduction Plot	\$58.00
<input type="checkbox"/>	069447-000 Trunk Support Swing Away / Detachable	\$158.00
<input type="checkbox"/>	069935 Positioning Belt	\$67.00
<input type="checkbox"/>	069444-ET Space Out Kit (2 cm)	\$5.00
<input type="checkbox"/>	064928 Std 70D Hanger	
	<input type="checkbox"/> 064928-001 - Left Silver	\$29.00
	<input type="checkbox"/> 064914 - Heel Loop (Single)	\$8.00
	<input type="checkbox"/> 064928-002 - Right Silver	\$29.00
	<input type="checkbox"/> 064914 - Heel Loop (Single)	\$8.00
<input type="checkbox"/>	069735 Hemiplegic Armrest (Left & Right Usable)	
	<input type="checkbox"/> 1 Unit	\$129.00
	<input type="checkbox"/> 2 Units	\$258.00
<input type="checkbox"/>	064280 Stump Pad <input type="checkbox"/> 064280-002 - Left <input type="checkbox"/> 064280-001 - Right	\$63.00
<input type="checkbox"/>	064914 Heel Loop	
	<input type="checkbox"/> Left Heel Loop	\$8.00
	<input type="checkbox"/> Right Heel Loop	\$8.00
<input type="checkbox"/>	069458 3 Wing Headrest	\$169.00

### ONTARIO FUNDING INFORMATION

WM000054	Breezy Relax <sup>2</sup> Wheelchair
SE0000117	Breezy Relax Comfort Back
SE0000118	Breezy Relax Comfort Cushion
SE0000119	Breezy Relax Headrest Pad