



April 1, 2017

ACCOUNT INFORMATION

Dealer Contact Name, Email Address, Phone Number, Fax Number, User Name or Number, User Weight, User Height, Seat Cushion Thickness, Knee to Heel Length, Bill to (Name), Bill to (Location), Ship to (Name), Address, City, Province, Postal Code, P.O. # or Quote

The order cannot be processed without the following:
- A user name or number, the client's measurements and weight
- Sunrise Medical is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Amysystems summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

STANDARD OPTIONS

The Quickie Xplore (HD) Hybrid Power Base, 6-wheel independent wheel suspension, an additional seat suspension system, tie-down hardware for strap systems, offboard charger, 8" front casters, and standard auto style buckle seat positioning strap, 5.9mph motors, Weight Capacity = 400lbs

QUICKIE XPLORE (HD) - 2 POWER WHEELCHAIR

- XPLOREHD2 QUICKIE XPLORE HD - 2 Power Wheelchair \$9,276

USER WEIGHT PACKAGES

- WCA4 400 lb Maximum user weight capacity \$536

Table with 2 columns: STANDARD FRAME COLORS and OPTIONNAL FRAME COLORS (May Increase lead times). Lists various color options like FCBLK, FCSBM, FCGMC, etc.

Note 1 Includes matching accessories with the Hydro Graphics option. The matching accessories option is not required with these colors

MATCHING ACCESSORIES TO THE FRAME COLOR

- FC***-M Matching Accessories to the Frame Color (Wheels, front access door, headrest mount and back shell when selected) \$242

TRANSPORT OPTION

- TTD4 2 4-Point wheelchair transport brackets N/C

Note 2 4-Point wheelchair transport brackets have not yet received WC-19 crash-testing approval. Use these transport brackets only to secure an unoccupied wheelchair during transport.

MOTOR PACKAGE

- MP65 6.5MPH 4 - Pole High-Speed Motor Package STD

ASSISTIVE DRIVING PACKAGE

- INTELLIDR 3, 4 Intelli-Drive \$1,475

Note 3 Intelli-Drive technology provides increased tracking capabilities with enhanced traction and steering precision. This is ideal for those with marginal hand function or switch control.
Note 4 Requires the R-net PM 120 Amp Expandable Controller (ECR120)

BATTERY CHARGER

- BCOB Offboard Battery Charger N/C
BCLES8OPT Lester 8 amp Offboard Battery Charger \$361

BATTERY TRAY AND FRAME SIZE

- FS24B** Battery Tray for Group 24 Batteries (Total base width = 25 3/8") N/C
 Please Select One Option: Front Caster Arms Flared out (STD) Straight (**Optional**)⁵

Note 5 The Straight Front Arm option is not available with any of the Center Mount Front Rigging Options due to Mechanical conflicts

BATTERIES

- BAT24PR** Gel Batteries - Group 24 (Each - On chair) \$830

DRIVE WHEEL

- | | | | |
|--------------------------|----------------|---|-------|
| <input type="checkbox"/> | DW14PT | 14" Cast Wheel Pneumatic Tire | STD |
| <input type="checkbox"/> | DW14PI | 14" Cast Wheel Pneumatic with Flat-Free inserts | \$242 |
| <input type="checkbox"/> | DW144PT | 6 14" x 4" Extra Wide Cast Wheel Pneumatic Tire | \$263 |
| <input type="checkbox"/> | DW144PI | 6 14" x 4" Extra Wide Cast Wheel Pneumatic with Flat-Free inserts | \$263 |
| <input type="checkbox"/> | DW12PT | 7 12 1/2" Cast Wheel Pneumatic Tire | N/C |
| <input type="checkbox"/> | DW12PI | 7 12 1/2" Cast Wheel Pneumatic with Flat-Free inserts | \$242 |

The 14" x 4" Extra Wide Wheel increases the overall width of the Base by 2 1/2"

Note 6 The option to Match the Wheels to the Shroud Color (FC...-M) is not available with the 14" x 4" Extra Wide wheels
 The 14" x 4" option is not available with VR2 electronics (ECV902A)

Note 7 12 1/2" Drive wheels reduce overall seat to floor height by 3/4", lowers ground clearance and speed, but increases motor torque

FRONT CASTERS

- | | | | |
|--------------------------|---------------|---|-------|
| <input type="checkbox"/> | FC8PTR | 8" Pneumatic Casters | STD |
| <input type="checkbox"/> | FC8PIR | 8" Pneumatic Casters with Flat-Free inserts | \$242 |
| <input type="checkbox"/> | FC9PT | 8 9" Pneumatic Casters | \$163 |
| <input type="checkbox"/> | FC9PI | 8 9" Pneumatic Casters with Flat-Free inserts | \$295 |

9" front casters increases seat to floor height by 1/2".

Note 8 9" casters are NOT compatible with the 12 1/2" Drive wheel configuration (DW12PT and DW12PI)

9" casters are NOT compatible with many front rigging options. Please see details in the front rigging sections on page #6

SEAT TO FLOOR HEIGHT (Approximate measurement to the top of seat pan ± 1/4")

Options for a Seat Only

- | | | | |
|--------------------------|-------|-----------------|---|
| <input type="checkbox"/> | 9, 11 | STFH1525 | 15 1/4" (Possible Seat angle of 0° only) |
| <input type="checkbox"/> | 11 | STFH16 | 16" (Possible Seat angle of 0° only) |
| <input type="checkbox"/> | 11 | STFH17 | 17" (Possible Seat angle of 0° and 3°) |
| <input type="checkbox"/> | | STFH18 | 18" (Possible Seat angle of 0° and 3°) |
| <input type="checkbox"/> | | STFH19 | 19" (Possible Seat angle of 0°, 3° and 6°) |
| <input type="checkbox"/> | | STFH20 | 20" (Possible Seat angle of 0°, 3°, 6°, 9°) |

Options for a Power Tilt or an HD Power Elevating Seat

- | | | | |
|--------------------------|---------|-----------------|--|
| <input type="checkbox"/> | 9,10,11 | STFH1575 | 15 3/4" (Possible Seat angle of 0° only) |
| <input type="checkbox"/> | 10, 11 | STFH165 | 16 1/2" (Possible Seat angle of 0° only) |
| <input type="checkbox"/> | 11 | STFH175 | 17 1/2" (Possible Seat angle of 0° only) |
| <input type="checkbox"/> | | STFH18 | 18" (Possible Seat angle of 0° only) |
| <input type="checkbox"/> | | STFH185 | 18 1/2" (Possible Seat angle of 0° only) |
| <input type="checkbox"/> | | STFH19 | 19" (Possible Seat angle of 0° only) |

Options for a Combo - Power Tilt with an HD Power Elevating Seat

- | | | | |
|--------------------------|---|-----------------|--|
| <input type="checkbox"/> | 9 | STFH1875 | 18 3/4" (Possible Seat angle of 0° only) |
| <input type="checkbox"/> | | STFH195 | 19 1/2" (Possible Seat angle of 0° only) |

Note 9 Requires the 12 1/2" Drive wheels (DW12PT or DW12PI). This lowers ground clearance and speed, but increases motor torque.

Note 10 These seat to floor heights are not available with the Heavy-Duty 11" Power Elevating Seat (SAPESXHD or ADDPESXHD)

Note 11 The Heavy Duty Swing-Away footrests (SOAAF) are not available with this seat to floor height

SEAT TO FLOOR HEIGHT SPECIFICATION (If not standard)

- SOMHSS** Specific seat to floor height requirement other than standard. RFQ
 Custom upcharges may apply. Required at: _____" (± 1/4" Range)

PRE-SET SEAT FRAME ANGLE

- | | | | |
|--------------------------|--------------|---|-----|
| <input type="checkbox"/> | SFE-0 | 0° Pre-set seat frame seat angle (Must be selected when using power options) | STD |
| <input type="checkbox"/> | SFE-3 | 3° Pre-set seat frame seat angle (Only available with a seat configuration without power options) | N/C |
| <input type="checkbox"/> | SFE-6 | 6° Pre-set seat frame seat angle (Only available with a seat configuration without power options) | N/C |
| <input type="checkbox"/> | SFE-9 | 9° Pre-set seat frame seat angle (Only available with a seat configuration without power options) | N/C |

LIGHTING PACKAGE

- EOLP** 12 LED Light Package (Includes front beams, side flashers, rear tail light with flashers & brake light) \$1,066
Note 12 Requires the upgraded Joystick with Light Operation or the OMNI (ECVJ6KL, ECELJL, ECRJCL, ECRJCTPL or EDIB)
 The Light Package is not compatible with the EPSEL option.

NON-EXPANDABLE ELECTRONICS - PG DRIVES VR2 CONTROLLER AND JOYSTICKS

<input type="checkbox"/>	ECV902A	VR2 90 Amp Drive Controller	STD
<input type="checkbox"/>	ECVJ4K	13 VR2 Standard 4-Key Joystick	STD
<input type="checkbox"/>	ECVJ6K	13 VR2 6-Key Joystick (Allows for a single power option to operate through this joystick)	N/C
<input type="checkbox"/>	ECVJ9KL	13 VR2 9-Key Joystick with Lighting Operation (Allows for a single power function to operate through this joystick)	\$158

Note 13 Requires the VR2 90A Drive Control (ECV902A).

VR2 Electronics can only accommodate a single power seat function. The only choices are Power Tilt or Elevating seat

EXPANDABLE ELECTRONICS - PG DRIVES RNET CONTROLLERS

<input type="checkbox"/>	ECEL90	R-net EL 90 Amp Expandable Controller	\$530
<input type="checkbox"/>	ECR120	R-net PM 120 Amp Expandable Controller	\$1,066
<input type="checkbox"/>	ECRNH	14 Harness for Expandable Controller	N/C

Note 14 Must be ordered with the expandable controller ECEL90 or ECR120 to interface with Joysticks or OMNI screen

EXPANDABLE ELECTRONICS - PG DRIVES RNET JOYSTICKS

<input type="checkbox"/>	ECO	Omit Joystick (For use with an Omni Display)	N/C
<input type="checkbox"/>	ECELJ	15 R-net (LED) Joystick	STD
<input type="checkbox"/>	ECRJJC	15 R-net Joystick with Color screen	\$158

JOYSTICKS WITH LIGHTING OPERATION

<input type="checkbox"/>	ECELJL	15 R-net (LED) Joystick with Lighting Operation	\$158
<input type="checkbox"/>	ECRJCL	15 R-net Joystick with Color screen and Lighting Operation	\$210
<input type="checkbox"/>	ECRJCTPL	15 R-net Joystick with Color screen w/ On & Off / Mode switch & speed pot and Lighting Operation	\$315

ADVANCED JOYSTICK

<input type="checkbox"/>	ECRJ35 New	15 R-net Advanced Joystick (3.5" Color screen). Lighting Operation, paddle switches, shortcut keys and stereo jacks	\$1,133
--------------------------	------------	---	---------

Note 15 This option requires the expandable controller ECEL90 or ECR120

JOYSTICK MOUNT

<input type="checkbox"/>	JSMF	Fixed Joystick Mount	N/C
<input type="checkbox"/>	JMSMA	Swing Away Joystick Mount - Height adjustable, Angle adjustable and Multi-Axis	\$388
		<input type="checkbox"/> Set-up with standard installation to swing-away sideways <input type="checkbox"/> Set-up at 35° to swing-away sideways and downwards	
<input type="checkbox"/>	JSMCM	Midline Mounted Swing Away Joystick Mount	\$388

JOYSTICK MOUNT SUPPORT

<input type="checkbox"/>	JSMFH	Joystick mount support Bracket	STD
--------------------------	-------	--------------------------------	-----

JOYSTICK HAND MOUNT

<input type="checkbox"/>	HML	Install Joystick on Left hand side	<input type="checkbox"/>	HMR	Install Joystick on Right hand side
--------------------------	-----	------------------------------------	--------------------------	-----	-------------------------------------

JOYSTICK HANDLE

<input type="checkbox"/>	JSH	Standard Joystick Gimble	STD
<input type="checkbox"/>	JBAT	Stick Handle 3" in Length	\$95
<input type="checkbox"/>	JBH	Ball Handle	\$95
<input type="checkbox"/>	JFBH	Foam Ball Handle	\$95
<input type="checkbox"/>	JTH	T-Handle	\$95
<input type="checkbox"/>	JGP	Goalpost	\$95

ATTENDANT DRIVE CONTROL OPTION

<input type="checkbox"/>	EADC	Attendant Hand Control for R-net EL & PM Controls	\$1,066
--------------------------	------	---	---------

ATTENDANT JOYSTICK HAND MOUNT

<input type="checkbox"/>	EADCL	Install Attendant Joystick on Rear Left hand side	<input type="checkbox"/>	EADCR	Install Attendant Joystick on Rear Right hand side
--------------------------	-------	---	--------------------------	-------	--

ENHANCED DISPLAY & DIGITAL INTERFACE BOX (Included with specialty controls)

<input type="checkbox"/>	EDIB	16 Omni (Enhanced Display & Digital Interface Box with Infrared)	\$2,033
--------------------------	------	--	---------

Note 16 Allows the switch operation of alternative devices through driver controls.
Requires one of the R-net fully Expandable Controllers (ECEL90 or ECR120)

ENHANCED DISPLAY MOUNT

<input type="checkbox"/>	EDSAM	User Swing Away Mount	<input type="checkbox"/>	EDML (Install Display on LH side)	<input type="checkbox"/>	EDMR (Install Display on RH side)	N/C
--------------------------	-------	-----------------------	--------------------------	-----------------------------------	--------------------------	-----------------------------------	-----

ENVIRONMENTAL CONTROL MODULE

<input type="checkbox"/>	ECU	17	Environment control module - (Output Module)	\$1,045
Note 17 Requires one of the fully expandable R-net controllers (ECEL90 or ECR120)				

BLUETOOTH MOUSE EMULATOR (Includes USB Dongle)

<input type="checkbox"/>	EDBT	18	Bluetooth Mouse Emulator (USB Dongle included)	\$1,066
<input type="checkbox"/>			Windows and MAC OS	
<input type="checkbox"/>			Android Devices	
<input type="checkbox"/>			Iphone, Ipad and Ipod	
Note 18 Requires one of the fully expandable R-net controllers (ECEL90 or ECR120) Not limited to only one selection. Any two can be installed on the same chair, but only one can be used at a time.				

POWER OPTIONS & SEATING

<input type="checkbox"/>	SAABA		Amyseat Rehab Seat - Seat Frame only	\$830
<input type="checkbox"/>	SAPTXXHD	19	Heavy-Duty CG Power Tilt	\$7,390
<input type="checkbox"/>	SAPRXHD	19, 20	Heavy-Duty Power Recline with Integrated Shear Reduction (ISR)	\$8,920
<input type="checkbox"/>	SAPTPRXHD	19, 20	Heavy-Duty CG Power Tilt and Power Recline with Integrated Shear Reduction (ISR)	\$13,255
<input type="checkbox"/>	SAPESXHD	19	Heavy-Duty 11" Power Elevating Seat Module for transfer assist	\$5,353
<input type="checkbox"/>	ADDPESXHD		Module - Add a Heavy-Duty 11" Power Elevating Seat to either SAPTXHD, SAPRXHD or SAPTPRXHD	\$4,282

Note 19 These options include the fully adjustable Amyseat (SAABA).

Note 20 It is recommended to use the Contour Back Support option with the power recline (STDSB or STDSB/CUSTOM)

ELECTRONICS - SINGLE OR DUAL (1 OR 2 ACTUATORS) POWER SEAT FUNCTION CONTROLLERS

<input type="checkbox"/>	EPSV	21, 22	Tilt or Elevating seat activator through 6-key VR2 Hand Control (ECVJ6K or ECVJ9KL)	N/C
<input type="checkbox"/>	EPSEL	21, 22	Tilt or Elevating seat activator through RNET 90 Amp Expandable Controller (ECEL90)	N/C
<input type="checkbox"/>	SAC		Single (1 Actuator) Power Seat Function Controller to Interface through R-net Electronics	N/C
<input type="checkbox"/>	DAC		Dual (2 Actuators) Power Seat Function Controller to Interface through R-net Electronics	\$530

With the SAC or DAC, please select one activation and position from the list below

<input type="checkbox"/>	SWTJSS	23, 24	Single actuator activation through R-net Joystick or Specialty Controls	STD
<input type="checkbox"/>	SWTJSD	23, 25	Dual actuator activation through R-net Joystick or Specialty Controls	STD
<input type="checkbox"/>	SWDUALACT		Dual Activation (Allows to operate all seat functions with external switches AND through the joystick)	\$330
<input type="checkbox"/>	SWTOGGLE	24 or 25	Single Toggle Switch (Single Mode: DAC - Dual Mode: SAC) <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation	N/C
<input type="checkbox"/>	SWPUSH1	24	Single Push Button Switch (Single Mode: SAC) <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation	N/C
<input type="checkbox"/>	SWPUSH2	24 or 25	Dual Push Button Switch (Single Mode: DAC - Dual Mode: SAC) <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation	N/C

Note 21 Allows for the activation of a single power option through the VR2 or R-net 90 Amp Controllers. **Only**
the Power Tilt or Elevating seat options are available with this feature.

Note 22 These options do not require the SAC controller as they connect via the chair's power module directly.

Note 23 Requires the R-net 90 Amp (ECEL90) or the 120 Amp (ECR120) Fully Expandable Controllers

Note 24 Requires the Single (1 Actuator) Power Seat Function Controller (SAC)

Note 25 Requires the Dual (2 Actuators) Power Seat Function Controller (DAC)

ELECTRONICS - MULTIPLE (3 to 6 ACTUATORS) POWER SEAT FUNCTION CONTROLLER

<input type="checkbox"/>	MAC		Multiple (3 to 6 Actuators) Power Seat Function Controller to Interface through R-net Electronics	\$2,564
--------------------------	-----	--	---	---------

With the MAC, please select one activation and position from the list below

<input type="checkbox"/>	SWTJSM	26, 27	Multiple Actuator activation through R-net Joystick or Specialty Controls	STD
<input type="checkbox"/>	SWDUALACT		Dual Activation (Allows to operate all seat functions with external switches AND through the joystick)	\$330
<input type="checkbox"/>	SWMULTI	27	4-Way Toggle Switch (Up to 4 Activations) - <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation	N/C
<input type="checkbox"/>	SWMULTIT	27	Triple Toggle Switch (Up to 5 Activations) - <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation	N/C
<input type="checkbox"/>	SWPUSH4	27	4-Way Push Button Switch - <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation	N/C
<input type="checkbox"/>	SWPUSH5	27	5-Way Push Button Switch - <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation	N/C

Note 26 Requires the Multiple (3 to 6 Actuators) Power Seat Function Controller (MAC)

Note 27 Requires the R-net 90 Amp (ECEL90) or 120 Amp (ECR120) Fully Expandable Controller, and the multiple Actuator Controller (MAC)

SWITCH ACTIVATION HAND MOUNT

Switch Installation:	<input type="checkbox"/>	TSML (Install Switch on LH side)	<input type="checkbox"/>	TSMR (Install Switch on RH side)	N/C
----------------------	--------------------------	----------------------------------	--------------------------	----------------------------------	-----

PRE-SET SEAT ANGLE (For Non-Reclining Backs)

SOBAA Please Select one: 86° 90° 98° 106°

N/C

STANDARD ADJUSTABLE SEAT WIDTH RANGE BETWEEN 19" and 26" STD

- SOWN19-24-22 (Width adjustment from 19" to 24" with a maximum Depth of 22")
- SOWN24-26-22 (Width adjustment from 24" to 26" with a maximum Depth of 22")

- | | |
|--|--|
| <input type="checkbox"/> 19" (SOWN19L) | <input type="checkbox"/> 24" (SOWN24) |
| <input type="checkbox"/> 20" (SOWN20L) | <input type="checkbox"/> 24" (SOWN24HDL) |
| <input type="checkbox"/> 21" (SOWN21) | <input type="checkbox"/> 25" (SOWN25HDL) |
| <input type="checkbox"/> 22" (SOWN22) | <input type="checkbox"/> 26" (SOWN26HDL) |
| <input type="checkbox"/> 23" (SOWN23) | |

ADJUSTABLE SEAT DEPTH RANGE BETWEEN 17" and 24"

- | | | | |
|---------------------------------------|---------------------------------------|---------------------------------------|---------------------------------------|
| <input type="checkbox"/> 17" (SODN17) | <input type="checkbox"/> 18" (SODN18) | <input type="checkbox"/> 19" (SODN19) | <input type="checkbox"/> 20" (SODN20) |
| <input type="checkbox"/> 21" (SODN21) | <input type="checkbox"/> 22" (SODN22) | <input type="checkbox"/> 23" (SODN23) | <input type="checkbox"/> 24" (SODN24) |

BACK CANE HEIGHT (The back cane height measurement is taken from the seat pan to the top of the push handles)

- | | | |
|--|--|--|
| <input type="checkbox"/> 16" (SOBCH16) | <input type="checkbox"/> 20" (SOBCH20) | <input type="checkbox"/> 24" (SOBCH24) |
| <input type="checkbox"/> 18" (SOBCH18) | <input type="checkbox"/> 22" (SOBCH22) ²⁸ | <input type="checkbox"/> 26" (SOBCH26) |

Note 28 22" is the minimum back cane height required with a power recline.

BACK SUPPORT

- | | | |
|--|--|---------|
| <input type="checkbox"/> STDSB ²⁹ | Contour back support with a maximum cushion height of 20" (19" and 20" Wide) | \$588 |
| <input type="checkbox"/> STDSB/CUSTOM | Contour back support (Non-Standard Sizes) Specify: Width: Height: | \$1,256 |

Note 29 It is recommended to use the Contour Back Support option with the power recline

³¹ **BACK SUPPORT TOTAL HEIGHT (MEASURED FROM THE SEAT PAN TO THE TOP OF THE BACK)**

- | | | | |
|---------------------------------------|---|---------------------------------------|---------------------------------------|
| <input type="checkbox"/> 16" (SOBH16) | <input type="checkbox"/> 20" (SOBH20) | <input type="checkbox"/> 24" (SOBH24) | <input type="checkbox"/> 28" (SOBH28) |
| <input type="checkbox"/> 18" (SOBH18) | <input type="checkbox"/> 22" (SOBH22) ³⁰ | <input type="checkbox"/> 26" (SOBH26) | |

Note 30 With a power recline, the minimum height is 22".

Note 31 The space between the bottom of the back and the seat pan is 4" with a non-reclining back and 6" with a reclining back

SEAT CUSHION

- | | | |
|---------------------------------------|---|---------|
| <input type="checkbox"/> SCPSC | Contour Seat Cushion with a maximum seat depth of 22" (15" to 20" Wide) | \$588 |
| <input type="checkbox"/> SCPSC/CUSTOM | Contour Seat Cushion (Non-Standard Sizes) Specify: Width: Depth: | \$1,281 |

CANTILEVER - HEIGHT ADJUSTABLE ARMRESTS - OPENS WITH A 5° ANGLE (For Non-Reclining Backs - Choose 2)

- | | | |
|-----------------------------------|---|-----|
| <input type="checkbox"/> SOAFBFLH | Flip-Back Cantilever - Full Length Left Side (Height Adjustable from 7" to 15" in 1" increments) | STD |
| <input type="checkbox"/> SOAFBFRH | Flip-Back Cantilever - Full Length Right Side (Height Adjustable from 7" to 15" in 1" increments) | STD |
| <input type="checkbox"/> SOAFBDLH | Flip-Back Cantilever - Desk Length Left Side (Height Adjustable from 7" to 15" in 1" increments) | N/C |
| <input type="checkbox"/> SOAFBDRH | Flip-Back Cantilever - Desk Length Right Side (Height Adjustable from 7" to 15" in 1" increments) | N/C |

DUAL POST - HEIGHT ADJUSTABLE ARMRESTS - INCLUDES TRANSFER GRIPS (For Non-Reclining Backs - Choose 2)

- | | | |
|-----------------------------------|--|-------|
| <input type="checkbox"/> SOADPFLH | Dual Post Height Adjustable - Full Length Left Side (Height Adjustable from 9 ½" - 13 ½") | \$190 |
| <input type="checkbox"/> SOADPFRH | Dual Post Height Adjustable - Full Length Right Side (Height Adjustable from 9 ½" - 13 ½") | \$190 |
| <input type="checkbox"/> SOADPDLH | Dual Post Height Adjustable - Desk Length Left Side (Height Adjustable from 9 ½" - 13 ½") | \$190 |
| <input type="checkbox"/> SOADPDRH | Dual Post Height Adjustable - Desk Length Right Side (Height Adjustable from 9 ½" - 13 ½") | \$190 |

RECLINING FLIP-UP, HEIGHT ADJUSTABLE ARMRESTS (For Reclining Backs - Choose 2)

- | | | |
|-----------------------------------|---|-----|
| <input type="checkbox"/> SOANIFLH | Reclining armrest - Full Length Left Side (Height Adjustable from 10 ½" to 14" in ½" increments) | N/C |
| <input type="checkbox"/> SOANIFRH | Reclining armrest - Full Length Right Side (Height Adjustable from 10 ½" to 14" in ½" increments) | N/C |
| <input type="checkbox"/> SOANIDLH | Reclining armrest - Desk Length Left Side (Height Adjustable from 10 ½" to 14" in ½" increments) | N/C |
| <input type="checkbox"/> SOANIDRH | Reclining armrest - Desk Length Right Side (Height Adjustable from 10 ½" to 14" in ½" increments) | N/C |

ARMREST HEIGHT ADJUSTMENT AT FACTORY

Please set the Armrest Height at: _____ " (The Armrest Height is taken from the seat pan to the top of the armpad. Not required for Dual Post)

ARMPADS

Foam Type Selection

<input type="checkbox"/>	SOFAPD	Desk Length Standard (10")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	STD
<input type="checkbox"/>	SOFAPF	Full Length Standard (14")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	STD
<input type="checkbox"/>	SOWAPD	33 Desk Length Waterfall (10")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$32
<input type="checkbox"/>	SOWAPF	33 Full Length Waterfall (14")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$32

Hybrid Gel and Foam Type Selection (The Gel always stays soft)

<input type="checkbox"/>	SOGELD new	Desk Length Hybrid Gel and Foam Flat armpad (10" x 2.5")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$32
<input type="checkbox"/>	SOGELF	Full Length Hybrid Gel and Foam Flat armpad (14" x 2.5")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$32
<input type="checkbox"/>	SOWGELD new	33 Desk Length Hybrid Gel and Foam Waterfall armpad (10")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$50
<input type="checkbox"/>	SOWGELF new	33 Full Length Hybrid Gel and Foam Waterfall armpad (14")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$50

Specialty Type Selection

<input type="checkbox"/>	SOCAR	34 Arm Trough with Hand Support (18" x 4 1/2")	-	Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$210
--------------------------	-------	--	---	--------------	---	--	-------

Note 33 The Waterfall armpads may interfere with flipping back the armrests when combined with Lateral Trunk supports

Note 34 Arm Trough Armpads may interfere with Armrests Flipping back.

FOOTRESTS

<input type="checkbox"/>	SOAAF	38 Heavy Duty Swing-Away footrests ----- Set at: <input type="checkbox"/> 60° ³⁹ <input type="checkbox"/> 70°					\$263
		Select Footrest Range:	<input type="checkbox"/> SOFLHD12 (12 1/2" - 15 1/2")	<input type="checkbox"/> SOFLHD14 (14 1/2" - 19 1/2")			
<input type="checkbox"/>	ADDPLHD	35, 38 Heavy Duty Power Elevating & Articulating Legrests	Select Function:	<input type="checkbox"/> Individual ³⁶	<input type="checkbox"/> Combined ³⁷		\$2,994
		Select Footrest Range:	<input type="checkbox"/> SOFLPW12 (12 1/2" - 15 1/2")	<input type="checkbox"/> SOFLPW14 (14 1/2" - 19 1/2")			

38 Please set the Footrest Length at: _____ " (The footrest length is taken from the seat pan to the heel)

Note 35 This option includes calf pads

Note 36 Individual legs count as two power functions.

Note 37 Combined legs count as one power function.

Note 38 The Footrest Length is taken from the top of the seat pan to the top of the footplate.

The 70° footrest angle is not compatible with the 9" front casters (FC9PT and FC9PI)

Note 39 The 70° footrest angle may interfere with the casters when selecting the straight front caster arms option

FOOTPLATES & ACCESSORIES

<input type="checkbox"/>	SOSFAFHD	Flip-up Footplates					STD
<input type="checkbox"/>	SOAFPHD	Angle Adjustable & Flip-up Footplates - Each					\$137
<input type="checkbox"/>	SOAHL	Heel Loops - Each					\$42
<input type="checkbox"/>	SOACPFR	Calf Pads for Adult Heavy-Duty Swing-Away footrests (SOAAF) with adjustable hardware - Each					\$137

CENTER MOUNTS

<input type="checkbox"/>	SOAAF90	42, 43, 44, 45, 46, 47 90° Fixed Center Mount Interface Bracket (Offers Manual adjustability from 95° to 70°)					\$242
		Select Foot Platform Range:	<input type="checkbox"/> SOFBA12 (12 1/2" - 16 1/2")	<input type="checkbox"/> SOFBA16 (16 1/2" - 19 1/2")			
<input type="checkbox"/>	ADDPFHD	40, 41, 42, 44, 45, 46 90° HD Power Center Mount Elevating & Articulating Foot Platform - *** Not available with the 12 1/2" Drive Wheel					\$2,994
		<input type="checkbox"/> SOFBPW12 (12 1/2" - 16 1/2")	<input type="checkbox"/> SOFBPW16 (16 1/2" - 20 1/2")				

45 Please set the Foot Platform Length at: _____ " (The foot Platform length is taken from the seat pan to the heel)

Note 40 The power center mount counts as one power function.

Note 41 This option includes calf pads

Note 42 Links and attaches the one piece center mount foot platform (SOFUFP) or the split footplates (SOAFPF) to the seat

Note 43 With the split footplates option (SOAFPF), this bracket offers independent height adjustability for each footplate

Note 44 Please Select the one piece center mount foot platform (SOFUFP) or the split footplates (SOAFPF)

Note 45 The Foot Platform Length is taken from the top of the seat pan to the top of the foot platform or split foot plates.

Note 46 Center Mount Front Riggings are not available with the Straight Front Arm option due to Mechanical conflicts

Note 47 This option is not compatible with the 9" Caster option

FOOT PLATFORMS & ACCESSORIES

<input type="checkbox"/>	SOFUFP	Flip-up one piece Foot Platform	N/C
<input type="checkbox"/>	SOAFPFB	Angle Adjustable & Flip-up Split Footplates for center mount - Each	\$137
<input type="checkbox"/>	SOACPCM	Calf Pads for 90° Fixed Center Mount (SOAAF90) with adjustable hardware - Each	\$137

SEAT POSITIONING

<input type="checkbox"/>	SOAUTO80-STD	80" Seat Positioning Strap	\$79
--------------------------	--------------	----------------------------	------

HEADREST

<input type="checkbox"/>	SONUAHR	Using headrest mounting hardware other than SOMAH located below. Specify: _____	
<input type="checkbox"/>	SOFHRS	Small Curved Headrest Pad 4" x 9"	\$215
<input type="checkbox"/>	SOFHRM	Medium Curved Headrest Pad 5" x 10"	\$215
<input type="checkbox"/>	SOFHRL	Large Curved Headrest Pad 6" x 11"	\$215
<input type="checkbox"/>	SOMAH	Multi-axis Removable Mounting Hardware	\$173
<input type="checkbox"/>	SOHRBC	Bracket to Mount on the Contour Back Support (Without Recline)	\$105
<input type="checkbox"/>	SOHRPR	Bracket to Mount on the Seat Frame (For Reclining Backs)	\$105

* Please select one of the headrest options with the mounting hardware and specify where the headrest will be installed

LATERAL TRUNK SUPPORTS

<input type="checkbox"/>	SOSLS3	Adjustable Laterals with Small Curved Pads (Each) - 3" x 5"	\$214
<input type="checkbox"/>	SOSLM4	Adjustable Laterals with Medium Curved Pads (Each) - 4" x 6"	\$214
<input type="checkbox"/>	SOSLL5	Adjustable Laterals with Large Curved Pads (Each) - 5" x 7"	\$214
<input type="checkbox"/>	SOSLHM	Adjustable Swing-Away Hardware Mounts - Each	\$316

** Please select one of the lateral pads and the adjustable & removable hardware.

HIP & THIGH SUPPORTS

<input type="checkbox"/>	SOHTSS10	Adjustable Hip & Thigh Support Pads (Each) - Small 10" x 4"	\$242
<input type="checkbox"/>	SOHTSM12	Adjustable Hip & Thigh Support Pads (Each) - Medium 12" x 5"	\$242
<input type="checkbox"/>	SOHTSL14	Adjustable Hip & Thigh Support Pads (Each) - Large 14" x 6"	\$242
<input type="checkbox"/>	SOHTSARH	Adjustable and Removable Hip & Thigh Support Hardware - Each	\$457

*** Please select one of the hip & thigh support pads and the adjustable & removable hardware.

ACCESSORIES

<input type="checkbox"/>	SOCH	Cup Holder (Self-Leveling)	\$95
<input type="checkbox"/>	SOPHP	Ram Mounts X-Grip Cell Phone Holder	\$163
<input type="checkbox"/>	SOTHS	Ram Mounts X-Grip Small Tablet Holder	\$184
<input type="checkbox"/>	SOTHL	Ram Mounts X-Grip Large Tablet Holder with Swivel	\$278

REFLECTOR KIT

<input type="checkbox"/>	REFLECT	Red & White Reflector Kit Stickers (They will not be installed - Included in a bag with the owner's manual)	N/C
--------------------------	---------	---	-----

SWITCHES AND OPTIONS

<input type="checkbox"/>	ASCPCP	PGDT RNET Dongle and PC programmer kit	\$615
<input type="checkbox"/>	ASCPVR2	PGDT VR2 Dongle and PC programmer kit	\$615
<input type="checkbox"/>	SSBS	Buddy Button Switch	\$115
<input type="checkbox"/>	SSES	Egg Switch-Black	\$115
<input type="checkbox"/>	SSML	Micro Light Switch	\$158

SPECIALTY DRIVER CONTROL OPTIONS

Therafin Sip-N-Puff Kit

<input type="checkbox"/>	SCSNPK	48, 49	Therafin Sip-N-Puff Kit		\$1,603
<input type="checkbox"/>	SCSNPBTk	49	Therafin Sip-N-Puff Breath Tube Kit		\$520
<input type="checkbox"/>	SCSRMSNP	49	Stealth Rock-n-Lock Swingaway Arm for Sip-N-Puff	<input type="checkbox"/> SWL (Swings on Left) <input type="checkbox"/> SWR (Swings on Right)	\$289
<input type="checkbox"/>	SCSRM	49	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/> SWL (Swings Left) <input type="checkbox"/> SWR (Swings Right)	\$520
<input type="checkbox"/>	SSES-STD	49, 51	Egg Switch for Changing Modes - Black	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$115
<input type="checkbox"/>	SSBS-STD	49	Buddy Button Switch (Stop Switch) - Red		\$115
<input type="checkbox"/>	SOMAH	49	Multi-axis Headrest Removable Mounting Hardware		\$173
<input type="checkbox"/>	SOFHRS		Small Curved Headrest Pad 4" x 9"		\$215
<input type="checkbox"/>	SOFHRM	50	Medium Curved Headrest Pad 5" x 10"		\$215
<input type="checkbox"/>	SOFHRL		Large Curved Headrest Pad 6" x 11"		\$215

Switch-It Proximity Head Array

<input type="checkbox"/>	SCSIHAA	48, 49	Adult - 3 Switch Proximity Head Array with 1/8" Mono Jack for mode/reverse		\$5,353
<input type="checkbox"/>	SCSRM	49	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/> SWL (Swings Left) <input type="checkbox"/> SWR (Swings Right)	\$520
<input type="checkbox"/>	SSES-STD	49, 51	Egg Switch for Changing Modes - Black	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$115
<input type="checkbox"/>	SSBS-STD	49	Buddy Button Switch (Stop Switch) - Red		\$115
<input type="checkbox"/>	SOMAH	49	Multi-axis Headrest Removable Mounting Hardware		\$173

Switch-It MicroPilot Light Touch Proportional Mini joystick for chin / lip control

<input type="checkbox"/>	SCSIMP	48, 49	MicroPilot Light touch mini Proportional Joystick		\$4,197
<input type="checkbox"/>	SCSRMMP	49	Stealth Rock-n-Lock Swingaway Mount for Chin M-P	<input type="checkbox"/> SWL (Swings Left) <input type="checkbox"/> SWR (Swings Right)	\$520
<input type="checkbox"/>	SCSRM	49	Stealth Rock-n-Lock Swingaway Mount for Switch	<input type="checkbox"/> SWL (Swings Left) <input type="checkbox"/> SWR (Swings Right)	\$520
<input type="checkbox"/>	SSES	49, 51	Egg Switch-Black	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$115
<input type="checkbox"/>	JCC		Chin Cup		\$55
<input type="checkbox"/>	SOMAH	49	Multi-axis Headrest Removable Mounting Hardware		\$173
<input type="checkbox"/>	SOFHRS		Small Curved Headrest Pad 4" x 9"		\$215
<input type="checkbox"/>	SOFHRM	50	Medium Curved Headrest Pad 5" x 10"		\$215
<input type="checkbox"/>	SOFHRL		Large Curved Headrest Pad 6" x 11"		\$215

Proportional Compact or Mini joystick for Midline Mount

<input type="checkbox"/>	SCPGMJS	48, 49	PG Drives JC200 Proportional Compact Remote Joystick for Midline Mount		\$998
<input type="checkbox"/>	SCGMJSM	52	Gatlin Midline Swing-away Mount for JC200	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$410
<input type="checkbox"/>	SCSIMP	48, 49	Switch-It MicroPilot Light touch mini Proportional Joystick		\$4,197
<input type="checkbox"/>	SCGMPPM	53	Gatlin Midline Swing-away Mount for MicroPilot	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$447
<input type="checkbox"/>	SCBTM	54	Bullet Hand Tray for MicroPilot	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$473

Proportional Compact or Mini joystick for Armrest Mount

<input type="checkbox"/>	SCPGMJS	48, 49	PG Drives JC200 Proportional Compact Remote Joystick for Armrest Mount		\$998
<input type="checkbox"/>	SCARMJSM	52	Swing-away and Multi-Axis armrest mount	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$410
<input type="checkbox"/>	SCSIMP	48, 49	Switch-It MicroPilot Light touch mini Proportional Joystick		\$4,197
<input type="checkbox"/>	SCARMPPM	53	Swing-away and Multi-Axis armrest mount	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$447
<input type="checkbox"/>	SCBTM	54	Bullet Hand Tray for MicroPilot	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$473

Proportional RIM Head Control Package

<input type="checkbox"/>	SCPHC	48, 49	Proportional RIM Head Control Package		\$2,621
<input type="checkbox"/>	SCSRM	49	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/> SWL (Swings Left) <input type="checkbox"/> SWR (Swings Right)	\$520
<input type="checkbox"/>	SSES-STD	49, 51	Egg Switch for Changing Modes - Black	<input type="checkbox"/> MLH (Mount on Left) <input type="checkbox"/> MRH (Mount on Right)	\$115
<input type="checkbox"/>	SOMAH	49	Multi-axis Headrest Removable Mounting Hardware		\$173

Note 48 Requires the Omni Enhanced Display & Digital Interface Box (EDIB)

Note 49 This option must be selected for the Kit

Note 50 Choose one (1) Headrest Pad

Note 51 The Mode Switch Type can be changed from this selection on Page #7

Note 52 Must be selected for use with the PG Drives JC200 Mini Joystick

Note 53 Must be selected for use with the Switch-It MicroPilot

Note 54 This option is only compatible with the Switch-It MicroPilot

55 CUSTOM SECTION

Note 55 Lead times will increase depending on the custom order pipeline at the time the purchase order is received
Please call for Quote and specifications for these options.

O² Holder

- | | | | | | | |
|--------------------------|---------|----|--|--|---|---------|
| <input type="checkbox"/> | SUPOXY | 56 | Custom Non-Articulating Oxygen Tank Holder | <input type="checkbox"/> MLH (Mount on Left) | <input type="checkbox"/> MRH (Mount on Right) | \$410 |
| <input type="checkbox"/> | SUPOXYA | 57 | Custom Articulating Oxygen Tank Holder | <input type="checkbox"/> MLH (Mount on Left) | <input type="checkbox"/> MRH (Mount on Right) | \$1,923 |
| | | | <input type="checkbox"/> D-Size Oxygen Tank (* 22" Tall x 4.3" Diameter) * The height includes the valve | | | |
| | | | <input type="checkbox"/> E-Size Oxygen Tank (* 31" Tall x 4.3" Diameter) * The height includes the valve | | | |

Note 56 Not available with a power tilt and/or power recline

Note 57 Required with a power tilt and/or power recline

VENT TRAY (Please complete in Full)

- | | | | | | | |
|--------------------------|---|--------|--|--------------|--|---------|
| <input type="checkbox"/> | SOVENT | | Custom Vent Tray | | | \$2,574 |
| | | 58 | <input type="checkbox"/> Articulating Vent Box (Required with Power Tilt and/or Recline) | Model: _____ | | |
| | | 58 | <input type="checkbox"/> Fixed Vent Box | Model: _____ | | |
| <input type="checkbox"/> | SOVTBE | 58, 59 | <input type="checkbox"/> Extended Battery Compartment | | | N/C |
| <input type="checkbox"/> | The Vent is being powered by 1 U1 Battery <input type="checkbox"/> Chair Powered by 2 Group 22NF Batteries <input type="checkbox"/> Chair Powered by 2 Group 24 Batteries | | | | | |
| <input type="checkbox"/> | The Vent is being powered by 1 Group 22NF Battery <input type="checkbox"/> Chair Powered by 2 Group 22NF Batteries | | | | | |
| <input type="checkbox"/> | The Vent is being powered by 1 Group 24 Battery <input type="checkbox"/> Chair Powered by 2 Group 22NF Batteries | | | | | |

Note 58 Amysystems does not supply the battery, charger or the wiring for the Vent.

These items are normally provided by the Vent Manufacturer.

Note 59 This option is not available with the Straight Front Caster Arms

OTHER CUSTOM OPTIONS (Request For Quote)

- | | | | | | | |
|--------------------------|------------|--------|---|--|--|-----|
| <input type="checkbox"/> | PAT | 60 | Power Anterior Tilt | | | RFQ |
| <input type="checkbox"/> | CUSTOM | 60, 61 | Power Flip-up Foot Platform (Please call CSR for availability with configuration and user weight capacity) | | | RFQ |
| <input type="checkbox"/> | CUSTOM | | Manual Wheel Locks | | | RFQ |
| <input type="checkbox"/> | CUSTOM | | Non Standard Frame Colors (Not all custom colors are available, please contact CSR for availability) | | | RFQ |
| <input type="checkbox"/> | CUSTOM | | Transfer Handles (Pair) | | | RFQ |
| <input type="checkbox"/> | CUSTOM | 60 | Endomorphic Back Canes (Pair) <input type="checkbox"/> 1" per side towards the inside = 2" in <input type="checkbox"/> 2" per side towards the inside = 4" in | | | RFQ |
| <input type="checkbox"/> | CUSTOM | 60 | Independent Depth Adjustability built into the Footrest Receiver that is attached to the Seat Pan | | | RFQ |
| <input type="checkbox"/> | CUSTOM # 1 | 60 | Other custom options. Specify in notes section bellow. | | | RFQ |
| <input type="checkbox"/> | CUSTOM # 2 | 60 | Other custom options. Specify in notes section bellow. | | | RFQ |

Note 60 Please specify details and specifications for these options.

Note 61 This option is not available with the Straight Front Caster Arms

NOTES / SPECIAL REQUESTS / CUSTOM SPECIFICATIONS
