



Quickie QM-7 Series Power Recline Seating

July 2018

EIPW11



Account: _____	Ship To: _____
Date: _____	Name: _____
PO#: _____	Address: _____
Buyer: _____	City, Province, Postal Code: _____
Marked For: _____	Phone: _____

✓ Option #	Definition	Price	✓ Option #	Definition	Price
------------	------------	-------	------------	------------	-------

1. QUICKIE QM-7 Series Power Recline Seating (Select One)

<u>Tilt AND Recline</u>			
QUICKIE QM-710™ MPC			
<input type="checkbox"/>	Quickie QM-710™ MPC Tilt and Recline		\$7,795.00
172PC69	Base Package		\$7,795.00
QUICKIE QM-715™ HD MPC			
<input type="checkbox"/>	Quickie QM-715™ HD MPC Tilt and Recline		\$8,365.00
172PC104	Base Package		\$8,365.00
QUICKIE QM-720™ MPD			
<input type="checkbox"/>	Quickie QM-720™ MPD Tilt and Recline		\$8,630.00
172PC87	Base Package		\$8,630.00

<u>Recline ONLY</u>			
QUICKIE QM-710™ SC			
<input type="checkbox"/>	Quickie QM-710™ SC for ASAP II Recline ONLY		\$7,795.00
172PC68	Base Package		\$7,795.00
QUICKIE QM-715™ HD SC			
<input type="checkbox"/>	Quickie QM-715™ HD SC for ASAP II Recline ONLY		\$8,365.00
172PC80	Base Package		\$8,365.00
QUICKIE QM-720™ SD			
<input type="checkbox"/>	Quickie QM-720™ SD for ASAP II Recline ONLY		\$8,630.00
172PC86	Base Package		\$8,630.00

MPC Seating Specifications

Typical Chair Weight	W/ Battery 375 lbs.	W/O Battery 285 lbs.
Transit	WC-19 occupied approved	
Seat Type	3.7 Power Recline. NA with QM-715	
Seat size range	Width 14"-22"	Depth 16"-22"
STFH range	MPC 16" - 19.3"	
Power Seating	50° CG tilt and Power recline w/ power shear reduction / 12" Seat elevator	

Actuators Required by Power Seating Function

Power Recline with Power Shear Reduction:	2 actuators
Combination : Recline, Tilt:	3 actuators
Combination : Recline, Seat elevate:	3 actuators
Combination : Recline, Power Legrests:	4 actuators
Combination : Recline, Tilt, Power Legrests:	5 actuators
Combination : Recline, Seat elevate, Power Legrests:	5 actuators
Combination : Recline, Tilt, Power Legrests, Seat elevate:	6 actuators



QM-710™ MPC
Shown w/ optional Seat Elevator

QUICKIE QM-7™ Base Specifications

Drive Wheel Position	Midwheel Drive	Drive Wheels	14" Drive Wheels w/ flat free insert (Black)
User weight Cap	300 lbs. - 350 lbs.	Tires/Casters	6" Front / 6" rear (Black)
Base size	Width 25" Length 36" Turn Radius 20"	Charger	8A Dual Mode Charger (Compact) CTE
Drive Control	R-net PM120 (Non-Expandable and Expandable)	Battery Size	Group 24 / Group 34.
Driving Range	Up to 13.5 miles	Shroud Colors	12 standard colors, 4 upgrade prints.
Drive Motors	4-pole 6 mph, Optional 5 mph and 8 mph	Ground Clearance	4"

2. BASE CHAIR OPTIONS

Speed Package (4-pole 300 lbs.)			
<input type="checkbox"/>	172PF85	QM-710 Standard 6 mph motor package	NC
<input type="checkbox"/>	172PF91	QM-710 5 mph High Torque motor package	\$557.00
<input type="checkbox"/>	172PF90	QM-710 8 mph Community motor package	\$1,170.00
<input type="checkbox"/>	172PF84	QM-715 HD Standard 5 mph motor package	NC
<input type="checkbox"/>	172PF89	QM-720 Standard 8 mph motor package	NC

Suretrac Drive Tracking			
<input type="checkbox"/>	172MI01	Suretrac Module and Motor Encoders	\$1,565.00

Shroud Color			
<input type="checkbox"/>	172EZ01	Candy Apple Red	NC
<input type="checkbox"/>	172EZ02	High Gloss Black	NC
<input type="checkbox"/>	172EZ03	Candy Blue	NC
<input type="checkbox"/>	172EZ04	Silver	NC
<input type="checkbox"/>	172EZ05	Copper	NC
<input type="checkbox"/>	172EZ08	Green Apple	NC
<input type="checkbox"/>	172EZ09	Yellow	NC
<input type="checkbox"/>	172EZ10	Pearl Pink	NC
<input type="checkbox"/>	172EZ11	Candy Purple	NC
<input type="checkbox"/>	172EZ12	Matte Black	NC
<input type="checkbox"/>	172EZ13	Midnight Blue	NC
<input type="checkbox"/>	172EZ15	Pearl White	NC
<input type="checkbox"/>	172EZ06	Digital Camo	\$557.00
<input type="checkbox"/>	172EZ07	Woodland Camo	\$557.00
<input type="checkbox"/>	172EZ14	Carbon Fiber Print	\$557.00
<input type="checkbox"/>	172EZ16	American Flag	\$557.00
<input type="checkbox"/>	172EZ17	Pink Camo	\$557.00
<input type="checkbox"/>	172EZ18	Grunge Skulls	\$557.00

User weight			
<input type="checkbox"/>	172WG06	Less than 150 lbs.	NC
<input type="checkbox"/>	172WG13	151-200 lbs.	NC
<input type="checkbox"/>	172WG14	201-250 lbs.	NC
<input type="checkbox"/>	172WG03	251-300 lbs.	NC
<input type="checkbox"/>	172WG15	301-350 lbs. Transit Package NA	NC

Drive Wheel			
<input type="checkbox"/>	172DW59	14" Cast Wheel Solid Dual Density Inserts Black Knobby Tires and Black Caster	Std

Transit Option			
<input type="checkbox"/>	172TO25	4-Point Transit Tie-down Package <i>Meets ISO-7176 for Occupant Use</i>	NC
<input type="checkbox"/>	172WN01	4-Point Transit Tie-down Package <i>Meets WC-19 for Occupant Use</i>	NC
<input type="checkbox"/>	172TO40	Omit Transit Tie-down	NC

WC-19 Transit Lap Belt			
<input type="checkbox"/>	172WN02	WC-19 Transit Lap Belt	\$222.00

Wheel locks			
<input type="checkbox"/>	172WL42	Wheel locks	\$150.00

Ventilator Options			
<input type="checkbox"/>	172VT4	Fixed Vent Tray (Pulmonetic LTV Models)	\$562.00

Pre-Set Seat Frame Angle <i>(STFH increases 1" per increment change)</i>			
<input type="checkbox"/>	172FA10	0°	NC

Seat Height (Use STFH Matrix below)			
<i>(Measured in the front) - Obtain Quote to verify est. STFH</i>			
<input type="checkbox"/>	172HT44	Position "0"	NC
<input type="checkbox"/>	172HT49	+5" STFH	NC
<input type="checkbox"/>	172HT45	+1" STFH	NC
<input type="checkbox"/>	172HT46	+2" STFH	NC
<input type="checkbox"/>	172HT47	+3" STFH	NC

LED Light Package			
<input type="checkbox"/>	172LI04	LED Light Package: Includes headlights, tail lights, turn indicators and flashers.	\$1,131.00

STFH Matrix					
Seating Type and Power Option Combinations	Seat-to-Floor Height at 0 tilt				
	0	+5"	+1"	+2"	+3"
SC Power Recline only			16.1	17.1	18.1
SC Power Recline and 9" Lift Seat		16.2	16.7	17.7	18.7
MPC Power Recline and Tilt.		16.3	17.3	18.3	18.3
MPC Power Recline with Tilt & 12" Lift combo		16**	16.5	17.5	18.5

STFH measured with batteries but no user. Actual STFH may vary by .5"
 *Adding seat angle will increase STFH by 1" for each incremental change.
 **Ultra low STFH may cause ground clearance of Power Centermount to drop below 2".

3. BATTERIES and CHARGERS

Batteries			
<input type="checkbox"/>	172AG21	GP22 Flat Top Batteries Pair w/strap	\$626.00
<input type="checkbox"/>	172AG04	GP24 MK Gel Flat Top Batts Pair w/strap <i>Recommended for Max Daily Range</i>	\$881.00
<input type="checkbox"/>	172AG22	GP24 Brightway AGM Flat Top Batteries Pair w/strap.	\$844.00
<input type="checkbox"/>	172AG11	GP34 MK Gel Flat Top Batteries (Pair) w/straps, <i>Required for lower STFH.</i>	\$881.00
<input type="checkbox"/>	172AG23	GP34 Brightway AGM Flat Top Batts (Pair) w/straps, <i>Required for lower STFH.</i>	\$796.00

Battery Charger			
<input type="checkbox"/>	172BJ19	8A Dual Mode (Compact) CTE	NC
<input type="checkbox"/>	172BJ29	10A Dual Mode Power First Upgrade	NC

4. ELECTRONICS

Drive Control Options			
Expandable Controllers			
<input type="checkbox"/>	172ER06	R-net Expandable Drive Control (120A)	\$1,131.00
	172EX05	Expandable 120A Controller - R-net	\$1,131.00
	172EX06	Harness for 120A Expandable Controller	NC

R-net Joystick			
<input type="checkbox"/>	172J70	Std LED Joystick for R-net Controls	NC
<input type="checkbox"/>	172J108	CJSM1: Joystick Standard: w/Stereo Jacks & Color Display	NC
<input type="checkbox"/>	172J114	CJSM2: Joystick Advanced: 3.5" Color Display, Paddle Switches, Shortcut Keys, IR, BT, and Stereo Jack.	NC
<input type="checkbox"/>	172J95	Omit JS. For use with OMNI	NC

Joystick Mount (Fixed, Swing-Away, Retractable, Midline)			
<input type="checkbox"/>	172JM10	Fixed Mount	NC
<input type="checkbox"/>	172JM11	Swing Away Mount	\$412.00
<input type="checkbox"/>	172JM15	Retractable Joystick mount w/45° stops	\$765.00
<input type="checkbox"/>	172JM14	Mid Line Mount	\$412.00

Joystick Mount (Location)			
<i>Height Adjustable.</i>			
<input type="checkbox"/>	172HM1	Right Side Mount	NC
<input type="checkbox"/>	172HM2	Left Side Mount	NC

Joystick Handle (Optional)			
<input type="checkbox"/>	172JH1	Ball Handle	\$101.00
<input type="checkbox"/>	172JH2	Foam Ball Handle	\$101.00
<input type="checkbox"/>	172JH4	T-Handle	\$101.00
<input type="checkbox"/>	172JH6	Goalpost	\$101.00

Attendant Control			
<input type="checkbox"/>	172A04	Attendant Control for R-net	\$1,131.00

Attendant Mounting location			
<input type="checkbox"/>	172DN1	Mounting location: Right Back Cane	NC
<input type="checkbox"/>	172DN2	Mounting Location: Left Back Cane	NC

R-net OMNI and Output Options			
<i>(For expandable controls, choose all that apply)</i>			
<input type="checkbox"/>	172ON6	R-net OMNI Input Control Module and display with built in Infrared and armrest mount. Choose mount below. (required for Specialty Input Devices.)	\$2,201.00
<input type="checkbox"/>	172ON7	R-net Output Module (ECM)	\$1,131.00
<input type="checkbox"/>	172ON8	R-net Bluetooth Mouse-mover for Windows / Android	\$1,131.00
<input type="checkbox"/>	172ON9	R-net Bluetooth for Apple Devices	\$1,131.00

OMNI Mount			
<i>(If used with Joystick must be mounted opposite side)</i>			
<input type="checkbox"/>	172DS04	Swing Away Left Side OMNI Mount	NC
<input type="checkbox"/>	172DS03	Swing Away Right Side OMNI Mount	NC

Switch inputs for Joysticks or OMNI			
<i>Drive or Mode / Profile Control</i>			
<input type="checkbox"/>	172SW22	Stereo to 2 Mono Splitter	\$79.00
<input type="checkbox"/>	172SW1	Buddy Button Switch	\$101.00
<input type="checkbox"/>	172SW4	Micro Light Switch	\$101.00
<input type="checkbox"/>	172SW24	Egg Switch-Black	\$127.00
<input type="checkbox"/>	172SW21	Dual toggle switch and bracket	\$227.00
<input type="checkbox"/>	172SW35	Switch-It Mini 1 button on Mono plug	\$284.00
<input type="checkbox"/>	172SW36	Switch-It Mini 2 button on 2 Mono plugs	\$284.00
<input type="checkbox"/>	172SW37	Switch-It Mini 2 button on 1 Stereo plug	\$284.00
<input type="checkbox"/>	172SW38	Switch-It Mini 4 button on DB9 connector	\$284.00

Switch Mount Bracket			
<input type="checkbox"/>	172EM01	Single Switch Mount Bracket	\$262.00

R-net Electronic Programmer			
<input type="checkbox"/>	172PR5	R-net PC Programmer with Dongle	\$652.00

5. Specialty Driver Control Options

Requires expandable R-net controls and OMNI unless otherwise noted.

5a. Proportional Control Packages

Head Control Proportional Packages			
<input type="checkbox"/> Switch-It Dual Pro Proportional Head Control w/Whitmyer and Facial Pads (Adult) \$6,939.00			
172ID72		Switch-It Dual Pro Head Control w/facial pad	\$6,513.00
172IJ04		SA Boom Mount with Egg Switch Control	\$426.00
<input type="checkbox"/> Switch-It Dual Pro Proportional Head Control w/Spot Pads Package \$6,939.00			
172ID70		Switch-It Dual Pro Head Control w/Spot pads	\$6,513.00
172IJ04		SA Boom Mount with Egg Switch Control	\$426.00
<input type="checkbox"/> Switch-It Dual Pro Proportional Head Control w/ Facial Pads and Gooseneck Mounted Proximity Switch Pkg \$6,939.00			
172ID72		Switch-It Dual Pro Head Control w/facial pad	\$6,513.00
172IJ11		SA Gooseneck Mount with Proximity Switch Control	\$426.00
<input type="checkbox"/> Switch-It Dual Pro Proportional Head Control w/Spot Pads and Gooseneck Mounted Proximity Switch Pkg \$6,939.00			
172ID70		Switch-It Dual Pro Head Control w/Spot pads	\$6,513.00
172IJ11		SA Gooseneck Mount with Proximity Switch Control	\$426.00

Hand Control Proportional Joystick Packages			
<input type="checkbox"/> Switch-It Versa-Guide Compact Joystick Midline with Gatlin Package \$2,513.00			
172ID71		Switch-It Versa-Guide Compact Proportional Joystick	\$2,163.00
172IJ07		S/A Gatlin Mount	\$350.00
<input type="checkbox"/> Switch-It Versa-Guide Compact Joystick Midline and prox flex switches with Gatlin Package \$3,201.00			
172ID71		Switch-It Versa-Guide Compact Proportional Joystick	\$2,163.00
172IM08		Dual Flex Proximity Switches (requires no force)	\$688.00
172IJ07		S/A Gatlin Mount	\$350.00
<input type="checkbox"/> Switch-It Micro-Pilot Joystick with Space Disc Package \$4,251.00			
172ID62		Micro-Pilot Proportional Joystick	\$3,749.00
172IJ08		S/A SpaceDisc mount	\$502.00
<input type="checkbox"/> Switch-It Micro-Guide Joystick with SpaceDisc Package \$5,126.00			
172ID63		Micro-Guide Proportional Joystick	\$4,624.00
172IJ08		S/A SpaceDisc mount	\$502.00

Chin Control Proportional Joystick Packages			
<input type="checkbox"/> Switch-It Versa-Guide Chin Control Compact Joystick and Egg Switch Package \$2,857.00			
172ID71		Switch-It Versa-Guide Compact Proportional Joystick	\$2,163.00
172IM07		S/A Mount off backrest	\$268.00
172IJ04		SA Boom Mount with egg switch	\$426.00
<input type="checkbox"/> Switch-It Versa-Guide Chin Control Compact Joystick and Proximity Flex Switches Package \$3,119.00			
172ID71		Switch-It Versa-Guide Compact Proportional Joystick	\$2,163.00
172IM08		Dual Flex Proximity Switches (requires no force)	\$688.00
172IM07		S/A Mount off backrest	\$268.00
<input type="checkbox"/> Switch-It Versa-Guide Bib Chin Control Compact Joystick with Proximity Flex switches Package Std 6" \$3,654.00			
172ID71		Switch-It Versa-Guide Compact Proportional Joystick	\$2,163.00
172IM08		Dual Flex Proximity Switches (requires no force)	\$688.00
172IM05		Bib mount Small	\$803.00
<input type="checkbox"/> Switch-It Versa-Guide Bib Chin Control Compact Joystick with Proximity Flex Switches Package Large 8" \$3,725.00			
172ID71		Switch-It Versa-Guide Compact Proportional Joystick	\$2,163.00
172IM08		Dual Flex Proximity Switches (requires no force)	\$688.00
172IM06		Bib mount Large	\$874.00
<input type="checkbox"/> Switch-It Versa-Guide Bib Chin Control Compact Joystick with Egg Switch Package Std 6" \$3,392.00			
172ID71		Switch-It Versa-Guide Compact Proportional Joystick	\$2,163.00
172IM05		Bib mount Small	\$803.00
172IJ04		SA Boom Mount with egg switch	\$426.00
<input type="checkbox"/> Switch-It Versa-Guide Bib Chin Control Compact Joystick with Egg Switch Package Large 8" \$3,463.00			
172ID71		Switch-It Versa-Guide Compact Proportional Joystick	\$2,163.00
172IM06		Bib mount Large	\$874.00
172IJ04		SA Boom Mount with egg switch	\$426.00
<input type="checkbox"/> Switch-It Micro-Pilot Bib Chin Control with Egg Switch Package Std 6" \$4,978.00			
172ID62		Micro-Pilot Proportional Joystick	\$3,749.00
172IM05		Bib mount Small	\$803.00
172IJ04		SA Boom Mount with egg switch	\$426.00
<input type="checkbox"/> Switch-It Micro-Pilot Bib Chin Control with Egg Switch Package Large 8" \$5,049.00			
172ID62		Micro-Pilot Proportional Joystick	\$3,749.00
172IM06		Bib mount Large	\$874.00
172IJ04		SA Boom Mount with egg switch	\$426.00
<input type="checkbox"/> Switch-It Micro-Guide Bib Chin Control with Egg Switch Package Std 6" \$5,853.00			
172ID63		Micro-Guide Proportional Joystick	\$4,624.00
172IM05		Bib mount Small	\$803.00
172IJ04		SA Boom Mount with egg switch	\$426.00
<input type="checkbox"/> Switch-It Micro-Guide Bib Chin Control with Egg Switch Package Large 8" \$5,924.00			
172ID63		Micro-Guide Proportional Joystick	\$4,624.00
172IM06		Bib mount Large	\$874.00
172IJ04		SA Boom Mount with egg switch	\$426.00

5. Specialty Driver Control Options (cont.)

5b. Switch Control

Requires Headrest and/or Headrest Mount

3 Switch Proximity Head Arrays Choose headrest mount below			
<input type="checkbox"/> Switch-It Head Array Package 3 Switch Curved (Long) \$5,679.00			
172II13		Switch-It Proximity Head Array Long	\$5,253.00
172IJ09		SA Boom Mount with Egg Switch Control	\$426.00
<input type="checkbox"/> Switch-It Head Array Package 3 switch curved (Stubby) \$6,037.00			
172II14		Switch-It Proximity Head Array Stubby	\$5,611.00
172IJ09		SA Boom Mount with Egg Switch Control	\$426.00

Whitmyer Headrest Mount for Specialty Controls			
<input type="checkbox"/>	172MH23	Whitmyer Linx Headrest Mount, Quick Release	\$515.00



Whitmyer Headrest Pad for Specialty controls			
<input type="checkbox"/>	172HE79	Plush Pad 8" for Specialty controls	\$328.00
<input type="checkbox"/>	172HE80	Plush Pad 10" for Specialty controls	\$328.00
<input type="checkbox"/>	172HE81	Plush Pad 14" for Specialty controls	\$328.00
<input type="checkbox"/>	172HE90	Adjusta-A-Plush Headrest Pad Narrow	\$328.00
<input type="checkbox"/>	172HE91	Adjusta-A-Plush Headrest Pad Standard	\$328.00

5 Switch Proximity Head Array (No Egg Switch Needed) Choose "Headrest Mount for Special Controls"			
<input type="checkbox"/>	172II04	Switch-It Head Array 5 switch flat	\$5,679.00

Extra Swing Away Head Boom Egg Switch and Mount			
<input type="checkbox"/>	172IJ05	Extra S/A Boom Mount with egg switch #1	\$122.00
<input type="checkbox"/>	172IJ06	Extra S/A Boom Mount with egg switch #2	\$122.00

Sip n' Puff Package Choose headrest and HR mount below			
<input type="checkbox"/> Quickie Sip n' Puff System \$3,100.00			
172ID34		Sip n' Puff Kit	\$2,189.00
172IH1		Sip n' Puff Tubing and SA Hardware	\$485.00
172IJ10		S/A Boom Mount with egg switch	\$426.00

Single Switches			
<input type="checkbox"/>	172IB02	Buddy Button Switch	\$114.00
<input type="checkbox"/>	172IB05	Micro Light Switch	\$101.00
<input type="checkbox"/>	172IB06	Egg Switch-Black	\$142.00
<input type="checkbox"/>	172IB07	Switch-It Mini 1 button on Mono plug	\$284.00

6. SEAT TYPE

<input type="checkbox"/>	172ST57	Power Recline	NC
--------------------------	---------	---------------	----

Recline Seat Frame size			
<input type="checkbox"/>	172FU26	Power Recline Seat Frame (G) 14x16 to 22x20 (Not width Adjustable)	NC
<input type="checkbox"/>	172FU27	Power Recline Seat Frame (H) 18x18 to 22"x22" (Not width Adjustable)	NC

7. SEAT SIZE

Seat Size Matrix. *Seat type is chosen in section 5.*

G/H = Recline Seat (Width not adjustable)
 Note: Power Recline Seat Depth measured from solid surface of back to front of seat pan. Please accommodate for backrest cushion thickness.

	Recline Seat Width					
	14"	16"	18"	20"	22"	
Seat Depth (Adjustable)	16"	G	G	G	G	G
	17"	G	G	G	G	G
	18"	G	G	G,H	G,H	G,H
	19"	G	G	G,H	G,H	G,H
	20"	G	G	G,H	G,H	G,H
	21"			H	H	H
	22"			H	H	H

Seat Width				
<input type="checkbox"/>	172R23	14"	Seat width	NC
<input type="checkbox"/>	172R5	16"	Seat width	NC
<input type="checkbox"/>	172R7	18"	Seat width	NC
<input type="checkbox"/>	172R9	20"	Seat width	NC
<input type="checkbox"/>	172R25	22"	Seat width	NC

Seat Depth				
<input type="checkbox"/>	172RF5	16"	Seat depth	NC
<input type="checkbox"/>	172RF6	17"	Seat depth	NC
<input type="checkbox"/>	172RF7	18"	Seat depth	NC
<input type="checkbox"/>	172RF8	19"	Seat depth	NC
<input type="checkbox"/>	172RF24	20"	Seat depth	NC
<input type="checkbox"/>	172RF26	21"	Seat depth	NC
<input type="checkbox"/>	172RF25	22"	Seat depth	NC

8. SEAT CUSHIONS

JAY Cushion Bundling			
<input type="checkbox"/>	172JC14	JAY Cushion On Chair Order <i>Choose Cushion in section 17c.</i>	NC

9. POWERED SEATING

Power Seating Options			
<input type="checkbox"/>	172PA24	Recline Only w/Power Shear Reduction (2 actuators)	\$9,090.00
<input type="checkbox"/>	172PA25	Tilt and Recline w/Power Shear Reduction (3 actuators)	\$13,075.00
<input type="checkbox"/>	172PA16	QM-715 HD Recline Only w/Power Shear Reduction Only (350 lbs.. weight limit) (2 actuators)	\$7,390.00
<input type="checkbox"/>	172PA17	QM-715 HD Tilt and Recline w/Power Shear Reduction (350 lbs.. weight limit) (3 actuators)	\$13,075.00

Powered Seat Elevator- (AKA: Power Lift)			
<input type="checkbox"/>	172PA26	9" Power Seat Elevate for use without tilt (<i>NA with QM-715 HD</i>)	\$3,456.00
<input type="checkbox"/>	172PA27	12" Power Seat Elevate for use with tilt (<i>Must be ordered with SPOT tilt, NA with QM-715 HD</i>)	\$3,663.00

Power Seating Control

Powered Seating Thru-drive Control Method.			
Thru-Drive 2+ Actuator			
<input type="checkbox"/>	172EO29	Thru Drive 2+ Actuators: (For operating Tilt, recline, lift, power legs) Includes Ctrl+5	\$3,809.00



Ctrl+5 Switchbox Options			
Ctrl+5 allows direct control to 5 power seating functions (10 total switches). Switches can also be used with assignable buttons for greater flexibility. Box includes 3-Amp USB port for charging portable devices			
<input type="checkbox"/>	172QG01	Ctrl+5 Switchbox mounted behind joystick	NC
<input type="checkbox"/>	172QG02	Ctrl+5 Switchbox mounted opposite joystick	NC
<input type="checkbox"/>	172QG03	Ctrl+5 Switchbox mounted on back cane	NC
<input type="checkbox"/>	172QG04	Ctrl+5 Switchbox mounted on seat rail	NC
<input type="checkbox"/>	172QG05	Omit Ctrl+5	(\$656.00)

Powered Seating Thru-drive Control Options			
Assignable Buttons			
"Assignable Button" feature allows for customization of joystick buttons and/or Ctrl+5 switches allowing access to various wheelchair and seating functions. When enabled, chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle.			
<input type="checkbox"/>	172PM01	Enable Assignable buttons	NC

Ctrl+5 Switch Type			
<input type="checkbox"/>	172EY03	Push buttons	NC
<input type="checkbox"/>	172EY04	Toggle switches	\$328.00

11. BACKREST AND BACK OPTIONS for 3.7 POWER RECLINE

Backrest Frame Type			
<input type="checkbox"/> 172BT89		Power Recline (Power Recline Seat only)	NC

Power Recline Backrest Preset Home Angle NC			
<input type="checkbox"/> 90° 172QK1		<input type="checkbox"/> 95° 172QK2 default	<input type="checkbox"/> 100° 172QK3

Back Type (Solid Recline)			
<input type="checkbox"/> 172BU83		Standard Curve back with 1.5" Foam cushion (Reduces Seat Depth)	NC
<input type="checkbox"/> 172JB03		JAY J3 for Power Recline w/1.5" Foam Choose type below. (Reduces Seat Depth)	NC

Back Height - Required for Standard Curved Recline Backs and Jay 3 recline backs.			
---	--	--	--

(Limited combinations of widths) Measured from pan to top of solid back

<input type="checkbox"/> 172B42		22"	NC
<input type="checkbox"/> 172B44		24"	NC
<input type="checkbox"/> 172B140		26"	NC
<input type="checkbox"/> 172B141		27"	NC
<input type="checkbox"/> 172B142		28"	NC

Push Handle Options			
<input type="checkbox"/> 172BP55		Fixed "Cane" Handles	NC
<input type="checkbox"/> 172BP56		Height Adjustable Push Handles w/2 utility hooks	\$57.00
<input type="checkbox"/> 172BP57		Omit push handles	NC

11a. JAY 3 BACK for 3.7 POWER RECLINE OPTIONS

JAY J3 Contour.			
<input type="checkbox"/> 172JP01		J3 PA Back (2.2" of Contour Depth) Available 14"-20" Wide (J3PAUTZYXX)	\$852.00
<input type="checkbox"/> 172JP02		J3 PA Back (2.2" of Contour Depth) Available 22" Wide (J3PAUTZ2FX)	\$1,812.00
<input type="checkbox"/> 172JP03		J3 PL Back (3.0" of Contour Depth) Available 14"-20" Wide (J3PLUTZYXX)	\$920.00
<input type="checkbox"/> 172JP04		J3 PDL Back (6.0" of Contour Depth) Available 14"-20" Wide (J3PDLUTZYXX)	\$987.00
<input type="checkbox"/> 172JP05		J3 PDC Back (6.0" of Contour Depth) Available 14"-20" Wide (J3PDCUTZYXX)	\$997.00
<input type="checkbox"/> 172JP06		J3 PD Back (6.0" of Contour Depth) Available 14"-20" Wide (J3PDUTZYXX)	\$997.00

JAY J3 Back Height. Backs are height adjustable			
<input type="checkbox"/> 172JR06		MT for 22" tall back. (Mid Thoracic Short 15" shell)	NC
<input type="checkbox"/> 172JR01		UT for 24" tall back. (Upper Thoracic Short 16.5" shell)	NC
<input type="checkbox"/> 172JR02		UT for 26" tall back. (Upper Thoracic Med 18" shell)	NC
<input type="checkbox"/> 172JR03		UT for 27" tall back. (Upper Thoracic Tall 19.5" shell)	NC
<input type="checkbox"/> 172JR04		SH for 28" tall back and widths 14-20" (Shoulder Height 21" shell,)	NC
<input type="checkbox"/> 172JR05		SH for 26-28" tall back and 22" wide seat (Shoulder Height 18" shell,)	NC

SPINE-ALIGN MINI FIT KIT 2

Includes: 5 wedges, 2 lat shims, 2 post/lat pelvic, 1 lg lumbar, 1 flat shim, and 1 curved

<input type="checkbox"/> 172JS01		SPINE-ALIGN MINI FIT KIT 2	\$202.00
----------------------------------	--	----------------------------	----------

JAY J3 Width


Back can be ordered 2" narrower than seat pan. Armrest will be adjusted to back width ordered. Future growth can be accomplished by replacing back with wider versions and adjusting arms outward.

<input type="checkbox"/> 172JQ01		Back width 14" (Use w/ 14-16" wide seat)	NC
<input type="checkbox"/> 172JQ02		Back width 16" (Use w/ 16-18" wide seat)	NC
<input type="checkbox"/> 172JQ03		Back width 18" (Use w/ 18-20" wide seat)	NC
<input type="checkbox"/> 172JQ04		Back width 20" (Use w/ 20-22" wide seat)	NC
<input type="checkbox"/> 172JQ05		Back width 22" (Use with 22" wide seat)	NC


12. HEADREST OPTIONS

Choose One Each

Whitmyer Headrest Mount			
<input type="checkbox"/> 172MH10		Axys Mount & Bracket	\$334.00
<input type="checkbox"/> 172MH11		Cobra Mount and Bracket	\$613.00
<input type="checkbox"/> 172MH1		Linx Mount and Bracket	\$251.00



Whitmyer Headrest Pad			
<input type="checkbox"/> 172HE65		Plush Head Support Pad 8"	\$251.00
<input type="checkbox"/> 172HE66		Plush Head Support Pad 10"	\$251.00
<input type="checkbox"/> 172HE78		Plush Head Support Pad 14"	\$251.00
<input type="checkbox"/> 172HE82		Adjust-A-Plush Head Support Pad Narrow	\$387.00
<input type="checkbox"/> 172HE83		Adjust-A-Plush Head Support Pad Standard	\$387.00



13. POSTURAL SUPPORTS

LATERAL THORACIC SUPPORT (Swing Away)

Available on all Solid Backs.

<input type="checkbox"/> LAT 4X4 S/A MOUNTS W/ 1" OFFSET - PR \$434.00		
172ZV14	LATERAL 4X4 W/ 1" OFFSET - PAIR	\$134.00
172ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$300.00
<input type="checkbox"/> LAT 4X5 S/A MOUNTS W/ 1" OFFSET - PR \$434.00		
172ZV15	LATERAL 4X5 W/ 1" OFFSET - PAIR	\$134.00
172ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$300.00
<input type="checkbox"/> LAT 4X6 S/A MOUNTS W/ 1" OFFSET - PR \$434.00		
172ZV16	LATERAL 4X6 W/ 1" OFFSET - PAIR	\$134.00
172ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$300.00
<input type="checkbox"/> LAT 6X4 S/A MOUNTS W/ 1" OFFSET - PR \$434.00		
172ZV17	LATERAL 6X4 W/ 1" OFFSET - PAIR	\$134.00
172ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$300.00
<input type="checkbox"/> LAT 6X5 S/A MOUNTS W/ 1" OFFSET - PR \$434.00		
172ZV18	LATERAL 6X5 W/ 1" OFFSET - PAIR	\$134.00
172ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$300.00
<input type="checkbox"/> LAT 6X6 S/A MOUNTS W/ 1" OFFSET - PR \$434.00		
172ZV19	LATERAL 6X6 W/ 1" OFFSET - PAIR	\$134.00
172ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$300.00

CHEST STRAP

<input type="checkbox"/> 172ZR03	Padded Chest Strap	\$300.00
----------------------------------	--------------------	----------

Swing Away Knee Adductor Supports

Available only with Fixed and Power Center Mount Legrests.

Not available with Dual Post Arm Rests or 2-point Locking Cantilever Arms

Knee Adductor Supports - LEFT (choose Pad and Arm) \$528.00		
Pad Size		
<input type="checkbox"/> 172KB5	Knee Adductor Pad SM- Left	\$228.00
<input type="checkbox"/> 172KB7	Knee Adductor Pad LG- Left	\$228.00
S/A Arm Size		
<input type="checkbox"/> 172XX04	Swing Away Short - Left	\$300.00
<input type="checkbox"/> 172XX06	Swing Away Long - Left	\$300.00
Knee Adductor Supports - RIGHT (choose Pad and Arm) \$528.00		
Pad Size		
<input type="checkbox"/> 172KB6	Knee Adductor Pad SM - Right	\$228.00
<input type="checkbox"/> 172KB8	Knee Adductor Pad LG - Right	\$228.00
S/A Arm Size		
<input type="checkbox"/> 172XX05	Swing Away Short - Right	\$300.00
<input type="checkbox"/> 172XX07	Swing Away Long - Right	\$300.00

Long Thigh Support w/Removable Mount

Available with Cantilever Arms, NA with Manual Recline Back, Dual Post, Single Post or Cantilever w/ 2 pt Locking arms.

<input type="checkbox"/> SPOT Long Thigh Support - Left \$434.00		
172TV2	Long Thigh Support - Left	\$134.00
172TX1	Removable Mount - Left	\$300.00
<input type="checkbox"/> SPOT Long Thigh Support - Right \$434.00		
172TV1	Long Thigh Support - Right	\$134.00
172TX2	Removable Mount - Right	\$300.00

14. ARMREST AND ARMPAD OPTIONS

Armrest for Power Recline Seating.			
(Armrest width can be adjusted 2" narrower than seat)			
<input type="checkbox"/> 172EG26		STD Cantilever/Reclining Arm 10-14" Height Range	NC
<input type="checkbox"/> 172EG27		Cantilever/Reclining Arm w/ 2 Point Locking Post. (10"-12" Ht Range.)	\$334.00
<input type="checkbox"/> 172EG28		Cantilever/Reclining Arm w/ 2 Point Locking Post. (11"-13" Ht Range.)	\$334.00

Reclining Armpad Options			
(Armpad can be angled inward 15 degrees)			
Standard Reclining Armpad			
Left Armpad NC		Right Armpad NC	
<input type="checkbox"/> 172AE05	St Reclining Arm Pad (16")	<input type="checkbox"/> 172AB05	
Comfort Reclining Armpad (Foam)			
Left Armpad (\$104.00)		Right Armpad (\$104.00)	
<input type="checkbox"/> 172AE22	Comfort 2" wide Reclining Arm Pad	<input type="checkbox"/> 172AB22	
<input type="checkbox"/> 172AE23	Comfort 4" wide Reclining Arm Pad	<input type="checkbox"/> 172AB23	
ON-CHAIR Comfort Reclining Armpad (Gel)			
Left Gel Armpad (\$197.00)		Right Gel Armpad (\$197.00)	
<input type="checkbox"/> 172AE25	Comfort 2" wide Reclining Arm Pad	<input type="checkbox"/> 172AB25	
<input type="checkbox"/> 172AE26	Comfort 4" wide Reclining Arm Pad	<input type="checkbox"/> 172AB26	
ON-CHAIR AES Trough/Elbow Block Reclining Armpads (Gel)			
Left Gel Armpad (\$274.00)		Right Gel Armpad (\$274.00)	
<input type="checkbox"/> 172AE27	Trough / Elbow Block Medium (Gel) (3.5" x 13") Reclining Arm Pad	<input type="checkbox"/> 172AB27	
<input type="checkbox"/> 172AE28	Trough / Elbow Block Large (Gel) (4" x 14.5") Reclining Arm Pad	<input type="checkbox"/> 172AB28	

15. HANGERS AND FOOTRESTS

Swing-away Hangers

Non-Powered Swing Away Legrests			
<i>Sold as Pairs</i>			
<input type="checkbox"/>	172H2	70° Swing Away (B)	NC
<input type="checkbox"/>	172H138	Manual ELR S/A (B)	\$374.00

Power Elevating Legrests (2 Actuators)			
<input type="checkbox"/>	172H126	Swing Away Power ELR w/ Articulation Pair (B footplates)	\$3,177.00

Footplates For Swing Away Hangers			
<input type="checkbox"/>	B - Adult Angle Adjust w/Heel Loop for use with 70° hanger and ELR (Pair)		\$217.00
	172F282	Adult Angle Adjustable Footplate	\$217.00
	172HL10	Heel Loop	NC
<input type="checkbox"/>	B - Adult Angle Adjust w/Heel Loop for use with 90° hanger (Pair)		\$217.00
	172F285	Adult Angle Adjustable Footplate	\$217.00
	172HL10	Heel Loop	NC
<input type="checkbox"/>	C - Kids Angle Adjust w/Leg Strap for use with 70° hanger (Pair)		\$217.00
	172F283	Kids Angle Adjustable Footplate	\$217.00
	172LS13	Leg Strap	NC
<input type="checkbox"/>	C - Kids Angle Adjust w/Leg Strap for use with 80° and 90° hanger(Pair)		\$217.00
	172F286	Kids Angle Adjustable Footplate	\$217.00
	172LS13	Leg Strap	NC
<input type="checkbox"/>	172F279	Omit Footplate and extension tube	NC

Omit Hangers			
<input type="checkbox"/>	172H144	Omit Hangers	NC

Hanger Extensions (See range listed in matrix below)			
<i>NOTE: Not all extension tube will be available; depending on the seat height, some extension tube options may not offer enough ground clearance.</i>			
<input type="checkbox"/>	172E16	High Mount Hanger	\$90.00
<input type="checkbox"/>	172E1	Short	NC
<input type="checkbox"/>	172E2	Medium	NC
<input type="checkbox"/>	172E5	Long	NC

Footrest Accessories			
<input type="checkbox"/>	1725A01	Toe Loop (Pair) Only avail on (opt B).	\$49.00
<input type="checkbox"/>	172LS14	Extra Leg Strap	\$57.00

Lower Leg Length Matrix (measured seat pan to the top of the footplate surface)				
Hanger Type	Standard	Short	Medium	Long
70° Swing Away		13.5"- 16"	16"-18.5"	18"-20.5"
Manual ELR S/A		12.5"-17"	14.5"-19"	17"-21"
Power ELR S/A		12"-14"	14"-17"	16.5"-19"
Fixed Cntrmnt Recline	3"-19"			
Power Cntrmnt	3"-19"			


Centermount Footplates

Non-Powered Centermount Legrests			
<input type="checkbox"/>	172H128	Fixed Centermount for Power Recline (F)	NC

Power Centermount (2 Actuators)			
<input type="checkbox"/>	172H127	Power Centermount ELR w/7.5" Power Articulation (F) Centermount footplates plus Custom Options Below	\$3,169.00

Power Centermount Footplate Extension Option (For Thru-Drive Operation)			
<input type="checkbox"/>	172IE01	Enable 7.5" Powered footplate extension for transfer assist. (Electronics programmed to allow footplates to be extended <u>with or without</u> legrest being raised)	\$1,639.00
<input type="checkbox"/>	172IE02	Disable Footplate extension for transfers. (Electronics programmed to <u>ONLY</u> extend footplates when legs are raised)	NC

Centermount Footplates			
<input type="checkbox"/>	F - Centermount Dual Footplate and Heel Block (Pair) (Individually Height (3"-19") and Angle Adjustable)		\$262.00
	172F315	CMT Angle Adjustable Footplates	\$262.00
	172HL14	Padded and Adjustable Heel Block	NC
<input type="checkbox"/>	I - Angle adj Centermount Footplate pair		\$218.00
	172F125	Angle Adjustable Footplates	\$218.00

Custom Centermount Options 			
<input type="checkbox"/>	172DI26	Extra Large Single Footplate for Centermount Footrest. Size 16"W x 14"D	\$328.00
<input type="checkbox"/>	172DI27	Gel Foot Pads for XL Single Centermount Footplate.	\$328.00
<input type="checkbox"/>	172DI28	Extra Large Pair Curved Calf Pads for Center Mount Legrests. Fully adjustable	\$328.00

16. ACCESSORIES

Wheelchair Tray Table			
<input type="checkbox"/>	1727A01	Tray Table SM Tuffak 12"-14" Pwr	\$351.00
<input type="checkbox"/>	1727A02	Tray Table MED Tuffak 15"-17" Pwr	\$351.00
<input type="checkbox"/>	1727A03	Tray Table LRG Tuffak 18"-20" Pwr	\$351.00

Backpack			
<input type="checkbox"/>	172A1	Black Utility	\$103.00

Positioning Belt - Select for any style seat			
<input type="checkbox"/>	172EJ01	Auto Style Buckle	NC
<input type="checkbox"/>	172EJ02	Belt 2" Aircraft Buckle	\$173.00
<input type="checkbox"/>	172EJ03	Belt 2" Aircraft Buckle Padded	\$188.00

Accessories			
<input type="checkbox"/>	1728A09	Ram Phone Mount on Round JS tube	\$71.00
<input type="checkbox"/>	1728A10	Ram Tablet Mount on Round JS tube	\$284.00
<input type="checkbox"/>	1728A11	Ram Cup holder	\$71.00
<input type="checkbox"/>	1728A12	Switch-It Remote Stop Switch	\$1,273.00

17c. JAY CUSHIONS

Seat width must match selection in section 7. Seat depth can be +/- 2" of selection in section 7

JAY Fusion Cushion				
<input type="checkbox"/>	172JU01	JFUSION1616	JAY Fusion 16" x 16"	\$976.00
<input type="checkbox"/>	172JU02	JFUSION1618	JAY Fusion 16" x 18"	\$976.00
<input type="checkbox"/>	172JU03	JFUSION1816	JAY Fusion 18" x 16"	\$976.00
<input type="checkbox"/>	172JU04	JFUSION1818	JAY Fusion 18" x 18"	\$976.00
<input type="checkbox"/>	172JU05	JFUSION1820	JAY Fusion 18" x 20"	\$976.00
<input type="checkbox"/>	172JU06	JFUSION2018	JAY Fusion 20" x 18"	\$976.00
<input type="checkbox"/>	172JU07	JFUSION2020	JAY Fusion 20" x 20"	\$976.00

JAY X2 Cushion				
<input type="checkbox"/>	172JU22	JX21616	JAY X2 16" x 16"	\$881.00
<input type="checkbox"/>	172JU23	JX21618	JAY X2 16" x 18"	\$881.00
<input type="checkbox"/>	172JU24	JX21816	JAY X2 18" x 16"	\$881.00
<input type="checkbox"/>	172JU25	JX21818	JAY X2 18" x 18"	\$881.00
<input type="checkbox"/>	172JU26	JX21820	JAY X2 18" x 20"	\$881.00
<input type="checkbox"/>	172JU27	JX22020	JAY X2 20" x 20"	\$881.00

JAY Union Cushion				
<input type="checkbox"/>	172JU15	JU1616	JAY UNION 16" x 16"	\$714.00
<input type="checkbox"/>	172JU16	JU1618	JAY UNION 16" x 18"	\$714.00
<input type="checkbox"/>	172JU17	JU1816	JAY UNION 18" x 16"	\$714.00
<input type="checkbox"/>	172JU18	JU1818	JAY UNION 18" x 18"	\$714.00
<input type="checkbox"/>	172JU19	JU1820	JAY UNION 18" x 20"	\$714.00
<input type="checkbox"/>	172JU20	JU2018	JAY UNION 20" x 18"	\$753.00
<input type="checkbox"/>	172JU21	JU2020	JAY UNION 20" x 20"	\$753.00

JAY Ion Cushion				
<input type="checkbox"/>	172JU08	JION1616	JAY ION 16" x 16"	\$613.00
<input type="checkbox"/>	172JU09	JION1618	JAY ION 16" x 18"	\$613.00
<input type="checkbox"/>	172JU10	JION1816	JAY ION 18" x 16"	\$613.00
<input type="checkbox"/>	172JU11	JION1818	JAY ION 18" x 18"	\$613.00
<input type="checkbox"/>	172JU12	JION1820	JAY ION 18" x 20"	\$613.00
<input type="checkbox"/>	172JU13	JION2018	JAY ION 20" x 18"	\$652.00
<input type="checkbox"/>	172JU14	JION2020	JAY ION 20" x 20"	\$652.00

JAY Custom Cushion

*To order any **JAY cushions not listed above**, please see the JAY price list. Cushion name or part number must be filled out w/ desired width and depth*

<input type="checkbox"/>	172JC20	JAY Bundled Cushion	Cushion Name or Part #	<input style="width: 100px;" type="text"/>	Width:	<input style="width: 50px;" type="text"/>	Depth:	<input style="width: 50px;" type="text"/>	Price List
--------------------------	---------	---------------------	------------------------	--	--------	---	--------	---	------------

18. CUSTOM OPTIONS

Looking for more options? Contact our custom "Built-4-Me" department.

To Inquire about any unique customization for Quickie products, contact:
Phone: (800) 263-3390 Fax: (800) 561-5834



Below are common upgrades that utilize **Ready Made Quotes**. These options can be quickly quoted to add a custom touch to your QM-7 power chair.

<input type="checkbox"/>	O2 Holder on Power Recline	Call for quote
--------------------------	----------------------------	----------------

Some Restrictions Apply

