



Zippie Xperience 2

July 2018

EIPW23



ACCOUNT INFORMATION

Dealer Contact Name:	_____	P.O. # or Quote:	_____
Email Address:	_____	User Name or Number:	_____
Phone Number:	_____	Fax Number:	_____
Bill to (Name):	_____	User Weight:	_____
Bill to (Location):	_____	User Height:	_____
Ship to (Name):	_____	Seat Cushion Thickness:	_____
Address:	_____	Knee to Heel Length: Left:	_____
City:	_____	Right:	_____
Province:	_____	Postal Code:	_____

The order cannot be processed without the following:

- A user name or number, the client's measurements and weight
- Sunrise Medical is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Sunrise Medical summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

STANDARD OPTIONS

Zippie Power Base, 6-wheel independent wheel suspension, an additional seat suspension system, tie-down hardware for strap systems, 8A offboard charger, 6" front casters on MWD or 8" on RWD, and standard auto style buckle seat positioning strap, 6.5 mph motors, **Weight Capacity = 300 lbs.**

(MID-WHEEL CONFIGURATION) ZIPPIE XPERIENCE 2 POWER WHEELCHAIR

<input type="checkbox"/> XPERPED2	ZIPPIE XPERIENCE 2 PED Power Wheelchair	\$8,545.00
-----------------------------------	---	------------

USER WEIGHT PACKAGES

<input type="checkbox"/> WCA3	300 lbs. Maximum user weight capacity	STD
-------------------------------	---------------------------------------	-----

STANDARD FRAME COLORS

<input type="checkbox"/> FCBLK	Black	STD
<input type="checkbox"/> FCSBM	Matte Black	STD
<input type="checkbox"/> FCGMC	Gun Metal (Charcoal)	STD
<input type="checkbox"/> FCSVM	Silver Metallic	STD
<input type="checkbox"/> FCCBU	Candy Blue	STD
<input type="checkbox"/> FCCRE	Candy Red	STD
<input type="checkbox"/> FCMPL	Midnight Purple	STD
<input type="checkbox"/> FCRZB	Razzleberry	STD
<input type="checkbox"/> FCMWT	Matte White	STD

OPTIONAL FRAME COLORS (May Increase lead times)

<input type="checkbox"/> FCMYG	Military Green (New)	N/C
<input type="checkbox"/> FCBOR	Burnt Orange Metallic	N/C
<input type="checkbox"/> FCBAO	Black Opal Metallic	N/C
<input type="checkbox"/> FCOPM	Sparkle Orange	N/C
<input type="checkbox"/> FCMYT	Military Tan (New)	N/C
<input type="checkbox"/> FCPKP	Pearl Pink	N/C
<input type="checkbox"/> FCHDD 1	Desert Digital Camo Hydro Graphics	\$340.00
<input type="checkbox"/> FCHWL 1	Woodland Camo Hydro Graphics	\$340.00
<input type="checkbox"/> FCHCF 1	Carbon Fiber Hydro Graphics (Matte Black base)	\$340.00

Note 1 Includes matching accessories with the Hydro Graphics option. The matching accessories option is not required with these colours.

MATCHING ACCESSORIES TO THE FRAME COLOR

<input type="checkbox"/> FC**-M	Matching Accessories to the Frame Color (Wheels, front access door, headrest mount and back shell when selected)	\$256.00
---------------------------------	--	----------

This option is included at no charge with the Hydro Graphics Frame Color option.

TRANSPORT OPTION

<input type="checkbox"/> TTD4	2	4-Point wheelchair transport brackets	STD
<input type="checkbox"/> WC19-ISO	New	WC-19 and ISO 7176-19 Occupied Transit Package	N/C

To confirm if your overall wheelchair configuration is WC-19 and ISO 7176-19 compliant, please refer to this website: <http://wc-transportation-safety.umtri.umich.edu/crash-tested-product-lists/wheelchairs>
Other restrictions may apply, please contact CSR for details. This option is only available on new power wheelchairs as an "on-chair" option. It is not available to install on a delivered wheelchair

Note 2 4-Point wheelchair transport brackets have not yet received WC-19 crash-testing approval.
Use these transport brackets only to secure an unoccupied wheelchair during transport.

MOTOR PACKAGE

<input type="checkbox"/> MP65	6.5 MPH 4-Pole Motor Package / with COG Release	STD
<input type="checkbox"/> MP85	8.5 MPH 4-Pole High Speed Motor Package (For Rear-Wheel configuration only-Not available with VR2 electronics ECV902A)	N/C

ASSISTIVE DRIVING PACKAGE

<input type="checkbox"/> INTELLIDR	3, 4	Intelli-Drive	\$1,694.00
------------------------------------	------	---------------	------------

Note 3 Intelli-Drive technology provides increased tracking capabilities with enhanced traction and steering precision. This is ideal for those with marginal hand function or switch control.
Note 4 Requires the R-net PM 120 Amp Expandable Controller (ECR120)

▼ BATTERY CHARGER			
<input type="checkbox"/>	BCOB	8A Offboard Battery Charger	N/C
<input type="checkbox"/>	BCLES8OPT	Lester 8 amp Offboard Battery Charger	\$383.00

▼ BATTERY TRAY AND FRAME SIZE			
<input type="checkbox"/>	FS22B	Battery Tray for Group 22NF Batteries (Total frame width = 23 3/8")	STD
<input type="checkbox"/>	FS2424B	CUSTOM Battery Tray for Group 24 Batteries with a Compact 24" frame width (Mid-Wheel Drive Configuration only)	N/C
Please Select One Option for the Hybrid RWD configuration: <input type="checkbox"/> Front Caster Arms <input type="checkbox"/> Flared out (STD) <input type="checkbox"/> Straight (Optional)			
Note 5 The 23 3/8" compact frame is not compatible with most of the Centermount Front Rigging Options due to Mechanical conflicts.			
Note 6 The 90 degree Power Centermount Elevating & Articulating Foot Platform (ADDPF) is not available with this option.			
Note 7 The Straight Front Arm option is not available with any of the Centermount Front Rigging Options due to Mechanical conflicts.			

▼ BATTERIES			
<input type="checkbox"/>	BAT22NFPR	GEL Batteries – Group 22NF – (Each-on chair)	\$658.00
<input type="checkbox"/>	BAT24PR	GEL Batteries – Group 24 – (Each-on chair)	\$881.00

▼ CHOOSE TIRE COLOR (APPLIES TO ALL TIRES OR CASTERS)			
<input type="checkbox"/>	BLACK (New)	<input type="checkbox"/>	GREY

▼ DRIVE WHEEL			
<input type="checkbox"/>	DW14PT	14" Cast Wheel Pneumatic Tire	STD
<input type="checkbox"/>	DW14PI	14" Cast Wheel Pneumatic with Flat-Free inserts	\$256.00
<input type="checkbox"/>	DW144PT	14" x 4" Extra Wide Cast Wheel Pneumatic Tire	\$279.00
<input type="checkbox"/>	DW144PI	14" x 4" Extra Wide Cast Wheel Pneumatic with Flat-Free inserts	\$279.00
<input type="checkbox"/>	DW12PT	12 1/2" Cast Wheel Pneumatic Tire	N/C
<input type="checkbox"/>	DW12PI	12 1/2" Cast Wheel Pneumatic with Flat-Free inserts	\$256.00
Note 8 The 14" x 4" Extra Wide Wheel increases the overall width of the Base by 2 1/2" The option to match the wheels to the shroud colour (FC...M) is not available with the 14" x 4" Extra Wide Wheels. The 14" x 4" option is not available with the VR2 electronics (ECV902A)			
Note 9 12 1/2" Drive wheels reduce overall seat to floor height by 3/4", lowers ground clearance and speed, but increases motor torque			

▼ FRONT CASTERS (Mid Wheel Configuration)			
<input type="checkbox"/>	FC6S	6" Solid Casters	STD
<input type="checkbox"/>	FC8PTM	8" Pneumatic Casters	\$173.00
<input type="checkbox"/>	FC8PIM	8" Pneumatic Casters with Flat-Free inserts	\$256.00
Note 10 The 8" casters are NOT compatible with the 12 1/2" Drive Wheel configuration (DW12PT and DW12PI)			

▼ FRONT CASTERS (Rear-Wheel Hybrid Configuration)			
<input type="checkbox"/>	FC8PTR	8" Pneumatic Casters	STD
<input type="checkbox"/>	FC8PIR	8" Pneumatic Casters with Flat-Free inserts	\$256
<input type="checkbox"/>	FCcharge with the Hydro Graphics F	9" Pneumatic Casters	\$173
<input type="checkbox"/>	FC9PI	9" Pneumatic Casters with Flat-Free inserts	\$204
Note 11 The 9" casters will increase the minimum seat to floor height by 1/2" and are NOT compatible with the 12 1/2" Drive wheel configuration (DW12PT and DW12PI) The 9" casters are NOT compatible with many front rigging options. Please see details in the front rigging sections			

▼ SEAT TO FLOOR HEIGHT (Approximate measurement to the top of seat pan ± 1/4")

► Options for a Seat Only									
<input type="checkbox"/>	STFH155	12, 13, 15	15 1/2"	(Possible Seat angle of 0° only)	<input type="checkbox"/>	STFH1825	13	18 1/4"	(Possible Seat angle of 0° and 3°)
<input type="checkbox"/>	STFH1625	13, 15	16 1/4"	(Possible Seat angle of 0° only)	<input type="checkbox"/>	STFH1925	13	19 1/4"	(Possible Seat angle of 0°, 3° and 6°)
<input type="checkbox"/>	STFH1725	13, 15	17 1/4"	(Possible Seat angle of 0° and 3°)	<input type="checkbox"/>	STFH2025	13	20 1/4"	(Possible Seat angle of 0°, 3°, 6°, 9°)

► Options for a Power Tilt or a STD Power Elevating Seat									
<input type="checkbox"/>	STFH155	12, 14, 15	15 1/2"	(Possible Seat angle of 0° only)	<input type="checkbox"/>	STFH1775	14	17 3/4"	(Possible Seat angle of 0° only)
<input type="checkbox"/>	STFH1625	14, 15	16 1/4"	(Possible Seat angle of 0° only)	<input type="checkbox"/>	STFH1825	14	18 1/4"	(Possible Seat angle of 0° only)
<input type="checkbox"/>	STFH1725	14	17 1/4"	(Possible Seat angle of 0° only)	<input type="checkbox"/>	STFH1875	14	18 3/4"	(Possible Seat angle of 0° only)

► Options for a Combo - Power Tilt with a STD Power Elevating Seat									
<input type="checkbox"/>	STFH175	12, 14	17 1/2"	(Possible Seat angle of 0° only)	<input type="checkbox"/>	STFH1825	14	18 1/4"	(Possible Seat angle of 0° only)

- Note 12 Requires the 12 1/2" Drive wheels (DW12PT or DW12PI). This lowers ground clearance and speed, but increases motor torque.
 Note 13 Seat to Floor Heights are taken with the Mid-Wheel Drive Configurations. Please remove 1/4" for the Rear-Wheel Drive Height.
 Note 14 Seat to Floor Heights are taken with the Mid-Wheel Drive Configurations. Please Add 1/4" for the Rear-Wheel Drive Height.
 Note 15 The Adult Heavy Duty Swing-Away footrests (SOAAF) are not available with this seat to floor height on the RWD configuration only.

▼ PRE-SET SEAT FRAME ANGLE

<input type="checkbox"/>	SFE-0	0° Pre-set seat frame seat angle (Must be selected when using power options)	STD
<input type="checkbox"/>	SFE-3	3° Pre-set seat frame seat angle (Only available with a seat configuration without power options)	N/C
<input type="checkbox"/>	SFE-6	6° Pre-set seat frame seat angle (Only available with a seat configuration without power options)	N/C
<input type="checkbox"/>	SFE-9	9° Pre-set seat frame seat angle (Only available with a seat configuration without power options)	N/C
Note 16 Please refer to the special notes in the above seat to floor height section to confirm seat angle possibilities. The Light Package is not compatible with the EPSEL option.			

▼ LIGHTING PACKAGE

<input type="checkbox"/>	EOLP	LED Light Package (Includes front beams, side flashers, rear tail light with flashers & brake light)	\$1,131.00
Note 17 Requires the upgraded Joystick with Light Operation or the OMNI (ECVJ6KL, ECELJL, ECRJCL, ECRJCTPL, ECRJC35B or EDIB).			

▼ NON-EXPANDABLE ELECTRONICS - PG DRIVES VR2 CONTROLLER AND JOYSTICKS

<input type="checkbox"/>	ECV902A	VR2 90 Amp Drive Controller	STD
<input type="checkbox"/>	ECVJ4K	18 VR2 Standard 4-Key Joystick	STD
<input type="checkbox"/>	ECVJ6K	18 VR2 6-Key Joystick (Allows for a single power option to operate through this joystick)	N/C
<input type="checkbox"/>	ECVJ9KL	18 VR2 9-Key Joystick with Lighting Operation (Allows for a single power function to operate through this joystick)	\$168.00

Note 18 Requires the VR2 90A Drive Control (ECV902A)
VR2 Electronics can only accommodate a single power seat function. The only choices are Power Tilt or Elevating seat

▼ EXPANDABLE ELECTRONICS - PG DRIVES R-NET CONTROLLERS

<input type="checkbox"/>	ECEL90	R-net EL 90 Amp Expandable Controller	\$562.00
<input type="checkbox"/>	ECR120	R-net PM 120 Amp Expandable Controller	\$1,131.00
<input type="checkbox"/>	ECRNH	19 Harness for Expandable Controller	N/C

Note 19 Must be ordered with the expandable controller ECEL90 or ECR120 to interface with Joysticks or OMNI screen

▼ EXPANDABLE ELECTRONICS - PG DRIVES R-NET JOYSTICKS

<input type="checkbox"/>	ECO	Omit Joystick (For use with an OMNI Display)	N/C
<input type="checkbox"/>	ECELJ	20 R-net (LED) Joystick	STD
<input type="checkbox"/>	ECRJJC	20 R-net Joystick with Color screen	\$168.00
► JOYSTICKS WITH LIGHTING OPERATION			
<input type="checkbox"/>	ECELJL	20 R-net (LED) Joystick with Lighting Operation	\$168.00
<input type="checkbox"/>	ECRJCL	20 R-net Joystick with Color screen and Lighting Operation	\$222.00
<input type="checkbox"/>	ECRJCTPL	20 R-net Joystick with Color screen w/ On & Off/ Mode switch & speed pot and Lighting Operation	\$334.00
► ADVANCED JOYSTICK			
<input type="checkbox"/>	ECRJJC35	20 R-net Advanced Joystick (3.5" Screen). Lighting Operation, paddle switches, shortcut keys and stereo jacks	\$1,202.00
<input type="checkbox"/>	ECRJJC35B	20 R-net Advanced Bluetooth Joystick (3.5" Screen). Lighting Operation, paddle switches, shortcut keys and stereo jacks	\$1,202.00

Note 20 This option requires the expandable controller ECEL90 or ECR120

▼ JOYSTICK MOUNT

<input type="checkbox"/>	JSMF	Fixed Joystick Mount	N/C
<input type="checkbox"/>	JMSA	Swing Away Joystick Mount - Height adjustable, Angle adjustable and Multi-Axis	\$412.00
		<input type="checkbox"/> Set-up with standard installation to swing-away sideways	<input type="checkbox"/> Set-up at 35° to swing-away sideways and downwards
<input type="checkbox"/>	JSMCM	Midline Mounted Swing Away Joystick Mount	\$412.00

▼ JOYSTICK MOUNT SUPPORT

<input type="checkbox"/>	JSMFH	Joystick mount support Bracket	STD
--------------------------	-------	--------------------------------	-----

▼ JOYSTICK HAND MOUNT

<input type="checkbox"/>	HML	Install Joystick on Left hand side	<input type="checkbox"/>	HMR	Install Joystick on Right hand side
--------------------------	-----	------------------------------------	--------------------------	-----	-------------------------------------

▼ JOYSTICK HANDLE

<input type="checkbox"/>	JSH	Standard Joystick Gimble	STD
<input type="checkbox"/>	JBAT	Stick Handle 3" in Length	\$101.00
<input type="checkbox"/>	JBH	Ball Handle	\$101.00
<input type="checkbox"/>	JFBH	Foam Ball Handle	\$101.00
<input type="checkbox"/>	JTH	T-Handle	\$101.00
<input type="checkbox"/>	JMUSH	2.5" Mushroom Handle	\$122.00
<input type="checkbox"/>	JGP	Goalpost Handle	\$101.00

▼ ATTENDANT DRIVE CONTROL OPTION

<input type="checkbox"/>	EADC	Attendant Hand Control for R-net EL & PM Controls	\$1,131.00
--------------------------	------	---	------------

▼ ATTENDANT JOYSTICK HAND MOUNT

<input type="checkbox"/>	EADCL	Install Attendant Joystick on Rear Left hand side	<input type="checkbox"/>	EADCR	Install Attendant Joystick on Rear Right hand side
--------------------------	-------	---	--------------------------	-------	--

▼ ENHANCED DISPLAY & DIGITAL INTERFACE BOX (Included with specialty controls)

<input type="checkbox"/>	EDIB	21 OMNI (Enhanced Display & Digital Interface Box with Infrared)	\$2,157.00
--------------------------	------	--	------------

Note 21 Allows the switch operation of alternative devices through driver controls. Requires one of the R-net fully Expandable Controllers (ECEL90 or ECR120)

▼ ENHANCED DISPLAY MOUNT

<input type="checkbox"/>	EDSAM	User Swing Away Mount	<input type="checkbox"/>	EDML (Install Display on LH side)	<input type="checkbox"/>	EDMR (Install Display on RH side)	N/C
--------------------------	-------	-----------------------	--------------------------	-----------------------------------	--------------------------	-----------------------------------	-----

▼ ENVIRONMENTAL CONTROL MODULE

<input type="checkbox"/> ECU	22	Environment control module - (Output Module)	\$1,131.00
		Note 22 Requires one of the fully expandable R-net controllers (ECEL90 or ECR120)	

▼ BLUETOOTH MOUSE EMULATOR (Includes USB Dongle)

<input type="checkbox"/> EDBT	23	Bluetooth Mouse Emulator (USB Dongle included)	\$1,131.00
		Note 23 Requires one of the fully expandable R-net controllers (ECEL90 or ECR120) . Not required if the ECRJC35B joystick has been selected Not limited to only one selection. Any two can be installed on the same chair, but only one can be used at a time.	
		<input type="checkbox"/> Windows and MAC OS	
		<input type="checkbox"/> Android Devices	
		<input type="checkbox"/> iPhone, iPad and iPod	

▼ POWER OPTIONS & SEATING

<input type="checkbox"/> SAABA		Amyseat Rehab Seat - Seat Frame only	\$881.00
<input type="checkbox"/> SAPTX	24	CG Power Tilt	\$7,385.00
<input type="checkbox"/> SAPRX	24, 25, 26	Power Recline with Integrated Shear Reduction (ISR)	\$9,090.00
<input type="checkbox"/> SAPTPRX	24, 25, 26	CG Power Tilt and Power Recline with Integrated Shear Reduction (ISR)	\$13,070.00
<input type="checkbox"/> SAPEX	24	11" Power Elevating Seat Module for transfer assist	\$3,663.00
<input type="checkbox"/> ADDPESX		Module - Add an 11" Power Elevating Seat to either SAPTX, SAPRX or SAPTPRX	\$3,404.00
<input type="checkbox"/> ADDMR	25, 26	Module - Manual Recline with Integrated Shear Reduction (ISR)	\$1,131.00
		Note 24 These options include the fully adjustable Amyseat (SAABA)	
		Note 25 This option is not available with the Pediatric Seat Frame (SOWN11-16)	
		Note 26 It is recommended to use the Contour Back Support option with the power recline (STDSB or STDSB/CUSTOM).	

▼ ELECTRONICS - SINGLE OR DUAL (1 OR 2 ACTUATORS) POWER SEAT FUNCTION CONTROLLERS

<input type="checkbox"/> EPSV	27, 28	Tilt or Elevating seat activation through 6-key VR2 Hand Control (ECVJ6K or ECVJ9KL)	N/C
<input type="checkbox"/> EPSEL	27, 28	Tilt or Elevating seat activation through RNET 90 Amp Expandable Controller (ECEL90)	N/C
<input type="checkbox"/> SAC		Single (1) Power Seat Function Controller to Interface through R-net Electronics	N/C
<input type="checkbox"/> DAC		Dual (2) Power Seat Function Controller to Interface through R-net Electronics	\$562.00
With the SAC or DAC, please select one activation and position from the list below			
<input type="checkbox"/> SWTJSS	29, 30	Single actuator activation through R-net Joystick or Specialty Controls EPSR2A	\$1,131.00
<input type="checkbox"/> SWTJSD	29, 31	Dual actuator activation through R-net Joystick or Specialty Controls EPSRA	\$1,131.00
<input type="checkbox"/> SWDUALACT		Dual Activation (Allows to operate all seat functions with external switches and through the joystick)	\$350.00
<input type="checkbox"/> SWTOGGLE	30 or 31	Single Toggle Switch (Single Mode: DAC - Dual Mode: SAC)	<input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation N/C
<input type="checkbox"/> SWPUSH1	30	Single Push Button Switch (Single Mode: SAC)	<input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation N/C
<input type="checkbox"/> SWPUSH2	30 or 31	Dual Push Button Switch (Single Mode: DAC - Dual Mode: SAC)	<input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation N/C
Note 27 Allows for the activation of a single power option through the VR2 or R-net 90 Amp Controllers. Only the Power Tilt or Elevating seat options are available with this feature.			
Note 28 These options do not require the SAC controller as they connect via the chair's power module directly			
Note 29 Requires the R-net 90 Amp (ECEL90) or the 120 Amp (ECR120) Fully Expandable Controllers			
Note 30 Requires the Single (1 Actuator) Power Seat Function (SAC)			
Note 31 Requires the Dual (2 Actuators) Power Seat Function (DAC)			

▼ ELECTRONICS - MULTIPLE (3 to 6 ACTUATORS) POWER SEAT FUNCTION CONTROLLER

<input type="checkbox"/> MAC		Multiple (3 to 6 Actuators) Power Seat Function Controller to Interface through R-net Electronics	\$2,720.00
With the MAC, please select one activation and position from the list below			
<input type="checkbox"/> SWTJSM	32, 33	Multiple Actuator activation through R-net Joystick or Specialty Controls MAC-EPSR	STD
<input type="checkbox"/> SWDUALACT		Dual Activation (Allows to operate all seat functions with external switches and through the joystick)	\$350.00
<input type="checkbox"/> SWMULTI	32	4-Way Toggle Switch (Up to 4 Activations)	<input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation N/C
<input type="checkbox"/> SWMULTIT	32	Triple Toggle Switch (Up to 5 Activations)	<input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation N/C
<input type="checkbox"/> SWPUSH4	32	4-Way Push Button Switch	<input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation N/C
<input type="checkbox"/> SWPUSH5	32	5-Way Push Button Switch	<input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation N/C
Note 32 Requires the Multiple (3 to 6 Actuators) Power Seat Function Controller (MAC)			
Note 33 Requires the R-net 90 Amp (ECEL90) or 120 Amp (ECR120) Fully Expandable Controller, and the multiple Actuator Controller (MAC)			

▼ SWITCH ACTIVATION HAND MOUNT

Switch Installation:	<input type="checkbox"/> TSML (Install Switch on LH side)	<input type="checkbox"/> TSMR (Install Switch on RH side)	N/C
----------------------	---	---	-----

▼ PEDIATRIC ADJUSTABLE SEAT WIDTH RANGE BETWEEN 11" AND 16"			
<input type="checkbox"/> SOWN11-16	Width adjustment from 11" to 16" with a maximum Depth of 16"		(Seat to Floor Height increases by 1/4")
<input type="checkbox"/> 11" (SOWN11)	<input type="checkbox"/> 13" (SOWN13)	<input type="checkbox"/> 15" (SOWN15)	
<input type="checkbox"/> 12" (SOWN12)	<input type="checkbox"/> 14" (SOWN14)	<input type="checkbox"/> 16" (SOWN16)	

▼ PEDIATRIC ADJUSTABLE SEAT DEPTH RANGE BETWEEN 11" AND 16"			
<input type="checkbox"/> 11" (SODN11)	<input type="checkbox"/> 13" (SODN13)	<input type="checkbox"/> 15" (SODN15)	
<input type="checkbox"/> 12" (SODN12)	<input type="checkbox"/> 14" (SODN14)	<input type="checkbox"/> 16" (SODN16)	

▼ STANDARD ADJUSTABLE SEAT WIDTH RANGE BETWEEN 15" AND 20"			
<input type="checkbox"/> SOWN15-20-19	(Width adjustment from 15" to 20" with a maximum Depth of 19")		N/C
<input type="checkbox"/> 15" (SOWN15)	<input type="checkbox"/> 17" (SOWN17)	<input type="checkbox"/> 19" (SOWN19S)	
<input type="checkbox"/> 16" (SOWN16)	<input type="checkbox"/> 18" (SOWN18)	<input type="checkbox"/> 20" (SOWN20S)	

▼ ADJUSTABLE SEAT DEPTH RANGE BETWEEN 15" AND 19"			
<input type="checkbox"/> 15" (SODN15)	<input type="checkbox"/> 17" (SODN17)	<input type="checkbox"/> 19" (SODN19)	
<input type="checkbox"/> 16" (SODN16)	<input type="checkbox"/> 18" (SODN18)		

▼ BACK CANE HEIGHT (The back cane height measurement is taken from the seat pan to the top of the push handles)			
<input type="checkbox"/> 16" (SOBCH16)	<input type="checkbox"/> 20" (SOBCH20)	<input type="checkbox"/> 24" (SOBCH24)	
<input type="checkbox"/> 18" (SOBCH18)	<input type="checkbox"/> 22" (SOBCH22) ³⁴	<input type="checkbox"/> 26" (SOBCH26)	

Note 34 22" is the minimum back cane height required with a power or manual recline.

▼ PRE-SET SEAT ANGLE (For Non-Reclining Backs)			
<input type="checkbox"/> SOBAA	Please Select one:	<input type="checkbox"/> 86° <input type="checkbox"/> 90° <input type="checkbox"/> 94° <input type="checkbox"/> 98° <input type="checkbox"/> 102° <input type="checkbox"/> 106° <input type="checkbox"/> 110° <input type="checkbox"/> 114° <input type="checkbox"/> 118° <input type="checkbox"/> 122°	N/C

▼ BACK SUPPORT			
<input type="checkbox"/> STDSB	³⁵	Standard contour back support 11" to 20" Wide with a maximum shell height of 22"	\$624.00
<input type="checkbox"/> STDSB/CUSTOM		Standard contour back support Non Standard Sizes up to 24" wide Width: _____ Depth: _____	\$1,333.00
Note 35 It is recommended to use the Back Support option with the power or manual recline.			

▼ BACK SUPPORT TOTAL HEIGHT (MEASURED FROM THE SEAT PAN TO THE TOP OF THE BACK)			
<input type="checkbox"/> 16" (SOBH16)	<input type="checkbox"/> 24" (SOBH24)		
<input type="checkbox"/> 18" (SOBH18)	<input type="checkbox"/> 26" (SOBH26)		
<input type="checkbox"/> 20" (SOBH20)	<input type="checkbox"/> 28" (SOBH28)		
<input type="checkbox"/> 22" (SOBH22) ³⁶			
Note 36 22" is the minimum back height required with a power or manual recline.			
Note 37 The space between the bottom of the back and the seat pan is 4" with a non-reclining back and 6" with a reclining back			

▼ SEAT CUSHIONS			
► Contour Seat Cushion			
<input type="checkbox"/> SCPCS	Contour Seat Cushion (11" to 20" Wide, with maximum 19" seat depth)		\$624.00
<input type="checkbox"/> SCPCS/CUSTOM	Contour Seat Cushion (Non Standard Sizes) Width: _____ Depth: _____		\$1,359.00

▼ CANTILEVER - HEIGHT ADJUSTABLE ARMRESTS - OPENS WITH A 5° ANGLE (For Non-Reclining Backs - Choose 2)			
<input type="checkbox"/> SOAFBFLH	Flip-Back Cantilever - Full Length Left Side (Height Adjustable from 7" to 15" with 1" increments)		STD
<input type="checkbox"/> SOAFBFRH	Flip-Back Cantilever - Full Length Right Side (Height Adjustable from 7" to 15" with 1" increments)		STD
<input type="checkbox"/> SOAFBDLH	Flip-Back Cantilever - Desk Length Left Side (Height Adjustable from 7" to 15" with 1" increments)		N/C
<input type="checkbox"/> SOAFBDRH	Flip-Back Cantilever - Desk Length Right Side (Height Adjustable from 7" to 15" with 1" increments)		N/C

▼ DUAL POST - HEIGHT ADJUSTABLE ARMRESTS - INCLUDES TRANSFER GRIPS (For Non-Reclining Backs - Choose 2)			
<input type="checkbox"/> SOADPPLH	Pediatric Dual Post Height Adjustable - Left Side (Height Adjustable from 7 1/4" - 10 1/4")		\$202.00
<input type="checkbox"/> SOADPPRH	Pediatric Dual Post Height Adjustable - Right Side (Height Adjustable from 7 1/4" - 10 1/4")		\$202.00
<input type="checkbox"/> SOADPFLH	Dual Post Height Adjustable - Full Length Left Side (Height Adjustable from 10 1/4" - 15 1/4")		\$202.00
<input type="checkbox"/> SOADPFRH	Dual Post Height Adjustable - Full Length Right Side (Height Adjustable from 10 1/4" - 15 1/4")		\$202.00
<input type="checkbox"/> SOADPDLH	Dual Post Height Adjustable - Desk Length Left Side (Height Adjustable from 10 1/4" - 15 1/4")		\$202.00
<input type="checkbox"/> SOADPDRH	Dual Post Height Adjustable - Desk Length Right Side (Height Adjustable from 10 1/4" - 15 1/4")		\$202.00

▼ RECLINING FLIP-UP, HEIGHT ADJUSTABLE ARMRESTS (For Reclining Backs - Choose 2)			
<input type="checkbox"/> SOANIFLH	Reclining armrest - Full Length Left Side (Height Adjustable from 10 1/2" to 14" with 1/2" increments)		N/C
<input type="checkbox"/> SOANIFRH	Reclining armrest - Full Length Right Side (Height Adjustable from 10 1/2" to 14" with 1/2" increments)		N/C
<input type="checkbox"/> SOANIDLH	Reclining armrest - Desk Length Left Side (Height Adjustable from 10 1/2" to 14" with 1/2" increments)		N/C
<input type="checkbox"/> SOANIDRH	Reclining armrest - Desk Length Right Side (Height Adjustable from 10 1/2" to 14" with 1/2" increments)		N/C

▼ ARMREST HEIGHT ADJUSTMENT AT FACTORY

Please set the Armrest Height at: _____ " (The Armrest Height is taken from the seat pan to the top of the armpad. Not required for Dual Post)

▼ ARMPADS

► Foam Type Selection

<input type="checkbox"/>	SOFAPD	Desk Length Standard (10")	Select Side: <input type="checkbox"/> Left Hand Side <input type="checkbox"/> Right Hand Side	STD
<input type="checkbox"/>	SOFAPF	Full Length Standard (14")	Select Side: <input type="checkbox"/> Left Hand Side <input type="checkbox"/> Right Hand Side	STD
<input type="checkbox"/>	SOWAPD	39 Desk Length Waterfall (10")	Select Side: <input type="checkbox"/> Left Hand Side <input type="checkbox"/> Right Hand Side	\$34.00
<input type="checkbox"/>	SOWAPF	39 Full Length Waterfall (14")	Select Side: <input type="checkbox"/> Left Hand Side <input type="checkbox"/> Right Hand Side	\$34.00

► Hybrid Gel and Foam Type Selection (The Gel always stays soft)

<input type="checkbox"/>	SOGELD	Desk Length Hybrid Gel and Foam Flat armpad (10" x 2.5")	Select Side: <input type="checkbox"/> Left Hand Side <input type="checkbox"/> Right Hand Side	\$34.00
<input type="checkbox"/>	SOGELF	Full Length Hybrid Gel and Foam Flat armpad (14" x 2.5")	Select Side: <input type="checkbox"/> Left Hand Side <input type="checkbox"/> Right Hand Side	\$34.00
<input type="checkbox"/>	SOWGELD	39 Desk Length Hybrid Gel and Foam Waterfall armpad (10")	Select Side: <input type="checkbox"/> Left Hand Side <input type="checkbox"/> Right Hand Side	\$54.00
<input type="checkbox"/>	SOWGELF	39 Full Length Hybrid Gel and Foam Waterfall armpad (14")	Select Side: <input type="checkbox"/> Left Hand Side <input type="checkbox"/> Right Hand Side	\$54.00

► Specialty Type Selection

<input type="checkbox"/>	SOCAR	40 Arm Trough with Hand Support (18" x 4 1/2")	Select Side: <input type="checkbox"/> Left Hand Side <input type="checkbox"/> Right Hand Side	\$222.00
<p>Note 39 The Waterfall armpads may interfere with flipping back the armrests when combined with Lateral Trunk supports</p> <p>Note 40 Arm Trough Armpads may interfere with Armrests Flipping back.</p>				

▼ FOOTRESTS

<input type="checkbox"/>	SOPAFAF	46 Pediatric Swing-Away Footrests	Set at: <input type="checkbox"/> 60° <input type="checkbox"/> 70° <input type="checkbox"/> 80° <input type="checkbox"/> 90°	\$334.00
			Select Footrest Range: <input type="checkbox"/> SOFLP5 (5 1/2" - 7 1/2") <input type="checkbox"/> SOFLP8 (8 1/2" - 10 1/2") <input type="checkbox"/> SOFLP10 (10 1/2" - 15 1/2")	
<input type="checkbox"/>	SOSAFAF	45, 46 Adult Swing-Away Footrests	Set at: <input type="checkbox"/> 60° <input type="checkbox"/> 70°	\$279.00
			Select Footrest Range: <input type="checkbox"/> SOFLA10 (10 1/2" - 14 1/2") <input type="checkbox"/> SOFLA14 (14 1/2" - 19 1/2")	
<input type="checkbox"/>	SOAAAF	45, 46 Adult Heavy Duty Swing-Away footrests	Set at: <input type="checkbox"/> 60° <input type="checkbox"/> 70°	\$390.00
			Select Footrest Range: <input type="checkbox"/> SOFLHD12 (12 1/2" - 16 1/2") <input type="checkbox"/> SOFLHD14 (14 1/2" - 19 1/2")	
<input type="checkbox"/>	ADDPL	42, 45, 46 Adult Power Elevating & Articulating Legrests	Please Select Function: <input type="checkbox"/> Individual <input type="checkbox"/> Combined	\$3,177.00
			Select Footrest Range: <input type="checkbox"/> SOFLPW12 (12 1/2" - 16 1/2") <input type="checkbox"/> SOFLPW14 (14 1/2" - 19 1/2")	

46 Please set the Footrest Length at: _____ "

- Note 41 Depending on the seat configuration, this option may not be available
- Note 42 This option includes calf pads
- Note 43 Individual legs count as two power functions.
- Note 44 Combined legs count as one power function.
- Note 45 This option is not available with the pediatric seat frame (SOWN11-16)
- Note 46 The Footrest Length is taken from the top of the seat pan to the top of the footplate. This footrest angle is not compatible with the 9" front casters (FC9PT and FC9PI)
- Note 47 This Footrest angle may interfere with the caster when selecting the straight front caster arm option.

▼ FOOTPLATES & ACCESSORIES

<input type="checkbox"/>	SOSFAF	Flip-up Footplates	STD
<input type="checkbox"/>	SOAFP	Angle Adjustable & Flip-up Footplates - Each	\$145.00
<input type="checkbox"/>	SOAHL	Heel Loops - Each	\$44.00
<input type="checkbox"/>	SOACPFR	Calf Pads for Adult Heavy-Duty Swing-Away footrests (SOAAF) with adjustable hardware - Each	\$145.00

▼ CENTER MOUNTS

<input type="checkbox"/>	SOAAF90	50, 51, 52, 55, 56 90° Fixed Center Mount Interface Bracket (Offers Manual adjustability from 95° to 70°)	\$256.00
			Select Foot Platform Range: <input type="checkbox"/> SOFBA12 (12 1/2" - 16 1/2") <input type="checkbox"/> SOFBA16 (16 1/2" - 19 1/2")
<input type="checkbox"/>	ADDPF	48, 49, 50, 52, 53, 54, 55 90° Power Center Mount Elevating & Articulating Foot Platform	\$3,177.00
			Select Foot Platform Range: <input type="checkbox"/> SOFBPW12 (12 1/2" - 16 1/2") <input type="checkbox"/> SOFBPW16 (16 1/2" - 20 1/2")

55 Please set the Foot Platform Length at: _____ "

- Note 48 The power center mount counts as one power function
- Note 49 This option includes calf pads
- Note 50 Links and attaches the one piece center mount foot platform (SOFUFP) or the split footplates (SOAFPFB) to the seat
- Note 51 With the split footplates option (SOAFPFB), this bracket offers independent height adjustability for each footplate
- Note 52 Please select the one piece center mount foot platform (SOFUFP) or the split footplates (SOAFPFB)
- Note 53 This option is not available with the pediatric seat frame (SOWN11-16)
- Note 54 This power center mount option is not available with the straight front caster arms (RWD) and Group 24/24 Frame (FS2424B-MWD)
- Note 55 The foot platform length is taken from the top of the seat pan to the top of the foot platform or split footplates.
- Note 56 This option is not compatible with the 9" Caster option on Rear Hybrid Configuration

▼ FOOT PLATFORMS & ACCESSORIES

<input type="checkbox"/>	SOFUFP	Flip-up one piece Foot Platform	N/C
<input type="checkbox"/>	SOAFPFB	Angle Adjustable & Flip-up Split Footplates for center mount - Each	\$145.00
<input type="checkbox"/>	SOACPCM	57 Calf Pads for 90° Fixed Center Mount (SOAAF90) with adjustable hardware - Each	\$145.00
<p>Note 57 Not available with the pediatric length (SOFBP7 - 7 1/2" - 12 1/2")</p>			

▼ SEAT POSITIONING (NOT COMPATIBLE WITH THE WC19 ----- ISO OPTION)

<input type="checkbox"/>	SOAUTO	55" Seat Positioning Strap	STD
--------------------------	--------	----------------------------	-----

▼ HEADREST

<input type="checkbox"/>	SONUHR	Using headrest mounting hardware other than SOMAH located below. Specify: _____	N/C
<input type="checkbox"/>	SOFHRS	Small Curved Headrest Pad 4" x 9"	\$228.00
<input type="checkbox"/>	SOFHRM	Medium Curved Headrest Pad 5" x 10"	\$228.00
<input type="checkbox"/>	SOFHRL	Large Curved Headrest Pad 6" x 11"	\$228.00
<input type="checkbox"/>	SOMAH	Multi-axis Removable Mounting Hardware	\$183.00
<input type="checkbox"/>	SOHRBC	Bracket to Mount on the Contour Back Support (Without Recline)	\$111.00
<input type="checkbox"/>	SOHRPR	Bracket to Mount on the Seat Frame (For Reclining Backs)	\$111.00

Please select one of the headrest options with the mounting hardware and specify where the headrest will be installed

▼ LATERAL TRUNK SUPPORTS

<input type="checkbox"/>	SOSLS3-C	Adjustable Laterals with Small Curved Pads (Each) - 3" x 5"	\$675.00
<input type="checkbox"/>	SOSLM4-C	Adjustable Laterals with Medium Curved Pads (Each) - 4" x 6"	\$675.00
<input type="checkbox"/>	SOSLL5-C	Adjustable Laterals with Large Curved Pads (Each) - 5" x 7"	\$675.00

▼ HIP & THIGH SUPPORTS

<input type="checkbox"/>	SOHTSS10	Adjustable Hip & Thigh Support Pads (Each) - Small 10" x 4"	\$256.00
<input type="checkbox"/>	SOHTSM12	Adjustable Hip & Thigh Support Pads (Each) - Medium 12" x 5"	\$256.00
<input type="checkbox"/>	SOHTSL14	Adjustable Hip & Thigh Support Pads (Each) - Large 14" x 6"	\$256.00
<input type="checkbox"/>	SOHTSARH	Adjustable and Removable Hip & Thigh Support Hardware - Each	\$485.00

Please select one of the hip & thigh support pads and the adjustable & removable hardware.

▼ ACCESSORIES

<input type="checkbox"/>	SOCH	Cup Holder (Self-Leveling)	\$101.00
<input type="checkbox"/>	SOCPH	Ram Mounts X-Grip Cell Phone Holder	\$173.00
<input type="checkbox"/>	SOTHS	Ram Mounts X-Grip Small Tablet Holder	\$196.00
<input type="checkbox"/>	SOTHL	Ram Mounts X-Grip Large Tablet Holder with Swivel	\$295.00

▼ REFLECTOR KIT

<input type="checkbox"/>	REFLECT	Red & White Reflector Kit Stickers (They will not be installed - Included in a bag with the owner's manual)	N/C
--------------------------	---------	---	-----

▼ SWITCHES AND OPTIONS

<input type="checkbox"/>	ASCPCP	PGDT R-NET Dongle and PC programmer kit	\$652.00
<input type="checkbox"/>	ASCPVR2	PGDT VR2 Dongle and PC programmer kit	\$652.00
<input type="checkbox"/>	SSBS	Buddy Button Switch	\$122.00
<input type="checkbox"/>	SSES	Egg Switch-Black	\$122.00
<input type="checkbox"/>	SSML	Ultra Light Switch	\$168.00

▼ SPECIALTY DRIVER CONTROL OPTIONS

▼ Therafin Sip-N-Puff Kit

<input type="checkbox"/>	SCSNPK	58, 59	Therafin Sip-N-Puff Kit		\$1,701.00
<input type="checkbox"/>	SCSNPBTk	59	Therafin Sip-N-Puff Breath Tube Kit		\$552.00
<input type="checkbox"/>	SCSRMSNP	59	Stealth Rock-n-Lock Swingaway Arm for Sip-N-Puff	<input type="checkbox"/> SWL (Swings on L) <input type="checkbox"/> SWR (Swings on Right)	\$307.00
<input type="checkbox"/>	SCSRM	59	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/> SWL (Swings on L) <input type="checkbox"/> SWR (Swings on Right)	\$552.00
<input type="checkbox"/>	SSES-STD	59, 61	Egg Switch for Changing Modes - Black	<input type="checkbox"/> MLH (Mount on L) <input type="checkbox"/> MRH (Mount on Right)	\$122.00
<input type="checkbox"/>	SSBS-STD	59	Buddy Button Switch (Stop Switch) - Red		\$122.00
<input type="checkbox"/>	SOMAH	59	Multi-axis Headrest Removable Mounting Hardware		\$183.00
<input type="checkbox"/>	SOFHRS	60	Small Curved Headrest Pad 4" x 9"		\$228.00
<input type="checkbox"/>	SOFHRM	60	Medium Curved Headrest Pad 5" x 10"		\$228.00
<input type="checkbox"/>	SOFHRL	60	Large Curved Headrest Pad 6" x 11"		\$228.00

▼ Switch-It Proximity Head Array

<input type="checkbox"/>	SCSIHAA	58, 59	Adult - 3 Switch Proximity Head Array with 1/8" Mono Jack for mode/reverse		\$5,679.00
<input type="checkbox"/>	SCSRM	59	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/> SWL (Swings on L) <input type="checkbox"/> SWR (Swings on Right)	\$552.00
<input type="checkbox"/>	SSES-STD	59, 61	Egg Switch for Changing Modes - Black	<input type="checkbox"/> MLH (Mount on L) <input type="checkbox"/> MRH (Mount on Right)	\$122.00
<input type="checkbox"/>	SSBS-STD	59	Buddy Button Switch (Stop Switch) - Red		\$122.00
<input type="checkbox"/>	SOMAH	59	Multi-axis Headrest Removable Mounting Hardware		\$183.00

▼ Proportional RIM Head Control Package

<input type="checkbox"/>	SCPHC	58, 59	Proportional RIM Head Control Package		\$2,781.00
<input type="checkbox"/>	SCSRM	59	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/> SWL (Swings on L) <input type="checkbox"/> SWR (Swings on Right)	\$552.00
<input type="checkbox"/>	SSES-STD	59, 61	Egg Switch for Changing Modes - Black	<input type="checkbox"/> MLH (Mount on L) <input type="checkbox"/> MRH (Mount on Right)	\$122.00
<input type="checkbox"/>	SOMAH	59	Multi-axis Headrest Removable Mounting Hardware		\$183.00

▼ Switch-It MicroPilot Light Touch Proportional Mini joystick for chin / lip control

<input type="checkbox"/>	SCSIMP	58, 59	MicroPilot Light touch mini Proportional Joystick		\$4,453.00
<input type="checkbox"/>	SCSRMMP	59	Stealth Rock-n-Lock Swingaway Arm for Chin MP	<input type="checkbox"/> SWL (Swings on L) <input type="checkbox"/> SWR (Swings on Right)	\$552.00
<input type="checkbox"/>	SCSRM	59	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/> SWL (Swings on L) <input type="checkbox"/> SWR (Swings on Right)	\$552.00
<input type="checkbox"/>	SSES-STD	59, 61	Egg Switch for Changing Modes - Black	<input type="checkbox"/> MLH (Mount on L) <input type="checkbox"/> MRH (Mount on Right)	\$122.00
<input type="checkbox"/>	JCC		Chin Cup		\$57.00
<input type="checkbox"/>	SOMAH	59	Multi-axis Headrest Removable Mounting Hardware		\$183.00
<input type="checkbox"/>	SOFHRS	60	Small Curved Headrest Pad 4" x 9"		\$228.00
<input type="checkbox"/>	SOFHRM	60	Medium Curved Headrest Pad 5" x 10"		\$228.00
<input type="checkbox"/>	SOFHRL	60	Large Curved Headrest Pad 6" x 11"		\$228.00

▼ Proportional Compact or Mini joystick for Midline Mount

<input type="checkbox"/>	SCPGMJS	58, 59	PG Drives JC200 Proportional Compact Remote Joystick for Midline Mount		\$1,051.00
<input type="checkbox"/>	SCGMJSM	62	Gatlin Midline Swing-away Mount for JC200	<input type="checkbox"/> MLH (Mount on L) <input type="checkbox"/> MRH (Mount on Right)	\$435.00
<input type="checkbox"/>	SCSIMP	58, 59	Switch-It MicroPilot Light touch mini Proportional Joystick		\$4,453.00
<input type="checkbox"/>	SCGMPPM	63	Gatlin Midline Swing-away Mount for MicroPilot	<input type="checkbox"/> MLH (Mount on L) <input type="checkbox"/> MRH (Mount on Right)	\$474.00
<input type="checkbox"/>	SCBTM	64	Bullet Hand Tray for MicroPilot	<input type="checkbox"/> MLH (Mount on L) <input type="checkbox"/> MRH (Mount on Right)	\$502.00

▼ Proportional Compact or Mini joystick for Armrest Mount

<input type="checkbox"/>	SCPGMJS	58, 59	PG Drives JC200 Proportional Compact Remote Joystick for Armrest Mount		\$1,059.00
<input type="checkbox"/>	SCARMJSM	62	Swing-away and Multi-Axis armrest mount	<input type="checkbox"/> MLH (Mount on L) <input type="checkbox"/> MRH (Mount on Right)	\$435.00
<input type="checkbox"/>	SCSIMP	58, 59	Switch-It MicroPilot Light touch mini Proportional Joystick		\$4,453.00
<input type="checkbox"/>	SCARMPM	63	Swing-away and Multi-Axis armrest mount	<input type="checkbox"/> MLH (Mount on L) <input type="checkbox"/> MRH (Mount on Right)	\$474.00
<input type="checkbox"/>	SCBTM	64	Bullet Hand Tray for MicroPilot	<input type="checkbox"/> MLH (Mount on L) <input type="checkbox"/> MRH (Mount on Right)	\$502.00

Note 58 Requires the OMNI Enhanced Display and Digital Interface Box (EDIB)

Note 59 This option must be selected for the kit

Note 60 Choose one (1) Headrest Pad

Note 61 The mode switch type can be changed from this selection on Page #7

Note 62 Must be selected for use with the PG Drives JC200 Mini Joystick

Note 63 Must be selected for use with the Switch-It MicroPilot

Note 64 This option is only compatible with the Switch-It MicroPilot

▼ CUSTOM SECTION

Note 65 Lead times will increase depending on the custom order pipeline at the time the purchase order is received. Please call for Quote and specifications for these options.

▼ O² Holder						
<input type="checkbox"/>	SUPOXY	66	Custom Non-Articulating Oxygen Tank Holder	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	\$435.00
<input type="checkbox"/>	SUPOXYA	67	Custom Articulating Oxygen Tank Holder	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	\$2,040.00
			<input type="checkbox"/> D-Size Oxygen Tank (* 22" Tall x 4.3" Diameter) * The height includes the valve			
			<input type="checkbox"/> E-Size Oxygen Tank (* 31" Tall x 4.3" Diameter) * The height includes the valve			
		Note 66	Not available with a power tilt and/or power recline			
		Note 67	Required with a power tilt and/or power recline			

▼ VENT TRAY (Please complete in Full)						
<input type="checkbox"/>	SOVENT		Custom Vent Tray			
<input type="checkbox"/>	SOVENT-A	68	Articulating Vent Box (Required with Power Tilt and/or Recline)	Model:		\$1,957.00
<input type="checkbox"/>	SOVENT-F	68	Fixed Vent Box	Model:		\$491.00
<input type="checkbox"/>	SOVTBE	68, 69, 70	Extended Battery Compartment			\$600.00
			<input type="checkbox"/> The Vent is being powered by 1 U1 Battery and the chair by 2 Group 22NF batteries			
			<input type="checkbox"/> The Vent is being powered by 1 U1 Battery and the chair by 2 Group 24 batteries			
			<input type="checkbox"/> The Vent is being powered by 1 Group 22NF Battery and the chair Powered by 2 Group 22NF Batteries			
			<input type="checkbox"/> The Vent is being powered by 1 Group 24 Battery and the chair Powered by 2 Group 22NF Batteries			
		Note 68	Amysystems does not supply the battery, charger or the wiring for the Vent. These items are normally provided by the Vent Manufacturer.			
		Note 69	This option is not available with the straight front caster arms.			
		Note 70	The Vent Tray option is not available with the CUSTOM Battery Tray for Group 24 Batteries with the 24" frame width (FS2424B)			

▼ OTHER CUSTOM OPTIONS (Request For Quote)						
Please specify details and specifications for these options.						
<input type="checkbox"/>	PAT	71	Power Anterior Tilt			RFQ
<input type="checkbox"/>	CUSTOM	71	Power Lateral Tilt 15° <input type="checkbox"/> Raise the Left Side only <input type="checkbox"/> Raise the Right Side only <input type="checkbox"/> Raise the Left & Right Side			RFQ
<input type="checkbox"/>	CUSTOM	71, 73	Power Elevating Stretcher (Please call CSR for availability with configuration and user weight capacity)			RFQ
<input type="checkbox"/>	CUSTOM	71, 73	Power Flip-up Foot Platform (Please call CSR for availability with configuration and user weight capacity)			RFQ
<input type="checkbox"/>	CUSTOM	71, 73	Power Vertical Lifting Foot Platform (Please call CSR for availability with configuration and user weight capacity)			RFQ
<input type="checkbox"/>	CUSTOM	71, 72, 74	Power Stander (Requires Power Tilt, Power Recline and 90° Power elevating & articulating center mount foot platform)			RFQ
<input type="checkbox"/>	CUSTOM	71	Independent Depth Adjustability built into the Footrest Receiver that is attached to the Seat Pan			RFQ
<input type="checkbox"/>	CUSTOM		Transfer Handles (Pair)			RFQ
<input type="checkbox"/>	CUSTOM		Manual Wheel Locks			RFQ
<input type="checkbox"/>	CUSTOM		Non Standard Frame Colors (Not all custom colors are available, please contact CSR for availability)			RFQ
<input type="checkbox"/>	CUSTOM # 1	71	Other custom options. Specify in notes section below.			RFQ
<input type="checkbox"/>	CUSTOM # 2	71	Other custom options. Specify in notes section below.			RFQ
		Note 71	Please specify details and specifications for these options.			
		Note 72	Please call CSR to verify if the power standing option is available with the required configuration for your user.			
		Note 73	This option is not available with the Straight Front Caster Arms on Rear Hybrid Configuration.			
			option is included at no charge with the H ₃ The power Stander is only available with the mid-wheel drive base configuration.			

▼ NOTES / SPECIAL REQUESTS / CUSTOM SPECIFICATIONS
