



July 01, 2015

ACCOUNT #: _____ SHIP TO: _____
 DATE: _____ NAME: _____
 P.O. #: _____ ADDRESS: _____
 CITY: _____ PROVINCE: _____
 POSTAL CODE: _____ PHONE #: _____
 BUYER: _____ MARKED FOR: _____



Product Info ▶

QUICKIE XPLORE SPECIFICATIONS

XPLORE WPA3QU270 Xplore™ Rear Wheel Drive \$8,486.00 XPLOREPS WPA3QU270 Xplore™ For Power Seating Options \$9,338.00 Kid: WPK3QU070

USER WEIGHT PACKAGES

- WCA3 WEPF 300 lb Standard User Weight Capacity STD
 WCA4 WEPG 400 lb Heavy - Duty Package (Re-inforced Frame) \$520.00

OPTIONAL MOTOR PACKAGE

- MP80-OPT 4 Pole 8.0 MPH High-Speed Motors (Max. user weight 250 lbs - Not available with VR2 electronics ECV902A) N/C

ASSISTIVE DRIVING PACKAGE

- INTELLIDR 1, 2 Intelli-Drive \$1,550.00
 Note 1 Intelli-Drive technology provides increased tracking capabilities with enhanced traction and steering precision. This is ideal for those with marginal hand function or switch control.
 Note 2 Requires the R-net PM 120 Amp Expandable Controller (ECR120)

TRANSPORT OPTION

- TTD4 3 4-Point wheelchair transport brackets 4-Point wheelchair transport brackets have not yet received WC-19 crash-testing approval. STD
 Note 3 Use these transport brackets only to secure an unoccupied wheelchair during transport.

BATTERY TRAY AND FRAME SIZE

- FS22B Battery Tray for Group 22NF Batteries (Total frame width = 24") STD
 FS24B Battery Tray for Group 24 Batteries (Total frame width = 25 3/8") STD
 Front Caster Arms Flared out (STD) Straight (Optional - Not available with the power Center Mount (ADDPF and ADPPFHD) N/C

BATTERIES

- BAT22NFPR MK Gel Batteries - Group 22NF *Sold as a pair Each \$301.00
 BAT24PR MK Gel Batteries - Group 24 *Sold as a pair Each \$403.00

BATTERY CHARGER

- BCOB Offboard Battery Charger N/C
 BCLES8OPT Lester 8 amp Offboard Battery Charger N/C

DRIVE WHEEL

- DW14PT 14" Cast Wheel Pneumatic Tire STD
 DW14PI 14" Cast Wheel Pneumatic with Flat-Free inserts \$235.00
 DW144PT 5, 7 14" x 4" Extra Wide Cast Wheel Pneumatic Tire **NEW!** \$255.00
 DW144PI 5, 7 14" x 4" Extra Wide Cast Wheel Pneumatic with Flat-Free inserts **NEW!** \$255.00
 DW12PT 6 12 1/2" Cast Wheel Pneumatic Tire STD
 DW12PI 6 12 1/2" Cast Wheel Pneumatic with Flat-Free inserts \$235.00

Note 5 The 14" x 4" Extra Wide Wheel increases the overall width of the Base by 2 1/2"
 The option to Match the Wheels to the Shroud Color (FC...-M) is not available with the 14" x 4" Extra Wide wheels
 Note 6 12 1/2" Drive wheels reduce overall seat to floor height by 3/4", lowers ground clearance and speed, but increases motor torque
 Note 7 The 14" x 4" option is not available with VR2 electronics (ECV902A)

FRONT CASTERS

- FC8PTR 6" Solid Casters STD
 FC8PIR 8 8" Pneumatic Casters \$235.00
 FC9PT 8 8" Pneumatic Casters with Flat-Free inserts \$158.00
 FC9PI 8 8" Pneumatic Casters with Flat-Free inserts \$286.00

Note 8 9" front casters increases seat to floor height by 1/2"
 The 9" casters are NOT compatible with the 12 1/2" Drive wheel configuration (DW12PT and DW12PI)

SEAT TO FLOOR HEIGHT (Approximate measurement to the top of seat pan ± 1/4")

Options for a Seat Only

- 9 STFH1525 15 1/4" (Possible Seat angle of 0° only)
 STFH16 16" (Possible Seat angle of 0° only)
 STFH17 17" (Possible Seat angle of 0° and 3°)
 STFH18 18" (Possible Seat angle of 0° and 3°)
 STFH19 19" (Possible Seat angle of 0°, 3° and 6°)
 STFH20 20" (Possible Seat angle of 0°, 3°, 6°, 9°)

Note 9 Requires the 12 1/2" Drive wheels (DW12PT or DW12PI). This lowers ground clearance and speed, but increases motor torque.
 Note 10 The Heavy-Duty 11" Power Elevating Seat increases seat to floor height by approximately 1 1/2" (ADDPESXHD)

Options for a Power Tilt or a STD Power Elevating Seat

- 9 STFH16 16" (Possible Seat angle of 0° only)
 10 STFH1675 16 3/4" (Possible Seat angle of 0° only)
 10 STFH1775 17 3/4" (Possible Seat angle of 0° only)
 STFH1825 18 1/4" (Possible Seat angle of 0° only)
 STFH1875 18 3/4" (Possible Seat angle of 0° only)
 STFH1925 19 1/4" (Possible Seat angle of 0° only)

Options for a Combo - Power Tilt with a STD Power Elevating Seat

- 9, 10 STFH1825 18 1/4" (Possible Seat angle of 0° only)
 10 STFH19 19" (Possible Seat angle of 0° only)

PRE-SET SEAT FRAME ANGLE

<input type="checkbox"/>	SFE-0	0° Pre-set seat frame seat angle (Must be selected when using power options)	STD
<input type="checkbox"/>	SFE-3	3° Pre-set seat frame seat angle (Only available with a seat configuration without power options)	STD
<input type="checkbox"/>	SFE-6	6° Pre-set seat frame seat angle (Only available with a seat configuration without power options)	STD
<input type="checkbox"/>	SFE-9	9° Pre-set seat frame seat angle (Only available with a seat configuration without power options)	STD

LIGHTING PACKAGE

<input type="checkbox"/>	EOLP	11 LED Light Package (Includes front beams, side flashers, rear tail light with flashers & brake light)	\$1,035.00
--------------------------	------	---	------------

Note 11 Requires the upgraded Joystick with Light Operation or the OMNI (ECVJ6KL, ECELJL, ECRJCL, ECRJCTPL or EDIB)

PG DRIVES 90 AMP STANDARD CONTROLLER & JOYSTICK

<input type="checkbox"/>	ECV902A	VR2 90 Amp Drive Controller	STD
<input type="checkbox"/>	ECVJ4K	12 VR2 Standard 4-Key Joystick	STD
<input type="checkbox"/>	ECVJ6K	12 VR2 6-Key Joystick (Allows for a single power option to operate through this joystick)	STD
<input type="checkbox"/>	ECVJ9KL	12 VR2 9-Key Joystick with Lighting Operation (Allows for a single power function to operate through this joystick)	\$153.00

Note 12 Requires the VR2 90A Drive Control (ECV902A)

PG DRIVES 90 AMP FULLY EXPANDABLE CONTROLLER & JOYSTICK FOR POWER SEATING & SPECIALTY CONTROLS

<input type="checkbox"/>	ECEL90	R-net EL 90 Amp Expandable Controller	\$515.00
<input type="checkbox"/>	ECRNH	13, 14 Harness for Expandable Controller ()	
<input type="checkbox"/>	ECELJ	13 R-net EL (LED) Multiple Drive Joystick	N/C
<input type="checkbox"/>	ECELJL	13 R-net EL (LED) Multiple Drive Joystick with Lighting Operation	\$148.00

Note 13 Requires the fully expandable EL 90 amp controller (ECEL90)
Note 14 Must be ordered with the expandable controller ECEL90 to interface with Joystick or OMNI screen

PG DRIVES 120 AMP FULLY EXPANDABLE CONTROLLER & JOYSTICK FOR POWER SEATING & SPECIALTY CONTROLS

<input type="checkbox"/>	ECR120	R-net PM 120 Amp Expandable Controller	\$1,035.00
<input type="checkbox"/>	ECRNH	15, 16 Harness for Expandable Controller	
<input type="checkbox"/>	ECRJJC	15 R-net Multiple Drive Standard Joystick with Color screen	\$153.00
<input type="checkbox"/>	ECRJCL	15 R-net Multiple Drive Joystick with Color screen and Lighting Operation	\$204.00
<input type="checkbox"/>	ECRJCTP	15 R-net Multiple Drive Joystick with Color screen w/ On & Off / Mode switch & speed pot	\$255.00
<input type="checkbox"/>	ECRJCTT	15 R-net Multiple Drive Joystick with Color screen w/ On & Off / Mode switch & profile	\$306.00
<input type="checkbox"/>	ECRJCTPL	15 R-net Multiple Drive Joystick with Color screen w/ On & Off / Mode switch & speed pot and Lighting Operation	\$250.00

Note 15 Requires the fully expandable R-net PM 120 amp controller (ECR120)
Note 16 Must be ordered with the expandable controller ECR120 to interface with Joystick or OMNI screen

JOYSTICK MOUNT

<input type="checkbox"/>	JSMF	Fixed Joystick Mount	N/C
<input type="checkbox"/>	JMSA	Swing Away Joystick Mount - Height adjustable, Angle adjustable and Multi-Axis	\$377.00
<input type="checkbox"/>	JSMCM	Center Mounted Swing Away Joystick Mount	\$377.00

JOYSTICK MOUNT SUPPORT

<input type="checkbox"/>	JSMFH	Fixed joystick mount support	STD
<input type="checkbox"/>	JSMAH	Adjustable height joystick mount support (Not required with the Swing Away Joystick Mount - JMSA)	\$133.00

JOYSTICK HAND MOUNT

<input type="checkbox"/>	HML	Install Joystick on Left hand side	N/C	<input type="checkbox"/>	HMR	Install Joystick on Right hand side	N/C
--------------------------	-----	------------------------------------	-----	--------------------------	-----	-------------------------------------	-----

JOYSTICK HANDLE

<input type="checkbox"/>	JSH	Stick Handle	STD
<input type="checkbox"/>	JBH	Ball Handle	\$92.00
<input type="checkbox"/>	JFBH	Foam Ball Handle	\$92.00
<input type="checkbox"/>	JTH	T-Handle	\$92.00
<input type="checkbox"/>	JGP	Goalpost	\$92.00

ATTENDANT DRIVE CONTROL OPTION

<input type="checkbox"/>	EADC	Attendant Hand Control for R-net EL & PM Controls	\$1,035.00
--------------------------	------	---	------------

ATTENDANT JOYSTICK HAND MOUNT

<input type="checkbox"/>	EADCL	Install Attendant Joystick on Rear Left hand side	N/C	<input type="checkbox"/>	EADCR	Install Attendant Joystick on Rear Right hand side	N/C
--------------------------	-------	---	-----	--------------------------	-------	--	-----

ENHANCED DISPLAY & DIGITAL INTERFACE BOX (Included with specialty controls)

<input type="checkbox"/>	EDIB	17 Omni (Enhanced Display & Digital Interface Box with Infrared)	\$1,974.00
--------------------------	------	--	------------

Note 17 Allows the switch operation of alternative devices through driver controls.
Requires the R-net PM (ECR120) Fully Expandable Controller

ENHANCED DISPLAY MOUNT

<input type="checkbox"/>	EDSAM	User Swing Away Mount	<input type="checkbox"/>	EDML	(Install Display on LH side)	<input type="checkbox"/>	EDMR	(Install Display on RH side)	N/C
--------------------------	-------	-----------------------	--------------------------	------	------------------------------	--------------------------	------	------------------------------	-----

INPUT / OUTPUT MODULE (ENVIRONMENTAL CONTROL MODULE)

<input type="checkbox"/>	IOM	18	IOM Input / Output Module (environment control module)	\$1,035.00
Note 18 Requires the R-net PM (ECR120) Fully Expandable Controller				

BLUETOOTH MOUSE EMULATOR (Includes USB Dongle)

<input type="checkbox"/>	EDBT	19	Bluetooth Mouse Emulator (USB Dongle included)	\$1,035.00
Note 19 Requires the fully expandable R-net PM 120 amp controller (ECR120)				

POWER OPTIONS & SEATING

<input type="checkbox"/>	SAABA	24	Eclipse Seat - Frame Only *Comes STD with XPLOREPS (Xplore - Power Seating Options)*	\$806.00
<input type="checkbox"/>	SAABA-COR3	21, 27	Rehab Seat with Corbeau Sport Back - Includes a manual Recline Release Lever <input type="checkbox"/> Lever LH <input type="checkbox"/> Lever RH NEW!	\$2,244.00
<input type="checkbox"/>	SAPTX SEC4P00	20, 22	CG Power Tilt	\$6,758.00
<input type="checkbox"/>	SAPTXHD	20, 24	Heavy-Duty CG Power Tilt NEW!	\$7,538.00
<input type="checkbox"/>	SAPRX SEC4P00	20,22,25	Power Recline with Integrated Shear Reduction (ISR)	\$8,318.00
<input type="checkbox"/>	SAPRXHD	20,24,25	Heavy-Duty Power Recline with Integrated Shear Reduction (ISR) NEW!	\$9,098.00
<input type="checkbox"/>	SAPTPRX SEC4	20,22,25	CG Power Tilt and Power Recline with Integrated Shear Reduction (ISR)	\$11,960.00
<input type="checkbox"/>	SAPTPRXHD	20,24,25	Heavy-Duty CG Power Tilt and Power Recline with Integrated Shear Reduction (ISR) NEW!	\$13,520.00
<input type="checkbox"/>	SAPESX	20, 22	11" Power Elevating Seat Module for transfer assist	\$3,351.00
<input type="checkbox"/>	SAPESXHD	20,23,26	Heavy-Duty 11" Power Elevating Seat Module for transfer assist NEW!	\$5,197.00
<input type="checkbox"/>	ADDPESX	22	Module - Add an 11" Power Elevating Seat to either SAPTX, SAPRX or SAPTPRX	\$3,116.00
<input type="checkbox"/>	ADDPESXHD	23, 26	Module - Add a Heavy-Duty 11" Power Elevating Seat to either SAPTXHD, SAPRXHD or SAPTPRXHD	\$4,157.00
<input type="checkbox"/>	ADDMR	22, 25	Module - Manual Recline with Integrated Shear Reduction (ISR)	\$1,035.00

Note 20 These options include the fully adjustable Amyseat (SAABA).

Note 21 The maximum user weight capacity with this option is 250 lbs

Note 22 The maximum user weight capacity with this option is 300 lbs

Note 23 The maximum user weight capacity with this option is 400 lbs

Note 24 The maximum user weight capacity with this option is 450 lbs

Note 25 It is recommended to use the Contour Back Support option with the power or manual reclines (STDSB or STDSB/CUSTOM)

Note 26 The Heavy-Duty 11" Power Elevating Seat increases seat to floor height by approximately 1 1/2" (ADDPESXHD)

The Corbeau Back is not available with the Power Recline options, but includes a manual Recline. However, the 50° Power Tilt,

Note 27 the 11" Power elevating seat module and the power elevating legs or center mount options are compatible.

The Corbeau Sport Back option only comes in one Size, but offers a seat width range from 16" to 20"

ELECTRONICS - SINGLE OR DUAL (1 OR 2 ACTUATORS) POWER SEAT FUNCTION CONTROLLERS

<input type="checkbox"/>	SAC	Single (1 Actuator) Power Seat Function Controller to Interface through R-net Electronics	N/C
<input type="checkbox"/>	DAC	Dual (2 Actuators) Power Seat Function Controller to Interface through R-net Electronics	\$515.00

With the SAC or DAC, please select one activation and position from the list below

<input type="checkbox"/>	SWTJSS	28, 29	Single actuator activation through R-net Joystick or Specialty Controls	\$1,035.00
<input type="checkbox"/>	SWTJSD	28, 30	Dual actuator activation through R-net Joystick or Specialty Controls	\$1,035.00
<input type="checkbox"/>	SWTOGGLE	29 or 30	Single Toggle Switch (Single Mode: DAC - Dual Mode: SAC) <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation	N/C
<input type="checkbox"/>	SWPUSH1	29	Single Push Button Switch (Single Mode: SAC) <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation	N/C
<input type="checkbox"/>	SWPUSH2	29 or 30	Dual Push Button Switch (Single Mode: DAC - Dual Mode: SAC) <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation	N/C
<input type="checkbox"/>	EPSV	31	Tilt or Elevating seat activation through 6-key VR2 Hand Control	N/C

Note 28 Requires the R-net EL (ECEL90) or the PM (ECR120) Fully Expandable Controllers

Note 29 Requires the Single (1 Actuator) Power Seat Function Controller (SAC)

Note 30 Requires the Dual (2 Actuators) Power Seat Function Controller (DAC)

Note 31 Allows for the activation of a power tilt or the 11" power elevating seat through the 6-key VR2 joystick (ECVJ6K)

ELECTRONICS - MULTIPLE (3 to 6 ACTUATORS) POWER SEAT FUNCTION CONTROLLER

<input type="checkbox"/>	MAC SEC4P50	Multiple (3 to 6 Actuators) Power Seat Function Controller to Interface through R-net Electronics	\$2,489.00
--------------------------	-------------	---	------------

With the MAC, please select one activation and position from the list below

<input type="checkbox"/>	SWTJSM	32, 33	Multiple Actuator activation through R-net Joystick or Specialty Controls	N/C
<input type="checkbox"/>	SWMULTI	32	4-Way Toggle Switch (Up to 4 Activations) - <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation	N/C
<input type="checkbox"/>	SWMULTIT	32	Triple Toggle Switch (Up to 5 Activations) - <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation	N/C
<input type="checkbox"/>	SWPUSH4	32	4-Way Push Button Switch - <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation	N/C
<input type="checkbox"/>	SWPUSH5	32	5-Way Push Button Switch - <input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation	N/C

Note 32 Requires the Multiple (3 to 6 Actuators) Power Seat Function Controller (MAC)

Note 33 Requires the R-net EL (ECEL90) or the PM (ECR120) Fully Expandable Controller, and the multiple Actuator Controller (MAC)

SWITCH ACTIVATION HAND MOUNT

Switch Installation:	<input type="checkbox"/> TSML (Install Switch on LH side)	<input type="checkbox"/> TSMR (Install Switch on RH side)	N/C
----------------------	---	---	-----

PRE-SET SEAT ANGLE (For Non-Reclining Backs)

SOBAA Please Select one: 86° 90° 94° 98° 102° 106° 110° 114° 118° 122° N/C

STANDARD ADJUSTABLE SEAT WIDTH RANGE BETWEEN 15" AND 20" WEPA STD

- SOWN15-20-19 (Width adjustment from 15" to 20" with a maximum Depth of 19")
- SOWN15-20-22 (Width adjustment from 15" to 20" with a maximum Depth of 22")³⁴
- 15" (SOWN15) 18" (SOWN18)³⁴
- 16" (SOWN16)³⁴ 19" (SOWN19S)³⁴
- 17" (SOWN17)³⁴ 20" (SOWN20S)³⁴

Note 34 For the Rehab Seat with Corbeau Sport Back (SAABA-COR3), please choose the long seat option SOWN15-20-22
The only possible Seat width options for the Corbeau Sport Seat (SAABA-COR3) are from 16" to 20" wide

LARGE ADJUSTABLE SEAT WIDTH RANGE BETWEEN 19" AND 24" WEPA \$449.00

- SOWN19-24-22 (Width adjustment from 19" to 24" and Depth range of 18" to 22")
- 19" (SOWN19L) 22" (SOWN22)
- 20" (SOWN20L) 23" (SOWN23)
- 21" (SOWN21) 24" (SOWN24)

ADJUSTABLE SEAT DEPTH RANGE BETWEEN 15" AND 22" WEPB

- 15" (SODN15) 16" (SODN16) 17" (SODN17) 18" (SODN18) STD
- 19" (SODN19) 20" (SODN20) 21" (SODN21) 22" (SODN22) \$306.00

BACK CANE HEIGHT³⁷ (This only considers the useful straight portion, and you must consider +3" to the top of the push handles for the total height) N/C

- 12" (SOBH12)³⁵ 15" (SOBH15) 18" (SOBH18) 21" (SOBH21)
- 13" (SOBH13) 16" (SOBH16) 19" (SOBH19) 22" (SOBH22)
- 14" (SOBH14) 17" (SOBH17)³⁶ 20" (SOBH20)

Note 35 This back cane height is not available with power options and some specialty controls.

Note 36 17" is the minimum back cane height required with a power or manual recline.

Note 37 The Back Cane option is not compatible for use with the Rehab Seat with Corbeau Sport Back (SAABA-COR3)

BACK SUPPORT

- STDSB [SEB4P0010](#)³⁸ Contour back support 15" to 24" Wide with a maximum shell height of 22" \$571.00
- STDSBO [SEB4P0](#) Omission of the Contour Back Support (CREDIT)
- STDSB/CUSTOM Contour Back Support for Non-Standard Sizes **Please Specify** Width: _____ Height: _____ \$1,244.00

Note 38 It is recommended to use the Contour Back Support option with the power or manual recline options. Not required for Corbeau

BACK SUPPORT HEIGHT (SHELL)

- 12" (SOBCH12) 16" (SOBCH16)³⁹ 20" (SOBCH20)
- 14" (SOBCH14) 18" (SOBCH18) 22" (SOBCH22)

Note 39 With a power or manual recline, the minimum back support shell height is 16".

The space between the bottom of the back and the seat pan is 3" with a fixed back and 6" with a recline

CANTILEVER - HEIGHT ADJUSTABLE ARMRESTS (For Non-Reclining Backs - Choose 2)

- SOACORFLH Corbeau Flip-Back Cantilever - Full Length Left Side (Height Adjustable from 9 ½" - 13 ½") N/C
- SOACORFRH Corbeau Flip-Back Cantilever - Full Length Right Side (Height Adjustable from 9 ½" - 13 ½") N/C
- SOACORDLH Corbeau Flip-Back Cantilever - Desk Length Left Side (Height Adjustable from 9 ½" - 13 ½") N/C
- SOACORDRH Corbeau Flip-Back Cantilever - Desk Length Right Side (Height Adjustable from 9 ½" - 13 ½") N/C
- SOAFBFLH Flip-Back Cantilever - Full Length Left Side (Height Adjustable from 11" to 15" in 1" increments) STD
- SOAFBFRH Flip-Back Cantilever - Full Length Right Side (Height Adjustable from 11" to 15" in 1" increments) STD
- SOAFBDLH Flip-Back Cantilever - Desk Length Left Side (Height Adjustable from 11" to 15" in 1" increments) N/C
- SOAFBDRH Flip-Back Cantilever - Desk Length Right Side (Height Adjustable from 11" to 15" in 1" increments) N/C

DUAL POST - HEIGHT ADJUSTABLE ARMRESTS (For Non-Reclining Backs - Choose 2)

- SOADPFLH Dual Post Height Adjustable - Full Length Left Side (Height Adjustable from 9 ½" - 13 ½") \$184.00
- SOADPFRH Dual Post Height Adjustable - Full Length Right Side (Height Adjustable from 9 ½" - 13 ½") \$184.00
- SOADPDLH Dual Post Height Adjustable - Desk Length Left Side (Height Adjustable from 9 ½" - 13 ½") \$184.00
- SOADPDRH Dual Post Height Adjustable - Desk Length Right Side (Height Adjustable from 9 ½" - 13 ½") \$184.00

RECLINING FLIP-UP, HEIGHT ADJUSTABLE ARMRESTS (For Reclining Backs - Choose 2)

- SOANIFLH Reclining armrest - Full Length Left Side (Height Adjustable from 10 ½" to 14" in ½" increments) N/C
- SOANIFRH Reclining armrest - Full Length Right Side (Height Adjustable from 10 ½" to 14" in ½" increments) N/C
- SOANIDLH Reclining armrest - Desk Length Left Side (Height Adjustable from 10 ½" to 14" in ½" increments) N/C
- SOANIDRH Reclining armrest - Desk Length Right Side (Height Adjustable from 10 ½" to 14" in ½" increments) N/C

ARMREST HEIGHT ADJUSTMENT AT FACTORY

Please set the Armrest Height at: _____ " (The Armrest Height is taken from the seat pan to the top of the armpad. Not required for Dual Post)

ARMPADS

<input type="checkbox"/>	SOWAPD	Desk Length Waterfall (10") - Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$31.00
<input type="checkbox"/>	SOWAPF	Full Length Waterfall (14") - Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$31.00
<input type="checkbox"/>	SOFAPD	Desk Length Standard (10") - Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	STD
<input type="checkbox"/>	SOFAPF	Full Length Standard (14") - Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	STD
<input type="checkbox"/>	SOCAR	40 Arm Trough with Hand Support (18" x 4 1/2") - Select Side:	<input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$204.00

Note 40 Arm Trough Armpads may interfere with Armrests Flipping back.

FOOTRESTS

<input type="checkbox"/>	SOSAFAF	45, 46 Adult Swing-Away Footrests ----- Set at: <input type="checkbox"/> 60° <input type="checkbox"/> 70° <input type="checkbox"/> 80°			\$255.00
		Select Footrest Range: <input type="checkbox"/> SOFLA10 (10 1/2" - 15 1/2")	<input type="checkbox"/> SOFLA14 (14 1/2" - 19 1/2")		
<input type="checkbox"/>	SOAAF	45, 47 Adult Heavy Duty Swing-Away footrests ----- Set at: <input type="checkbox"/> 60° <input type="checkbox"/> 70° <input type="checkbox"/> 80°			\$357.00
		Select Footrest Range: <input type="checkbox"/> SOFLHD12 (12 1/2" - 15 1/2")	<input type="checkbox"/> SOFLHD14 (14 1/2" - 19 1/2")		
<input type="checkbox"/>	ADDPL	42, 45, 46 Adult Power Elevating & Articulating Legrests	Please Select Function: <input type="checkbox"/> Individual ⁴³	<input type="checkbox"/> Combined ⁴⁴	\$2,907.00
	SEF4P000	Select Footrest Range: <input type="checkbox"/> SOFLPW12 (12 1/2" - 15 1/2")	<input type="checkbox"/> SOFLPW14 (14 1/2" - 19 1/2")		
<input type="checkbox"/>	ADDPLHD	42, 45, 47 Heavy Duty Power Elevating & Articulating Legrests	Select Function: <input type="checkbox"/> Individual ⁴³	<input type="checkbox"/> Combined ⁴⁴	\$3,111.00
		Select Footrest Range: <input type="checkbox"/> SOFLPW12 (12 1/2" - 15 1/2")	<input type="checkbox"/> SOFLPW14 (14 1/2" - 19 1/2")		

⁴⁵ Please set the Footrest Length at: _____ "

Note 41 Depending on the seat configuration, this option may not be available.

Note 42 This option includes calf pads

Note 43 Individual legs count as two power functions.

Note 44 Combined legs count as one power function.

Note 45 The Footrest Length is taken from the top of the seat pan to the top of the footplate.

Note 46 The Standard user weight capacity is 300 Lbs.

Note 47 The Standard user weight capacity is compatible for 450 Lbs.

FOOTPLATES & ACCESSORIES

<input type="checkbox"/>	SOSFAF	Flip-up Footplates			STD
<input type="checkbox"/>	SOAFP	WEPJ Angle Adjustable & Flip-up Footplates			\$133.00
<input type="checkbox"/>	SOAHL	Heel Loops - Pair			\$41.00
<input type="checkbox"/>	SOACPFR	Calf Pads for Adult Heavy-Duty Swing-Away footrests (SOAAF) with adjustable hardware - Pair			\$133.00

CENTER MOUNTS

<input type="checkbox"/>	SOAAF90	50,51,52 90° Fixed Center Mount Interface Bracket			\$235.00
		Select Foot Platform Range: <input type="checkbox"/> SOFBA12 (12 1/2" - 16 1/2")	<input type="checkbox"/> SOFBA16 (16 1/2" - 20 1/2")		
<input type="checkbox"/>	ADDPF	48, 49, 50, 52, 53 90° Power Center Mount Elevating & Articulating Foot Platform - *** Not available with the 12 1/2" Drive Wheel			\$2,907.00
	SEF4P000	<input type="checkbox"/> SOFBPW12 (12 1/2" - 16 1/2")	<input type="checkbox"/> SOFBPW16 (16 1/2" - 20 1/2")		
<input type="checkbox"/>	ADDPFHD	48, 49, 50, 52, 54 90° HD Power Center Mount Elevating & Articulating Foot Platform - *** Not available with the 12 1/2" Drive Wheel			\$3,111.00
		<input type="checkbox"/> SOFBPW12 (12 1/2" - 16 1/2")	<input type="checkbox"/> SOFBPW16 (16 1/2" - 20 1/2")		

⁵⁰ Please set the Foot Platform Length at: _____ "

Note 48 The power center mount counts as one power function.

Note 49 This option includes calf pads

Note 50 The Foot Platform Length is taken from the top of the seat pan to the top of the foot platform or split foot plates.

Note 51 Links and attaches the one piece center mount foot platform (SOFUFP) or the split footplates (SOAFPFB) to the seat

Note 52 Please Select the one piece center mount foot platform (SOFUFP) or the split footplates (SOAFPFB)

Note 53 The Standard user weight capacity is 300 Lbs.

Note 54 The Standard user weight capacity is compatible for 450 Lbs.

FOOT PLATFORMS & ACCESSORIES

<input type="checkbox"/>	SOFUFP	WEPJ Flip-up one piece Foot Platform			N/C
<input type="checkbox"/>	SOAFPFB	WEPJ Angle Adjustable & Flip-up Split Footplates for center mount			\$133.00
<input type="checkbox"/>	SOACPCM	Calf Pads for 90° Fixed Center Mount (SOAAF90) with adjustable hardware - Pair			\$133.00



SEAT POSITIONING

<input type="checkbox"/>	SOAUTO	55" Seat Positioning Strap	STD
<input type="checkbox"/>	SOAUTO80	80" Seat Positioning Strap	\$77.00

HEADREST ⁵⁵

<input type="checkbox"/>	SONUAHR	Using headrest mounting hardware other than SOSAH or SOMAH located below. Specify: _____	
<input type="checkbox"/>	SOFHRS	SEH4P0010 Small Curved Headrest Pad 4" x 9"	\$209.00
<input type="checkbox"/>	SOFHRM	SEH4P0010 Medium Curved Headrest Pad 5" x 10"	\$209.00
<input type="checkbox"/>	SOFHRL	SEH4P0010 Large Curved Headrest Pad 6" x 11"	\$209.00
<input type="checkbox"/>	SOSAH	Single-Axis Removable Mounting Hardware	\$148.00
<input type="checkbox"/>	SOMAH	Multi-axis Removable Mounting Hardware	\$168.00
<input type="checkbox"/>	SOHRBC	Bracket to Mount on the Contour Back Support (Without Recline for SOSAH or SOMAH)	\$102.00
<input type="checkbox"/>	SOHRSB	Bracket to Mount on the Seat Frame Crossbar (Without Recline for SOSAH or SOMAH)	\$102.00
<input type="checkbox"/>	SOHRPR	Bracket to Mount on the Seat Frame (With Power Recline for SOSAH or SOMAH)	\$102.00

Please select one of the headrest options with the mounting hardware and specify where the headrest will be installed (The frame or back cushion)

Note 55 The Headrest option is not compatible with the Rehab Seat with Corbeau Sport Back option (SAABA-COR3)

LATERAL TRUNK SUPPORTS ⁵⁶

<input type="checkbox"/>	SOSLS3	Adjustable Laterals with Small Curved Pads (Pair) - 3" x 5"	\$617.00
<input type="checkbox"/>	SOSLM4	Adjustable Laterals with Medium Curved Pads (Pair) - 4" x 6"	\$617.00
<input type="checkbox"/>	SOSLL5	Adjustable Laterals with Large Curved Pads (Pair) - 5" x 7"	\$617.00
<input type="checkbox"/>	SOSLHM	Adjustable Swing-Away Hardware Mounts - Each	\$306.00

Please select one of the lateral pads and the adjustable & removable hardware.

Note 56 The Lateral Trunk Support option is not compatible with the Rehab Seat with Corbeau Sport Back option (SAABA-COR3)

HIP & THIGH SUPPORTS

<input type="checkbox"/>	SOHTSS10	Adjustable Hip & Thigh Support Pads (Pair) - Small 10" x 4"	\$235.00
<input type="checkbox"/>	SOHTSM12	Adjustable Hip & Thigh Support Pads (Pair) - Medium 12" x 5"	\$235.00
<input type="checkbox"/>	SOHTSL14	Adjustable Hip & Thigh Support Pads (Pair) - Large 14" x 6"	\$235.00
<input type="checkbox"/>	SOHTSARH	Adjustable and Removable Hip & Thigh Support Hardware - Pair	\$444.00

Please select one of the hip & thigh support pads and the adjustable & removable hardware.

ACCESSORIES

<input type="checkbox"/>	SOCH	Cup Holder NEW!	\$92.00
<input type="checkbox"/>	SOCPH	Ram Mounts X-Grip Cell Phone Holder NEW!	\$158.00
<input type="checkbox"/>	SOTHS	Ram Mounts X-Grip Small Tablet Holder NEW!	\$179.00
<input type="checkbox"/>	SOTHL	Ram Mounts X-Grip Large Tablet Holder NEW!	\$270.00

SWITCHES AND OPTIONS

<input type="checkbox"/>	ASCDDT	Hand Held DTT Programmer Pad	\$1,244.00
<input type="checkbox"/>	ASCPCP	P&G Dongle and PC compatible	\$597.00
<input type="checkbox"/>	SSBS	Buddy Button Switch	\$112.00
<input type="checkbox"/>	SSES	Egg Switch-Black	\$112.00
<input type="checkbox"/>	SSML	Micro Light Switch	\$153.00



SPECIALTY DRIVER CONTROL OPTIONS

Therafin Sip-N-Puff Kit ⁶⁴

<input type="checkbox"/>	SCSNPK	57, 58	Therafin Sip-N-Puff Kit			\$1,556.00
<input type="checkbox"/>	SCSNPBTK	58	Therafin Sip-N-Puff Breath Tube Kit			\$505.00
<input type="checkbox"/>	SCSRMSNP	58	Stealth Rock-n-Lock Swingaway Arm for Sip-N-Puff	<input type="checkbox"/> SWL (Swings on Left)	<input type="checkbox"/> SWR (Swings on Right)	\$281.00
<input type="checkbox"/>	SCSRM	58	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/> SWL (Swings Left)	<input type="checkbox"/> SWR (Swings Right)	\$459.00
<input type="checkbox"/>	SSES-STD	58, 60	Egg Switch for Changing Modes - Black	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	\$112.00
<input type="checkbox"/>	SSBS-STD	58	Buddy Button Switch (Stop Switch) - Red			\$112.00
<input type="checkbox"/>	SOMAH	58	Multi-axis Headrest Removable Mounting Hardware			\$168.00
<input type="checkbox"/>	SOFHRS		Small Curved Headrest Pad 4" x 9"			\$209.00
<input type="checkbox"/>	SOFHRM	59	Medium Curved Headrest Pad 5" x 10"			\$209.00
<input type="checkbox"/>	SOFHRL		Large Curved Headrest Pad 6" x 11"			\$209.00

Switch-It Proximity Head Array ⁶⁴

<input type="checkbox"/>	SCSIHAA	57, 58	Adult - 3 Switch Proximity Head Array with 1/8" Mono Jack for mode/reverse			\$5,197.00
<input type="checkbox"/>	SCSRM	58	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/> SWL (Swings Left)	<input type="checkbox"/> SWR (Swings Right)	\$459.00
<input type="checkbox"/>	SSES-STD	58, 60	Egg Switch for Changing Modes - Black	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	\$112.00
<input type="checkbox"/>	SSBS-STD	58	Buddy Button Switch (Stop Switch) - Red			\$112.00
<input type="checkbox"/>	SOMAH	58	Multi-axis Headrest Removable Mounting Hardware			\$168.00

Switch-It Proximity Head Array & Sip-N-Puff Combo ⁶⁴

<input type="checkbox"/>	SCSICMBOA	57, 58	Adult - 4 Switch Combo Head Array with two proximity sensors for Left & Right and Low pressure Sip-N-Puff for forward and reverse and 1/8" jack for mode			\$5,197.00
<input type="checkbox"/>	SCSNPBTK	58	Therafin Sip-N-Puff Breath Tube Kit			\$505.00
<input type="checkbox"/>	SCSRMSNP	58	Stealth Rock-n-Lock Swingaway Arm for Sip-N-Puff	<input type="checkbox"/> SWL (Swings on Left)	<input type="checkbox"/> SWR (Swings on Right)	\$281.00
<input type="checkbox"/>	SCSRM	58	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/> SWL (Swings Left)	<input type="checkbox"/> SWR (Swings Right)	\$459.00
<input type="checkbox"/>	SSES-STD	58, 60	Egg Switch for Changing Modes - Black	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	\$112.00
<input type="checkbox"/>	SSBS-STD	58	Buddy Button Switch (Stop Switch) - Red			\$112.00
<input type="checkbox"/>	SOMAH	58	Multi-axis Headrest Removable Mounting Hardware			\$168.00

Switch-It MicroPilot Light Touch Proportional Mini joystick for chin / lip control ⁶⁴

<input type="checkbox"/>	SCSIMP	57, 58	MicroPilot Light touch mini Proportional Joystick			\$4,075.00
<input type="checkbox"/>	SCSRMMP	58	Stealth Rock-n-Lock Swingaway Mount for Chin M-P	<input type="checkbox"/> SWL (Swings Left)	<input type="checkbox"/> SWR (Swings Right)	\$505.00
<input type="checkbox"/>	SCSRM	58	Stealth Rock-n-Lock Swingaway Mount for Switch	<input type="checkbox"/> SWL (Swings Left)	<input type="checkbox"/> SWR (Swings Right)	\$459.00
<input type="checkbox"/>	SSES	58, 60	Egg Switch-Black	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	\$112.00
<input type="checkbox"/>	JCC		Chin Cup			\$51.00
<input type="checkbox"/>	SOMAH	58	Multi-axis Headrest Removable Mounting Hardware			\$168.00
<input type="checkbox"/>	SOFHRS		Small Curved Headrest Pad 4" x 9"			\$209.00
<input type="checkbox"/>	SOFHRM	59	Medium Curved Headrest Pad 5" x 10"			\$209.00
<input type="checkbox"/>	SOFHRL		Large Curved Headrest Pad 6" x 11"			\$209.00

Proportional Mini joystick for Midline Mount

<input type="checkbox"/>	SCPGMJS	57, 58	PG Drives JC200 Proportional Compact Remote Joystick for Midline Mount			\$969.00
<input type="checkbox"/>	SCGMJSM	61	Gatlin Midline Swing-away Mount for JC200	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	\$398.00
<input type="checkbox"/>	SCSIMP	57, 58	Switch-It MicroPilot Light touch mini Proportional Joystick			\$4,075.00
<input type="checkbox"/>	SCGMPPM	62	Gatlin Midline Swing-away Mount for MicroPilot	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	\$434.00
<input type="checkbox"/>	SCBTM	63	Bullet Hand Tray for MicroPilot	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	\$459.00

Proportional Mini joystick for Armrest Mount

<input type="checkbox"/>	SCPGMJS	57, 58	PG Drives JC200 Proportional Compact Remote Joystick for Armrest Mount			\$969.00
<input type="checkbox"/>	SCARMJSM	61	Swing-away armrest mount for JC200	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	\$398.00
<input type="checkbox"/>	SCSIMP	57, 58	Switch-It MicroPilot Light touch mini Proportional Joystick			\$4,075.00
<input type="checkbox"/>	SCARMPM	62	Swing-away and Multi-Axis armrest mount	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	\$434.00
<input type="checkbox"/>	SCBTM	63	Bullet Hand Tray for MicroPilot	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	\$459.00

Proportional RIM Head Control Package ⁶⁴

<input type="checkbox"/>	SCPHC	57, 58	Proportional RIM Head Control Package			\$2,545.00
<input type="checkbox"/>	SCSRM	58	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/> SWL (Swings Left)	<input type="checkbox"/> SWR (Swings Right)	\$459.00
<input type="checkbox"/>	SSES-STD	58, 60	Egg Switch for Changing Modes - Black	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	\$112.00
<input type="checkbox"/>	SOMAH	58	Multi-axis Headrest Removable Mounting Hardware			\$168.00

- Note 57 Requires the Omni Enhanced Display & Digital Interface Box (EDIB)
- Note 58 This option must be selected for the Kit
- Note 59 Choose one (1) Headrest Pad
- Note 60 The Mode Switch Type can be changed from this selection on Page #6
- Note 61 Must be selected for use with the PG Drives JC200 Mini Joystick
- Note 62 Must be selected for use with the Switch-It MicroPilot
- Note 63 This option is only compatible with the Switch-It MicroPilot
- Note 64 This option is not compatible for use with the Rehab Seat with Corbeau Sport Back (SAABA-COR3)



CUSTOM SECTION ⁶⁵

Note 65 Lead times will increase depending on the custom order pipeline at the time the purchase order is received
Please call for Quote and specifications for these options.

O² Holder

- SUPOXY **NEW!** 66 Custom Non-Articulating Oxygen Tank Holder MLH (Mount on Left) MRH (Mount on Right) \$398.00
- SUPOXYA **NEW!** 66 Custom Articulating Oxygen Tank Holder MLH (Mount on Left) MRH (Mount on Right) \$1,867.00
 - D-Size Oxygen Tank (* 22" Tall x 4.3" Diameter) * The height includes the valve
 - E-Size Oxygen Tank (* 31" Tall x 4.3" Diameter) * The height includes the valve

Note 66 The O² Holder option is not compatible for use with the Rehab Seat with Corbeau Sport Back (SAABA-COR3)

VENT TRAY (Please complete in Full)

- SOVENT **NEW!** 67, 68, 69 Custom Vent Tray - Articulating \$2,499.00
 - Articulating Vent Box (Bi-Pap) Model: _____
 - Laptop Style (LTV Style) Model: _____
- SOVTBE 67, 68, 69 Extended Battery Compartment N/C
- The Vent is being powered by 1 U1 Battery Chair Powered by 2 Group 22NF Batteries Chair Powered by 2 Group 24 Batteries
- The Vent is being powered by 1 Group 22NF Battery Chair Powered by 2 Group 22NF Batteries Chair Powered by 2 Group 24 Batteries
- The Vent is being powered by 1 Group 24 Battery Chair Powered by 2 Group 22NF Batteries

Note 67 Amysystems does not supply the battery charger or the wiring for the Vent.
These items are normally provided by the Vent Manufacturer.

Note 68 The Vent Tray option is not compatible for use with the Rehab Seat with Corbeau Sport Back (SAABA-COR3)

Note 69 The Vent Tray option is not available with the CUSTOM Battery Tray for Group 24 Batteries with a 24" frame width (FS2424B)

OTHER CUSTOM OPTIONS (Request For Quote)

- PAT 70, 71 Anterior Tilt RFQ
- CUSTOM 70, 71 15° Lateral Tilt Raise the Left Side only Raise the Right Side only Raise the Left & Right Side RFQ
- CUSTOM 70, 71 Elevating Stretcher (Please call CSR for availability with configuration and user weight capacity) RFQ
- CUSTOM 70 Power Flip-up Foot Platform (Please call CSR for availability with configuration and user weight capacity) RFQ
- CUSTOM 70 Vertical Lifting Foot Platform (Please call CSR for availability with configuration and user weight capacity) RFQ
- CUSTOM 70, 71, 72, 73 Stander (Requires Power Tilt, Power Recline and 90° Power elevating & articulating center mount foot platform) RFQ
- CUSTOM Manual Wheel Locks \$383.00
- CUSTOM Non Standard Frame Colors (Not all custom colors are available, please contact CSR for availability) \$765.00
- CUSTOM Red & White Reflector Kit Stickers (They will not be installed - Included in a bag with the owner's manual) N/C
- CUSTOM # 1 70 Other custom options. Specify in notes section below. RFQ
- CUSTOM # 2 70 Other custom options. Specify in notes section below. RFQ

Note 70 Please specify details and specifications for these options.

Note 71 The Maximum User Weight Capacity is 250 Lbs.

Note 72 Please call CSR to verify if the power standing option is available with the required configuration for your user.

Note 73 This option is not available with the Rehab Seat and Corbeau Sport Back (SAABA-COR3)

SEAT TO FLOOR HEIGHT SPECIFICATION (If not standard)

- SOMHSS Specific seat to floor height requirement other than standard. RFQ
Custom upcharges may apply. Required at: _____ " (± 1/4" Range)

FRAME COLOURS

STANDARD FRAME COLOURS

- FCBLK Black N/C
- FCCBU Blue Pearl N/C
- FCCRE Red Pearl N/C
- FCBAO Black Opal Metallic N/C
- FCGMC Gun Metal (Charcoal Grey) **NEW!** N/C
- FCSVM Silver Metallic **NEW!** N/C
- FCSBM Stealth Black Matte **NEW!** N/C

OPTIONAL FRAME COLOURS

- FCTXG Toxic Green N/C
- FCBOR Burnt Orange Metallic N/C
- FCMLP Midnight Purple N/C
- FCRZB Razzleberry N/C
- FCOPM Orange Pearl Metallic N/C
- FCYLW Yellow N/C

MATCHING ACCESSORIES TO THE SHROUD COLOUR

- FC***-M Matching Accessories to the Shroud Colour (Wheels, front access door & headrest mount when selected) \$235.00