

Account: _____
 Date: _____
 PO#: _____
 Buyer: _____
 Marked For: _____

Ship To: _____
 Name: _____
 Address: _____
 City, Province, Postal Code: _____
 Phone: _____

Option #	Definition	Price	Option #	Definition	Price
----------	------------	-------	----------	------------	-------

1. QUICKIE QM-7 Series Multi-Power Seating (Select One)

QUICKIE QM-710™ MPC for SEDEO Ergo

NEW!

<input type="checkbox"/>	Quickie QM-710™ MPC Tilt/Recline For SEDEO Ergo				\$8,060.00
	172PC114	Base Package: Includes SEDEO Ergo seating, LED lights and indicators,			\$8,060.00

QUICKIE QM-715™ HD MPC for SEDEO Ergo

NEW!

<input type="checkbox"/>	Quickie QM-715™ HD MPC Tilt/Recline For SEDEO Ergo				\$8,575.00
	172PC115	Base Package: Includes SEDEO Ergo seating, LED lights and indicators for user weight 300-350lbs,			\$8,575.00

QUICKIE QM-720™ MPD for SEDEO Ergo

NEW!

<input type="checkbox"/>	Quickie QM-720™ MPD Tilt/Recline				\$9,013.00
	172PC116	Base Package: Includes SEDEO Ergo seating, LED lights and indicators, 8mph Max speed.			\$9,013.00

MPC Seating Specifications

Typical Chair Weight	W/battery	375lbs	W/O Battery	285lbs
Transit	ISO Transit			
Seat Type	SEDEO Ergo			
Seat size range	Width	16"-22"	Depth	16"-22"
STFH range	MPC	16.4" - 19.4"		
Power Seating	50° CG tilt and Power recline w/ 7" Mechanical shear reduction / 12" Seat elevator			



1.1 QM-710™ MPC

QUICKIE QM-7 Series Base Specifications

Drive Wheel Position	Midwheel Drive	Drive Wheels	14" Drive Wheels w/ flat free insert (Black)
User weight Cap	300-350 lb	Tires/Casters	6" Front / 6" rear (Black)
Base size	Width 25" Length 36" Turn Radius 20"	Charger	8A Dual Mode Charger (Compact) CTE
Drive Control	R-Net PM120 (Expandable)	Battery Size	Group 24 / Group 34.
Driving Range	Up to 13.5 miles	Color Accents	Matte Black with 6 color accents
Drive Motors	4-pole 6 mph, Optional 5mph and 8mph	Ground Clearance	4"



Lights and Indicators

2. BASE CHAIR OPTIONS

Speed Package (4-pole 300lb)

<input type="checkbox"/>	172PF85	-	QM-710 Standard 6 mph motor package	N/C
<input type="checkbox"/>	172PF91	-	QM-710 5 mph High Torque motor pkg	\$541.00
<input type="checkbox"/>	172PF90	-	QM-710 8 mph Community motor package	\$1,136.00
<input type="checkbox"/>	172PF84	-	QM-715 HD Standard 5 mph motor pkg	N/C
<input type="checkbox"/>	172PF89	-	QM-720 Standard 8 mph motor package	N/C

Suretrac Drive Tracking **NEW!**

<input type="checkbox"/>	172MI01	-	Suretrac Module and Motor Encoders	\$1,519.00
--------------------------	---------	---	------------------------------------	------------

Color Accents **NEW!**

<input type="checkbox"/>	172EZ19	-	Matte Black Shroud w/ White Accents	N/C
<input type="checkbox"/>	172EZ20	-	Matte Black Shroud w/ Apple Green Acc	N/C
<input type="checkbox"/>	172EZ21	-	Matte Black Shroud w/ Blue Accents	N/C
<input type="checkbox"/>	172EZ22	-	Matte Black Shroud w/ Black Accents	N/C
<input type="checkbox"/>	172EZ23	-	Matte Black Shroud w/ Orange Accents	N/C
<input type="checkbox"/>	172EZ24	-	Matte Black Shroud w/ Red Accents	N/C

User weight

<input type="checkbox"/>	172WG06	-	Less than 150lbs	N/C
<input type="checkbox"/>	172WG13	-	151-200lb	N/C
<input type="checkbox"/>	172WG14	-	201-250lb	N/C
<input type="checkbox"/>	172WG03	-	251-300lb NA with 16" Seat depth and seat elevator	N/C
<input type="checkbox"/>	172WG15	-	301-350lb	N/C

Drive Wheel

<input type="checkbox"/>	172DW59	-	14" Cast Wheel Solid Dual Density Inserts Black Knobby Tires and Black Caster	N/C
--------------------------	---------	---	--	-----

Transit Option

<input type="checkbox"/>	172TO25	-	ISO Transit Package (NA with QM-715)	\$384.00
<input type="checkbox"/>	172TO40	-	Omit Transit	NC

Wheel locks

<input type="checkbox"/>	172WL42	-	Wheel locks	\$146.00
--------------------------	---------	---	-------------	----------

Ventilator Option **NEW!**

<input type="checkbox"/>	172VT10	-	Fixed Vent Tray (Pulmonetic LTV Models)	\$546.00
--------------------------	---------	---	---	----------

O2 Bottle Option **NEW!**

<input type="checkbox"/>	172A702	-	O2 Bottle mount to backrest D and E-style	\$350.00
--------------------------	---------	---	---	----------

Pre-Set Seat Frame Angle

<input type="checkbox"/>	172FA10	-	0°	N/C
--------------------------	---------	---	----	-----

Seat Height (Use STFH Matrix below)

(Measured in the front) - Obtain Quote to verify est. STFH

<input type="checkbox"/>	172HT44	-	Position "0"	N/C
<input type="checkbox"/>	172HT49	-	+5" STFH	Requires Group 34 Batteries
<input type="checkbox"/>	172HT45	-	+1" STFH	N/C
<input type="checkbox"/>	172HT46	-	+2" STFH	N/C
<input type="checkbox"/>	172HT47	-	+3" STFH	N/C

LED lights

<input type="checkbox"/>	172LI06	-	Include LED lights and indicators	NEW! N/C
<input type="checkbox"/>	172LI02	-	Omit LED lights and indicators	N/C

STFH Matrix					
Seating Type and Power Option Combinations	Seat-to-Floor Height at 0 tilt				
	0	+5"	+1"	+2"	+3"
Ultra-Low STFH Seat Frame "A" <i>(Upcharge see Section 6)</i> Reduces STFH by 1" from standard. Raises recline pivot point by 2" Armrest and backrest also raised 2"	MPC Power Recline and Tilt.	16.5"	17"	18"	19"
	MPC Power Recline and HD Tilt.	16.8"	17.3"	18.3"	19.3"
MPC Power Recline with Tilt and 12" Lift combo	16.4"	16.9"	17.4"	18.4"	19.4"
Standard STFH Seat Frame "B" Seat pan is in line with recline pivot for standard STFH	MPC Power Recline and Tilt.		*17.5"	18"	19"
	MPC Power Recline and HD Tilt.		*17.8"	18.3"	19.3"
	MPC Power Recline with Tilt and 12" Lift combo	17.4"	17.9"	18.4"	19.4"

* Seat Depth Limited to 19" with the STFH configuration
STFH measured with batteries but no user. Actual STFH may vary by .5"

3. BATTERIES and CHARGERS

Batteries

<input type="checkbox"/>	172AG04	-	GP24 MK Gel Flat Top Batts Pair w/strap Recommended for Max Daily Range	\$855.00
<input type="checkbox"/>	172AG22	-	GP24 Brightway AGM Flat Top Batteries Pair w/strap.	\$819.00
<input type="checkbox"/>	172AG11	-	GP34 MK Gel Flat Top Batteries (Pair) w/straps. <i>Required for lower STFH.</i>	\$855.00
<input type="checkbox"/>	172AG23	-	GP34 Brightway AGM Flat Top Batts (Pair) w/straps. <i>Required for lower STFH.</i>	\$773.00

Battery Charger

<input type="checkbox"/>	172BJ19	-	8A Dual Mode (Compact) CTE	N/C
--------------------------	---------	---	----------------------------	-----



Color Accents

4. ELECTRONICS

Expandable Controllers			
<input type="checkbox"/>	172ER08	R-Net Expandable Drive Control (120A)	\$3,781.0
	172EX08	Expandable 120A Controller - R-net Includes Memory Seating Positions	
	172EX07	Harness for 120A Expandable Controller	

Memory Seating NEW!			
"Memory Seating" feature allows for 6 custom seating positions that are accessed either thru the joystick, Ctrl+S, Joystick Shortcut keys or programmable jacks on OMNI and joystick. When enabled, chair will be pre-programmed with: Seated, Supine, Pressure Relief, Feet Above Heart Transfer and C-ME. Settings are dealer settable using mono switch (not included).			
<input type="checkbox"/>	172OG01	- Include Memory Seating	N/C
<input type="checkbox"/>	172OG02	- Omit Memory Seating	N/C

R-net Joystick			
<input type="checkbox"/>	172J113	CJSM: Joystick w/Stereo Jacks & Color Display with Light and indicator controls.	N/C
<input type="checkbox"/>	172J112	- CJSM2: Joystick Advanced: 3.5" Color Display, Paddle Switches, Programmable Shortcut Keys, ICONS, Stereo Jack and IR control	\$1,167.00
<input type="checkbox"/>	172J114	- CJSM2: Joystick Advanced: 3.5" Color Display, Paddle Switches, Programmable Shortcut Keys, ICONS, Stereo Jack, IR and dual (IOS/Android) Bluetooth control	\$1,167.00
<input type="checkbox"/>	172J95	- Omit JS. For use with OMNI	N/C

Joystick Mount (Fixed, Retractable, Mid)			
<input type="checkbox"/>	172JM16	Fixed Mount for Ergo	N/C
<input type="checkbox"/>	172JM17	Retractable Joystick mount for Ergo	\$743.00
<input type="checkbox"/>	172JM14	Mid Line Mount	\$400.00

Joystick Mount (Location)			
<i>Height Adjustable.</i>			
<input type="checkbox"/>	172HM1	- Right Side Mount	N/C
<input type="checkbox"/>	172HM2	- Left Side Mount	N/C

Joystick Handle (Optional)			
<input type="checkbox"/>	172JH1	Ball Handle	\$98.00
<input type="checkbox"/>	172JH2	Foam Ball Handle	\$98.00
<input type="checkbox"/>	172JH4	T-Handle	\$98.00
<input type="checkbox"/>	172JH6	Goalpost	\$98.00

Electronic Programmers			
<input type="checkbox"/>	172PR5	R-net PC Programmer with Dongle	\$633.00

R-net OMNI and Output Options			
<i>(For expandable controls, choose all that apply)</i>			
<input type="checkbox"/>	172ON6	- Rnet OMNI Input Control Module and display with built in Infrared and armrest mount. Choose mount below. (required for Specialty Input Devices.)	\$2,137.00
<input type="checkbox"/>	172ON7	Rnet Output Module (ECM)	\$1,098.00
<input type="checkbox"/>	172ON8	- Rnet Bluetooth Mouse-mover for Windows / Android	\$1,098.00
<input type="checkbox"/>	172ON9	- Rnet Bluetooth for Apple Devices	\$1,098.00

Attendant Control			
<input type="checkbox"/>	172OA04	Attendant Control for R-net	\$1,098.00

Attendant Mounting location			
<input type="checkbox"/>	172DN1	Mounting location: Right Back Cane	N/C
<input type="checkbox"/>	172DN2	Mounting Location: Left Back Cane	N/C

OMNI Mount			
<i>(If used with Joystick must be mounted opposite side)</i>			
<input type="checkbox"/>	172DS04	Swing Away Left Side OMNI Mount	NC
<input type="checkbox"/>	172DS03	Swing Away Right Side OMNI Mount	NC

Switch inputs for Joysticks or OMNI			
<i>Switches can be used with "Assignable Button" feature to control Power, Horn, Mode, Speeds, Profile, Seating, IR, Bluetooth, Lighting, Output Module, Memory Settings and Pressure Relief Reminder Mute.</i>			
<input type="checkbox"/>	172SW22	Stereo to 2 Mono Splitter	\$77.00
<input type="checkbox"/>	172SW6	Ribbon Switch	\$123.00
<input type="checkbox"/>	172SW1	Buddy Button Switch	\$98.00
<input type="checkbox"/>	172SW2	Disc Switch	\$208.00
<input type="checkbox"/>	172SW9	Wobble Switch	\$276.00
<input type="checkbox"/>	172SW4	Micro Light Switch	\$103.00
<input type="checkbox"/>	172SW24	Egg Switch-Black	\$98.00
<input type="checkbox"/>	172SW21	Dual toggle switch and bracket	\$183.00
<input type="checkbox"/>	172SW35	Switch-it Mini 1 button on Mono plug	\$276.00
<input type="checkbox"/>	172SW36	Switch-it Mini 2 button on 2 Mono plugs	\$276.00
<input type="checkbox"/>	172SW37	Switch-it Mini 2 button on 1 Stereo plug	\$276.00
<input type="checkbox"/>	172SW38	Switch-it Mini 4 button on DB9 connector	\$276.00

Switch Mount Bracket			
<input type="checkbox"/>	172EM01	Single Switch Mount Bracket	\$211.00

Default Memory Positions



M1: Seated

M2: Supine

M3: Pressure Relief

M4: Feet above Heart

M5: Transfer

M6: C-ME

5. Specialty Driver Control Options

Requires expandable Rnet controls and OMNI unless otherwise noted.

5a. Proportional Control Packages

Head Control Proportional Packages		
□ Switch-It Dual Pro Proportional Head Control with Whitmyer and Facial Pads (Adult) \$6,737.00		
172ID72	Switch-It Dual Pro Head Control w/facial pad (Adult)	\$6,323.00
172J09	SA Boom Mount with Egg Switch Control	\$414.00
□ Switch-It Dual Pro Proportional Head Control with Spot Pads Package \$6,737.00		
172ID70	Switch-It Dual Pro Head Control w/Spot pads	\$6,323.00
172J09	SA Boom Mount with Egg Switch Control	\$414.00

Hand Control Proportional Joystick Packages

□ Switch-It Versa-Guide Compact Joystick Midline with Gatlin Package \$2,100.00		
172ID71	Switch-It Versa-Guide Compact Proportional Joystick	\$1,650.00
172J07	S/A Gatlin Mount	\$340.00
□ Switch-It Versa-Guide Compact Joystick Midline and prox flex switches with Gatlin Package \$3,108.00		
172ID71	Switch-It Versa-Guide Compact Proportional Joystick	\$2,100.00
172IM08	Dual Flex Proximity Switches (requires no force)	\$668.00
172J07	S/A Gatlin Mount	\$340.00
□ Switch-It Micro-Pilot Joystick with Space Disc Package \$4,127.00		
172ID62	Micro-Pilot Proportional Joystick	\$3,640.00
172J08	S/A SpaceDisc mount	\$487.00
□ Switch-It Micro-Guide Joystick with SpaceDisc Package \$4,489.00		
172ID63	Micro-Guide Proportional Joystick	\$4,489.00
172J08	S/A SpaceDisc mount	
□ ASL Proportional Micro Extremity Joystick (ASL130) Midline with Gatlin Package \$4,547.00		
172ID59	Micro Extremity Joystick	\$4,207.00
172J07	S/A Gatlin Mount	\$340.00
□ ASL Mushroom Joystick (ALS135) Midline with Gatlin Package \$3,454.00		
172ID60	Mushroom Joystick	\$3,114.00
172J07	S/A Gatlin Mount	\$340.00
□ ASL Compact Joystick w Integrated buttons (ASL134) Midline with Gatlin Package \$1,576.00		
172ID61	ASL Compact Joystick	\$1,236.00
172J07	S/A Gatlin Mount	\$340.00
□ ASL Micro Mini Joystick w Hand Pad (ASL136) Midline with Gatlin Package \$4,132.00		
172ID68	Micro Mini Joystick	\$3,792.00
172J07	S/A Gatlin Mount	\$340.00



Dual Pro Head Control

Chin Control Proportional Joystick Packages		
□ Switch-It Versa-Guide Chin Control Compact Joystick and Egg Switch Package \$2,478.00		
172ID71	Switch-It Versa-Guide Compact Proportional Joystick	\$2,100.00
172IM07	S/A Mount off backrest	\$260.00
172J04	SA Boom Mount with egg switch	\$118.00
□ Switch-It Versa-Guide Chin Control Compact Joystick and Proximity Flex Switches Package \$3,028.00		
172ID71	Switch-It Versa-Guide Compact Proportional Joystick	\$2,100.00
172IM08	Dual Flex Proximity Switches (requires no force)	\$668.00
172IM07	S/A Mount off backrest	\$260.00
□ Switch-It Versa-Guide Bib Chin Control Compact Joystick with Proximity Flex switches Package Std 6" \$3,548.00		
172ID71	Switch-It Versa-Guide Compact Proportional Joystick	\$2,100.00
172IM08	Dual Flex Proximity Switches (requires no force)	\$668.00
172IM05	Bib mount Small	\$780.00
□ Switch-It Versa-Guide Bib Chin Control Compact Joystick with Proximity Flex Switches Package Large 8" \$3,617.00		
172ID71	Switch-It Versa-Guide Compact Proportional Joystick	\$2,100.00
172IM08	Dual Flex Proximity Switches (requires no force)	\$668.00
172IM06	Bib mount Large	\$849.00
□ Switch-It Versa-Guide Bib Chin Control Compact Joystick with Egg Switch Package Std 6" \$2,998.00		
172ID71	Switch-It Versa-Guide Compact Proportional Joystick	\$2,100.00
172IM05	Bib mount Small	\$780.00
172J04	SA Boom Mount with egg switch	\$118.00
□ Switch-It Versa-Guide Bib Chin Control Compact Joystick with Egg Switch Package Large 8" \$3,067.00		
172ID71	Switch-It Versa-Guide Compact Proportional Joystick	\$2,100.00
172IM06	Bib mount Large	\$849.00
172J04	SA Boom Mount with egg switch	\$118.00
□ Switch-It Micro-Pilot Bib Chin Control with Egg Switch Package Std 6" \$4,538.00		
172ID62	Micro-Pilot Proportional Joystick	\$3,640.00
172IM05	Bib mount Small	\$780.00
172J04	SA Boom Mount with egg switch	\$118.00
□ Switch-It Micro-Pilot Bib Chin Control with Egg Switch Package Large 8" \$4,607.00		
172ID62	Micro-Pilot Proportional Joystick	\$3,640.00
172IM06	Bib mount Large	\$849.00
172J04	SA Boom Mount with egg switch	\$118.00
□ Switch-It Micro-Guide Bib Chin Control with Egg Switch Package Std 6" \$5,387.00		
172ID63	Micro-Guide Proportional Joystick	\$4,489.00
172IM05	Bib mount Small	\$780.00
172J04	SA Boom Mount with egg switch	\$118.00
□ Switch-It Micro-Guide Bib Chin Control with Egg Switch Package Large 8" \$5,456.00		
172ID63	Micro-Guide Proportional Joystick	\$4,489.00
172IM06	Bib mount Large	\$849.00
172J04	SA Boom Mount with egg switch	\$118.00
□ ASL Proportional Micro Extremity Joystick Chin / Lip with Egg Switch Control Package \$4,585.00		
172ID59	Micro Extremity Joystick	\$4,207.00
172IM07	S/A Mount off backrest	\$260.00
172J04	SA Boom Mount with egg switch	\$118.00

5. Specialty Driver Control Options (cont.)

5b. Switch Control

Requires Headrest and/or Headrest Mount

3 Switch Proximity Head Arrays Choose headrest mount below

<input type="checkbox"/> ASL Proximity Head Array Package (Adult) \$5,928.00		
172II02	ASL Proximity Head Array (Adult)	\$5,514.00
172IJ09	SA Boom Mount with Egg Switch Control	\$414.00
<input type="checkbox"/> ASL Proximity Head Array Package (Ped) \$5,928.00		
172II03	ASL Proximity Head Array (Ped)	\$5,514.00
172IJ09	SA Boom Mount with Egg Switch Control	\$414.00
<input type="checkbox"/> Switch-It Head Array Package 3 Switch Curved (Long) \$5,928.00		
172II13	Switch-It Proximity Head Array Long	\$5,514.00
172IJ09	SA Boom Mount with Egg Switch Control	\$414.00
<input type="checkbox"/> Switch-It Head Array Package 3 switch curved (Stubby) \$5,862.00		
172II14	Switch-It Proximity Head Array Stubby	\$5,448.00
172IJ09	SA Boom Mount with Egg Switch Control	\$414.00

5 Switch Proximity Head Array (No Egg Switch Needed) Choose "Headrest Mount for Special Controls"

<input type="checkbox"/>	172II04	E2330	Switch-It Head Array 5 switch flat	\$5,514.00
--------------------------	---------	-------	------------------------------------	------------

Multiple Switches.


<input type="checkbox"/>	172SW31	Switch-It 3 Naked proximity switches on DB9 connector. (requires OMNI)	\$3,130.00
<input type="checkbox"/>	172SW32	Switch-It 4 Naked proximity switches on DB9 connector. (requires OMNI)	\$3,830.00
<input type="checkbox"/>	172SW33	Switch-It 5 Naked proximity switches on DB9 connector. (requires OMNI)	\$4,270.00
<input type="checkbox"/>	172SW34	ASL (502) 5-Switch adaptor box on DB9 connector. (for use with 5 mono switches, not included)	\$4,270.00
<input type="checkbox"/>	172SW23	4-way toggle switch button DB9 Connector. Includes bracket	\$690.00

Sip N Puff Package Choose headrest and HR mount below

<input type="checkbox"/> Quickie Sip n' Puff System \$865.00		
172ID34	Sip N' Puff Kit	\$276.00
172IH1	Sip n' Puff Tubing and SA Hardware	\$471.00
172IJ10	SA Boom Mount with egg switch	\$118.00

Headrest Mount for Specialty Controls

NEW!

<input type="checkbox"/>	172MH23	Whitmyer Linx Headrest Mount, Quick Release	 \$500.00
--------------------------	---------	---	--

Whitmyer Headrest Pad for Specialty controls

<input type="checkbox"/>	172HE79	Plush Pad 8" for Specialty controls	\$318.00
<input type="checkbox"/>	172HE80	Plush Pad 10" for Specialty controls	\$318.00
<input type="checkbox"/>	172HE81	Plush Pad 14" for Specialty controls	\$318.00
<input type="checkbox"/>	172HE90	Adjusta-A-Plush Headrest Pad Narrow	\$318.00
<input type="checkbox"/>	172HE91	Adjusta-A-Plush Headrest Pad Standard	\$318.00

Extra Swing Away Head Boom Egg Switch and Mount

<input type="checkbox"/>	172IJ05	Extra SA Boom Mount with egg switch #1	\$118.00
<input type="checkbox"/>	172IJ06	Extra SA Boom Mount with egg switch #2	\$118.00



Single Switches

<input type="checkbox"/>	172IB01	Ribbon Switch	\$123.00
<input type="checkbox"/>	172IB02	Buddy Button Switch	\$111.00
<input type="checkbox"/>	172IB03	Disc Switch	\$208.00
<input type="checkbox"/>	172IB04	Wobble Switch	\$649.00
<input type="checkbox"/>	172IB05	Micro Light Switch	\$98.00
<input type="checkbox"/>	172IB06	Egg Switch-Black	\$138.00
<input type="checkbox"/>	172IB07	Switch-It Mini 1 button on Mono plug	\$276.00

6. SEAT TYPE

<input type="checkbox"/>	172ST63	-	SEDEO Ergo Recline	N/C	<input type="checkbox"/>	172FU30	-	Ultra-Low STFH Seat Frame "A"	\$361.00
<input type="checkbox"/>					<input type="checkbox"/>	172FU31	-	Standard Seat Frame "B"	N/C

7. SEAT SIZE



Critical Measurements	
A: Seat Pan Width (Section 7):	Measured on the seat pan from inside track rail to inside track rail, Infinite adjustable from 16"-22" wide.
B: Seat Cushion Width (Section 8):	When ordered with ERGO thigh supports, the seat cushion will be 2" narrower than selected seat width as ERGO thigh supports consume 2" of effective seat width. JAY Bundled cushion will be shipped actual selected size so be sure to order appropriately if combined with ERGO thigh supports.
C: ERGO Thigh Support (Bolster - Section 9) Width:	ERGO thigh supports measure 3.5" wide each, plus they sit on top of the outer track rails and consume 2" of effective seat width. ERGO thigh supports and track rail unlatch, swing away and lift off for lateral transfers.
D: Seat Depth (Section 7):	Measured from back shell surface to front of seat pan. Most backs cushions consume 1" of effective seat depth. Seat Pan has Infinite adjustable from 16"-22" deep.

Seat Width (Pan Width)			
<input type="checkbox"/>	172R5	16"	Seat width (N/C)
<input type="checkbox"/>	172R6	17"	Seat width (N/C)
<input type="checkbox"/>	172R7	18"	Seat width (N/C)
<input type="checkbox"/>	172R8	19"	Seat width (N/C)
<input type="checkbox"/>	172R9	20"	Seat width (N/C)
<input type="checkbox"/>	172R52	21"	Seat width (N/C)
<input type="checkbox"/>	172R25	22"	Seat width (N/C)

Seat Depth (Section 7)			
<input type="checkbox"/>	172RF5	16"	Seat depth NA above 250lbs (N/C)
<input type="checkbox"/>	172RF6	17"	Seat depth (N/C)
<input type="checkbox"/>	172RF7	18"	Seat depth (N/C)
<input type="checkbox"/>	172RF8	19"	Seat depth (N/C)
<input type="checkbox"/>	172RF24	20"	Seat depth (N/C)
<input type="checkbox"/>	172RF26	21"	Seat depth (N/C)
<input type="checkbox"/>	172RF25	22"	Seat depth (N/C)

Some No-go's, see STFH Matrix

8. SEAT CUSHIONS

JAY Ion Cushion for SEDEO ERGO NEW!			
<input type="checkbox"/>	172CU45	JAY ION cushion with 3DX cover	\$458.00
<input type="checkbox"/>	172CU46	JAY ION cushion with 3DX and Black Stretch	\$458.00
<input type="checkbox"/>	172CU47	JAY ION cushion w/ Premium leatherette cover	\$458.00

9. POSTURAL SUPPORTS

LATERAL THORACIC SUPPORT (Swing Away)			
<i>Available on all Solid Backs.</i>			
<input type="checkbox"/>	LAT 4X4 S/A MOUNTS W/ 1" OFFSET - PR		\$421.00
	172ZV14	LATERAL 4X4 W/ 1" OFFSET - PAIR	\$130.00
	172ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$291.00
<input type="checkbox"/>	LAT 4X5 S/A MOUNTS W/ 1" OFFSET - PR		\$421.00
	172ZV15	LATERAL 4X5 W/ 1" OFFSET - PAIR	\$130.00
	172ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$291.00
<input type="checkbox"/>	LAT 4X6 S/A MOUNTS W/ 1" OFFSET - PR		\$421.00
	172ZV16	LATERAL 4X6 W/ 1" OFFSET - PAIR	\$130.00
	172ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$291.00
<input type="checkbox"/>	LAT 6X4 S/A MOUNTS W/ 1" OFFSET - PR		\$421.00
	172ZV17	LATERAL 6X4 W/ 1" OFFSET - PAIR	\$130.00
	172ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$291.00
<input type="checkbox"/>	LAT 6X5 S/A MOUNTS W/ 1" OFFSET - PR		\$421.00
	172ZV18	LATERAL 6X5 W/ 1" OFFSET - PAIR	\$130.00
	172ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$291.00
<input type="checkbox"/>	LAT 6X6 S/A MOUNTS W/ 1" OFFSET - PR		\$421.00
	172ZV19	LATERAL 6X6 W/ 1" OFFSET - PAIR	\$130.00
	172ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$291.00

CHEST STRAP			
<input type="checkbox"/>	172ZR03	E0960	Padded Chest Strap \$291.00

Ergo Thigh Supports <i>(Upholstered in 3DX unless back is ordered in Leatherette).</i>			
<input type="checkbox"/>	Ergo Thigh Support - Left (Seat Bolster)		\$345.00
	172TV5	Ergo Long Thigh Support - Left	\$108.00
	172TX3	Swing-away Mount - Left	\$237.00
<input type="checkbox"/>	Ergo Thigh Support - Right (Seat Bolster)		\$345.00
	172TV6	Ergo Long Thigh Support - Right	\$108.00
	172TX4	Swing-away Mount - Right	\$237.00

Long Thigh Support w/Removable Mount			
<input type="checkbox"/>	SPOT Long Thigh Support - Left		\$421.00
	172TV2	Long Thigh Support - Left	\$130.00
	172TX1	Removable Mount - Left	\$291.00
<input type="checkbox"/>	SPOT Long Thigh Support - Right		\$421.00
	172TV1	Long Thigh Support - Right	\$130.00
	172TX2	Removable Mount - Right	\$291.00

Swing Away Knee Adductor Supports			
Knee Adductor Supports - LEFT (choose Pad and Arm)			
Pad Size			
<input type="checkbox"/>	172KB5	Knee Adductor Pad SM- Left	\$221.00
<input type="checkbox"/>	172KB7	Knee Adductor Pad LG- Left	\$221.00
S/A Arm Size			
<input type="checkbox"/>	172XX04	Swing Away Short - Left	\$291.00
<input type="checkbox"/>	172XX06	Swing Away Long - Left	\$291.00
Knee Adductor Supports - RIGHT (choose Pad and Arm)			
Pad Size			
<input type="checkbox"/>	172KB6	Knee Adductor Pad SM- Right	\$221.00
<input type="checkbox"/>	172KB8	Knee Adductor Pad LG- Right	\$221.00
S/A Arm Size			
<input type="checkbox"/>	172XX05	Swing Away Short - Right	\$291.00
<input type="checkbox"/>	172XX07	Swing Away Long - Right	\$291.00

10. POWERED SEATING

Power Seating Options

NEW!

<input type="checkbox"/>	172PA28		50 Degree Tilt and Ergo Recline w/Shear Reduction and optional Biometric Repositioning (2 actuators)	\$12,824.00
<input type="checkbox"/>	172xxx		5° anterior tilt (reduces max posterior tilt to 45°)	\$5,514.00
<input type="checkbox"/>	172PA29		QM-715 HD 50 Degree Tilt and Ergo Recline w/Shear Reduction and optional Biometric Repositioning (350 weight limit) (2 actuators)	\$12,824.00

Biometric Repositioning Recline

NEW!

<input type="checkbox"/>	172OG03	-	Include Biometric Reposition recline. (Requires Memory positions.)	N/C
<input type="checkbox"/>	172OG04		Omit Biometric Reposition recline.	N/C

Powered Seat Elevator

(NA with QM-715 HD)

NEW!

<input type="checkbox"/>	172PA30		12" Power Seat Elevate (Optional C-ME)	\$3,116.00
C-ME (3mph drive speed with full elevate)				
<input type="checkbox"/>	172OG05	-	Include C-ME with Seat elevator (NA with QM-720 or QM-715HD)	N/C
<input type="checkbox"/>	172OG06	-	Omit C-ME with Seat elevator	N/C



Power Seating Control

Powered Seating Thru-drive Control Method.

Thru-Drive 2+ Actuator

<input type="checkbox"/>	172EO31		Thru Drive 2+ Actuators: (For operating tilt, recline, lift, power legs) Includes Ctrl+5 and optional PRAR	\$2,601.00
--------------------------	---------	--	--	------------

Powered Seating Thru-drive Control Options

Assignable Buttons

"Assignable Button" feature allows for customization of joystick buttons and/or Ctrl+5 Switches allowing access to various wheelchair and seating functions. When enabled, chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle.

<input type="checkbox"/>	172PM01	-	Enable Assignable buttons	N/C
<input type="checkbox"/>	172PM02		Disable Assignable buttons	N/C

PRAR Pressure Relief Alert and Recording

NEW!

PRAR Pressure Relief Alert and Record feature provides programmable alerts to be set into the Rnet electronics that will prompt the user when it is time to perform a pressure relief. The met system will record how often the pressure relief is completed or ignored. Requires Memory Seating function.

<input type="checkbox"/>	172OG07	-	Include PRAR	N/C
<input type="checkbox"/>	172OG08		Omit PRAR	N/C

Ctrl+5 Switchbox Options

Ctrl+5 allows direct control to 5 power seating functions (10 total switches).

Switches can also be used with assignable buttons and Memory seating for greater flexibility. Box includes 3-Amp USB port for charging portable devices.

#1 Ctrl+5 Switch Location

<input type="checkbox"/>	172QG01	-	Ctrl+5 Switchbox mounted behind joystick	N/C
<input type="checkbox"/>	172QG02		Ctrl+5 Switchbox mounted opposite joystick	N/C
<input type="checkbox"/>	172QG03		Ctrl+5 Switchbox mounted on back cane	N/C

#1 Ctrl+5 Switch Type

<input type="checkbox"/>	172EY03	-	Push buttons	N/C
<input type="checkbox"/>	172EY04		Toggle switches	\$318.00

A 2nd Ctrl+5 can be added to provide additional control for Memory seating or environmental controls. It can also be mounted as an attendant control all functions of #1 are available.

#2 Ctrl+5 Switch Location

NEW!

<input type="checkbox"/>	172QG06	-	Ctrl+5 Switchbox mounted behind joystick	\$1,293.00
<input type="checkbox"/>	172QG07		Ctrl+5 Switchbox mounted opposite joystick	\$1,293.00
<input type="checkbox"/>	172QG08		Ctrl+5 Switchbox mounted on back cane	\$1,293.00

#2 Ctrl+5 Switch Type

<input type="checkbox"/>	172EY05	-	Push buttons	N/C
<input type="checkbox"/>	172EY06		Toggle switches	\$258.00



Icon on Ctrl+5 flashes Blue when it is time to do a pressure relief and / or

Joystick emits an Alert Tune



Ctrl+5 Push Button option



Ctrl+5 Toggle Switch option



11. BACKREST AND BACK OPTIONS FOR ERGO RECLINE

Backrest Frame Type			
<input type="checkbox"/>	172BT89	-	Power Recline (Power Recline Seat only)

N/C

Back Upholstery Type			
<input type="checkbox"/>	172BU47	-	Ergo Back Cushion: choose options below
<input type="checkbox"/>	172JB03	-	Jay J3 for Power Recline w/1.5" Foam Choose type below.
<input type="checkbox"/>	172BU14	-	Omit Back. Includes bracket for aftermarket backs

N/C

N/C

N/C

Push Handle Options			
<input type="checkbox"/>	172BP55	-	Fixed "Cane" Handles. NA with O2 holder
<input type="checkbox"/>	172BP56	-	Height Adjustable Push Handles. w/2 utility hooks
<input type="checkbox"/>	172BP57	-	Omit push handles

N/C

\$55.00

N/C

Power Recline Backrest Preset Home Angle N/C		
<input type="checkbox"/>	90° 172QK1	<input type="checkbox"/>
<input type="checkbox"/>	95° 172QK2 default	<input type="checkbox"/>
<input type="checkbox"/>	100° 172QK3	

N/C

11a. SEDEO ERGO BACK CUSHION AND OPTIONS

NEW!

SEDEO Ergo Back Type SEDEO Ergo Back coding: (xx=height, YY=width, zz=Upholstery) Back cushion reduces effective seat depth by 1"			
<input type="checkbox"/>	172OD21	Ergo Back Cushion: ERGOxxYYzz For 16" to 20" seat width	\$922.00
<input type="checkbox"/>	172OD22	Ergo Back Cushion: ERGOxxYYzz For 22" seat width	\$922.00

SEDEO Ergo Back Width (YY)			
<i>Back can be ordered a different width than seat pan. Armrest will be adjusted to back width ordered.</i>			
<input type="checkbox"/>	172BN18	-	16" Back width
<input type="checkbox"/>	172BN19	-	18" Back width
<input type="checkbox"/>	172BN20	-	20" Back width
<input type="checkbox"/>	172BN21	-	22" Back width

N/C

N/C

N/C

N/C

SEDEO Ergo Back Height (xx) (measured from seat pan to top of back shell)					
<input type="checkbox"/>	172B152	-	22"	NA with Ultra Low STFH Seat Frame	N/C
<input type="checkbox"/>	172B153	-	24"		N/C
<input type="checkbox"/>	172B154	-	26"		N/C
<input type="checkbox"/>	172B158	-	28"	NA with Standard Seat Frame	N/C





3DX Option



SEDEO Ergo Back Upholstery options (zz)					
<input type="checkbox"/>	172ZD76	-	sf	3DX™ spacer fabric (sf)	N/C
<input type="checkbox"/>	172ZD77	-	xs	3DX™ spacer fabric and Lycra® stretch outer layer (xs)	N/C
<input type="checkbox"/>	172ZD78	-	pl	Premium Leatherette with perforated inserts (pl)	N/C




Premium Leatherette Option

12. HEADREST OPTIONS						
Headrest Mount		Whitmyer Headrest Pad				
<input type="checkbox"/>	172MH18	Ergo Multi-position headrest mount, Quick release	<input type="checkbox"/>	172HE65	Whitmyer Plush Pad 8"	\$244.00
<input type="checkbox"/>	172MH22	Whitmyer Linx Headrest Mount, Quick Release 	<input type="checkbox"/>	172HE66	Whitmyer Plush Pad 10"	\$244.00
			<input type="checkbox"/>	172HE78	Whitmyer Plush Pad 14"	\$244.00
			<input type="checkbox"/>	172HE82	Whitmyer Adjusta-A-Plush Pad Narrow	\$376.00
			<input type="checkbox"/>	172HE83	Whitmyer Adjusta-A-Plush Pad Standard	\$376.00

13. ARMREST AND ARMPAD OPTIONS							
Armrest for Power Recline Seating. (Armrest width is adjustable from 16" to 22")			Reclining Arm pad Options (Arm pad can be angled inward 15 degrees)				
<input type="checkbox"/>	172EG29	 Cantilever / Reclining Arm Standard STFH 10"-14" Arm Height Range Ultra-low STFH 12"-16" Arm Height Range	N/C				
			Standard Waterfall Reclining Arm pad (Upholstered in 3DX unless back is ordered in Leatherette).				
			Left Arm Pad (N/C) Right Arm pad (N/C)				
			<input type="checkbox"/>	172AE29	St Reclining Arm Pad (16")	<input type="checkbox"/>	172AB29
			Comfort Reclining Arm pad (Foam)				
			Left Armpad (\$101.00) Right Armpad (101.00)				
			<input type="checkbox"/>	172AE22	Comfort 2" wide Reclining Arm Pad	<input type="checkbox"/>	172AB22
			<input type="checkbox"/>	172AE23	Comfort 4" wide Reclining Arm Pad	<input type="checkbox"/>	172AB23
			Comfort Reclining Arm pad (Gel)				
			Left Gel Armpad (\$191.00) Right Gel Armpad (191.00)				
			<input type="checkbox"/>	172AE25	Comfort 2" wide Reclining Arm Pad	<input type="checkbox"/>	172AB25
			<input type="checkbox"/>	172AE26	Comfort 4" wide Reclining Arm Pad	<input type="checkbox"/>	172AB26
			AES Trough/Elbow Block Reclining Arm pads (Gel)				
			Left Gel Armpad (\$266.00) Right Gel Armpad (\$266.00)				
			<input type="checkbox"/>	172AE27	Trough / Elbow Block Medium (Gel) (3.5" x 13") Reclining Arm Pad	<input type="checkbox"/>	172AB27
			<input type="checkbox"/>	172AE28	Trough / Elbow Block Large (Gel) (4" x 14.5") Reclining Arm Pad	<input type="checkbox"/>	172AB28

14. CENTERMOUNT FOOTREST				
Power Elevating Legrests		Footplates		
<input type="checkbox"/>	172H127	Power Centermount ELR w/7.5" Power Articulation (F) Centermount footplates plus Custom Options Below	\$3,077.00	
Power Centermount Footplate Extension Option (For Thru-Drive Operation)				
<input type="checkbox"/>	172E01	 7.5" Powered footplate extension for transfer assist. (Electronics programmed to allow footplates to be extended <u>with or without</u> legrest being raised)	\$530.00	
<input type="checkbox"/>	172E02	 Footplate extension for transfers. (Electronics programmed to <u>ONLY</u> extend footplates when legs are raised.)	N/C	
Fixed Center Mount Legrests				
<input type="checkbox"/>	172H128	Fixed Center Mount Power Recline (Choose individual angle adjustable footplates and custom options below) (F)	N/C	

15. ACCESSORIES				
Wheelchair Tray Table		Positioning Belt - Select for any style seat		
<input type="checkbox"/>	1727A06	Swing-Away Tray 	\$341.00	
<input type="checkbox"/>	172EJ01	Auto Style Buckle	N/C	
<input type="checkbox"/>	172EJ02	Belt 2" Aircraft Buckle	\$168.00	
<input type="checkbox"/>	172EJ03	Belt 2" Aircraft Buckle Padded	\$183.00	
Backpack				
<input type="checkbox"/>	172A1	Black Utility	\$100.00	
Accessories				
<input type="checkbox"/>	1728A09	Ram Phone Mount	\$69.00	
<input type="checkbox"/>	1728A10	Ram Tablet Mount	\$276.00	
<input type="checkbox"/>	1728A11	Ram Cup holder	\$69.00	
<input type="checkbox"/>	1728A12	Switch-it Remote Stop Switch	\$1,236.00	

16. CUSTOM OPTIONS



Looking for more options? Contact our custom "built 4 me" department.

To inquire about any unique customization for Quicke products, contact:

Phone: (800) 263-3390 or Fax: (800) 561-5834

ORDER INFORMATION: Customer Service: 800-263-3390 Fax: 800-561-5834 www.sunrisemedical.ca
Specifications are subject to change without notice. Copyright © 2017 Sunrise Medical Canada Inc.100338 REV 2