



Quickie QM-7 Series SEDEO Ergo

EIPW11, April 2017

Account: _____
 Date: _____
 PO#: _____
 Buyer: _____
 Marked For: _____

Ship To: _____
 Name: _____
 Address: _____
 City, State, PC: _____
 Phone: _____

✓ Option #	Definition	Price	✓ Option #	Definition	Price
------------	------------	-------	------------	------------	-------

1. QUICKIE QM-7 Series Multi-Power Seating (Select One)

QUICKIE QM-710™ MPC for SEDEO Ergo (Group 3 K0861/K0011)

NEW!

<input type="checkbox"/>	Quickie QM-710™ MPC Tilt/Recline For SEDEO Ergo				\$7,825.00
	172PC114	Base Package: Includes SEDEO Ergo seating, LED lights and indicators,			\$7,825.00

QUICKIE QM-715™ HD MPC for SEDEO Ergo (Group 3 K0862/K0011)

NEW!

<input type="checkbox"/>	Quickie QM-715™ HD MPC Tilt/Recline For SEDEO Ergo				\$8,325.00
	172PC115	Base Package: Includes SEDEO Ergo seating, LED lights and indicators for user weight 350lbs,			\$8,325.00

QUICKIE QM-720™ MPD for SEDEO Ergo (Group 4 K0884/K0014)

NEW!

<input type="checkbox"/>	Quickie QM-720™ MPD Tilt/Recline				\$8,570.00
	172PC116	Base Package: Includes SEDEO Ergo seating, LED lights and indicators, 8mph Max speed.			\$8,570.00

MPC Seating Specifications

HGPC Code	K0861 / K0862 / K0884		
Typical Chair Weight	W/battery	375lbs	W/O Battery 285lbs
Transit	ISO Transit		
Seat Type	SEDEO Ergo		
Seat size range	Width	16"-22"	Depth 16"-22"
STFH range	MPC	16.4" - 19.4"	
Power Seating	50° CG tilt and Power recline w/ 7" Mechanical sheer reduction / 12" Seat elevator		



1.1 QM-710™ MPC

QUICKIE QM-7 Series Base Specifications

Drive Wheel Position	Midwheel Drive	Drive Wheels	14" Drive Wheels w/ flat free insert (Black)
User weight Cap	300-350 lb	Tires/Casters	6" Front / 6" rear (Black)
Base size	Width 25" Length 36" Turn Radius 20"	Charger	8A Dual Mode Charger (Compact) CTE
Drive Control	R-Net PM120 (Expandable)	Battery Size	Group 24 / Group 34.
Driving Range	Up to 13.5 miles	Color Accents	Matte Black with 6 color accents
Drive Motors	4-pole 6 mph, Optional 5mph and 8mph	Ground Clearance	4"



Lights and Indicators

2. BASE CHAIR OPTIONS

Speed Package (4-pole)			
<input type="checkbox"/>	172PF85	QM-710 Standard 6 mph motor package	N/C
<input type="checkbox"/>	172PF91	QM-710 5 mph High Torque motor pkg	\$525.00
<input type="checkbox"/>	172PF90	QM-710 8 mph Community motor package	\$1,103.00
<input type="checkbox"/>	172PF84	QM-715 HD Standard 5 mph motor pkg	N/C
<input type="checkbox"/>	172PF89	QM-720 Standard 8 mph motor package	N/C

Suretrac Drive Tracking NEW!			
<input type="checkbox"/>	172MI01	Suretrac Module and Motor Encoders	\$1,475.00

Color Accents NEW!			
<input type="checkbox"/>	172EZ19	Matte Black Shroud w/ White Accents	N/C
<input type="checkbox"/>	172EZ20	Matte Black Shroud w/ Apple Green Acc	N/C
<input type="checkbox"/>	172EZ21	Matte Black Shroud w/ Blue Accents	N/C
<input type="checkbox"/>	172EZ22	Matte Black Shroud w/ Black Accents	N/C
<input type="checkbox"/>	172EZ23	Matte Black Shroud w/ Orange Accents	N/C
<input type="checkbox"/>	172EZ24	Matte Black Shroud w/ Red Accents	N/C

User weight			
<input type="checkbox"/>	172WG06	Less than 150lbs	N/C
<input type="checkbox"/>	172WG13	151-200lb	N/C
<input type="checkbox"/>	172WG14	201-250lb	N/C
<input type="checkbox"/>	172WG03	251-300lb NA with 16" Seat depth and seat elevator	N/C
<input type="checkbox"/>	172WG15	301-350lb	N/C

Drive Wheel			
<input type="checkbox"/>	172DW59	14" Cast Wheel Solid Dual Density Inserts Black Knobby Tires and Black Caster	N/C

Transit Option			
<input type="checkbox"/>	172T025	ISO Transit Package (NA with QM-715)	\$373.00
<input type="checkbox"/>	172T040	Omit Transit	N/C

Wheel locks			
<input type="checkbox"/>	172WL42	Wheel locks	\$142.00

Ventilator Option NEW!			
<input type="checkbox"/>	172VT10	Fixed Vent Tray (Pulmonetic LTV Models)	\$525.00

O2 Bottle Option NEW!			
<input type="checkbox"/>	172A702	O2 Bottle mount to backrest D and E-style	\$340.00

Pre-Set Seat Frame Angle			
<input type="checkbox"/>	172FA10	0°	N/C

Seat Height (Use STFH Matrix below)				
<i>(Measured in the front) - Obtain Quote to verify est. STFH</i>				
<input type="checkbox"/>	172HT44	Position "0"	Requires Group 34 Batteries	N/C
<input type="checkbox"/>	172HT49	+.5" STFH		N/C
<input type="checkbox"/>	172HT45	+1" STFH		N/C
<input type="checkbox"/>	172HT46	+2" STFH		N/C
<input type="checkbox"/>	172HT47	+3" STFH		N/C

LED lights			
<input type="checkbox"/>	172LI06	Include LED lights and indicators NEW!	N/C
<input type="checkbox"/>	172LI02	Omit LED lights and indicators	N/C

	STFH Matrix					
	Seating Type and Power Option Combinations	Seat-to-Floor Height at 0 tilt				
		0	+5"	+1"	+2"	+3"
Ultra-Low STFH Seat Frame "A" <i>(Upcharge see Section 6)</i> Reduces STFH by 1" from standard. Raises recline pivot point by 2" Armrest and backrest also raised 2"	MPC Power Recline and Tilt.		16.5"	17"	18"	19"
	MPC Power Recline and HD Tilt.		16.8"	17.3"	18.3"	19.3"
	MPC Power Recline with Tilt and 12" Lift combo	16.4"	16.9"	17.4"	18.4"	19.4"
Standard STFH Seat Frame "B" Seat pan is in line with recline pivot for standard STFH	MPC Power Recline and Tilt.		*17.5"	18"	19"	
	MPC Power Recline and HD Tilt.		*17.8"	18.3"	19.3"	
	MPC Power Recline with Tilt and 12" Lift combo	17.4"	17.9"	18.4"	19.4"	
* Seat Depth Limited to 19" with the STFH configuration STFH measured with batteries but no user. Actual STFH may vary by .5"						

3. BATTERIES and CHARGERS

Batteries			
<input type="checkbox"/>	172AG04	GP24 MK Gel Flat Top Batts Pair w/strap Recommended for Max Daily Range	\$830.00
<input type="checkbox"/>	172AG22	GP24 Brightway AGM Flat Top Batteries Pair w/strap.	\$795.00
<input type="checkbox"/>	172AG11	GP34 MK Gel Flat Top Batteries (Pair) w/straps, Required for lower STFH.	\$830.00
<input type="checkbox"/>	172AG23	GP34 Brightway AGM Flat Top Batts (Pair) w/straps, Required for lower STFH.	\$750.00

Battery Charger			
<input type="checkbox"/>	172BJ19	8A Dual Mode (Compact) CTE	N/C



ORDERING INFORMATION: PHONE: 1.800.263.3390 FAX: 1.800.561.5834 WEBSITE: www.sunrisemedical.ca

Specifications are subject to change without notice. Copyright © 2016 100337 REV 2

4. ELECTRONICS

Expandable Controllers

<input type="checkbox"/>	172ER08	R-Net Expandable Drive Control (120A)	\$1,066.0
	172EX08	Expandable 120A Controller - R-net Includes Memory Seating Positions	N/C
	172EX07	Harness for 120A Expandable Controller	N/C

Memory Seating

NEW!

"Memory Seating" feature allows for 6 custom seating positions that are accessed either thru the joystick, Ctrl+5, Joystick Shortcut keys or programmable jacks on OMNI and joystick. When enabled, chair will be pre-programmed with:

Seated, Supine, Pressure Relief, Feet Above Heart Transfer and C-ME.
Settings are dealer settable using mono switch (not included).

<input type="checkbox"/>	172OG01	Include Memory Seating	N/C
<input type="checkbox"/>	172OG02	Omit Memory Seating	N/C

R-net Joystick

<input type="checkbox"/>	172J113	CJSM: Joystick w/Stereo Jacks & Color Display with Light and indicator controls.	N/C
<input type="checkbox"/>	172J112	CJSM2: Joystick Advanced: 3.5" Color Display, Paddle Switches, Programmable Shortcut Keys, ICONs, Stereo Jack and IR control	\$1,133.00
<input type="checkbox"/>	172J114	CJSM2: Joystick Advanced: 3.5" Color Display, Paddle Switches, Programmable Shortcut Keys, ICONs, Stereo Jack, IR and dual (IOS/Android) Bluetooth control	\$1,133.00
<input type="checkbox"/>	172J95	Omit JS. For use with OMNI	N/C

Joystick Mount (Fixed, Retractable, Mid)

<input type="checkbox"/>	172JM16	Fixed Mount for Ergo	N/C
<input type="checkbox"/>	172JM17	Retractable Joystick mount for Ergo	\$460.00
<input type="checkbox"/>	172JM14	Mid Line Mount	\$388.00

Joystick Mount (Location)

Height Adjustable.

<input type="checkbox"/>	172HM1	Right Side Mount	N/C
<input type="checkbox"/>	172HM2	Left Side Mount	N/C

Joystick Handle (Optional)

<input type="checkbox"/>	172JH1	Ball Handle	\$95.00
<input type="checkbox"/>	172JH2	Foam Ball Handle	\$95.00
<input type="checkbox"/>	172JH4	T-Handle	\$95.00
<input type="checkbox"/>	172JH6	Goalpost	\$95.00

Electronic Programmers

<input type="checkbox"/>	172PR5	R-net PC Programmer with Dongle	\$615.00
--------------------------	--------	---------------------------------	----------

R-net OMNI and Output Options

(For expandable controls, choose all that apply)

<input type="checkbox"/>	172ON6	Rnet OMNI Input Control Module and display with built in Infrared and armrest mount. Choose mount below. (required for Specialty Input Devices.)	\$2,075.00
<input type="checkbox"/>	172ON7	Rnet Output Module (ECM)	\$1,066.00
<input type="checkbox"/>	172ON8	Rnet Bluetooth Mouse-mover for Windows / Android	\$1,066.00
<input type="checkbox"/>	172ON9	Rnet Bluetooth for Apple Devices	\$1,066.00

Attendant Control

<input type="checkbox"/>	172OA04	Attendant Control for R-net	\$1,066.00
--------------------------	---------	-----------------------------	------------

Attendant Mounting location

<input type="checkbox"/>	172DN1	Mounting location: Right Back Cane	N/C
<input type="checkbox"/>	172DN2	Mounting Location: Left Back Cane	N/C

OMNI Mount

(If used with Joystick must be mounted opposite side)

<input type="checkbox"/>	172DS04	Swing Away Left Side OMNI Mount	N/C
<input type="checkbox"/>	172DS03	Swing Away Right Side OMNI Mount	N/C

Switch inputs for Joysticks or OMNI

Switches can be used with "Assignable Button" feature to control Power, Horn, Mode, Speeds, Profile, Seating, IR, Bluetooth, Lighting, Output Module, Memory Settings and Pressure Relief Reminder Mute.

<input type="checkbox"/>	172SW22	Stereo to 2 Mono Splitter	\$75.00
<input type="checkbox"/>	172SW6	Ribbon Switch	\$119.00
<input type="checkbox"/>	172SW1	Buddy Button Switch	\$95.00
<input type="checkbox"/>	172SW2	Disc Switch	\$202.00
<input type="checkbox"/>	172SW9	Wobble Switch	\$630.00
<input type="checkbox"/>	172SW4	Micro Light Switch	\$95.00
<input type="checkbox"/>	172SW24	Egg Switch-Black	\$119.00
<input type="checkbox"/>	172SW21	Dual toggle switch and bracket	\$214.00
<input type="checkbox"/>	172SW35	Switch-it Mini 1 button on Mono plug	\$268.00
<input type="checkbox"/>	172SW36	Switch-it Mini 2 button on 2 Mono plugs	\$268.00
<input type="checkbox"/>	172SW37	Switch-it Mini 2 button on 1 Stereo plug	\$268.00
<input type="checkbox"/>	172SW38	Switch-it Mini 4 button on DB9 connector	\$268.00

Switch Mount Bracket

<input type="checkbox"/>	172EM01	Single Switch Mount Bracket	\$247.00
--------------------------	---------	-----------------------------	----------

Default Memory Positions



M1: Seated



M2: Supine



M3: Pressure Relief



M4: Feet above Heart



M5: Transfer



M6: C-ME

5. Specialty Driver Control Options

Requires expandable Rnet controls and OMNI unless otherwise noted.

5a. Proportional Control Packages

Head Control Proportional Packages			
□ Switch-it Dual Pro Proportional Head Control with Whitmyer and Facial Pads (Adult)			\$6,541.00
172ID72		Switch-it Dual Pro Head Control w/facial pad (Adult)	\$6,139.00
172IJ09		SA Boom Mount with Egg Switch Control	\$402.00
□ Switch-it Dual Pro Proportional Head Control with Spot Pads Package			\$6,541.00
172ID70		Switch-it Dual Pro Head Control w/Spot pads	\$6,139.00
172IJ09		SA Boom Mount with Egg Switch Control	\$402.00

Hand Control Proportional Joystick Packages			
□ Switch-it Versa-Guide Compact Joystick Midline with Gatlin Package			\$2,369.00
172ID71		Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
172IJ07		S/A Gatlin Mount	\$330.00
□ Switch-it Versa-Guide Compact Joystick Midline and prox flex switches with Gatlin Package			\$3,018.00
172ID71		Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
172IM08		Dual Flex Proximity Switches (requires no force)	\$649.00
172IJ07		S/A Gatlin Mount	\$330.00
□ Switch-it Micro-Pilot Joystick with Space Disc Package			\$4,007.00
172ID62		Micro-Pilot Proportional Joystick	\$3,534.00
172IJ08		S/A SpaceDisc mount	\$473.00
□ Switch-it Micro-Guide Joystick with SpaceDisc Package			\$4,831.00
172ID63		Micro-Guide Proportional Joystick	\$4,358.00
172IJ08		S/A SpaceDisc mount	\$473.00
□ ASL Proportional Micro Extremity Joystick (ASL130) Midline with Gatlin Package			\$4,414.00
172ID59		Micro Extremity Joystick	\$4,084.00
172IJ07		S/A Gatlin Mount	\$335.00
□ ASL Mushroom Joystick (ALS135) Midline with Gatlin Package			\$3,353.00
172ID60		Mushroom Joystick	\$3,023.00
172IJ07		S/A Gatlin Mount	\$330.00
□ ASL Compact Joystick w Integrated buttons (ASL134) Midline with Gatlin Package			\$1,530.00
172ID61		ASL Compact Joystick	\$1,200.00
172IJ07		S/A Gatlin Mount	\$330.00
□ ASL Micro Mini Joystick w Hand Pad (ASL136) Midline with Gatlin Package			\$4,012.00
172ID68		Micro Mini Joystick	\$3,682.00
172IJ07		S/A Gatlin Mount	\$330.00



Dual Pro Head Control

Chin Control Proportional Joystick Packages			
□ Switch-it Versa-Guide Chin Control Compact Joystick and Egg Switch Package			\$2,406.00
172ID71		Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
172IM07		S/A Mount off backrest	\$252.00
172IJ04		SA Boom Mount with egg switch	\$115.00
□ Switch-it Versa-Guide Chin Control Compact Joystick and Proximity Flex Switches Package			\$2,940.00
172ID71		Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
172IM08		Dual Flex Proximity Switches (requires no force)	\$649.00
172IM07		S/A Mount off backrest	\$252.00
□ Switch-It Versa-Guide Bib Chin Control Compact Joystick with Proximity Flex switches Package Std 6"			\$3,445.00
172ID71		Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
172IM08		Dual Flex Proximity Switches (requires no force)	\$649.00
172IM05		Bib mount Small	\$757.00
□ Switch-It Versa-Guide Bib Chin Control Compact Joystick with Proximity Flex Switches Package Large 8"			\$3,512.00
172ID71		Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
172IM08		Dual Flex Proximity Switches (requires no force)	\$649.00
172IM06		Bib mount Large	\$824.00
□ Switch-It Versa-Guide Bib Chin Control Compact Joystick with Egg Switch Package Std 6"			\$2,911.00
172ID71		Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
172IM05		Bib mount Small	\$757.00
172IJ04		SA Boom Mount with egg switch	\$115.00
□ Switch-It Versa-Guide Bib Chin Control Compact Joystick with Egg Switch Package Large 8"			\$2,978.00
172ID71		Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
172IM06		Bib mount Large	\$824.00
172IJ04		SA Boom Mount with egg switch	\$115.00
□ Switch-It Micro-Pilot Bib Chin Control with Egg Switch Package Std 6"			\$4,406.00
172ID62		Micro-Pilot Proportional Joystick	\$3,534.00
172IM05		Bib mount Small	\$757.00
172IJ04		SA Boom Mount with egg switch	\$115.00
□ Switch-It Micro-Pilot Bib Chin Control with Egg Switch Package Large 8"			\$4,473.00
172ID62		Micro-Pilot Proportional Joystick	\$3,534.00
172IM06		Bib mount Large	\$824.00
172IJ04		SA Boom Mount with egg switch	\$115.00
□ Switch-It Micro-Guide Bib Chin Control with Egg Switch Package Std 6"			\$5,230.00
172ID63		Micro-Guide Proportional Joystick	\$4,358.00
172IM05		Bib mount Small	\$757.00
172IJ04		SA Boom Mount with egg switch	\$115.00
□ Switch-It Micro-Guide Bib Chin Control with Egg Switch Package Large 8"			\$5,297.00
172ID63		Micro-Guide Proportional Joystick	\$4,358.00
172IM06		Bib mount Large	\$824.00
172IJ04		SA Boom Mount with egg switch	\$115.00
□ ASL Proportional Micro Extremity Joystick Chin / Lip with Egg Switch Control Package			\$4,498.00
172ID59		Micro Extremity Joystick	\$4,084.00
172IM07		S/A Mount off backrest	\$299.00
172IJ04		SA Boom Mount with egg switch	\$115.00

5. Specialty Driver Control Options (cont.)

5b. Switch Control

Requires Headrest and/or Headrest Mount

3 Switch Proximity Head Arrays Choose headrest mount below

<input type="checkbox"/> ASL Proximity Head Array Package (Adult) \$5,755.00		
172II02	ASL Proximity Head Array (Adult)	\$5,353.00
172IJ09	SA Boom Mount with Egg Switch Control	\$402.00
<input type="checkbox"/> ASL Proximity Head Array Package (Ped) \$5,755.00		
172II03	ASL Proximity Head Array (Ped)	\$5,353.00
172IJ09	SA Boom Mount with Egg Switch Control	\$402.00
<input type="checkbox"/> Switch-It Head Array Package 3 Switch Curved (Long) \$5,755.00		
172II13	Switch-it Proximity Head Array Long	\$5,353.00
172IJ09	SA Boom Mount with Egg Switch Control	\$402.00
<input type="checkbox"/> Switch-It Head Array Package 3 switch curved (Stubby) \$5,691.00		
172II14	Switch-it Proximity Head Array Stubby	\$5,289.00
172IJ09	SA Boom Mount with Egg Switch Control	\$402.00

5 Switch Proximity Head Array (No Egg Switch Needed) Choose "Headrest Mount for Special Controls"

<input type="checkbox"/> 172II04	Switch-It Head Array 5 switch flat	\$4,295.00
----------------------------------	------------------------------------	------------

Multiple Switches.


<input type="checkbox"/> 172SW31	Switch-it 3 Naked proximity switches on DB9 connector. (requires OMNI)	\$3,039.00
<input type="checkbox"/> 172SW32	Switch-it 4 Naked proximity switches on DB9 connector. (requires OMNI)	\$3,718.00
<input type="checkbox"/> 172SW33	Switch-it 5 Naked proximity switches on DB9 connector. (requires OMNI)	\$4,146.00
<input type="checkbox"/> 172SW34	ASL (502) 5-Switch adaptor box on DB9 connector. (for use with 5 mono switches, not included)	\$4,146.00
<input type="checkbox"/> 172SW23	4-way toggle switch button DB9 Connector. Includes bracket	\$670.00

Sip N Puff Package Choose headrest and HR mount below

<input type="checkbox"/> Quickie Sip n' Puff System \$2,635.00		
172ID34	Sip N' Puff Kit	\$2,063.00
172IH1	Sip n' Puff Tubing and SA Hardware	\$457.00
172IJ10	SA Boom Mount with egg switch	\$115.00

Headrest Mount for Specialty Controls

NEW!

<input type="checkbox"/> 172MH23	Whitmyer Linx Headrest Mount, Quick Release	 ON-CHAIR	\$485.00
----------------------------------	---	---	----------

Whitmyer Headrest Pad for Specialty controls

<input type="checkbox"/> 172HE79	Plush Pad 8" for Specialty controls	\$309.00
<input type="checkbox"/> 172HE80	Plush Pad 10" for Specialty controls	\$309.00
<input type="checkbox"/> 172HE81	Plush Pad 14" for Specialty controls	\$309.00
<input type="checkbox"/> 172HE90	Adjusta-A-Plush Headrest Pad Narrow	\$309.00
<input type="checkbox"/> 172HE91	Adjusta-A-Plush Headrest Pad Standard	\$309.00

Extra Swing Away Head Boom Egg Switch and Mount

172IJ05	Extra SA Boom Mount with egg switch #1	\$115.00
172IJ06	Extra SA Boom Mount with egg switch #2	\$115.00

Single Switches



<input type="checkbox"/> 172IB01	Ribbon Switch	\$119.00
<input type="checkbox"/> 172IB02	Buddy Button Switch	\$108.00
<input type="checkbox"/> 172IB03	Disc Switch	\$202.00
<input type="checkbox"/> 172IB04	Wobble Switch	\$630.00
<input type="checkbox"/> 172IB05	Micro Light Switch	\$95.00
<input type="checkbox"/> 172IB06	Egg Switch-Black	\$134.00
<input type="checkbox"/> 172IB07	Switch-it Mini 1 button on Mono plug	\$268.00

6. SEAT TYPE

<input type="checkbox"/>	172ST63	SEDEO Ergo Recline	N/C
--------------------------	---------	--------------------	-----

<input type="checkbox"/>	172FU30	Ultra-Low STFH Seat Frame "A"	\$350.00
<input type="checkbox"/>	172FU31	Standard Seat Frame "B"	N/C

7. SEAT SIZE



Critical Measurements

A: Seat Pan Width (Section 7):

Measured on the seat pan from inside track rail to inside track rail, Infinite adjustable from 16"-22" wide.

B: Seat Cushion Width (Section 8):

When ordered with ERGO thigh supports, the seat cushion will be 2" narrower than selected seat width as ERGO thigh supports consume 2" of effective seat width. JAY Bundled cushion will be shipped actual selected size so be sure to order appropriately if combined with ERGO thigh supports.

C: ERGO Thigh Support (Bolster - Section 9) Width:

ERGO thigh supports measure 3.5" wide each, plus they sit on top of the outer track rails and consume 2" of effective seat width. ERGO thigh supports and track rail unlatch, swing away and lift off for lateral transfers.

D: Seat Depth (Section 7):

Measured from back shell surface to front of seat pan. Most backs cushions consume 1" of effective seat depth. Seat Pan has Infinite adjustable from 16"-22" deep.

Seat Width (Pan Width)

<input type="checkbox"/>	172R5	16"	Seat width	(N/C)
<input type="checkbox"/>	172R6	17"	Seat width	(N/C)
<input type="checkbox"/>	172R7	18"	Seat width	(N/C)
<input type="checkbox"/>	172R8	19"	Seat width	(N/C)
<input type="checkbox"/>	172R9	20"	Seat width	(N/C)
<input type="checkbox"/>	172R52	21"	Seat width	(N/C)
<input type="checkbox"/>	172R25	22"	Seat width	(N/C)

Seat Depth (Section 7)

<input type="checkbox"/>	172RF5	16"	Seat depth	NA above 250lbs	(N/C)
<input type="checkbox"/>	172RF6	17"	Seat depth		(N/C)
<input type="checkbox"/>	172RF7	18"	Seat depth		(N/C)
<input type="checkbox"/>	172RF8	19"	Seat depth		(N/C)
<input type="checkbox"/>	172RF24	20"	Seat depth	see STFH Matrix	(N/C)
<input type="checkbox"/>	172RF26	21"	Seat depth		(N/C)
<input type="checkbox"/>	172RF25	22"	Seat depth		(N/C)

8. SEAT CUSHIONS

JAY Ion Cushion for SEDEO ERGO

NEW!

<input type="checkbox"/>	172CU45	JAY ION cushion with 3DX cover	\$445.00
<input type="checkbox"/>	172CU46	JAY ION cushion with 3DX and Black Stretch	\$445.00
<input type="checkbox"/>	172CU47	JAY ION cushion w/ Premium leatherette cover	\$445.00

JAY Cushion Bundling

<input type="checkbox"/>	172JC14	JAY Cushion <u>On Chair Order</u>	NC
--------------------------	---------	-----------------------------------	----

Note: see measurement "B" and "C" above. When ordered with ERGO thigh guide, order bundled cushion 2" narrower than effective seat width.

JAY X2 Cushion

<input type="checkbox"/>	172JU22	JX21616	JAY X2 16" x 16"	\$830.00
<input type="checkbox"/>	172JU23	JX21618	JAY X2 16" x 18"	\$830.00
<input type="checkbox"/>	172JU24	JX21816	JAY X2 18" x 16"	\$830.00
<input type="checkbox"/>	172JU25	JX21818	JAY X2 18" x 18"	\$830.00
<input type="checkbox"/>	172JU26	JX21820	JAY X2 18" x 20"	\$830.00
<input type="checkbox"/>	172JU27	JX22020	JAY X2 20" x 20"	\$830.00

JAY Fusion Cushion

<input type="checkbox"/>	172JU01	JFUSION1616	JAY Fusion 16" x 16"	\$920.00
<input type="checkbox"/>	172JU02	JFUSION1618	JAY Fusion 16" x 18"	\$920.00
<input type="checkbox"/>	172JU03	JFUSION1816	JAY Fusion 18" x 16"	\$920.00
<input type="checkbox"/>	172JU04	JFUSION1818	JAY Fusion 18" x 18"	\$920.00
<input type="checkbox"/>	172JU05	JFUSION1820	JAY Fusion 18" x 20"	\$920.00
<input type="checkbox"/>	172JU06	JFUSION2018	JAY Fusion 20" x 18"	\$920.00
<input type="checkbox"/>	172JU07	JFUSION2020	JAY Fusion 20" x 20"	\$920.00

JAY Union Cushion

<input type="checkbox"/>	172JU15	JU1616	JAY UNION 16" x 16"	\$673.00
<input type="checkbox"/>	172JU16	JU1618	JAY UNION 16" x 18"	\$673.00
<input type="checkbox"/>	172JU17	JU1816	JAY UNION 18" x 16"	\$673.00
<input type="checkbox"/>	172JU18	JU1818	JAY UNION 18" x 18"	\$673.00
<input type="checkbox"/>	172JU19	JU1820	JAY UNION 18" x 20"	\$673.00
<input type="checkbox"/>	172JU20	JU2018	JAY UNION 20" x 18"	\$673.00
<input type="checkbox"/>	172JU21	JU2020	JAY UNION 20" x 20"	\$673.00

JAY Custom Cushion

To order any **JAY cushions not listed above**, please see the JAY price list. Cushion name or part number must be filled out w/ desired width and depth

<input type="checkbox"/>	172JC20	JAY Bundled Cushion	Cushion Name or Part #	Width:	Depth:	Price List
--------------------------	---------	---------------------	------------------------	--------	--------	------------

9. POSTURAL SUPPORTS

LATERAL THORACIC SUPPORT (Swing Away)

Available on all Solid Backs.

<input type="checkbox"/> LAT 4X4 S/A MOUNTS W/ 1" OFFSET - PR \$805.00		
172ZV14	LATERAL 4X4 W/ 1" OFFSET - PAIR	\$260.00
172ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$545.00
<input type="checkbox"/> LAT 4X5 S/A MOUNTS W/ 1" OFFSET - PR \$805.00		
172ZV15	LATERAL 4X5 W/ 1" OFFSET - PAIR	\$260.00
172ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$545.00
<input type="checkbox"/> LAT 4X6 S/A MOUNTS W/ 1" OFFSET - PR \$805.00		
172ZV16	LATERAL 4X6 W/ 1" OFFSET - PAIR	\$260.00
172ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$545.00
<input type="checkbox"/> LAT 6X4 S/A MOUNTS W/ 1" OFFSET - PR \$805.00		
172ZV17	LATERAL 6X4 W/ 1" OFFSET - PAIR	\$260.00
172ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$545.00
<input type="checkbox"/> LAT 6X5 S/A MOUNTS W/ 1" OFFSET - PR \$805.00		
172ZV18	LATERAL 6X5 W/ 1" OFFSET - PAIR	\$260.00
172ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$545.00
<input type="checkbox"/> LAT 6X6 S/A MOUNTS W/ 1" OFFSET - PR \$805.00		
172ZV19	LATERAL 6X6 W/ 1" OFFSET - PAIR	\$260.00
172ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$545.00

CHEST STRAP

<input type="checkbox"/> 172ZR03	Padded Chest Strap	\$230.00
----------------------------------	--------------------	----------

Ergo Thigh Supports

(Upholstered in 3DX unless back is ordered in Leatherette).

<input type="checkbox"/> Ergo Thigh Support - Left (Seat Bolster) NEW! \$335.00		
172TV5	Ergo Long Thigh Support - Left	\$105.00
172TX3	Swing-away Mount - Left	\$230.00
<input type="checkbox"/> Ergo Thigh Support - Right (Seat Bolster) NEW! \$335.00		
172TV6	Ergo Long Thigh Support - Right	\$105.00
172TX4	Swing-away Mount - Right	\$230.00

Long Thigh Support w/Removable Mount

<input type="checkbox"/> SPOT Long Thigh Support - Left \$409.00		
172TV2	Long Thigh Support - Left	\$126.00
172TX1	Removable Mount - Left	\$283.00
<input type="checkbox"/> SPOT Long Thigh Support - Right \$409.00		
172TV1	Long Thigh Support - Right	\$126.00
172TX2	Removable Mount - Right	\$283.00

Swing Away Knee Adductor Supports

Knee Adductor Supports - LEFT (choose Pad and Arm) **\$489.00**

Pad Size		
<input type="checkbox"/> 172KB5	Knee Adductor Pad SM- Left	\$215.00
<input type="checkbox"/> 172KB7	Knee Adductor Pad LG- Left	\$215.00
S/A Arm Size		
<input type="checkbox"/> 172XX04	Swing Away Short - Left	\$283.00
<input type="checkbox"/> 172XX06	Swing Away Long - Left	\$283.00

Knee Adductor Supports - RIGHT (choose Pad and Arm) **\$489.00**

Pad Size		
<input type="checkbox"/> 172KB6	Knee Adductor Pad SM- Right	\$215.00
<input type="checkbox"/> 172KB8	Knee Adductor Pad LG- Right	\$215.00
S/A Arm Size		
<input type="checkbox"/> 172XX05	Swing Away Short - Right	\$283.00
<input type="checkbox"/> 172XX07	Swing Away Long - Right	\$283.00

10. POWERED SEATING

Power Seating Options **NEW!**

<input type="checkbox"/> 172PA28	50 Degree Tilt and Ergo Recline w/Shear Reduction and optional Biometric Repositioning (2 actuators)	\$12,450.00
<input type="checkbox"/> 172PA29	QM-715 HD 50 Degree Tilt and Ergo Recline w/Shear Reduction and optional Biometric Repositioning (350 weight limit) (2 actuators)	\$12,450.00

Biometric Repositioning Recline **NEW!**

<input type="checkbox"/> 172OG03	Include Biometric Reposition recline. (Requires Memory positions.)	N/C
<input type="checkbox"/> 172OG04	Omit Biometric Reposition recline.	N/C


Powered Seat Elevator **NEW!**

(NA with QM-715 HD)

<input type="checkbox"/> 172PA30	12" Power Seat Elevate (Optional C-ME)	\$3,025.00
----------------------------------	---	------------

C-ME (3mph drive speed with full elevate)

<input type="checkbox"/> 172OG05	Include C-ME with Seat elevator (NA with QM-720 or QM-715HD)	N/C
<input type="checkbox"/> 172OG06	Omit C-ME with Seat elevator	N/C



Biometric Repositioning

Power Seating Control

Powered Seating Thru-drive Control Method.

Thru-Drive 2+ Actuator

<input type="checkbox"/>	172EO31	Thru Drive 2+ Actuators: (For operating Tilt, recline, lift, power legs) Includes Ctrl+5 and optional PRAR	\$2,525.00
--------------------------	---------	--	------------

Powered Seating Thru-drive Control Options

Assignable Buttons

"Assignable Button" feature allows for customization of joystick buttons and/or Ctrl+5 Switches allowing access to various wheelchair and seating functions. When enabled, chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle.

<input type="checkbox"/>	172PM01	Enable Assignable buttons	N/C
<input type="checkbox"/>	172PM02	Disable Assignable buttons	N/C

PRAR Pressure Relief Alert and Recording

NEW!

PRAR Pressure Relief Alert and Record feature provides programmable alerts to be set into the Rnet electronics that will prompt the user when it is time to perform a pressure relief. The net system will record how often the pressure relief is completed or ignored. Requires Memory Seating function.

<input type="checkbox"/>	172OG07	Include PRAR	N/C
<input type="checkbox"/>	172OG08	Omit PRAR	N/C



Icon on Ctrl+5 flashes Blue when it is time to do a pressure relief and / or

Joystick emits an Alert Tune



Ctrl+5 Switchbox Options

Ctrl+5 allows direct control to 5 power seating functions (10 total switches). Switches can also be used with assignable buttons and Memory seating for greater flexibility. Box includes 3-Amp USB port for charging portable devices

#1 Ctrl+5 Switch Location

<input type="checkbox"/>	172QG01	Ctrl+5 Switchbox mounted behind joystick	N/C
<input type="checkbox"/>	172QG02	Ctrl+5 Switchbox mounted opposite joystick	N/C
<input type="checkbox"/>	172QG03	Ctrl+5 Switchbox mounted on back cane	N/C

#1 Ctrl+5 Switch Type

<input type="checkbox"/>	172EY03	Push buttons	N/C
<input type="checkbox"/>	172EY04	Toggle switches	\$250.00

A 2nd Ctrl+5 can be added to provide additional control for Memory seating or environmental controls. It can also be mounted as an attendant control all functions of #1 are available.

#2 Ctrl+5 Switch Location

NEW!

<input type="checkbox"/>	172QG06	Ctrl+5 Switchbox mounted behind joystick	\$1,255.00
<input type="checkbox"/>	172QG07	Ctrl+5 Switchbox mounted opposite joystick	\$1,255.00
<input type="checkbox"/>	172QG08	Ctrl+5 Switchbox mounted on back cane	\$1,255.00

#2 Ctrl+5 Switch Type

<input type="checkbox"/>	172EY05	Push buttons	N/C
<input type="checkbox"/>	172EY06	Toggle switches	\$250.00



Ctrl+5 Push Button option



Ctrl+5 Toggle Switch option

11. BACKREST AND BACK OPTIONS FOR ERGO RECLINE

Backrest Frame Type			
<input type="checkbox"/>	172BT89	Power Recline (Power Recline Seat only)	N/C

Back Upholstery Type			
<input type="checkbox"/>	172BU47	Ergo Back Cushion: choose options below	N/C
<input type="checkbox"/>	172JB03	Jay J3 for Power Recline w/1.5" Foam Choose type below.	N/C
<input type="checkbox"/>	172BU14	Omit Back. Includes bracket for aftermarket backs	N/C

Push Handle Options			
<input type="checkbox"/>	172BP55	Fixed "Cane" Handles. NA with O2 holder	N/C
<input type="checkbox"/>	172BP56	Height Adjustable Push Handles. w/2 utility hooks	\$53.00
<input type="checkbox"/>	172BP57	Omit push handles	N/C

Power Recline Backrest Preset Home Angle N/C					
<input type="checkbox"/>	90° 172QK1	<input type="checkbox"/>	95° 172QK2 default	<input type="checkbox"/>	100° 172QK3

11a. SEDEO ERGO BACK CUSHION AND OPTIONS



SEDEO Ergo Back Type SEDEO Ergo Back coding: (xx=height, YY=width, zz=Upholstery) Back cushion reduces effective seat depth by 1"			
<input type="checkbox"/>	172OD21	Ergo Back Cushion ERGOxxYYzz : For 16" to 20" seat width	\$895.00
<input type="checkbox"/>	172OD22	Ergo Back Cushion: ERGOxxYYzz For 22" seat width	\$895.00

SEDEO Ergo Back Height (xx) (measured from seat pan to top of back shell)				
<input type="checkbox"/>	172B152	22"	NA with Ultra Low STFH Seat Frame	N/C
<input type="checkbox"/>	172B153	24"		N/C
<input type="checkbox"/>	172B154	26"		N/C
<input type="checkbox"/>	172B158	28"	NA with Standard Seat Frame	N/C

SEDEO Ergo Back Width (YY)				
<i>Back can be ordered a different width than seat pan. Armrest will be adjusted to back width ordered.</i>				
<input type="checkbox"/>	172BN18	16"	Back width	N/C
<input type="checkbox"/>	172BN19	18"	Back width	N/C
<input type="checkbox"/>	172BN20	20"	Back width	N/C
<input type="checkbox"/>	172BN21	22"	Back width	N/C

SEDEO Ergo Back Upholstery options (zz)				
<input type="checkbox"/>	172ZD76	sf	3DX™ spacer fabric (sf)	N/C
<input type="checkbox"/>	172ZD77	xs	3DX™ spacer fabric and Lycra® stretch outer layer (xs)	N/C
<input type="checkbox"/>	172ZD78	pl	Premium Leatherette with perforated inserts (pl)	N/C



3DX Option



Premium Leatherette Option

11b. JAY J3 BACK for ERGO POWER RECLINE OPTIONS

JAY J3 Back Contour. Back cushion reduces effective seat depth by 1"			
<input type="checkbox"/>	172JP01	J3 PA Back (2.2" of Contour Depth) Available 16"-20" Wide (J3PAUTZYXX)	\$803.00
<input type="checkbox"/>	172JP03	J3 PL Back (3.0" of Contour Depth) Available 16"-20" Wide (J3PLUTZYXX)	\$867.00
<input type="checkbox"/>	172JP04	J3 PDL Back (6.0" of Contour Depth) Available 16"-20" Wide (J3PDLUTZYXX)	\$930.00
<input type="checkbox"/>	172JP05	J3 PDC Back (6.0" of Contour Depth) Available 16"-20" Wide (J3PDCUTZYXX)	\$940.00
<input type="checkbox"/>	172JP06	J3 PD Back (6.0" of Contour Depth) Available 16"-20" Wide (J3PDUTZYXX)	\$940.00


JAY J3 Back Height Setting. Backs are height adjustable <i>Measured from pan to top of solid back.</i>			
<input type="checkbox"/>	172B42	22"	N/C
<input type="checkbox"/>	172B44	24"	N/C
<input type="checkbox"/>	172B140	26"	N/C
<input type="checkbox"/>	172B141	27"	N/C
<input type="checkbox"/>	172B142	28"	N/C

SPINE-ALIGN MINI FIT KIT 2 <i>Includes: 5 wedges, 2 lat shims, 2 post/lat pelvic, 1 lg lumbar, 1 flat shim, and 1 curved</i>				
<input type="checkbox"/>	172JS01	K0108	SPINE-ALIGN MINI FIT KIT 2	\$101.00

JAY J3 Back Shell Options and Height Ranges. Jay Shells are infinitely height adjustable							
	Choose Frame Type in Section 6	Frame A		Frame B			
		min	max	min	max		
<input type="checkbox"/>	172JR06	Mid Thoracic Short 15" shell	21"	24"	19"	22"	N/C
<input type="checkbox"/>	172JR01	Upper Thoracic Short 16.5" shell	22"	25"	20"	23"	N/C
<input type="checkbox"/>	172JR02	Upper Thoracic Med 18" shell	23"	26"	21"	24"	N/C
<input type="checkbox"/>	172JR03	Upper Thoracic Tall 19.5" shell	24"	27"	22"	25"	N/C
<input type="checkbox"/>	172JR04	Shoulder Height 21" shell	25"	28"	23"	26"	N/C

JAY J3 Back Width <i>Back can be ordered a different width than seat pan. Armrest will be adjusted to back width ordered.</i>			
<input type="checkbox"/>	172JQ01	Back width 14"	N/C
<input type="checkbox"/>	172JQ02	Back width 16"	N/C
<input type="checkbox"/>	172JQ03	Back width 18"	N/C
<input type="checkbox"/>	172JQ04	Back width 20"	N/C

12. HEADREST OPTIONS

Headrest Mount			
<input type="checkbox"/>	172MH18	Ergo Multi-position headrest mount, Quick release	\$300.00
<input type="checkbox"/>	172MH22	Whimyer Linx Headrest Mount, Quick Release 	\$485.00

Whitmyer Headrest Pad			
<input type="checkbox"/>	172HE65	Whitmyer Plush Pad 8"	\$237.00
<input type="checkbox"/>	172HE66	Whitmyer Plush Pad 10"	\$237.00
<input type="checkbox"/>	172HE78	Whitmyer Plush Pad 14"	\$237.00
<input type="checkbox"/>	172HE82	Whitmyer Adjusta-A-Plush Pad Narrow	\$365.00
<input type="checkbox"/>	172HE83	Whitmyer Adjusta-A-Plush Pad Standard	\$365.00

13. ARMREST AND ARMPAD OPTIONS

Armrest for Power Recline Seating. <i>(Armrest width is adjustable from 16" to 22")</i>			
<input type="checkbox"/>	172EG29	NEW! Cantilever / Reclining Arm Standard STFH 10"-14" Arm Height Range Ultra-low STFH 12"-16" Arm Height Range	N/C


Reclining Arm pad Options <i>(Arm pad can be angled inward 15 degrees)</i>			
Standard Waterfall Reclining Arm pad (Upholstered in 3DX unless back is ordered in Leatherette).			
Left Arm Pad (N/C)		Right Arm pad (N/C)	
<input type="checkbox"/>	172AE29	St Reclining Arm Pad (16")	<input type="checkbox"/>
Comfort Reclining Arm pad (Foam)			
Left Arm pad (\$98.00)		Right Arm pad (98.00)	
<input type="checkbox"/>	172AE22	Comfort 2" wide Reclining Arm Pad	<input type="checkbox"/>
<input type="checkbox"/>	172AE23	Comfort 4" wide Reclining Arm Pad	<input type="checkbox"/>
Comfort Reclining Arm pad (Gel)			
Left Gel Arm pad (\$185.00)		Right Gel Arm pad (185.00)	
<input type="checkbox"/>	172AE25	Comfort 2" wide Reclining Arm Pad	<input type="checkbox"/>
<input type="checkbox"/>	172AE26	Comfort 4" wide Reclining Arm Pad	<input type="checkbox"/>
AES Trough/Elbow Block Reclining Arm pads (Gel)			
Left Gel Arm pad (\$258.00)		Right Gel Arm pad (\$258.00)	
<input type="checkbox"/>	172AE27	Trough / Elbow Block Medium (Gel) (3.5" x 13") Reclining Arm Pad	<input type="checkbox"/>
<input type="checkbox"/>	172AE28	Trough / Elbow Block Large (Gel) (4" x 14.5") Reclining Arm Pad	<input type="checkbox"/>

14. Centermount Footrest options

Power Elevating Legrests			
<input type="checkbox"/>	172H127	Power Centermount ELR w/7.5" Power Articulation (F) Centermount footplates plus Custom Options Below	\$2,987.00
Power Centermount Footplate Extension Option (For Thru-Drive Operation)			
<input type="checkbox"/>	172IE01	Enable 7.5" Powered footplate extension for transfer assist. (Electronics programmed to allow footplates to be extended with or without legrest being raised)	\$515.00
<input type="checkbox"/>	172IE02	Disable Footplate extension for transfers. (Electronics programmed to ONLY extend footplates when legs are raised.)	N/C

Footplates			
<input type="checkbox"/>		F - Centermount Single/Dual Footplate <i>(Single Footplate can be unlocked to form individual footplates)</i>	\$247.00
	172F315	CMT Angle Adjustable Footplates	\$247.00

Fixed Center Mount Legrests			
<input type="checkbox"/>	172H128	Fixed Center Mount Power Recline (Choose individual angle adjustable footplates and custom options below) (F)	N/C

Custom Centermount Options 			
<input type="checkbox"/>	172DI26	Extra Large Single Footplate for Centermount Footrest. Size 16"W x 14"D	\$309.00
<input type="checkbox"/>	172DI27	Gel Foot Pads for XL Single Centermount Footplate.	\$309.00
<input type="checkbox"/>	172DI28	Extra Large Pair Curved Calf Pads for Center Mount Legrests. Fully adjustable	\$309.00

15. ACCESSORIES

Wheelchair Tray Table NEW!			
<input type="checkbox"/>	1727A06	Swing-Away Tray	\$331.00

Backpack			
<input type="checkbox"/>	172A1	Black Utility	\$97.00

Positioning Belt - Select for any style seat			
<input type="checkbox"/>	172EJ01	Auto Style Buckle	N/C
<input type="checkbox"/>	172EJ02	Belt 2" Aircraft Buckle	\$163.00
<input type="checkbox"/>	172EJ03	Belt 2" Aircraft Buckle Padded	\$178.00

Accessories			
<input type="checkbox"/>	1728A09	Ram Phone Mount	\$67.00
<input type="checkbox"/>	1728A10	Ram Tablet Mount	\$268.00
<input type="checkbox"/>	1728A11	Ram Cup holder	\$67.00
<input type="checkbox"/>	1728A12	Switch-it Remote Stop Switch	\$1,200.00

16. CUSTOM OPTIONS



Looking for more options? Contact our custom "built 4 me" department.

To Inquire about any unique customization for Quickie products, contact:
Phone: (800) 333-4000, Option #3 or Fax: (800) 333-9011