



# Quickie QM-7 Series Non-Powered Seating

March 2017

EIPW11

Account: _____	Ship To: _____
Date: _____	Name: _____
PO#: _____	Address: _____
Buyer: _____	City, Province, Postal Code: _____
Marked For: _____	Phone: _____

✓	Option #	Definition	Price	✓	Option #	Definition	Price
---	----------	------------	-------	---	----------	------------	-------

## 1. QUICKIE QM-7 Series Non-Powered Seating

### QUICKIE QM-710™ BC

<input type="checkbox"/>	Quickie QM-710™ BC w/ ASAP II Seat		
172PC67	Base Package		\$7,349.00

### QUICKIE QM-715™ HD BC

<input type="checkbox"/>	Quickie QM-715™ HD BC w/ ASAP II Seat		
172PC79	Base Package		\$7,885.00

### QUICKIE QM-720™ BD

<input type="checkbox"/>	Quickie QM-720™ BD w/ ASAP II Seat		
172PC81	Base Package		\$8,137.00



### BC Seating Specifications

<b>Typical Chair Weight</b>	<b>W/battery</b>	350lbs	<b>W/O Battery</b>	260lbs
<b>Transit</b>	ISO occupied approved, NA with QM-715 HD			
<b>Seat Type</b>	ASAP II			
<b>Seat size range</b>	<b>Width</b>	12"-24"	<b>Depth</b>	12"-22"
<b>STFH range</b>	<b>ASAP II</b>	16.3" -21"		
<b>Power Seat Elevator</b>	Power 9" Seat elevator, NA with QM-715 HD			

### Actuators Required by Power Seating Function

<b>12" Seat Elevate:</b>	1 actuator
<b>Power Leg Function (Centermount or Swing away):</b>	2 actuators
<b>Combination : Lift and Power Elevating Legrest. (Centermount or Swing away):</b>	3 actuators

### QUICKIE QM-7 Base Specifications

<b>Drive Wheel Position</b>	Midwheel Drive	<b>Drive Wheels</b>	14" Drive Wheels w/ flat free insert (Black or Grey)
<b>User weight Cap</b>	300 lb / 400LB	<b>Tires/Casters</b>	6" Front / 6" rear (Black or Grey)
<b>Base size</b>	<b>Width</b> 25" <b>Length</b> 36" <b>Turn Radius</b> 20"	<b>Charger</b>	8A Dual Mode Charger (Compact) CTE
<b>Drive Control</b>	R-Net PM120 (Non-Expandable and Expandable)	<b>Battery Size</b>	Group 24 / Group 34.
<b>Driving Range</b>	Up to 13.5 miles	<b>Shroud Colors</b>	12 standard colors, 6 upgrade prints.
<b>Drive Motors</b>	4-pole 6 mph, Optional 5mph and 8mph	<b>Ground Clearance</b>	4"

## 2. BASE CHAIR OPTIONS

Speed Package (4-pole)			
<input type="checkbox"/>	172PF85	<b>QM-710 Standard</b> 6 mph motor package	N/C
<input type="checkbox"/>	172PF91	QM-710 5 mph High Torque motor package	\$525.00
<input type="checkbox"/>	172PF90	QM-710 8 mph Community motor package	\$1,103.00
<input type="checkbox"/>	172PF84	<b>QM-715 HD Standard</b> 5 mph motor package	N/C
<input type="checkbox"/>	172PF89	<b>QM-720 Standard</b> 8 mph motor package	N/C

Shroud Color			
<input type="checkbox"/>	172EZ01	Candy Apple Red	N/C
<input type="checkbox"/>	172EZ02	High Gloss Black	N/C
<input type="checkbox"/>	172EZ03	Candy Blue	N/C
<input type="checkbox"/>	172EZ04	Silver	N/C
<input type="checkbox"/>	172EZ05	Copper	N/C
<input type="checkbox"/>	172EZ08	Green Apple	N/C
<input type="checkbox"/>	172EZ09	Yellow	N/C
<input type="checkbox"/>	172EZ10	Pearl Pink	N/C
<input type="checkbox"/>	172EZ11	Candy Purple	N/C
<input type="checkbox"/>	172EZ12	Matte Black	N/C
<input type="checkbox"/>	172EZ13	Midnight Blue	N/C
<input type="checkbox"/>	172EZ15	Pearl White	N/C
<input type="checkbox"/>	172EZ06	Digital Camo.	\$525.00
<input type="checkbox"/>	172EZ07	Woodland Camo.	\$525.00
<input type="checkbox"/>	172EZ14	Carbon Fiber Print	\$525.00
<input type="checkbox"/>	172EZ16	American Flag.	\$525.00
<input type="checkbox"/>	172EZ17	Pink Camo	\$525.00
<input type="checkbox"/>	172EZ18	Grunge Skulls	\$525.00

User weight			
<input type="checkbox"/>	172WG06	Less than 150lbs	N/C
<input type="checkbox"/>	172WG13	151-200lb	N/C
<input type="checkbox"/>	172WG14	201-250lb	N/C
<input type="checkbox"/>	172WG03	251-300lb	N/C
<input type="checkbox"/>	172WG15	301-350lb	N/C
<input type="checkbox"/>	172WG16	351-400lb	N/C
		Transit Package NA	N/C

Drive Wheel			
<input type="checkbox"/>	172DW59	14" Cast Wheel Solid Dual Density Inserts <b>Black Knobby Tires</b> and <b>Black Caster</b>	Std

Suretrac Drive Tracking <span style="background-color: orange; color: white; padding: 2px;">NEW!</span>			
<input type="checkbox"/>	172MI01	Suretrac Module and Motor Encoders	\$1,475.00


Transit Option			
<input type="checkbox"/>	172TO25	4-Point Transit Tie-down Package <i>Meets ISO-7176 for Occupant Use (NA with QM715 H)</i>	\$373.00
<input type="checkbox"/>	172TO40	Omit Transit Tie-down	NC

Wheel locks			
<input type="checkbox"/>	172WL42	Wheel locks	\$142.00

Ventilator Options			
<input type="checkbox"/>	172VT4	Fixed Vent Tray (Pulmonetic LTV Models)	\$525.00

Pre-Set Seat Frame Angle (STFH increases 1" per increment change)			
<input type="checkbox"/>	172FA10	0°	N/C
<input type="checkbox"/>	172FA3	3°	Not available with Tilt or Lift
<input type="checkbox"/>	172FA4	6°	N/C

Seat Height (Use STFH Matrix below)			
<input type="checkbox"/>	172HT44	Position "0"	Requires Group 34 Batteries
<input type="checkbox"/>	172HT49	+5" STFH	N/C
<input type="checkbox"/>	172HT45	+1" STFH	N/C
<input type="checkbox"/>	172HT46	+2" STFH	N/C
<input type="checkbox"/>	172HT47	+3" STFH	N/C

LED Light Package			
<input type="checkbox"/>	172LI05	 LED Light Package: Includes headlights, tail lights, turn indicators and flashers.	\$1,066.00

STFH Matrix					
Seating Type and Power Option Combinations	Seat-to-Floor Height				
	0	+5"	+1"	+2"	+3"
BC ASAP II No power seating.*	16.3	16.8	17.3	18.3	19.3
BC ASAP II 9" Lift only.	17	17.5	18	19	20

\*STFH measured with batteries but no user. Actual STFH may vary by .5"  
 \*\*Adding seat angle will increase STFH by 1" for each incremental change.  
 \*\*Ultra low STFH may cause ground clearance of Power Centermount to drop below 2".

## 3. BATTERIES and CHARGERS

Batteries			
<input type="checkbox"/>	172AG04	MK Gr 24 Flat Top Gel Cell (Pair) w/straps <i>Recommended for Max Daily Range</i>	\$830.00
<input type="checkbox"/>	172AG11	MK Gr 34 Flat Top Gel Cell (Pair) w/straps, <i>Required for lower STFH.</i>	\$830.00

Battery Charger			
<input type="checkbox"/>	172BJ19	8A Dual Mode (Compact) CTE	N/C

## 4. ELECTRONICS

### Drive Control Options

#### Non-Expandable Controllers

<input type="checkbox"/>	172ER05	R-Net Non-Expandable Drive Control (120A)	N/C
--------------------------	---------	---	-----

#### Expandable Controllers

<input type="checkbox"/>	172ER06	R-Net Expandable Drive Control (120A)	\$1,066.00
	172EX05	Expandable 120A Controller - R-net	\$1,066.00
	172EX06	Harness for 120A Expandable Controller	N/C

### R-net Joystick

<input type="checkbox"/>	172J70	Std LED Joystick for Rnet Controls	N/C
<input type="checkbox"/>	172J108	Joystick w/Stereo Jacks & Color Display	N/C
<input type="checkbox"/>	172J88	Joystick W/Toggle, Spd Pot, Stereo Jack & Color Display	\$158.00
<input type="checkbox"/>	172J112	Joystick Advanced: 3.5" Color Display, Paddle Switches, Shortcut Keys and Stereo Jack.	\$1,133.00
<input type="checkbox"/>	172J95	Omit JS. For use with OMNI	N/C

### Joystick Mount (Fixed, Swing Away, Retractable, Mid)

<input type="checkbox"/>	172JM10	Fixed Mount	N/C
<input type="checkbox"/>	172JM11	Swing Away Mount	\$388.00
<input type="checkbox"/>	172JM15	Retractable Joystick mount w/45° stops	\$721.00
<input type="checkbox"/>	172JM14	Mid Line Mount	\$388.00

### Joystick Mount (Location)

*Height Adjustable.*

<input type="checkbox"/>	172HM1	Right Side Mount	N/C
<input type="checkbox"/>	172HM2	Left Side Mount	N/C

### Joystick Handle (Optional)

<input type="checkbox"/>	172JH1	Ball Handle	\$95.00
<input type="checkbox"/>	172JH2	Foam Ball Handle	\$95.00
<input type="checkbox"/>	172JH4	T-Handle	\$95.00
<input type="checkbox"/>	172JH6	Goalpost	\$95.00

### Attendant Control

<input type="checkbox"/>	1720A04	Attendant Control for R-net	\$1,066.00
--------------------------	---------	-----------------------------	------------

### Attendant Mounting location

<input type="checkbox"/>	172DN1	Mounting location: Right Back Cane	N/C
<input type="checkbox"/>	172DN2	Mounting Location: Left Back Cane	N/C

### R-net OMNI and Output Options

*(For expandable controls, choose all that apply)*

<input type="checkbox"/>	172ON6	Rnet OMNI Input Control Module and display with built in Infrared and armrest mount. Choose mount below. (required for Specialty Input Devices.)	\$2,075.00
<input type="checkbox"/>	172ON7	Rnet Output Module (ECM)	\$1,066.00
<input type="checkbox"/>	172ON8	Rnet Bluetooth Mouse-mover for Windows / Android	\$1,066.00
<input type="checkbox"/>	172ON9	Rnet Bluetooth for Apple Devices	\$1,066.00

### OMNI Mount

*(If used with Joystick must be mounted opposite side)*

<input type="checkbox"/>	172DS04	Swing Away Left Side OMNI Mount	N/C
<input type="checkbox"/>	172DS03	Swing Away Right Side OMNI Mount	N/C

### Switch inputs for Joysticks or OMNI

*Drive or Mode / Profile Control*

<input type="checkbox"/>	172SW22	Stereo to 2 Mono Splitter	\$75.00
<input type="checkbox"/>	172SW6	Ribbon Switch	\$119.00
<input type="checkbox"/>	172SW1	Buddy Button Switch	\$95.00
<input type="checkbox"/>	172SW2	Disc Switch	\$202.00
<input type="checkbox"/>	172SW9	Wobble Switch	\$630.00
<input type="checkbox"/>	172SW4	Micro Light Switch	\$95.00
<input type="checkbox"/>	172SW24	Egg Switch-Black	\$119.00
<input type="checkbox"/>	172SW21	Dual toggle switch and bracket	\$214.00
<input type="checkbox"/>	172SW35	Switch-it Mini 1 button on Mono plug	\$268.00
<input type="checkbox"/>	172SW36	Switch-it Mini 2 button on 2 Mono plugs	\$268.00
<input type="checkbox"/>	172SW37	Switch-it Mini 2 button on 1 Stereo plug	\$268.00
<input type="checkbox"/>	172SW38	Switch-it Mini 4 button on DB9 connector	\$268.00

### Single Switch Mount Bracket

<input type="checkbox"/>	172EM01	Single Switch Mount Bracket	\$247.00
--------------------------	---------	-----------------------------	----------

### Electronic Programmers

<input type="checkbox"/>	172PR5	R-net PC Programmer with Dongle	\$615.00
--------------------------	--------	---------------------------------	----------

## 5. Specialty Driver Control Options

Requires expandable Rnet controls and OMNI unless otherwise noted. Not available with Upholstered backs.

### 5a. Proportional Control Packages

Head Control Proportional Packages		
<input type="checkbox"/> <b>Switch-it Dual Pro Proportional Head Control with Whitmyer and Facial Pads (Adult)</b> <b>\$6,541.00</b>		
172ID72	Switch-it Dual Pro Head Control w/facial pad (Adult)	\$6,139.00
172J09	SA Boom Mount with Egg Switch Control	\$402.00
<input type="checkbox"/> <b>Switch-it Dual Pro Proportional Head Control with Spot Pads Package</b> <b>\$6,541.00</b>		
172ID70	Switch-it Dual Pro Head Control w/Spot pads	\$6,139.00
172J09	SA Boom Mount with Egg Switch Control	\$402.00

Chin Control Proportional Joystick Packages		
<input type="checkbox"/> <b>Switch-it Versa-Guide Chin Control Compact Joystick and Egg Switch Package</b> <b>\$2,154.00</b>		
172ID71	Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
172IM07	S/A Mount off backrest	NC
172J04	SA Boom Mount with egg switch	\$115.00
<input type="checkbox"/> <b>Switch-it Versa-Guide Chin Control Compact Joystick and Proximity Flex Switches Package</b> <b>\$2,940.00</b>		
172ID71	Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
172IM08	Dual Flex Proximity Switches (requires no force)	\$649.00
172IM07	S/A Mount off backrest	\$252.00
<input type="checkbox"/> <b>Switch-It Versa-Guide Bib Chin Control Compact Joystick with Proximity Flex switches Package Std 6"</b> <b>\$2,688.00</b>		
172ID71	Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
172IM08	Dual Flex Proximity Switches (requires no force)	\$649.00
172IM05	Bib mount Small	\$757.00
<input type="checkbox"/> <b>Switch-It Versa-Guide Bib Chin Control Compact Joystick with Proximity Flex Switches Package Large 8"</b> <b>\$3,512.00</b>		
172ID71	Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
172IM08	Dual Flex Proximity Switches (requires no force)	\$649.00
172IM06	Bib mount Large	\$824.00
<input type="checkbox"/> <b>Switch-It Versa-Guide Bib Chin Control Compact Joystick with Egg Switch Package Std 6"</b> <b>\$2,911.00</b>		
172ID71	Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
172IM05	Bib mount Small	\$757.00
172J04	SA Boom Mount with egg switch	\$115.00
<input type="checkbox"/> <b>Switch-It Versa-Guide Bib Chin Control Compact Joystick with Egg Switch Package Large 8"</b> <b>\$2,978.00</b>		
172ID71	Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
172IM06	Bib mount Large	\$824.00
172J04	SA Boom Mount with egg switch	\$115.00
<input type="checkbox"/> <b>Switch-It Micro-Pilot Bib Chin Control with Egg Switch Package Std 6"</b> <b>\$4,406.00</b>		
172ID62	Micro-Pilot Proportional Joystick	\$3,534.00
172IM05	Bib mount Small	\$757.00
172J04	SA Boom Mount with egg switch	\$115.00
<input type="checkbox"/> <b>Switch-It Micro-Pilot Bib Chin Control with Egg Switch Package Large 8"</b> <b>\$4,473.00</b>		
172ID62	Micro-Pilot Proportional Joystick	\$3,534.00
172IM06	Bib mount Large	\$824.00
172J04	SA Boom Mount with egg switch	\$115.00
<input type="checkbox"/> <b>Switch-It Micro-Guide Bib Chin Control with Egg Switch Package Std 6"</b> <b>\$5,230.00</b>		
172ID63	Micro-Guide Proportional Joystick	\$4,358.00
172IM05	Bib mount Small	\$757.00
172J04	SA Boom Mount with egg switch	\$115.00
<input type="checkbox"/> <b>Switch-It Micro-Guide Bib Chin Control with Egg Switch Package Large 8"</b> <b>\$5,297.00</b>		
172ID63	Micro-Guide Proportional Joystick	\$4,358.00
172IM06	Bib mount Large	\$824.00
172J04	SA Boom Mount with egg switch	\$115.00
<input type="checkbox"/> <b>ASL Proportional Micro Extremity Joystick Chin / Lip with Egg Switch Control Package</b> <b>\$4,451.00</b>		
172ID59	Micro Extremity Joystick	\$4,084.00
172IM07	S/A Mount off backrest	\$252.00
172J04	SA Boom Mount with egg switch	\$115.00

Hand Control Proportional Joystick Packages		
<input type="checkbox"/> <b>Switch-it Versa-Guide Compact Joystick Midline with Gatlin Package</b> <b>\$2,369.00</b>		
172ID71	Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
172J07	S/A Gatlin Mount	\$330.00
<input type="checkbox"/> <b>Switch-it Versa-Guide Compact Joystick Midline and prox flex switches with Gatlin Package</b> <b>\$3,018.00</b>		
172ID71	Switch-it Versa-Guide Compact Proportional Joystick	\$2,039.00
172IM08	Dual Flex Proximity Switches (requires no force)	\$649.00
172J07	S/A Gatlin Mount	\$330.00
<input type="checkbox"/> <b>Switch-it Micro-Pilot Joystick with Space Disc Package</b> <b>\$4,007.00</b>		
172ID62	Micro-Pilot Proportional Joystick	\$3,534.00
172J08	S/A SpaceDisc mount	\$473.00
<input type="checkbox"/> <b>Switch-it Micro-Guide Joystick with SpaceDisc Package</b> <b>\$4,831.00</b>		
172ID63	Micro-Guide Proportional Joystick	\$4,358.00
172J08	S/A SpaceDisc mount	\$473.00
<input type="checkbox"/> <b>ASL Proportional Micro Extremity Joystick (ASL130) Midline with Gatlin Package</b> <b>\$4,414.00</b>		
172ID59	Micro Extremity Joystick	\$4,084.00
172J07	S/A Gatlin Mount	\$330.00
<input type="checkbox"/> <b>ASL Mushroom Joystick (ALS135) Midline with Gatlin Package</b> <b>\$3,353.00</b>		
172ID60	Mushroom Joystick	\$3,023.00
172J07	S/A Gatlin Mount	\$330.00
<input type="checkbox"/> <b>ASL Compact Joystick w Integrated buttons (ASL134) Midline with Gatlin Package</b> <b>\$1,530.00</b>		
172ID61	ASL Compact Joystick	\$1,200.00
172J07	S/A Gatlin Mount	\$330.00
<input type="checkbox"/> <b>ASL Micro Mini Joystick w Hand Pad (ASL136) Midline with Gatlin Package</b> <b>\$4,012.00</b>		
172ID68	Micro Mini Joystick	\$3,682.00
172J07	S/A Gatlin Mount	\$330.00

## 5. Specialty Driver Control Options (cont.)

### 5b. Switch Control

*Requires Headrest and/or Headrest Mount*

<b>3 Switch Proximity Head Arrays</b> Choose headrest mount below			
<input type="checkbox"/> <b>ASL Proximity Head Array Package (Adult)</b> <b>\$5,755.00</b>			
172II02		ASL Proximity Head Array (Adult)	\$5,353.00
172J09		SA Boom Mount with Egg Switch Control	\$402.00
<input type="checkbox"/> <b>ASL Proximity Head Array Package (Ped)</b> <b>\$5,755.00</b>			
172II03		ASL Proximity Head Array (Ped)	\$5,353.00
172J09		SA Boom Mount with Egg Switch Control	\$402.00
<input type="checkbox"/> <b>Switch-It Head Array Package 3 Switch Curved (Long)</b> <b>\$5,755.00</b>			
172II13		Switch-it Proximity Head Array Long	\$5,353.00
172J09		SA Boom Mount with Egg Switch Control	\$402.00
<input type="checkbox"/> <b>Switch-It Head Array Package 3 switch curved (Stubby)</b> <b>\$5,691.00</b>			
172II14		Switch-it Proximity Head Array Stubby	\$5,289.00
172J09		SA Boom Mount with Egg Switch Control	\$402.00

<b>5 Switch Proximity Head Array (No Egg Switch Needed)</b> Choose "Headrest Mount for Special Controls"			
<input type="checkbox"/>	172II04	Switch-It Head Array 5 switch flat	\$5,353.00

<b>Multiple Switches.</b>			
<input type="checkbox"/>	172SW31	Switch-it 3 Naked proximity switches on DB9 connector. (requires OMNI)	\$3,039.00
<input type="checkbox"/>	172SW32	Switch-it 4 Naked proximity switches on DB9 connector. (requires OMNI)	\$3,718.00
<input type="checkbox"/>	172SW33	Switch-it 5 Naked proximity switches on DB9 connector. (requires OMNI)	\$4,146.00
<input type="checkbox"/>	172SW34	ASL (502) 5-Switch adaptor box on DB9 connector. (for use with 5 mono switches, not included)	\$4,146.00
<input type="checkbox"/>	172SW23	4-way toggle switch button DB9 Connector. Includes bracket	\$670.00

<b>Sip N Puff Package</b> Choose headrest and HR mount below			
<input type="checkbox"/> <b>Quickie Sip n' Puff System</b> <b>\$2,635.00</b>			
172ID34		Sip N' Puff Kit	\$2,063.00
172IH1		Sip n' Puff Tubing and SA Hardware	\$457.00
172J10		SA Boom Mount with egg switch	\$115.00

<b>Headrest Mount for Specialty Controls</b>			
<input type="checkbox"/>	172MH12	Axys Mount & Bracket for Specialty Controls	\$309.00
<input type="checkbox"/>	172MH13	Cobra Mount & Bracket for Specialty Controls	\$386.00
<input type="checkbox"/>	172MH14	Linx Mount & Bracket for Specialty Controls	\$309.00

<b>Headrest Pad for Specialty controls</b>			
<input type="checkbox"/>	172HE79	Plush Pad 8" for Specialty controls	\$309.00
<input type="checkbox"/>	172HE80	Plush Pad 10" for Specialty controls	\$309.00
<input type="checkbox"/>	172HE81	Plush Pad 14" for Specialty controls	\$309.00

<b>Extra Swing Away Head Boom Egg Switch and Mount</b>			
	172J05	Extra SA Boom Mount with egg switch #1	\$115.00
	172J06	Extra SA Boom Mount with egg switch #2	\$115.00

<b>Single Switches</b>			
<input type="checkbox"/>	172IB01	Ribbon Switch	\$119.00
<input type="checkbox"/>	172IB02	Buddy Button Switch	\$108.00
<input type="checkbox"/>	172IB03	Disc Switch	\$202.00
<input type="checkbox"/>	172IB04	Wobble Switch	\$630.00
<input type="checkbox"/>	172IB05	Micro Light Switch	\$95.00
<input type="checkbox"/>	172IB06	Egg Switch-Black	\$134.00
<input type="checkbox"/>	172IB07	Switch-it Mini 1 button on Mono plug	\$268.00

## 6. SEAT TYPE

<input type="checkbox"/> 172ST56	<b>Quickie A.S.A.P II</b> (Adjustable Seating and Positioning) Choose Seat Frame Size Platform.	\$830.00
----------------------------------	---	----------

### ASAP II Frame size

(each frame is width and depth adjustable within given range. )

<input type="checkbox"/> 172FU17	Adjustable 12x12 to 16X19 ( A )	N/C
<input type="checkbox"/> 172FU18	Adjustable 16x12 to 20x20 ( B )	N/C
<input type="checkbox"/> 172FU19	Adjustable 18x16 to 22x22 ( C )	N/C

## 7. SEAT SIZE

### Seat Size Matrix. Seat type is chosen in section 5.

**A = ASAP II Seat Frame**  
(Adjustable from 12x12 to 16X19)

**B = ASAP II Seat Frame**  
(Adjustable from 16X12 to 20X20)

**C = ASAP II Seat Frame**  
(Adjustable from 18X16 to 22X22)

		ASAP II Seat Width									
		12"	13"	14"	15"	16"	17"	18"	19"	20"	22"
Seat Depth	12"	A	A	A	A	A,B	B	B	B	B	
	13"	A	A	A	A	A,B	B	B	B	B	
	14"	A	A	A	A	A,B	B	B	B	B	
	15"	A	A	A	A	A,B	B	B	B	B	
	16"	A	A	A	A	A,B	B	B,C	B,C	B,C	C
	17"	A	A	A	A	A,B	B	B,C	B,C	B,C	C
	18"	A	A	A	A	A,B	B	B,C	B,C	B,C	C
	19"	A	A	A	A	A,B	B	B,C	B,C	B,C	C
	20"					B	B	B,C	B,C	B,C	C
	21"							C	C	C	C
	22"							C	C	C	C

Seat Width			
<input type="checkbox"/> 172R21	12"	Seat width	(N/C)
<input type="checkbox"/> 172R22	13"	Seat width	(N/C)
<input type="checkbox"/> 172R23	14"	Seat width	(N/C)
<input type="checkbox"/> 172R24	15"	Seat width	(N/C)
<input type="checkbox"/> 172R5	16"	Seat width	(N/C)
<input type="checkbox"/> 172R6	17"	Seat width	(N/C)
<input type="checkbox"/> 172R7	18"	Seat width	(N/C)
<input type="checkbox"/> 172R8	19"	Seat width	(N/C)
<input type="checkbox"/> 172R9	20"	Seat width	(N/C)
<input type="checkbox"/> 172R25	22"	Seat width	(N/C)

Seat Depth			
<input type="checkbox"/> 172RF20	12"	Seat depth	(N/C)
<input type="checkbox"/> 172RF21	13"	Seat depth	(N/C)
<input type="checkbox"/> 172RF22	14"	Seat depth	(N/C)
<input type="checkbox"/> 172RF23	15"	Seat depth	(N/C)
<input type="checkbox"/> 172RF5	16"	Seat depth	(N/C)
<input type="checkbox"/> 172RF6	17"	Seat depth	(N/C)
<input type="checkbox"/> 172RF7	18"	Seat depth	(N/C)
<input type="checkbox"/> 172RF8	19"	Seat depth	(N/C)
<input type="checkbox"/> 172RF24	20"	Seat depth	(N/C)
<input type="checkbox"/> 172RF26	21"	Seat depth	(N/C)
<input type="checkbox"/> 172RF25	22"	Seat depth	(N/C)

24" Seat Width			
<input type="checkbox"/> 172VU16	24"	Seat width (fixed width, depth is adjustable 16" to 22" deep) <i>N/A with Power or manual recline back or sling back upholstery.</i>	\$525.00

## 8. SEAT CUSHIONS

JAY Cushion Bundling			
<input type="checkbox"/> 172JC14		Jay Cushion <u>On Chair Order</u> <u>See order Form for codes</u>	See Order Form

## 9. POWERED SEATING

### Powered Seat Elevator- (AKA: Power Lift)

<input type="checkbox"/>	172PA27	12" Power Seat Elevate for use with tilt ( <i>Must be ordered with SPOT tilt, NA with QM-715 HD</i> )(1 Actuator)	\$3,452.00
--------------------------	---------	---	------------

### Power Seating Control

#### Powered Seating Switch Control Method

(For use with power seating actuators. Switch activation is stand-alone and does not require expandable or thru-drive electronics.)

#### Single Actuator Switch Control (Lift)

<input type="checkbox"/>	172EO10	Switch Control 1 actuator with <u>Dual Toggle</u>	N/C
<input type="checkbox"/>	172EO11	Switch Control 1 actuator with <u>Buddy Button</u>	N/C
<input type="checkbox"/>	172EO19	Switch Control 1 actuator with <u>Switch-it Mini 1 button</u>	N/C
<input type="checkbox"/>	172EO20	Switch Control 1 actuator with <u>Switch-it Mini 2 button</u>	N/C

#### 2-3 Actuator Switch Control (Lift/Legs, Lift/Tilt or Legs only)

<input type="checkbox"/>	172EO21	Switch Control 2-3 actuator with <u>Dual Toggle</u>	\$750.00
<input type="checkbox"/>	172EO12	Switch Control 2-3 actuator with <u>4-way toggle</u>	\$750.00
<input type="checkbox"/>	172EO23	Switch Control 2-3 actuator with <u>Switch-it Mini 2 button</u>	\$750.00
<input type="checkbox"/>	172EO22	Switch Control 2-3 actuator with <u>Switch-it Mini 4 button</u>	\$750.00

#### Switch Mount

*(for use with switch control powered seating)*

<input type="checkbox"/>	172EY01	Switch Mount Right	N/C
<input type="checkbox"/>	172EY02	Switch Mount Left	N/C

#### Powered Seating Thru-drive Control Method.

##### Thru-Drive 1 Actuator

<input type="checkbox"/>	172EO28	Thru Drive 1 Actuator (For operating Lift thru Joystick.) Includes Ctrl+5	\$1,751.00
--------------------------	---------	---	------------

##### Thru-Drive 2+ Actuator

<input type="checkbox"/>	172EO29	Thru Drive 2+ Actuators: (For operating lift, power legs) Includes Ctrl+5	\$2,525.00
--------------------------	---------	---	------------

#### Powered Seating Thru-drive Control Options

##### Assignable Buttons

"Assignable Button" feature allows for customization of joystick buttons and/or Ctrl+5 Switches allowing access to various wheelchair and seating functions. When enabled, chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle.

<input type="checkbox"/>	172PM01	Enable Assignable buttons	N/C
--------------------------	---------	---------------------------	-----

##### Ctrl+5 Switchbox Options

Ctrl+5 allows direct control to 5 power seating functions (10 total switches). Switches can also be used with assignable buttons for greater flexibility. Box includes 3-Amp USB port for charging portable devices



<input type="checkbox"/>	172QG01	Ctrl+5 Switchbox mounted behind joystick	N/C
<input type="checkbox"/>	172QG02	Ctrl+5 Switchbox mounted opposite joystick	N/C
<input type="checkbox"/>	172QG03	Ctrl+5 Switchbox mounted on back cane	N/C
<input type="checkbox"/>	172QG05	Omit Ctrl+5	-\$618.00

#### Ctrl+5 Switch Type

<input type="checkbox"/>	172EY03	Push buttons	N/C
<input type="checkbox"/>	172EY04	Toggle switches	\$309.00

## 10. BACKREST AND BACK OPTIONS for ASAP II

#### Backrest Frame Type

<input type="checkbox"/>	172BT8	Angle Adjustable ( <i>ASAP II Seat only</i> )	N/C
<input type="checkbox"/>	172BT98	Manual Recline 40° ( <i>ASAP II Seat only</i> ) <i>not width adjustable</i>	\$536.00

#### Manual Recline Adjustment Method Option

*Back Cane Height 23".*

<input type="checkbox"/>	172MN01	Tool	N/C
<input type="checkbox"/>	172MN02	Knob (NA with Cantilever or dual post arms.)	N/C

#### Angle Adjustable Backrest Options

*(NA with Manual recline)*

<input type="checkbox"/>	172OD8	Non-Folding Backrest	N/C
<input type="checkbox"/>	172OD9	Folding Backrest	\$137.00

#### Cane Angle

<input type="checkbox"/>	172BP40	Straight Back Post	N/C
<input type="checkbox"/>	172BP41	8° Bend Back Post ( <i>Avail. only w/ Upholstery or J3 Back</i> )	\$126.00

#### Towel Bar Depth

<input type="checkbox"/>	172QT03	Shallow (NA with Upholstery backs.)	N/C
<input type="checkbox"/>	172QT04	Deep (Extends 4" from back canes for deep backs)	N/C

#### Back Cane

*(Upholstery Height Measured pan to top of upholstery.)*

<input type="checkbox"/>	Low (15.5") 172B1	<input type="checkbox"/>	Med (17.5")172B5	<input type="checkbox"/>	Tall (19.5") 172B9
--------------------------	-------------------	--------------------------	------------------	--------------------------	--------------------

#### Back Type (Upholstery or Solid)

<input type="checkbox"/>	172BU1	Upholstery ( <i>Not avail with tilt or recline</i> )	N/C
<input type="checkbox"/>	172BU70	Tension adj Upholstery Ballistic Nylon. (with 1" width adjustment.) ( <i>Not avail with tilt or recline</i> )	\$394.00
<input type="checkbox"/>	172BU71	Tension adj Upholstery 3DX Vented. (with 1" width adjustment) ( <i>Not avail with tilt or recline</i> )	\$415.00
<input type="checkbox"/>	172BU76	S.P.O.T. Curved	N/C
<input type="checkbox"/>	172BU77	Jay SPO Posterior Positioning ( <i>ASAP II only</i> )	\$598.00
<input type="checkbox"/>	172BU78	ART SPO Custom Back ( <i>ASAP II only</i> )	See order form
<input type="checkbox"/>	172BU14	Omit Back.	-\$15.00

#### Jay 3 Option (select Omit above)

<input type="checkbox"/>	172JB08	Jay J3 <u>On Chair Order</u> <i>Choose back, headrest, laterals</i>	See Order Form
--------------------------	---------	--	----------------

#### Back Height - Required for SPOT Curved, SPO Backs.




*(Limited combinations of widths) Measured from pan to top of solid back.*

<input type="checkbox"/>	172B132	14"		N/C
<input type="checkbox"/>	172B36	16"		N/C
<input type="checkbox"/>	172B38	18"		N/C
<input type="checkbox"/>	172B40	20"		N/C
<input type="checkbox"/>	172B42	22"		N/C
<input type="checkbox"/>	172B44	24"		N/C
<input type="checkbox"/>	172B140	26"		N/C

## 12. HEADREST OPTIONS


Choose One Each.

Headrest Mount			
<input type="checkbox"/>	172MH10	Axys Mount & Bracket	\$315.00
<input type="checkbox"/>	172MH11	Cobra Mount and Bracket	\$578.00
<input type="checkbox"/>	172MH1	Linx Mount and Bracket	\$237.00



**NEW!**

Headrest Pad			
<input type="checkbox"/>	172HE65	Plush Head Support Pad 8"	\$237.00
<input type="checkbox"/>	172HE66	Plush Head Support Pad 10"	\$237.00
<input type="checkbox"/>	172HE78	Plush Head Support Pad 14"	\$237.00
<input type="checkbox"/>	172HE82	Adjusta Plush Head Support Pad Narrow	\$365.00
<input type="checkbox"/>	172HE83	Adjusta Plush Head Support Pad std	\$365.00



## 13. POSTURAL SUPPORTS

### LATERAL THORACIC SUPPORT (Swing Away)

*Available on all Solid Backs.*

<input type="checkbox"/>	<b>LAT 4X4 S/A MOUNTS W/ 1" OFFSET - PR</b>		<b>\$409.00</b>
	172ZV14	LATERAL 4X4 W/ 1" OFFSET - PAIR	\$126.00
	172ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$283.00
<input type="checkbox"/>	<b>LAT 4X5 S/A MOUNTS W/ 1" OFFSET - PR</b>		<b>\$409.00</b>
	172ZV15	LATERAL 4X5 W/ 1" OFFSET - PAIR	\$126.00
	172ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$283.00
<input type="checkbox"/>	<b>LAT 4X6 S/A MOUNTS W/ 1" OFFSET - PR</b>		<b>\$409.00</b>
	172ZV16	LATERAL 4X6 W/ 1" OFFSET - PAIR	\$126.00
	172ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$283.00
<input type="checkbox"/>	<b>LAT 6X4 S/A MOUNTS W/ 1" OFFSET - PR</b>		<b>\$409.00</b>
	172ZV17	LATERAL 6X4 W/ 1" OFFSET - PAIR	\$126.00
	172ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$283.00
<input type="checkbox"/>	<b>LAT 6X5 S/A MOUNTS W/ 1" OFFSET - PR</b>		<b>\$409.00</b>
	172ZV18	LATERAL 6X5 W/ 1" OFFSET - PAIR	\$126.00
	172ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$283.00
<input type="checkbox"/>	<b>LAT 6X6 S/A MOUNTS W/ 1" OFFSET - PR</b>		<b>\$409.00</b>
	172ZV19	LATERAL 6X6 W/ 1" OFFSET - PAIR	\$126.00
	172ZS07	LATERAL SWING AWAY MOUNT - PAIR	\$283.00

### CHEST STRAP

<input type="checkbox"/>	172ZR03	Padded Chest Strap	\$283.00
--------------------------	---------	--------------------	----------

### Swing Away Knee Adductor Supports

*Available only with Fixed and Power Center Mount Legrests.*

*Not available with Dual Post Arm Rests or 2-point Locking Cantilever Arms*

Knee Adductor Supports - LEFT (choose Pad and Arm)			
<b>Pad Size</b>			
<input type="checkbox"/>	172KB5	Knee Adductor Pad SM- Left	\$215.00
<input type="checkbox"/>	172KB7	Knee Adductor Pad LG- Left	\$215.00
<b>S/A Arm Size</b>			
<input type="checkbox"/>	172XX04	Swing Away Short - Left	\$283.00
<input type="checkbox"/>	172XX06	Swing Away Long - Left	\$283.00
Knee Adductor Supports - RIGHT (choose Pad and Arm)			
<b>Pad Size</b>			
<input type="checkbox"/>	172KB6	Knee Adductor Pad SM- Right	\$215.00
<input type="checkbox"/>	172KB8	Knee Adductor Pad LG- Right	\$215.00
<b>S/A Arm Size</b>			
<input type="checkbox"/>	172XX05	Swing Away Short - Right	\$283.00
<input type="checkbox"/>	172XX07	Swing Away Long - Right	\$283.00

### Long Thigh Support w/Removable Mount

*Available with Cantilever Arms, NA with Manual Recline Back, Dual Post, Single Post or Cantilever w/ 2 pt Locking arms.*

<input type="checkbox"/>	<b>SPOT Long Thigh Support - Left</b>		<b>\$409.00</b>
	172TV2	Long Thigh Support - Left	\$126.00
	172TX1	Removable Mount - Left	\$283.00
<input type="checkbox"/>	<b>SPOT Long Thigh Support - Right</b>		<b>\$409.00</b>
	172TV1	Long Thigh Support - Right	\$126.00
	172TX2	Removable Mount - Right	\$283.00

## 14. ARMREST AND ARMPAD OPTIONS

Cantilever Armrest Mounting			
<b>Cantilever Armrests</b>			
<i>(Adjustable length arm pad included (choose Pad setting below, except for Low Arm Option) NA with Manual Recline Back</i>			
<input type="checkbox"/>	172EG12	<b>STD Cantilever Arm 9"-12" Ht Range</b>	Std
<input type="checkbox"/>	172EG13	<b>Low Cantilever Arm 7"-10" Ht Range</b> (Choose ASAPII pad type and size below)	NC
<input type="checkbox"/>	172EG10	<b>Standard Cantilever Arm w/ 2 Point Locking Post.</b> (9"-10" Ht Range.) NA w/ seat depths below 16"	\$289.00
<input type="checkbox"/>	172EG11	<b>Standard Cantilever Arm w/ 2 Point Locking Post.</b> (11"-12" Ht Range.) NA w/ seat depths below 16"	\$289.00

Cantilever Armpad				
Pad Setting Options for <b>Standard Cantilever Arm</b>				
<b>Left Arm Pad (Std)</b>		<b>Right Arm Pad (Std)</b>		
<input type="checkbox"/>	172AE11	Full Length Cantilever (14")	<input type="checkbox"/>	172AB11
<input type="checkbox"/>	172AE10	Desk Length Cantilever (10")	<input type="checkbox"/>	172AB10

Single - Dual Post Armrest Mounting			
<b>Single Post Armrests</b>			
<input type="checkbox"/>	172EG14	Single Post Height Adjustable - Std.	\$289.00
<input type="checkbox"/>	172EG15	Single Post Height Adjustable - Low	\$289.00

Dual Post Armrests (NA with Infinite Seat Depth)			
<input type="checkbox"/>	172EG09	Shrt Dual Post Flip Back HT Adj 10"-15"	\$378.00
<input type="checkbox"/>	172EG23	Lng Dual Post Flip Back HT Adj 10"-15"	\$378.00

Pad Type / Size Options for <b>Dual Post / Single Post / Low Cantilever</b>				
<b>Left Arm Pad (N/C)</b>		<b>Right Arm Pad (N/C)</b>		
<input type="checkbox"/>	172AE01	Desk Length Waterfall (10")	<input type="checkbox"/>	172AB01
<input type="checkbox"/>	172AE02	Full Length Waterfall (14") <i>NA w/Low Cantilever arm</i>	<input type="checkbox"/>	172AB02
<input type="checkbox"/>	172AE03	Desk Length Standard (10")	<input type="checkbox"/>	172AB03
<input type="checkbox"/>	172AE04	Full Length Standard (14") <i>NA w/Low Cantilever arm</i>	<input type="checkbox"/>	172AB04



## 15. HANGERS AND FOOTRESTS

### Swing-away Hangers

#### Non-Powered Swing Away Legrests

*Sold as Pairs*

<input type="checkbox"/>	172H2	70° Swing Away (B)	N/C
<input type="checkbox"/>	172H6	90° Swing Away (B) (5.5-9")	N/C
<input type="checkbox"/>	172H138	Manual ELR S/A (B) (Avail. with Fixed Centermount, choose below)	\$352.00

#### Power Elevating Legrests (2 Actuators)

<input type="checkbox"/>	172H126	Swing Away Power ELR w/ Articulation Pair (B footplates) (N/A for 12" or 13" wide)	\$2,994.00
--------------------------	---------	---	------------

#### Footplates

<input type="checkbox"/>		<b>B - Adult Angle Adjust w/Heel Loop for use with 70° hanger and ELR (Pair)</b>	<b>\$205.00</b>
	172F282	Adult Angle Adjustable Footplate	\$205.00
	172HL10	Heel Loop	N/C
<input type="checkbox"/>		<b>B - Adult Angle Adjust w/Heel Loop for use with 90° hanger (Pair)</b>	<b>\$205.00</b>
	172F285	Adult Angle Adjustable Footplate	\$205.00
	172HL10	Heel Loop	N/C
<input type="checkbox"/>		<b>C - Kids Angle Adjust w/Leg Strap for use with 70° hanger (Pair)</b>	<b>\$205.00</b>
	172F283	Kids Angle Adjustable Footplate	\$205.00
	172LS14	Leg Strap	N/C
<input type="checkbox"/>		<b>C - Kids Angle Adjust w/Leg Strap for use with 80° and 90° hanger(Pair)</b>	<b>\$205.00</b>
	172F286	Kids Angle Adjustable Footplate	\$205.00
	172LS14	Leg Strap	N/C
<input type="checkbox"/>	172F279	Omit Footplate and extension tube	N/C

#### Omit Hangers.

<input type="checkbox"/>	172H144	Omit Hangers	N/C
--------------------------	---------	--------------	-----

#### Hanger Extensions. (See range listed in matrix below)

*NOTE: Not all extension tube will be available; depending on the seat height, some extension tube options may not offer enough ground clearance.*

<input type="checkbox"/>	172E16	High Mount Hanger	\$84.00
<input type="checkbox"/>	172E1	Short	N/C
<input type="checkbox"/>	172E2	Medium	N/C
<input type="checkbox"/>	172E5	Long	N/C

#### Footrest Accessories

<input type="checkbox"/>	1725A01	Toe Loop (Pair) Only avail on (opt B).	\$47.00
<input type="checkbox"/>	172LS14	Extra Leg Strap	\$32.00

### Centermount Legrests

#### Non-Powered Center Mount Legrests

##### Fixed Center Mount Legrests (includes footplate)

**Plus options below**

<input type="checkbox"/>	172H134	Fixed Center ASAP II Mount Short	N/C
<input type="checkbox"/>	172H135	Fixed Center ASAP II Mount Long	N/C

#### Angle Adjustable Center Mount Legrests for ASAP II

<input type="checkbox"/>	172H104	Angle adjustable Center mount Med (I) <i>Independent leg length and angle adjustable</i>	Built-4-Me Only
<input type="checkbox"/>	172H108	Angle adjustable Center mount Large (I) <i>Independent leg length and angle adjustable</i>	Built-4-Me Only

#### Power Centermount (2 Actuators)

<input type="checkbox"/>	172H127	Power Centermount ELR w/7.5" Power Articulation (F) Centermount footplates plus Custom Options Below	\$2,987.00
--------------------------	---------	--	------------

#### Power Centermount Footplate Extension Option (For Thru-Drive Operation)

<input type="checkbox"/>	172IE01	<b>Enable</b> 7.5" Powered footplate extension for transfer assist. (Electronics programmed to allow footplates to be extended <b>with or without</b> legrest being raised)	\$515.00
<input type="checkbox"/>	172IE02	<b>Disable</b> Footplate extension for transfers. (Electronics programmed to <b>ONLY</b> extend footplates when legs are raised.)	N/C

#### Centermount Footplates

<input type="checkbox"/>		<b>F - Centermount Dual Footplate (Pair)</b> (Individually Height (3"-19") and Angle Adjustable)	<b>\$247.00</b>
	172F315	CMT Angle Adjustable Footplates	\$247.00
<input type="checkbox"/>		<b>I - Angle adj Centermount Footplate pair</b>	<i>Built-4-Me Only</i>
	172F125	Angle Adjustable Footplates	<i>Built-4-Me Only</i>

#### Custom Centermount Options



<input type="checkbox"/>	172DI26	Extra Large Single Footplate for Centermount Footrest. Size 16"W x 14"D	\$309.00
<input type="checkbox"/>	172DI27	Gell Foot Pads for XL Single Centermount Footplate.	\$309.00
<input type="checkbox"/>	172DI28	Extra Large Pair Curved Calf Pads for Center Mount Legrests. Fully adjustable	\$309.00

Lower Leg Length Matrix (measured seat pan to the top of the footplate surface)						Footplate Matrix			
Hanger Type	Standard	High Mount	Short	Medium	Long	B	F	I	J
70° Swing Away		7"-14"	13.5"-16"	16"-18.5"	18"-20.5"	X			
90° Swing Away		5.5-9" Footplate mounts to front of hanger.				X			
Manual ELR S/A			12.5"-17"	14.5"-19"	17"-21"	X			
Power ELR S/A			12"-14"	14"-17"	16.5"-19"	X			
Fixed Cntrmnt ASAP			12.5"-15.5"		14.5"-19.5"			Single Footplate	
Angle Adj Cntrmnt				3"-13"	3"-15"			X	X
Power Cntrmnt	3"-19"						X		

**16. ACCESSORIES**

Wheelchair Tray Table			
<input type="checkbox"/>	1727A01	TRAY TABLE SM TUFFAK 12"-14" PWR	\$331.00
<input type="checkbox"/>	1727A02	TRAY TABLE MED TUFFAK 15"-17" PWR	\$331.00
<input type="checkbox"/>	1727A03	TRAY TABLE LRG TUFFAK 18"-20" PWR	\$331.00

Positioning Belt - Select for any style seat			
<input type="checkbox"/>	172EJ01	Auto Style Buckle	N/C
<input type="checkbox"/>	172EJ02	Belt 2" Aircraft Buckle	\$163.00
<input type="checkbox"/>	172EJ03	Belt 2" Aircraft Buckle Padded	\$178.00

Backpack			
<input type="checkbox"/>	172A1	Black Utility	\$97.00

Accessories			
<input type="checkbox"/>	1728A09	Ram Phone Mount on Round JS tube	\$67.00
<input type="checkbox"/>	1728A10	Ram Tablet Mount on Round JS tube	\$268.00
<input type="checkbox"/>	1728A11	Ram Cup holder	\$67.00
<input type="checkbox"/>	1728A12	Switch-it Remote Stop Switch	\$1,200.00

**17. ACCESSORIES**



Save time! Wheelchair and seating will be delivered together.  
Hardware pre-installed.

Cushions			
<input type="checkbox"/>	Jay Cushion:		
	Width:		Depth:
	Notes:		

Backs			
<input type="checkbox"/>	Jay Back:		
	Width:		Height:
	Notes:		

Headrest	
<input type="checkbox"/>	System:
	Pads:
	Notes:
	Hardware:
	Notes:

Jay Positioning Components	
<input type="checkbox"/>	Anterior Trunk Supports:
	Size:
	Hand and Arm Supports:
	Size:
	Pelvic Positioning:
	Size:
	Foot Positioning:
	Size:
	Notes: