



Quickie Xplore2 HD

July 2018

EIPW24



ACCOUNT INFORMATION

Dealer Contact Name:	_____	P.O. # or Quote:	_____
Email Address:	_____	User Name or Number:	_____
Phone Number:	_____ Fax Number: _____	User Weight:	_____
Bill to (Name):	_____	User Height:	_____
Bill to (Location):	_____	Seat Cushion Thickness:	_____
Ship to (Name):	_____	Knee to Heel Length: Left: _____ Right: _____	
Address:	_____		
City:	_____		
Province:	_____ Postal Code: _____		

The order cannot be processed without the following:
 - A user name or number, the client's measurements and weight
 - Sunrise Medical is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
 - Order confirmations will be provided by Sunrise Medical summarizing the items selected. Please review this carefully.
 - Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

STANDARD OPTIONS

Quickie Xplore 2 STD Power Base, 6-wheel independent wheel suspension, an additional seat suspension system, tie-down hardware for strap systems, 8A offboard charger, 8" front casters, and standard auto style buckle seat positioning strap, 6.5 mph motors, Weight Capacity = 400 lbs.

QUICKIE Xplore 2 STD POWER WHEELCHAIR

XPLOREHD2... QUICKIE Xplore 2 HD Power Wheelchair \$9,275.00

USER WEIGHT PACKAGES

WCA34 400 lbs. Maximum user weight capacity STD

STANDARD FRAME COLORS

<input type="checkbox"/> FCBLK	Black	STD
<input type="checkbox"/> FCSBM	Matte Black	STD
<input type="checkbox"/> FCGMC	Gun Metal (Charcoal)	STD
<input type="checkbox"/> FCSVM	Silver Metallic	STD
<input type="checkbox"/> FCCBU	Candy Blue	STD
<input type="checkbox"/> FCCRE	Candy Red	STD
<input type="checkbox"/> FCMPL	Midnight Purple	STD
<input type="checkbox"/> FCRZB	Razzeberry	STD
<input type="checkbox"/> FCMWT	Matte White	STD

OPTIONAL FRAME COLORS (May Increase lead times)

<input type="checkbox"/> FCMYG	Military Green (New)	N/C
<input type="checkbox"/> FCBOR	Burnt Orange Metallic	N/C
<input type="checkbox"/> FCBAO	Black Opal Metallic	N/C
<input type="checkbox"/> FCOPM	Sparkle Orange	N/C
<input type="checkbox"/> FCMYT	Military Tan (New)	N/C
<input type="checkbox"/> FCPKP	Pearl Pink	N/C
<input type="checkbox"/> FCHDD 1	Desert Digital Camo Hydro Graphics	\$340.00
<input type="checkbox"/> FCHWL 1	Woodland Camo Hydro Graphics	\$340.00
<input type="checkbox"/> FCHCF 1	Carbon Fiber Hydro Graphics (Matte Black base)	\$340.00

Note 1 This options is included at no charge with the Hydro Graphics Frame Color option.

MATCHING ACCESSORIES TO THE FRAME COLOR

FC***-M Matching Accessories to the Frame Color (Wheels, front access door, headrest mount and back shell when selected) \$256.00

TRANSPORT OPTION

TTD4 2 4-Point wheelchair transport brackets STD
 WC19-ISO New WC-19 and ISO 7176-19 Occupied Transit Package N/C

To confirm if your overall wheelchair configuration is WC-19 and ISO 7176-19 compliant, please refer to this website: <http://wc-transportation-safety.umtri.umich.edu/crash-tested-product-lists/wheelchairs>
 Other restrictions may apply, please contact CSR for details. This option is only available on new power wheelchairs as an "on-chair" option. It is not available to install on a delivered wheelchair

Note 2 4-Point wheelchair transport brackets have not yet received WC-19 crash-testing approval.
 Use these transport brackets only to secure an unoccupied wheelchair during transport.

MOTOR PACKAGE

MP65 6.5 MPH 4-Pole Motor Package / with COG Release STD

ASSISTIVE DRIVING PACKAGE

INTELLIDR Intelli-Drive \$1,694.00

Note 3 Intelli-Drive technology provides increased tracking capabilities with enhanced traction and steering precision.
 This is ideal for those with marginal hand function or switch control.
 Note 4 Requires the R-net PM 120 Amp Expandable Controller (ECR120)

▼ BATTERY CHARGER			
<input type="checkbox"/>	BCOB	8A Offboard Battery Charger	N/C
<input type="checkbox"/>	BCLES8OPT	Lester 8 amp Offboard Charger	\$383.00

▼ BATTERY TRAY AND FRAME SIZE			
<input type="checkbox"/>	FS24B	Battery Tray for Group 24 Batteries (Total frame width = 25 3/8")	N/C
Please Select One Option: <input type="checkbox"/> Front Caster Arms <input type="checkbox"/> Flared out (STD) <input type="checkbox"/> Straight (Optional)			
Note 5 The Straight Front Arm option is not available with any of the Centermount Front Rigging Options due to Mechanical conflicts			

▼ BATTERIES			
<input type="checkbox"/>	BAT24PR	Gel Batteries-Group 24 (Each-On chair)	\$881.00

▼ CHOOSE TIRE COLOR (APPLIES TO ALL TIRES OR CASTERS)			
<input type="checkbox"/>	BLACK (New)	<input type="checkbox"/>	GREY

▼ DRIVE WHEEL			
<input type="checkbox"/>	DW14PT	14" Cast Wheel Pneumatic Tire	STD
<input type="checkbox"/>	DW14PI	14" Cast Wheel Pneumatic with Flat-Free Inserts	\$256.00
<input type="checkbox"/>	DW144PT	6 14" x 4" Extra Wide Cast Wheel Pneumatic Tire	\$279.00
<input type="checkbox"/>	DW144PI	6 14" x 4" Extra Wide Cast Wheel Pneumatic with Flat-Free Inserts	\$279.00
<input type="checkbox"/>	DW12PT	7 12 1/2" Cast Wheel Pneumatic Tire	N/C
<input type="checkbox"/>	DW12PI	7 12 1/2" Cast Wheel Pneumatic with Flat-Free Inserts	\$256.00
Note 6 The 14" x 4" Extra Wide Wheel increases the overall width of the Base by 2 1/2"			
The option to match the wheels to the shroud colour (FC...M) is not available with the 14" x 4" Extra Wide Wheels.			
The 14" x 4" option is not available with the VR2 electronics (ECV902A)			
Note 7 12 1/2" Drive wheels reduce overall seat to floor height by 1/4", lowers ground clearance and speed, but increases motor torque			

▼ FRONT CASTERS			
<input type="checkbox"/>	FC8PTR	8" Pneumatic Casters	STD
<input type="checkbox"/>	FC8PIR	8" Pneumatic Casters with Flat-Free Inserts	\$256.00
<input type="checkbox"/>	FC9PT	8 9" Pneumatic Casters	\$173.00
<input type="checkbox"/>	FC9PI	8 9" Pneumatic Casters with Flat-Free Inserts	\$204.00
Note 8 The 9" casters will increase the minimum seat to floor height by 1/2" and are NOT compatible with the 12 1/2" Drive wheel configuration (DW12PT and DW12PI)			
The 9" casters are NOT compatible with many front rigging options. Please see details in the front rigging section on pg #6.			

▼ SEAT TO FLOOR HEIGHT (Approximate measurement to the top of seat pan ± 1/4")					
► Options for a Seat Only					
<input type="checkbox"/>	STFH1525	15 1/2" (Possible Seat angle of 0° only)	<input type="checkbox"/>	STFH18	18" (Possible Seat angle of 0° and 3°)
<input type="checkbox"/>	STFH16	16" (Possible Seat angle of 0° only)	<input type="checkbox"/>	STFH19	19" (Possible Seat angle of 0°, 3° and 6°)
<input type="checkbox"/>	STFH17	17" (Possible Seat angle of 0° and 3°)	<input type="checkbox"/>	STFH20	20" (Possible Seat angle of 0°, 3°, 6°, 9°)
► Options for a Power Tilt or a STD Power Elevating Seat					
<input type="checkbox"/>	STFH1575	15 3/4" (Possible Seat angle of 0° only)	<input type="checkbox"/>	STFH18	18" (Possible Seat angle of 0° only)
<input type="checkbox"/>	STFH165	16 1/2" (Possible Seat angle of 0° only)	<input type="checkbox"/>	STFH185	18 1/2" (Possible Seat angle of 0° only)
<input type="checkbox"/>	STFH175	17 1/2" (Possible Seat angle of 0° only)	<input type="checkbox"/>	STFH19	19" (Possible Seat angle of 0° only)
► Options for a Combo - Power Tilt with a STD Power Elevating Seat					
<input type="checkbox"/>	STFH1875	18 3/4" (Possible Seat angle of 0° only)			
<input type="checkbox"/>	STFH195	19 1/2" (Possible Seat angle of 0° only)			
Note 9 Requires the 12 1/2" Drive wheels (DW12PT or DW12PI). This lowers ground clearance and speed, but increases motor torque.					
Note 10 These seat to floor height measurement are not available with the Heavy-Duty 11" Power Elevating Seat (SAPEXHD or ADDEPXHD)					
Note 11 The Heavy Duty Swing-Away footrests (SOAAF) are not available with this seat to floor height.					

▼ SEAT TO FLOOR HEIGHT SPECIFICATION (If not standard)			
<input type="checkbox"/>	SOMHSS	Specific seat to floor height requirement other than standard. Custom upcharges may apply. Required at: _____" (± 1/4" Range)	RFQ

▼ PRE-SET SEAT FRAME ANGLE			
<input type="checkbox"/>	SFE-0	0° Pre-set seat frame seat angle (Must be selected when using power options)	STD
<input type="checkbox"/>	SFE-3	3° Pre-set seat frame seat angle (Only available with a seat configuration without power options)	N/C
<input type="checkbox"/>	SFE-6	6° Pre-set seat frame seat angle (Only available with a seat configuration without power options)	N/C
<input type="checkbox"/>	SFE-9	9° Pre-set seat frame seat angle (Only available with a seat configuration without power options)	N/C

▼ LIGHTING PACKAGE			
<input type="checkbox"/>	EOLP	12 LED Light Package (Includes front beams, side flashers, rear tail light with flashers & brake light)	\$1,131.00
Note 12 Requires the upgraded Joystick with Light Operation or the OMNI (ECVJ6KL, ECELJL, ECRJCL, ECRJCTPL, ECRJC35B or EDIB). The Light Package is not compatible with the EPSEL option.			

▼ NON-EXPANDABLE ELECTRONICS - PG DRIVES VR2 CONTROLLER AND JOYSTICKS

<input type="checkbox"/>	ECV902A	VR2 90 Amp Drive Controller	STD
<input type="checkbox"/>	ECVJ4K	13 VR2 Standard 4-Key Joystick	STD
<input type="checkbox"/>	ECVJ6K	13 VR2 6-Key Joystick (Allows for a single power option to operate through this joystick)	N/C
<input type="checkbox"/>	ECVJ9KL	13 VR2 9-Key Joystick with Lighting Operation (Allows for a single power function to operate through this joystick)	\$168.00
Note 13 Requires the VR2 90A Drive Control (ECV902A)			
VR2 Electronics can only accommodate a single power seat function. The only choices are Power Tilt or Elevating seat			

▼ EXPANDABLE ELECTRONICS - PG DRIVES R-NET CONTROLLERS

<input type="checkbox"/>	ECEL90	R-net EL 90 Amp Expandable Controller	\$562.00
<input type="checkbox"/>	ECR120	R-net PM 120 Amp Expandable Controller	\$1,131.00
<input type="checkbox"/>	ECRNH	14 Harness for Expandable Controller	N/C
Note 14 Must be ordered with the expandable controller ECEL90 or ECR120 to interface with Joysticks or OMNI screen			

▼ EXPANDABLE ELECTRONICS - PG DRIVES R-NET JOYSTICKS

<input type="checkbox"/>	ECO	Omit Joystick (For use with an Omni Display)	N/C
<input type="checkbox"/>	ECELJ	R-net (LED) Joystick	STD
<input type="checkbox"/>	ECRJCL	R-net Joystick with Color screen	\$168.00
▶ JOYSTICKS WITH LIGHTING OPERATION			
<input type="checkbox"/>	ECELJL	R-net (LED) Joystick with Lighting Operation	\$168.00
<input type="checkbox"/>	ECRJCL	R-net Joystick with Color screen and Lighting Operation	\$222.00
<input type="checkbox"/>	ECRJCTPL	R-net Joystick with Color screen w/ On & Off/ Mode switch & speed pt and Lighting Operation.	
▶ ADVANCED JOYSTICK			
<input type="checkbox"/>	ECRJ35	15 R-net Advanced Joystick (3.5" Screen). Lighting Operation, paddle switches, shortcut keys and stereo jacks	\$1,202.00
<input type="checkbox"/>	ECRJ35B	15 R-net Advanced Bluetooth Joystick (3.5" Screen). Lighting Operation, paddle switches, shortcut keys and stereo jacks	\$1,202.00
Note 15 This option requires the expandable controller ECEL90 or ECR120			

▼ JOYSTICK MOUNT

<input type="checkbox"/>	JSMF	Fixed Joystick Mount	N/C
<input type="checkbox"/>	JMSA	Swing Away Joystick Mount - Height adjustable, Angle adjustable and Multi-Axis	\$412.00
		<input type="checkbox"/> Set-up with standard installation to swing-away sideways	<input type="checkbox"/> Set-up at 35° to swing-away sideways and downwards
<input type="checkbox"/>	JSMCM	Midline Mounted Swing Away Joystick Mount	\$412.00

▼ JOYSTICK MOUNT SUPPORT

<input type="checkbox"/>	JSMFH	Joystick mount support Bracket	STD
--------------------------	-------	--------------------------------	-----

▼ JOYSTICK HAND MOUNT

<input type="checkbox"/>	HML	Install Joystick on Left hand side	<input type="checkbox"/>	HMR	Install Joystick on Right hand side
--------------------------	-----	------------------------------------	--------------------------	-----	-------------------------------------

▼ JOYSTICK HANDLE

<input type="checkbox"/>	JSH	Standard Joystick Gimble	STD
<input type="checkbox"/>	JBAT	Stick Handle 3" in Length	\$101.00
<input type="checkbox"/>	JBH	Ball Handle	\$101.00
<input type="checkbox"/>	JFBH	Foam Ball Handle	\$101.00
<input type="checkbox"/>	JTH	T-Handle	\$101.00
<input type="checkbox"/>	JMUSH	New 2.5" Mushroom Handle	\$122.00
<input type="checkbox"/>	JGP	Goalpost Handle	\$101.00

▼ ATTENDANT DRIVE CONTROL OPTION

<input type="checkbox"/>	EADC	Attendant Hand Control for R-net EL & PM Controls	\$1,131.00
--------------------------	------	---	------------

▼ ATTENDANT JOYSTICK HAND MOUNT

<input type="checkbox"/>	EADCL	Install Attendant Joystick on Rear Left hand side	<input type="checkbox"/>	EADCR	Install Attendant Joystick on Rear Right hand side
--------------------------	-------	---	--------------------------	-------	--

▼ ENHANCED DISPLAY & DIGITAL INTERFACE BOX (Included with specialty controls)

<input type="checkbox"/>	EDIB	16 OMNI (Enhanced Display & Digital Interface Box with Infrared)	\$2,157.00
Note 16 Allows the switch operation of alternative devices through driver controls. Requires one of the R-net fully Expandable Controllers (ECEL90 or ECR120)			

▼ ENHANCED DISPLAY MOUNT

<input type="checkbox"/>	EDSAM	User Swing Away Mount	<input type="checkbox"/>	EDML (Install Display on LH side)	<input type="checkbox"/>	EDMR (Install Display on RH side)	N/C
--------------------------	-------	-----------------------	--------------------------	-----------------------------------	--------------------------	-----------------------------------	-----

▼ ENVIRONMENTAL CONTROL MODULE

<input type="checkbox"/> ECU	17	Environment control module - (Output Module)	\$1,131.00
	Note 17	Requires one of the fully expandable R-net controllers (ECEL90 or ECR120)	

▼ BLUETOOTH MOUSE EMULATOR (Includes USB Dongle)

<input type="checkbox"/> EDBT	18	Bluetooth Mouse Emulator (USB Dongle included)	\$1,131.00
	Note 18	Requires one of the fully expandable R-net controllers (ECEL90 or ECR120) . Not required if the ECRJC35B joystick has been selected Not limited to only one selection. Any two can be installed on the same chair, but only one can be used at a time.	
		<input type="checkbox"/> Windows and MAC OS	<input type="checkbox"/> Android Devices
			<input type="checkbox"/> iPhone, iPad and iPod

▼ POWER OPTIONS & SEATING

<input type="checkbox"/> SAABA		Amyseat Rehab Seat - Seat Frame only	\$881.00
<input type="checkbox"/> SAPTXHD	19	Heavy-Duty CG Power Tilt	\$8,235.00
<input type="checkbox"/> SAPRXHD	19, 20	Heavy-Duty Power Recline with Integrated Shear Reduction (ISR)	\$9,940.00
<input type="checkbox"/> SAPTPRXHD	19, 20	Heavy-Duty CG Power Tilt and Power Recline with Integrated Shear Reduction (ISR)	\$14,775.00
<input type="checkbox"/> SAPHXHD	19	Heavy-Duty 11" Power Elevating Seat Module for transfer assist	\$5,679.00
<input type="checkbox"/> ADDPESXHD		Module - Add a Heavy-Duty 11" Power Elevating Seat to either SAPTXHD, SAPRXHD or SAPTPRXHD	\$4,542.00
	Note 19	These options include the fully adjustable Amyseat (SAABA)	
	Note 20	It is recommended to use the Contour Back Support option with the power reclines (STDSB or STDSB/CUSTOM)	

▼ ELECTRONICS - SINGLE OR DUAL (1 OR 2 ACTUATORS) POWER SEAT FUNCTION CONTROLLERS

<input type="checkbox"/> EPSV	21, 22	Tilt or Elevating seat activation through 6-key VR2 Hand Control (ECVJ6K or ECVJ9KL)	N/C
<input type="checkbox"/> EPSEL	21, 22	Tilt or Elevating seat activation through R-NET 90 Amp Expandable Controller (ECEL90)	N/C
<input type="checkbox"/> SAC		Single (1) Power Seat Function Controller to Interface through R-net Electronics	N/C
<input type="checkbox"/> DAC		Dual (2) Power Seat Function Controller to Interface through R-net Electronics	\$562.00
With the SAC or DAC, please select one activation and position from the list below			
<input type="checkbox"/> SWTJSS	23, 24	Single actuator activation through R-net Joystick or Specialty Controls EPSR2A	\$1,131.00
<input type="checkbox"/> SWTJSD	23, 25	Dual actuator activation through R-net Joystick or Specialty Controls EPSRA	\$1,131.00
<input type="checkbox"/> SWDUALACT		Dual Activation (Allows to operate all seat functions with external switches and through the joystick)	\$350.00
<input type="checkbox"/> SWTOGGLE	24 or 25	Single Toggle Switch (Single Mode: DAC - Dual Mode: SAC)	<input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation
<input type="checkbox"/> SWPUSH1	24	Single Push Button Switch (Single Mode: SAC)	<input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation
<input type="checkbox"/> SWPUSH2	24 or 25	Dual Push Button Switch (Single Mode: DAC - Dual Mode: SAC)	<input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation
	Note 21	Allows for the activation of a single power option through the VR2 or R-net 90 Amp Controllers.	
		Only the Power Tilt or Elevating seat options are available with this feature.	
	Note 22	These options do not require the SAC controller as they connect via the chair's power module directly	
	Note 23	Requires the R-net 90 Amp (ECEL90) or the 120 Amp (ECR120) Fully Expandable Controllers	
	Note 24	Requires the Single (1 Actuator) Power Seat Function (SAC)	
	Note 25	Requires the Dual (2 Actuators) Power Seat Function (DAC)	

▼ ELECTRONICS - MULTIPLE (3 to 6 ACTUATORS) POWER SEAT FUNCTION CONTROLLER

<input type="checkbox"/> MAC		Multiple (3 to 6 Actuators) Power Seat Function Controller to Interface through R-net Electronics	\$2,720.00
With the MAC, please select one activation and position from the list below			
<input type="checkbox"/> SWTJSM	26, 27	Multiple Actuator activation through R-net Joystick or Specialty Controls MAC-EPSR	STD
<input type="checkbox"/> SWDUALACT		Dual Activation (Allows to operate all seat functions with external switches and through the joystick)	\$350.00
<input type="checkbox"/> SWMULTI	27	4-Way Toggle Switch (Up to 4 Activations)	<input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation
<input type="checkbox"/> SWMULTIT	27	Triple Toggle Switch (Up to 5 Activations)	<input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation
<input type="checkbox"/> SWPUSH4	27	4-Way Push Button Switch	<input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation
<input type="checkbox"/> SWPUSH5	27	5-Way Push Button Switch	<input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation
	Note 26	Requires the Multiple (3 to 6 Actuators) Power Seat Function Controller (MAC)	
	Note 27	Requires the R-net 90 Amp (ECEL90) or 120 Amp (ECR120) Fully Expandable Controller, and the multiple Actuator Controller (MAC)	

▼ SWITCH ACTIVATION HAND MOUNT

Switch Installation:	<input type="checkbox"/> TSML (Install Switch on LH side)	<input type="checkbox"/> TSMR (Install Switch on RH side)	N/C
----------------------	---	---	-----

▼ LARGE ADJUSTABLE SEAT WIDTH RANGE BETWEEN 19" AND 26"			
<input type="checkbox"/> SOWN19-24-22	Width adjustment from 19" to 24" with a maximum Depth 22"		\$490.00
<input type="checkbox"/> SOWN24-26-22	Width adjustment from 24" to 26" with a maximum Depth 22"		\$490.00
<input type="checkbox"/> 19" (SOWN19L)	<input type="checkbox"/> 22" (SOWN22)	<input type="checkbox"/> 24" (SOWN24HDL)	
<input type="checkbox"/> 20" (SOWN20L)	<input type="checkbox"/> 23" (SOWN23)	<input type="checkbox"/> 25" (SOWN25HDL)	
<input type="checkbox"/> 21" (SOWN21)	<input type="checkbox"/> 24" (SOWN24)	<input type="checkbox"/> 26" (SOWN26HDL)	
▼ ADJUSTABLE SEAT DEPTH RANGE BETWEEN 17" AND 22"			
<input type="checkbox"/> 17" (SODN17)	<input type="checkbox"/> 18" (SODN18)		STD
<input type="checkbox"/> 19" (SODN19)	<input type="checkbox"/> 20" (SODN20)	<input type="checkbox"/> 21" (SODN21)	<input type="checkbox"/> 22" (SODN22) \$334.00
▼ BACK CANE HEIGHT (The back cane height measurement is taken from the seat pan to the top of the push handles)			
<input type="checkbox"/> 16" (SOBCH16)	<input type="checkbox"/> 22" (SOBCH22) ²⁸		
<input type="checkbox"/> 18" (SOBCH18)	<input type="checkbox"/> 24" (SOBCH24)		
<input type="checkbox"/> 20" (SOBCH20)	<input type="checkbox"/> 26" (SOBCH26)		
Note 28 22" is the minimum back cane height required with a power or manual recline.			
▼ PRE-SET SEAT ANGLE (For Non-Reclining Backs)			
<input type="checkbox"/> SOBAA	Please Select one:	<input type="checkbox"/> 82° <input type="checkbox"/> 90° <input type="checkbox"/> 98° <input type="checkbox"/> 106°	N/C
▼ BACK SUPPORT			
<input type="checkbox"/> STDSB	²⁹ Contour back support 15" to 20" Wide with a maximum shell height of 22"		\$624.00
<input type="checkbox"/> STDSB/CUSTOM	Contour back support Non Standard Sizes up to 24" wide	Width: _____ Depth: _____	\$1,333.00
Note 29 It is recommended to use the Back Support option with the power recline.			
▼ BACK SUPPORT TOTAL HEIGHT (MEASURED FROM THE SEAT PAN TO THE TOP OF THE BACK)			
<input type="checkbox"/> 16" (SOBH16)	<input type="checkbox"/> 24" (SOBH24)		
<input type="checkbox"/> 18" (SOBH18)	<input type="checkbox"/> 26" (SOBH26)		
<input type="checkbox"/> 20" (SOBH20)	<input type="checkbox"/> 28" (SOBH28)		
<input type="checkbox"/> 22" (SOBH22) ³⁰			
Note 30 22" is the minimum back height required with a power recline.			
Note 31 The space between the bottom of the back and the seat pan is 4" with a non-reclining back and 6" with a reclining back.			
▼ SEAT CUSHIONS			
► Contour Seat Cushion			
<input type="checkbox"/> SCPSC	Contour Seat Cushion (15" to 20" Wide, with maximum 22" seat depth)		\$624.00
<input type="checkbox"/> SCPSC/CUSTOM	Contour Seat Cushion (Non Standard Sizes) Width: _____ Depth: _____		\$1,359.00
▼ CANTILEVER - HEIGHT ADJUSTABLE ARMRESTS - OPENS WITH A 5° ANGLE (For Non-Reclining Backs - Choose 2)			
<input type="checkbox"/> SOAFBFLH	Flip-Back Cantilever - Full Length Left Side (Height Adjustable from 7" to 15" with 1" increments)		STD
<input type="checkbox"/> SOAFBFRH	Flip-Back Cantilever - Full Length Right Side (Height Adjustable from 7" to 15" with 1" increments)		STD
<input type="checkbox"/> SOAFBDLH	Flip-Back Cantilever - Desk Length Left Side (Height Adjustable from 7" to 15" with 1" increments)		N/C
<input type="checkbox"/> SOAFBDRH	Flip-Back Cantilever - Desk Length Right Side (Height Adjustable from 7" to 15" with 1" increments)		N/C
▼ DUAL POST - HEIGHT ADJUSTABLE ARMRESTS - INCLUDES TRANSFER GRIPS (For Non-Reclining Backs - Choose 2)			
<input type="checkbox"/> SOADPFLH	Dual Post Height Adjustable - Full Length Left Side (Height Adjustable from 10 ¼" - 15 ¼")		\$202.00
<input type="checkbox"/> SOADPFRH	Dual Post Height Adjustable - Full Length Right Side (Height Adjustable from 10 ¼" - 15 ¼")		\$202.00
<input type="checkbox"/> SOADPDHL	Dual Post Height Adjustable - Desk Length Left Side (Height Adjustable from 10 ¼" - 15 ¼")		\$202.00
<input type="checkbox"/> SOADPDRH	Dual Post Height Adjustable - Desk Length Right Side (Height Adjustable from 10 ¼" - 15 ¼")		\$202.00
▼ RECLINING FLIP-UP, HEIGHT ADJUSTABLE ARMRESTS (For Reclining Backs - Choose 2)			
<input type="checkbox"/> SOANIFLH	Reclining armrest - Full Length Left Side (Height Adjustable from 10 ½" to 14" with ½" increments)		N/C
<input type="checkbox"/> SOANIFRH	Reclining armrest - Full Length Right Side (Height Adjustable from 10 ½" to 14" with ½" increments)		N/C
<input type="checkbox"/> SOANIDLH	Reclining armrest - Desk Length Left Side (Height Adjustable from 10 ½" to 14" with ½" increments)		N/C
<input type="checkbox"/> SOANIDRH	Reclining armrest - Desk Length Right Side (Height Adjustable from 10 ½" to 14" with ½" increments)		N/C

▼ ARMREST HEIGHT ADJUSTMENT AT FACTORY

Please set the Armrest Height at: _____ " (The Armrest Height is taken from the seat pan to the top of the armpad. Not required for Dual Post)

▼ ARMPADS

► Foam Type Selection

<input type="checkbox"/>	SOFAPD	Desk Length Standard (10")	Select Side: <input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	STD
<input type="checkbox"/>	SOFAPF	Full Length Standard (14")	Select Side: <input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	STD
<input type="checkbox"/>	SOWAPD	33 Desk Length Waterfall (10")	Select Side: <input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$34.00
<input type="checkbox"/>	SOWAPF	33 Full Length Waterfall (14")	Select Side: <input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$34.00

► Hybrid Gel and Foam Type Selection (The Gel always stays soft)

<input type="checkbox"/>	SOGELD	Desk Length Hybrid Gel and Foam Flat armpad (10" x 2.5")	Select Side: <input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$34.00
<input type="checkbox"/>	SOGELF	Full Length Hybrid Gel and Foam Flat armpad (14" x 2.5")	Select Side: <input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$34.00
<input type="checkbox"/>	SOWGELD	33 Desk Length Hybrid Gel and Foam Waterfall armpad (10")	Select Side: <input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$54.00
<input type="checkbox"/>	SOWGELF	33 Full Length Hybrid Gel and Foam Waterfall armpad (14")	Select Side: <input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$54.00

► Specialty Type Selection

<input type="checkbox"/>	SOCAR	34 Arm Trough with Hand Support (18" x 4 1/2")	Select Side: <input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	\$222.00
	Note 33	The Waterfall armpads may interfere with flipping back the armrests when combined with Lateral Trunk supports			
	Note 34	Arm Trough Armpads may interfere with Armrests Flipping back.			

▼ FOOTRESTS

<input type="checkbox"/>	SOAAF	38 Adult Heavy Duty Swing-Away footrests	Set at: <input type="checkbox"/> 60°	<input type="checkbox"/> 70° ³⁹	\$390.00
			Select Footrest Range: <input type="checkbox"/> SOFLHD12 (12 1/2" - 15 1/2")	<input type="checkbox"/> SOFLHD14 (14 1/2" - 19 1/2")	
<input type="checkbox"/>	ADDPLHD	35, 38 Adult Power Elevating & Articulating Legrests	Please Select Function: <input type="checkbox"/> Individual ³⁶	<input type="checkbox"/> Combined ³⁷	\$3,399.00
			Select Footrest Range: <input type="checkbox"/> SOFLPW12 (12 1/2" - 15 1/2")	<input type="checkbox"/> SOFLPW14 (14 1/2" - 19 1/2")	

38 Please set the Footrest Length at: _____ "

Note 35 This options includes calf pads
Note 36 Individual legs count as two power functins.
Note 37 Combined legs count as one power function.
Note 38 The Footrest Length is taken from the top of the seat pan to the top of the footplate.
Note 39 The 70 degree footrest angle is not compatible with the 9" front casters (FC9PT and FC9PI)
 The 70 degree footrest angle may interfere with the casters when selecting the straight front caster arms option.

▼ FOOTPLATES & ACCESSORIES

<input type="checkbox"/>	SOSFAFHD	Flip-up Footplates		STD
<input type="checkbox"/>	SOAFPHD	Angle Adjustable & Flip-up Footplates - Each		\$145.00
<input type="checkbox"/>	SOAHL	Heel Loops - Each		\$44.00
<input type="checkbox"/>	SOACPFRR	Calf Pads for Adult Heavy-Duty Swing-Away footrests (SOAAF) with adjustable hardware - Each		\$145.00

▼ CENTER MOUNTS

<input type="checkbox"/>	SOAAF90	42, 43, 44, 45, 46, 47 90° Fixed Center Mount Interface Bracket (Offers Manual adjustability from 95° to 70°)	Select Foot Platform Range: <input type="checkbox"/> SOFBA12 (12 1/2" - 16 1/2")	<input type="checkbox"/> SOFBA16 (16 1/2" - 19 1/2")	\$256.00
<input type="checkbox"/>	ADDPFHD	40, 41, 42, 44, 45, 46 90° Power Center Mount Elevating & Articulating Foot Platform	Select Foot Platform Range: <input type="checkbox"/> SOFBPW12 (12 1/2"- 16 1/2")	<input type="checkbox"/> SOFBPW16 (16 1/2"- 20 1/2")	\$3,399.00

45 Please set the Foot Platform Length at: _____ "

Note 40 The power center mount counts as one power function
Note 41 This option includes calf pads
Note 42 Links and attaches the one piece center mount foot platform (SOFUFP) or the split footplates (SOAFPF) to the seat
Note 43 With the split footplates option (SOAFPF), this bracket offers independent height adjustability for each footplate.
Note 44 Please select the one piece center mount foot platform (SOFUFP) or the split footplates (SOAFPF)
Note 45 The Foot Platform length is taken from the top of the seat pan to the top of the foot platform or split footplates.
Note 46 Center Mount Front Riggings are not available with the Straight Front Arm option due to Mechanical conflicts
Note 47 This option is not compatible with the 9" Caster option

▼ FOOT PLATFORMS & ACCESSORIES

<input type="checkbox"/>	SOFUFP	Flip-up one piece Foot Platform		N/C
<input type="checkbox"/>	SOAFPF	Angle Adjustable & Flip-up Split Footplates for center mount - Each		\$145.00
<input type="checkbox"/>	SOACPCM	Calf Pads for 90° Fixed Center Mount (SOAAF90) with adjustable hardware - Each		\$145.00

▼ SEAT POSITIONING (NOT COMPATIBLE WITH THE WC19 ----- ISO OPTION)			
<input type="checkbox"/>	SOAUTO80	80" Seat Positioning Strap	\$83.00
▼ HEADREST			
<input type="checkbox"/>	SONUAHR	Using headrest mounting hardware other than SOMAH located below. Specify: _____	N/C
<input type="checkbox"/>	SOFHRS	Small Curved Headrest Pad 4" x 9"	\$228.00
<input type="checkbox"/>	SOFHRM	Medium Curved Headrest Pad 5" x 10"	\$228.00
<input type="checkbox"/>	SOFHRL	Large Curved Headrest Pad 6" x 11"	\$228.00
<input type="checkbox"/>	SOMAH	Multi-axis Removable Mounting Hardware	\$183.00
<input type="checkbox"/>	SOHRBC	Bracket to Mount on the Contour Back Support (Without Recline)	\$111.00
<input type="checkbox"/>	SOHRPR	Bracket to Mount on the Seat Frame (For Reclining Backs)	\$111.00
* Please select one of the headrest options with the mounting hardware and specify where the headrest will be installed			
▼ LATERAL TRUNK SUPPORTS			
<input type="checkbox"/>	SOSLS3-C	Adjustable Laterals with Small Curved Pads (Each) - 3" x 5"	\$675.00
<input type="checkbox"/>	SOSLM-C	Adjustable Laterals with Medium Curved Pads (Each) - 4" x 6"	\$675.00
<input type="checkbox"/>	SOSLL5-C	Adjustable Laterals with Large Curved Pads (Each) - 5" x 7"	\$675.00
▼ HIP & THIGH SUPPORTS			
<input type="checkbox"/>	SOHTSS10	Adjustable Hip & Thigh Support Pads (Each) - Small 10" x 4"	\$256.00
<input type="checkbox"/>	SOHTSM12	Adjustable Hip & Thigh Support Pads (Each) - Medium 12" x 5"	\$256.00
<input type="checkbox"/>	SOHTSL14	Adjustable Hip & Thigh Support Pads (Each) - Large 14" x 6"	\$256.00
<input type="checkbox"/>	SOHTSARH	Adjustable and Removable Hip & Thigh Support Hardware - Each	\$485.00
** Please select one of the hip & thigh support pads and the adjustable & removable hardware.			
▼ ACCESSORIES			
<input type="checkbox"/>	SOCH	Cup Holder (Self-Leveling)	\$101.00
<input type="checkbox"/>	SOPHP	Ram Mounts X-Grip Cell Phone Holder	\$173.00
<input type="checkbox"/>	SOTHS	Ram Mounts X-Grip Small Tablet Holder	\$196.00
<input type="checkbox"/>	SOTHL	Ram Mounts X-Grip Large Tablet Holder with Swivel	\$295.00
▼ REFLECTOR KIT			
<input type="checkbox"/>	REFLECT	Red & White Reflector Kit Stickers (They will not be installed - Included in a bag with the owner's manual)	N/C
▼ SWITCHES AND OPTIONS			
<input type="checkbox"/>	ASCPCP	PGDT R-NET Dongle and PC programmer kit	\$652.00
<input type="checkbox"/>	ASCPVR2	PGDT VR2 Dongle and PC programmer kit	\$652.00
<input type="checkbox"/>	SSBS-STD	Buddy Button Switch	\$122.00
<input type="checkbox"/>	SSES-STD	Egg Switch-Black	\$122.00
<input type="checkbox"/>	SSML	Ultra Light Switch	\$168.00

▼ SPECIALTY DRIVER CONTROL OPTIONS

▼ Therafin Sip n' Puff Kit

<input type="checkbox"/>	SCSNPK	48, 49	Therafin Sip n' Puff Kit				\$ 1,701.00	
<input type="checkbox"/>	SCSNPBTk	49	Therafin Sip n' Puff Breath Tube Kit				\$ 552.00	
<input type="checkbox"/>	SCSRMSNP	49	Stealth Rock-n-Lock Swingaway Arm for Sip-N-Puff	<input type="checkbox"/>	SWL (Swings on Left)	<input type="checkbox"/>	SWR (Swings on Right)	\$ 307.00
<input type="checkbox"/>	SCSRM	49	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/>	SWL (Swings on Left)	<input type="checkbox"/>	SWR (Swings on Right)	\$552.00
<input type="checkbox"/>	SSES-STD	49, 51	Egg Switch for Changing Modes - Black	<input type="checkbox"/>	MLH (Mount on Left)	<input type="checkbox"/>	MRH (Mount on Right)	\$122.00
<input type="checkbox"/>	SSBS-STD	49	Buddy Button Switch (Stop Switch) - Red				\$122.00	
<input type="checkbox"/>	SOMAH	49	Multi-axis Headrest Removable Mounting Hardware				\$183.00	
<input type="checkbox"/>	SOFHRS	50	Small Curved Headrest Pad 4" x 9"				\$228.00	
<input type="checkbox"/>	SOFHRM	50	Medium Curved Headrest Pad 5" x 10"				\$228.00	
<input type="checkbox"/>	SOFHRL	50	Large Curved Headrest Pad 6" x 11"				\$228.00	

▼ Switch-It Proximity Head Array

<input type="checkbox"/>	SCSIHAA	48, 49	Adult - 3 Switch Proximity Head Array with 1/8" Mono Jack for mode/reverse				\$5,679.00	
<input type="checkbox"/>	SCSRM	49	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/>	SWL (Swings on Left)	<input type="checkbox"/>	SWR (Swings on Right)	\$552.00
<input type="checkbox"/>	SSES-STD	49, 51	Egg Switch for Changing Modes - Black	<input type="checkbox"/>	MLH (Mount on Left)	<input type="checkbox"/>	MRH (Mount on Right)	\$122.00
<input type="checkbox"/>	SSBS-STD	49	Buddy Button Switch (Stop Switch) - Red				\$122.00	
<input type="checkbox"/>	SOMAH	49	Multi-axis Headrest Removable Mounting Hardware				\$183.00	

▼ Proportional RIM Head Control Package

<input type="checkbox"/>	SCPHC	48, 49	Proportional RIM Head Control Package				\$2,781.00	
<input type="checkbox"/>	SCSRM	49	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/>	SWL (Swings on Left)	<input type="checkbox"/>	SWR (Swings on Right)	\$552.00
<input type="checkbox"/>	SSES-STD	49, 51	Egg Switch for Changing Modes - Black	<input type="checkbox"/>	MLH (Mount on Left)	<input type="checkbox"/>	MRH (Mount on Right)	\$122.00
<input type="checkbox"/>	SOMAH	49	Multi-axis Headrest Removable Mounting Hardware				\$183.00	

▼ Switch-It MicroPilot Light Touch Proportional Mini joystick for chin / lip control

<input type="checkbox"/>	SCSIMP	48, 49	MicroPilot Light touch mini Proportional Joystick				\$4,453.00	
<input type="checkbox"/>	SCSRMMP	49	Stealth Rock-n-Lock Swingaway Arm for Chin MP	<input type="checkbox"/>	SWL (Swings on Left)	<input type="checkbox"/>	SWR (Swings on Right)	\$552.00
<input type="checkbox"/>	SCSRM	49	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/>	SWL (Swings on Left)	<input type="checkbox"/>	SWR (Swings on Right)	\$552.00
<input type="checkbox"/>	SSES-STD	49, 51	Egg Switch for Changing Modes - Black	<input type="checkbox"/>	MLH (Mount on Left)	<input type="checkbox"/>	MRH (Mount on Right)	\$122.00
<input type="checkbox"/>	JCC		Chin Cup				\$57.00	
<input type="checkbox"/>	SOMAH	49	Multi-axis Headrest Removable Mounting Hardware				\$183.00	
<input type="checkbox"/>	SOFHRS	50	Small Curved Headrest Pad 4" x 9"				\$228.00	
<input type="checkbox"/>	SOFHRM	50	Medium Curved Headrest Pad 5" x 10"				\$228.00	
<input type="checkbox"/>	SOFHRL	50	Large Curved Headrest Pad 6" x 11"				\$228.00	

▼ Proportional Compact or Mini joystick for Midline Mount

<input type="checkbox"/>	SCPGMJS	48, 49	PG Drives JC200 Proportional Compact Remote Joystick for Midline Mount				\$1,059.00	
<input type="checkbox"/>	SCGMJSM	52	Gatlin Midline Swing-away Mount for JC200	<input type="checkbox"/>	MLH (Mount on Left)	<input type="checkbox"/>	MRH (Mount on Right)	\$435.00
<input type="checkbox"/>	SCSIMP	48, 49	Switch-It MicroPilot Light touch mini Proportional Joystick				\$4,453.00	
<input type="checkbox"/>	SCGMPM	53	Gatlin Midline Swing-away Mount for MicroPilot	<input type="checkbox"/>	MLH (Mount on Left)	<input type="checkbox"/>	MRH (Mount on Right)	\$474.00
<input type="checkbox"/>	SCBTM	54	Bullet Hand Tray for MicroPilot	<input type="checkbox"/>	MLH (Mount on Left)	<input type="checkbox"/>	MRH (Mount on Right)	\$502.00

▼ Proportional Compact or Mini joystick for Armrest Mount

<input type="checkbox"/>	SCPGMJS	48, 49	PG Drives JC200 Proportional Compact Remote Joystick for Armrest Mount				\$1,059.00	
<input type="checkbox"/>	SCARMJSM	52	Swing-away and Multi-Axis armrest mount	<input type="checkbox"/>	MLH (Mount on Left)	<input type="checkbox"/>	MRH (Mount on Right)	\$435.00
<input type="checkbox"/>	SCSIMP	48, 49	Switch-It MicroPilot Light touch mini Proportional Joystick				\$4,453.00	
<input type="checkbox"/>	SCARMPM	53	Swing-away and Multi-Axis armrest mount	<input type="checkbox"/>	MLH (Mount on Left)	<input type="checkbox"/>	MRH (Mount on Right)	\$474.00
<input type="checkbox"/>	SCBTM	54	Bullet Hand Tray for MicroPilot	<input type="checkbox"/>	MLH (Mount on Left)	<input type="checkbox"/>	MRH (Mount on Right)	\$502.00

Note 48 Requires the OMNI Enhanced Display and Digital Interface Box (EDIB)

Note 49 This option must be selected for the kit

Note 50 Choose one (1) Headrest Pad

Note 51 The mode switch type can be changed from this selection on Page #6

Note 52 Must be selected for use with the PG Drives JC200 Mini Joystick

Note 53 Must be selected for use with the Switch-It MicroPilot

Note 54 This option is only compatible with the Switch-It MicroPilot

▼ CUSTOM SECTION

Note 55 Lead times will increase depending on the custom order pipeline at the time the purchase order is received. Please call for Quote and specifications for these options.

▼ O² Holder

<input type="checkbox"/> SUPOXY	Custom Non-Articulating Oxygen Tank Holder	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	\$435.00
<input type="checkbox"/> SUPOXYA	Custom Articulating Oxygen Tank Holder	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	\$2,040.00
	<input type="checkbox"/> D-Size Oxygen Tank (* 22" Tall x 4.3" Diameter) * The height includes the valve			
	<input type="checkbox"/> E-Size Oxygen Tank (* 31" Tall x 4.3" Diameter) * The height includes the valve			

Note 56 Not available with a power tilt and/or power recline

Note 57 Required with a power tilt and/or power recline

▼ VENT TRAY (Please complete in Full)

<input type="checkbox"/> SOVENT	Custom Vent Tray			
<input type="checkbox"/> SOVENT-A 58	Articulating Vent Box (Required with Power Tilt and/or Recline)	Model:		\$1,957.00
<input type="checkbox"/> SOVENT-F 58	Fixed Vent Box	Model:		\$491.00
<input type="checkbox"/> SOVTBE 58, 59	Extended Battery Compartment (K0108)			\$600.00

- The Vent is being powered by 1 U1 Battery and the chair by 2 Group 22NF batteries
- The Vent is being powered by 1 U1 Battery and the chair by 2 Group 24 batteries
- The Vent is being powered by 1 Group 22NF Battery and the chair Powered by 2 Group 22NF Batteries
- The Vent is being powered by 1 Group 24 Battery and the chair Powered by 2 Group 22NF Batteries

Note 58 Sunrise Medical does not supply the battery, charger or the wiring for the Vent.

These items are normally provided by the Vent Manufacturer.

Note 59 This option is not available with the Straight Front Caster Arms.

▼ OTHER CUSTOM OPTIONS (Request For Quote)

Please specify details and specifications for these options.

<input type="checkbox"/> CUSTOM 60	Power Anterior Tilt			RFQ
<input type="checkbox"/> CUSTOM 60, 61	Power Flip-up Foot Platform (Please call CSR for availability with configuration and user weight capacity)			RFQ
<input type="checkbox"/> CUSTOM 60	Independent Depth Adjustability built into the Footrest Receiver that is attached to the Seat Pan			RFQ
<input type="checkbox"/> CUSTOM 60	Endomorphic Back Canes (Pair) <input type="checkbox"/> 1" per side towards the inside = 2" in <input type="checkbox"/> 2" per side towards the inside = 4" in			RFQ
<input type="checkbox"/> CUSTOM	Transfer Handles (Pair)			RFQ
<input type="checkbox"/> CUSTOM	Manual Wheel Locks			RFQ
<input type="checkbox"/> CUSTOM	Non Standard Frame Colors (Not all custom colors are available, please contact CSR for availability)			RFQ
<input type="checkbox"/> CUSTOM # 1 60	Other custom options. Specify in notes section below.			RFQ
<input type="checkbox"/> CUSTOM # 2 60	Other custom options. Specify in notes section below.			RFQ

Note 60 Please specify details and specifications for these options

Note 61 This option is not available with the Straight Front Caster Arms

▼ NOTES / SPECIAL REQUESTS / CUSTOM SPECIFICATIONS
