



AADL #W848

**Quickie M6™  
EIF1**

For Alberta Aids to Daily Living only



July 01, 2015

ACCOUNT #: \_\_\_\_\_ SHIP TO: \_\_\_\_\_  
 DATE: \_\_\_\_\_ NAME: \_\_\_\_\_  
 P.O. #: \_\_\_\_\_ ADDRESS: \_\_\_\_\_  
 CITY: \_\_\_\_\_ PROVINCE: \_\_\_\_\_  
 POSTAL CODE: \_\_\_\_\_ PHONE #: \_\_\_\_\_  
 BUYER: \_\_\_\_\_ MARKED FOR: \_\_\_\_\_



Product Info ►

**1. FRAME STYLE & SETUP****Base Model**

<input type="checkbox"/>	115WM21	Quickie M6	\$4,065.00
--------------------------	---------	------------	------------

**Frame Color**

<input type="checkbox"/>	115C2	Black	NC
--------------------------	-------	-------	----

<input type="checkbox"/>	115C38	Blue Opal	NC
--------------------------	--------	-----------	----

**Frame Width**

<input type="checkbox"/>	115W13	22"	NC
<input type="checkbox"/>	115W14	24"	NC
<input type="checkbox"/>	115W15	26"	NC
<input type="checkbox"/>	115W16	28"	NC
<input type="checkbox"/>	115W17	30"	NC

**Frame Length**

<input type="checkbox"/>	115FL9	Regular (18"-19" depth)	NC
<input type="checkbox"/>	115FL18	Long (20"-22" depth)	NC

**Seat Height**

<input type="checkbox"/>	115HT65	17" <i>Includes low adjustable axle plate and shorter fork. Only available with 5" caster</i>	\$474.00
<input type="checkbox"/>	115HT11	18.5" <i>Only available w/ 5" casters</i>	NC
<input type="checkbox"/>	115HT14	20"	NC

**Seat Depth**

<input type="checkbox"/>	115HT11	18"	NC
<input type="checkbox"/>	115HT12	19"	NC
<input type="checkbox"/>	115SD13	20" <i>(Long frame only)</i>	NC
<input type="checkbox"/>	115SD14	21" <i>(Long frame only)</i>	\$342.00
<input type="checkbox"/>	115SD33	22" <i>(Long frame only)</i>	\$342.00

**Backrest**

<input type="checkbox"/>	115BT34	8° w/ Push Handles <i>Adjustable (16-3/4" to 20-1/8")</i>	NC
--------------------------	---------	--	----

**Back Height Preset**

<input type="checkbox"/>	115B91	16-3/4"	NC
<input type="checkbox"/>	115B96	17-7/8"	NC
<input type="checkbox"/>	115B39	19"	NC
<input type="checkbox"/>	115B98	20-1/8"	NC

**Armrest**

<input type="checkbox"/>	115AA33	Dual Post Height Adj - Desk 10"	\$321.00
<input type="checkbox"/>	115AA34	Dual Post Height Adj - Full 14"	\$321.00

**2. FRONT SEAT HEIGHT SET UP**

**Casters**

<input type="checkbox"/>	115CW6	5" Soft Roll	\$107.00
<input type="checkbox"/>	115CW36	8" x 2" Pneumatic	NC
<input type="checkbox"/>	115CW37	8" x 2" Pneumatic w/Airless Insert	\$332.00

**Front Hangers**

<input type="checkbox"/>	115H1	60° Swing Away	NC
<input type="checkbox"/>	115H146	70° Heavy-Duty Ft. Mount Lift Off <i>Only available with 5" caster and ang adj footplates. NA w/ extension tubes</i>	\$179.00
<input type="checkbox"/>	115H132	80° Heavy-Duty Ft. Mount Lift Off <i>Only available with 5" caster and ang adj footplates. NA w/ extension tubes</i>	\$179.00
<input type="checkbox"/>	115H41	ELR - Heavy-Duty Reinforced <i>Only available with aluminum footplates</i>	\$882.00

**Footplates**

*Heel loop or leg strap not included, See Footplate Options Section*

<input type="checkbox"/>	115F1	Composite	NC
<input type="checkbox"/>	115F25	Adult Angle Adj	\$143.00
<input type="checkbox"/>	115F5	Aluminum	NC

**Footplate Options**

<input type="checkbox"/>	115HL1	Heel Loop	\$36.00
--------------------------	--------	-----------	---------

**Extension Tubes**

<input type="checkbox"/>	115E1	Short (13"-15.5")	NC
<input type="checkbox"/>	115E2	Medium (15.5"-18")	NC
<input type="checkbox"/>	115E5	Long (18"-20.5") <i>NA w/ 5" caster</i>	NC

**Lower Leg Length Range Matrix**

Hanger	Ext Tube	Footplate		
		Composite	Alum Angle Adj	Aluminum
60°	Short	11.50" - 13.25"	11.75" - 14.00"	11.25" - 13.50"
	Medium	13.75" - 16.00"	14.00" - 16.25"	13.50" - 14.00"
	Long	16.00" - 16.50"	16.00" - 17.75"	NA
70°	No Ext Tube	NA	9.75" - 13.50"	NA
80°	No Ext Tube	NA	10.75" - 14.75"	NA
ELR	Short	NA	NA	13.00" - 17.00"
	Medium	NA	NA	14.25" - 17.00"

*All values are ± 0.25" and based on flat seat heights and do not account for the minimum 2" ground clearance or footplate thickness. Full range not available with all seat heights. ELR values based on 70° hanger angle*

**Seat Height Matrix**

Caster	Front/Rear
5"	17" ( <i>Different Axle &amp; Fork</i> )
5"	18.5"
8"	20"

**3. REAR SEAT HEIGHT SET UP**

**Adjustable Axle Plate**

<input type="checkbox"/>	115AL2	Standard (for 18.5", 20" Seat Ht.)	NC
<input type="checkbox"/>	115AL28	Low (only for 17" Seat Ht.)	NC

**Axle**

<input type="checkbox"/>	115AX1	Quick Release	NC
--------------------------	--------	---------------	----

**Rear Anti-Tip**

<input type="checkbox"/>	115AT1	Rear Anti-Tip	\$102.00
<input type="checkbox"/>	115AT3	Omit Rear Anti-Tip	NC

**Wheel Rim/Size**

<input type="checkbox"/>	115RW2	24" Spoke Wheel	NC
--------------------------	--------	-----------------	----

**Rear Wheel Tire**

<input type="checkbox"/>	115RT1	Pneumatic Tire	NC
<input type="checkbox"/>	115RT2	Pneumatic Airless Insert Tire *	\$209.00

**Handrim**

<input type="checkbox"/>	115HR1	Aluminum Anodized	NC
<input type="checkbox"/>	115HR4	Plastic Coated	\$133.00
<input type="checkbox"/>	115HR7	Omit Handrim	-\$24.00

**Wheel Lock**

<input type="checkbox"/>	115WL1	High Mount Push Wheel Lock	NC
<input type="checkbox"/>	115WL12	Do Not Mount Wheel Lock	NC
<input type="checkbox"/>	115WL10	Omit Wheel Lock	NC

**Extension Handle**

<input type="checkbox"/>	115EH1	6" Pair	\$61.00
--------------------------	--------	---------	---------

\* Options marked by an asterisk are not covered by AADL and are the client's financial responsibility.