



# Quickie P222 SE EIPW2

For Alberta Aids to Daily Living only



July 01, 2015

ACCOUNT #: \_\_\_\_\_ SHIP TO: \_\_\_\_\_  
 DATE: \_\_\_\_\_ NAME: \_\_\_\_\_  
 P.O. #: \_\_\_\_\_ ADDRESS: \_\_\_\_\_  
 CITY: \_\_\_\_\_ PROVINCE: \_\_\_\_\_  
 POSTAL CODE: \_\_\_\_\_ PHONE #: \_\_\_\_\_  
 BUYER: \_\_\_\_\_ MARKED FOR: \_\_\_\_\_



Product Info ▶

## QUICKIE P222 SE GROUP 4 CATEGORY

### QUICKIE P222 SE

- |                          |         |   |            |
|--------------------------|---------|---|------------|
| <input type="checkbox"/> | 116PC19 | Quickie P222 SE BD<br>(Solid Seat - No Power)               | \$8,007.00 |
| <input type="checkbox"/> | 116PC1  | Quickie P222 SE SD<br>(Solid seat/single power seat option) | \$8,007.00 |
| <input type="checkbox"/> | 116PC3  | Quickie P222 SE<br>(Base Only)                              | \$8,007.00 |

### PERFORMANCE PACKAGE

- |                          |        |   |    |
|--------------------------|--------|---|----|
| <input type="checkbox"/> | 116PF8 | 8.5 (8.5 mph / 350 LB Capacity)<br>(Max. Wt. cap. For CG Tilt 250lbs) | NC |
|--------------------------|--------|---|----|

### POWER CG TILT (53° TILT)

See electronic section for control options

- |                          |        |  |            |
|--------------------------|--------|--|------------|
| <input type="checkbox"/> | 116ST2 | P222 SE Power CG Tilt<br>Max. Wt. capacity for CG Tilt 250lbs.<br>Not available with base only | \$6,110.00 |
|--------------------------|--------|--|------------|

### TILT WEIGHT CAPACITY

- |                          |        |                               |    |
|--------------------------|--------|-------------------------------|----|
| <input type="checkbox"/> | 116LB3 | Up to 175 lb Weight Capacity  | NC |
| <input type="checkbox"/> | 116LB4 | 176 to 250 lb Weight Capacity | NC |

### LED LIGHT PACKAGE (R-NET ONLY)

- |                          |         |                                  |            |
|--------------------------|---------|----------------------------------|------------|
| <input type="checkbox"/> | 116LI04 | LED Light Package (NA Base only) | \$1,035.00 |
|--------------------------|---------|----------------------------------|------------|

## FRAME COLOUR OPTIONS

- |                          |         |  |    |                               |         |   |          |
|--------------------------|---------|--|----|-------------------------------|---------|---|----------|
| <input type="checkbox"/> | 116C140 | Aztec Gold                             | NC | <input type="checkbox"/>      | 116C103 | Rootbeer  | NC       |
| <input type="checkbox"/> | 116C2   | Black (Not available as base w/ Zebra) | NC | <input type="checkbox"/>      | 116C114 | Sparkle Silver  | NC       |
| <input type="checkbox"/> | 116C141 | Black Cherry                           | NC | <input type="checkbox"/>      | 116C145 | Sunrise Orange  | NC       |
| <input type="checkbox"/> | 116C10  | Black Opal                             | NC | <input type="checkbox"/>      | 116C113 | Titanium Colour Paint   | NC       |
| <input type="checkbox"/> | 116C38  | Blue Opal                              | NC | <input type="checkbox"/>      | 116C8   | Yellow  | NC       |
| <input type="checkbox"/> | 116C36  | Candy Blue                             | NC | <b>STANDARD MATTE COLOURS</b> |         |   |          |
| <input type="checkbox"/> | 116C14  | Candy Purple                           | NC | <input type="checkbox"/>      | 116C116 | Matte Black   | NC       |
| <input type="checkbox"/> | 116C35  | Candy Red                              | NC | <input type="checkbox"/>      | 116C151 | Matte Black Cherry  | NC       |
| <input type="checkbox"/> | 116C142 | Electric Blue                          | NC | <input type="checkbox"/>      | 116C152 | Matte Purple  | NC       |
| <input type="checkbox"/> | 116C143 | Evergreen                              | NC | <input type="checkbox"/>      | 116C150 | Matte Electric Blue   | NC       |
| <input type="checkbox"/> | 116C126 | Glow                                   | NC | <input type="checkbox"/>      | 116C149 | Matte Evergreen   | NC       |
| <input type="checkbox"/> | 116C144 | Green Apple                            | NC | <b>CUSTOM COLOURS</b>         |         |   |          |
| <input type="checkbox"/> | 116C147 | Hot Sparkle Pink                       | NC | <input type="checkbox"/>      | 116C42  | Stars & Stripes   | \$311.00 |
| <input type="checkbox"/> | 116C148 | Mauve Pink                             | NC | <input type="checkbox"/>      | 116C43  | Mossy Oak Camouflage  | \$311.00 |
|                          |         |  |    | <input type="checkbox"/>      | 116C112 | Desert Camouflage   | \$311.00 |
|                          |         |  |    | <input type="checkbox"/>      | 116C108 | Zebra (Choose above base color)<br>Std with White unless above base paint is<br>chosen NA with Black. | \$311.00 |

## BATTERIES & CHARGERS

### BATTERIES (GEL CELLS)

- |                          |         |   |          |
|--------------------------|---------|---|----------|
| <input type="checkbox"/> | 116AG03 | GP24 STD Post Pair w/ Straps  | \$806.00 |
| <input type="checkbox"/> | 116AG05 | Omit Batteries (Includes harness & straps<br>for standard post batteries) | NC       |

### BATTERY CHARGER

- |                          |         |                            |          |
|--------------------------|---------|----------------------------|----------|
| <input type="checkbox"/> | 116BJ19 | 8A Dual Mode (Compact) CTE | NC       |
| <input type="checkbox"/> | 116BJ3  | Omit Charger               | -\$41.00 |

**SEATING**

**SEAT WIDTH**

<input type="checkbox"/>	116S5	14" (not avail. w/ 19" & 20" seat depths)	NC
<input type="checkbox"/>	116S6	15" (not avail. w/ 19" & 20" seat depths)	NC
<input type="checkbox"/>	116S7	16" (not avail. w/19" & 20" seat depths and dual post arms )	NC
<input type="checkbox"/>	116S8	17"	NC
<input type="checkbox"/>	116S9	18"	NC
<input type="checkbox"/>	116S10	19"	NC
<input type="checkbox"/>	116S11	20"	\$418.00
<input type="checkbox"/>	116S27	22" (not avail. w/CG Power Tilt; Base width is 26.5" w/ 12"x4" tires)	\$418.00
<input type="checkbox"/>	116S28	24" (not avail. w/CG Power Tilt; Base width is 26.5" w/ 12"x4" tires)	\$520.00

**SEAT HEIGHT**

<input type="checkbox"/>	116HT10	18" (not avail. w/long ext. tubes; add 1.75" for CG Tilt)	NC
<input type="checkbox"/>	116HT12	19" ( add 1.75" for CG Tilt)	NC
<input type="checkbox"/>	116HT14	20" ( add 1.75" for CG Tilt)	NC

**SEAT BASE TYPE**

<input type="checkbox"/>	116SL2	Seat Sling Parapak (n/a w/CG Power Tilt)	NC
<input type="checkbox"/>	116SL5	Seat Pan Aluminum	\$168.00
<input type="checkbox"/>	116SL6	Omit Seat Sling	-\$16.00

**SEAT DEPTH**

<input type="checkbox"/>	116SD7	14"	NC
<input type="checkbox"/>	116SD8	15"	NC
<input type="checkbox"/>	116SD9	16"	NC
<input type="checkbox"/>	116SD10	17"	NC
<input type="checkbox"/>	116SD11	18"	NC
<input type="checkbox"/>	116SD12	19"	NC
<input type="checkbox"/>	116SD73	20"	NC

**SEAT CUSHIONS**

<input type="checkbox"/>	116CU1	2" Cushion	\$158.00
--------------------------	--------	------------	----------

**JAY Cushions (Optional)**

To order any JAY cushions on-chair, please see the JAY price list.

Cushion name or part number must be filled out with desired width and depth

<input type="checkbox"/>	116JC20	JAY Bundled Cushion	Price List
		Cushion Name or Part #	
		Width: _____	Depth: _____

**BACKS**

**BACKREST TYPE**

<input type="checkbox"/>	116BT74	Folding Angle Adj. Backrest	\$153.00
<input type="checkbox"/>	116BT3	Non-Folding Angle Adj. Backrest	NC

**BACK HEIGHT (UPHOLSTERY)**

<input type="checkbox"/>	116B4	Low 15.5"	\$112.00
<input type="checkbox"/>	116B8	Medium 17.5"	NC
<input type="checkbox"/>	116B13	Tall 19.5"	\$112.00

**BACK UPHOLSTERY**

<input type="checkbox"/>	116BU27	Parapak Upholstery (Select back height)	NC
<input type="checkbox"/>	116BU26	Jay Solid Backrest (NC with Power CG Tilt) (Select backheight; Decrease seat depth by 1")	\$515.00
<input type="checkbox"/>	116BU14	Omit Back Upholstery	-\$16.00

**JAY SOLID BACKREST BACK HEIGHT**

<input type="checkbox"/>	116B36	16"	NC
<input type="checkbox"/>	116B37	17"	NC
<input type="checkbox"/>	116B38	18"	NC
<input type="checkbox"/>	116B39	19"	NC
<input type="checkbox"/>	116B40	20"	NC
<input type="checkbox"/>	116B41	21"	NC
<input type="checkbox"/>	116B42	22"	NC
<input type="checkbox"/>	116B43	23"	NC
<input type="checkbox"/>	116B44	24"	NC

**BACK POST TYPE**

<input type="checkbox"/>	116BP1	Straight Back Post	NC
<input type="checkbox"/>	116BP41	8 Bend Back Post - Not avail. w/ Jay Solid backrest	NC

**JAY 3 OPTION**

<input type="checkbox"/>	116JB01	Jay J3 (Select OMIT Back Uph above) See J3 OF	
--------------------------	---------	---	--

**POSTURAL SUPPORTS**

**HEADREST (INCLUDED ON SPECIFIED SPECIALTY CONTROLS)**

*JAY Solid Backrest Only*

<input type="checkbox"/>	<b>Headrest (8" x 4.5")</b>	<b>\$337.00</b>
	116HE65 Curved Headrest 8" x 4.5"	
	116MH1 Removable Mount	
<input type="checkbox"/>	<b>Headrest (9" x 4.5")</b>	<b>\$337.00</b>
	116HE66 Curved Headrest 9" x 4.5"	
	116MH1 Removable Mount	

**PAD STYLE**

*JAY Solid Backrest Only*

<input type="checkbox"/>	<b>Lateral Pads Curved</b>	<b>\$490.00</b>
	116KE1 Lateral Pad Curved Pair	
	116LM1 Swing Away Mount Pair	
<input type="checkbox"/>	<b>Lateral Pads Straight</b>	<b>\$490.00</b>
	116KE2 Lateral Pad Straight Pair	
	116LM1 Swing Away Mount Pair	

**PAD SIZE**

<input type="checkbox"/>	116KD1 Lateral Support 3" x 6" (pr)	NC
<input type="checkbox"/>	116KD2 Lateral Support 4" x 5" (pr)	NC
<input type="checkbox"/>	116KD3 Lateral Support 4" x 6" (pr)	NC
<input type="checkbox"/>	116KD4 Omit Lateral Support	-\$26.00

**ARMRESTS**

**ARMREST (PAIR)**

*Single Post Height Adj. Armrests are N/C on Power CG Tilt*

<input type="checkbox"/>	116AA131 Dual Post Flip Back Short (14-16in)	NC
<input type="checkbox"/>	116AA132 Dual Post Flip Back Long (16-20in)	NC
<input type="checkbox"/>	116AA51 Dual Post Hgt Adj - Short (14-16 in)	\$321.00
<input type="checkbox"/>	116AA52 Dual Post Hgt Adj - Long (16-20 in)	\$321.00
<input type="checkbox"/>	116AA8 Single Post Height Adjustable	\$286.00
<input type="checkbox"/>	116AA14 Omit Armrests and mounts	NC

**ARMPADS (PAIR)**

<input type="checkbox"/>	116AM7 Desk Length Waterfall (10") (n/a. w/ sgl. post Armrest w/Jay Solid Backrest)	NC
<input type="checkbox"/>	116AM8 Full Length Waterfall (14") (n/a. w/ sgl. post Armrest w/Jay Solid Backrest)	NC
<input type="checkbox"/>	116AM1 Desk Length Standard (10")	NC
<input type="checkbox"/>	116AM2 Full Length Standard (14")	NC

**FOOTRESTS**

**HANGER (PAIR)**

<input type="checkbox"/>	116H12 70° Fixed Straight (only avail. w/ solid 1-piece footplate.; n/a. w/ CG Power Tilt)	NC
<input type="checkbox"/>	116H2 70° Swing Away Hanger (Not. Avail w/ solid 1-piece footplate.)	NC
<input type="checkbox"/>	116H13 80° Fixed Tapered (only avail. w/ solid 1-piece footplate.; n/a. w/ CG Power Tilt)	NC
<input type="checkbox"/>	116H37 90° Swing Away (Only available with kids angle adjustable footplates)	NC
<input type="checkbox"/>	116H16 Elevating Hanger (Req. aluminum footplate or angle adj.)	\$342.00
<input type="checkbox"/>	116H32 Omit Hangers	-\$16.00

**FOOTREST (PAIR)**

<input type="checkbox"/>	<b>Composite w/ Heel Loops</b>	<b>NC</b>
	Not Avail. w/ fixed hangers	
	116F1 Composite Footplate	
	116HL1 Heel Loop	
<input type="checkbox"/>	<b>Adult Angle Adjust. w/ heel loop</b>	<b>\$148.00</b>
	Not Avail. w/ fixed hangers	
	116F282 Adult Angle Adj. Footplate	
	116HL10 Heel Loops	
<input type="checkbox"/>	<b>Kids Angle Adjust. w/ heel loop</b>	<b>\$148.00</b>
	Only avail. w/ 90° hangers	
	116F283 Kids Angle Adj. Footplate	
	116LS13 Leg Straps	
	116F287 Aluminum (Only avail. with ELRs)	
<input type="checkbox"/>	<b>116F60 Solid One Piece</b>	<b>NC</b>
	not avail. w/Power CG Tilt	
<input type="checkbox"/>	<b>Two-piece Locking w/Heel Loop</b>	<b>NC</b>
	116F58 Two-piece Locking Footplate	
	116HL10 Heel Loop	
<input type="checkbox"/>	<b>Multi-Positional w/Heel Loop</b>	<b>\$209.00</b>
	Not Avail. w/ fixed hangers	
	116F47 Multi-Positional Footplate	
	116HL10 Heel Loop	
<input type="checkbox"/>	<b>116F6 Omit Foot Plate w/Heel Loop</b>	<b>-\$16.00</b>

**EXTENSION TUBE**

*Not available with 90° hangers.*

<input type="checkbox"/>	116E1 Short (13" - 15")	NC
<input type="checkbox"/>	116E2 Medium (15" - 18") Not available with multi-position footplate.	NC
<input type="checkbox"/>	116E5 Long (18" - 20") Not available with 70° or 80° hangers.	NC

**FOOTREST ACCESSORIES**

<input type="checkbox"/>	116LS3 Extra Leg Strap, Black	\$51.00
--------------------------	-------------------------------	---------

**ELECTRONICS**

**DRIVE CONTROL OPTIONS**

<input type="checkbox"/>	116ER04	R-Net Non-Expandable Drive Control (120A)	NC
<input type="checkbox"/>	<b>R-Net Expandable Drive Control Upgrade (120A)</b>		<b>\$1,035.00</b>
	116EX03	Expandable 120A Controller - R-Net	
	116EX04	Harnessing for 120A Expandable Controller	

**R-NET JOYSTICK (FOR EXPANDABLE ELECTRONICS)**

<input type="checkbox"/>	116J70	Std LED Joystick for Expandable Ctrl	NC
<input type="checkbox"/>	116J93	Joystick w/stereo Jacks and Monochrome Display	\$245.00
<input type="checkbox"/>	116J88	Joystick W/Toggle, Spd Pot, Stereo Jacks & Color Display	\$245.00
<input type="checkbox"/>	116J95	Omit JS for use with OMNI	NC

**JOYSTICK MOUNT (FIXED OR SWING AWAY)**

<input type="checkbox"/>	116JM10	Fixed Mount	NC
<input type="checkbox"/>	116JM11	Swing Away Mount	\$377.00

**JOYSTICK MOUNT (LOCATION)**

*Height Adjustable.*

<input type="checkbox"/>	116HM1	Right Side Mount	NC
<input type="checkbox"/>	116HM2	Left Side Mount	NC

**JOYSTICK HANDLE (OPTIONAL)**

<input type="checkbox"/>	116JH1	Ball Handle	\$92.00
<input type="checkbox"/>	116JH2	Foam Ball Handle	\$92.00
<input type="checkbox"/>	116JH4	T-Handle	\$92.00
<input type="checkbox"/>	116JH6	Goalpost	\$92.00

**R-NET OMNI & OUPUT OPTIONS**

*(For expandable controls, choose all that apply)*

<input type="checkbox"/>	116ON1	Rnet OMNI Input Control Module and display with built in Infrared and armrest mount. Choose mount below. (required for Specialty Input Devices.)	\$1,933.00
--------------------------	--------	--	------------

**OMNI MOUNT**

*(If used with Joystick must be mounted opposite side)*

<input type="checkbox"/>	116DS03	Swing Away Left Side OMNI Mount	\$240.00
<input type="checkbox"/>	116DS04	Swing Away Right Side OMNI Mount	\$240.00

**SWITCH INPUTS FOR JOYSTICKS OR OMNI**

*Drive or Mode / Profile Control*

<input type="checkbox"/>	116SW22	Stereo to 2 Mono Splitter	\$80.00
<input type="checkbox"/>	116SW6	Ribbon Switch	\$61.00
<input type="checkbox"/>	116SW1	Buddy Button Switch	\$71.00
<input type="checkbox"/>	116SW2	Disc Switch	\$158.00
<input type="checkbox"/>	116SW9	Wobble Swich	\$678.00
<input type="checkbox"/>	116SW4	Micro Light Switch	\$97.00
<input type="checkbox"/>	116SW24	Egg Switch-Black	\$117.00
<input type="checkbox"/>	116SW21	Dual toggle switch only	\$173.00
<input type="checkbox"/>	116SW23	4-way toggle switch w/5th button DB9 Connector.	\$224.00

**SWITCH MOUNT BRACKET**

<input type="checkbox"/>	116EM01	Single Switch Mount Bracket	\$199.00
<input type="checkbox"/>	116EM02	Dual toggle switch mount bracket	\$199.00
<input type="checkbox"/>	116EM03	4-way toggle switch mount bracket.	\$199.00

**SPECIALTY DRIVER CONTROL OPTIONS**

*Requires expandable controls and OMNI. All packages must be added to the specialty control on-chair package. Not available with Upholstered backs.*

<input type="checkbox"/>	<b>Quickie Sip n' Puff System</b>		<b>\$3,641.00</b>
	116ID34	Sip N' Puff Kit	
	116II01	Headrest	
	116IJ01	Headrest Removable Mount	
	116IH1	Sip n' Puff Tubing and Hardware	
<input type="checkbox"/>	<b>ASL Head Array Package (Adult)</b>		<b>\$5,197.00</b>
	116II02	ASL Proximity Head Array (Adult)	
	116IJ01	Headrest Removable Mount	
<input type="checkbox"/>	<b>ASL Head Array Package (Pediatric)</b>		<b>\$5,197.00</b>
	116II03	ASL Proximity Head Array (Ped)	
	116IJ01	Headrest Removable Mount	
<input type="checkbox"/>	<b>Switch-It Head Array Package (Adult)</b>		<b>\$5,197.00</b>
	116II04	Switch-it Proximity Head Array	
	116IJ01	Headrest Removable Mount	
<input type="checkbox"/>	<b>Switch-It Head Array Package (Pediatric)</b>		<b>\$5,197.00</b>
	116II05	Switch-it Proximity Head Array	
	116IJ01	Headrest Removable Mount	
<input type="checkbox"/>	<b>ASL Proportional Mini Joystick Chin / Lip Control Package</b>		<b>\$5,197.00</b>
	116ID36	Mini Joystick	
	116II01	Headrest	
	116IJ01	Headrest Removable Mount	
<input type="checkbox"/>	<b>ASL Proportional Mini Joystick Midline with Gatlin Package</b>		<b>\$5,197.00</b>
	116ID37	Mini joystick	
	116IJ03	ASL Proportional S/A Joystick Mount	
<input type="checkbox"/>	<b>Quickie Proportional Chin Control Package</b>		<b>\$5,197.00</b>
	116ID45	Proportional Joystick	
	116IJ03	ASL Proportional S/A Joystick Mount	

**ATTENDANT CONTROL**

<input type="checkbox"/>	1160A02	Attendant Control for R-net	\$1,035.00
--------------------------	---------	-----------------------------	------------

**ATTENDANT MOUNTING LOCATION**

<input type="checkbox"/>	116DN1	Mounting Location: Right Back Cane	NC
<input type="checkbox"/>	116DN2	Mounting Location: Left Back Cane	NC

**ELECTRONIC PROGRAMMERS**

<input type="checkbox"/>	116PR5	R-Net PC Programmer with Dongle	\$603.00
<input type="checkbox"/>	116PR6	DTT R-Net cables included	\$1,244.00



**POWERED SEATING**

**POWERED SEATING CONTROL METHOD (TILT)**

*For use with power seating actuators. Switch activation is stand-alone and does not require expandable or thru-drive electronics.*

<input type="checkbox"/>	116EO10	Switch Control 1 actuator w/dual toggle and mount. (For operating Tilt)	NC
<input type="checkbox"/>	116EO11	Switch Control 1 actuator w/buddy button & mount. (For operating Tilt)	NC

**SWITCH MOUNT**

*(For use with switch control powered seating)*

<input type="checkbox"/>	116EY01	Switch Mount Right	NC
<input type="checkbox"/>	116EY02	Switch Mount Left	NC

**POWERED SEATING THRU-DRIVE CONTROL METHOD**

"Assignable Button" feature listed below (available on Rnet Systems) allows for direct "shortcut" access to various seating functions via joystick buttons. Functions accessed by depressing the button for more than 1 second. By selecting the thru-drive w/assignable button the chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle.

<input type="checkbox"/>	116EO15	Thru Drive 1 Actuator w/out assignable buttons. For use with R-Net (For operating Tilt thru Joystick.)	\$1,015.00
<input type="checkbox"/>	116EO13	Thru Drive 1 Actuator with assignable buttons. (For operating Tilt thru Joystick and/or Joystick. Requires R-Net.)	\$1,015.00

**WHEELS**

**DRIVE WHEEL TYPE**

<input type="checkbox"/>	116DW5	12" x 4" Pneumatic Tire	NC
<input type="checkbox"/>	116DW6	12" x 4" Pneumatic w/Ins.	\$173.00
<input type="checkbox"/>	116DW27	12" x 3" Pneumatic Tire	NC
<input type="checkbox"/>	116DW28	12" x 3" Pneumatic w/Ins.	\$173.00

*The 12" x 3" wheel size reduces the overall width specification by approx. 2" and also alters the overall appearance of the chair.*

**DRIVE WHEEL TIRE**

<input type="checkbox"/>	116DT1	Knobby Tire	NC
--------------------------	--------	-------------	----

**WHEEL LOCK**

<input type="checkbox"/>	116WL26	Push to Lock	\$82.00
<input type="checkbox"/>	116WC1	Do Not Mount	NC

**WHEEL LOCK EXTENSION HANDLE**

*Wheel-lock extensions handles are not available with 22" & 24" seats*

<input type="checkbox"/>	116EH1	Extension Handle 6"	\$61.00
--------------------------	--------	---------------------	---------

**CASTERS**

**FRONT CASTER WHEEL**

<input type="checkbox"/>	116CW17	8" Pneumatic	NC
<input type="checkbox"/>	116CW20	8" Airless with Insert	\$143.00
<input type="checkbox"/>	116CW29	9" Pneumatic (not avail. With Suspension Fork)	\$102.00
<input type="checkbox"/>	116CW28	9" Pneumatic w/Insert (not avail. With Suspension Fork)	\$158.00

**CASTER FORK**

*Suspension Fork Options are only available with 8" casters.*

<input type="checkbox"/>	116CF31	Solid Caster Fork	NC
<input type="checkbox"/>	116CF35	Suspension Fork Med. 150-250 lbs. (Only avail. w/ 8" front casters)	\$306.00
<input type="checkbox"/>	116CF36	Suspension Fork Heavy. 251-350 lbs. (Only avail. w/ 8" front casters)	\$306.00

**MANUAL CONVERSION KIT**

116CV2 Manual Wheelchair Conversion Kit \$959.00

**AXLE**

116AX1 Quick Release NC  
 116AX11 Quad Quick Release \$153.00

**WHEEL RIM/SIZE**

116RW1 24" Mag NC  
 116RW2 24" Spoke (not avail w/14" seat width) NC  
 116RW7 Omit Manual Rear Wheel NC

**REAR WHEEL TIRE**

116RT1 Pneumatic NC  
 116RT2 Pneumatic Airless Insert \$194.00  
 116RT3 Full Polyurethane NC  
 116RT4 Low Polyurethane NC  
 116RT6 High Pressure Clincher (not available with mag wheels) \$87.00

**HAND RIM**

116HR1 Aluminum Anodized NC  
 116HR4 Plastic Coated \$133.00  
 116HR7 Omit handrim -\$20.00

**PROJECTIONS**

116HR8 8 Oblique \$122.00

**WHEEL LOCKS**

116WC2 Push to Lock Wheel Lock \$82.00  
 116WC1 Do not mount NC  
 116ET1 6" Extension Handles \$51.00

**ACCESSORIES**

**BACKPACK**

116A1 Black Utility Backpack \$94.00  
 116A2 Black/Charcoal Kids \$87.00

**SEAT POUCH**

116A8 Black/Charcoal \$71.00

**POSITIONING BELT**

116A13 2" Aircraft Buckle Padded \$173.00  
 116A12 2" Aircraft Buckle \$117.00  
 116A142 1.5" Auto Buckle \$82.00

**WHEEL CHAIR TRAY TABLE**

*Requires standard armpads to mount.*

116A31 Medium Tuffak (15"-17") Requires std. armpads to mount tray \$321.00  
 116A32 Large Tuffak (18"-20") Requires std. armpads to mount tray \$321.00  
 116A33 X-Large Tuffak (21"-22") Requires std. armpads to mount tray \$321.00

**MISCELLANEOUS**

116A82 Caddy \$77.00  
 116A114 Dual Post Armrest Pouch \$24.00

**ADDITIONAL SPECIFICATIONS**

**Weight:** 115 lbs (without batteries)  
**Overall Width:** Using 12" x 4" Drive Tire: 24.5" using 14"-20" Seat Width  
 26.5" using 22" and 24" Seat Width  
**Overall Length:** 31.5" (Without Footrests)  
 Deduct 2" if using 12" x 3" Tire  
**Range:** As far as 30 miles (With Group 24 batteries depending on conditions and user environment)  
**Top Speed:** Up to 8.5 mph