



AADL #W135

Zippie Kidz EIKZ

For Alberta Aids to Daily Living only



July 01, 2015

ACCOUNT #: _____ SHIP TO: _____
 DATE: _____ NAME: _____
 P.O. #: _____ ADDRESS: _____
 CITY: _____ PROVINCE: _____
 POSTAL CODE: _____ PHONE #: _____
 BUYER: _____ MARKED FOR: _____



Product Info ►

ZIPPIE KIDZ

013WM36 Zippie KIDZ \$2,810.00

1. FRAME STYLE

<input type="checkbox"/> 013FC3 Standard Configuration STD	<input type="checkbox"/> 013CV1 Reverse Configuration \$87.00
<input type="checkbox"/> 013FC4 Reverse Configuration (Lg. wheels in front & casters in back. Not avail. w/ 20" rear wheels or 5" casters. NC	Conversion Kit (Not avail. w/ 20" rear wheels or 5" casters. Used to reverse configur chair w/Lg. wheels in front & casters in back. NC

2. FRAME COLOUR

STANDARD COLOURS

<input type="checkbox"/> 013C2 Black NC
<input type="checkbox"/> 013C10 Black Opal NC
<input type="checkbox"/> 013C38 Blue Opal NC
<input type="checkbox"/> 013C36 Candy Blue NC

3. FRAME WIDTH

(Seat Width 1/2" Narrower)

<input type="checkbox"/> 013W1 10" STD	<input type="checkbox"/> 013W4 13" NC
<input type="checkbox"/> 013W2 11" NC	<input type="checkbox"/> 013W5 14" NC
<input type="checkbox"/> 013W3 12" NC	

4. SEAT UPHOLSTERY DEPTH

12" rear wheels on regular frame are not available with high-pull wheel locks.

REGULAR FRAME

<input type="checkbox"/> 013SD1 8" Seat Depth STD
<input type="checkbox"/> 013SD2 9" Seat Depth NC
<input type="checkbox"/> 013SD3 10" Seat Depth NC

LONG FRAME

<input type="checkbox"/> 013SD4 11" Seat Depth STD
<input type="checkbox"/> 013SD5 12" Seat Depth NC
<input type="checkbox"/> 013SD6 13" Seat Depth NC
<input type="checkbox"/> 013SD7 14" Seat Depth NC



5. SEATING

<input type="checkbox"/>	013SL1	Seat Upholstery	STD
<input type="checkbox"/>	013SL6	Omit Seat Upholstery	-\$18.00
<input type="checkbox"/>	013CU1	2" Seat Cushion	\$158.00

FORWARD-FOLDING

<input type="checkbox"/>	013B101	Regular Height 8"-10"	\$102.00
<input type="checkbox"/>	013B102	Tall Height 10"-12"	\$102.00

NON-FOLDING

<input type="checkbox"/>	013B103	Regular Height 8"-10"	\$102.00
<input type="checkbox"/>	013B104	Tall Height 10"-12"	NC

6. BACKREST

<input type="checkbox"/>	013SR1	Stroller Handle	\$138.00
--------------------------	--------	-----------------	----------

<input type="checkbox"/>	013BU14	Omit Back Upholstery	-\$18.00
<input type="checkbox"/>	013BT14	Omit Backrest & Back Upholstery	-\$43.00

7. FOOTREST

<input type="checkbox"/>	013F54	Regular 8 1/2" to 10 3/4"	STD
<input type="checkbox"/>	013F46	Long 11" to 13 1/4"	NC

<input type="checkbox"/>	013LS5	Leg Strap	NC
<input type="checkbox"/>	013F38	High Mount	\$82.00

8. ARMRESTS (T-POST)

Adjust in 1/4" increments

<input type="checkbox"/>	013AA12	Kidz Std 4 1/2" to 5 3/4"	STD
<input type="checkbox"/>	013AM14	Foam Pad	NC
<input type="checkbox"/>	013AM1	Desk Pad	NC

<input type="checkbox"/>	013AA13	Kidz Std 5 1/2" to 6 3/4"	NC
<input type="checkbox"/>	013AM14	Foam Pad	STD
<input type="checkbox"/>	013AM1	Desk Pad	NC
<input type="checkbox"/>	013AA14	Omit Armrest	NC

9. CASTERS

CASTER WHEEL

<input type="checkbox"/>	013CW1	3" Micro	NC
<input type="checkbox"/>	013CW30	3" Micro Light up	\$61.00
<input type="checkbox"/>	013CW5	5" Low-Profile Polyurethane	STD

CASTER OPTIONS

<input type="checkbox"/>	013CF4	3/4" Ext. Fork Stem Bolt	NC
<input type="checkbox"/>	013CW12	Omit Caster Wheels	NC

REAR WHEEL SIZE

<input type="checkbox"/>	013RW13	12" Mag (Avail. w/ 3" casters only.) 12" rear wheels on reg. frame are not avail. With high pull wheel locks.	NC
<input type="checkbox"/>	013RW14	16" Mag (Not available with full poly tires)	STD
<input type="checkbox"/>	013RW3	20" Mag Wheel	NC
<input type="checkbox"/>	013RW74	20" Lite Spoke	\$270.00

REAR WHEEL TIRE

<input type="checkbox"/>	013RT1	Pneumatic	NC
<input type="checkbox"/>	013RT2	Pneumatic Airless Insert	\$209.00
<input type="checkbox"/>	013RT3	Full Polyurethane	NC



AXLES

AXLES

<input type="checkbox"/>	013AX1	Quick Release Axles	STD
<input type="checkbox"/>	013AX11	Quad Quick Release Axles (Only avail. w/20" rear wheels. Not avail. w/16" rear wheel rev. config.)	\$163.00

HANDRIM

**Not available with 12" or 16" rear wheels.*

<input type="checkbox"/>	013HR1	Aluminum Anodized	NC
<input type="checkbox"/>	013HR4	Plastic Coated	\$133.00
<input type="checkbox"/>	013HR7	Omit Handrims	-\$12.00

WHEEL LOCK

<input type="checkbox"/>	013WL9	High Push	STD
<input type="checkbox"/>	013WL8	High Pull (12" rear wheels on reg. frame are not avail. w/ high pull wheel locks. Not avail. w/16" rear wheel rev. config.)	NC
<input type="checkbox"/>	013WL10	Omit	NC

EXTENSION HANDLE

<input type="checkbox"/>	013EH1	6" Ext. Handles (Avail. w/ high- push wheel locks only.)	\$61.00
--------------------------	--------	---	---------

ANTI-TIP (STANDARD CONFIGURATION)

<input type="checkbox"/>	013AT1	Rear Anti-Tip Tubes	\$102.00
<input type="checkbox"/>	013AO1	Front Anti-Tip Tubes	\$133.00

ANTI-TIP (REVERSE CONFIGURATION)

<input type="checkbox"/>	013AT19	Rear Anti-Tip Tubes	\$122.00
<input type="checkbox"/>	013AO1	Front Anti-Tip Tubes	\$122.00

ACCESSORIES

UPPER EXTREMITY SUPPORT TRAY

<input type="checkbox"/>	013A29	X-Sm (Suggested 10"-12" Width)	\$321.00
<input type="checkbox"/>	013A30	Sm (Suggested 13"-14" Width)	\$321.00

SPOKE GUARD

(Not available with 12" & 16" rear wheels)

<input type="checkbox"/>	013A23	Black	\$133.00
<input type="checkbox"/>	013A303	Clear Zippie	\$133.00

BACK PACK (KIDS BACKPACK INCLUDED)

<input type="checkbox"/>	013A1	Black Utility Backpack	\$94.00
<input type="checkbox"/>	013A2	Black Kids	\$87.00