


Customer				Ship To			
Account				Name			
Account #		Date		Address			
☐ Quote	☐ Order	PO #		City		Province	
Buyer				PC	Phone		
Mark For				Email / Fax			

Technical Data Overview			
Product Width	14" to 20"	Base	Pre-contoured closed cell foam
Product Depths	14" to 20"	Insert	Multizone elastomer insert with 2 standard pressure adjustment pods
Product Height	2"	Cover	Double spacer mesh outer cover with durable side bands and non-slip base
Product Weight	2 lbs.	JAY Your Way Available	Yes
User Weight Capacity	330 lbs.		

Measuring and Sizing Guide



X = Width

Measurement made from outside edge to outside edge of the cushion.

Example: JLP1618 | X=16 (Width)

Y = Depth

Measurement from the front lower edge of the cushion to the back lower edge of the cushion.

Example: JLP1618 | Y=18 (Depth)

		Width							
		14	15	16	17	18	19	20	
Depth	14	●	○	○					
	15	●	●	●	○	○			
	16	●	●	●	●	●	●	○	
	17	●	●	●	●	●	●	●	
	18	○	●	●	●	●	●	●	
	19	○	○	○	●	●	●	●	
	20		○	○	○	●	●	●	

● Standard Offering

○ Available Through JAY Your Way

JAY LitePro Cushion									
✓	OPTION		DESCRIPTION	RETAIL	✓	OPTION		DESCRIPTION	RETAIL
☐	JLP1414		JAY LitePro 14" x 14"	\$925	☐	JLP1719		JAY LitePro 17" x 19"	\$925
☐	JLP1415		JAY LitePro 14" x 15"	\$925	☐	JLP1816		JAY LitePro 18" x 16"	\$925
☐	JLP1416		JAY LitePro 14" x 16"	\$925	☐	JLP1817		JAY LitePro 18" x 17"	\$925
☐	JLP1417		JAY LitePro 14" x 17"	\$925	☐	JLP1818		JAY LitePro 18" x 18"	\$925
☐	JLP1515		JAY LitePro 15" x 15"	\$925	☐	JLP1819		JAY LitePro 18" x 19"	\$925
☐	JLP1516		JAY LitePro 15" x 16"	\$925	☐	JLP1820		JAY LitePro 18" x 20"	\$925
☐	JLP1517		JAY LitePro 15" x 17"	\$925	☐	JLP1916		JAY LitePro 19" x 16"	\$925
☐	JLP1518		JAY LitePro 15" x 18"	\$925	☐	JLP1917		JAY LitePro 19" x 17"	\$925
☐	JLP1615		JAY LitePro 16" x 15"	\$925	☐	JLP1918		JAY LitePro 19" x 18"	\$925
☐	JLP1616		JAY LitePro 16" x 16"	\$925	☐	JLP1919		JAY LitePro 19" x 19"	\$925
☐	JLP1617		JAY LitePro 16" x 17"	\$925	☐	JLP1920		JAY LitePro 19" x 20"	\$925
☐	JLP1618		JAY LitePro 16" x 18"	\$925	☐	JLP2017		JAY LitePro 20" x 17"	\$925
☐	JLP1716		JAY LitePro 17" x 16"	\$925	☐	JLP2018		JAY LitePro 20" x 18"	\$925
☐	JLP1717		JAY LitePro 17" x 17"	\$925	☐	JLP2019		JAY LitePro 20" x 19"	\$925
☐	JLP1718		JAY LitePro 17" x 18"	\$925	☐	JLP2020		JAY LitePro 20" x 20"	\$925

JAY Pressure Adjustment Pods

Standard pods (pair) included in base cushion selection above - only select JLPSTD if extra required. Firm pods (pair) not included in base cushion selection above - select JLPFIRM if required.

✓	OPTION		DESCRIPTION	RETAIL	✓	OPTION		DESCRIPTION	RETAIL
☐	JLPSTD		JAY LitePro Standard Pods (Pair)	\$75	☐	JLPFIRM		JAY LitePro Firm Pods (Pair)	\$75

JAY Your Way Modifications

Add an "M" in front of the item number to order a modification through JAY Your Way.

OPTION	DESCRIPTION	✓		✓		RETAIL	✓	OPTION	DESCRIPTION	RETAIL
IW	Increase Width	☐	+1"	☐	+2"	N/C	☐	CVR	Reverse Dartex Cover	N/C
IL	Increase Length	☐	+1"	☐	+2"	N/C				

JAY LitePro Product Images

