

Design Template

OCTANE



EIRGK03



Rider Information

Name: _____

Phone: _____

Email: _____

Height: _____ **Weight:** _____

Current Chair: _____ **Current Cushion:** _____

Ability/Diagnosis: _____

RGK Contact Authorization:

RGK would like to send you exclusive information, news and promotional offers on RGK.

I would like to be contacted by RGK:

Signature: _____ **Date:** _____

*We do not pass your personal information onto any third party.

Account Information

Account Name : _____

Account #: _____ **PO#:** _____

Contact Name : _____

Phone: _____

Email: _____

Ship to: **Account** ☐ **Rider** ☐

Name: _____


Address: _____

RGK Genius: _____


Frame Styles

- | | | | |
|--------------------------|----------|--|---------|
| ☐ | 703XZ010 | Octane (Max User Weight 110kg/242lbs.) | \$8,887 |
| ☐ | 703XA070 | Open design | STD |
| ☐ | 703XA080 | Dual-Tube design/Closed frame design | NCO |
| ☐ | 703XW009 | Quote Only (No CAD for Approval) | STD |
| ☐ | 703XW011 | Order Only (No CAD for Approval) | STD |
| ☐ | 703XW012 | Quote CAD drawing for approval (1st CAD is NCO) | \$550 |
| ☐ | 703XW010 | Order CAD drawing for approval (1st CAD is NCO) | \$300 |

CAD approved by: ☐ 703XW013 Customer ☐ 703XW014 RGK Genius



Open Design



Dual-Tube Design

Measurements

Please Note: All RGK chairs are hand built and have a manufacturing tolerance of +/- 5mm (.25"). All measurements can be done in 5mm (.25") increments.

Please select either Millimeters or Inches.

Measurements: ☐ **Millimeters (mm.)** ☒ **Inches (in.)**

703XY010 Rear Seat Width Measured from inside of side guard to inside of opposite side guard. Or outside seat tube to outside of opposite seat tube near backrest. (Min. 250mm/9.84" and Max. 500mm/19.69")

Current Chair: mm. in. New Chair: mm. in.

703XY020 Seat Depth/Seat Length Measured from front of back post to the front edge of seat sling. Frame bend will start .5" in front of seat sling unless Extra Frame Length requested. (Min. 250mm/9.84" and Max. 500mm/19.69")

Current Chair: mm. in. New Chair: mm. in.

☐ **703XY030 Extra Frame Length** Measured from front edge of seat sling to beginning of Front frame bend. (Seat length + Extra seat length)



Current Chair: mm. in. New Chair: mm. in.

703XY050 Front Frame Height Measured from floor to apex of seat tube. (Min. 380mm/14.96" and Max. 550mm/21.65")

Current Chair: mm. in. New Chair: mm. in.

703XY060 Rear Frame Height Measured from floor to top of seat tube at back edge of seat sling. (Min. 360mm/14.17" and Max. 550mm/21.65")

Current Chair: mm. in. New Chair: mm. in.



Design Template

(New Options in Blue) **OCTANE**

703XY070 Occupied Frame Length/Frame Length Measured from the front of the backrest tube to the front edge of the footrest. (Min. 390mm/15.35" and Max. 750mm/29.53")

Current Chair: New Chair:

703XY072 Caster Stem Position Measured from the front edge of the footrest to the centerline of the caster barrel. **703XY073** RGK to Optimize. ☐

Current Chair: New Chair:

☐ **703XA020 Ergonomic Seat** Measured from front of back post to beginning of upward bend of seat tube. \$125

Current Chair: New Chair:

☐ **703XY080 Front Seat Taper/Seat Taper Design** Measured from inside of front frame tube at the front frame bend to inside of opposite front frame tube at the front frame bend.

Current Chair: New Chair:

☐ **703XY091 Footrest Width/Parallel Taper** Measured from inside of front frame tube at footrest adjustment to inside of opposite front frame tube at footrest adjustment. (Tip: when ordered together 703XY091 is usually 20mm (1") less than 703XY080).

Current Chair: New Chair:

703XY115 Front Footrest Height Measured from ground up to the top of the Front footrest tube or cover.

Current Chair: New Chair:

703XY120 Rear Footrest Height Measured from ground up to the top of the Rear footrest tube or cover. Standard is 20mm less than Front Footrest Height.





Current Chair: New Chair:

703XY125 Footrest Length Measured from front edge of footrest to back edge of footrest. (Max: 140mm for certified weight)

Current Chair: New Chair:

☐ **703XB200 High Mount Footrest/High Level Footrest** Measured from ground up to the top of the Front footrest tube or cover. \$207

Current Chair: New Chair:



Design Template

(New Options in Blue) **OCTANE**

703XY330 **Stem Bolt/Caster Pin Width** Measured from centerline of caster pin to the centerline of the opposite caster pin. 703XY331 RGK to Optimize: ☐

Current Chair: New Chair:

703XY231 **Backrest Height** Measured from top of seat tube at rear of frame to top of back-post.

Current Chair: New Chair:

703XY360 **Backrest Angle** Measured as the angle from top front of the backrest tube to the floor. (Max. -10° to +10°) ☐ -5° ☐ -3° ☐ 0° ☐ 3° ☐ 5°

Current Chair: New Chair:

☐ 703XG110 **Inward Taper** Measured from inside backrest tube to inside of opposite backrest tube at the top of the backrest. \$150

Current Chair: New Chair:

☐ 703XG100 **Ergonomic Backrest** Measured from top of seat tube at rear of frame to the start of the rearward bend of back-post. \$150

Current Chair: New Chair:

☐ 703XG120 **Outward Flare** Measured from inside backrest tube to inside of opposite backrest tube at the top of the backrest. \$150






Current Chair: New Chair:

703XY390 **Camber Angle** Place a straight edge across the outside of rear wheel and use an angle finder to show the degree of camber.

Current Chair:	☐ 0°	New Chair:	☐ 703XY396	0°
	☐ 1°		☐ 703XY397	1°
	☐ 2°		☐ 703XY398	2°
	☐ 3°		☐ 703XY399	3°
	☐ 5°		☐ 703XY400	5°

703XH165 **Center of Gravity** Measured from front of back post at the seat tube to the centerline of the rear axle.

Current Chair: New Chair:



Design Template



(New Options in Blue) **OCTANE**

703XY460 **Rear Wheel Gap** Measured from outer most point of chair (seat back post, side guard or armrest) to inside of rear tire.

Current Chair:	☐ .6"(15)	New Chair:	☐ 703XY464 .6"(15)
	☐ .78"(20)		☐ 703XY465 .78"(20)
	☐ 1"(25)		☐ 703XY466 1"(25)

☐ 703XR060 **Welded Calf Strap Mounting Loops** Measured from the ground to the top of where the calf strap is required. 703XR061 RGK to Optimize: ☐ \$154

Current Chair:			New Chair:		
----------------	----------------------	----------------------	------------	----------------------	----------------------



Footrest Design

☐	703XY100	Round Design	STD
☐	703XY105	Square Design	NCO

Footrest Styles

☐	703XB020	Tubular Titanium Footrest	STD
☐	703XB070	Angle adjustable aluminum footrest	\$213
☐	703XB080	Tubular footrest welded into frame (No height adjustment)	\$474
☐	703XB030	Flip up plastic footrest (Dual-tube frame design only)	\$363
☐	703XB150	Heavy duty tubular footrest	\$185

Footrest Cover

☐	793XB050	Lightweight 3k Carbon Fiber Footrest cover	\$255
☐	703XB155	Lightweight Spinerger CX6 Cast Carbon Fiber Footrest cover	\$252
☐	703XB040	Plastic Cover (Not available with 703XB070 & 703XB030)	\$101
☐	703XB170	Aluminum cover (Not available with 703XB070 & 703XB030)	\$101
☐	703XB165	Rubber grip tape (3M Safety walk. Not available with 703XB030)	\$89

Front Wheels

☐	703XC060	Aluminum Center 4"x1.3" (Black Hub only)	STD
☐	703XC110	Aluminum Center 4"x1" (Silver hub only)	\$239
☐	703XC135	Aluminum Center 4"x1.25"	\$239
	☐	703XC1131 Sliver hub & black tire	
	☐	703XC1132 Black hub & black tire	
	☐	703XC1133 Purple hub & black tire	
	☐	703XC1134 Green hub & black tire	
	☐	703XC1135 Red hub & black tire	
	☐	703XC1136 Blue hub & black tire	
☐	703XC155	4" x 1.25" Spinerger CX6 Cast Carbon fiber front wheel	\$125
	☐	703XC156 Cerakoted Spinerger CX6 Cast Carbon fiber front wheel (Hub only. See Cerakot color options on page 8)	\$242
☐	703XC070	Soft roll 5"x1.3"	NCO
☐	703XC100	Pneumatic Aluminum Center 5"x1"	\$255
☐	703XC120	Aluminum Center 5"x1"	\$239
☐	703XC140	Aluminum Center 5"x1.4" (Sliver hub & black tire)	\$239

Design Template

(New Options in Blue) **OCTANE**

☐	703XC145	Aluminum Center 6"x1.4" (Sliver hub & black tire)	\$239
☐	703XC150	Plastic Center 6"x1.4"soft roll (Black hub & grey tire)	\$239
Forks			
☐	703XC192	High Trail fork	NCO
☐	703XC300	Omobic Single Sided fork	\$250
☐	703XC210	Froglegs suspension fork	\$652
☐	703XC211	Black fork	STD
Wheel Lock			
☐	703XD020	Underseat Scissor	NCO
☐	703XD030	Compact Scissor	\$274
☐	703XD035	Compact Scissor with Direct Mount	\$363
☐	703XD050	Push To Lock	NCO
☐	703XD085	Eagle (Indicate color)	\$482
☐	703XD086	Eagle with Direct Mount (Indicate color)	\$523
☐	703XD088	Black	☐ 703XD089 Purple
☐	703XD090	Red	☐ 703XD091 Blue
☐	703XD094	Omobic all aluminum	\$149
☐	703XD095	Omobic all aluminum with Direct Mount	\$178
☐	703XD100	Push Lox Dynamic Hub Locking system	\$1,518
☐	703XD010	Omit Wheel Lock	NCO
Side Guards			
☐	703XE030	Screw on aluminum	\$337
☐	703XE080	Lift out aluminum	\$548
☐	703XE130	Fold flat aluminum	\$690
☐	703XE275	Screw on Spinergy CX6 Cast Carbon fiber	\$751
☐	703XE100	Lift Out 3k Carbon fiber	\$789
☐	703XE106	Lift Out Spinergy CX6 Cast Carbon fiber	\$1,045
☐	703XE150	Fold flat 3k carbon fiber	\$1,084
☐	703XE155	Fold flat Spinergy CX6 Cast Carbon fiber	\$1,084
☐	703XE070	Screw on ABS plastic	\$172
☐	703XE120	Lift out ABS plastic	\$502
☐	703XE170	Fold flat ABS plastic	\$855
☐	703XEM05	25mm above tire (Follow tire)	NCO
☐	703XEM06	Individual above tire (Follow tire - max. 25mm)	NCO
☐	703XEM02	Level with tire	STD
☐	703XEM07	Side guard length	STD
☐	703XEM08	Side guard height	STD
☐	703XE040	Screw on aluminum with fender	\$287
☐	703XE090	Lift out aluminum with fender	\$596
☐	703XE140	Fold flat aluminum with fender	\$728
☐	703XE060	Screw on 3k carbon fiber with fender	\$981
☐	703XE110	Lift out 3k carbon fiber with fender	\$1,166

Design Template

(New Options in Blue) **OCTANE**

☐	703XE160	Fold flat 3k carbon fiber with fender		\$1,285
☐	703XEM03	10mm above tire		NCO
☐	703XEM01	Individual above tire (Follow tire - max. 25mm)		NCO
☐	703XEM09	Center of tire		STD
☐	703XEM10	Outside of tire		NCO
☐	703XE190	Omit side guards		NCO
Armrest				
☐	703XF010	Lift out armrest with plastic sideguard	Height (Top of seatrail to top of armrest)	\$610
☐	703XF020	Lift out armrest with aluminum sideguard	Height (Top of seatrail to top of armrest)	\$481
☐	703XF130	Swing way tubular armrest with welded in receiver (Requires 703XH110 Square back adjustable center gravity)		\$715
☐		Height (Top of seatrail to top of armrest)	Tubular armrest length	
☐	703XF140	Omit armrest		NCO
Backrest Height				
☐	703XY230	Fixed height		STD
☐	703XY245	Adjustable height (Range: +/- 25mm)		\$255
Backrest Design				
☐	703XG010	Fixed backrest (non-folding)		STD
☐	703XG020	Folding backrest with 2 position snap lock and angle adjustable (Range: +/- 5°)		\$713
Push Handles				
☐	703XG070	Screw in push handles (Not available with 703XY245)		\$232
☐	703XG080	Flip down push handles (Not available with 703XY245)		\$454
☐	703XG085	Flip down push handles height adjustable (Not available with 703XY245)		\$742
☐	703XG090	Stroller push handles height adjustable (Measure height from floor to the top of the push handle)		\$454
☐		703XG098 RGK to Optimize: ☐		
☐	703XG097	Integrated push handles		\$696
☐	703XG030	Omit push handles		NCO
Backrest Options				
☐	703XG040	Heavy duty backrest upgrade (Not available with 703XY245)		\$213
☐	703XG130	Extra deep brace/rigidizer bar (Recommended for JAY backs)		\$118
Center of Gravity				
☐	703XH110	Square back adjustable of gravity (Range: 25mm - 150mm)		STD
☐	703XH120	Round back adjustable of gravity (Range: 60mm - 150mm)		NCO
☐	703XH100	Single post fixed center of gravity (Range: 25mm - 150mm)		\$749
☐	703XH090	V-Style fixed center of gravity (Range: 25mm - 150mm)		\$749
Axle Pin Options				
☐	703XK227	Standard axle pin		STD
☐	703XK220	Titanium 130mm quick release axle pin		\$165
☐	703XK240	Quad release axle pin		\$248

Design Template

(New Options in Blue) **OCTANE**

Anti-Tip

☐	703XJ010	Swing-away anti tip LEFT side (Only available with square back adjustable center of gravity. Not available with 703XH140)	\$172
☐	703XJ020	Swing-away anti tip RIGHT side (Only available with square back adjustable center of gravity. Not available with 703XH140)	\$172
☐	703XJ100	Single sport anti-tip (Centered mounted)	\$519

Frame and Footrest Finishes

Non painted finishes

☐	703XN030	Bead Blast frame	STD
☐	703XN050	Polished frame	\$1,372

Painted finishes

☐	703XN500	CeraKote frame	\$550
☐	703XN600	CeraKote components (Component color will match frame if not otherwise specified. Components include: Forks, Wheel Lock clamps, and Camber Tube clamps.)	\$550
☐	703XN800	CeraKote Spinergy wheels (See Rear Wheel Type on page 8. Hubs and Rims are CeraKoted on Rear Wheels, Hubs on Front Wheels.)	\$787

Frame	Components	Spinergy	Colors	Frame	Components	Spinergy	Colors
☐ 703XN501	☐ 703XN601	☐ 703XN801	Midnight Black	☐ 703XN509	☐ 703XN609	☐ 703XN809	Rose Gold
☐ 703XN502	☐ 703XN602	☐ 703XN802	Java Bronze	☐ 703XN510	☐ 703XN610	☐ 703XN810	Retro Blue
☐ 703XN503	☐ 703XN603	☐ 703XN803	Pure Titanium	☐ 703XN511	☐ 703XN611	☐ 703XN811	Hot Pink
☐ 703XN504	☐ 703XN604	☐ 703XN804	Fire Red	☐ 703XN512	☐ 703XN612	☐ 703XN812	Wicked Green
☐ 703XN505	☐ 703XN605	☐ 703XN805	Prime Blue	☐ 703XN513	☐ 703XN613	☐ 703XN813	Bright White
☐ 703XN506	☐ 703XN606	☐ 703XN806	Sahara Gold	☐ 703XN514	☐ 703XN614	☐ 703XN814	Tangerine Orange
☐ 703XN507	☐ 703XN607	☐ 703XN807	Black Cherry	☐ 703XN515	☐ 703XN615	☐ 703XN815	Sunbeam Yellow
☐ 703XN508	☐ 703XN608	☐ 703XN808	Fresh Green	☐ 703XN516	☐ 703XN616	☐ 703XN816	Matte Navy

Accent Color

☐	703XN251	Black	STD
☐	703XN253	Red	NCO
☐	703XN255	Purple	NCO
☐	703XN252	Silver	NCO
☐	703XN254	Blue	NCO
☐	703XN250	Individual Accent Color	\$328

Seat Upholstery

☐	703XP040	Seat Sling Adjustable Strap	STD
☐	703XR050	Carbon Fiber Seat Plate	\$494

Back Upholstery

☐	703XP035	Sub4 Airtech Adjustable Strap	STD
☐	703XP050	Airtech Adjustable Strap	NCO
	☐ 703XP070	Backrest upgrade to leather corners (Only with Airtech upholstery 703XP050)	\$125
	☐ 703XP110	Backrest zip pocket (Only with Airtech upholstery 703XP050)	\$131
	Upholstery Stitching Color		
☐	703XP550	Black	STD
☐	703XP551	Silver	NCO
☐	703XP552	Red	NCO
☐	703XP553	Orange	NCO
☐	703XP554	Blue	NCO

Design Template

(New Options in Blue) **OCTANE**

☐	703XP555	Purple			NCO
☐	703XP556	Individual	Pantone color number		NCO
		RGK Backrest Logo Color			
☐	703XP557	White (Standard if nothing selected)			STD
☐	703XP558	Red			STD
☐	703XP559	Black			NCO
☐	703XP560	Blue			NCO
☐	703XP561	Purple			NCO
☐	703XP570	Omit back upholstery			NCO
☐	703XP562	Individual	Pantone color number		NCO
Rear Wheel Size					
703XY420	☐	24" (540) Rear wheel size	☐	703XY430 25" (559) Rear wheel size	STD
Rear Wheel Type					
☐	703XK160	Spinergy ES Wire (20 spokes. CeraKote option available on Rim only. Not available with Spinergy anodized color options.)			STD
☐	703XK015	Proton (Not available with Spinergy anodized color options.)			STD
☐	703XK320	Spinergy XLX (12 spoke, cross radial spokes.)			\$1,099
☐	703XK050	Spinergy LX (12 spoke)			\$1,030
☐	703XK060	Spinergy Spox (18 spoke. Not available with Spinergy anodized color options or spoke colors, purple, orange, pink and green.)			\$685
☐	703XK070	Spinergy BLXL (18 bladed spoke.)			\$1,099
☐	703XK150	Spinergy CLX (18 bladed spokes, carbon fiber rim. Spinergy color option available on hub only. Not available with silver hub)			\$3,233
☐	703XK155	RGK EXCLUSIVE Spinergy CX6 Cast Carbon fiber (Cast Carbon Fiber Hub, Matte Black Carbon Fiber Rim and Black Bladed Spokes with 2 White Bladed Spokes Around Valve Stem. Not available with Spinergy anodized color options.)			\$3,700
		Spinergy spoke colors			
☐	703XK650	Black			STD
☐	703XK651	Red			NCO
☐	703XK652	Purple			NCO
☐	703XK653	White			NCO
☐	703XK654	Blue			NCO
☐	703XK655	Orange			NCO
☐	703XK656	Yellow			NCO
☐	703XK657	Pink			NCO
☐	703XK658	Green			NCO
☐	703XK180	Individual spoke design			\$219
		Spinergy anodized colors (Rims and Hubs are anodized not spokes, on Spinergy BLXL (703XK075), Hubs only on 703XK150)			
☐	703XK190	Black			STD
☐	703XK202	Purple			\$421
☐	703XK203	Blue			\$421
☐	703XK201	Red			\$421
☐	703XK204	Gold			\$421
☐	703XK205	Individual hub/rims			\$585

Supply Without Rear Wheels

☐	703K040	Omit Rear Wheels			NCO
		703XK041 Wheel Type:		703XK042 Tire Type:	

Handrims Style

☐	703XY485	Tab Length Wide			STD
☐	703XY480	Tab Length Narrow			NCO
		☐ 703XL090 Remove excess mounting tabs (Recommended for small wheel gaps)			\$150
☐	703XL010	Aluminum Hard Anodized 6 Tab			STD
☐	703XL160	Elipse Ergonomic 6 Tab			\$650
☐	703XL020	QGrip 6 Tab			\$412
☐	703XL040	Titanium 6 Tab			\$459
☐	703XL035	Natural-Fit LT	☐ 703XL036 Standard thumb ring	☐ 703XL037 Super Grip thumb ring	☐ 703XL038 Omit thumb ring
☐	703XL150	Carbolife Gekko			\$765
☐	703XL220	Carbolife Gekko S			\$819
					\$654
☐	703XL170	Carbolife Curve L			\$716
☐	703XL230	Carbolife Curve			\$569
☐	703XL180	Carbolife Curve L Tetra			\$912
☐	703XL240	Carbolife Curve Tetra			\$727
☐	703XL190	Carbolife Quadro			\$912
☐	703XL200	Spinergy TC			\$758
☐	703XL210	Spinergy TC Max			\$890

Tires

☐	703XM060	Schwalbe One (Lightest option)			\$137
☐	703XM150	IRC EXERACER			\$194
☐	703XM010	Kenda	☐ 703XM011 Grey	☐ 703XM012 Blue	☐ 703XM013 Red
					STD
☐	703XM020	Schwalbe Marathon Plus			\$101
☐	703XM040	Shox Solid Tire			\$165
☐	703XM070	Schwalbe Rocket Ron mountain bike tire			\$69
☐	703XM050	Kenda K885 Mountain Bike Tire			\$41

Inner Tubes

☐	703XY500	Schraeder Valve Inner Tube			STD
--------------------------	----------	-----------------------------------	--	--	-----

Cushions

☐	703XP080	Cushion foam			\$125
		Density	☐ 703XP001 Soft	☐ 703XP002 Medium	
		Thickness	☐ 703XP003 50mm	☐ 703XP004 75mm	☐ 703XP005 100mm

Cushion Covers

☐	703XP085	Nylon Cushion Cover (Black)			STD
☐	703XP090	Airlite Cushion Cover (Black)			\$101

Jay Products - Please Include Jay Order Form

- ☐ 703XP900 Jay Cushion Refer to Jay Order Form
- ☐ 703XP910 Jay Backrest Refer to Jay Order Form

Empulse Products - Please Include Empulse Order Form

- ☐ 703EP100 ☒ Empulse Add-On Power-Assist Devices (IN11- USEWITHRGK) ☐ Ship with chair or ☐ Ship seperate from chair (uninstalled)
- ☐ F35 - Pull Device
- ☐ M90 - Wheel Add-On
- ☐ R90 - Push-Assist

Accessories

- ☐ 703XR080 Frontwheel Assembly \$1,085
- ☐ 703XR240 Reverse Velcro on seating ☐ 703XR241 Double sided tape ☐ 703XR242 Sewn in \$138
- ☐ 703XP120 Zip Black Frame Protectors Neoprene Fasten (Black) \$87
- ☐ 703XP420 Velcro Black Frame Protectors Nylon Fasten (Black) \$83
- ☐ 703XR010 Black Seat Pouch Velcro Attach (Black) \$61
- ☐ 703XR200 Black Rucksack with RGK Logo (Black) \$156
- ☐ 703XR230 Spinergy Wheel Bag \$207
- ☐ 703XR210 Recon Mini Rider-Carry Air Compressor \$235
- ☐ 703XR090 Tool Kit STD
- ☐ 703XS060 Calf Strap STD

Ordering Information

Estimated MSRP*

0

Quote

Order

(*Price does not include any one off requested options)

Specification Notes

