

GEMINO 60 WALKER GEMINO 60 M WALKER

April 2025



Product sold with backrest in Canada.

ORDER FORM

Invoice address or customer number:

Customer number:

Name:

Street:

Town:

Postal code:

Country:

E-mail:

Tel:

Delivery address:

☐ Deliver to Invoice address

Customer number:

Name:

Street:

Town:

Postal code:

Country:

E-mail:

Tel:

Order information

☐ Order

☐ Quote

Order date:


Your reference:

Reimbursement code: ADP #: AA0001082AA1

Customer Service: cscanada@sunmed.com

Base price \$ 1,326.00

DIMENSIONS



	Gemino 60	Gemino 60 M
A Total height	98.5-115.5 cm	89.5-98.5 cm
B Total length	74 cm	74 cm
C Total width	64 cm	64 cm
D Width between push handles	32 cm	32 cm
E Seat height	62 cm	53 cm
Height folded	109 cm	101 cm
Length folded	74 cm	74 cm
Width folded	35.5 cm	35.5 cm
Turning circle	95 cm	95 cm
F Front wheels	28 x 4.2 cm	28 x 4.2 cm
G Rear wheels	23 x 4.2 cm	23 x 4.2 cm
Total weight without basket	11.8 kg	11.6 kg
H Weight basket, standard	0.3 kg	0.25 kg
Max. user weight	130 kg	130 kg
Max. weight in basket	10 kg	10 kg
Recommended user height	150-200 cm	135-170 cm

SIZES AND COLOURS		Price \$	Quantity	Total
Gemino 60 Walker				
☐ HC-7160020	Gemino 60 Walker, silver grey	1,326.00		
Gemino 60 M Walker				
☐ HC-7160220	Gemino 60 M Walker, silver grey	1,326.00		

ACCESSORIES		Price \$	Quantity	Total
Other accessories				
☐ HC-7109011	Crutch holder	106.00		
☐ HC-7109140	Table tray	106.00		
☐ HC-7109060	Umbrella, black	264.00		
☐ HC-7109050	Basket for oxygen cylinder - Gemino 60 (not for Gemino 60 M)	264.00		
☐ HC-7169080	Slowdown brake, set	252.00		
☐ HC-7109095	Light	129.00		
☐ HC-7109100	Infusion holder	264.00		

GEMINO 60 WALKER / GEMINO 60 M WALKER			Price \$
Base price			
Option price			
Total price			

GEMINO ACCESSORIES

GEMINO

OTHER ACCESSORIES



HC-7109011 **\$106.00**
Crutch holder



HC-7109140 **\$106.00**
Table tray



HC-7109095 **\$129.00**
Light



HC-7109060 **\$264.00**
Umbrella, black



HC-7109050 & HC-7109051 **\$264.00**
Basket for oxygen cylinder*



HC-7109100 **\$264.00**
Infusion holder



Sunrise Medical
1000 Creditstone Road, Unit 2
Concord, ON
L4K 4P8, Canada
Tel. +1.800.263.3390
E-mail: cscanada@sunmed.com
www.sunrisemedical.ca

*Not available for Gemino 60 M/60 M Walker.