

GEMINO 60 WALKER
GEMINO 60 M WALKER

February 2024



Product sold with backrest in Canada.

ORDER FORM

Invoice address or customer number:
 Customer number:
 Name:
 Street:
 Town: Postal code:
 Country:
 E-mail:
 Tel:


Delivery address: Deliver to Invoice address
 Customer number:
 Name:
 Street:
 Town: Postal code:
 Country:
 E-mail:
 Tel:

Order information
 Order Quote
 Order date:
 Your reference:
 Reimbursement code: ADP #: AA0001082AA1

Customer Service: cscanada@sunmed.com

Base price \$ 1,268.00

DIMENSIONS



	Gemino 60	Gemino 60 M
A Total height	98.5-115.5 cm	89.5-98.5 cm
B Total length	74 cm	74 cm
C Total width	64 cm	64 cm
D Width between push handles	32 cm	32 cm
E Seat height	62 cm	53 cm
Height folded	109 cm	101 cm
Length folded	74 cm	74 cm
Width folded	35.5 cm	35.5 cm
Turning circle	95 cm	95 cm
F Front wheels	28 x 4.2 cm	28 x 4.2 cm
G Rear wheels	23 x 4.2 cm	23 x 4.2 cm
Total weight without basket	11.8 kg	11.6 kg
H Weight basket, standard	0.3 kg	0.25 kg
Max. user weight	130 kg	130 kg
Max. weight in basket	10 kg	10 kg
Recommended user height	150-200 cm	135-170 cm

SIZES AND COLOURS		Price \$	Quantity	Total
Gemino 60 Walker				
<input type="checkbox"/> HC-7160020	Gemino 60 Walker, silver grey	1,268.00		
Gemino 60 M Walker				
<input type="checkbox"/> HC-7160220	Gemino 60 M Walker, silver grey	1,268.00		

ACCESSORIES		Price \$	Quantity	Total
Other accessories				
<input type="checkbox"/> HC-7109011	Crutch holder	97.00		
<input type="checkbox"/> HC-7109140	Table tray	97.00		
<input type="checkbox"/> HC-7109060	Umbrella, black	240.00		
<input type="checkbox"/> HC-7109050	Basket for oxygen cylinder - Gemino 60 (not for Gemino 60 M)	240.00		
<input type="checkbox"/> HC-7169080	Slowdown brake, set	426.00		
<input type="checkbox"/> HC-7109095	Light	117.00		
<input type="checkbox"/> HC-7109100	Infusion holder	240.00		

GEMINO 60 WALKER / GEMINO 60 M WALKER		Price \$
Base price		
Option price		
Total price		



Sunrise Medical
 237 Romina Drive, Unit 3
 Concord, ON
 L4K 4V3, Canada
 Tel. +1.800.263.3390
 Fax: +1.800.561.5834
 E-mail: cscanada@sunmed.com
www.sunrisemedical.ca



OTHER ACCESSORIES



// 7109011 \$97.00
Crutch holder



// 7109140 \$97.00
Table tray



// 7109095 \$117.00
Light



// 7109060 \$240.00
Umbrella, black



// 7109050 & 7109051 \$240.00
Basket for oxygen cylinder*



// 7109100 \$240.00
Infusion holder

*Not available for Gemino 60 M/60 M Walker.