



JAY® X2 Cushion



February 2024

Mark For: _____
 Date: _____
 Dealer Acct #: _____
 Dealer: _____
 Dealer Contact: _____
 Dealer Address: _____
 Dealer City: _____ Prov: _____ PC: _____
 Dealer Phone: () _____ Fax: () _____
 Confirmation Email: _____
 Confirm Via: Fax Email

Submitting for: Quote Order
 PO#: _____
ADDITIONAL SHIPPING INFORMATION
 Ship To: _____
 Attention: _____
 Address: _____
 Address: _____
 Ship To City: _____ Prov: _____ PC: _____
 Ship To Phone: () _____ Fax: () _____

JAY X2 Cushion

JAY X2 CUSHION Less than 22" in Width (300lb weight capacity)

✓	Option #	HCPCS	Definition	Actual Width	Actual Depth	Price
<input type="checkbox"/>	JX21414	E2622	JAY X2 Cushion 14" X 14"	13.75"	14"	\$930.00
<input type="checkbox"/>	JX21416	E2622	JAY X2 Cushion 14" X 16"	13.75"	16"	\$930.00
<input type="checkbox"/>	JX21515	E2622	JAY X2 Cushion 15" X 15"	15.25"	15"	\$930.00
<input type="checkbox"/>	JX21516	E2622	JAY X2 Cushion 15" X 16"	15.25"	16"	\$930.00
<input type="checkbox"/>	JX21616	E2622	JAY X2 Cushion 16" X 16"	15.75"	16"	\$930.00
<input type="checkbox"/>	JX21617	E2622	JAY X2 Cushion 16" X 17"	15.75"	17"	\$930.00
<input type="checkbox"/>	JX21618	E2622	JAY X2 Cushion 16" X 18"	15.75"	18"	\$930.00
<input type="checkbox"/>	JX21620	E2622	JAY X2 Cushion 16" X 20"	15.75"	20"	\$930.00
<input type="checkbox"/>	JX21716	E2622	JAY X2 Cushion 17" X 16"	17.25"	16"	\$930.00
<input type="checkbox"/>	JX21717	E2622	JAY X2 Cushion 17" X 17"	17.25"	17"	\$930.00
<input type="checkbox"/>	JX21816	E2622	JAY X2 Cushion 18" X 16"	17.75"	16"	\$930.00
<input type="checkbox"/>	JX21818	E2622	JAY X2 Cushion 18" X 18"	17.75"	18"	\$930.00
<input type="checkbox"/>	JX21820	E2622	JAY X2 Cushion 18" X 20"	17.75"	20.25"	\$930.00
<input type="checkbox"/>	JX22020	E2622	JAY X2 Cushion 20" X 20"	19.75"	20"	\$930.00

JAY X2 CUSHION with Extra Cover - Less than 22" in Width (300lb weight capacity)

<input type="checkbox"/>	JX21414-2	E2622	JAY X2 Cushion 14" X 14"	13.75"	14"	\$988.00
<input type="checkbox"/>	JX21416-2	E2622	JAY X2 Cushion 14" X 16"	13.75"	16"	\$988.00
<input type="checkbox"/>	JX21515-2	E2622	JAY X2 Cushion 15" X 15"	15.25"	15"	\$988.00
<input type="checkbox"/>	JX21516-2	E2622	JAY X2 Cushion 15" X 16"	15.25"	16"	\$988.00
<input type="checkbox"/>	JX21616-2	E2622	JAY X2 Cushion 16" X 16"	15.75"	16"	\$988.00
<input type="checkbox"/>	JX21617-2	E2622	JAY X2 Cushion 16" X 17"	15.75"	17"	\$988.00
<input type="checkbox"/>	JX21618-2	E2622	JAY X2 Cushion 16" X 18"	15.75"	18"	\$988.00
<input type="checkbox"/>	JX21620-2	E2622	JAY X2 Cushion 16" X 20"	15.75"	20"	\$988.00
<input type="checkbox"/>	JX21716-2	E2622	JAY X2 Cushion 17" X 16"	17.25"	16"	\$988.00
<input type="checkbox"/>	JX21717-2	E2622	JAY X2 Cushion 17" X 17"	17.25"	17"	\$988.00
<input type="checkbox"/>	JX21816-2	E2622	JAY X2 Cushion 18" X 16"	17.75"	16"	\$988.00
<input type="checkbox"/>	JX21818-2	E2622	JAY X2 Cushion 18" X 18"	17.75"	18"	\$988.00
<input type="checkbox"/>	JX21820-2	E2622	JAY X2 Cushion 18" X 20"	17.75"	20.25"	\$988.00
<input type="checkbox"/>	JX22020-2	E2622	JAY X2 Cushion 20" X 20"	19.75"	20"	\$988.00

Please visit www.sunparts.us for spare parts information.

JAY YOUR WAY MODIFICATIONS

Note: To order sizes offered via JAY Your Way add a "M" to the front of the part number and follow the prompts.

Active Fluid Bladder

AFB Active Fluid Bladder NC
 Notes: 40% Larger fluid bladder that provides fluid support to the outside edges of the cushion.

Fluid Volume Change

OF Fluid Volume Overfill NC
 5% 10% 15% 20% 25% 30% 35% 40%

UF Fluid Volume Underfill NC
 -5% -10% -15% -20%

Obliquity Fluid Volume Change

OOFR Obliquity Overfill - Right Side NC
 5% 10% 15% 20% 25% 30% 35% 40%

OOFL Obliquity Overfill - Left Side NC
 5% 10% 15% 20% 25% 30% 35% 40%

OUFRR Obliquity Underfill - Right Side NC
 -5% -10% -15% -20%

OUFLR Obliquity Underfill - Left Side NC
 -5% -10% -15% -20%

Cover Options

CVR Cover Options NC
 X-static w/ Stretch Incontinent Microclimatic


JAY X2 Size Matrix

		WIDTH											
		12	13	14	15	16	17	18	19	20	21	22	
DEPTH	12	N/A	JYW	JYW	JYW	JYW	JYW	JYW	N/A	N/A	N/A	N/A	
	13	N/A	JYW	JYW	JYW	JYW	JYW	JYW	N/A	N/A	N/A	N/A	
	14	N/A	JYW		JYW	JYW	JYW	JYW	N/A	N/A	N/A	N/A	
	15	N/A	JYW	JYW		JYW	JYW	JYW	N/A	N/A	N/A	N/A	
	16	N/A	JYW						JYW	JYW	N/A	N/A	
	17	N/A	N/A	N/A	JYW				JYW	JYW	JYW	N/A	N/A
	18	N/A	N/A	N/A	JYW				JYW	JYW	N/A	N/A	
	19	N/A	N/A	N/A	JYW	JYW	JYW	JYW	JYW	JYW	N/A	N/A	
	20	N/A	N/A	N/A	JYW	JYW	JYW	JYW	JYW		N/A	N/A	
	21	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	
	22	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	

= Standard Offering
 JYW = Offered Via Jay Your Way
 N/A = Not available

JAY X2 Measuring Guide

Measurements made from completed cushion with covers.



Width

X = Width

Note: Measurement made from rear bottom edge to bottom edge of the cushion.
Ex: JX21618 X=16"

Depth

Y = Depth

Note: Measurement made from front edge of rear cutout and front lower of cushion.
Ex: JX21618 Y = 18"

Size Change

RL Reduce Length NC
 -1" -2" -3" -4"

RW Reduce Width NC
 -1"

Positioning Components

APC Add Positioning Components NC

Notes: Add 1" tall Lateral, Medial, and Pelvic Positioning Pieces.

Don't see the option you want? NO PROBLEM! Write your request below and we can evaluate it and provide a quote.

Notes:



Sunrise Medical Canada Inc
237 Romina Drive, Unit 3 • Concord, ON L4K 4V3 • Canada