



ORDER FORM

Everyday Activity Seat

Activity chair for children or young adults with moderate needs who still require pelvic stability during the day.

Customer			
Account			
Account #		Date	
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #	
Buyer			
Mark For			

Shipping Address			
Name			
Address			
City		Prov	
PC		Phone	
Email / Fax			

Everyday Activity Seat

Size 1 Weight Capacity: 66 lbs | Size 2 Weight Capacity: 110 lbs | Size 3 Weight Capacity: 220 lbs

✓ OPTION	DESCRIPTION	RETAIL
Base Frame		\$
Selection Required.		
<input type="checkbox"/> 151-6102	Size 1 Seat Base, Backrest and Hi-Low Chassis with Gas Spring	\$3,485
<input type="checkbox"/> 151-6202	Size 2 Seat Base, Backrest and Hi-Low Chassis with Gas Spring	\$3,588
<input type="checkbox"/> 151-6302	Size 3 Seat Base, Backrest and Hi-Low Hydraulic Assist Chassis	\$7,073
Chair Cover		\$
Selection Required. Must match base frame size.		
<input type="checkbox"/> 151-6104-03VIN	Size 1 Vinyl Seat and Back Covers - Blue	\$403
<input type="checkbox"/> 151-6206-03VIN	Size 2 Vinyl Seat and Back Covers - Blue	\$460
<input type="checkbox"/> 151-6304-03VIN	Size 3 Vinyl Seat and Back Covers - Blue	\$762
Lateral Support		\$
Must match base frame size.		
<input type="checkbox"/> 151-6106-07VIN	Size 1 Rigid Chest Laterals - Black Vinyl (Pair)	\$428
<input type="checkbox"/> 151-6210-07VIN	Size 2 Rigid Chest Laterals - Black Vinyl (Pair)	\$428
<input type="checkbox"/> 151-6306-07VIN	Size 3 Rigid Chest Laterals - Black Vinyl (Pair)	\$576
<input type="checkbox"/> 151-6108-07VIN	Size 1 Swing-Away Chest Laterals - Black Vinyl (Pair)	\$529
<input type="checkbox"/> 151-6212-07VIN	Size 2 Swing-Away Chest Laterals - Black Vinyl (Pair)	\$529
<input type="checkbox"/> 151-6308-07VIN	Size 3 Swing-Away Chest Laterals - Black Vinyl (Pair)	\$708
<input type="checkbox"/> 151-6218-07VIN	Size 1 and 2 Seat Protraction Pads - Black Vinyl (Pair)	\$519
<input type="checkbox"/> 151-6314-07VIN	Size 3 Seat Protraction Pads - Black Vinyl (Pair)	\$555
<input type="checkbox"/> 137-808-05	Chest Support - Size 1 & 2 (Requires rigid or swing-away laterals)	\$201
<input type="checkbox"/> 133-625	Chest Support - Size 3 (Requires rigid or swing-away laterals)	\$201
Harness		\$
<input type="checkbox"/> 155-N001-L600	Butterfly Harness (Small)	\$409
<input type="checkbox"/> 155-N001-L601	Butterfly Harness (Medium)	\$409
<input type="checkbox"/> 155-N001-L602	Butterfly Harness (Large)	\$409
Upper Extremity Positioning		\$
Must match base frame size. Tray requires armrests.		
<input type="checkbox"/> 151-6117	Size 1 Armrests (Pair)	\$466
<input type="checkbox"/> 151-6226	Size 2 Armrests (Pair)	\$466
<input type="checkbox"/> 151-3747	Size 3 Armrests (Pair)	\$466
<input type="checkbox"/> 151-1644	Size 1 Tray - Plastic Gray	\$393
<input type="checkbox"/> 151-2644	Size 2 Tray - Plastic Gray	\$393
<input type="checkbox"/> 151-3767	Size 3 Tray - Plastic Black	\$393

Everyday Activity Seat

✓ OPTION	DESCRIPTION	RETAIL
Head Supports		\$
<input type="checkbox"/> 151-6208-07VIN	Size 1 and 2 Flat Headrest Support and Cushion - Black Vinyl Must match base frame size	\$364
<input type="checkbox"/> 151-6209-07VIN	Size 1 and 2 Laterals for Flat Headrest - Black Vinyl (Pair) Must match base frame size. Includes width adjustable lateral pads & covers. Requires Flat Headrest Support and Cushion.	\$309
<input type="checkbox"/> 151-3751-05	Size 3 Contoured Headrest - Black Neoprene Must match base frame size	\$535
<input type="checkbox"/> 151-1756	Size 1 and 2 Contoured Headrest Adapter Bracket	\$84
<input type="checkbox"/> 120-859	Size 1 and 2 Contoured Headrest Support Hardware	\$390
<input type="checkbox"/> 137-800-05	Size 1 and 2 Contoured Headrest Cushion - Black	\$185
Foot Positioning		\$
Footplates are optional. Footplate selection required with Sandals		
<input type="checkbox"/> 151-1745	Size 1 One Piece Footplate - Metal	\$523
<input type="checkbox"/> 151-2745	Size 2 One Piece Footplate - Metal	\$523
<input type="checkbox"/> 151-3745	Size 3 One Piece Footplate - Metal	\$523
<input type="checkbox"/> 152-1600	Small Sandals w/ Straps (Pair)	\$269
<input type="checkbox"/> 152-2600	Medium Sandals w/ Straps (Pair)	\$287
<input type="checkbox"/> 152-3600	Large Sandals w/ Straps (Pair)	\$310
<input type="checkbox"/> 151-1760	QTY Size 1 and 2 Sandal Raiser 1" (Each) Requires Sandal Selection, must match Sandal size.	\$62
<input type="checkbox"/> 151-3760	QTY Size 3 Sandal Raiser 1" (Each) Requires Sandal Selection, must match Sandal size.	\$69
Pelvic Fixation		\$
<input type="checkbox"/> 151-1660	Size 1 Pelvic Belt 4-Point	\$162
<input type="checkbox"/> 151-2660	Size 2 Pelvic Belt 4-Point	\$174
<input type="checkbox"/> 151-3660	Size 3 Pelvic Belt 4-Point	\$185
<input type="checkbox"/> 151-6115-07VIN	Size 1 Side Pads - Black Vinyl (Pair)	\$566
<input type="checkbox"/> 151-6222-07VIN	Size 2 Side Pads - Black Vinyl (Pair)	\$672
<input type="checkbox"/> 151-6317-07VIN	Size 3 Side Pads - Black Vinyl (Pair)	\$857
<input type="checkbox"/> 151-6224-07VIN	Size 2 Adduction Side Pad - Black Vinyl	\$319
<input type="checkbox"/> 151-6320-07VIN	Size 3 Adduction Side Pad - Black Vinyl	\$319
<input type="checkbox"/> 151-6120-07VIN	Size 1 Wedge Shape Pommel - Black Vinyl	\$241
<input type="checkbox"/> 151-6229-07VIN	Size 2 Wedge Shape Pommel - Black Vinyl	\$297
<input type="checkbox"/> 151-6334-07VIN	Size 3 Wedge Shape Pommel - Black Vinyl	\$303

Technical Specifications - Imperial

Frame Sizing Specification	Size 2	Size 3
Max Load	110 lbs	220 lbs
Weight (Seat and Hi-Low Chassis)	75 lbs	108 lbs
Seat Height (Hi-Low Chassis)	12.6" - 23.6"	17.7" - 23.6"
Seat Width	10.2" - 15"	13.4" - 19.7"
Seat Depth	9" - 15"	13" - 18.9"
Lateral Distance (Chest)	6.7" - 9.4"	9.4" - 14.6"
Back Height	15.7" - 20.5"	17.7" - 24.8"
Back Angle	30° Supine / 10° Prone	30° Supine / 10° Prone
Armrest Height	7.9" - 11.8"	7" - 12.7"
Lower Leg Length	9.8" - 17.7"	9.8" - 17.7"
Lower Leg Angle	-10° to +50°	-10° to +50°
Footplate Angle	-15° to +10°	-15° to +10°
Tilt-In-Space	30° Supine / 10° Prone	30° Supine / 10° Prone

Technical Specifications - Metric

Frame Sizing Specification	Size 2	Size 3
Max Load	50 kgs	100 kgs
Weight (Seat and Hi-Low Chassis)	34 kgs	49 kgs
Seat Height (Hi-Low Chassis)	320 mm - 600 mm	450 mm - 600 mm
Seat Width	260 mm - 380 mm	340 mm - 500 mm
Seat Depth	230 mm - 380 mm	330 mm - 480 mm
Lateral Distance (Chest)	170 mm - 240 mm	240 mm - 370 mm
Back Height	400 mm - 520 mm	450 mm - 630 mm
Back Angle	30° Supine / 10° Prone	30° Supine / 10° Prone
Armrest Height	200 mm - 300 mm	180 mm - 325 mm
Lower Leg Length	250 mm - 450 mm	250 mm - 450 mm
Lower Leg Angle	-10° to +50°	-10° to +50°
Footplate Angle	-15° to +10°	-15° to +10°
Tilt-In-Space	30° Supine / 10° Prone	30° Supine / 10° Prone



Sunrise Medical Canada Inc.
237 Romina Drive, Unit 3 • Concord, ON L4K 4V3 • Canada

Customer Service: (800) 263-3390 E-mail: cscanada@sunmed.com www.sunrisemedical.ca
 Specifications subject to change without notice. © 2024 Sunrise Medical Canada Inc. MK-129846-CAN Rev 11

