



ORDER FORM

Everyday Activity Seat EILK03

Activity chair for children or young adults with moderate needs who still require pelvic stability during the day.

Customer			
Account			
Account #		Date	
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #	
Buyer			
Mark For			

Shipping Address			
Name			
Address			
City		Prov	
PC		Phone	
Email / Fax			

Everyday Activity Seat

Size 1 Weight Capacity: 66 lbs | Size 2 Weight Capacity: 110 lbs | Size 3 Weight Capacity: 220 lbs

✓ OPTION	DESCRIPTION	RETAIL
Base Frame		\$
Selection Required.		
<input type="checkbox"/> 151BF04	Size 1 Seat Base, Backrest and Hi-Low Chassis with Gas Spring	\$3,485
<input type="checkbox"/> 151BF05	Size 2 Seat Base, Backrest and Hi-Low Chassis with Gas Spring	\$3,588
<input type="checkbox"/> 151BF06	Size 3 Seat Base, Backrest and Hi-Low Hydraulic Assist Chassis	\$7,073
Chair Cover		\$
Selection Required. Must match base frame size.		
<input type="checkbox"/> 151IU16	Size 1 Vinyl Seat and Back Covers - Blue	\$403
<input type="checkbox"/> 151IU17	Size 2 Vinyl Seat and Back Covers - Blue	\$460
<input type="checkbox"/> 151IU18	Size 3 Vinyl Seat and Back Covers - Blue	\$762
Lateral Support		\$
Must match base frame size.		
<input type="checkbox"/> 151ZY19	Size 1 Rigid Chest Laterals - Black Vinyl (Pair)	\$428
<input type="checkbox"/> 151ZY20	Size 2 Rigid Chest Laterals - Black Vinyl (Pair)	\$428
<input type="checkbox"/> 151ZY21	Size 3 Rigid Chest Laterals - Black Vinyl (Pair)	\$576
<input type="checkbox"/> 151ZY22	Size 1 Swing-Away Chest Laterals - Black Vinyl (Pair)	\$529
<input type="checkbox"/> 151ZY23	Size 2 Swing-Away Chest Laterals - Black Vinyl (Pair)	\$529
<input type="checkbox"/> 151ZY24	Size 3 Swing-Away Chest Laterals - Black Vinyl (Pair)	\$708
<input type="checkbox"/> 151KE47	Size 1 and 2 Seat Protraction Pads - Black Vinyl (Pair)	\$519
<input type="checkbox"/> 151KE48	Size 3 Seat Protraction Pads - Black Vinyl (Pair)	\$555
<input type="checkbox"/> 151ZY25	Chest Support - Size 1 & 2 (Requires rigid or swing-away laterals)	\$201
<input type="checkbox"/> 151ZY26	Chest Support - Size 3 (Requires rigid or swing-away laterals)	\$201
Harness		\$
<input type="checkbox"/> 151NH105	Butterfly Harness (Small)	\$409
<input type="checkbox"/> 151NH106	Butterfly Harness (Medium)	\$409
<input type="checkbox"/> 151NH107	Butterfly Harness (Large)	\$409
Upper Extremity Positioning		\$
Must match base frame size. Tray requires armrests.		
<input type="checkbox"/> 151AA167	Size 1 Armrests (Pair)	\$360
<input type="checkbox"/> 151AA168	Size 2 Armrests (Pair)	\$360
<input type="checkbox"/> 151AA169	Size 3 Armrests (Pair)	\$466
<input type="checkbox"/> 151A217	Size 1 Tray - Plastic Gray	\$393
<input type="checkbox"/> 151A218	Size 2 Tray - Plastic Gray	\$393
<input type="checkbox"/> 151A219	Size 3 Tray - Plastic Black	\$390



Flip Away Laterals



Flip Away Laterals and Chest Support



Protraction Pads



Armrest

Everyday Activity Seat

✓ OPTION	DESCRIPTION	RETAIL
Head Supports		\$
<input type="checkbox"/> 151Y053	Size 1 and 2 Flat Headrest Support and Cushion - Black Vinyl Must match base frame size	\$364
<input type="checkbox"/> 151Y052	Size 1 and 2 Laterals for Flat Headrest - Black Vinyl (Pair) Must match base frame size. Includes width adjustable lateral pads & covers. Requires Flat Headrest Support and Cushion.	\$309
<input type="checkbox"/> 151Y055	Size 3 Contoured Headrest - Black Neoprene Must match base frame size	\$535
<input type="checkbox"/> 151Y056	Size 1 and 2 Contoured Headrest Adapter Bracket	\$84
<input type="checkbox"/> 151AW24	Size 1 and 2 Contoured Headrest Support Hardware	\$390
<input type="checkbox"/> 151Y027	Size 1 and 2 Contoured Headrest Cushion - Black	\$185
Foot Positioning		\$
<input type="checkbox"/> 151F130	Size 1 One Piece Footplate - Metal	\$523
<input type="checkbox"/> 151F131	Size 2 One Piece Footplate - Metal	\$523
<input type="checkbox"/> 151F132	Size 3 One Piece Footplate - Metal	\$523
<input type="checkbox"/> 1518A107	Small Sandals w/ Straps (Pair)	\$269
<input type="checkbox"/> 1518A108	Medium Sandals w/ Straps (Pair)	\$287
<input type="checkbox"/> 1518A109	Large Sandals w/ Straps (Pair)	\$310
<input type="checkbox"/> 1519A25	Size 1 and 2 Sandal Raiser 1" (Each) Requires Sandal Selection, must match Sandal size.	\$62
<input type="checkbox"/> 1519A26	Size 3 Sandal Raiser 1" (Each) Requires Sandal Selection, must match Sandal size.	\$69
Pelvic Fixation		\$
<input type="checkbox"/> 151ZZ36	Size 1 Pelvic Belt 4-Point	\$162
<input type="checkbox"/> 151ZZ37	Size 2 Pelvic Belt 4-Point	\$174
<input type="checkbox"/> 151ZZ38	Size 3 Pelvic Belt 4-Point	\$185
<input type="checkbox"/> 151KE49	Size 1 Side Pads - Black Vinyl (Pair)	\$566
<input type="checkbox"/> 151KE50	Size 2 Side Pads - Black Vinyl (Pair)	\$672
<input type="checkbox"/> 151KE51	Size 3 Side Pads - Black Vinyl (Pair)	\$857
<input type="checkbox"/> 151KE52	Size 2 Adduction Side Pad - Black Vinyl	\$319
<input type="checkbox"/> 151KE53	Size 3 Adduction Side Pad - Black Vinyl	\$319
<input type="checkbox"/> 1518A111	Size 1 Wedge Shape Pommel - Black Vinyl	\$241
<input type="checkbox"/> 1518A112	Size 2 Wedge Shape Pommel - Black Vinyl	\$297
<input type="checkbox"/> 1518A113	Size 3 Wedge Shape Pommel - Black Vinyl	\$303



Flat Headrest



Contoured Headrest



Footplate



Sandals



Pelvic Belt 4-Point



Wedge Shape Pommel



Side Pads

Technical Specifications - Imperial

Frame Sizing Specification	Size 1	Size 2	Size 3
Max Load	66 lbs	110 lbs	220 lbs
Weight (Seat and Hi-Low Chassis)	65 lbs	75 lbs	108 lbs
Seat Height (Hi-Low Chassis)	9.8" - 23.6"	12.6" - 23.6"	17.7" - 23.6"
Seat Width	6.3" - 11"	10.2" - 15"	13.4" - 19.7"
Seat Depth	5.9" - 11"	9" - 15"	13" - 18.9"
Lateral Distance (Chest)	5.9" - 8.7"	6.7" - 9.4"	9.4" - 14.6"
Back Height	11.8" - 15.7"	15.7" - 20.5"	17.7" - 24.8"
Back Angle	30° Supine / 10° Prone	30° Supine / 10° Prone	30° Supine / 10° Prone
Armrest Height	5.9" - 9.8"	7.9" - 11.8"	7" - 12.7"
Lower Leg Length	5.9" - 9.8"	9.8" - 17.7"	9.8" - 17.7"
Lower Leg Angle	-10° to +50°	-10° to +50°	-10° to +50°
Footplate Angle	-15° to +10°	-15° to +10°	-15° to +10°
Tilt-In-Space	30° Supine / 10° Prone	30° Supine / 10° Prone	30° Supine / 10° Prone

Technical Specifications - Metric

Frame Sizing Specification	Size 1	Size 2	Size 3
Max Load	30 kg	50 kg	100 kg
Weight (Seat and Hi-Low Chassis)	29.5 kg	34 kg	49 kg
Seat Height (Hi-Low Chassis)	250 mm - 600 mm	320 mm - 600 mm	450 mm - 600 mm
Seat Width	160 mm - 280 mm	260 mm - 380 mm	340 mm - 500 mm
Seat Depth	150 mm - 280 mm	230 mm - 380 mm	330 mm - 480 mm
Lateral Distance (Chest)	150 mm - 220 mm	170 mm - 240 mm	240 mm - 370 mm
Back Height	300 mm - 400 mm	400 mm - 520 mm	450 mm - 630 mm
Back Angle	30° Supine / 10° Prone	30° Supine / 10° Prone	30° Supine / 10° Prone
Armrest Height	150 mm - 250 mm	200 mm - 300 mm	180 mm - 325 mm
Lower Leg Length	150 mm - 250 mm	250 mm - 450 mm	250 mm - 450 mm
Lower Leg Angle	-10° to +50°	-10° to +50°	-10° to +50°
Footplate Angle	-15° to +10°	-15° to +10°	-15° to +10°
Tilt-In-Space	30° Supine / 10° Prone	30° Supine / 10° Prone	30° Supine / 10° Prone

Sunrise Medical Canada Inc.
 1000 Creditstone Road, Unit 2
 Concord, Ont. L4K 4P8



MK-129846-CAN © 07/2024

