



Serial No:



ORDER FORM			
Customer Dealer Information		Shipping Address	Identical To Billing Address
Dealer		Client	
Client		Address	
Address			
E-Mail:		E-Mail:	
Phone:		Phone:	
Order Information		User Information	
Purchase Order No.:		Current Wheelchair:	
Special Requests (INTERNAL USE ONLY):		User Disability:	
		User Weight:	kg User height
Build4me			

Product Style

GSAXZ010

Grandslam

\$7,862.00

Max User Weight 125kg (275lbs)

PLEASE NOTE: ALL RGK CHAIRS ARE HAND BUILT AND HAVE A MANUFACTURING MEASUREMENT TOLERANCE OF +/- 5MM

CAD Approval

XW009	Quote only (No CAD For Approval)	
XW010	Order CAD Drawing For Approval (Includes 2 CAD Amendments)	\$300.00
XW011	Order (No CAD Required)	
XW012	Quote & CAD For Approval (Includes 2 CAD Amendments)	\$550.00

XY010	Seat Width (Inside to inside of sideguards)	
	mm (Max - 500mm)	
XY020	Seat Length	
	mm (Max - 500mm)	
XY030	Extra Seat Length (Provides additional leg support & hand grip)	
	mm	
XY050	Front frame height	(Min - 360mm) (Max - 630mm)
	mm	
XY060	Rear frame height	(Min - 360mm) (Max - 630mm)
	mm	
XY070	Frame Length	(Min = XY020 + XY030 + 140mm)
	mm	
XA020	Ergonomic seat	(Standard = 120mm)
	Specify Length	
	mm	\$125.00

Front Frame and Seat

Front Frame Options

XY080	Seat taper design	XA130	V-Seat design	\$316.00
Specify taper width		Specify taper width		
mm		mm		
		NOTE: Inside to inside wheel (XY470) required		
XA136	Parallel taper design (Not with seat taper XY080)			\$93.00
Specify taper width				
mm				
XA131	V-seat and taper design			\$341.00
Specify inside width at front of sideguard (A)	Specify taper width (B)			
mm	mm			
	NOTE: Inside to inside wheel (XY470) required			

Knee Support Options

XA200	Knee support ratchet system (Tubes on the outside)	No Cost
Specify length		
mm	Optimise (Smallest)	
XA202	Plastic shin pads (Includes 5mm hard padding)	No Cost
XA203	Carbon fibre shin pads (Includes 5mm hard padding)	\$762.00
XA204	Extra padding + cover for shin pads	\$186.00
Specify padding thickness	Specify padding density	
mm	Med Hard	

Seat Options

XR040	Aluminium seat plate (riveted to frame)	\$199.00
Specify length of plate	(Std XY020 - 20mm)	
mm		
XR050	Carbon fibre seat plate (riveted to frame)	\$494.00
Specify length of plate	(Std XY020 - 20mm)	
mm		

Footrest

GSAXB080 Fixed position footrest **Std**

GSAXB130 Adjustable position footrest (50mm range) **\$316.00**

Front footrest height

mm (min = 35mm fixed/40mm adjustable)

Rear footrest height

mm (min = 35mm fixed/40mm adjustable)

Footrest position

(front of footplate to castor pot)

mm

Footrest size

Footrest width

Not ABS

mm (Std=200mm)

Footrest length

mm (Std=250mm)

Fixed Footrest Plate Options

XB180 Aluminium footplate (2.5mm) **No Cost**

XB190 Carbon fibre footplate **\$577.00**

XB100 ABS Moulded foot tray **\$188.00**

Specify length

Max = 380mm (Full)
mm Std = 330mm

GSAXB140 Spacer kit **\$131.00**

Spacer kit

(Allows for 25mm increase in height)

Adjustable Footrest Options

XB180 Aluminium footplate (2.5mm) **No Cost**

XB100 ABS moulded foot tray **\$188.00**

Specify length

Max = 380mm (Full)
mm Std = 330mm

Footrest Plate Options

XB240 Aluminium raised edges **\$194.00**

(Cost per edge)

Front

Specify height

Specify Depth

Back

Left

(Max = 60mm)

Right

XB250 Carbon Fibre raised edge **\$1,391.00**

(1 side only)

Front

Specify height

Back

(Max = 45mm)

Castor pin width

XY330 Castor pin width

mm **MAX 320mm**

Optimise measurements

Camber Bar, Sideguards, Backrest and Anti-Tip

GSAXY510	Fixed centre of gravity	(Note: Dimension is taken from front of backrest to centre of axle)	Std
	Specify measurement		
	mm		
GSAXY515	Backrest COG adjustment		\$408.00
	(not with welded sides) (+/- 20mm)		
	Specify measurement		
	mm		
GSAXY560	Adjustable centre of gravity		NC
	Specify measurement		
	mm		

Specify Camber Angle

XY390	Camber angle	GSAXJ040	Single fixed anti tip	No Cost
	16 18 20 22		Specify distance behind rear wheel	
			mm (Std=100mm)	

Sideguards

GSAXE240	Welded aluminium sport	\$1,084.00
	(Not with backrest COG Adj - GSAXY515)	
XE030	Screw on aluminium	\$337.00
XE050	Screw on carbon fibre	\$748.00

Sideguard measurements

Sideguard length	25mm above tyre
mm	(Follow tyre)
Sideguard height	Individual above tyre
mm	(Follow tyre)
	mm
	Level with tyre
	(Follow tyre)

Backrest

XY230	Fixed backrest height	Std	XY245	Adjustable backrest height	\$255.00
	Specify height			(+/-25mm)	
	mm			Specify height	
XY360	Backrest angle		XG100	Ergonomic backrest	\$150.00
-5 -3 0 +3 +5				Vertical height	
				mm	
Individual	Deg			Rearward position	
(Max -10/+10 Deg)				mm	
XG110	Inward taper		XG120	Outward flare	\$150.00
	(Inside of backrest posts)			(Inside of backrest posts)	
	Specify width or per side			Specify width or per side	
	mm			mm	

Castor Options

XC160	Sport 61mm Castor (Black)	XC181	Bones Reds bearings (All castor bearings)	
XC170	Sport 72mm Castor (Black)	XC182	Bones Swiss Ceramic bearings (All castor bearings)	\$619.00

Rear Wheels

Rear Wheel Size

XY420	24" Wheels	XY430	25" Wheels	XY440	26" Wheels	XY450	700c Wheels
-------	------------	-------	------------	-------	------------	-------	-------------

Rear Wheel Gap

XY460	Rear Wheel Gap	XY470	Inside to inside of wheel
15	20	25	

Individual

mm

mm

Rear Wheel Options

XK030	Equalizer 36/48 Spoke (Black hub & spokes)				Std
XK080	Spinergy SLX 24 Spoke				\$1,036.00
XK100	Spinergy XSLX 24 Spoke				\$1,092.00
XK040	Supply without wheels and axle pins	Wheel type:	Tyre type:		NC

Spinergy Colour Options

Black Spokes	White Spokes	Yellow Spokes	Red Spokes	Green Spokes
Blue Spokes	Pink Spokes	Purple Spokes	Orange Spokes	
XK180	Individual spoke design (Specify colour)	Inside: Outside:		\$219.00

Spinergy Anodise Options

XK190	Black hub (Std option)	XK200	Silver hub	\$131.00
XK201	Red hub/rims			\$421.00
XK202	Purple hub/rims			\$421.00
XK203	Blue hub/rims			\$421.00
XK204	Gold hub/rims			\$421.00

Spare Wheel

GSAXK500	Spare wheel spinergy (Built up + axle pin)	\$2,019.00
XK520	Spare wheel upcharge - Titanium rims	\$634.00
XK522	Spare wheel upcharge - Remove excess tabs	\$123.00
XK524	Spare wheel upcharge - Individual Spoke Colour	\$209.00

Handrim Position/Options

XY480	Narrow set handrim	
XY485	Wide set handrim	
XL090	Remove excess mounting tabs (recommended for small wheel gaps)	\$150.00
XL080	Super narrow set handrim (only with XL010 & XL070 and includes XL090)	\$294.00

Handrim Style

XL010	Aluminium hard anodised 6 tab	Std
XL040	Titanium 6 tab (Not with 700C)	\$459.00
XL070	Aluminium hard anodised 12 tab (24" not available)	\$150.00
XL020	Tetra Grip 6 tab	\$412.00
XL130	Aluminium hard anodised 6 tab painted (Select Colour)	\$166.00
XL100	Para grip 6 tab (700C not available)	\$412.00

Tyre Options

XM070	Kenda Kaliente	Red (Std)	Blue	Std
XM150	IRC Exeracer Pro	(Red Only)		\$194.00
XM080	Tufo MS3 tubeless	(Blue Only)		\$83.00
XM095	Kenda Wide	(Grey Only)		\$80.00

Inner Tube Valve Options

XY500	Schraeder inner tube valve	No cost
XY490	Presta inner tube valve (Standard if no option selected)	Std

Finish Options

XN010 Powder coated frame No Cost

Frame	Gloss	Frame	Matt	Frame	Special
	Black gloss 103		Black matt 105		Chrome mirror 300
	White gloss 100		White matt 101		Sparkle silver 170
	Grey gloss 132		Grey matt 131		Metallic black matt 104
	Cream gloss 135		Light grey matt 133		Metallic grey matt 130
	Red gloss 108		Navy blue matt 152		Black sandpaper 201
	Yellow gloss 120				Grey sandpaper 202
	Forest green gloss 143				White sandpaper 200
	Green gloss 107				
	Dark blue gloss 151				
	Royal blue gloss 154				
	Blue gloss 174				
	Duckegg Blue gloss 141				
	Pink gloss 188				
	Orange gloss 180				
	Purple gloss 157				

XN150 Powder coated components No Cost
(Frame colour if not specified)
Specify component colour code from above list

GSAXN030	Brushed Frame (All parts black unless XN160 chosen)	\$582.00
XN160	Brushed Components	\$125.00
GSAXN050	Polished Frame (All parts black unless XN170 chosen)	\$685.00
XN170	Polished Components	\$219.00

Accent Colour

Note: Accent colour includes forks only

Black accent	Silver accent
Red accent	Blue accent
Purple accent	



Upholstery

Seat upholstery

XP040	Seat sling strap adjustable (Std)	No seat sling (Seat Plate)	Std
-------	-----------------------------------	----------------------------	-----

Backrest upholstery

XP055	Airtech sport backrest strap adjustable (Std)	XP070	Backrest Leather corners	\$125.00
-------	---	-------	--------------------------	----------

Backrest Stitching options

Silver	Black	Blue	Individual:
Red	Orange	Purple	

RGK Logo colour

White (std)	Green	Blue	Individual:
Red	Orange	Purple	

Model Logo colour

White (std)	Green	Blue	Individual:
Red	Orange	Purple	

XP170	Individual backrest embroidery (Add details on page 12)	\$247.00
-------	--	----------

Cushion

XP080	Cushion Foam (Select density and size below)	\$125.00
-------	---	----------

Soft (CMS040)		Medium (CMS050)		
50mm	75mm	100mm	Other	mm

XP340	Hard Plastazote cushion upgrade (Indicate max height)	\$118.00
	mm	

Cushion Cover Options

XP085	Nylon cushion cover (Black)	No Cost
XP380	Neoprene cushion cover (Black)	\$131.00
XP090	Air Lite cushion cover (Black)	\$101.00

Straps

XS050	100mm Pro adjustable velcro strap	\$178.00
XS055	50mm Pro adjustable velcro strap	\$139.00
XS057	50mm Pro knee strap (measurement from tube to tube over knees)	\$258.00

Strap 1	Type	Measurement		
Strap 2	Type	Measurement	Specify strap length	Distance from back post
Strap 3	Type	Measurement	mm	mm

XS020	Sport ratchet strap (Click strap) (11" x 2" pad)	Quantity:	1	2	\$144.00
XS021	Sport ratchet strap (Single)	Quantity:	1	2	\$185.00
XS025	Pelvi Loc 2DA ratchet strap (For thighs & knees)	Quantity:	1	2	\$418.00
	Select size:				
	Medium (9")	Large (11")			

XS030	Pelvi Loc 3DA ratchet strap (For hips)	Quantity:	1	2	\$432.00
	Select size:				
	XS (8")	S (12")	M (14")	L (17")	

XS035	M2 Ratchet strap (For thighs and knees)	Quantity:	1	2	\$291.00
	6" 8" 10" 12"				

XS036	M2 Ratchet strap (Single) (For Feet)	Quantity:	1	2	\$308.00
	6" 8" 10" 12"				

Strap 1	Type	Strap 2	Type	Strap 3	Type
					Specify distance from back post
					mm

	Specify distance from back post	Specify distance from back post			
	mm	mm			
XS040	Fixed mounting points (One set per click strap)	x1	x2	x3	\$165.00
XS045	Clamped mounting points (One set per click strap)	x1	x2	x3	\$99.00

Frame options

XR060	Calf strap mounting loops	\$159.00
XR065	Strap location loops	\$122.00

Accessories

XR030	Triple wheel bag (Black)			\$207.00
XP130	Frame protectors full length			\$131.00
XS100	Body point waist strap	Small	Large	\$106.00
XS070	Foot strap (select if required)			No Cost
XS060	Calf strap (select if required)			No Cost
XP400	Back bar protector			No Cost
XR090	Tool kit			No Cost



1000 Creditstone Road, Unit 2, Concord, ON L4K 4P8

Customer Service: (800) 263-3390 E-Mail: cscanada@sunmed.com www.sunrisemedical.ca
Specifications are subject to change without notice. © 9/2026 Sunrise Medical Canada Inc.
CN-100063-CAN REV (7).