



### Order Form

### EIPW32


|             |  |
|-------------|--|
| Account:    |  |
| Date:       |  |
| PO#:        |  |
| Buyer Name: |  |
| Marked For: |  |

|                                 |  |
|---------------------------------|--|
| Ship To:                        |  |
| Name:                           |  |
| Address:                        |  |
| City, Province,<br>Postal Code: |  |
| Phone:                          |  |

Sunrise Medical Canada Inc.  
 237 Romina Drive, Unit 3  
 Concord, ON L4K 4V3  
 1 (800) 263-3390

# Client Measurements

|     |                         |      |
|-----|-------------------------|------|
| a)  | Top of head to seat     | Inch |
| b)  | Top of shoulder to seat | Inch |
| c)  | Underarm to seat        | Inch |
| d)  | Elbow to seat           | Inch |
| e)  | Top of leg to seat      | Inch |
| f)  | Shoulder width          | Inch |
| g)  | Chest width             | Inch |
| i)  | Hip width               | Inch |
| j)  | Actual sitting depth    | Inch |
| m)  | Foot length             | Inch |
| (n) | Knee to foot            | Inch |



## User Weight

Notes:

## Chair Specifications

|  |   |                         |  |
|--|---|-------------------------|--|
| <b>Weights starting at:<br/>(With Batteries)</b> | <i>NPC: 342 lbs.<br/>SPC: 367 lbs.<br/>MPC: 381 lbs.</i>  | <b>Seat Type</b>        | <i>SEDEO PRO</i>                       |
| <b>Power Seating</b>                             | <i>50° of CG Tilt and 150° Power Recline w/ 4" mechanical shear reduction / 12" powered seat elevator</i> | <b>Seat Size Range</b>  | <i>Width 16"-22"<br/>Depth 14"-22"</i> |
|  |   | <b>STF Height Range</b> | <i>16.5"-19"</i>                       |

## Hybrid Drive Base Specifications

|                               |  |                         |   |
|-------------------------------|--|-------------------------|---|
| <b>Drive Wheel Position</b>   | <i>Hybrid Drive</i>  | <b>Drive Wheels</b>     | <i>14" and 13" drive wheels w/flat free inserts<br/>14" and 13" w/pneumatic inserts</i> |
| <b>User Weight Capacity</b>   | <i>300 lbs. / 350 lbs.HD</i>                                 | <b>Tires/Casters</b>    | <i>9" Front / 4" Rear (Black)</i>   |
| <b>Minimum Turning Radius</b> | <i>25.5" (excludes footplates)</i>                           | <b>Charger</b>          | <i>10A Dual Mode (Compact)</i>  |
| <b>Base Size</b>              | <i>Width 24.5", Length 36"</i>                               | <b>Battery Size</b>     | <i>Group 24 / Group 34</i>  |
| <b>Drive Control</b>          | <i>VR-2 90, R-net 90A, R-net 120A</i>                        | <b>Color Accents</b>    | <i>7 Colour Accents</i>   |
| <b>Driving Range</b>          | <i>Up to 17.5 miles (w/ Gp24 batteries and 5 mph motors)</i> | <b>Ground Clearance</b> | <i>3"</i>   |
| <b>Drive Motors</b>           | <i>4-pole 6mph, 5mph HD, 8mph</i>                            |                         |   |

# Base

## Standard Duty Hybrid Drive Base = End Item Number EIPW32 (AADL Contract Price)

|  |          |   |            |
|--|----------|---|------------|
| <input type="checkbox"/>   | 320PC248 | Quickie Q500 H NPC No Power Option for SEDEO PRO Seat     | \$8,995.00 |
| <i>Included in Base Pricing : Q500 H Base, VR-2 90amp Electronics for No Power Seat Functions, Six Wheel Independent Suspension System, 14" Drive Wheels with Dual Density Inserts, 9" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5".</i> |          |   |            |
| <input type="checkbox"/>   | 320PC249 | Quickie Q500 H SPC Single Power Option for Sedeo PRO Seat | \$8,995.00 |
| <i>Included in Base Pricing : Q500 H Base, VR-2 90amp Electronics for No Power Seat Functions, Six Wheel Independent Suspension System, 14" Drive Wheels with Dual Density Inserts, 9" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5".</i> |          |   |            |
| <input type="checkbox"/>   | 320PC250 | Quickie Q500 H MPC Tilt/Recline For SEDEO PRO Seat        | \$8,995.00 |
| <i>Included in Base Pricing : Q500 H Base, VR-2 90amp Electronics for No Power Seat Functions, Six Wheel Independent Suspension System, 14" Drive Wheels with Dual Density Inserts, 9" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5".</i> |          |   |            |

## Heavy Duty Hybrid Drive Base = End Item Number EIPW32 (AADL Contract Price)

|  |          |  |             |
|--|----------|--|-------------|
| <input type="checkbox"/>   | 320PC251 | Quickie Q500 H NPC HD No Power Option for SEDEO PRO Seat     | \$10,000.00 |
| <i>Included in Base Pricing : Q500 H Base, VR-2 90amp Electronics for No Power Seat Functions, Six Wheel Independent Suspension System, 14" Drive Wheels with Dual Density Inserts, 9" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 5mph, Weight Capacity = 350lbs., STFH = 16.5".</i> |          |  |             |
| <input type="checkbox"/>   | 320PC252 | Quickie Q500 H SPC HD Single Power Option for SEDEO PRO Seat | \$10,000.00 |
| <i>Included in Base Pricing : Q500 H Base, VR-2 90amp Electronics for No Power Seat Functions, Six Wheel Independent Suspension System, 14" Drive Wheels with Dual Density Inserts, 9" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 5mph, Weight Capacity = 350lbs., STFH = 16.5".</i> |          |  |             |
| <input type="checkbox"/>   | 320PC253 | Quickie Q500 H MPC HD Tilt/Recline For SEDEO PRO Seat        | \$10,000.00 |
| <i>Included in Base Pricing : Q500 H Base, VR-2 90amp Electronics for No Power Seat Functions, Six Wheel Independent Suspension System, 14" Drive Wheels with Dual Density Inserts, 9" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 5mph, Weight Capacity = 350lbs., STFH = 16.5".</i> |          |  |             |

## High Performance 8mph Standard Duty Hybrid Drive Base = End Item Number EIPW32 (Client Cost)

|  |          |   |             |
|--|----------|---|-------------|
| <input type="checkbox"/>   | 320PC254 | Quickie Q500 H NPD No Power Option for SEDEO PRO Seat     | \$10,657.00 |
| <i>Included in Base Pricing : Q500 H Base, R-net 120amp Electronics for No Power Seat Functions, Six Wheel Independent Suspension System, 14" Drive Wheels with Dual Density Inserts, 9" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 8mph, Weight Capacity = 300lbs., STFH = 16.5".</i> |          |   |             |
| <input type="checkbox"/>   | 320PC255 | Quickie Q500 H SPD Single Power Option for SEDEO PRO Seat | \$10,657.00 |
| <i>Included in Base Pricing : Q500 H Base, R-net 120amp Electronics for No Power Seat Functions, Six Wheel Independent Suspension System, 14" Drive Wheels with Dual Density Inserts, 9" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 8mph, Weight Capacity = 300lbs., STFH = 16.5".</i> |          |   |             |
| <input type="checkbox"/>   | 320PC256 | Quickie Q500 H MPD Tilt/Recline For SEDEO PRO Seat        | \$10,657.00 |
| <i>Included in Base Pricing : Q500 H Base, R-net 120amp Electronics for No Power Seat Functions, Six Wheel Independent Suspension System, 14" Drive Wheels with Dual Density Inserts, 9" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 8mph, Weight Capacity = 300lbs., STFH = 16.5".</i> |          |   |             |

# Base

## User Weight

|                          |         |                |    |
|--------------------------|---------|----------------|----|
| <input type="checkbox"/> | 320WG15 | 301 - 350 lbs. | NC |
|--------------------------|---------|----------------|----|

## Motor Package

|                          |          |   |            |
|--------------------------|----------|---|------------|
| <input type="checkbox"/> | 320PF98  | 4-pole 8 kph / 5 mph <i>(Included with HD base)</i>                                       | NC         |
| <input type="checkbox"/> | 320PF99  | 4-pole 10 kph / 6 mph <i>(Included with Standard base)</i>                                | NC         |
| <input type="checkbox"/> | 320PF101 | 4-pole 12.5 kph / 8 mph <i>(Upgrade for Standard base. NA with HD base)</i> (Client Cost) | \$1,625.00 |

## Tracking

|                          |         |   |            |
|--------------------------|---------|---|------------|
| <input type="checkbox"/> | 320MI01 | SureTrac Module & Motor Encoders <i>(Improves Drive Tracking)</i> | \$1,565.00 |
|--------------------------|---------|---|------------|

## Drive Wheels

|                          |         |   |          |
|--------------------------|---------|---|----------|
| <input type="checkbox"/> | 320DW58 | Drive Wheel, 14", Dual Density, Knobby, Black | NC       |
| <input type="checkbox"/> | 320DW67 | Drive Wheel, 14", Pneumatic, Knobby, Black    | \$384.00 |
| <input type="checkbox"/> | 320DW69 | Drive Wheel, 13", Dual Density, Knobby, Black | NC       |
| <input type="checkbox"/> | 320DW76 | Drive Wheel, 13", Pneumatic, Knobby, Black    | \$384.00 |

## Caster Wheels

|                          |         |                                    |    |
|--------------------------|---------|------------------------------------|----|
| <input type="checkbox"/> | 320CW57 | Caster Wheels 9", Solid, Black     | NC |
| <input type="checkbox"/> | 320CW56 | Caster Wheels 9", Pneumatic, Black | NC |

## Wheel Locks

|                          |         |             |          |
|--------------------------|---------|-------------|----------|
| <input type="checkbox"/> | 320WL42 | Wheel Locks | \$130.00 |
|--------------------------|---------|-------------|----------|

## Color Accents

|                          |         |             |    |
|--------------------------|---------|-------------|----|
| <input type="checkbox"/> | 320EZ20 | Apple Green | NC |
| <input type="checkbox"/> | 320EZ22 | Black       | NC |
| <input type="checkbox"/> | 320EZ21 | Blue        | NC |
| <input type="checkbox"/> | 320EZ23 | Orange      | NC |
| <input type="checkbox"/> | 320EZ24 | Red         | NC |
| <input type="checkbox"/> | 320EZ25 | Silver      | NC |
| <input type="checkbox"/> | 320EZ19 | White       | NC |

## ISO Transit Tie Downs

|                          |         |                     |    |
|--------------------------|---------|---------------------|----|
| <input type="checkbox"/> | 320TO50 | ISO Transit Package | NC |
|--------------------------|---------|---------------------|----|

## LED Lights (Client Cost)

|                          |         |                           |            |
|--------------------------|---------|---------------------------|------------|
| <input type="checkbox"/> | 320LI04 | LED Lights and Indicators | \$1,131.00 |
| <input type="checkbox"/> | 320LI02 | Omit LED Lights           | NC         |

## Seat Heights

Note: Available Seat-to-Floor heights are dependant on 2 options: 1. Battery size and 2. Powered Seating.

Use the matrix below to determine available ranges for your desired options. Measured from front of seat pan to floor with 0° frame angle.

|                          | Seat Height |          |
|--------------------------|-------------|----------|
| <input type="checkbox"/> | 320HT103    | 16.5" NC |
| <input type="checkbox"/> | 320HT55     | 17" NC   |
| <input type="checkbox"/> | 320HT104    | 17.5" NC |
| <input type="checkbox"/> | 320HT56     | 18" NC   |
| <input type="checkbox"/> | 320HT105    | 18.5" NC |
| <input type="checkbox"/> | 320HT57     | 19" NC   |

## Batteries & Chargers (AADL Contract Pricing)

### Group 24 Battery Options *(Recommended for max daily range)*

|                          |         |  |          |
|--------------------------|---------|--|----------|
| <input type="checkbox"/> | 320AG25 | GP24 Flat Top Pair w/Strap MK Gel        | \$844.00 |
| <input type="checkbox"/> | 320AG26 | GP24 Flat Top Pair w/Strap Brightway AGM | \$881.00 |

### Group 34 Battery Options *(Required for lower seat-to-floor heights)*

|                          |         |   |          |
|--------------------------|---------|---|----------|
| <input type="checkbox"/> | 320AG27 | GP34 Flat Top Pair w/Straps MK Gel        | \$918.00 |
| <input type="checkbox"/> | 320AG29 | GP34 Flat Top Pair w/Straps Brightway AGM | \$796.00 |

### Omit Batteries

|                          |         |  |    |
|--------------------------|---------|--|----|
| <input type="checkbox"/> | 320AG24 | Battery Harnesses Only <i>(Omit Batteries)</i> | NC |
|--------------------------|---------|--|----|

### Chargers

|                          |         |                                |          |
|--------------------------|---------|--------------------------------|----------|
| <input type="checkbox"/> | 320BJ30 | 10A Dual Mode <i>(Compact)</i> | \$150.00 |
| <input type="checkbox"/> | 320BJ19 | 8A Dual Mode <i>(Compact)</i>  | NC       |

# Powered Seating

## Power Seating Options

### Seat Elevator (Client Cost)

|                          |         |  |            |
|--------------------------|---------|--|------------|
| <input type="checkbox"/> | 320PI07 | Seat Elevator, Powered, 300 mm / 12" (N/A w/ HD or 8 mph motors)   | \$3,209.00 |
| <input type="checkbox"/> | 320PI08 | C-Me Seat Elevator, Powered, 300 mm / 12" (Not available with HD model or 8 mph motors)<br>C-Me speed up to 3.0 mph. | \$3,209.00 |

### Seat Tilt (AADL W993 for Power Tilt)

|                          |         |                       |            |
|--------------------------|---------|-----------------------|------------|
| <input type="checkbox"/> | 320PA35 | Posterior 50° Tilt    | \$7,995.00 |
| <input type="checkbox"/> | 320PA36 | Posterior 50° HD Tilt | \$8,235.00 |

### Seat Recline (AADL W993 for Power Recline)

|                          |         |   |            |
|--------------------------|---------|---|------------|
| <input type="checkbox"/> | 320PA37 | 150° Powered Recline w/ 4" Shear Reduction    | \$7,385.00 |
| <input type="checkbox"/> | 320PA39 | 150° HD Powered Recline w/ 4" Shear Reduction | \$8,235.00 |

### Seat Tilt/Recline (AADL W993 for Power Tilt and Recline)

|                          |         |  |             |
|--------------------------|---------|--|-------------|
| <input type="checkbox"/> | 320PA41 | Posterior 50° Tilt and 150° Recline Powered, w/ 4" Shear Reduction       | \$12,252.00 |
| <input type="checkbox"/> | 320PA43 | Posterior 50° HD Tilt and 150° HD Recline Powered, w/ 4" Shear Reduction | \$13,282.00 |

## Power Recline Backrest Preset Home Angle

|                          |        |                        |    |
|--------------------------|--------|------------------------|----|
| <input type="checkbox"/> | 320QK1 | 90°                    | NC |
| <input type="checkbox"/> | 320QK2 | 95° (Default Position) | NC |
| <input type="checkbox"/> | 320QK3 | 100°                   | NC |

## Pre Set Seat Frame Angle (No Power Seating Options, Increases Seat to Floor Height)

|                          |         |                      |    |
|--------------------------|---------|----------------------|----|
| <input type="checkbox"/> | 320FA10 | 0° (Adds 0" to STFH) | NC |
| <input type="checkbox"/> | 320FA3  | 3° (Adds 1" to STFH) | NC |
| <input type="checkbox"/> | 320FA4  | 6° (Adds 2" to STFH) | NC |
| <input type="checkbox"/> | 320FA5  | 9° (Adds 3" to STFH) | NC |

## Backrest Frame Type

|                          |         |  |          |
|--------------------------|---------|--|----------|
| <input type="checkbox"/> | 320BF08 | Captain Seat Frame (Includes Back Cushion)   | NC       |
| <input type="checkbox"/> | 320BF10 | Tool Angle Adjustable Back, 90°-120°   | NC       |
| <input type="checkbox"/> | 320BF11 | User Adjustable Manual Semi-Recline Back, 90°-130°, Right (Recline lever mounted on opposite side of joystick) | \$760.00 |
| <input type="checkbox"/> | 320BF12 | User Adjustable Manual Semi-Recline Back, 90°-130°, Left (Recline lever mounted on opposite side of joystick)  | \$760.00 |

## Backrest Frame Option

|                          |         |                          |          |
|--------------------------|---------|--------------------------|----------|
| <input type="checkbox"/> | 320BF13 | Flip Down, Quick Release | \$165.00 |
|--------------------------|---------|--------------------------|----------|

## Back Cushion Type

|                          |          |   |          |
|--------------------------|----------|---|----------|
| <input type="checkbox"/> | 320BT134 | SEDEO Back Cushion (Choose Seating in Section A)  | NC       |
| <input type="checkbox"/> | 320BT136 | JAY J3 Interface Only (Allows Custom JAY Backs)   | NC       |
| <input type="checkbox"/> | 320BT137 | SEDEO PRO Aftermarket Back Adapter (Height Adjustable Round Tubes for Traditional Mounting) | \$585.00 |
| <input type="checkbox"/> | 320BT133 | OMIT (For Custom Mounting. No Headrest Mount)   | NC       |

## A - SEDEO Backrest Dimensions

### Back Width

|                          |         |             |    |
|--------------------------|---------|-------------|----|
| <input type="checkbox"/> | 320BN28 | 400mm / 16" | NC |
| <input type="checkbox"/> | 320BN29 | 460mm / 18" | NC |
| <input type="checkbox"/> | 320BN30 | 500mm / 20" | NC |
| <input type="checkbox"/> | 320BN31 | 560mm / 22" | NC |

### Back Shell Height

|                          |         | Tool Angle Adjustable | 130° Recline | 150° Recline | Installed<br>Back<br>Height |         |    |
|--------------------------|---------|-----------------------|--------------|--------------|-----------------------------|---------|----|
| <input type="checkbox"/> | 320B255 | 460mm / 18"           | 22"-23"      | 22"-24"      |                             | 23"-25" | NC |
| <input type="checkbox"/> | 320B256 | 560mm / 22"           | 24"-27"      | 25"-27"      |                             | 26"-29" | NC |
| <input type="checkbox"/> | 320B257 | 660mm / 26"           | 28"-30"      | 28"-30"      |                             | 30"     | NC |

### Installed Back Height - (From Pan to Top of Back)

|                          |         |     |    |
|--------------------------|---------|-----|----|
| <input type="checkbox"/> | 320JR32 | 22" | NC |
| <input type="checkbox"/> | 320JR33 | 23" | NC |
| <input type="checkbox"/> | 320JR34 | 24" | NC |
| <input type="checkbox"/> | 320JR35 | 25" | NC |
| <input type="checkbox"/> | 320JR36 | 26" | NC |
| <input type="checkbox"/> | 320JR37 | 27" | NC |
| <input type="checkbox"/> | 320JR38 | 28" | NC |
| <input type="checkbox"/> | 320JR39 | 29" | NC |
| <input type="checkbox"/> | 320JR40 | 30" | NC |

# Seating

## Headrests (Requires Pad, Upholstery and Mount) (Seating Benefits)

### Headrest Pad

|                          |         |   |          |
|--------------------------|---------|---|----------|
| <input type="checkbox"/> | 320HE36 | SEDEO Pad Small With Cover (LxH) 240 x 150 mm / 10" x 6"<br><i>(Only Available with SEDEO Multi-Position Headrest Mount and SEDEO Upholstery)</i> | \$210.00 |
| <input type="checkbox"/> | 320HE37 | SEDEO Pad Flexible, (LxH) 400 x 110 mm / 16" x 4.5"<br><i>(Only Available with SEDEO Multi-Position Headrest Mount and SEDEO Upholstery)</i>      | \$340.00 |
| <input type="checkbox"/> | 320HE38 | Whitmyer Plush Headrest Pad 203 mm / 8"   | \$328.00 |
| <input type="checkbox"/> | 320HE39 | Whitmyer Plush Headrest Pad 254 mm / 10"  | \$328.00 |
| <input type="checkbox"/> | 320HE44 | Whitmyer Plush Headrest Pad 356 mm / 14"  | \$328.00 |
| <input type="checkbox"/> | 320HE45 | Whitmyer Adjust-A-Plush Headrest Pad Small (Narrow) (Adjustable from 6" to 14" Wide x 3" Tall)  | \$387.00 |
| <input type="checkbox"/> | 320HE46 | Whitmyer Adjust-A-Plush Headrest Pad Medium (Adjustable from 7" to 16" Wide x 4" Tall)  | \$387.00 |
| <input type="checkbox"/> | 320HE47 | Whitmyer Adjust-A-Plush Headrest Pad Large (Standard) (Adjustable from 8" to 18" Wide x 5" Tall)  | \$387.00 |
| <input type="checkbox"/> | 320HE48 | OMIT Headrest Pad   | NC       |

### Headrest Pad Upholstery

|                          |         |                            |    |
|--------------------------|---------|----------------------------|----|
| <input type="checkbox"/> | 320NU10 | SEDEO Upholstery           | NC |
| <input type="checkbox"/> | 320NU04 | Syncra Premium Leatherette | NC |
| <input type="checkbox"/> | 320NU05 | Syncra Stretch             | NC |
| <input type="checkbox"/> | 320NU06 | Syncra 3DX Microclimatic   | NC |
| <input type="checkbox"/> | 320NU11 | Whitmyer Lycra             | NC |

### Headrest Mount

|                          |         |   |          |
|--------------------------|---------|---|----------|
| <input type="checkbox"/> | 320MH32 | SEDEO Multi-Position Headrest Mount with Quick Release <i>(Only available with SEDEO Pad)</i> | \$515.00 |
| <input type="checkbox"/> | 320MH33 | Whitmyer Linx Headrest Mount with Quick Release   | \$215.00 |
| <input type="checkbox"/> | 320MH34 | Whitmyer Linx-2 Headrest Mount with Quick Release   | \$260.00 |
| <input type="checkbox"/> | 320MH35 | Whitmyer AXYS Headrest Mount with Quick Release   | \$215.00 |
| <input type="checkbox"/> | 320MH36 | Whitmyer AXYS Offset Headrest Mount with Quick Release  | \$215.00 |
| <input type="checkbox"/> | 320MH37 | Whitmyer Cobra FB Headrest Mount with Quick Release   | \$215.00 |
| <input type="checkbox"/> | 320MH17 | OMIT Headrest Mount   | NC       |

### Push Handles

|                          |         |   |         |
|--------------------------|---------|---|---------|
| <input type="checkbox"/> | 320PH27 | Height Adjustable Trailing Push Handles | \$57.00 |
|--------------------------|---------|---|---------|

### Seat Frame Rails

|                          |         |  |    |
|--------------------------|---------|--|----|
| <input type="checkbox"/> | 320SN14 | Seat Rail Short <i>(For Effective Seat Depth 14" - 18") (No charge with power seating)</i> | NC |
| <input type="checkbox"/> | 320SN15 | Seat Rail Long <i>(For Effective Seat Depth 17" - 22") (No charge with power seating)</i>  | NC |

### Seat Cushion Type

|                          |         |       |    |
|--------------------------|---------|-------|----|
| <input type="checkbox"/> | 320ST79 | SEDEO | NC |
| <input type="checkbox"/> | 320ST75 | Omit  | NC |

### Seat Size

#### Seat Pan Width (Odd sizes not available on SEDEO cushions)

|                          |        |     |    |
|--------------------------|--------|-----|----|
| <input type="checkbox"/> | 320R53 | 16" | NC |
| <input type="checkbox"/> | 320R54 | 17" | NC |
| <input type="checkbox"/> | 320R55 | 18" | NC |
| <input type="checkbox"/> | 320R56 | 19" | NC |
| <input type="checkbox"/> | 320R57 | 20" | NC |
| <input type="checkbox"/> | 320R58 | 21" | NC |
| <input type="checkbox"/> | 320R59 | 22" | NC |

#### Seat Pan Depth (Front of Seat Pan to front of back cushion)

|                          |         |     |    |
|--------------------------|---------|-----|----|
| <input type="checkbox"/> | 320RF52 | 14" | NC |
| <input type="checkbox"/> | 320RF53 | 15" | NC |
| <input type="checkbox"/> | 320RF51 | 16" | NC |
| <input type="checkbox"/> | 320RF45 | 17" | NC |
| <input type="checkbox"/> | 320RF46 | 18" | NC |
| <input type="checkbox"/> | 320RF47 | 19" | NC |
| <input type="checkbox"/> | 320RF48 | 20" | NC |
| <input type="checkbox"/> | 320RF49 | 21" | NC |
| <input type="checkbox"/> | 320RF50 | 22" | NC |



# Positioning Supports & Accessories

## Knee and Thigh Swing Away Support Pads (Seating Benefits)

### Lateral Knee Support Pads and Hardware (Must Pick Pad and Mount)

|                          |         |  |          |
|--------------------------|---------|--|----------|
| <input type="checkbox"/> | 320XY05 | Lateral Knee Support Pad, Right (4" round) | \$228.00 |
| <input type="checkbox"/> | 320XY06 | Lateral Knee Support Pad, Left (4" round)  | \$228.00 |
| <input type="checkbox"/> | 320XX08 | Lateral Swing Away Bracket Knee - Right    | \$300.00 |
| <input type="checkbox"/> | 320XX09 | Lateral Swing Away Bracket Knee - Left     | \$300.00 |

### Lateral Thigh Support Pads and Hardware (Must pick pad and mount)

|                          |         |  |          |
|--------------------------|---------|--|----------|
| <input type="checkbox"/> | 320KF22 | Lateral Thigh Support Pad Short, Right (4.8" H x 4.8" L / 120mm x 120mm) | \$134.00 |
| <input type="checkbox"/> | 320KF23 | Lateral Thigh Support Pad Short, Left (4.8" H x 4.8" L / 120mm x 120mm)  | \$134.00 |
| <input type="checkbox"/> | 320KF24 | Lateral Thigh Support Pad Long, Right (4.8" H x 7" L / 120 mm x 180mm)   | \$134.00 |
| <input type="checkbox"/> | 320KF25 | Lateral Thigh Support Pad Long, Left (4.8" H x 7" L / 120mm x 180mm)     | \$134.00 |
| <input type="checkbox"/> | 320KF26 | Lateral Thigh Removable Mount - Right                                    | \$300.00 |
| <input type="checkbox"/> | 320KF27 | Lateral Thigh Removable Mount - Left                                     | \$300.00 |

## JAY Lateral Thigh Supports (Must Pick Pad, Upholstery and Mount) (Seating Benefits)

### JAY Lateral Thigh Support Pads

|                          |         |  |          |
|--------------------------|---------|--|----------|
| <input type="checkbox"/> | 320KF16 | Lateral Thigh Support Pad Standard, Right (8.5" x 3.25") | \$134.00 |
| <input type="checkbox"/> | 320KF17 | Lateral Thigh Support Pad Standard, Left (8.5" x 3.25")  | \$134.00 |
| <input type="checkbox"/> | 320KF18 | Lateral Thigh Support Pad Large, Right (12" x 3.25")     | \$134.00 |
| <input type="checkbox"/> | 320KF19 | Lateral Thigh Support Pad Large, Left (12" x 3.25")      | \$134.00 |

### JAY Lateral Thigh Support Pad Upholstery

|                          |         |                            |    |
|--------------------------|---------|----------------------------|----|
| <input type="checkbox"/> | 320NU07 | Syncra Premium Leatherette | NC |
| <input type="checkbox"/> | 320NU08 | Syncra Stretch             | NC |
| <input type="checkbox"/> | 320NU09 | Syncra 3DX Microclimatic   | NC |

### JAY Lateral Thigh Removable Mount

|                          |         |                                       |          |
|--------------------------|---------|---------------------------------------|----------|
| <input type="checkbox"/> | 320KF20 | Lateral Thigh Removable Mount - Right | \$300.00 |
| <input type="checkbox"/> | 320KF21 | Lateral Thigh Removable Mount - Left  | \$300.00 |

## Armrest Structure

|                          |         |   |          |
|--------------------------|---------|---|----------|
| <input type="checkbox"/> | 320EG35 | Reclining, Flip Back, Long (Height Adjustable 9" to 15")                                  | NC       |
| <input type="checkbox"/> | 320EG36 | Reclining, Flip Back with 2-point post locking arm (Height Adjustable 11" to 15" x .625") | \$325.00 |
| <input type="checkbox"/> | 320EG32 | Single Post, Flip Back (N/A w/ thigh support)(Height Adjustable 10" to 15")               | NC       |
| <input type="checkbox"/> | 320EG33 | Single Post, Removable (N/A w/ thigh support)(Height Adjustable 10" to 15")               | NC       |

## Arm Pads

### SEDEO Armrest Pad

|                          |         |  |    |
|--------------------------|---------|--|----|
| <input type="checkbox"/> | 320AB43 | 300 mm / 12", Integral Foam, Right   | NC |
| <input type="checkbox"/> | 320AE43 | 300 mm / 12", Integral Foam, Left  | NC |
| <input type="checkbox"/> | 320AB44 | 400 mm / 16", Integral Foam, Right   | NC |
| <input type="checkbox"/> | 320AE44 | 400 mm / 16", Integral Foam, Left  | NC |
| <input type="checkbox"/> | 320AB45 | Wide, 450 mm / 18", Integral Foam, Right (Increases Armpad width by 1.5". May prevent flip-back) | NC |
| <input type="checkbox"/> | 320AE45 | Wide, 450 mm / 18", Integral Foam, Left (Increases Armpad width by 1.5". May prevent flip-back)  | NC |

### Comfort Arm Pad (Foam) (Seating Benefits)

|                          |         |               |          |
|--------------------------|---------|---------------|----------|
| <input type="checkbox"/> | 320AB46 | 2" Wide Right | \$104.00 |
| <input type="checkbox"/> | 320AE46 | 2" Wide Left  | \$104.00 |
| <input type="checkbox"/> | 320AB47 | 4" Wide Right | \$104.00 |
| <input type="checkbox"/> | 320AE47 | 4" Wide Left  | \$104.00 |

### Comfort Arm Pad (Gel) (Seating Benefits)

|                          |         |               |          |
|--------------------------|---------|---------------|----------|
| <input type="checkbox"/> | 320AB48 | 2" Wide Right | \$197.00 |
| <input type="checkbox"/> | 320AE48 | 2" Wide Left  | \$197.00 |
| <input type="checkbox"/> | 320AB49 | 4" Wide Right | \$197.00 |
| <input type="checkbox"/> | 320AE49 | 4" Wide Left  | \$197.00 |

### AES Trough/Elbow Block Arm pads (Gel) (Seating Benefits)

|                          |         |                  |          |
|--------------------------|---------|------------------|----------|
| <input type="checkbox"/> | 320AB50 | 3.5" x 13" Right | \$274.00 |
| <input type="checkbox"/> | 320AE50 | 3.5" x 13" Left  | \$274.00 |
| <input type="checkbox"/> | 320AB51 | 4" x 14.5" Right | \$274.00 |
| <input type="checkbox"/> | 320AE51 | 4" x 14.5" Left  | \$274.00 |

# Positioning Supports & Accessories

## Side Guards (Pair) (Client Cost)

|                          |         |  |          |
|--------------------------|---------|--|----------|
| <input type="checkbox"/> | 320SG43 | Side Guards (9"L x 5.5"H) (N/A with Reclining Armrests)  | \$228.00 |
| <input type="checkbox"/> | 320SG44 | Side Guards Integral Foam (11"L x 5.5"H) (N/A with Reclining Armrests)   | \$328.00 |
| <input type="checkbox"/> | 320SG45 | Side Guards Wedge Shaped, Padded (18"L x 5"H, 3" Thick in Front 1.5" in Rear) (Only w/ Single Post Removable Armrests) | \$328.00 |
| <input type="checkbox"/> | 320SG46 | Side Guards Mounted On Seat Rail Bracket (Only with Reclining Armrests) (N/A w/ Thigh Supports)                        | \$328.00 |

## SEDEO PRO Swing-Away Legrest (Pair) (320H179 is Client Cost)

|                          |         |   |            |
|--------------------------|---------|---|------------|
| <input type="checkbox"/> | 320H177 | Legrests, Swing Away, 70°, SEDEO PRO (14" to 21" Leg Length)  | NC         |
| <input type="checkbox"/> | 320H178 | Legrests, Manual Elevating, Swing Away, SEDEO PRO (Includes Padded Calf Pads) (14".5 to 21.5" Leg Length) | \$558.00   |
| <input type="checkbox"/> | 320H179 | Legrests, Power Elevating, Swing Away, SEDEO PRO (Includes Padded Calf Pads) (15.5" to 22.5" Leg Length)  | \$3,169.00 |
| <input type="checkbox"/> | 320H180 | Hanger Receivers Only, SEDEO PRO  | NC         |
| <input type="checkbox"/> | 320H176 | Omit Legrests   | NC         |

## Swing Away Legrest Footplates

|                          |         |                                     |          |
|--------------------------|---------|-------------------------------------|----------|
| <input type="checkbox"/> | 320F369 | Angle Adjustable Flip-Up Footplates | \$262.00 |
|--------------------------|---------|-------------------------------------|----------|

## Swing Away Legrest Calf Supports

|                          |         |   |          |
|--------------------------|---------|---|----------|
| <input type="checkbox"/> | 320D141 | Calf Support, Depth Adjustable, Integral Foam (7" Wide x 10" Tall, N/A 16-17" W)                | \$384.00 |
| <input type="checkbox"/> | 320D142 | Calf Support, Depth Adjustable, Padded (6.5" to 8" Wide, Depending On Chair Width, x 9.5" Tall) | \$472.00 |

## Swing Away Legrest Options

|                          |         |  |          |
|--------------------------|---------|--|----------|
| <input type="checkbox"/> | 320F364 | Footplates, Coupled                              | \$262.00 |
| <input type="checkbox"/> | 320D144 | Heel Loops                                       | NC       |
| <input type="checkbox"/> | 320D140 | Calf Strap (Only For 70° Swing-Away Legrests)    | \$130.00 |
| <input type="checkbox"/> | 320KG4  | Knee Pads, Integral Foam (2.5" Wide x 4.5" Tall) | \$260.00 |

## SEDEO Lower Extremity Swing Away Support

|                          |         |   |          |
|--------------------------|---------|---|----------|
| <input type="checkbox"/> | 320VD01 | Lower Extremity Support, Padded, Swing Away, SEDEO, Right | \$472.00 |
| <input type="checkbox"/> | 320VD02 | Lower Extremity Support, Padded, Swing Away, SEDEO, Left  | \$472.00 |

## Pelvic Positioning Belts

|                          |         |  |         |
|--------------------------|---------|--|---------|
| <input type="checkbox"/> | 320A837 | Positioning Belt, Auto Style Buckle, 50 mm 1600 mm / 2"        | NC      |
| <input type="checkbox"/> | 320A839 | Positioning Belt, Aircraft Buckle, 50 mm x 1300 mm / 2"        | NC      |
| <input type="checkbox"/> | 320A840 | Positioning Belt, Aircraft Buckle, Padded 50 mm x 1300 mm / 2" | \$95.00 |

## Accessories

### Ventilator Option

|                          |         |   |          |
|--------------------------|---------|---|----------|
| <input type="checkbox"/> | 320VT12 | Vent Tray Fixed (Pulmonetic LTV Models) | \$562.00 |
|--------------------------|---------|---|----------|

### O2 Tank Holder

|                          |         |   |          |
|--------------------------|---------|---|----------|
| <input type="checkbox"/> | 320A718 | O2 Tank Holder Mounted to Backrest (D and E Styles) | \$361.00 |
|--------------------------|---------|---|----------|

### Wheelchair Tray Table

|                          |         |                  |          |
|--------------------------|---------|------------------|----------|
| <input type="checkbox"/> | 3207A09 | Tray, Swing Away | \$351.00 |
|--------------------------|---------|------------------|----------|

### Backpack (Client Cost)

|                          |       |               |          |
|--------------------------|-------|---------------|----------|
| <input type="checkbox"/> | 320A1 | Black Utility | \$103.00 |
|--------------------------|-------|---------------|----------|

### Medial Thigh Support (Seating Benefits)

|                          |         |                            |          |
|--------------------------|---------|----------------------------|----------|
| <input type="checkbox"/> | 3208A80 | Medial Thigh Support       | \$224.00 |
| <input type="checkbox"/> | 3208A81 | Pommel Swing Away Hardware | \$330.00 |

### Misc. Accessories (Client Cost)

|                          |         |                                      |          |
|--------------------------|---------|--------------------------------------|----------|
| <input type="checkbox"/> | 3208A59 | Transfer Bar (Pair)                  | \$455.00 |
| <input type="checkbox"/> | 3208A52 | Ram Mount Phone                      | \$71.00  |
| <input type="checkbox"/> | 3208A53 | Ram Mount Tablet                     | \$284.00 |
| <input type="checkbox"/> | 3208A54 | Ram Mount Cup Holder                 | \$71.00  |
| <input type="checkbox"/> | 3208A58 | Rear View Mirror                     | \$71.00  |
| <input type="checkbox"/> | 320A611 | Crutch Holder                        | \$125.00 |
| <input type="checkbox"/> | 3208A73 | Bag For Armrest                      | \$71.00  |
| <input type="checkbox"/> | 3208A75 | Reflectors (Backrest shroud mounted) | \$71.00  |



# Joystick / Electronics

## Joystick (Remote)

|                          |         |   |            |
|--------------------------|---------|---|------------|
| <input type="checkbox"/> | 320J127 | VR-2 Joystick   | NC         |
| <input type="checkbox"/> | 320J121 | Std LED Joystick for R-net Controls   | NC         |
| <input type="checkbox"/> | 320J116 | Color LCD (CJSM1)   | NC         |
| <input type="checkbox"/> | 320J117 | Color LCD, Heavy Duty (Fixed Joystick Mount Only)   | \$1,495.00 |
| <input type="checkbox"/> | 320J118 | High Resolution Color LCD, (CJSM2+BT) (Client Cost)<br>3.5" Color Display, Paddle Switches, Programmable Shortcut Keys, ICONs, Stereo Jack, IR and dual (IOS/Android) Bluetooth control | \$1,167.00 |
| <input type="checkbox"/> | 320J119 | Omit Remote Control (Required with Arm Mounted Specialty Controls)  | -\$25.00   |

## Joystick Mounting

|                          |         |  |          |
|--------------------------|---------|--|----------|
| <input type="checkbox"/> | 320JM34 | Fixed on Armrest                               | NC       |
| <input type="checkbox"/> | 320JM35 | Parallel Swing Away                            | \$412.00 |
| <input type="checkbox"/> | 320JM36 | Link-It Mid Line Mount on Armrest (Swing Away) | \$412.00 |

## Joystick Mount (Location)

|                          |         |                    |    |
|--------------------------|---------|--------------------|----|
| <input type="checkbox"/> | 320HM16 | Right Hand armrest | NC |
| <input type="checkbox"/> | 320HM17 | Left Hand armrest  | NC |

## Display Language Selection (R-net)

|                          |       |         |    |
|--------------------------|-------|---------|----|
| <input type="checkbox"/> | 320L1 | English | NC |
| <input type="checkbox"/> | 320L3 | French  | NC |
| <input type="checkbox"/> | 320L7 | Spanish | NC |

## Joystick Handle (Optional, N/A with HD Joystick)

|                          |         |   |          |
|--------------------------|---------|---|----------|
| <input type="checkbox"/> | 320JH13 | Ball Handle                                 | \$101.00 |
| <input type="checkbox"/> | 320JH14 | Foam Ball Handle                            | \$101.00 |
| <input type="checkbox"/> | 320JH15 | T-Handle                                    | \$101.00 |
| <input type="checkbox"/> | 320JH16 | Goalpost                                    | \$101.00 |
| <input type="checkbox"/> | 320JH24 | Bodypoint - Rubber Dome Joystick Handle     | \$200.00 |
| <input type="checkbox"/> | 320JH25 | Bodypoint - U-Shaped Joystick Handle (4 in) | \$200.00 |
| <input type="checkbox"/> | 320JH26 | Bodypoint - Grooved Mushroom Joystick       | \$200.00 |

## Electronic Programmers

|                          |         |                                 |          |
|--------------------------|---------|---------------------------------|----------|
| <input type="checkbox"/> | 320PR10 | R-net PC Programmer with Dongle | \$652.00 |
| <input type="checkbox"/> | 320PR11 | VR-2 Hand Held Programmer       | \$449.00 |

## Power Modules

### Non-Expandable

|                          |         |  |    |
|--------------------------|---------|--|----|
| <input type="checkbox"/> | 320ER12 | Power Module, VR-2, 90A, Non-Expandable (Only available with 6mph motors)    | NC |
| <input type="checkbox"/> | 320ER18 | Power Module, R-net, 90A, Non-Expandable (Only available with 5mph and 8mph) | NC |
| <input type="checkbox"/> | 320ER11 | Power Module, R-net, 120A, Non-Expandable (Only available with Group 4)      | NC |

### Expandable

|                          |         |   |    |
|--------------------------|---------|---|----|
| <input type="checkbox"/> | 320ER13 | Power Module, R-net, 90A, Expandable (Required for Specialty Controls or 2+ Actuators)                                    | NC |
|                          | 320EX17 | Expandable 90A Controller - R-net   | NC |
|                          | 320EX18 | Harness for 90A Expandable Controller   | NC |
| <input type="checkbox"/> | 320ER10 | Power Module, R-net, 120A, Expandable (Only available with Group 4)<br>(Required for Speciality controls or 2+ actuators) | NC |
|                          | 320EX09 | Expandable 120A Controller - R-net  | NC |
|                          | 320EX10 | Harness for 120A Expandable Controller  | NC |

# Joystick / Electronics

## Seating Control

### Powered Seating Control

|                          |         |  |            |
|--------------------------|---------|--|------------|
| <input type="checkbox"/> | 320EO40 | Omit Thru-drive ( <i>Requires Switch Control</i> ) | NC         |
| <input type="checkbox"/> | 320EO37 | Thru-drive Control for 1 Actuator                  | \$1,332.00 |
| <input type="checkbox"/> | 320EO38 | Thru-drive Control for 2+ Actuators                | \$2,679.00 |

### Switch Control (*Available with Omit Thru-Drive Only*) (1 actuator for Tilt or Elevate Only)

|                          |         |  |    |
|--------------------------|---------|--|----|
| <input type="checkbox"/> | 320EO33 | Dual Toggle for 1 actuator   | NC |
| <input type="checkbox"/> | 320EO34 | Buddy Button for 1 actuator ( <i>Button Switches are not mounted</i> ) | NC |
| <input type="checkbox"/> | 320EO35 | Switch-It Mini 1 button for 1 actuator                                 | NC |
| <input type="checkbox"/> | 320EO36 | Switch-It Mini 2 button for 1 actuator                                 | NC |
| <input type="checkbox"/> | 320EO41 | Ctrl+5 Switch Box for 2 actuators or more                              | NC |

### Additional Power Seating Control (*Available with Thru-Drive Only*)

|                          |         |  |          |
|--------------------------|---------|--|----------|
| <input type="checkbox"/> | 320QG11 | Ctrl+5 Switch Box With Free Assignable Functions                             | \$995.00 |
| <input type="checkbox"/> | 320QG12 | Ctrl+5 Switch Box With Free Assignable Functions ( <i>Multi-Power only</i> ) | NC       |

### Dual Toggle Mounting Location

|                          |        |                     |    |
|--------------------------|--------|---------------------|----|
| <input type="checkbox"/> | 320TN1 | Toggle Mount, Right | NC |
| <input type="checkbox"/> | 320TN2 | Toggle Mount, Left  | NC |

### Ctrl+5 Mounting Location

|                          |         |                  |    |
|--------------------------|---------|------------------|----|
| <input type="checkbox"/> | 320EY19 | Right Armrest    | NC |
| <input type="checkbox"/> | 320EY20 | Left Armrest     | NC |
| <input type="checkbox"/> | 320EY21 | Back Cane, Right | NC |
| <input type="checkbox"/> | 320EY24 | Back Cane, Left  | NC |

### Ctrl+5 Switch Type

|                          |         |                 |          |
|--------------------------|---------|-----------------|----------|
| <input type="checkbox"/> | 320EY22 | Push Buttons    | NC       |
| <input type="checkbox"/> | 320EY23 | Toggle Switches | \$328.00 |

## Attendant Controls (Prior Approval Required)

|                          |         |  |            |
|--------------------------|---------|--|------------|
| <input type="checkbox"/> | 3200A05 | Dual Control Module for Attendant                        | \$1,131.00 |
| <input type="checkbox"/> | 3200A06 | R-net Color LCD ( <i>CJSM1</i> ), Attendant Control Only | \$1,131.00 |

## Attendant Mounting Location

|                          |        |                 |    |
|--------------------------|--------|-----------------|----|
| <input type="checkbox"/> | 320DN5 | Right Back Cane | NC |
| <input type="checkbox"/> | 320DN6 | Left Back Cane  | NC |

## R-net Programming

### Assignable Buttons

*Assignable Button feature allows for customization of joystick buttons and/or Ctrl+5 Switches allowing access to various wheelchair and seating functions. When enabled, chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle.*

|                          |         |                            |    |
|--------------------------|---------|----------------------------|----|
| <input type="checkbox"/> | 320PM12 | Include Assignable buttons | NC |
|--------------------------|---------|----------------------------|----|

# Specialty Input Devices (AADL W100 Grant)

## R-net Omni2 And Output Options (For Expandable Controls, Choose All That Apply) (AADL W100 Grant)

|                          |         |  |            |
|--------------------------|---------|--|------------|
| <input type="checkbox"/> | 320ON14 | R-Net Omni2 Only Input Control Module and 3.5" LCD Display w/ Built-In Infrared and Bluetooth Choose Mount Below                                       | \$2,201.00 |
| <input type="checkbox"/> | 320ON15 | R-Net Omni2 Bundled Input Control Module And 3.5" LCD Display w/ Built-In Infrared and Bluetooth. (Bundled Omni2 available w/ Switch-It Input Devices) | NC         |

## Omni2 Mount (If Used with Standard Joystick, SunDisk, or Bullet Hand Tray, Must Be Mounted Opposite Side) (AADL W100 Grant)

|                          |         |   |          |
|--------------------------|---------|---|----------|
| <input type="checkbox"/> | 320DS14 | Link-It with Parallel Swing Away Omni2 Mount w/ Rear Mounted Interface, Right | \$644.00 |
| <input type="checkbox"/> | 320DS15 | Link-It with Parallel Swing Away Omni2 Mount w/ Rear Mounted Interface, Left  | \$644.00 |

## Hand Control Packages (Requires Standard Joystick Omit / Only one hand control device can be ordered) (AADL W100 Grant)

### Proportional Input Devices

|                          |         |                         |            |
|--------------------------|---------|-------------------------|------------|
| <input type="checkbox"/> | 320EU14 | Switch-It MicroGuide    | \$4,284.00 |
| <input type="checkbox"/> | 320EU15 | Switch-It MicroPilot    | \$4,284.00 |
| <input type="checkbox"/> | 320EU17 | Switch-It VersaGuide    | \$2,360.00 |
| <input type="checkbox"/> | 320EU18 | Switch-It VersaGuide EZ | \$2,789.00 |

### Mounting Position (Required)

|                          |         |   |            |
|--------------------------|---------|---|------------|
| <input type="checkbox"/> | 320EM04 | Link-It to the Backrest (Not Available w/ Omit Back Cushion Type) | \$1,235.00 |
| <input type="checkbox"/> | 320EM05 | Link-It to the Fixed Joystick Mount                               | \$384.00   |
| <input type="checkbox"/> | 320EM06 | Link-It to the Parallel Swing Away Joystick Mount                 | \$644.00   |
| <input type="checkbox"/> | 320EM07 | Link-It to the Midline Swing Away                                 | \$1,365.00 |
| <input type="checkbox"/> | 320EM08 | Omit Mount  | NC         |

### Mounting Platforms (Optional) Not available with Link-It to the Backrest

|                          |         |   |          |
|--------------------------|---------|---|----------|
| <input type="checkbox"/> | 320MG16 | SunDisk Hand Tray (Only available w/ MicroSeries joysticks) | \$384.00 |
| <input type="checkbox"/> | 320MG15 | Bullet Hand Tray  | \$384.00 |

### Mounting Location (Required)

|                          |         |       |    |
|--------------------------|---------|-------|----|
| <input type="checkbox"/> | 320IM09 | Right | NC |
| <input type="checkbox"/> | 320IM10 | Left  | NC |

### Switch Options

See Switch Selections (Maximum 2)

## Head Control Packages (Select headrest mount on next page) (AADL W100 Grant)

Only one head control device can be ordered.

### Proportional Input Devices

|                          |         |  |            |
|--------------------------|---------|--|------------|
| <input type="checkbox"/> | 320EU25 | Switch-It DualPro Proportional Head Array - Facial Pads                              | \$6,513.00 |
| <input type="checkbox"/> | 320EU44 | Switch-It DualPro Proportional Head Array - Facial Pads - Dedicated Mode Outer Left  | \$6,513.00 |
| <input type="checkbox"/> | 320EU45 | Switch-It DualPro Proportional Head Array - Facial Pads - Dedicated Mode Outer Right | \$6,513.00 |
| <input type="checkbox"/> | 320EU46 | Switch-It DualPro Proportional Head Array - Facial Pads - Dedicated Mode & Reverse   | \$6,513.00 |
| <input type="checkbox"/> | 320EU26 | Switch-It DualPro Proportional Head Array - Spot Pads                                | \$6,513.00 |
| <input type="checkbox"/> | 320EU27 | Switch-It RIM Proportional Head Control  | \$6,513.00 |

### Switched Input Devices

|                          |         |   |            |
|--------------------------|---------|---|------------|
| <input type="checkbox"/> | 320EU28 | Switch-It 3-Switch Head Array - Short Curved Wings                          | \$5,679.00 |
| <input type="checkbox"/> | 320EU29 | Switch-It 3-Switch Head Array - Short Straight Wings                        | \$5,679.00 |
| <input type="checkbox"/> | 320EU30 | Switch-It 3-Switch Head Array - Long Curved Wings                           | \$5,679.00 |
| <input type="checkbox"/> | 320EU31 | Switch-It 3-Switch Head Array - Long Straight Wings                         | \$5,679.00 |
| <input type="checkbox"/> | 320EU42 | Switch-It 4-Switch Head Array - Long Straight Wings Only - Mode Outer Left  | \$5,679.00 |
| <input type="checkbox"/> | 320EU43 | Switch-It 4-Switch Head Array - Long Straight Wings Only - Mode Outer Right | \$5,679.00 |
| <input type="checkbox"/> | 320EU33 | Switch-It 5-Switch Head Array - Long Straight Wings Only - Mode Outer Left  | \$5,679.00 |
| <input type="checkbox"/> | 320EU47 | Switch-It 5-Switch Head Array - Long Straight Wings Only - Mode Outer Right | \$5,679.00 |

### Switch Options

See Switch Selections (maximum 2)

## Specialty Input Devices (AADL W100 Grant)

### Sip 'N Puff (AADL W100 Grant)

|                          |         |   |            |
|--------------------------|---------|---|------------|
| <input type="checkbox"/> | 320ID73 | Switch-It Sip 'N Puff Combo Head Array with Link-It Hardware <i>(Requires headrest mount)</i> | \$6,494.00 |
| <input type="checkbox"/> | 320ID74 | Switch-It OptiStop with Sip 'N Puff   | \$2,529.00 |

### Headrest Mounts *(Applies to Head/Sip 'N Puff control packages)* (Seating Benefits)

|                          |         |  |          |
|--------------------------|---------|--|----------|
| <input type="checkbox"/> | 320IJ13 | Whitmyer Linx Headrest Mount with Quick Release        | \$215.00 |
| <input type="checkbox"/> | 320IJ14 | Whitmyer Linx-2 Headrest Mount with Quick Release      | \$260.00 |
| <input type="checkbox"/> | 320IJ15 | Whitmyer AXYS Headrest Mount with Quick Release        | \$215.00 |
| <input type="checkbox"/> | 320IJ16 | Whitmyer AXYS Offset Headrest Mount with Quick Release | \$215.00 |
| <input type="checkbox"/> | 320IJ17 | Whitmyer Cobra FB Headrest Mount with Quick Release    | \$515.00 |

### Chin Control Packages (AADL W100 Grant)

#### Proportional Input Devices

|                          |         |   |            |
|--------------------------|---------|---|------------|
| <input type="checkbox"/> | 320EU34 | Switch-It MicroGuide <i>(Switches mounted to the headrest only)</i> | \$4,284.00 |
| <input type="checkbox"/> | 320EU35 | Switch-It MicroPilot <i>(Switches mounted to the headrest only)</i> | \$4,284.00 |
| <input type="checkbox"/> | 320EU36 | Switch-It VersaGuide  | \$2,360.00 |
| <input type="checkbox"/> | 320EU37 | Switch-It VersaGuide EZ   | \$2,789.00 |

#### Mounting Position

|                          |         |  |            |
|--------------------------|---------|--|------------|
| <input type="checkbox"/> | 320EM09 | Link-It To Chin Harness, Small <i>(150mm / 6")</i> | \$1,625.00 |
| <input type="checkbox"/> | 320EM10 | Link-It To Chin Harness, Large <i>(220mm / 8")</i> | \$1,625.00 |
| <input type="checkbox"/> | 320EM11 | Link-It To Headrest <i>(Swing Away)</i>            | \$1,625.00 |
| <input type="checkbox"/> | 320EM53 | Link-It To Manual Chin Boom - Right Side           | \$4,934.00 |
| <input type="checkbox"/> | 320EM54 | Link-It To Manual Chin Boom - Left Side            | \$4,934.00 |
| <input type="checkbox"/> | 320EM16 | Omit Mount   | NC         |

#### Input Device Mounting Location

|                          |         |       |    |
|--------------------------|---------|-------|----|
| <input type="checkbox"/> | 320IM13 | Right | NC |
| <input type="checkbox"/> | 320IM14 | Left  | NC |

#### Switch Options

|   |   |   |
|---|---|---|
| - | See Switch Selections <i>(Minimum 1 for the Power Chin Boom, Maximum 3)</i> | - |
|---|---|---|

### Safety Devices (AADL W100 Grant)

|                          |         |  |            |
|--------------------------|---------|--|------------|
| <input type="checkbox"/> | 320SA06 | Switch-It Remote Stop Switch - XLR Connector | \$1,216.00 |
|--------------------------|---------|--|------------|

# Specialty Input Devices Cont.

## Switched Driving Packages (AADL W100 Grant)

|                          |         |  |            |
|--------------------------|---------|--|------------|
| <input type="checkbox"/> | 320SW44 | Switch-It Cool Cube with DB9 connector for proximity / Proportional switch driving.<br><i>(Pro Spot Proportional Switch's are required if proportional driving is desired)</i><br><i>(For Non-Proportional driving, Single Switches can be used)</i> | \$3,645.00 |
|--------------------------|---------|--|------------|

## Switches

Single switches can be programmed to activate any function of the wheelchair including: On/Off, Mode, Profile, Speeds, Seating, Memory Positions, Horn, Lights & Indicators, IR functions, Bluetooth Mouse Clicks, Environmental Controls when plugged into assignable port

Switches can connect via: Joysticks (Standard, Micro series, Versa series, All-round series), Omni2, Cool Cube, Head Controls

Switch ordering instructions:

1. Select Switch type and quantity in section "A" using columns below. *(Max 6 switches. see limitations below)*
2. Select Switch Mounting location from Section "B" using corresponding switch column.
3. Select Switch Mounting Side from Section "C" using corresponding switch column.

### Section A - Switch Selections

#### Single Function Switches - Mechanical Switches

|  | Switch / Quantity        |                          |                          |                          |                          |                          | \$       |
|--|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|----------|
|  | Switch 1                 | Switch 2                 | Switch 3                 | Switch 4                 | Switch 5                 | Switch 6                 |          |
| 2.5" Single Switch - Black               | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | \$101.00 |
| 2.5" Single Switch - Blue                | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | \$101.00 |
| 2.5" Single Switch - Green               | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | \$101.00 |
| 2.5" Single Switch - Red                 | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | \$101.00 |
| 2.5" Single Switch - Yellow              | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | \$101.00 |
| Egg Switch - Black                       | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | \$101.00 |
| Egg Switch - Blue                        | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | \$101.00 |
| Egg Switch - Green                       | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | \$101.00 |
| Egg Switch - Red                         | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | \$101.00 |
| Egg Switch - Yellow                      | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | \$101.00 |
| MicroLite <i>(Includes color labels)</i> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | \$284.00 |
| Mini Button 1 - Blue <i>(Velcro)</i>     | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | \$284.00 |
| Mini Button 1 - Green <i>(Velcro)</i>    | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | \$284.00 |
| Mini Button 1 - Red <i>(Velcro)</i>      | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | \$284.00 |
| 1" Single Switch - Black                 | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | \$358.00 |
| 1" Single Switch - Blue                  | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | \$358.00 |
| 1" Single Switch - Green                 | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | \$358.00 |
| 1" Single Switch - Orange                | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | \$358.00 |
| 1" Single Switch - Red                   | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | \$358.00 |
| 1" Single Switch - Yellow                | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | \$358.00 |

#### Electronic Switches

|  |                          |                          |                          |                          |                          |                          |          |
|--|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|----------|
| Flex Switch <i>(Headrest mount only)</i> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | \$280.00 |
| Zero Touch Proximity Sensor              | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | \$254.00 |

#### Pro Spot - Proportional Switches

|                                       |                          |                          |                          |                          |                          |                          |          |
|---------------------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|----------|
| Pro Spot Proportional Switch - Red    | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | \$500.00 |
| Pro Spot Proportional Switch - Green  | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | \$500.00 |
| Pro Spot Proportional Switch - Blue   | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | \$500.00 |
| Pro Spot Proportional Switch - Yellow | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | \$500.00 |
| Pro Spot Proportional Switch - Orange | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | \$500.00 |

### Section B - Switch Mounting Position

|  |                          |                          |                          |                          |                          |                          |          |
|--|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|----------|
| Link-It to the Backrest <i>(Up to 1 per side)</i>  | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | \$845.00 |
| Link-It to the Headrest - Gooseneck & Rigid <i>(Up to 2 per side)</i>                          | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | \$845.00 |
| Link-It to the Fixed Joystick Mount <i>(Up to 2 per side)</i>                                  | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | \$644.00 |
| Link-It to the Parallel Swing Away Joystick Mount <i>(Up to 2 per side)</i>                    | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | \$644.00 |
| Link-It to the Seat Rail <i>(Up to 2 per side)</i><br><i>(N/A w/ Thigh &amp; Knee Support)</i> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | \$514.00 |
| Omit Mount <i>(Includes Velcro)</i>  | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | NC       |

### Section C - Switch Mounting Location

|       |                          |                          |                          |                          |                          |                          |    |
|-------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|----|
| Right | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | NC |
| Left  | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | NC |

### Joystick Integrated Switches (Pair)

|                          |         |   |  |            |
|--------------------------|---------|---|--|------------|
| <input type="checkbox"/> | 320JO03 | - | Proximity Switches on Flex Mount, Pair <i>(Used with VersaGuide/VersaGuide EZ)</i> | \$1,294.00 |
|--------------------------|---------|---|--|------------|

### Additional Switches (velcro)

|                          |         |  |          |
|--------------------------|---------|--|----------|
| <input type="checkbox"/> | 320SA08 | Mini 2 - Button <i>(Stereo Plug)</i>           | \$284.00 |
| <input type="checkbox"/> | 320SA09 | Mini 2 - Button <i>(Dual Mono Plug)</i>        | \$284.00 |
| <input type="checkbox"/> | 320SA10 | Mini 4 - Button In-Line <i>(DB9 connector)</i> | \$284.00 |