

February 2024

**Customer Contact Details**

Customer Name: _____	Ship to Customer Address: <input type="checkbox"/>
Advisor: _____	Shipping Name: _____
Address: _____	Shipping Address: _____
_____	_____
_____	_____
Phone: _____	Phone: _____
Email: _____	Email: _____

Additional Information

Notes: _____	Current Wheelchair: _____
_____	Customer Weight: _____
_____	Customer Height: _____
_____	Customer Disability: _____

Guide:

- Base product
- Required option
- Optional Extra
- Critical Information
- Information Needed

PRODUCT STYLE

FXTXZ010 OCTANE FX
\$13,488.00

XA070 Open design

Max User Weight 110kg (242lbs)

Yes, I am within the max user weight

CAD APPROVAL

<p><input type="checkbox"/> XW012 QUOTE CAD drawing for approval \$550 FLAT FEE (Includes two free revisions)</p> <p><input type="checkbox"/> XW009 QUOTE only no CAD \$0</p>	<p><input type="checkbox"/> XW010 ORDER CAD drawing for approval \$300 FLAT FEE (Includes two free revision)</p> <p><input type="checkbox"/> XW011 ORDER only no CAD \$0</p>
---	--

KEY MEASUREMENTS

<p><input type="checkbox"/> XY010 Seat width (Inside of sideguards)</p> <p> <input type="text" value="mm"/> (Max - 460mm) (Min - 250mm)</p> <p><input type="checkbox"/> XY020 Seat length</p> <p> <input type="text" value="mm"/> (Max - 460mm) (Min - 250mm)</p> <p><input type="checkbox"/> XY030 Extra seat length (Provides additional leg support & hand grip)</p> <p> <input type="text" value="mm"/></p> <p><input type="checkbox"/> XY050 Front frame height</p> <p> <input type="text" value="mm"/> (Max - 550mm) (Min - 480mm)</p> <p><input type="checkbox"/> XY060 Rear frame height</p> <p> <input type="text" value="mm"/> (Max - 550mm) (Min - 360mm)</p> <p><input type="checkbox"/> XY070 Frame length</p> <p> <input type="text" value="mm"/> (Max 625mm) (Min XY020 + XY030 + 140mm)</p> <p><input type="checkbox"/> XY072 Caster position</p> <p> <input type="text" value="mm"/></p> <p><input type="checkbox"/> XA020 Ergonomic seat \$119.00</p> <p> Specify length</p> <p><input type="text" value="mm"/> (Standard = 120mm)</p>	
---	--

Front Frame Options

XY080 Seat taper design

Specify width

mm

Optimize measurements
(XY010 Seat Width -100mm)

Footrest Design (Choose One)

XY100 Round design

Caster Size	Max Height	With XB200
3"	115mm	165mm
4"	125mm	175mm
5"	135mm	185mm

XB200 High level footrest (Only with XB010, XB080 & XB150)
\$198.00

XY105 Square design

Caster Size	Max Height	With XB200
3"	135mm	185mm
4"	145mm	195mm
5"	155mm	205mm

XB200 High level footrest (Only with XB010, XB080 & XB150)
\$198.00

Footrest Individual Measurements

XY115 Front footrest height

mm

XY120 Rear footrest height

mm

(Std = 20 lower than front height)

XY125 Footrest length

mm

(Std = 100mm)
(160mm and over requires heavy duty upgrade XB150)

Footrest Style and Cover

FTXB080 Tubular footrest welded w/ Aluminium plate
(Supplied with spacer for +20mm adjustment)

FTXB110 Tubular footrest welded w/ Carbon Fibre plate
\$321.00 (Supplied with spacer for +20mm adjustment)

FTXB120 Clamp-on aluminium plate footrest
(Angle/height adjustable)

Footrest Options

XB165 Rubber grip tape (3M safety walk)
\$85.00

CASTER WHEELS, FORKS & WHEEL LOCKS

Caster Wheels

75mm (3")

- XC150 Aluminium centre 3"x1.4"
\$229.00 (Silver hub and black urethane tire)

100mm (4")

- XC060 Soft roll 4"x1.3"
(Black)
- XC110 Aluminium centre 4"x1"
\$229.00 (Silver hub and black urethane tire)
- XC130 Aluminium centre 4"x1.4"
\$229.00 (choose colour below)



- | | | | |
|--|--|--|---|
| <input type="checkbox"/> Silver hub & black tire | <input type="checkbox"/> Black hub & blue tire | <input type="checkbox"/> Purple hub & black tire | <input type="checkbox"/> Silver hub & blue tire |
| <input type="checkbox"/> Black hub & black tire | <input type="checkbox"/> Red hub & black tire | <input type="checkbox"/> Gold hub & black tire | <input type="checkbox"/> Silver hub & red tire |
| <input type="checkbox"/> Black hub & red tire | <input type="checkbox"/> Blue hub & black tire | <input type="checkbox"/> Green hub & black tire | |

125mm (5")

- | | |
|--|--|
| <input type="checkbox"/> XC070 Soft roll 5"x1.3"
(Black) | <input type="checkbox"/> XC120 Aluminium centre 5"x1"
\$229.00 (Silver hub and black urethane tire) |
| <input type="checkbox"/> XC100 Pneumatic Aluminium centre 5"x1"
\$244.00 (Grey) | <input type="checkbox"/> XC140 Aluminium centre 5"x1.4"
\$229.00 (Silver hub and black urethane tire) |

150mm (6")

- | | |
|--|--|
| <input type="checkbox"/> XC080 Soft roll 6"x1.3" | <input type="checkbox"/> XC145 Aluminium centre 6"x1.4"
\$305.00 (Silver hub and black urethane tire) |
|--|--|

Caster Forks

- | | |
|--|--|
| <input type="checkbox"/> XC190 Low ride caster fork
(Std when no option selected) | <input type="checkbox"/> XC210 Froglegs suspension caster fork
\$624.00 |
|--|--|


Note: Colour follows accent



- Black
- Silver

Caster Pin Width

- XY330 Caster pin width
(XY010 plus 70mm)



Wheel Lock Options (Choose one)




- | | |
|--|---|
| <input type="checkbox"/> XD020 Underseat Scissor | <input type="checkbox"/> XD050 Push to lock |
| <input type="checkbox"/> XD030 Compact scissor
\$262.00 | <input type="checkbox"/> XD046 Out Front push to lock
\$403.00 |
| <input type="checkbox"/> XD035 Compact scissor with lightweight fixing
\$347.00 (Only with fixed centre of gravity) | |
| <input type="checkbox"/> XD045 Out Front Compact scissor
\$403.00 | <input type="checkbox"/> XD010 Omit wheel locks |

SIDEGUARDS AND ARMRESTS

Sideguard		Sideguard with Fender	
Aluminium Options		Aluminium Options	
<input type="checkbox"/> XE080	Lift out aluminium \$524.00	<input type="checkbox"/> XE090	Lift out aluminium with fender \$570.00
<input type="checkbox"/> FXTXE130	Fold flat aluminium \$730.00	<input type="checkbox"/> FXTXE140	Fold flat aluminium with fender \$842.00
Carbon Fibre Options		Carbon Fibre Options	
<input type="checkbox"/> XE100	Lift out carbon fibre \$755.00	<input type="checkbox"/> XE110	Lift out carbon fibre with fender \$1,115.00
<input type="checkbox"/> FXTXE150	Fold flat carbon fibre \$1,460.00	<input type="checkbox"/> FXTXE160	Fold flat carbon fibre with fender \$1,607.00
ABS Plastic Options		Sideguard Measurements	
<input type="checkbox"/> XE120	Lift out ABS plastic \$480.00	Fender height	
<input type="checkbox"/> FXTXE170	Fold flat ABS plastic \$730.00	<input type="checkbox"/> 10mm above tire <i>(Standard)</i>	
Armrest Options		<input type="checkbox"/> Individual above tire <i>(Min = 10mm)</i>	
<input type="checkbox"/> XF010	Lift out armrest with plastic sideguard \$583.00 <i>(Indicate fixed height)</i>	<input type="checkbox"/> <input type="text" value="mm"/>	
<input type="checkbox"/> XF020	Lift out armrest with aluminium sideguard \$460.00 <i>(Indicate fixed height)</i>	Fender position	
Sideguard Measurements		<input type="checkbox"/> Centre of tire <i>(Standard if no option selected)</i>	
Optimized measurements		<input type="checkbox"/> Outside of tire	
<input type="checkbox"/> 25mm above tire <i>(Follow tire)</i>		Individual measurements	
<input type="checkbox"/> Individual above tire <i>(Follow tire)</i>		<input type="checkbox"/> Sideguard length	
<input type="checkbox"/> Level with tire <i>(Follow tire)</i>		<input type="checkbox"/> Sideguard height	
<input type="checkbox"/> <input type="text" value="mm"/>		<input type="checkbox"/> <input type="text" value="mm"/>	

BACKREST	
Backrest Height	
<input type="checkbox"/> XY230 Fixed height backrest (Maximum 400mm) Specify height <div style="border: 1px solid black; width: 100px; height: 20px; margin-top: 5px;"></div> mm	<input type="checkbox"/> XY245 Adjustable height backrest \$244.00 (Maximum height 400mm) (+/-25mm) Specify height <div style="border: 1px solid black; width: 100px; height: 20px; margin-top: 5px;"></div> mm
Backrest Design	
<input checked="" type="checkbox"/> FXTXG060 Folding backrest with Q Lock (Black)	
Ergonomic Options	
<input checked="" type="checkbox"/> XY360 Backrest angle -5 -3 0 +3 +5 <div style="display: flex; justify-content: space-around; width: 100px;"> <div style="border: 1px solid black; width: 15px; height: 15px;"></div> <div style="border: 1px solid black; width: 15px; height: 15px;"></div> <div style="border: 1px solid black; width: 15px; height: 15px;"></div> <div style="border: 1px solid black; width: 15px; height: 15px;"></div> <div style="border: 1px solid black; width: 15px; height: 15px;"></div> </div> Individual _____ Deg (Max -7/+7 Deg)	<input type="checkbox"/> XG100 Ergonomic backrest \$144.00 Vertical height <div style="border: 1px solid black; width: 100px; height: 20px; margin-top: 5px;"></div> mm Rearward position <div style="border: 1px solid black; width: 100px; height: 20px; margin-top: 5px;"></div> mm
<input type="checkbox"/> XG110 Inward taper (inside of backrest posts) Specify width or per side <div style="border: 1px solid black; width: 100px; height: 20px; margin-top: 5px;"></div> mm	<input type="checkbox"/> XG120 Outward flare (inside of backrest posts) Specify width or per side <div style="border: 1px solid black; width: 100px; height: 20px; margin-top: 5px;"></div> mm
Push handle Options	
<input type="checkbox"/> XG070 Screw in push handles \$222.00 (N/A with XY245) <input type="checkbox"/> XG080 Flip down push handles \$434.00 (N/A with XY245) <input type="checkbox"/> XG085 Flip down push handles height adjustable \$709.00 (N/A with XY245) <input type="checkbox"/> XG097 Integrated push handle \$666.00 (quickie design)	

CAMBER BAR

Adjustable Centre of Gravity	Fixed Centre of Gravity
<p><input type="checkbox"/> XH110 Square back adjustable centre of gravity (25-150mm)</p> <p>(Max seat width 500mm)</p> <p>Specify measurement <input style="width: 100px;" type="text"/> mm</p>  <p style="text-align: center;">XH110</p> <p><input type="checkbox"/> XH120 Round back adjustable centre of gravity (60-150mm)</p> <p>Specify measurement <input style="width: 100px;" type="text"/> mm</p>  <p style="text-align: center;">XH120</p>	<p><input type="checkbox"/> FXTXH090 V-style fixed centre of gravity (25-125mm)</p> <p>\$716.00</p> <p>Specify measurement <input style="width: 100px;" type="text"/> mm</p>  <p style="text-align: center;">FXAXH090</p>
Specify Camber Angle	Specify Camber Angle
<p><input checked="" type="checkbox"/> XY390 Camber angle</p> <p><input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3 <input type="checkbox"/> 5</p>	<p><input checked="" type="checkbox"/> XY390 Camber angle</p> <p><input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3 <input type="checkbox"/> 5 <input style="width: 50px;" type="text"/> Individual</p>
Camber Bar Options	
<p><input type="checkbox"/> FXTXK227 Standard Axle pin (120 or 130mm)</p> <p><input type="checkbox"/> FXTXK226 100mm Titanium Axle pin (Fixed Axle only) \$446.00</p>	<p><input type="checkbox"/> XK220 Titanium 130mm quick release axle pins (pair) \$158.00 (N/A with XK010, XK070 or FXTXK226)</p> <p><input type="checkbox"/> XK240 Quad release axle nuts (pair) \$237.00 (N/A with FXTXK226)</p>
Anti-Tip Options	
<p><input type="checkbox"/> XJ010 Swingaway anti tip LEFT side \$165.00 (Only with square back frame. Not with XH040)</p> <p><input type="checkbox"/> XJ020 Swingaway anti tip RIGHT side \$165.00 (Only with square back frame. Not with XH040)</p>	<p><input type="checkbox"/> XJ100 Single sport anti tip \$497.00 (centre mounted)</p>

REAR WHEELS			
Rear Wheel Size			
<input type="checkbox"/> XY420 24" Wheels			<input type="checkbox"/> XY430 25" Wheels
Rear Wheel Gap			
<div style="display: flex; align-items: center;"> <div style="margin-right: 20px;"> <input checked="" type="checkbox"/> XY460 Rear wheel gap <div style="display: flex; gap: 10px; margin-top: 5px;"> <input type="checkbox"/> 15 <input type="checkbox"/> 20 <input type="checkbox"/> 25 </div> <div style="border: 1px solid black; width: 100px; height: 20px; margin-top: 5px; display: flex; align-items: center; justify-content: center;"> Individual mm </div> </div> <div style="margin-left: 20px;"> </div> </div>			
Rear Wheel Options			
<input type="checkbox"/> XK010 Fusion 16 spoke <i>(Silver Hub, black spokes)</i>	<input type="checkbox"/> XK320 Spinergy XLX 12 spoke, cross/radial spoked \$1,051.00		
<input type="checkbox"/> XK050 Spinergy LX 12 spoke \$985.00	<input type="checkbox"/> XK060 Spinergy SPOX 18 spoke \$655.00		
<input type="checkbox"/> XK150 Spinergy CLX 18 carbon blade spoke \$3,093.00 <i>(N/A in silver hub)</i>	<input type="checkbox"/> XK070 Spinergy BLXL 18 blade spoke \$1,051.00 <i>(N/A in silver hub)</i>		
Spinergy Colour Options:			
<input type="checkbox"/> Black spokes	<input type="checkbox"/> White spokes	<input type="checkbox"/> Yellow spokes	<input type="checkbox"/> XK180 Individual spoke design \$210.00 <i>(Specify colour)</i>
<input type="checkbox"/> Red spokes	<input type="checkbox"/> Blue spokes	<input type="checkbox"/> Pink spokes	<i>Inside:</i> _____
<input type="checkbox"/> Purple spokes	<input type="checkbox"/> Orange spokes	<input type="checkbox"/> Green spokes	<i>Outside:</i> _____
Spinergy Anodise Options			
<input type="checkbox"/> XK190 Black hub <i>(std option)</i>	<input type="checkbox"/> XK200 Silver hub \$125.00 <i>(n/a with BLXL/CLX)</i>	<input type="checkbox"/> XK201 Red hub/rims \$403.00 <i>(LX, XLX & BLXL only)</i>	<input type="checkbox"/> XK205 Individual hub/rims \$560.00 Select colour: _____
<input type="checkbox"/> XK202 Purple hub/rims \$403.00 <i>(LX, XLX & BLXL only)</i>	<input type="checkbox"/> XK203 Blue hub/rims \$403.00 <i>(LX, XLX & BLXL only)</i>	<input type="checkbox"/> XK204 Gold hub/rims \$531.00 <i>(LX, XLX & BLXL only)</i>	
Supply Without Wheels			
<input type="checkbox"/> XK040 Supply without Wheels	<input style="width: 150px;" type="text" value="Wheel Type:"/>		<input style="width: 150px;" type="text" value="Tire Type:"/>
Handrim Position / Options			
<input type="checkbox"/> XY480 Narrow set handrim			<input type="checkbox"/> XY485 Wide set handrim
<input type="checkbox"/> XL090 Remove excess mounting tabs \$144.00 <i>(recommended for small wheel gaps)</i>			
Handrim Styles			
<input type="checkbox"/> XL010 Aluminium hard anodized 6 tab <i>(Standard if no option selected)</i>			<input type="checkbox"/> XL040 Titanium 6 tab \$439.00
<input type="checkbox"/> XL020 Tetra grip 6 tab \$394.00			<input type="checkbox"/> XL050 Stainless steel 6 tab \$244.00
<input type="checkbox"/> XL100 Para grip 6 tab \$394.00			
Tire Options			
<input type="checkbox"/> XM010 Kenda	<input type="checkbox"/> Grey	<input type="checkbox"/> Blue	<input type="checkbox"/> Red
<input type="checkbox"/> XM130 Kenda black \$27.00			
<input type="checkbox"/> XM020 Schwalbe Marathon Plus \$97.00			<input type="checkbox"/> XM060 Schwalbe One <i>(24", 25" only)</i> \$131.00
<input type="checkbox"/> XM040 Shox solid airless (Black) \$158.00			<input type="checkbox"/> XM050 Kenda K885 mountain bike tire <i>(24", 25" only)</i> \$39.00
Inner Tube Valve Options			
<input type="checkbox"/> XY500 Schraeder inner tube valve <i>(Standard if no option selected)</i>			<input type="checkbox"/> XY490 Presta inner tube valve
Wheel Accessories			
<input type="checkbox"/> XK310 Offroad wheel set 24" rim with handrim, MTB tire and axle pins (Replaces 25" wheel with XM010) \$919.00			<input type="checkbox"/> XK300 Offroad wheel set 22" rim with handrim, MTB tire and axle pins (Replaces 24" wheel with XM010) \$919.00

FINISH OPTIONS			
Frame & Component Finish			
<input type="checkbox"/> XN010 Powder coated frame	<input type="checkbox"/> XN150 Powder coated components <i>(Frame colour if not specified)</i>		
Frame Comp.	Gloss	Frame Comp.	Matt
<input type="checkbox"/> <input type="checkbox"/> Black gloss 103 <input type="checkbox"/> <input type="checkbox"/> White gloss 100 <input type="checkbox"/> <input type="checkbox"/> Grey gloss 132 <input type="checkbox"/> <input type="checkbox"/> Cream gloss 135 <input type="checkbox"/> <input type="checkbox"/> Red gloss 108 <input type="checkbox"/> <input type="checkbox"/> Yellow gloss 120 <input type="checkbox"/> <input type="checkbox"/> Forest green gloss 143 <input type="checkbox"/> <input type="checkbox"/> Green gloss 107 <input type="checkbox"/> <input type="checkbox"/> Dark blue gloss 151 <input type="checkbox"/> <input type="checkbox"/> Royal blue gloss 154 <input type="checkbox"/> <input type="checkbox"/> Blue gloss 174 <input type="checkbox"/> <input type="checkbox"/> Duckegg blue gloss 141 <input type="checkbox"/> <input type="checkbox"/> Pink gloss 188 <input type="checkbox"/> <input type="checkbox"/> Orange gloss 180 <input type="checkbox"/> <input type="checkbox"/> Purple gloss 157 <input type="checkbox"/> FXTXN030 Beadblast frame \$617.00 <i>(All parts black unless XN160 chosen)</i> <input type="checkbox"/> FXTXN050 Polished frame \$1,533.00 <i>(All parts black unless XN170 chosen)</i>	<input type="checkbox"/> <input type="checkbox"/> Black matt 105 <input type="checkbox"/> <input type="checkbox"/> White matt 101 <input type="checkbox"/> <input type="checkbox"/> Grey matt 131 <input type="checkbox"/> <input type="checkbox"/> Light grey matt 133 <input type="checkbox"/> <input type="checkbox"/> Navy blue matt 152 <input type="checkbox"/> XN200 Individual paint colour (State Colour/RAL code) \$557.00 _____	<input type="checkbox"/> <input type="checkbox"/> Chrome Mirror 300 <input type="checkbox"/> <input type="checkbox"/> Sparkle silver 170 <input type="checkbox"/> <input type="checkbox"/> Metallic black matt 104 <input type="checkbox"/> <input type="checkbox"/> Metallic grey matt 130 <input type="checkbox"/> <input type="checkbox"/> Black sandpaper 201 <input type="checkbox"/> <input type="checkbox"/> Grey sandpaper 202 <input type="checkbox"/> <input type="checkbox"/> White sandpaper 200 <input type="checkbox"/> XN160 Beadblast components \$119.00 <input type="checkbox"/> XN170 Polished components \$210.00	
Accent Colour (Choose One)			
<i>Note: Accent colour includes forks only</i>			
<input type="checkbox"/> Black accent <input type="checkbox"/> Red accent <input type="checkbox"/> Purple accent	<input type="checkbox"/> Silver accent <input type="checkbox"/> Blue accent <input type="checkbox"/> Carbon accent	<input type="checkbox"/> RGK-XN250 Individual Accent Colour \$310.00 _____	
UPHOLSTERY			
Seat Upholstery			
<input type="checkbox"/> XP040 Seat sling strap adjustable <input type="checkbox"/> FXTXP045 Sub4 Seat sling strap adjustable \$110.00	<input type="checkbox"/> XR050 Carbon fibre seat plate \$473.00		
Backrest Upholstery			
<input type="checkbox"/> XP050 Airtech backrest strap adjustable <input type="checkbox"/> FXTXP035 Airtech Sub4 backrest strap adjustable \$230.00	<input type="checkbox"/> XP060 Backrest upgrade to full leather \$493.00 <i>(Black)</i> <input type="checkbox"/> XP070 Backrest upgrade to leather corners \$119.00 <i>(Black)</i>		
Backrest Upholstery Stitching Colour			
<input type="checkbox"/> Silver <input type="checkbox"/> Red	<input type="checkbox"/> Black <input type="checkbox"/> Light blue	<input type="checkbox"/> Blue <input type="checkbox"/> Purple	<input type="checkbox"/> Individual: _____
Backrest Upholstery Logo Colour			
<input type="checkbox"/> White (Std) <input type="checkbox"/> Red	<input type="checkbox"/> Black <input type="checkbox"/> Light blue	<input type="checkbox"/> Blue <input type="checkbox"/> Purple	<input type="checkbox"/> Individual: _____
Backrest Accessories			
<input type="checkbox"/> XP110 Backrest with zip pocket \$125.00 <i>(N/A with Sub4 backrest)</i>	<input type="checkbox"/> XP170 Individual backrest embroidery \$236.00 <i>(Provide .jpeg file and positioning)</i>		

CUSHION								
Cushion foam options								
<input type="checkbox"/> XP080 \$119.00	Cushion foam <i>(select density and size below)</i>	<table style="width: 100%; border: none;"> <tr> <td style="width: 33%; text-align: center; padding: 5px;"><input type="checkbox"/> Soft (CMS040)</td> <td style="width: 33%; text-align: center; padding: 5px;"><input type="checkbox"/> Medium (CMS050)</td> <td style="width: 33%; text-align: center; padding: 5px;"><input type="checkbox"/> Hard (9lb recon)</td> </tr> <tr> <td style="text-align: center; padding: 5px;"><input type="checkbox"/> 50mm</td> <td style="text-align: center; padding: 5px;"><input type="checkbox"/> 75mm</td> <td style="text-align: center; padding: 5px;"><input type="checkbox"/> 100mm</td> </tr> </table>	<input type="checkbox"/> Soft (CMS040)	<input type="checkbox"/> Medium (CMS050)	<input type="checkbox"/> Hard (9lb recon)	<input type="checkbox"/> 50mm	<input type="checkbox"/> 75mm	<input type="checkbox"/> 100mm
<input type="checkbox"/> Soft (CMS040)	<input type="checkbox"/> Medium (CMS050)	<input type="checkbox"/> Hard (9lb recon)						
<input type="checkbox"/> 50mm	<input type="checkbox"/> 75mm	<input type="checkbox"/> 100mm						
Cushion cover upgrades								
<input type="checkbox"/> XP081	Nylon cushion cover <i>(Black)</i>							
<input type="checkbox"/> XP090 \$97.00	Air Lite cushion cover <i>(Black)</i>							
<input type="checkbox"/> XP100 \$250.00	Leather cushion cover <i>(Black)</i>							
Jay Products								
<input type="checkbox"/> JAY Cushion <i>(Price dependant on size / option)</i>	<input type="checkbox"/> <i>Jay Order form attached</i>	<div style="border: 1px solid black; height: 80px; width: 100%;"></div> <p style="font-size: small; margin-top: 5px;">Notes:</p>						
<input type="checkbox"/> JAY Backrest <i>(Price dependant on size / option)</i>	<input type="checkbox"/> <i>Jay Order form attached</i>							
<input type="checkbox"/> JAY Whitmyer Headrest <i>(Price dependant on size / option)</i>	<input type="checkbox"/> <i>Jay Order form attached</i>							
ACCESSORIES								
<input type="checkbox"/> XR080 \$1,037.00	Frontwheel assembly <i>(Only with tubular style footrest)</i>							
<input type="checkbox"/> XP120 \$83.00	Frame protectors Neoprene zip fasten <i>(Black)</i>							
<input type="checkbox"/> XP420 \$80.00	Frame protectors Nylon velcro fasten <i>(Black)</i>							
<input type="checkbox"/> XR010 \$58.00	Seat pouch velcro attach <i>(Black)</i>							
<input type="checkbox"/> XR200 \$149.00	Rucksack with RGK logo (Black)	NEW!						
<input type="checkbox"/> XR100 \$329.00	FX carry bag							
<input checked="" type="checkbox"/> XR090	Tool kit	NEW!						
<input checked="" type="checkbox"/> XS060	Calf strap							

NOTES
