



EIRGK02



Rider Information

Name: _____

Phone: _____

Email: _____

Height: _____ **Weight:** _____

Current Chair: _____ **Current Cushion:** _____

Ability/Diagnosis: _____

RGK Contact Authorization:

RGK would like to send you exclusive information, news and promotional offers on RGK.

I would like to be contacted by RGK:

Signature: _____ **Date:** _____

*We do not pass your personal information onto any third party.

Account Information

Account Name : _____

Account #: _____ **PO#:** _____

Contact Name : _____

Phone: _____

Email: _____

Ship to: **Account** **Rider**

Name: _____

Address: _____


RGK Genius: _____

Frame Styles

702XZ010 **Octane FX** (Max User Weight 110kg/242lbs.)
702XA070 **Open design**

702XW009 **Quote Only** (No CAD for Approval)
702XW011 **Order Only** (No CAD for Approval)
702XW012 **Quote CAD drawing for approval** (1st CAD is NCO)
702XW010 **Order CAD drawing for approval** (1st CAD is NCO)

CAD approved by: 702XW013 Customer 702XW014 RGK Genius



Open Design

\$14,101
STD

STD
STD
\$550
\$300

Measurements

Please Note: All RGK chairs are hand built and have a manufacturing tolerance of +/- 5mm (.25"). All measurements can be done in 5mm (.25") increments.

Please select either Millimeters or Inches.

Measurements: ☐ **Millimeters (mm.)** ☐ **Inches (in.)**

702XY010 Rear Seat Width Measured from inside of side guard to inside of opposite side guard. Or outside seat tube to outside of opposite seat tube near backrest. (Min. 250mm/9.84" and Max. 460mm/18.11")

mm. in. mm. in.
Current Chair: New Chair:

702XY020 Seat Depth/Seat Length Measured from front of back post to the front edge of seat sling. Frame bend will start .5" in front of seat sling unless Extra Frame Length requested. (Min. 250mm/9.84" and Max. 460mm/18.11")

Current Chair: New Chair:

702XY030 Extra Frame Length Measured from front edge of seat sling to beginning of Front frame bend. (Max. 460mm/18.11" Seat length + Extra seat length)



Current Chair: New Chair:

702XY050 Front Frame Height Measured from floor to apex of seat tube. (Min. 480mm/18.89" and Max. 550mm/21.65")

Current Chair: New Chair:

702XY060 Rear Frame Height Measured from floor to top of seat tube at back edge of seat sling. (Min. 360mm/14.17" and Max. 550mm/21.65")

Current Chair: New Chair:



Design Template

702XY070 Occupied Frame Length/Frame Length Measured from the front of the backrest tube to the front edge of the footrest. (Min. 390mm/15.35" and Max. 625mm/24.6")

Current Chair: New Chair:

702XY072 Caster Stem Position Measured from the front edge of the footrest to the centerline of the caster barrel. **702XY073** RGK to Optimize:

Current Chair: New Chair:

702XA020 Ergonomic Seat Measured from front of back post to beginning of upward bend of seat tube. \$125

Current Chair: New Chair:

702XY080 Front Seat Taper/Seat Taper Design Measured from inside of front frame tube at the front frame bend to inside of opposite front frame tube at the front frame bend. **702XY081** RGK to Optimize.

Current Chair: New Chair:

702XY115 Front Footrest Height Measured from ground up to the top of the Front footrest tube or cover.

Current Chair: New Chair:

702XY120 Rear Footrest Height Measured from ground up to the top of the Rear footrest tube or cover. Standard is 20mm less than Front Footrest Height.

Current Chair: New Chair:

702XY125 Footrest Length Measured from front edge of footrest to back edge of footrest. (160mm and over requires heavy duty upgrade)


Current Chair: New Chair:

702XB200 High Mount Footrest/High Level Footrest Measured from ground up to the top of the Front footrest tube or cover. \$207

Current Chair: New Chair:

702XY330 Stem Bolt/Caster Pin Width Measured from centerline of caster pin to the centerline of the opposite caster pin. **702XY331** RGK to Optimize:

Current Chair: New Chair:



Design Template

702XY231 **Backrest Height** Measured from top of seat tube at rear of frame to top of back-post.

Current Chair: New Chair:

702XY360 **Backrest Angle** Measured as the angle from top front of the backrest tube to the floor. (Max. -10° to +10°) -5° -3° 0° 3° 5°

Current Chair: New Chair:

702XG110 **Inward Taper** Measured from inside backrest tube to inside of opposite backrest tube at the top of the backrest. \$150

Current Chair: New Chair:

702XG100 **Ergonomic Backrest** Measured from top of seat tube at rear of frame to the start of the rearward bend of back-post. \$150

Current Chair: New Chair:

702XG120 **Outward Flare** Measured from inside backrest tube to inside of opposite backrest tube at the top of the backrest. \$150

Current Chair: New Chair:

702XY390 **Camber Angle** Place a straight edge across the outside of rear wheel and use an angle finder to show the degree of camber.






Current Chair:	0°	New Chair:	702XY396 0°
	1°		702XY397 1°
	2°		702XY398 2°
	3°		702XY399 3°
	5°		702XY400 5°

702XH165 **Center of Gravity** Measured from front of back post at the seat tube to the centerline of the rear axle.

Current Chair: New Chair:

702XY460 **Rear Wheel Gap** Measured from outer most point of chair (seat back post, side guard or armrest) to inside of rear tire.

Current Chair:	.6"(15)	New Chair:	702XY464 .6"(15)
	.78"(20)		702XY465 .78"(20)
	1"(25)		702XY466 1"(25)



Footrest Design

702XY100	Round Design	STD
702XY105	Square Design	NCO

Footrest Styles

702XB080	Tubular footrest welded into frame aluminum plate (No height adjustment)	STD
702XB110	Tubular footrest welded into frame with carbon fiber footplate (No height adjustment)	\$336
702XB050	Lightweight 3k Carbon Fiber Footrest cover	STD
702XB155	Lightweight Spinergy CX6 Cast Carbon Fiber Footrest cover	NCO
702XB120	Clamp-on aluminum footrest plate (Angle and height adjustment)	NCO
702XB165	Rubber grip tape (3M Safety walk. Not available with 702XB030)	\$89
702XB170	Omit Footplate	NCO

Front Wheels

702XC060	Aluminum Center 4"x1.3" (Black hub only)	STD
702XC110	Aluminum Center 4"x1" (Silver hub only)	\$239
702XC135	Aluminum Center 4"x1.25"	\$320
702XC1131	Sliver hub & black tire	
702XC1132	Black hub & black tire	
702XC1133	Purple hub & black tire	
702XC1134	Green hub & black tire	
702XC1135	Red hub & black tire	
702XC1136	Blue hub & black tire	
702XC155	4" x 1.25" Spinergy CX6 Cast Carbon fiber front wheel	\$125
702XC156	Cerakoted Spinergy CX6 Cast Carbon fiber front wheel (Hub only. See Cerakot color options on page 6)	\$242
702XC070	Soft roll 5"x1.3"	NCO
702XC100	Pneumatic Aluminum Center 5"x1"	\$255
702XC120	Aluminum Center 5"x1"	\$239
702XC140	Aluminum Center 5"x1.4" (Sliver hub & black tire)	\$239
702XC145	Aluminum Center 6"x1.4" (Sliver hub & black tire)	\$239

Forks

702XC190	Low Ride Fork	STD
702XC192	High Trail fork	NCO
702XC300	Omobic Single Sided fork	\$248
702XC210	Froglegs suspension fork	\$652
702XC211	Black fork	STD

Wheel Lock

702XD020	Underseat Scissor	NCO
702XD030	Compact Scissor	\$274
702XD035	Compact Scissor with Direct Mount	\$363
702XD050	Push To Lock	NCO
702XD060	Omobic all aluminum	\$155
702XD095	Omobic all aluminum with Direct Mount	\$185
702XD100	Push Lox Dynamic Hub Locking system	\$1,518

Design Template

702XD085	Eagle (Indicate color)				\$482
702XD086	Eagle with Direct Mount (Indicate color)				\$523
	702XD088 Black	702XD089 Purple	702XD090 Red	702XD091 Blue	
702XD010	Omit Wheel Lock				NCO
Side Guards					
702XE080	Lift out aluminum				\$548
702XE130	Fold flat aluminum				\$690
702XE100	Lift Out 3k Carbon Fiber				\$789
702XE106	Lift Out Spinerger CX6 Cast Carbon fiber				\$1,045
702XE150	Fold flat 3k carbon fiber				\$1,084
702XE155	Fold flat Spinerger CX6 Cast Carbon fiber				\$1,084
702XE120	Lift out ABS plastic				\$502
702XE170	Fold flat ABS plastic				\$855
	702XEM05 25mm above tire (Follow tire)				NCO
	702XEM06 Individual above tire (Follow tire - max. 25mm)				NCO
	702XEM02 Level with tire				STD
	702XEM07 Side guard length				STD
	702XEM08 Side guard height				STD
702XE090	Lift out aluminum with fender				\$596
702XE140	Fold flat aluminum with fender				\$728
702XE110	Lift Out 3k Carbon Fiber with fender				\$1,166
702XE211	Lift Out Spinerger CX6 Cast Carbon fiber with fender				\$901
	702XEM03 10mm above tire				NCO
	702XEM01 Individual above tire (Follow tire - max. 25mm)				NCO
	702XEM09 Center of tire				STD
	702XEM10 Outside of tire				NCO
702XE195	Omit side guards				NCO
Armrest					
702XF010	Lift out armrest with plastic sideguard	Height (Top of seatrail to top of armrest)			\$610
702XF020	Lift out armrest with aluminum sideguard	Height (Top of seatrail to top of armrest)			\$481
702XF130	Swing way tubular armrest with welded in receiver (Requires 702XH110 Square back adjustable center gravity)				\$715
	Height (Top of seatrail to top of armrest)	Tubular armrest length			
702XF140	Omit armrest				NCO
Backrest Height					
702XY230	Fixed height				STD
702XY245	Adjustable height (Maximum height 15.75(400) Range: +/- 25mm)				\$255
Backrest Design					
702XG060	Folding backrest with Q Lock				STD
Push Handles					
702XG070	Screw in push handles (Not available with 702XY245)				\$232
702XG080	Flip down push handles (Not available with 702XY245)				\$454
702XG085	Flip down push handles height adjustable (Not available with 702XY245)				\$742

Design Template

702XG097	Integrated push handles								\$696
702XG030	Omit push handles								NCO
Center of Gravity									
702XH110	Square back adjustable of gravity (Range: 25mm - 150mm)								STD
702XH120	Round back adjustable of gravity (Range: 60mm - 150mm)								NCO
702XH090	V-Style fixed center of gravity (Range: 25mm - 125mm)								\$749
Axle Pin Options									
702XK227	Standard axle pin								STD
702XK226	Titanium 100mm quick release axle pin								\$625
702XK220	Titanium 130mm quick release axle pin								\$165
702XK240	Quad release axle pin								\$248
Anti-Tip									
702XJ010	Swing-away anti tip LEFT side (Only available with square back adjustable center of gravity. Not available with 702XH140)								\$172
702XJ020	Swing-away anti tip RIGHT side (Only available with square back adjustable center of gravity. Not available with 702XH140)								\$172
702XJ100	Single sport anti-tip (Centered mounted)								\$519
702XJ110	Omit Anti-Tip								NCO
Frame and Footrest Finishes									
Non painted finishes									
702XN030	Bead blast frame								STD
702XN050	Polished frame								\$1,372
Painted finishes									
702XN500	CeraKote frame								\$550
702XN600	CeraKote components (Component color will match frame if not otherwise specified. Components include: Forks, Wheel Lock clamps, and Camber Tube clamps.)								\$550
702XN800	CeraKote Spinergy wheels (See Rear Wheel Type on page 8. Hubs and Rims are CeraKoted on Rear Wheels, Hubs on Front Wheels.)								\$675
Frame	Components	Spinergy	Colors	Frame	Components	Spinergy	Colors		
702XN501	702XN601	702XN801	Midnight Black	702XN509	702XN609	702XN809	Rose Gold		
702XN502	702XN602	702XN802	Java Bronze	702XN510	702XN610	702XN810	Retro Blue		
702XN503	702XN603	702XN803	Pure Titanium	702XN511	702XN611	702XN811	Hot Pink		
702XN504	702XN604	702XN804	Fire Red	702XN512	702XN612	702XN812	Wicked Green		
702XN505	702XN605	702XN805	Prime Blue	702XN513	702XN613	702XN813	Bright White		
702XN506	702XN606	702XN806	Sahara Gold	702XN514	702XN614	702XN814	Tangerine Orange		
702XN507	702XN607	702XN807	Black Cherry	702XN515	702XN615	702XN815	Sunbeam Yellow		
702XN508	702XN608	702XN808	Fresh Green	702XN516	702XN616	702XN816	Matte Navy		
Accent Color									
702XN251	Black								STD
702XN253	Red								NCO
702XN255	Purple								NCO
702XN252	Silver								NCO
702XN254	Blue								NCO
702XN250	Individual Accent Color								\$328

Seat Upholstery

702XP040	Seat Sling Adjustable Strap	STD
702XR050	Carbon Fiber Seat Plate	\$494

Back Upholstery

702XP050	Airtech Adjustable Strap	NCO
702XP060	Backrest upgrade to black full leather	\$515
702XP070	Backrest upgrade to black leather corners	\$125
702XP110	Backrest zip pocket	\$131

RGK Backrest Logo Color

702XP557	White (Standard if nothing selected)	STD
702XP558	Red	STD
702XP559	Black	NCO
702XP560	Blue	NCO
702XP561	Purple	NCO
702XP562	Individual Pantone color number	NCO
702XP570	Omit back upholstery	NCO

Rear Wheel Size

702XY420	24" (540) Rear wheel size	702XY430	25" (559) Rear wheel size	STD
----------	---------------------------	----------	---------------------------	-----

Rear Wheel Type

702XK160	Spinergy ES Wire (20 spokes. CeraKote option available on Rim. Not available with Spinergy anodized color options.)	STD
702XK015	Proton (Not available with Spinergy anodized color options.)	STD
702XK320	Spinergy XLX (12 spoke, cross radial spokes.)	\$1,099
702XK050	Spinergy LX (12 spoke)	\$1,030
702XK060	Spinergy SPOX (18 spoke. Not available with Spinergy anodized color options or spoke colors, purple, orange, pink and green.)	\$685
702XK070	Spinergy BLXL (18 bladed spoke.)	\$1,099
702XK150	Spinergy CLX (18 bladed spokes, carbon fiber rim. Spinergy color option available on hub only. Not available with silver hub)	\$3,233
702XK330	RGK EXCLUSIVE Spinergy CX6 Cast Carbon Fiber (Cast Carbon Fiber Hub, Matte Black Carbon Fiber Rim and Black Bladed Spokes with 2 White Bladed Spokes Around Valve Stem. Not available with Spinergy anodized color options.)	\$1,915
	Spinergy spoke colors	
702XK650	Black	STD
702XK651	Red	NCO
702XK652	Purple	NCO
702XK653	White	NCO
702XK654	Blue	NCO
702XK655	Orange	NCO
702XK656	Yellow	NCO
702XK657	Pink	NCO
702XK658	Green	NCO
702XK180	Individual spoke design	\$219
	Spinergy anodized colors (Rims and Hubs are anodized not spokes, on Spinergy BLXL (XK075), Hubs only on 702XK150)	
702XK190	Black	STD
702XK202	Purple (LX, XLX & BLXL only)	\$421

Design Template

702XK203	Blue (LX, XLX & BLXL only)				\$421
702XK201	Red (LX, XLX & BLXL only)				\$421
702XK204	Gold (LX, XLX & BLXL only)				\$555
702XK205	Individual hub/rims				\$585
Supply Without Rear Wheels					
702XK040	Omit Rear Wheels				NCO
	702XK041 Wheel Type:		702XK042 Tire Type:		
Handrims Style					
702XY485	Tab Length Wide				STD
702XY480	Tab Length Narrow				NCO
	702XL090	Remove excess mounting tabs (Recommended for small wheel gaps XY460)			\$150
702XL010	Aluminum Hard Anodized 6 Tab				STD
702XL160	Elipse Ergonomic 6 Tab				\$650
702XL020	QGrip 6 Tab				\$412
702XL040	Titanium 6 Tab				\$459
702XL035	Natural-Fit LT	702XL036 Standard thumb ring	702XL037 Super Grip thumb ring	702XL038 Omit thumb ring	\$765
702XL150	Carbolife Gekko				\$819
702XL220	Carbolife Gekko S				\$654
702XL170	Carbolife Curve L				\$716
702XL230	Carbolife Curve				\$569
702XL180	Carbolife Curve L Tetra				\$912
702XL240	Carbolife Curve Tetra				\$727
702XL190	Carbolife Quadro				\$912
702XL200	Spinergy TC				\$758
702XL210	Spinergy TC Max				\$890
Tires					
702XM060	Schwalbe One (Lightest option)				\$137
702XM150	IRC EXERACER				\$194
702XM010	Kenda	702XM011 Grey	702XM012 Blue	702XM013 Red	STD
702XM020	Schwalbe Marathon Plus				\$101
702XM040	Shox Solid Tire				\$165
702XM070	Schwalbe Rocket Ron mountain bike tire				\$69
702XM050	Kenda K885 Mountain Bike Tire				\$41
Inner Tubes					
702XY500	Schraeder Valve Inner Tube				STD
Cushions					
702XP080	Cushion foam				\$125
	Density	702XP001 Soft	702XP002 Medium		
	Thickness	702XP003 50mm	702XP004 75mm	702XP005 100mm	
Cushion Covers					
702XP085	Nylon Cushion Cover (Black)				STD

Design Template

702XP090	Airlite Cushion Cover (Black)				\$101
Jay Products - Please Include Jay Order Form					
702XP900	Jay Cushion	Refer to Jay Order Form			
702XP910	Jay Backrest	Refer to Jay Order Form			
Empulse Products - Please Include Empulse Order Form					
702EP100	⊘ Empulse Add-On Power-Assist Devices (IN11- USEWITHRGK) F35 - Pull Device M90 - Wheel Add-On R90 - Push-Assist	Ship with chair	or	Ship seperate from chair (uninstalled)	
Accessories					
702XR080	Frontwheel Assembly (Only with 702XB010, 702XB080, 702XB150)				\$1,085
702XR240	Reverse Velcro on seating	702XR241 Double sided tape	702XR242 Sewn in		\$138
702XP120	Black Frame Protectors Neoprene Zip Fasten (Black)				\$87
702XP420	Black Frame Protectors Nylon Velcro Fasten (Black)				\$83
702XR010	Black Seat Pouch Velcro Attach (Black)				\$61
702XR200	Black Rucksack with RGK Logo (Black)				\$156
702XR100	FX carry bag				\$344
702XR230	Spinergy Wheel Bag				\$207
702XR210	Recon Mini Rider-Carry Air Compressor				\$235
702XR090	Tool Kit				STD
702XS060	Calf Strap				STD

Ordering Information

Estimated MSRP*
(*Price does not include any one off requested options)

Quote

Order

Specification Notes

