



Order Form

EIPW32

Account:	
Date:	
PO#:	
Buyer Name:	
Marked For:	

Ship To:	
Name:	
Address:	
City, Province, Postal Code:	
Phone:	


Sunrise Medical Canada Inc.
 237 Romina Drive, Unit 3
 Concord ON L4K 4V3
 1 (800) 263-3390



Customer Service: (800) 263-3390 Fax: (800) 561-5834 Email: cscanada@sunmed.com sunrisemedical.ca
 Specifications are subject to change without notice. © 2023 Sunrise Medical Canada Inc. MKc 129728c AADL, REV

Client Measurements

a)	Top of head to seat	Inch
b)	Top of shoulder to seat	Inch
c)	Underarm to seat	Inch
d)	Elbow to seat	Inch
e)	Top of leg to seat	Inch
f)	Shoulder width	Inch
g)	Chest width	Inch
i)	Hip width	Inch
j)	Actual sitting depth	Inch
m)	Foot length	Inch
(n)	Knee to foot	Inch



User Weight

Notes:

Chair Specifications

Power Seating	50° of CG Tilt / 12" powered seat elevator	Seat Size Range	Width 16"-22" Depth 14"-22"
Seat Type	SEDEO LITE	STF Height Range	16.5"-19"

Hybrid Drive Base Specifications

Drive Wheel Position	Hybrid Drive	Drive Wheels	14" drive wheels w/flat free inserts 14" w/pneumatic inserts
User Weight Capacity	300 lbs., 350 lbs HD	Tires/Casters	9" Front / 4" Rear (Black)
Minimum Turning Radius	25.5" (excludes footplates)	Charger	10A Dual Mode (Compact)
Base Size	Width 24.5", Length 36"	Battery Size	Group 24 / Group 34
Drive Control	VR-2 90, R-net 90A,	Color Accents	7 Colour Accents
Drive Motors	4-pole 6mph, & 5mph	Ground Clearance	3"

Base

Standard Duty Hybrid Drive Base = End Item Number EIPW32 (AADL Contract Price)

<input type="checkbox"/>	320PC248	Quickie Q500 H NPC No Power Option for SEDEO LITE Seat	\$8,995.00
<i>Included in Base Pricing : Q500 H Base, VR-2 90amp Electronics for No Power Seat Functions, Six Wheel Independent Suspension System, 14" Drive Wheels with Dual Density Inserts, 9" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5".</i>			

<input type="checkbox"/>	320PC249	Quickie Q500 H SPC Single Power Option for Sedeo LITE Seat	\$8,995.00
<i>Included in Base Pricing : Q500 H Base, VR-2 90amp Electronics for No Power Seat Functions, Six Wheel Independent Suspension System, 14" Drive Wheels with Dual Density Inserts, 9" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5".</i>			

Heavy Duty Hybrid Drive Base = End Item Number EIPW32 (AADL Contract Price)

<input type="checkbox"/>	320PC251	Quickie Q500 H NPC HD No Power Option for SEDEO LITE Seat	\$10,000.00
<i>Included in Base Pricing : Q500 H Base, VR-2 90amp Electronics for No Power Seat Functions, Six Wheel Independent Suspension System, 14" Drive Wheels with Dual Density Inserts, 9" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 5mph, Weight Capacity = 350lbs., STFH = 16.5".</i>			

<input type="checkbox"/>	320PC252	Quickie Q500 H SPC HD Single Power Option for SEDEO LITE Seat	\$10,000.00
<i>Included in Base Pricing : Q500 H Base, VR-2 90amp Electronics for No Power Seat Functions, Six Wheel Independent Suspension System, 14" Drive Wheels with Dual Density Inserts, 9" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 5mph, Weight Capacity = 350lbs., STFH = 16.5".</i>			

Base

User Weight

<input type="checkbox"/>	320WG23	Up to 300 lbs.	NC
<input type="checkbox"/>	320WG15	301 - 350 lbs.	NC

Motor Package

<input type="checkbox"/>	320PF98	4-pole 8 kph / 5 mph <i>(Included with HD base)</i>	NC
<input type="checkbox"/>	320PF99	4-pole 10 kph / 6 mph <i>(Included with Standard base)</i>	NC

Tracking

<input type="checkbox"/>	320MI01	SureTrac Module & Motor Encoders <i>(Improves Drive Tracking)</i>	\$1,565.00
--------------------------	---------	---	------------

Drive Wheels

<input type="checkbox"/>	320DW58	Drive Wheel, 14", Dual Density, Knobby, Black	NC
<input type="checkbox"/>	320DW67	Drive Wheel, 14", Pneumatic, Knobby, Black	\$384.00

Caster Wheels

<input type="checkbox"/>	320CW57	Caster Wheels 9", Solid, Black	NC
<input type="checkbox"/>	320CW56	Caster Wheels 9", Pneumatic, Black	NC

Wheel Locks

<input type="checkbox"/>	320WL42	Wheel Locks	\$130.00
--------------------------	---------	-------------	----------

Color Accents

<input type="checkbox"/>	320EZ20	Apple Green	NC
<input type="checkbox"/>	320EZ22	Black	NC
<input type="checkbox"/>	320EZ21	Blue	NC
<input type="checkbox"/>	320EZ23	Orange	NC
<input type="checkbox"/>	320EZ24	Red	NC
<input type="checkbox"/>	320EZ25	Silver	NC
<input type="checkbox"/>	320EZ19	White	NC

ISO Transit Tie Downs

<input type="checkbox"/>	320TO50	ISO Transit Package	NC
--------------------------	---------	---------------------	----

LED Lights **(Client Cost)**

<input type="checkbox"/>	320LI04	LED Lights and Indicators	\$1,131.00
<input type="checkbox"/>	320LI02	Omit LED Lights	NC

Seat Heights

Note: Available Seat-to-Floor heights are dependant on 2 options: 1. Battery size and 2. Powered Seating.

Use the matrix below to determine available ranges for your desired options. Measured from front of seat pan to floor with 0° frame angle.

	Seat Height	
<input type="checkbox"/>	320HT103	16.5" NC
<input type="checkbox"/>	320HT55	17" NC
<input type="checkbox"/>	320HT104	17.5" NC
<input type="checkbox"/>	320HT56	18" NC
<input type="checkbox"/>	320HT105	18.5" NC
<input type="checkbox"/>	320HT57	19" NC

Batteries & Chargers **(AADL Contract Pricing)**

Group 24 Battery Options *(Recommended for max daily range)*

<input type="checkbox"/>	320AG25	GP24 Flat Top Pair w/Strap MK Gel	\$844.00
<input type="checkbox"/>	320AG26	GP24 Flat Top Pair w/Strap Brightway AGM	\$881.00

Group 34 Battery Options *(Required for lower seat-to-floor heights)*

<input type="checkbox"/>	320AG27	GP34 Flat Top Pair w/Straps MK Gel	\$918.00
<input type="checkbox"/>	320AG29	GP34 Flat Top Pair w/Straps Brightway AGM	\$796.00

Omit Batteries

<input type="checkbox"/>	320AG24	Battery Harnesses Only <i>(Omit Batteries)</i>	NC
--------------------------	---------	--	----

Chargers

<input type="checkbox"/>	320BJ30	10A Dual Mode <i>(Compact)</i>	\$150.00
<input type="checkbox"/>	320BJ19	8A Dual Mode <i>(Compact)</i>	NC

Powered Seating

Power Seating Options

Seat Elevator (Client Cost)

<input type="checkbox"/>	320PI07	Seat Elevator, Powered, 300 mm / 12" (Not available with HD model)	\$3,209.00
<input type="checkbox"/>	320PI08	C-Me Seat Elevator, Powered, 300 mm / 12" (Not available with HD model) C-Me speed up to 3.0 mph.	\$3,209.00

Seat Tilt (AADL W993 for Power Tilt)

<input type="checkbox"/>	320PA35	Posterior 50° Tilt	\$7,995.00
--------------------------	---------	--------------------	------------

Pre Set Seat Frame Angle (No Power Seating Options, Increases Seat to Floor Height)

<input type="checkbox"/>	320FA10	0° (Adds 0" to STFH)	NC
<input type="checkbox"/>	320FA3	3° (Adds 1" to STFH)	NC
<input type="checkbox"/>	320FA4	6° (Adds 2" to STFH)	NC
<input type="checkbox"/>	320FA5	9° (Adds 3" to STFH)	NC

Backrest Frame Type

<input type="checkbox"/>	320BF10	Tool Angle Adjustable Back (SEDEO LITE: 86°-106°)	NC
--------------------------	---------	---	----

Backrest Frame Option

<input type="checkbox"/>	320BF13	Flip Down, Quick Release	\$165.00
<input type="checkbox"/>	320BF24	Removable Back Canes	NC

SEDEO Lite Backrest Frame

Back Canes Height

<input type="checkbox"/>		Back Canes, Low (21" from seat pan to top of push handle, 19" upholstery height)	NC
<input type="checkbox"/>		Back Canes, Med (23" from seat pan to top of push handle, 21" upholstery height)	NC

Back Cushion Type

<input type="checkbox"/>		Standard Back Upholstery (19" upholstery height low back, 21" upholstery height medium back)	NC
<input type="checkbox"/>	320BT133	Omit (For Custom Mounting. No Headrest Mount)	NC

Seat Size

Seat Pan Width

<input type="checkbox"/>	320R78	15"	NC
<input type="checkbox"/>	320R53	16"	NC
<input type="checkbox"/>	320R54	17"	NC
<input type="checkbox"/>	320R55	18"	NC
<input type="checkbox"/>	320R56	19"	NC
<input type="checkbox"/>	320R57	20"	NC

Seat Pan Depth

<input type="checkbox"/>	320RF52	14"	NC
<input type="checkbox"/>	320RF53	15"	NC
<input type="checkbox"/>	320RF51	16"	NC
<input type="checkbox"/>	320RF45	17"	NC
<input type="checkbox"/>	320RF46	18"	NC
<input type="checkbox"/>	320RF47	19"	NC
<input type="checkbox"/>	320RF48	20"	NC
<input type="checkbox"/>	320RF49	21"	NC
<input type="checkbox"/>	320RF50	22"	NC

Positioning Supports & Accessories

Knee and Thigh Swing Away Support Pads (Seating Benefits)

Lateral Knee Support Pads and Hardware (Must Pick Pad and Mount)

<input type="checkbox"/>	320XY05	Lateral Knee Support Pad, Right (4" round)	\$228.00
<input type="checkbox"/>	320XY06	Lateral Knee Support Pad, Left (4" round)	\$228.00
<input type="checkbox"/>	320XX08	Lateral Swing Away Bracket Knee - Right	\$300.00
<input type="checkbox"/>	320XX09	Lateral Swing Away Bracket Knee - Left	\$300.00

Lateral Thigh Support Pads and Hardware (Must pick pad and mount)

<input type="checkbox"/>	320KF22	Lateral Thigh Support Pad Short, Right (4.8" H x 4.8" L / 120mm x 120mm)	\$134.00
<input type="checkbox"/>	320KF23	Lateral Thigh Support Pad Short, Left (4.8" H x 4.8" L / 120mm x 120mm)	\$134.00
<input type="checkbox"/>	320KF24	Lateral Thigh Support Pad Long, Right (4.8" H x 7" L / 120 mm x 180mm)	\$134.00
<input type="checkbox"/>	320KF25	Lateral Thigh Support Pad Long, Left (4.8" H x 7" L / 120mm x 180mm)	\$134.00
<input type="checkbox"/>	320KF26	Lateral Thigh Removable Mount - Right	\$300.00
<input type="checkbox"/>	320KF27	Lateral Thigh Removable Mount - Left	\$300.00

Armrest Structure

<input type="checkbox"/>	320EG35	Reclining, Flip Back, Long (Height Adjustable 8" to 13")	NC
<input type="checkbox"/>	320EG36	Reclining, Flip Back with 2-point post locking arm (Height Adjustable 11" to 13" x .625")	\$325.00
<input type="checkbox"/>	320EG34	Single Post, SEDEO LITE, Removable (Not available with thigh support)(Height Adjustable 10" to 15")	NC

Arm Pads

SEDEO Armrest Pad

<input type="checkbox"/>	320AB43	300 mm / 12", Integral Foam, Right	NC
<input type="checkbox"/>	320AE43	300 mm / 12", Integral Foam, Left	NC
<input type="checkbox"/>	320AB44	400 mm / 16", Integral Foam, Right	NC
<input type="checkbox"/>	320AE44	400 mm / 16", Integral Foam, Left	NC
<input type="checkbox"/>	320AB61	SEDEO LITE, Right	NC
<input type="checkbox"/>	320AE61	SEDEO LITE, Left	NC

Comfort Arm Pad (Foam) (Seating Benefits)

<input type="checkbox"/>	320AB46	2" Wide Right	\$104.00
<input type="checkbox"/>	320AE46	2" Wide Left	\$104.00
<input type="checkbox"/>	320AB47	4" Wide Right	\$104.00
<input type="checkbox"/>	320AE47	4" Wide Left	\$104.00

Comfort Arm Pad (Gel) (Seating Benefits)

<input type="checkbox"/>	320AB48	2" Wide Right	\$197.00
<input type="checkbox"/>	320AE48	2" Wide Left	\$197.00
<input type="checkbox"/>	320AB49	4" Wide Right	\$197.00
<input type="checkbox"/>	320AE49	4" Wide Left	\$197.00

AES Trough/Elbow Block Arm pads (Gel) (Seating Benefits)

<input type="checkbox"/>	320AB50	3.5" x 13" Right	\$274.00
<input type="checkbox"/>	320AE50	3.5" x 13" Left	\$274.00
<input type="checkbox"/>	320AB51	4" x 14.5" Right	\$274.00
<input type="checkbox"/>	320AE51	4" x 14.5" Left	\$274.00

Positioning Supports & Accessories

Side Guards (Pair)

<input type="checkbox"/>	320SG46	Side Guards Mounted On Seat Rail Bracket (Only with Reclining Armrests) (N/A w/ Thigh Supports)	\$328.00
--------------------------	---------	---	----------

SEDEO PRO Swing-Away Legrest (Pair) (320H179 is Client Cost)

<input type="checkbox"/>	320H177	Legrests, Swing Away, 70°, SEDEO PRO (14" to 21" Leg Length)	NC
<input type="checkbox"/>	320H178	Legrests, Manual Elevating, Swing Away, SEDEO PRO (Includes Padded Calf Pads) (14".5 to 21.5" Leg Length)	\$558.00
<input type="checkbox"/>	320H179	Legrests, Power Elevating, Swing Away, SEDEO PRO (Includes Padded Calf Pads) (15.5" to 22.5" Leg Length)	\$3,169.00
<input type="checkbox"/>	320H180	Hanger Receivers Only, SEDEO PRO	NC
<input type="checkbox"/>	320H176	Omit Legrests	NC

Swing Away Legrest Footplates

<input type="checkbox"/>	320F369	Angle Adjustable Flip-Up Footplates	\$262.00
--------------------------	---------	-------------------------------------	----------

Swing Away Legrest Calf Supports

<input type="checkbox"/>	320DI41	Calf Support, Depth Adjustable, Integral Foam (7" Wide x 10" Tall, N/A 16-17" W)	\$384.00
<input type="checkbox"/>	320DI42	Calf Support, Depth Adjustable, Padded (6.5" to 8" Wide, Depending On Chair Width, x 9.5" Tall)	\$472.00

Swing Away Legrest Options

<input type="checkbox"/>	320F364	Footplates, Coupled	\$262.00
<input type="checkbox"/>	320DI44	Heel Loops	NC
<input type="checkbox"/>	320DI40	Calf Strap (Only For 70° Swing-Away Legrests)	\$130.00
<input type="checkbox"/>	320KG4	Knee Pads, Integral Foam (2.5" Wide x 4.5" Tall)	\$260.00

SEDEO Lower Extremity Swing Away Support

<input type="checkbox"/>	320VD01	Lower Extremity Support, Padded, Swing Away, SEDEO, Right	\$472.00
<input type="checkbox"/>	320VD02	Lower Extremity Support, Padded, Swing Away, SEDEO, Left	\$472.00

Pelvic Positioning Belts

<input type="checkbox"/>	320A837	Positioning Belt, Auto Style Buckle, 50 mm 1600 mm / 2"	NC
<input type="checkbox"/>	320A839	Positioning Belt, Aircraft Buckle, 50 mm x 1300 mm / 2"	NC
<input type="checkbox"/>	320A840	Positioning Belt, Aircraft Buckle, Padded 50 mm x 1300 mm / 2"	\$95.00

Accessories

Ventilator Option

<input type="checkbox"/>	320VT12	Vent Tray Fixed (Pulmonetic LTV Models)	\$562.00
--------------------------	---------	---	----------

O2 Tank Holder

<input type="checkbox"/>	320A718	O2 Tank Holder Mounted to Backrest (D and E Styles)	\$361.00
--------------------------	---------	---	----------

Wheelchair Tray Table

<input type="checkbox"/>	3207A09	Tray, Swing Away	\$351.00
--------------------------	---------	------------------	----------

Backpack (Client Cost)

<input type="checkbox"/>	320A1	Black Utility	\$103.00
--------------------------	-------	---------------	----------

Medial Thigh Support (Seating Benefits)

<input type="checkbox"/>	3208A80	Medial Thigh Support	\$224.00
<input type="checkbox"/>	3208A81	Pommel Swing Away Hardware	\$330.00

Misc. Accessories (Client Cost)

<input type="checkbox"/>	3208A59	Transfer Bar (Pair)	\$455.00
<input type="checkbox"/>	3208A52	Ram Mount Phone	\$71.00
<input type="checkbox"/>	3208A53	Ram Mount Tablet	\$284.00
<input type="checkbox"/>	3208A54	Ram Mount Cup Holder	\$71.00
<input type="checkbox"/>	3208A58	Rear View Mirror	\$71.00
<input type="checkbox"/>	320A611	Crutch Holder	\$125.00
<input type="checkbox"/>	3208A73	Bag For Armrest	\$71.00
<input type="checkbox"/>	3208A75	Reflectors (Backrest shroud mounted)	\$71.00

Joystick / Electronics

Joystick (Remote)

<input type="checkbox"/>	320J127	VR-2 Joystick	NC
<input type="checkbox"/>	320J121	Std LED Joystick for R-net Controls	NC
<input type="checkbox"/>	320J116	Color LCD (<i>CJSM1</i>)	NC
<input type="checkbox"/>	320J117	Color LCD, Heavy Duty (<i>Fixed Joystick Mount Only</i>)	\$1,495.00
<input type="checkbox"/>	320J118	High Resolution Color LCD, (CJSM2+BT) (Client Cost) 3.5" Color Display, Paddle Switches, Programmable Shortcut Keys, ICONs, Stereo Jack, IR and dual (IOS/Android) Bluetooth control	\$1,167.00
<input type="checkbox"/>	320J119	Omit Remote Control (<i>Required with Arm Mounted Specialty Controls</i>)	-\$25.00

Joystick Mounting

<input type="checkbox"/>	320JM34	Fixed on Armrest	NC
<input type="checkbox"/>	320JM35	Parallel Swing Away	\$412.00
<input type="checkbox"/>	320JM36	Link-It Mid Line Mount on Armrest (<i>Swing Away</i>)	\$412.00

Joystick Mount (Location)

<input type="checkbox"/>	320HM16	Right Hand armrest	NC
<input type="checkbox"/>	320HM17	Left Hand armrest	NC

Display Language Selection (*R-net*)

<input type="checkbox"/>	320L1	English	NC
<input type="checkbox"/>	320L3	French	NC
<input type="checkbox"/>	320L7	Spanish	NC

Joystick Handle (*Optional, N/A with HD Joystick*)

<input type="checkbox"/>	320JH13	Ball Handle	\$101.00
<input type="checkbox"/>	320JH14	Foam Ball Handle	\$101.00
<input type="checkbox"/>	320JH15	T-Handle	\$101.00
<input type="checkbox"/>	320JH16	Goalpost	\$101.00
<input type="checkbox"/>	320JH24	Bodypoint - Rubber Dome Joystick Handle	\$200.00
<input type="checkbox"/>	320JH25	Bodypoint - U-Shaped Joystick Handle (<i>4 in</i>)	\$200.00
<input type="checkbox"/>	320JH26	Bodypoint - Grooved Mushroom Joystick	\$200.00

Electronic Programmers

<input type="checkbox"/>	320PR10	R-net PC Programmer with Dongle	\$652.00
<input type="checkbox"/>	320PR11	VR-2 Hand Held Programmer	\$449.00

Power Modules

Non-Expandable

<input type="checkbox"/>	320ER12	Power Module, VR-2, 90A, Non-Expandable (<i>Only available with 6mph motors</i>)	NC
<input type="checkbox"/>	320ER18	Power Module, R-net, 90A, Non-Expandable (<i>Only available with 5mph and 8mph</i>)	NC

Expandable

<input type="checkbox"/>	320ER13	Power Module, R-net, 90A, Expandable (<i>Required for Specialty Controls or 2+ Actuators</i>)	NC
	320EX17	Expandable 90A Controller - R-net	NC
	320EX18	Harness for 90A Expandable Controller	NC

Joystick / Electronics

Seating Control

Powered Seating Control

<input type="checkbox"/>	320EO40	Omit Thru-drive <i>(Requires Switch Control)</i>	NC
<input type="checkbox"/>	320EO37	Thru-drive Control for 1 Actuator	\$1,332.00
<input type="checkbox"/>	320EO38	Thru-drive Control for 2+ Actuators	\$2,679.00

Switch Control *(Available with Omit Thru-Drive Only) (1 actuator for Tilt or Elevate Only)*

<input type="checkbox"/>	320EO33	Dual Toggle for 1 actuator	NC
<input type="checkbox"/>	320EO34	Buddy Button for 1 actuator <i>(Button Switches are not mounted)</i>	NC
<input type="checkbox"/>	320EO35	Switch-It Mini 1 button for 1 actuator	NC
<input type="checkbox"/>	320EO36	Switch-It Mini 2 button for 1 actuator	NC
<input type="checkbox"/>	320EO41	Ctrl+5 Switch Box for 2 actuators or more	NC

Additional Power Seating Control *(Available with Thru-Drive Only)*

<input type="checkbox"/>	320QG11	Ctrl+5 Switch Box With Free Assignable Functions	\$995.00
<input type="checkbox"/>	320QG12	Ctrl+5 Switch Box With Free Assignable Functions <i>(Multi-Power only)</i>	NC

Dual Toggle Mounting Location

<input type="checkbox"/>	320TN1	Toggle Mount, Right	NC
<input type="checkbox"/>	320TN2	Toggle Mount, Left	NC

Ctrl+5 Mounting Location

<input type="checkbox"/>	320EY19	Right Armrest	NC
<input type="checkbox"/>	320EY20	Left Armrest	NC
<input type="checkbox"/>	320EY21	Back Cane, Right	NC
<input type="checkbox"/>	320EY24	Back Cane, Left	NC

Ctrl+5 Switch Type

<input type="checkbox"/>	320EY22	Push Buttons	NC
<input type="checkbox"/>	320EY23	Toggle Switches	\$328.00

Attendant Controls *(Prior Approval Required)*

<input type="checkbox"/>	3200A05	Dual Control Module for Attendant	\$1,131.00
<input type="checkbox"/>	3200A06	R-net Color LCD <i>(CJSM1)</i> , Attendant Control Only	\$1,131.00

Attendant Mounting Location

<input type="checkbox"/>	320DN5	Right Back Cane	NC
<input type="checkbox"/>	320DN6	Left Back Cane	NC

Assignable Buttons

Assignable Button feature allows for customization of joystick buttons and/or Ctrl+5 Switches allowing access to various wheelchair and seating functions. When enabled, chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle.

<input type="checkbox"/>	320PM12	Include Assignable buttons	NC
--------------------------	---------	----------------------------	----