

**QUICKIE**<sup>®</sup>

**QRi**<sup>™</sup>

A hybrid of strength and innovation



## Live Without Limits

The Quickie® QRi™ is an ultra lightweight rigid chair with a strong and innovative frame design perfect for users on the go.

# QUICKIE® QRi



### Frame Design

- Ultra lightweight and durable aluminum frame with a 265 lb. weight capacity
- Open frame and fluid caster housing design for a stylish look

### Commuter Options

- WC-19 Transit Approved at 265 lbs.
- Small footprint for maneuverability and easy transport

### Individualized Fit

- Standard 1" front frame inset
- .5"-3" center of gravity axle plate adjustment
- 75°, 80°, and 85° front frame angle options



Transit option



Standard 1" frame inset



Stable, adjustable caster design



Broad selection of wheels and tires

SPECIFICATIONS	QRi
HCPCS Code	K0005
Frame Style	Rigid - Open
Overall Weight	22 lbs.
Weight Capacity	265 lbs.
Seat Widths	12" to 20"
Seat Depths	12" to 20"
Back Angle Adjustment	-8° to 17°
Front STFH	15" to 20"
Rear STFH	14" to 20"
Frame Angles	75°, 80°, 85°

