



ORDER FORM

February 2024

Freestander

Upright accessible standing frame
from Early Intervention to adult.

Customer			
Account			
Account #		Date	
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #	
Buyer			
Mark For			

Shipping Address			
Name			
Address			
City		Prov	
PC		Phone	
Email / Fax			

Freestander

Size 1 Weight Capacity: 132 lbs | Size 2 Weight Capacity: 176 lbs | Size 3 Weight Capacity: 198 lbs

✓ OPTION	DESCRIPTION	RETAIL
----------	-------------	--------

Freestander

\$

Selection Required. Includes Chest and Hip Pad, Flexible Lateral, PU Knee Supports, Footplate, Tray Cover and Activity Bowl. Sandal Size and positioning belt selections required

<input type="checkbox"/>	LFS/1-MW	Size 1 Freestander - Meadow	\$2,604
<input type="checkbox"/>	LFS/2-MB	Size 2 Freestander - Mondeo Blue	\$2,870
<input type="checkbox"/>	LFS/3-MB	Size 3 Freestander - Mondeo Blue	\$3,229

Positioning

\$

If Solid Pelvic Plate selected, pelvic band selection required.

<input type="checkbox"/>	LFS/1/B-MW	Size 1 Hip Belt - Meadow	\$277
<input type="checkbox"/>	LFS/2/B-MB	Size 2 Hip Belt - Mondeo Blue	\$287
<input type="checkbox"/>	LFS/3/B-MB	Size 3 Hip Belt - Mondeo Blue	\$311
<input type="checkbox"/>	LFS/1/L-MW	Size 1 Solid Pelvic Plate - Meadow	\$447
<input type="checkbox"/>	LFS/2/L-MB	Size 2 Solid Pelvic Plate - Mondeo Blue	\$505
<input type="checkbox"/>	LFS/3/L-MB	Size 3 Solid Pelvic Plate - Mondeo Blue	\$565
<input type="checkbox"/>	LPS/1/B-MW	Size 1 Pelvic Band - Meadow	\$261
<input type="checkbox"/>	LPS/2/B-MB	Size 2 Pelvic Band - Mondeo Blue	\$287
<input type="checkbox"/>	LPS/3/B-MB	Size 3 Pelvic Band - Mondeo Blue	\$311

Additional Knee Supports

\$

<input type="checkbox"/>	LPS/1/K-MW	Size 1 Upholstered Knee Cups - Meadow (Pair)	\$179
<input type="checkbox"/>	LPS/2/K-MB	Size 2 Upholstered Knee Cups - Mondeo Blue (Pair)	\$179
<input type="checkbox"/>	LPS/3/K-MB	Size 3 Upholstered Knee Cups - Mondeo Blue (Pair)	\$179
<input type="checkbox"/>	LPS/4/K-MB	Size 4 Upholstered Knee Cups - Mondeo Blue (Pair)	\$179

Foot Positioning

\$

Sandal selection required.

<input type="checkbox"/>	152-1600	Small Sandals with Straps (Pair)	\$269
<input type="checkbox"/>	152-2600	Medium Sandals with Straps (Pair)	\$287
<input type="checkbox"/>	152-3600	Large Sandals with Straps (Pair)	\$310
<input type="checkbox"/>	117-787	Sandal Raiser 1" (Each) Requires Sandal Selection.	\$105

Accessories

\$

<input type="checkbox"/>	LFS-9500	Freestander Care Kit (Includes spare nuts, bolts, screws and washers)	\$534
--------------------------	----------	---	-------

Technical Specifications - Imperial

Frame Sizing Specification	Size 1	Size 2	Size 3
Age (approximate)	1 1/2 - 5 years	4 - 10 years	9 - 18 years
User Body Height (Min / Max)	31.5" / 45.3"	43.3" / 57"	55.1" / 72.8"
Maximum User Weight	132 lbs	176 lbs	198 lbs
Chest Pad Height (Min / Max)	20.7" / 33.5"	26.8" / 43.3"	38.4" / 56.1"
Width Between Chest Supports (Min / Max)	6" / 10"	8" / 12"	10" / 14"
Activity Tray	16" x 21"	18" x 23"	20" x 25"

Technical Specifications - Metric

Frame Sizing Specification	Size 1	Size 2	Size 3
Age (approximate)	1 1/2 - 5 years	4 - 10 years	9 - 18 years
User Body Height (Min / Max)	800 mm / 1151 mm	1100 mm / 1448 mm	1400 mm / 1849 mm
Maximum User Weight	60 kgs	80 kgs	90 kgs
Chest Pad Height (Min / Max)	470 mm / 680 mm	680 mm / 975 mm	975 mm / 1350 mm
Width Between Chest Supports (Min / Max)	152 mm / 254 mm	203 mm / 305 mm	254 mm / 356 mm
Activity Tray	400 mm x 505 mm	450 mm x 575 mm	500 mm x 625 mm



Sunrise Medical Canada Inc.
237 Romina Drive, Unit 3 • Concord, ON L4K 4V3 • Canada

Customer Service: (800) 263-3390 • E-mail: cscanada@sunmed.com • www.sunrisemedical.ca

©2024 Sunrise Medical Canada Inc. MK-129857-CAN Rev (6)

