



Quickie 2 Lite

April 2025

EIQ2X

Mark For:

Date: _____

Dealer Acct #: _____

Dealer: _____

Dealer Contact: _____

Dealer Address: _____

Dealer City: _____ Prov: _____ PC: _____

Dealer Phone: () _____ Fax: () _____

Confirmation Email: _____

Confirm Via: ☐ Fax ☐ Email

Submitting for:

☐ Quote

☐ Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____

Attention: _____

Address: _____

Address: _____

Ship To City: _____ Prov: _____ PC: _____

Ship To Phone: () _____ Fax: () _____

✓	Option #	Definition	Price	✓	Option #	Definition	Price
---	----------	------------	-------	---	----------	------------	-------

1. FRAME

Base model

☐	002WM58	Quickie 2 Lite QLT	\$2,953.00
--------------------------	---------	--------------------	------------

Transit Option

☐	002TO49	Transit System WC19	NC
--------------------------	---------	---------------------	----

Frame Style (265 lbs. Weight Capacity)

See figure 1. Measure from front frame to floor

☐	002FS31	Swing Away Frame	NC
☐	002FS32	Fixed Front 70°	NC
☐	002FS33	Fixed Front 80°	NC

Frame Inset (Fixed Front Only)

See **FIGURE 2** Frame inset is the distance from the outside of the seat tube to the outside of the footrest. Footrest width = frame width minus 2.

☐	002FF1	0"	NC
☐	002FF2	1" N/A for 14" and 15" Wide using composite angle adj footrest	NC

2. SIZE / SEATING

Frame Width

See **FIGURE 2** Measured from outside of seat tube to the outside of opposite seat tube.

☐	002W5	14"	NC
☐	002W6	15"	NC
☐	002W7	16"	NC
☐	002W8	17"	NC
☐	002W9	18"	NC
☐	002W10	19"	NC
☐	002W11	20"	\$308.00

Seat Depth

See **FIGURE 1** and **Seat Squeeze Matrix** Measured from back canes to front edge of upholstery.

☐	002SD7	14" Front of seat sling 2" behind front edge of frame.	NC
☐	002SD8	15"	NC
☐	002SD9	16"	NC
☐	002SD10	17"	NC
☐	002SD11	18"	NC
☐	002SD72	19"	\$308.00
☐	002SD13	20"	\$308.00

Frame Depth

If growth is desired Seat Depth (*) may be ordered up to 2" smaller on swing away frames and up to 3" on fixed front frames.

☐	002FD5	14"	NC
☐	002FD6	15"	NC
☐	002FD7	16" (*)	NC
☐	002FD8	17" (*)	NC
☐	002FD9	18" (*)	NC
☐	002FD10	19" (*)	NC
☐	002FD36	20" (*)	NC

Seat Sling

☐	002SL1	Standard Seat Sling	NC
☐	002SL6	Omit Seat Sling	NC



JAY Cushions (Optional)

To order any JAY cushions on-chair, please see the JAY price list.

Cushion name or part number must be filled out with desired width and depth

☐	002JC20	JAY Bundled Cushion		
		Cushion Name or Part #		Price List
		Width: _____	Depth: _____	

2. SIZE / SEATING (CONT)



Seat Squeeze Matrix

Seat Depth	Max Seat Squeeze
14"-15"	2.0"
16"-19"	2.5"
20"	3.0"

3. CASTERS / SEAT HEIGHT

Front Seat Height

20" Max height for Fixed front. 20.5"-21" Only available in Swing Away

☐ 002SH42	14.5"	NC	☐ 002SH18	18.0"	NC
☐ 002SH12	15.0"	NC	☐ 002SH20	18.5"	NC
☐ 002SH13	15.5"	NC	☐ 002SH21	19.0"	NC
☐ 002SH14	16.0"	NC	☐ 002SH23	19.5"	NC
☐ 002SH15	16.5"	NC	☐ 002SH24	20.0"	NC
☐ 002SH16	17.0"	NC	☐ 002SH26	20.5"	NC
☐ 002SH17	17.5"	NC	☐ 002SH27	21.0"	NC


Caster Fork

☐ 002CF65	Multiple Position	NC
☐ 002CF69	Multiple Position Single Sided Fork	\$176.00
☐ 002CF53	Frog Legs Suspension (0-100 lbs)	\$574.00
☐ 002CF54	Frog Legs Suspension (101-190 lbs)	\$574.00
☐ 002CF55	Frog Legs Suspension (191-265 lbs)	\$574.00

Multi Position

Single side

Frog Leg



Caster Wheels

☐ 002CW1	3" Micro	NC
☐ 002CW3	4" Low Profile Polyurethane	NC
☐ 002CW31	4" Micro Lighted	\$110.00
☐ 002CW72	4" x 1" Aluminum Polyurethane	\$73.00
☐ 002CW44	4" x 1.25" Semi Pneumatic	NC
☐ 002CW64	4" x 1.5" Aluminum Soft Roll	\$110.00
☐ 002CW5	5" Low Profile Polyurethane	NC
☐ 002CW51	5" Micro Lighted	\$110.00
☐ 002CW6	5" Soft Roll (Grey) N/A w/ Frog Legs	\$98.00
☐ 002CW65	5" x 1" Aluminum Polyurethane	\$126.00
☐ 002CW7	5" x 1.5" Semi-Pneumatic	\$110.00
☐ 002CW66	5" x 1.5" Aluminum Soft Roll	NC
☐ 002CW34	6" Polyurethane	NC
☐ 002CW33	6" Pneumatic	\$110.00
☐ 002CW35	6" x 1.5" Semi Pneumatic	\$110.00
☐ 002CW68	6" x 1.5" Aluminum Soft Roll	NC

Caster Housing

☐ 002CP15	Integral Caster Housing	NC
☐ 002CP12	Hemi Integral Cstr Housing Subtract 2" from Standard Front STFH Matrix	NC

Front Seat Height Matrix

Integral Caster Housing (subtract 2" from FSH for Hemi Housing)												
Caster Size	*FSH - SA Front Frame	*FSH - Fixed Front Frame	Fork Hole (1 = Bottom Hole)									
			Multi Pos Low			Multi Pos Std			Frog Legs			
			Spacer Size:			Spacer Size:			Spacer Size:			
			0"	0.5"	1"	0"	0.5"	1"	0"	0.5"	1"	
3"	16.5	15.5	2									
	17.0	16.0	1	2								
	17.5	16.5		1	2	2						
	18.0	17.0			1	1	2		2			
	18.5	17.5					1	2	1	2		
	19.0	18.0						1		1	2	
4"	19.5	18.5										1
	20.0	19.0										
	17.0	16.0	2									
	17.5	16.5	1	2		3						
	18.0	17.0		1	2	2	3		3			
	18.5	17.5			1	1	2	3	2	3		
5"	19.0	18.0					1	2	3	2	3	
	19.5	18.5						1		1	2	
	20.0	19.0									1	
	20.5	19.5						1				1
	19.0	18.0				2						
	19.5	18.5				1	2					
6"	20.0	19.0					1	2	1			
	20.5	19.5						1		1		
	21.0	20.0									1	

*FSH values given are estimates. Actual seat to floor height may vary based on seat depth, squeeze & caster selection.

✓	Option #	Definition	Price	✓	Option #	Definition	Price
4. REAR WHEELS / SEAT HEIGHT							

Rear Seat Height						
13" & 13.5" N/A w/ transit						
☐	002RH11	13.0"	NC	☐	002RH18	16.5"
☐	002RH12	13.5"	NC	☐	002RH19	17.0"
☐	002RH13	14.0"	NC	☐	002RH26	17.5"
☐	002RH14	14.5"	NC	☐	002RH20	18.0"
☐	002RH15	15.0"	NC	☐	002RH27	18.5"
☐	002RH16	15.5"	NC	☐	002RH21	19.0"
☐	002RH17	16.0"	NC	☐	002RH37	19.5"

Axle			
☐	002AX1	Quick Release Steel	NC
☐	002AX11	Quad Quick Release Mountain Wheel only	\$182.00

Xtender Power Assist (265 weight capacity)			
☐	002XT6	22" Xtender Power Assist JWX-2 NICKEL	\$9,996.00
☐	002XT7	22" Xtender Power Assist JWX-2 LITHIUM	\$11,604.00
☐	002XT4	24" Xtender Power Assist JWX-2 NICKEL	\$9,996.00
☐	002XT5	24" Xtender Power Assist JWX-2 LITHIUM	\$11,604.00

Please select one of the following handrim options

☐	002XH5	Plastic Coated	NC
☐	002XH6	Stainless Steel	NC

Xtender - Other Options			
☐	002XA10	Spare NiMH Power Pack - NICKEL	\$1,925.00
☐	002XA11	Spare Li-ion Power Pack - LITHIUM	\$2,479.00

Xtender Wheel Options			
☐	002XA8	SPARE 22" Lite Spoke Pneu Alum	NC
☐	002XA9	SPARE 22" Lite Spoke Pneu Plastic	\$144.00
☐	002XA2	SPARE 24" Lite Spoke Pneu Alum	NC
☐	002XA4	SPARE 24" Lite Spoke Pneu Plastic	\$144.00

Rear Wheel Type			
-----------------	--	--	--

See availability matrix. Not applicable w/ Xtender Power Assist Option.

Ultra Lightweight Wheels only available w/ Alum Anodized Handrim

☐	002RW103	20" Mag Wheel 5 Spoke	NC
☐	002RW104	22" Mag Wheel 5 Spoke	NC
☐	002RW102	24" Mag Wheel 5 Spoke	NC
☐	002RW90	20" Lite Spoke Wheel	\$320.00
☐	002RW91	22" Lite Spoke Wheel	\$320.00
☐	002RW92	24" Lite Spoke Wheel	NC
☐	002RW93	25" Lite Spoke Wheel	NC
☐	002RW94	26" Lite Spoke Wheel	\$126.00
☐	002RW77	22" Ultra Lightweight Spoke Wheel	\$357.00
☐	002RW78	24" Ultra Lightweight Spoke Wheel	\$357.00
☐	002RW28	22" Spinerger SPOX Black	\$864.00
☐	002RW95	24" Spinerger SPOX Black	\$864.00
☐	002RW96	25" Spinerger SPOX Black	\$864.00
☐	002RW97	26" Spinerger SPOX Black	\$864.00
☐	002RW32	22" Spinerger LX Black	\$1,111.00
☐	002RW98	24" Spinerger LX Black	\$1,111.00
☐	002RW99	25" Spinerger LX Black	\$1,111.00
☐	002RW40	26" Spinerger LX Black	\$1,111.00
☐	002RW50	24" Mountain	\$393.00

Rear Wheel Type (Cont)			
* Omit Rear wheel will default chair configuration set up to a 24" wheel.			
☐	002RW25	Omit Rear Wheel and Axle *	NC
☐	002RW49	Omit Rear Wheel (Use w/24" Xtender)	NC
☐	002RW111	Omit Rear Wheel (Use w/22" Xtender)	NC

Rear Wheel Tire			
☐	002RT1	Pneumatic	NC
☐	002RT2	Pneumatic Airless Insert N/A with Xtender.	\$247.00
☐	002RT24	Pneumatic Airless Insert - Infinergy® Lightweight - E-TPU material	\$247.00
☐	002RT3	Full Poly	NC
☐	002RT38	SHOX Performance Poly	\$146.00
☐	002RT25	Primo Orion Poly	\$271.00
☐	002RT39	High-Pressure Clincher Black	\$105.00
☐	002RT27	Knobby V-TRAK	\$121.00
☐	002RT21	Schwalbe Marathon Plus Evolution	\$157.00
☐	002RT12	Mountain Tire Only available on 24" Mountain Rear Wheel	NC

Handrim			
---------	--	--	--

Natural Fits available in 22", 24", and 25". Ultra lightweight wheels ONLY available w/ alum ano handrim. Natural Fit w/ 5 spoke MAG only available in 24". Not all handrims available on all wheels

☐	002HR1	Aluminum Anodized	NC
☐	002HR4	Plastic Coated -N/A on every rim	\$146.00
☐	002HR22	Natural Fit - Standard Grip	\$374.00
☐	002HR25	Natural Fit - No Thumb Grip	\$374.00
☐	002HR23	Natural Fit - Super Grip	\$434.00
☐	002HR30	Natural Fit LT- Standard Grip	\$521.00
☐	002HR31	Natural Fit LT- Super Grip	\$592.00
☐	002HR32	Natural Fit LT - No Thumb Grip	\$557.00
☐	002HR26	The Surge - Oval w/Gription Strip	\$557.00
☐	002HR27	The Surge LT - Oval w/Gription Strip	\$592.00
☐	002HR28	Q-Grip Neoprene Coated	\$471.00
☐	002HR24	Spinerger Flex Rim Available on 24" & 25" Spinerger LX wheel	\$701.00
☐	002HR7	Omit Handrim Option may not be available on all wheels	NC
☐	002HR8	8 Oblique N/A w/ Xtender, Only available with 24" or 26" or Lite Spoke rims.	\$146.00

Rear Seat to Floor Height Matrix																
----------------------------------	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--

Note: These values are representative of a 16" seat depth with 1" of seat squeeze and standard axle plate. Due to variances caused by seat angle, rear seat height ranges may shift up or down by as much as 1/2". Call customer service or online configurator for more information.

Wheel/Tire	13.0	13.5	14.0	14.5	15.0	15.5	16.0	16.5	17.0	17.5	18.0	18.5	19.0	19.5
20" All	X	X	X	X	X	X	X	X	N/A	N/A	N/A	N/A	N/A	N/A
22" Poly	N/A	N/A	X	X	X	X	X	X	X	X	N/A	N/A	N/A	N/A
22" Pneu	N/A	N/A	N/A	X	X	X	X	X	X	X	X	N/A	N/A	N/A
24" Pneu	N/A	N/A	N/A	N/A	X	X	X	X	X	X	X	X	N/A	N/A
24" All but Pneu	N/A	N/A	N/A	N/A	X	X	X	X	X	X	X	X	N/A	N/A
25" SHOX	N/A	N/A	N/A	N/A	X	X	X	X	X	X	X	X	N/A	N/A
25" All but SHOX	N/A	N/A	N/A	N/A	N/A	X	X	X	X	X	X	X	X	N/A
26" All	N/A	N/A	N/A	N/A	N/A	N/A	X	X	X	X	X	X	X	X

✓	Option #	Definition	Price	✓	Option #	Definition	Price
4. REAR WHEELS / SEAT HEIGHT (cont)							

Wheel	5 Spoke Mag	Mag	Lite Spoke	Ultra Lt wt	Spinergy
Tire					
Pneumatic	20,22,24	20,22,24	20,22,24,26		24,26
Pneumatic w/Airless Insert	20,22,24	20,22,24	20,22,24,26		24,26
Full Poly Tire	20,22,24	20,22,24	20,22,24		
High Pressure Clincher black			22,24,26	22,24	22,24,26
Knobby V-TRAK			24		24
Schwalbe Marathon			24,25,26	22,24	24,25,26
SHOX Perform. Poly	20,22,24		24,25		24,25
20" = 451 mm 22" = 501 mm 24" = 540 mm 25" = 559 26" = 590					

Camber			
☐	002C11	0°	NC
☐	00217C1	2° (default)	NC
☐	002C13	4°	NC

Center of Gravity Preset			
☐	002CC2	1.0" (default)	NC
☐	002CC3	1.5"	NC
☐	002CC4	2.0"	NC
☐	002CC5	2.5"	NC
☐	002CC6	3.0"	NC
☐	002CC7	3.5"	NC

5. WHEEL LOCKS

Wheel Lock			
No Gos may exist depending on wheel size and position. Ergo & compact N/A w/ 13"-14" frame depth			
☐	002WL46	High Mount Push	NC
☐	002WL50	High Mount Pull	NC
☐	002WL75	Ergo Scissor	\$157.00
☐	002WL64	Compact	\$162.00
☐	002WL24	Unilateral Wheel Lock Right	\$320.00
☐	002WL23	Unilateral Wheel Lock Left	\$320.00
☐	002WL10	Omit	NC

Wheel Lock Option			
☐	002WL12	Do Not Mount	NC

Extension Handle			
N/A w/ ergo scissor or compact wheel locks			
☐	002EH7	Extension Handle (Pair)	\$73.00

Grade Aid			
Only available w/ high mount push or pull to lock wheel locks on pneumatic or pneumatic w/ airless tire options			
☐	002GA1	Grade Aid	\$157.00

6. HANGER & FOOTREST

Swing - Away Hangers			
Footplate mounts with Extension Tube using clamping hardware			
☐	002H145	60° Swing In/Out Hanger	NC
☐	002H139	70° Swing In/Out Hanger	NC
☐	002H141	80° Swing In/Out Hanger	NC
☐	002H32	Omit Hanger and Footplate	NC

Footrest Extension Tubes			
See Extension matrix for leg length range			
☐	002E1	Short	NC
☐	002E2	Medium	NC
☐	002E5	Long	NC

Z-finity Hangers/Footrest System			
(Maximum weight capacity for Z-Finity footrests is 165 lbs)			
Some interference may exist in certain configurations.			
☐	002H162	Z-finity Hanger	NC

Z-finity Footrest Mount			
☐	002DI22	Direct Hanger Mt/Fixed Frame - High	NC
☐	002DI23	Direct Hanger Mt/Fixed Frame - Low	NC
☐	002DI24	Extension Tube Mount	NC

Z-finity Footrest Pivot Arm			
☐	002DI20	Short 2"	NC
☐	002DI21	Long 3 1/4"	NC

Z-finity Footplate			
Some interference may exist in certain configurations.			
☐	002F320	Multi Angle Adjustable (Pair) (not avail. w/8"-9" widths)	\$247.00
☐	002F321	Angle Adjustable Platform Footboard	\$247.00
☐	002F6	Omit Footplate	NC

Z-finity ExtensionTube Length (When Choosing Ext. Tube Mount)			
☐	002E61	Hemi Short	NC
☐	002E1	Short	NC
☐	002E2	Medium	NC
☐	002E5	Long	NC

Extension Matrix

Frame Type	Footrest Angle	Inset		Footrest Type and Extension							
				Composite/ Aluminum			Platform		Tubular/ Autofold		
				Short	Medium	Long	Short	Long	Short	Long	
Fixed Front	70°	0°	Min	13.5	13.5	14.5	13.5	15.5	12	12.5	
			Max	17.5	19.5	20.5	16.5	20.5	16.5	18	
		1°	Min	13.5	15	16	13.5	17	12	13.5	
			Max	17.5	19.5	20.5	16.5	20.5	16.5	18	
	80°	0°	Min	14.5	14.5	15.5	14.5	16.5	13	13.5	
			Max	17.5	19.5	20.5	18	20.5	18	19.5	
		1°	Min	14.5	15	17	14.5	18	13	14.5	
			Max	17.5	19.5	20.5	18	20.5	18	19.5	
Swing Away	60°		Min	12	12	12.5	12	12.5			
			Max	13.5	15.5	17	13	17			
		70°	Min	12	12	14.5	12	14			
			Max	16.5	18.5	20	14	18.5			
	80°		Min	12	12	14	12	13.5			
			Max	15.5	17.5	19	13	17.5			
		70°		Min	12	12	14.5	12	14		
				Max	16.5	18.5	20	14	18.5		

✓	Option #	Definition	Price	✓	Option #	Definition	Price
6. HANGER & FOOTREST (cont)							
Footplates				Lower Leg Length Setting			
Heel loop or leg strap not included, See Footplate Options Section				Recommend 2" footrest clearance to ground surface			
☐	002F1	Composite (14"-20"w)	NC	☐	002LO07	12.0"	NC
☐	002F271	Composite Angle Adj-(16"-20"w)	\$176.00	☐	002LO08	12.5"	NC
☐	002F3	Aluminum Angle Adj (14"-20"w)	\$176.00	☐	002LO09	13.0"	NC
☐	002F48	Platform Flip Up	\$230.00	☐	002LO10	13.5"	NC
☐	002F245	Tubular Flip Up Fixed front frame only	\$150.00	☐	002LO11	14.0"	NC
☐	002F250	Auto Fold Angle Adj Fixed front frame only	\$298.00	☐	002LO12	14.5"	NC
Footplate Options				☐	002LO13	15.0"	NC
☐	002HL1	Heel Loop	NC	☐	002LO14	15.5"	NC
☐	002HL5	Leg Strap	NC	☐	002LO15	16.0"	NC
7. BACKREST / SEATING							
Back Type				JAY Back (Optional)			
☐	002BT10	Standard (9"-20")	NC	☐	002JB01	J3 Back (See J3 Order Form)	Price List
☐	002BT84	Angle Adjustable Back (3" & 4")	\$357.00				
Back Post Options				Angle Adjustable Back Setting			
Standard Back type options:				☐	002AH1	3" Pivot Point Setting	NC
☐	002PH6	Integral Push Handles	NC	☐	002AH2	4" Pivot Point Setting	NC
☐	002PH11	Integral Push Handles w/8° Medium & Tall only, N/A w/ Angle Adj Back	NC	Select desired setting			
☐	002PH2	Omit Push Handles	NC	☐	-14°	☐	-12°
☐	002PH12	Omit Push Handles w/8° Medium & Tall only, N/A w/ Angle Adj Back	NC	☐	-9°	☐	-7°
☐	002PH17	Folding Push Handle with Ht. Adj. Back N/A w/Angle Adj Back	\$310.00	☐	-1°	☐	0°
Angle Adj Back type options:				☐	+4°	☐	+5°
☐	002PH5	Fold-Down Push Handles N/A w/ standard Upholstery and J3 Back	\$168.00				
Back Height (Std or Angle Adjustable)				Backrest Height Matrix			
See Matrix							
☐	002B1	Low	NC				
☐	002B5	Medium	NC				
☐	002B9	Tall	NC				
Back Height Setting (Std or Angle Adjustable)							
17" max height for shipping. 21" only available w/ Ang Adj Back							
☐	002UX4	9"	NC	☐	002UX18	16"	NC
☐	002UX6	10"	NC	☐	002UX20	17"	NC
☐	002UX8	11"	NC	☐	002UX22	18"	NC
☐	002UX10	12"	NC	☐	002UX24	19"	NC
☐	002UX12	13"	NC	☐	002UX26	20"	NC
☐	002UX14	14"	NC	☐	002UX28	21" (*)	NC
☐	002UX16	15"	NC				
Backrest Upholstery							
Std EXO upholstery and Tension adjustable binding is always black.							
☐	002BU1	Standard Upholstery	NC				
☐	002BU95	Standard Upholstery EXO	NC				
Color Binding: ☐ Silver ☐ Orange ☐ Red ☐ Blue							
☐	002BU70	Backrest Tension Adjust Ballistic *Reference 186BU70 for billing	\$190.00				
☐	002BU89	Backrest Tension Adj EXO *Reference XXXBU89 for billing	\$449.00				
Color Binding: ☐ Silver ☐ Orange ☐ Red ☐ Blue							
☐	002BU71	Backrest Tension Adjust 3DX Vented *Reference 186BU71 for billing	\$471.00				
☐	002BU82	Jay GO Backrest	\$598.00				
☐	002BU14	Omit Back Upholstery	NC				

Back Type	Push Handles	Backrest Height	Avail. Ht. Range
Standard - (Fixed)	Integral	Low	9"-12"
		Medium	13"-16"
		Tall	17"-20"
	Fold Down	Medium	12"-15"
		Tall	16"-19"
		Low	9"-12"
Angle Adjustable (3" Pivot)	Omit	Medium	13"-16"
		Tall	17"-20"
		Low	13"-16"
	Integral	Medium	17"-20"
		Tall	N/A
		Fold	Medium
Angle Adjustable (4" Pivot)	Omit	Low	13"-16"
		Medium	17"-20"
		Tall	N/A
	Integral	Low	14"-17"
		Medium	18"-21"
		Tall	N/A
Fold	Medium	17"-20"	
	Low	14"-17"	
	Medium	18"-21"	
Folding Push Handle - (Stroller Handle w/ Knuckle)	N/A	Low	9"-12"
		Medium	13"-16"
		Tall	17"-20"

✓	Option #	Definition	Price	✓	Option #	Definition	Price
8. ARMREST / SIDE GUARDS							

Padded Swing Away			
<i>Some limitations when used w/ angle adj back in the rear most setting</i>			
☐	002AA18	Padded Swing-Away STD (9" to 13")	\$278.00
	Setting	☐ 9" ☐ 10" ☐ 11" ☐ 12" ☐ 13"	

Height Adjustable			
<i>4 Point (FS) J3 mounting bracket will cause interference with Flip Back Height Adj armrest.</i>			
☐	002AA11	Single Post Hgt. Adj. Std (8.5"-13.5")	\$339.00
☐	002AA29	Single Post Hgt. Adj. Low (7"-11")	\$339.00
☐	002AA65	Flip Back Height Adj <i>LITE</i> (8.5"-13.5")	\$381.00

Omit Armrest			
☐	002AA14	Omit Armrest	NC

Armpad for Height Adjustable			
<i>Ergo armpad only available with Single post Height adj & Flip back armrest</i>			
☐	002AM23	Desk Length - Classic (10")	NC
☐	002AM24	Full Length - Classic (14")	NC
☐	002AM27	Ergo Pad (9.5"-13.5")	NC

Side Guard			
☐	002SG33	Plastic Standard (8.5" High)	\$188.00
☐	002SG29	Plastic Low (6.5" High)	\$188.00

9. KEY OPTIONS & ACCESSORIES

Anti-Tip			
☐	002AT1	Rear Anti-Tip (NC w/ Xtender)	\$121.00

Backpack			
☐	002A1	Quickie Backpack Black	\$112.00
☐	002A2	Quickie Kids Backpack Black	\$105.00

Seat Pouch			
☐	0021A01	Seat Pouch Black /Charcoal	\$83.00

JAY Positioning			
☐	002K1	JAY Lower Extremity S/A support	NC
		☐ Left ☐ Right Position	
☐	002K2	Omit Hanger	NC
		☐ Left ☐ Right	
☐	002K7	Pad Size	NC
		☐ 6" x 8" ☐ 8" x 10"	

Additional Cover & Pack Accessories			
☐	0028A02	Gel Cover Hanger/Frame	\$162.00
☐	0028A01	Gel Cover Swing Away Arm	\$162.00

Positioning Belts			
☐	002A141	Belt 1.5" Auto Buckle w/Clamps	\$97.00
☐	002A117	Belt 2" Aircraft Buckle w/Clamps	\$139.00
☐	002A118	Belt 2" Aircraft Buckle Pad w/Clamp	\$204.00
☐	002A120	Positioning Belt Slide Clamps	\$61.00

Spoke Guards			
<i>Only available w/ Mountain, and Lite Spoke wheels</i>			
<i>Not available with 5 spoke MAG wheels</i>			
☐	0024A01	Quickie Black Spoke Guards	\$157.00
☐	0024A02	Quickie Clear Spoke Guards	\$157.00

Miscellaneous			
☐	0028A03	Caddy	\$89.00
☐	002IG2	Impact Guard Neoprene	\$61.00
☐	002LS14	Extra Leg Strap Black	\$37.00
☐	0028A21	Tool Kit	\$43.00
☐	002A999	Quickie Water Bottle & Holder	\$81.00

10. COLOR

Standard Frame Colors			
☐	002C140	Aztec Gold	NC
☐	002C2	Black (<i>Not avail as base w/Zebra</i>)	NC
☐	002C141	Black Cherry	NC
☐	002C10	Black Opal	NC
☐	002C38	Blue Opal	NC
☐	002C36	Candy Blue	NC
☐	002C14	Candy Purple	NC
☐	002C35	Candy Red	NC
☐	002C142	Electric Blue	NC
☐	002C143	Evergreen	NC
☐	002C44	Glossy White	NC
☐	002C126	Glow	NC
☐	002C144	Green Apple	NC
☐	002C147	Hot Sparkle Pink	NC
☐	002C148	Mauve Pink	NC
☐	002C103	Rootbeer	NC
☐	002C114	Sparkle Silver	NC
☐	002C145	Sunrise Orange	NC
☐	002C113	Titanium Color Paint	NC
☐	002C8	Yellow	NC

Matte Frame Colors			
☐	002C116	Matte Black	NC
☐	002C151	Matte Black Cherry	NC
☐	002C150	Matte Electric Blue	NC
☐	002C149	Matte Evergreen	NC
☐	002C152	Matte Purple	NC

Personalized Frame Color			
☐	002C112	Desert Camouflage	\$368.00
☐	002C43	Mossy Oak Camouflage	\$368.00
☐	002C42	Stars & Stripes	\$368.00
☐	002C108	Zebra (<i>Choose above base color</i>) Std with White unless above base paint is chosen N/A w/ Black	\$368.00

Anodized Color Package			
<i>Axle plate, Multi Forks and Aluminum Soft Roll casters only. N/A with Frog Legs.</i>			
☐	002EB07	Black - Caster only	NC
☐	002EB08	Blue	\$212.00
☐	002EB13	Gold	\$212.00
☐	002EB09	Lime Green	\$212.00
☐	002EB10	Purple	\$212.00
☐	002EB11	Red	\$212.00
☐	002EB12	Ti Grey	\$212.00



Sunrise Medical Canada Inc.

1000 Creditstone Road, Unit 2, Concord, ON L4K 4P8

Customer Service: 800-263-3390 Email: cscanada@sunmed.com www.sunrisemedical.ca © 4/2025

MK-100402-CAN Rev12